IN

cera**motion**[®]
One Touch

my first
ceramic
in paste

Kurzanweisung
Short instructions for use
Mode d'emploi condensé
Brevi istruzioni per l'uso
Instrucciones abreviadas

D
DENTAURUM

DE

EN

FR

IT

ES

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Technical data.

Classification CE 0483

ceraMotion® One Touch pastes are dental ceramics Class 1 (according to DIN EN ISO 6872:2015) for the finishing of frameworks made of zirconia and lithium disilicate.

Indications

Staining technique, cut-back technique, characterization and glazing of ceramic materials:

- Zirconia (fully anatomical, partial and full veneering with ceraMotion® Zr zirconia ceramic)
- Lithium disilicate (fully anatomical, partial and full veneering with ceraMotion® Zr zirconia ceramic)

Contraindications

- ceraMotion® One Touch is contraindicated for all applications that are not listed under the heading "Indications".
- ceraMotion® One Touch may not be used if there is a known intolerance to any of the ingredients.

ceraMotion®
One Touch

Physical data

	Test standard	Measurements	Werte
Coefficient of thermal expansion	ISO 6872 : 2015	CTE (25 °C/77 °F - Tg) [10 ⁻⁶ K ⁻¹]	2D pastes = 8.6 3D pastes = 9.1 Paste Glaze = 8.9
Flexural strength (3 point flex test)	ISO 6872 : 2015	> 50 MPa	> 115 MPa
Chemical solubility	ISO 6872 : 2015	< 100 µg/cm ²	≈ 30 µg/cm ²
Transformation temperature (Tg)	ISO 6872 : 2015	–	≈ 480 °C – 490 °C (896 °F - 914 °F)

Steps from beginning to end.

1

Basic

Framework preparation

- Lithium disilicate
Blast with Al₂O₃ 50µm/1–2 bar,
steam-clean (blasting abrasive,
REF 128-017-00).
- Zirconia
(Prepare according to
manufacturer's instructions).

2

Basic

Apply
ceraMotion® Zr Paste Glaze
The crown must be completely
covered by a thin layer of
ceraMotion® Zr Paste Glaze.

3

Individual

Apply 2D and 3D paste
You can use the enclosed brush
(3D Brush, REF 260-906-10)
to apply the 3D pastes.

4

Basic

First firing

Fire the finished restoration
at 730 °C/1346 °F-

5

Finished!

6

Individual

Corrections (optional)

- 2D pastes
- 3D pastes
- Zr Paste Glaze

7

Individual

Second firing

Fire the finished
restoration at
720 °C/1328 °F

Firing Chart.

The firing temperatures indicated are guidelines. Deviations are possible for different furnaces and the values listed below may need to be altered accordingly.

We recommend you run a trial firing to find the right temperature for your furnace. This is the only way to find the right firing procedure for your furnace.

 Please take note of the manufacturer's instructions when creating the framework.

Firing Chart

	Temperature at start (°C/°F)	Drying time (m)	Heat rate* (°C/°F per m)	Vacuum start (°C/°F)	Vacuum stop (°C/°F)	Firing temperature (°C/°F)	Hold-time** (m)
First firing cM 2D- and 3D-Paste cM Zr Paste Glaze	450 / 842	8	55 / 131	450 / 842	730 / 1346	730 / 1346	1 (without vacuum)
Second firing Correction cM 2D- and 3D-Paste cM Zr Paste Glaze	450 / 842	8	55 / 131	450 / 842	720 / 1328	720 / 1328	1 (without vacuum)

* Reducing the heat rate improves the firing quality in the case of large restorations.

** Extend the hold-time in the case of large restorations to compensate for the poor thermal conductivity of ZrO₂.

Processing tips.

- The ceraMotion® Stains Universal / colors in powder form can be used in addition. Please mix these either as usual with the ceraMotion® Stains Liquid or with the Diluting Liquid.

If you choose to use the ceraMotion® Universal Stains in powder form to characterize the restoration, you can mix this into the ceraMotion® Zr Paste Glaze. This will give you an individual 2D paste.

Please use the Diluting Liquid if necessary.

- If you would like to produce an individually colored 3D paste, mix the Universal Stains into the paste 3D neutral. Use the Refreshing Liquid if necessary.
- The required surface structure can be achieved by modeling the paste.
- The 3D pastes are used for the cut-back technique.

- To reduce the glaze level, turn the firing temperature down to 715 °C/1319 °F.
- To increase the glaze level, turn the firing temperature up to 750 °C/1382 °F max.
- You can influence the level of glazing individually by polishing the restoration mechanically.

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