




CE 0483

Triceram®

The ceramics for
titanium and zirconium oxide
Product Information and Instructions for Use

Titanium / Advantages at a glance

For contented patients

Titanium has a high resistance to corrosion

Titanium does not corrode. Following the establishment of a protective oxide layer, it is very slow to react. This is also referred to as „chemically inert“.

Titanium is bio-compatible

Titanium does not release ions and does not trigger pathological reactions in the body.

Titanium is the most bio-compatible metal available

It has proven to be very effective in human medicine and dentistry.

Titanium has low thermal conductivity

The patient can enjoy hot or cold food without experiencing dental pain. Gold has a thermal conductivity of 15 times higher than titanium, causing the patient to feel greater sensitivity to hot or cold.

Titanium has a neutral taste

Titanium’s passivated surface prevents a reaction between metal, saliva and food. The patient does not experience a metallic taste as with many other dental alloys.

Titanium has a multitude of indications

It is possible to use titanium for a variety of different indications, with the advantage of requiring just one metal for all restorations.

Titanium is light

This proves more comfortable for the patient. Gold is four times heavier, which is noticeable in large bridge framework.

Titanium has a high strength, without being too hard or too soft

Titanium can withstand the tremendous forces generated during mastication. Titanium is easy to process in the dental laboratory.

Illustrations: Photos: by Dentaurum
Page 20: Dentallabor Bischoff, Northeim

Contents

	Page
1. Technical Data	4
2. Dental Technical Instructions for Use on Titanium Framework.....	6
2.1. Framework design and conditioning	6
2.1.1. Basic factors	6
2.1.2. Trimming and grinding titanium framework.....	7
2.2. Ceramic build-up	8
2.2.1. Colour-combination table	8
2.2.2. Firing chart	8
2.2.3. Firing temperature control	9
2.2.4. Application of bonder.....	10
2.2.5. AESTHETIC UNIVERSAL Bonder	11
2.2.6. Applying the opaque	12
2.2.7. Use of shoulder material.....	13
2.2.8. Build-up with dentin and incisal material	15
2.2.8.1. Layer diagram	15
2.2.8.2. Note on materials and liquids	15
2.2.8.3. Layer build-up (standard) step-by-step	18
2.2.8.4. Individual layer build-up step-by-step.....	20
2.2.9. Use of gingival material	21
2.2.10. Use of correction material.....	21
2.2.11. Preparation of Triceram® ceramics	21
3. Dental Technical Instructions for Use on Zirconium Oxide Framework.....	22
4. Triceram® Product Components and Assortments.....	23
4.1. Assortments.....	23
4.2. Product components.....	25
5. Firing Charts for Titanium Framework.....	28
6. Firing Charts for Zirconium Oxide Framework.....	29

Guarantee

All data on products and processing are given to the best of our knowledge and technical experience, they are however without obligation. Our general terms for supply and payment apply exclusively. Our products undergo continual development. We therefore reserve the right to make changes at our own discretion.

1. Technical Data

CTE and Firing Temperatures

Pure titanium requires a ceramics with special characteristics for bonding.

On the one hand the ceramics must have a CTE (coefficient of thermal expansion) which is compatible with titanium ($9.6 \times 10^{-6} + K^{-1}$ 25 - 500°C) and on the other hand the ceramic firing must not heat the titanium over 800° C. This is due to intensive chemical reactions which occur on titanium when this temperature is exceeded.

Triceram® is a monophase, bonding ceramics with a unique capability of meeting these specific demands. The firing temperatures do not exceed 795° C and the CTE is compatible (see table: CTE values of individual Triceram® materials and harmonises ideally with titanium throughout firing). See Thermal expansion of titanium and Triceram®.



Flexural Strength and Resistance to Chemical Corrosion


Even More: Triceram® successfully fulfils and surpasses all minimum requirements according to the international standards (EN ISO 6872, EN ISO 9693) and has far lower chemical solubility levels than the maximum values permitted (see table: Solubility of individual Triceram® materials).

Bonding Strength

The bonding strength between ceramic and underlying framework played an important role in development. Triceram® exceeds all international standard demands EN ISO 9693 (25 MPa) by far, with results of 42 MPa (see graph: Bonding strength characterization). Aesthetical freedom in marginal design and at the same time maximal bonding strength is achieved through the application of the white coloured Triceram® bonder. New aesthetical limits and even greater bond strength can be achieved through the application of AESTHETIC UNIVERSAL bonder.

Reliability

Throughout the whole technological development of Triceram® bonding ceramics, the choice of raw materials with their various characteristics and ability to create high quality aesthetics, was of great importance. The production of CARMEN®, the bonding ceramics from Esprident for conventional alloys, has shown that the avoidance of natural raw materials and use of well proven production methods such as the thermo-coloration technique, all aid in achieving reliability.




Thermal expansion of titanium and Triceram®

Coefficient of thermal expansion (CTE 25-400 °C):

Triceram® bonder:	$9,2 \times 10^{-6} K^{-1}$
Triceram® opaque:	$8,8 \times 10^{-6} K^{-1}$
Triceram® dentin:	$8,6 \times 10^{-6} K^{-1}$
Triceram® incisal:	$8,5 \times 10^{-6} K^{-1}$

CTE values of individual Triceram® materials



Flexural strength

Solubility:


The maximum permitted solubility in the hydrolysis test all to ISO 9693 is 100 µg/cm²

Figures obtained for Triceram® titanium ceramics:

Opaque:	55 µg/cm²
Dentin:	31 µg/cm²
Incisal:	31 µg/cm²

The solubility is therefore well below the specified limit.

Solubility of individual Triceram® materials



Bonding strength characterization

2. Dental Technical Instructions for Use on Titanium Framework

2.1. Framework design and conditioning

2.1.1. Basic factors


Designing and conditioning the framework correctly to suit titanium is an indispensable factor in the production of durable, aesthetical, ceramic veneered restorations. A sound knowledge and understanding of titanium's characteristics is essential to achieve successful results.

This information can be acquired from manufacturers and suppliers of titanium casting machines, for example the Rematitan®-System from Dentaaurum.

Please follow the manufacturer's instructions for use regarding the wax-up, sprue construction, investment and casting.

There may be processing differences according to which manufacturer the material is supplied from. The following instructions for designing and conditioning the titanium framework are based upon the Rematitan®-system by Dentaaurum. Many aspects however can be applied to titanium framework produced through other systems.


Titanium frameworks can be made by casting, milling or eroding, as well as by using pre-fabricated elements (e.g. implant prosthetic components). However, the surface of the titanium framework must be entirely free of the α -case layer.



2.1.2. Trimming and grinding titanium framework

Important! Do not overheat the cast framework during grinding. Avoid red heat and cool well with water. Do not tilt the separating disc.

Recommended separating discs:	
① ST-separating disc for thick sprues	(REF 130-100-00)
② ST-mini separating disc (rapid grinder or handpiece)	(REF 130-110-00)
③ TX-separating disc for thin sprues	(REF 130-512-00)



Grinding in preparation for ceramic veneering:

Notes in general: Surface processing must be recognised as the most important factor in the manufacturing process of titanium. Depending on the casting system, the investment material and the wall thickness, the α -case layer formed by cast titanium may vary in thickness. This is enriched with oxygen and other substances and has a negative effect on the ceramic coat. More or less intensive grinding of the surface is therefore necessary, depending on the casting system, the investment material and the wall thickness.

Recommended rule of thumb:	
• Light grinding of crowns with thin to normal wall thickness.	
• Intense grinding of crowns with thicker walls or solid pontics .	

Use carbide cutters **only** and work in one direction. Especially suitable are the cross-cut carbide cutters specially designed for titanium. The α -case layer typical of titanium must be removed completely.

Recommended cutters of the Dentaaurum product range:

④ carbide cutter, mini	(REF 123-610-00)
⑤ carbide cutter, midi	(REF 123-611-00)
⑥ carbide cutter, maxi	(REF 123-612-00)
⑦ carbide cutter, maxi plus	(REF 123-613-00)
Maximum speed: 10 000 min ⁻¹ , low pressure.	



Due to titanium's low thermal conductivity, excessive pressure during grinding may result in overheating. This in turn prevents the production of metal filings, subsequently causing overlapping.

Never use tools for processing titanium with any other type of metal or material. Do not use grinding points or diamond burs for areas to be covered with dental ceramics (risk of gas bubbles in the ceramics!) Ensure that the cutter does not become clogged with metal chips. If necessary, clean it with a wire brush. Carefully sandblast the surface of the framework with a microblaster using aluminium oxide (125-250 μ m) and 2-3 bars air pressure. Recommended jet angle 45°.

At lower pressure, the conditioning of the surface is inadequate. At higher pressure there is a risk of surface contamination.

After having sandblasted and cleaned the metal framework (using a steam cleaner), do not touch it with your fingers.

Never immerse titanium frameworks into hydrofluoric acid as this severely corrodes titanium.

2.2. Ceramic build-up

2.2.1. Colour-combination table

Colour	A1	A2	A3	A3,5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
Opaque	O A1	O A2	O A3	O A3,5	O A4	O B1	O B2	O B3	O B4	O C1	O C2	O C3	O C4	O D2	O D3	O D4
Opaque dentin	OD A1	OD A2	OD A3	OD A3,5	OD A4	OD B1	OD B2	OD B3	OD B4	OD C1	OD C2	OD C3	OD C4	OD D2	OD D3	OD D4
Dentin	D A1	D A2	D A3	D A3,5	D A4	D B1	D B2	D B3	D B4	D C1	D C2	D C3	D C4	D D2	D D3	D D4
Incisal	IT 57	IT 57	IT 59	IT 59	IT 60	IT 57	IT 59	IT 59	IT 60	IT 60	IT 59	IT 59	IT 60	IT 60	IT 59	IT 59
Shoulder*	1/2 A	2/3 A	3/4 A	1/1 A	1/1 A	1/2 B	2/3 B	1/1 B	1/1 B	2/3 C	3/4 C	3/4 C	1/1 C	1/3 A 1/3 D	1/2 A 1/4 D	1/1 D

* Portion of corresponding shoulder material, remainder SM-modifier white.

2.2.2. Firing chart

	Base temperature*	Drying time	Heat rate	Vacuum start	Vacuum end	Final temperature	Holding time	Cooling time
Powder bonder bake	500 °C 932 °F	4 min	65 °C/min 149 °F/min	500 °C 932 °F	795 °C 1463 °F	795 °C 1463 °F	1 min under vacuum	0 min
Paste bonder bake	500 °C 932 °F	6 min	65 °C/min 149 °F/min	500 °C 932 °F	795 °C 1463 °F	795 °C 1463 °F	1 min under vacuum	0 min
Opaque bake 1+2	500 °C 932 °F	4 min	65 °C/min 149 °F/min	500 °C 932 °F	795 °C 1463 °F	795 °C 1463 °F	1 min under vacuum	0 min
Shoulder bake	500 °C 932 °F	6 min	55 °C/min 131 °F/min	500 °C 932 °F	785 °C 1445 °F	785 °C 1445 °F	1 min under vacuum	0 min
Dentin bake (with liquid LV+ univ.)**	500 °C 932 °F	6 min	55 °C/min 131 °F/min	500 °C 932 °F	755 °C 1391 °F	755 °C 1391 °F	1 min under vacuum	0 min
Correction bake (with liquid LV+ univ.)**	500 °C 932 °F	4 min	55 °C/min 131 °F/min	500 °C 932 °F	755 °C 1391 °F	755 °C 1391 °F	1 min under vacuum	0 min
Glaze bake	500 °C 932 °F	2 min	55 °C/min 131 °F/min	****	****	755 °C*** 1391 °F	1 min****	0 min
Correction material	500 °C 932 °F	4 min	55 °C/min 131 °F/min	500 °C 932 °F	715 °C 1319 °F	715 °C 1319 °F	1 min under vacuum	0 min
Gingival material	500 °C 932 °F	6 min	55 °C/min 131 °F/min	500 °C 932 °F	755 °C 1391 °F	755 °C 1391 °F	1 min under vacuum	0 min

* Positioning of objects only at correct standby temperatures.

** When using modelling liquid MV universal, observe different bake sequence (see page 16).

*** By reducing (adjusting) the firing temperature, variation in glaze can be achieved.

**** The glaze bake can be carried out under or without vacuum.The grade of glaze can be increased by prolonging the holding time.

The data shown are valid for fine-silver calibrated furnaces. These may vary from one furnace to another or from one type of furnace to another. Correction of the firing temperature is urgently recommended, as the furnaces are generally calibrated for higher temperatures. It is recommended to use a furnace specially for titanium ceramics.

When different alloys are fired in the same furnace, it is necessary to make frequent purge firings to the furnace according to the manufacturer's instructions. This prevents the metal framework from becoming contaminated with other alloy constituents and other negative effects of impurity.

2.2.3. Firing temperature control

In order to adjust the firing temperature of your furnace, we recommend a firing test, as this is the only way to judge the correct bake sequence.

For that purpose, use NT-material (neutral transparent) mixed with modelling liquid LV universal and fire at the following temperatures:

- Base temperature 500 °C/932 °F
- Drying time 6 minutes
- Heating rate 55 °C/min. / 131 °F/min.
- Vacuum start 500 °C/932 °F
- Vacuum end at final temperature 755 °C/1391 °F
- Holding time 1 minute under vacuum

Put the bake specimen on platinum foil, not on firing cotton, as this may cause dulling.



Fig. 1: Bake specimen, correct.



Fig. 1a: Bake specimen, temperature too low.

The temperature of the furnace is correct when the specimen is clear and translucent with sharp edges (fig. 1). If the final temperature is too high, the specimen has a glossy appearance and the edges are not sharply defined. If the final temperature is too low, the specimen appears milky white (fig. 1a). According to your result, please raise or lower the final temperature in steps of 10 °C/50 °F and fire a new specimen.

Furnace cleaning

The ceramic furnace must be cleaned regularly in order to remove contamination from the inner surfaces of the muffle. We recommend:

- Frequent cleansing firing with carbon fibre chips (REF 260-317-00)
 - Cleansing of firing tray with pins
 - Base temperature: 600 °C/1112 °F
 - Drying time: 1 min
 - Heating rate: 100 – 120 °C/212 – 248 °F per min
 - Final temperature: 1050 °C/1922 °F
 - Holding time: 10 min.
- Firing program without vacuum. Please follow furnace manufacturer's instructions!

Always keep the furnace closed. Close the furnace after use or switch to night modus to prevent uptake of humidity.

2.2.4. Application of bonder

After sandblasting, a resting time of 10 minutes in air is necessary to allow the surface oxide to form. Do not wait for longer than 30 minutes before applying the bonder.

	Base temperature	Drying time	Heat rate	Vacuum start	Vacuum end	Final temperature	Holding time	Cooling time
Powder bonder bake	500 °C 932 °F	4 min	65 °C/min 149 °F/min	500 °C 932 °F	795 °C 1463 °F	795 °C 1463 °F	1 min under vacuum	0 min
Paste bonder bake	500 °C 932 °F	6 min	65 °C/min 149 °F/min	500 °C 932 °F	795 °C 1463 °F	795 °C 1463 °F	1 min under vacuum	0 min

When using the paste bonder, please be sure to mix thoroughly in the pot before application. If necessary add universal paste liquide, REF 299-170-40, to create a creamy consistency. Caution: Please observe the longer drying times in the firing chart.

Using a glass spatula the bonder powder is mixed with unversal liquid B.O.L., REF 299-180-40, to a creamy consistency (fig. 2).

Apply bonder thinly on the framework by glass spatula or brush, as if wetting the surface. Avoid the formation of puddles or drops on framework surface. The actual masking of the framework is achieved after having applied the opaque.

The bonder is not the opaque!

After baking, the bonder has an even shiny surface.



Fig. 2:
Mixing of the bonder to a creamy consistency with the glass spatula.



Fig. 3a:
3-part bridge on model. Metal framework fully conditioned. 1 crown shortened for application of shoulder material.



Fig. 3b:
3-part bridge with bonder applied (before bonder bake).



Fig. 3c:
3-part bridge (after bonder bake).

2.2.5. AESTHETIC UNIVERSAL Bonder

Application

Note: When using AESTHETIC UNIVERSAL bonder on a titanium framework, Triceram® titanium bonder should no longer be applied!

1. Prepare the metal framework according to chapter 2.1. Sandblast the metal surface and clean thoroughly. Warm the titanium framework gently using the heat from an open ceramic furnace.

Storage information:
During long periods of storage the gold content inside the syringe can separate from the paste. For this reason we recommend turning the syringe from time to time. Should a dispersion of the mixture be visible, please proceed as follows:

- Remove the lid.
- Withdraw the plunger as far as possible and stand the syringe upright with the plunger at the base.
- The paste will sink toward the base.
- Push the plunger back to its original position until the paste can be seen in the syringe head again. Repeat this procedure two or three times.
- If the AESTHETIC UNIVERSAL bonder becomes difficult to spread it is possible to add a small amount of AESTHETIC UNIVERSAL bonder liquid, REF 299-200-10. Place the paste as required in a small bowl and mix with a small amount of liquid.

Important: Use the liquid sparingly! Once in use, the object must be dried for a longer period of time, otherwise danger of bubbling.

2. Apply the AESTHETIC UNIVERSAL bonder in a thin layer to the titanium framework using a clean brush (Esprident brush no. 4, REF 260-905-00). Use gentle pressure and light rotation of the brush on the surface to cover the titanium framework uniformly.
3. Apply a thin covering layer of AESTHETIC UNIVERSAL bonder. Apply AESTHETIC UNIVERSAL bonder uniformly to ensure that no air becomes trapped.
4. After application, dry the covered surface for at least 5 min. in front of an opened ceramic furnace at the standby temperature of 400 °C / 752 °F. The shiny paste will turn matt. The object must not exceed 150 °C / 302 °F. If necessary, a second layer of AESTHETIC UNIVERSAL bonder may be applied on top of the dried bonder. Again dry thoroughly.

Important: It is essential to dry the AESTHETIC UNIVERSAL bonder thoroughly before firing. Otherwise the remaining bonder liquid within the superficially dried paste will break free later during ceramic firing. These inclusions can cause bubbles to occur. During drying, the objects should reach a temperature of between 100- 120 °C/212- 248 °F irrespective of the furnace program, but never exceed 150 °C/302 °F. If the open furnace does not generate enough heat, raise the standby temperature. The pre-heating furnace can also be used as an alternative (increase heat slowly to 130°C/266°F and hold for 5 minutes) or with a hot air blower (set at 120 °C/248 °F for 5 minutes).

5. After firing, an even distribution of ceramic particles should be visible on the titanium surface (rough surface). If necessary, apply a second layer
6. Then apply the opaque according to chapter 2.2.6. Only one firing is generally required.

Availability:

Set
(syringe 4g, brush,
liquid and small bowl)
REF 292-200-15

Syringe
(4g)
REF 292-200-10

Brush No. 4
REF 260-905-00

Liquid
REF 299-200-10

Firing chart

	Base temperature	Drying time	Heat rate	Vacuum start	Vacuum end	Final temperature	Holding time	Cooling time
AESTHETIC UNIVERSAL bonder	400 °C 752 °F	min. 5 min	55 °C/min 131 °F/min	400 °C 752 °F	800 °C 1472 °F	800 °C 1472 °F	1 min under vacuum	0 min

2.2.6. Applying the opaque

Using a glass spatula, mix the opaque with the universal liquid B.O.L. to a creamy consistency.

	Base temperature	Drying time	Heat rate	Vacuum start	Vacuum end	Final temperature	Holding time	Cooling time
Opaque bake 1+2	500 °C 932 °F	4 min	65 °C/min 149 °F/min	500 °C 932 °F	795 °C 1463 °F	795 °C 1463 °F	1 min under vacuum	0 min



Fig. 4a:
3-part bridge before first opaque bake.

For the first bake, apply the opaque thinly but evenly to the bonder covered framework (fig. 4a). Use the rounded tip of the glass spatula. After baking, the opaque appears silky and shiny (fig. 4b).

For the second opaque bake, the opaque is also applied in a thin even layer, taking care to cover the entire surface. The result should appear dense and shiny (fig. 4c).



Fig. 4b:
3-part bridge after first opaque bake.



Fig. 4c:
3-part bridge after second opaque bake.

2.2.7. Use of shoulder material

The Triceram® system comprises four different shoulder materials which are divided into the colour categories A-B-C-D.

All colour variations from A1 to D4 can be achieved by mixing the SM-modifier "white" in accordance with the table shown opposite.

The transparency of the shoulder and the effect of depth can be increased by adding the SM-modifier "transparent". For mixing use the shoulder material liquid universal.

To achieve the best results when lifting off the applied shoulder material and avoid any discoloration in the shoulder zone, it is recommended to work with the isolation SM-Isokit.

	A	B	C	D	Blanc
A1	1/2	–	–	–	1/2
A2	2/3	–	–	–	1/3
A3	3/4	–	–	–	1/4
A3,5	1/1	–	–	–	–
A4	1/1	–	–	–	–
B1	–	1/2	–	–	1/2
B2	–	2/3	–	–	1/3
B3	–	1/1	–	–	–
B4	–	1/1	–	–	–
C1	–	–	2/3	–	1/3
C2	–	–	3/4	–	1/4
C3	–	–	3/4	–	1/4
C4	–	–	1/1	–	–
D2	1/3	–	–	1/3	1/3
D3	1/2	–	–	1/4	1/4
D4	–	–	–	1/1	–



SM-Isokit
REF 260-324-01

25 ml

Isolates plaster from ceramics and ceramics from ceramics. Smoothes out and seals at the same time.

The separating effect is fantastic. Optimal results are achieved on dies and basal areas which have not been previously treated with hardener, sealant or die spacer. The surface must be dry. The isolation should be applied twice.

Application of shoulder material

Isolate the die with SM-Isoprotector in the area where the shoulder material is to be applied. The shoulder material, which has been selected to match the tooth shade, can be applied immediately.

Mix the shoulder material with corresponding shoulder material liquid universal. Apply the shoulder material from the framework towards the preparation line. Dry the material (hair-dryer or tissue). Then remove crown from the die and fire according to indicated program.

	Base temperature	Drying time	Heat rate	Vacuum start	Vacuum end	Final temperature	Holding time	Cooling time
Shoulder bake	500 °C 932 °F	6 min	55 °C/min 131 °F/min	500 °C 932 °F	785 °C 1445 °F	785 °C 1445 °F	1 min under vacuum	0 min

To compensate for shrinkage of the shoulder material which is to be expected due to the reduction of the crown (according to the preparation), a second bake is required. Re-separate the areas for correction with SM-Isoprotector and add the missing material. Bake as for first shoulder bake.

After the second bake, remount the crown onto the die and adjust the marginal edge with adequate finishing burs. (Observe the speed of the handpiece max. 15 000 min⁻¹). Continue the process in the usual way.



Fig. 5a:
Crown shortened after opaque bake. Shoulder material applied before shoulder bake.



Fig. 5b:
Crown after first shoulder bake.



Fig. 5c:
Crown before second shoulder bake. Shoulder material applied.



Fig. 5d:
Crown after fitting.

2.2.8. Build-up with dentin and incisal material

2.2.8.1. Layer diagram

Selection of the ceramic materials according to tooth shade. For build-up of the shade refer to the layer diagram opposite.

Note on layer build-up: To avoid visible colour demarcations, use opaque dentin (OD). It also prevents the visibility of the opaqued titanium framework. OD materials can also be used to control brightness (by mixing dentin and opaque dentin). At layer thicknesses from 0.9 mm on, opaque dentin is not required.

2.2.8.2. Note on materials and liquids


All the materials – opaque dentin, dentin, incisal, neutral transparent, transpa-effect, gingival and correction materials – are mixed using the modelling liquid LV universal. The ceramic materials should be kept as evenly moist as possible. Any additional moistening is done using modelling liquid LV universal.

Opaque dentin OD

- diffusing effect
- three-dimensional effect of depth (when completed)
- complementing of actual tooth shade (even in an unfavourable space situation)
- intensifies depth reflections

Dentin D

- 16 dentin colours
- easy handling
- stability of shape and edges
- bake stability
- low shrinkage
- colour stability
- prismatic effect
- natural fluorescence



Incisal transparent IT

Classification of the incisal materials:

A1	IT 57	B4	IT 60
A2	IT 57	C1	IT 60
A3	IT 59	C2	IT 59
A3,5	IT 59	C3	IT 59
A4	IT 60	C4	IT 60
B1	IT 57	D2	IT 60
B2	IT 59	D3	IT 59
B3	IT 59	D4	IT 59

With IT 58 it is possible to achieve additional grey effects in the area of the posteriors.

Incisal opalescent IO

- for incisal alternative layer
- opalescent bluish effect in reflected light, slightly orange effect in permeated light
- intensifies natural refraction of light
- application in incisal area

Transpa-effect materials TE

In order to increase the lively effect in the incisal area, or to obtain individual or characteristic effects, TE-materials can also be used with Triceram®.

- incisal alternative layer
- application in incisal area
- natural refraction of light

Stains Universal ST

To obtain modifications in shade. They can be laid in and mixed with the materials. Stains can be mixed up to 10 mass% with the ceramic materials (SM, OD, D, IT, IO, TE, NT + G). If they are mixed with Stains Universal liquid, these materials are also suitable for painting the surface of the ceramic material.

This material has highly individual characteristics and special depth of colour.

Body ST A	ST 7 dark pink
Body ST B	ST 8 purple
Body ST C	ST 9 blue
ST 0 neutral	ST 10 grey
ST 1 white	ST 11 olive-green
ST 2 vanilla	ST 12 olive-yellow
ST 3 yellow	ST 13 medium brown
ST 45 orange plus	ST 14 red brown
ST 6 pink	ST 15 black

Note:

Stains Universal ST can also be used in combination with the CARMEN® system (for mixing and painting).

The liquids:

Stains Universal liquid

- for mixing with Stains Universal ST
- for the glaze bake

Shoulder material liquid universal

- for mixing shoulder materials
- increases the stability of shape and edges

Universal liquid B.O.L.

- for mixing powder bonder, opaque and liner

Universal paste liquid

- for mixing paste bonder

Modelling liquid MV universal

- for increasing plasticity
- especially in dry warm climates
- for prolonging the carving ability

Modelling liquid LV Standard universal and liquid LV+ unsiversal

LV Standard universal

- universal standard liquid for mixing and re-moistening of previously mixed materials

LV+ universal

- the material shrinkage has been reduced to a minimum
- new additives keep the mixing liquid moist for longer. Ceramic remains moist and easy to use throughout the entire build-up

Important:

When using modelling liquid MV universal, note the differences in bake sequence from modelling liquid LV universal (see firing chart). Do not mix the universal liquids LV and MV as this may cause bubbles!

Note:

Universal liquids can also be used in combination with the CARMEN® system.

Bake sequence for dentin and correction bakes, see firing chart:

For standard modelling liquid LV universal

	Base temperature	Drying time	Holding time at base temp.	Heat rate	Vacuum start	Vacuum end	Final temperature	Holding time	Cooling time
Dentin bake (with liquid LV+ universal)	500 °C 932 °F	6 min	0 min	55 °C/min 131 °F/min	500 °C 932 °F	755 °C 1391 °F	755 °C 1391 °F	1 min under vacuum	0 min
Correction bake (with liquid LV+ universal)	500 °C 932 °F	4 min	0 min	55 °C/min 131 °F/min	500 °C 932 °F	755 °C 1391 °F	755 °C 1391 °F	1 min under vacuum	0 min

For modelling liquid MV universal

Dentin bake (with liquid MV univ.)	500 °C 932 °F	4 min	4 min	55 °C/min 131 °F/min	500 °C 932 °F	755 °C 1391 °F	755 °C 1391 °F	1 min under vacuum	0 min
Correction bake (with liquid MV univ.)	500 °C 932 °F	3 min	3 min	55 °C/min 131 °F/min	500 °C 932 °F	755 °C 1391 °F	755 °C 1391 °F	1 min under vacuum	0 min



Triceram®

Incisal opalescent material

For reinforcing natural opalescence

2.2.8.3. Layer build-up (standard) step-by-step



1. Opaque dentin (OD) is applied to the entire labial surface. The OD projects well beyond the incisal edge of the titanium framework (increase in light dispersion).



2. Build-up of the complete anatomical tooth form with dentin.



3. Cutting back the dentin in incisal area to make room for completion with incisal transparent material (IT).



4. Layering of incisal area with IT material.



5. When building up a bridge, separate the ceramics between the teeth, down to the opaque layer, using a scalpel (labially as well as palatinally).



6. Correction of shape after dentin bake (before correction bake).



7. After the correction bake, make the necessary corrections to shape (for suitable separating discs see 4, page 23).



8. Grind the bridge lightly using a diamond bur (suitable diamonds see 4 page 23).



9. The glaze bake is carried out without glaze material. Special colour effects can be achieved by applying stains to the ceramic surface. Mix the Stains with the stains universal liquid.



10. Finished work after glaze bake.

Bake sequence for dentin and correction bakes, see firing chart:

For standard modelling liquid LV universal

	Base temperature	Drying time	Holding time at base temp.	Heat rate	Vacuum start	Vacuum end	Final temperature	Holding time	Cooling time
Dentin bake (with liquid LV+ universal)	500 °C 932 °F	6 min	0 min	55 °C/min 131 °F/min	500 °C 932 °F	755 °C 1391 °F	755 °C 1391 °F	1 min under vacuum	0 min
Correction bake (with liquid LV+ universal)	500 °C 932 °F	4 min	0 min	55 °C/min 131 °F/min	500 °C 932 °F	755 °C 1391 °F	755 °C 1391 °F	1 min under vacuum	0 min

For modelling liquid MV universal

Dentin bake (with liquid MV univ.)	500 °C 932 °F	4 min	4 min	55 °C/min 131 °F/min	500 °C 932 °F	755 °C 1391 °F	755 °C 1391 °F	1 min under vacuum	0 min
Correction bake (with liquid MV univ.)	500 °C 932 °F	3 min	3 min	55 °C/min 131 °F/min	500 °C 932 °F	755 °C 1391 °F	755 °C 1391 °F	1 min under vacuum	0 min

Bake sequence for glaze bake, see firing chart

	500 °C 932 °F	2 min		55 °C/min 131 °F/min	--	--	755 °C 1391 °F	1 min **	0 min
--	------------------	-------	--	-------------------------	----	----	-------------------	----------	-------

* By reducing (adjusting) the firing temperature, variation in glaze can be achieved.
** The glaze bake can be carried out under or without vacuum. The grade of glaze can be increased by prolonging the holding time.

2.2.8.4. Individual layer build-up step-by-step



1. Colour modifications with stains on baked opaque.



2. Application of opaque dentin. More individual effects can be achieved by carefully mixing OD and stains together.



3. An individual build-up using dentin, dentin and stains mixed together and finished off with incisal and transpa effect materials..



4. Shape corrections and application of effect or transpa materials as required.



5. After the correction bake and final shaping, stains can be used to achieve special effects on the surface.



6. Finished bridge after glaze bake.

2.2.9. Use of gingival material

Gingival materials G (light) and GD (dark)

- build-up and reconstruction of gingival and papilla situation
- compensation of mucous-membrane deficiencies
- for masking parts of titanium framework with gingival coloured opaque wherever necessary

	Base temperature	Drying time	Heat rate	Vacuum start	Vacuum end	Final temperature	Holding time	Cooling time
Gingival material	500 °C 932 °F	6 min	55 °C/min 131 °F/min	500 °C 932 °F	755 °C 1391 °F	755 °C 1391 °F	1 min under vacuum	0 min

2.2.10. Use of correction material

Use for correcting minor defects (contact points, occlusion area).

In order to obtain perfect results from an aesthetic point of view, it is recommended to mix the correction material with dentin D or incisal material IT (proportion as required - 50:50).

	Base temperature	Drying time	Heat rate	Vacuum start	Vacuum end	Final temperature	Holding time	Cooling time
Gingival material	500 °C 932 °F	6 min	55 °C/min 131 °F/min	500 °C 932 °F	715 °C 1319 °F	715 °C 1319 °F	1 min under vacuum	0 min

2.2.11. Preparation of Triceram® ceramics

To achieve good surface density and homogeneity, it is recommended to grind over the surface with super-fine sintered diamonds or sintered diamond discs. Finishing, structuring and polishing (before glaze bake) are carried out with diadurite-impregnated special fibre polishers and micro-fine diamond-impregnated felt polishers (for interdental and occlusal polishing). With Triceram® glaze materials are not required.

3. Dental Technical Instructions for Use on Zirconium Oxide Framework

Zirconium oxide has different thermal characteristics to titanium. For this reason, in order to achieve optimal firing results of Triceram® on zirconium oxide, the firing temperature must be raised and the holding time lengthened. The same basic prerequisites for metal alloy framework, also apply for zirconium oxide. It is essential that the size of the framework is a smaller version of the final anatomical shape desired. A successful restoration relies on an even thickness of bonding ceramics. Please note the firing recommendations in the table below.

1. Preparing the framework
Observe instruction for use of your zirconium oxide provider.

2. Applying the Triceram® Fluo liner
Apply the Fluo liner according to the selected tooth shade (see colour-combination table) and mix with universal liquid B.O.L., REF 299-180-40, to a creamy consistency. For the first bake, apply the liner thinly to the framework. For the second liner bake, the liner is also applied in a thin even layer, taking care to cover the entire surface. After the second bake, the liner has a silky, shiny appearance

3. Ceramic build-up
(see page 18)
Please note that all firing temperatures must be increased by 5 °C - 10 °C compared to the firing recommendations for titanium, (see firing table below).
The holding time for the dentin bakes must be prolonged by 1.5 – 2 minutes.
Different furnaces may require individual variations in temperature. Tips for furnace control see page 9 Triceram® Instructions for Use.

Colour-combination table

tooth shade	FL 1	FL 2	FL 3	FL 4	FL 5	FL 6
A1	1/3	–	–	–	2/3	–
A2	2/3	–	–	–	1/3	–
A3	1	–	–	–	–	–
A3,5	3/4	–	–	–	–	1/4
A4	1/2	–	–	–	–	1/2
B1	–	1/3	–	–	2/3	–
B2	–	2/3	–	–	1/3	–
B3	–	3/4	–	–	1/4	–
B4	–	1/1	–	–	–	–
C1	–	–	1/3	–	2/3	–
C2	–	–	2/3	–	1/3	–
C3	–	–	3/4	–	1/4	–
C4	–	–	1	–	–	–
D2	–	–	–	2/3	1/3	–
D3	–	–	–	3/4	1/4	–
D4	–	1/3	–	2/3	–	–

Recommended Firing Temperatures for Zirconium Oxide Framework

	Base temperature	Drying time	Heat rate	Vacuum start	Vacuum end	Final temperature	Holding time	Cooling time
Fluo Liner 1+2 bake	500 °C 932 °F	4 min	65 °C/min 149 °F/min	500 °C 932 °F	800 °C 1472 °F	800 °C 1472 °F	1 min under vacuum	0 min
Shoulder bake	500 °C 932 °F	6 min	55 °C/min 131 °F/min	500 °C 932 °F	790 °C 1454 °F	790 °C 1454 °F	1 min under vacuum	0 min
1 Dentin bake	500 °C 932 °F	6 min	55 °C/min 131 °F/min	500 °C 932 °F	760 °C 1400 °F	760 °C 1400 °F	1,5 - 2 min under vacuum	0 min
2. Dentin bake	500 °C 932 °F	4 min	55 °C/min 131 °F/min	500 °C 932 °F	760 °C 1400 °F	760 °C 1400 °F	1,5 - 2 min under vacuum	0 min
Glaze bake	500 °C 932 °F	2 min	55 °C/min 131 °F/min	--	--	760 °C 1400 °F	1 min.*	0 min

* The glaze bake can be carried out under or without vacuum. The grade of glaze can be increased by prolonging the holding time.
Instead of using the Triceram® liner you can also use Triceram® bonder and Triceram® opaque.

4. Triceram® Product Components and Assortments

4.1. Assortments

Triceram® Starter set

contains:

1 Powder bonder	3 g
1 Paste bonder	3 g
1 Fluo Liner 1	20 g
1 Opaque O A3	15 g
1 Opaque Dentin OD A3	15 g
1 Dentin D A3	15 g
1 Incisal transparent IT 59	15 g
1 Neutral transparent NT	15 g
1 Universal B.O.L. liquid	20 ml
1 Modelling liquid LV+ universal	20 ml
1 Modelling liquid LV standard	20 ml
1 Glass spatula for applying bonder and opaque	
1 Opaque brush	

290-113-02	1 piece
------------	---------

CE 0483



Triceram® Master set Ti (1)

contains:

1 Powder bonder	15 g
17 Opaque O A1-D4 + G	15 g each
16 Opaque Dentin OD A1-D4	15 g each
16 Dentin D A1-D4	15 g each
4 Incisal transparent IT 57-60	15 g each
5 Incisal opalescent IO 1-5	15 g each
1 Neutral transparent NT	15 g
1 Correction CM	3 g
6 Transpa-effect TE 1-6	15 g each
4 Shoulder SM A-D	15 g each
2 Shoulder modifiers SM W, SM T	15 g each
2 Gingival G, GD	15 g each
15 Stains universal ST 0-15	3 g each
1 Stains universal liquid	20 ml
1 Shoulder material liquid univ.	20 ml
1 Universal B.O.L. liquid	100 ml
1 Modelling liquid MV univ.	20 ml
1 Modelling liquid LV+ univ.	100 ml
1 Modelling liquid LV standard	20 ml
1 SM-Isofit	20 ml
1 Glass spatula for applying bonder and opaque	
1 Opaque brush No. 6	
1 Shade guide I Basic materials	
1 Shade guide II Effect materials	
1 Shade guide III Effect materials	
1 Shade guide for Stains Universal	

290-112-01	1 piece
------------	---------

CE 0483





290-112-11	1 piece
------------	---------

CE 0483

Triceram® Master set Ti (2)

Contents identical to REF 290-112-01, with paste bonder instead of powder bonder.



290-112-03	1 piece
------------	---------

CE 0483

Triceram® Master set Zr (3)

Contents identical to REF 290-112-01, with
6 Fluo Liner FL 1-FL 6 20 g each
4 Intensive Fluo Dentin 15 g each
3 Chroma Dentin 15 g each

Without powder bonder, paste bonder and opaque A1-D4.



296-000-01	1 piece
------------	---------

CE 0483

Stains Universal set

contains:

15 Stains Universal 0-15 3 g each
1 Stains Universal liquid 20 ml
1 Shade guide for Stains Universal



Aesthetic Universal Bonder set	292-200-15	4 g
Aesthetic Universal Bonder	292-200-10	4 g

CE 0483

Aesthetic Universal Bonder

- stronger bond between the metal and ceramics
- optimum reproduction of the shade, especially in the marginal area of the crown



Powder bonder	292-100-10	3 g
Powder bonder	292-100-30	15 g
Powder bonder	292-100-60	40 g
Paste bonder	292-110-05	3 g
Paste bonder	292-110-10	8 g

CE 0483

Bonder

- creates a connection between the titanium and the ceramic (bond)
- recommended liquid: Universal paste liquid

4.2. Product components

FL Fluo Liner

Fluo Liner 1	292-514-40	20 g
Fluo Liner 2	292-524-40	20 g
Fluo Liner 3	292-534-40	20 g
Fluo Liner 4	292-544-40	20 g
Fluo Liner 5	292-554-40	20 g
Fluo Liner 6	292-564-40	20 g

O Opaque

Opaque A1	292-211-30	15 g
Opaque A1	292-211-60	40 g
Opaque A2	292-212-30	15 g
Opaque A2	292-212-60	40 g
Opaque A3	292-213-30	15 g
Opaque A3	292-213-60	40 g
Opaque A3.5	292-214-30	15 g
Opaque A3.5	292-214-60	40 g
Opaque A4	292-215-30	15 g
Opaque A4	292-215-60	40 g
Opaque B1	292-221-30	15 g
Opaque B1	292-221-60	40 g
Opaque B2	292-222-30	15 g
Opaque B2	292-222-60	40 g
Opaque B3	292-223-30	15 g
Opaque B3	292-223-60	40 g
Opaque B4	292-224-30	15 g
Opaque B4	292-224-60	40 g
Opaque C1	292-231-30	15 g
Opaque C1	292-231-60	40 g
Opaque C2	292-232-30	15 g
Opaque C2	292-232-60	40 g
Opaque C3	292-233-30	15 g
Opaque C3	292-233-60	40 g
Opaque C4	292-234-30	15 g
Opaque C4	292-234-60	40 g
Opaque D2	292-242-30	15 g
Opaque D2	292-242-60	40 g

Opaque Gingiva

Opaque D3	292-243-30	15 g
Opaque D3	292-243-60	40 g
Opaque D4	292-244-30	15 g
Opaque D4	292-244-60	40 g
Opaque Gingiva	292-250-30	15 g
Chroma Opaque CO 1 Bleach	292-301-30	15 g

SM Shoulder material

Shoulder material A	293-110-30	15 g
Shoulder material B	293-120-30	15 g
Shoulder material C	293-130-30	15 g
Shoulder material D	293-140-30	15 g

SM Modifier

SM Modifier white	293-150-30	15 g
SM Modifier transparent	293-199-30	15 g

OD Opaque Dentine

Opaque Dentine A1	294-111-30	15 g
Opaque Dentine A1	294-111-60	40 g
Opaque Dentine A2	294-112-30	15 g
Opaque Dentine A2	294-112-60	40 g
Opaque Dentine A3	294-113-30	15 g
Opaque Dentine A3	294-113-60	40 g
Opaque Dentine A3.5	294-114-30	15 g
Opaque Dentine A3.5	294-114-60	40 g
Opaque Dentine A4	294-115-30	15 g
Opaque Dentine A4	294-115-60	40 g
Opaque Dentine B1	294-121-30	15 g
Opaque Dentine B1	294-121-60	40 g
Opaque Dentine B2	294-122-30	15 g
Opaque Dentine B2	294-122-60	40 g
Opaque Dentine B3	294-123-30	15 g
Opaque Dentine B3	294-123-60	40 g
Opaque Dentine B4	294-124-30	15 g
Opaque Dentine B4	294-124-60	40 g
Opaque Dentine C1	294-131-30	15 g
Opaque Dentine C1	294-131-60	40 g

OD Opaque Dentine

Opaque Dentine C2	294-132-30	15 g
Opaque Dentine C2	294-132-60	40 g
Opaque Dentine C3	294-133-30	15 g
Opaque Dentine C3	294-133-60	40 g
Opaque Dentine C4	294-134-30	15 g
Opaque Dentine C4	294-134-60	40 g
Opaque Dentine D2	294-142-30	15 g
Opaque Dentine D2	294-142-60	40 g
Opaque Dentine D3	294-143-30	15 g
Opaque Dentine D3	294-143-60	40 g
Opaque Dentine D4	294-144-30	15 g
Opaque Dentine D4	294-144-60	40 g

D Dentine

Dentine A1	294-211-30	15 g
Dentine A1	294-211-60	40 g
Dentine A2	294-212-30	15 g
Dentine A2	294-212-60	40 g
Dentine A3	294-213-30	15 g
Dentine A3	294-213-60	40 g
Dentine A3.5	294-214-30	15 g
Dentine A3.5	294-214-60	40 g
Dentine A4	294-215-30	15 g
Dentine A4	294-215-60	40 g
Dentine B1	294-221-30	15 g
Dentine B1	294-221-60	40 g
Dentine B2	294-222-30	15 g
Dentine B2	294-222-60	40 g
Dentine B3	294-223-30	15 g
Dentine B3	294-223-60	40 g
Dentine B4	294-224-30	15 g
Dentine B4	294-224-60	40 g
Dentine C1	294-231-30	15 g
Dentine C1	294-231-60	40 g
Dentine C2	294-232-30	15 g
Dentine C2	294-232-60	40 g
Dentine C3	294-233-30	15 g
Dentine C3	294-233-60	40 g
Dentine C4	294-234-30	15 g
Dentine C4	294-234-60	40 g
Dentine D2	294-242-30	15 g
Dentine D2	294-242-60	40 g
Dentine D3	294-243-30	15 g
Dentine D3	294-243-60	40 g
Dentine D4	294-244-30	15 g
Dentine D4	294-244-60	40 g

Bleach

Chroma Opaque CO 1 Bleach	292-301-30	15 g
Chroma Dentine CD 1 Bleach	294-251-30	15 g
Chroma Dentine CD 2 Bleach	294-252-30	15 g
Chroma Incisal CI 1 Bleach	295-701-30	15 g

CD Chroma Dentine

Chroma Dentine A	294-410-30	15 g
Chroma Dentine B	294-420-30	15 g
Chroma Dentine C	294-430-30	15 g

IFD Intensive Fluo Dentine

Intensive Fluo Dentine creme	294-751-30	15 g
Intensive Fluo Dentine yellow	294-760-30	15 g
Intensive Fluo Dentine orange	294-770-30	15 g
Intensive Fluo Dentine orange plus	294-773-30	15 g

IT Incisal transparent

Incisal 57	295-150-30	15 g
Incisal 57	295-150-60	40 g
Incisal 58	295-154-30	15 g
Incisal 58	295-154-60	40 g
Incisal 59	295-160-30	15 g
Incisal 59	295-160-60	40 g
Incisal 60	295-170-30	15 g
Incisal 60	295-170-60	40 g

IO Incisal opalescence

Incisal opal.1 ice white	295-210-30	15 g
Incisal opal.2 polar white	295-215-30	15 g
Incisal opal.3 snow white	295-250-30	15 g
Incisal opal.4 sky blue	295-270-30	15 g
Incisal opal.5 amber	295-280-30	15 g
Incisal opalescence 57	295-300-30	15 g
Incisal opalescence 58	295-305-30	15 g
Incisal opalescence 59	295-310-30	15 g
Incisal opalescence 60	295-315-30	15 g

NT Neutral transparent

Neutral transparent	295-299-30	15 g
Neutral transparent	295-299-60	40 g

TE Transpa Effect

Transpa Effect 1 white	295-350-30	15 g
Transpa Effect 2 yellow	295-360-30	15 g
Transpa Effect 3 pink	295-385-30	15 g
Transpa Effect 4 homey	295-370-30	15 g
Transpa Effect 5 blue	295-380-30	15 g
Transpa Effect 6 brown	295-367-30	15 g

ST Stains universal

Body Stains A	296-110-10	3 g
Body Stains B	296-120-10	3 g
Body Stains C	296-130-10	3 g
Stains 0 neutral	296-199-10	3 g
Stains 1 white	296-150-10	3 g
Stains 2 vanilla	296-162-10	3 g
Stains 3 yellow	296-160-10	3 g
Stains 45 orange plus	296-171-10	3 g
Stains 6 pink	296-185-10	3 g
Stains 7 dark pink	296-187-10	3 g
Stains 8 purple	296-189-10	3 g
Stains 9 blue	296-180-10	3 g
Stains 10 grey	296-155-10	3 g
Stains 11 olive green	296-193-10	3 g
Stains 12 olive yellow	296-192-10	3 g
Stains 13 medium brown	296-167-10	3 g
Stains 14 reddish brown	296-168-10	3 g
Stains 15 black	296-198-10	3 g

CM Correction

Correction material	297-199-10	3 g
---------------------	------------	-----

G Gingiva

Gingiva G (light)	298-110-30	15 g
Gingiva GD (dark)	298-120-30	15 g

Liquids

Stains universal	299-110-40	20ml
Shoulder material univ.	299-120-40	20ml
Modelling MV univ.	299-150-40	20ml
Modelling MV univ.	299-150-80	100ml
Modelling liquid LV standard univ.	299-160-41	20ml
Modelling liquid LV standard univ.	299-160-81	100ml
Modelling liquid LV standard univ.	299-160-91	500ml
Modelling liquid LV+ univ.	299-160-40	20ml
Modelling liquid LV+ univ.	299-160-80	100ml
Modelling liquid LV+ univ.	299-160-99	500ml
Universal paste liquid	299-170-40	20ml
Universal paste liquid	299-170-80	100ml
Universal B.O.L. liquid	299-180-40	20ml
Universal B.O.L. liquid	299-180-80	100ml
SM-Isofit	260-324-01	20ml
AESTHETIC Bonder liquid	299-200-10	5ml

5. Firing charts – for Titanium Framework

Austromat 3001	
Powder bonder bake	C500 T60 T180 · L9 V9 TO65 · C795 T60 V0 CO LO T2 C500
Paste bonder bake	C500 T60 T300 · L9 V9 TO65 · C795 T60 V0 CO LO T2 C500
Opaque bake 1 + 2	C500 T60 T180 · L9 V9 TO65 · C795 T60 V0 CO LO T2 C500
Shoulder bake	C500 T60 T360 · L9 V9 TO55 · C785 T60 V0 CO LO T2 C500
Dentin bake 1 (with liquid LV universal) **	C500 T60 T360 · L9 V9 TO55 · C755 T60 V0 CO LO T2 C500
Dentin bake 2 (with liquid LV universal) **	C500 T60 T240 · L9 V9 TO55 · C755 T60 V0 CO LO T2 C500
Glaze bake	C500 T60 T120 · L9 TO55 · C755 T60 CO LO T2 C500
Correction material	C500 T60 T240 · L9 V9 TO55 · C715 T60 V0 CO LO T2 C500
Gingival material	C500 T60 T360 · L9 V9 TO55 · C755 T60 V0 CO LO T2 C500

Vacumat 100/200/50						
	Final temperature	Predrying time (min.)	Heat rate (min.)	Holding time	Vacuum (min.)	Cooling
Powder bonder bake	795°C 1463°F	4	5	1 min. under vacuum	6	—
Paste bonder bake	795°C 1463°F	6	5	1 min. under vacuum	6	—
Opaque bake 1 + 2	795°C 1463°F	4	5	1 min. under vacuum	6	—
Shoulder bake	785°C 1445°F	6	5	1 min. under vacuum	6	—
Dentin bake 1 (with liquid LV universal) **	755°C 1391°F	6	5	1 min. under vacuum	6	—
Dentin bake 2 (with liquid LV universal) **	755°C 1391°F	4	5	1 min. under vacuum	6	—
Glaze bake	755°C 1391°F	2	5	1 min. ***	***	—
Correction material	715°C 1319°F	6	5	1 min. under vacuum	6	—
Gingival material	755°C 1391°F	6	5	1 min. under vacuum	6	—

Multimat MCII (Mach1/Mach2)								
	Preheating temperature	Drying (min.)	Preheat- ing (min.)	Vacuum (min.)	Firing time (min.)	Firing tempera- ture	Heat rate/min.	Vacuum
Powder bonder bake	500°C 932°F	3	1	1,0	2,0	815°C 1499°F	65°C 149°F	50 122
Paste bonder bake	500°C 932°F	5	1	1,0	2,0	815°C 1499°F	65°C 149°F	50 122
Opaque bake 1 + 2	500°C 932°F	3	1	1,0	2,0	815°C 1499°F	65°C 149°F	50 122
Shoulder bake	500°C 932°F	4	2	1,0	2,0	805°C 1481°F	55°C 131°F	50 122
Dentin bake 1 (with liquid LV universal) **	500°C 932°F	5	2	1,0	2,0	775°C 1427°F	55°C 131°F	50 122
Dentin bake 2 (with liquid LV universal) **	500°C 932°F	4	2	1,0	2,0	775°C 1427°F	55°C 131°F	50 122
Glaze bake	500°C 932°F	2	2	***	1,5	775°C 1427°F	55°C 131°F	***
Correction material	500°C 932°F	4	2	1,0	2,0	735°C 1355°F	55°C 131°F	50 122
Gingival material	500°C 932°F	4	2	1,0	2,0	775°C 1427°F	55°C 131°F	50 122

** When using modelling liquid MV universal, observe different bake sequence.

*** The glaze bake can be carried out under or without vacuum. The grade of glaze can be increased by prolonging the holding time.

6. Firing charts – for Zirconium Oxide Framework

Austromat 3001	
Liner bake 1 + 2	C500 T60 T180 · L9 V9 TO65 · C800 T60 V0 CO LO T2 C500
Shoulder bake	C500 T60 T360 · L9 V9 TO65 · C790 T60 V0 CO LO T2 C500
Dentin bake 1 (with liquid LV universal) **	C500 T60 T360 · L9 V9 TO55 · C760 T90-120 V0 CO LO T2 C500
Dentin bake 2 (with liquid LV universal) **	C500 T60 T240 · L9 V9 TO55 · C760 T90-120 V0 CO LO T2 C500
Glaze bake ***	C500 T60 T120 · L9 TO55 · C760 T60 CO LO T2 C500
	The holding time for dentin and correction bake should be between 90 and 120 seconds.


Vacumat 100/200/50						
	Final temperature	Pre-drying time (min.)	Heat rate (min.)	Holding time	Vacuum (min.)	Cooling
Liner bake 1 + 2	800°C 1472°F	4	5	1 min. under vacuum	6	—
Shoulder bake	790°C 1454°F	6	5	1 min. under vacuum	6	—
Dentin bake 1 (with liquid LV universal) **	760°C 1400°F	6	5	1,5-2 min. under vacuum	6,5-7,0	—
Dentin bake 2 (with liquid LV universal) **	760°C 1400°F	4	5	1,5-2 min. under vacuum	6,5-7,0	—
Glaze bake	760°C 1400°F	2	5	1 min. ***	***	—


Multimat MCII (Mach1/Mach2)								
	Pre-heating temperature	Drying (min.)	Pre- heating (min.)	Vacuum (min.)	Firing time (min.)	Firing temperature	Heat rate/min	Vacuum
Liner bake 1 + 2	500°C 932°F	3	1	1,0	2,0	820°C 1508°F	65°C 149°F	50 122
Shoulder bake	500°C 932°F	4	2	1,0	2,0	810°C 1490°F	55°C 131°F	50 122
Dentin bake 1 (with liquid LV universal) **	500°C 932°F	4	2	1,0	2,5-3,0	780°C 1436°F	55°C 131°F	50 122
Dentin bake 2 (with liquid LV universal) **	500°C 932°F	4	2	1,0	2,5-3,0	780°C 1436°F	55°C 131°F	50 122
Glaze bake	500°C 932°F	2	2	***	1,5	780°C 1436°F	55°C 131°F	***

Please note:

The values given are only guidelines and should be adjusted individually depending on the respective situation, age and manufacturer of the furnace.

The firing charts are based on newly fine-silver calibrated furnaces.

	291-000-04	1 piece	I Basic materials Triceram® Shade guide
	291-000-05	1 piece	II Effect materials Triceram® Shade guide
	291-000-06	1 piece	III Effect materials Triceram® Shade guide
	291-000-03	1 piece	Shade guide for Stains universal ø 35 mm
	260-315-00	1 piece	Shade sample former For easy preparation of bakesamples.
	260-302-00	1 piece	Glass slab for stains and materials 280 x 160 x 10 mm
	260-303-00	5 pieces	Fleece replacement for glass slab

Porcelain slab for stains 295 x 70 x 40 mm	260-304-00	1 piece	
SM-Isofit To smooth out and seal atthe same time. To isolate plaster from ceramics and ceramics from ceramics.	260-324-01	20 ml	
Teeth Light Small lamp for demonstration of light dispersion in the ceramic work.	260-323-00	1 piece	
Agate spatula	260-300-05	1 piece	
Portion spoon	260-300-06	1 piece	
Glass instrument for bonder and opaque	291-000-00	1 piece	
Firing tray with holders	260-300-03	1 piece	
Holders for firing tray	260-300-13	6 pieces	
Carbon fibre chips For purge of ceramic furnaces.	260-317-00	10 pieces	

Dentaurum Group

Germany | Benelux | España | France | Italia | Switzerland | Australia | Canada | USA
and in more than 130 countries of the world.



➔ For more information on our products and services, please visit www.dentaurum.de

Date of information: 07/13

Subject to modifications



D
DENTAURUM