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
**REF 302-602-00**

**Dear customer**

Thank you for choosing a quality product from Dentaurum.

It is essential to read these instructions carefully and adhere to them to ensure safe, efficient use and ensure that you and your patients gain full benefit.

Instructions for use cannot describe every eventuality and possible application. In case of questions or ideas, please contact your local representative.

As our products are regularly upgraded, we recommend that you always carefully read the current Instructions for use supplied with the product and stored in the Internet at [www.dentaurum.com](http://www.dentaurum.com), even though you may frequently use the same product.

#### 1. Manufacturer

Dentaurum GmbH & Co. KG · Turnstr. 31 · 75228 Ispringen · Germany

#### 2. General product description

The tomas®-RPE eyelet (Fig. 1) can be used to fabricate skeletally anchored appliances for rapid palatal expansion (RPE). The appliance can be held in place in the palate with two tomas®-pin (hybrid RPE appliance with occlusal support, Fig. 2) or four tomas®-pin (skeletal RPE appliance / MARPE with occlusal support, Fig. 3). The tomas®-RPE eyelet serves to couple the tomas®-pin directly with the RPE screw. Each pack contains two supports to aid placement of the tomas®-RPE eyelet on the working cast.

#### 3. Precautions

Before use, operators should ensure that they have carefully read and followed the Instructions for use supplied. We recommend operators complete a relevant tomas® training course prior to use, as the Instructions for use can only cover some of the many areas of application.

Only original components may be used with the tomas®-RPE eyelet according to the Instructions for use. The tomas®-RPE eyelet was especially developed for the tomas®-pin and is not suitable for use with other mini implants. The tomas®-RPE eyelet should be used preferably with the tomas®-pin EP. It is, however, also possible to use it with the tomas®-pin SD.

RPE appliances that are intended to be fixed in place skeletally with the help of the tomas®-RPE eyelet must have an occlusal support (preferably in the area of the first premolars, Fig. 2 and 3).

#### 4. Application

##### 4.1 Planning of the position of the tomas®-pin (dentist)

An impression should be taken of the maxilla and the mandible and models made to aid planning and fabrication of the appliance. The models and further diagnostics are used to plan the position of the tomas®-pin.

For the skeletally anchored RPE appliance to function safely, it is essential that the tomas®-pin is firmly anchored in the bone (primary stability). For a hybrid RPE appliance, the position of the two tomas®-pin EP should each be three to five millimeters paramedian at the level of the first premolars (Fig. 4). Four tomas®-pin EP are required for a purely skeletally anchored RPE appliance (MARPE). The two anterior tomas®-pin EP should be placed in the same way as for the hybrid RPE appliance. The two posterior tomas®-pin EP should each be placed three to five millimeters paramedian at the level of the second molars (Fig. 5). The position of the pins and their insertion path should be communicated to the dental technician in a suitable manner.

There are several aspects to consider when selecting the length of the tomas®-pin. The portion of the tomas®-pin in the bone must be at least as long as the portion of the pin not in the bone. It is important to consider not only the mucosal thickness in the direction of insertion, but also that, when using the tomas®-RPE eyelet, the conical neck of the tomas®-pin is not in the mucosa. This results in a height of at least 1.5 mm. A further aspect is that the tomas®-pin used for a MARPE should be anchored in the posterior section of the palate bicortically. For this reason, the anterior tomas®-pin EP should be 10 mm long and the posterior tomas®-pin EP should be at least 8 mm long.

##### 4.2 Fabrication of the RPE appliance (dental technician)

The tomas®-RPE eyelet has been optimized for application with the tomas®-pin and RPE screws with retention legs of 1.5 mm diameter. These retention legs fit exactly on the circumferential groove (Ø 1.5 mm) on the tomas®-RPE eyelet (Fig. 1). The active element of the RPE appliance is the RPE screw placed centrally in the palate. The RPE appliance should be fabricated in accordance with the manufacturer's instructions. The RPE appliance must have an occlusal support (in the area of the first premolars).

The marking of the position of the tomas®-pin serves as orientation for aligning the anterior retention legs of the RPE screw. These should be bent in the direction of the markings. The retention legs should preferably be fitted to the tomas®-RPE eyelet from the mesial direction. If this is not possible, they can also be fitted laterally.

A small hole (Ø 2.3 mm) is drilled into the model at the positions marked by the dentist for the anterior tomas®-pin. The tomas®-RPE eyelet is fitted on the support. The large opening (Ø 3.0 mm) should face the grip on the holder (Fig. 1 and 6). The holder with the tomas®-RPE eyelet is then adjusted to the direction of insertion as stipulated by the dentist and the holder is fixed in position with wax. Final adjustments are then made to the anterior retention legs of the RPE screw on the tomas®-RPE eyelet. The legs should grasp the circumference of the eyelet as much as possible. The retention legs of the RPE screw are welded to the tomas®-RPE eyelet. Solder should not be used for reasons of biocompatibility and stability.

The posterior retention legs of the RPEs®-pin are connected to the first molars. The skeletally anchored RPE appliance with the tomas®-RPE eyelet must have occlusal support – preferably in the area of the first premolars (Fig. 7 and Fig. 8). Should the RPE appliance be fixed on four tomas®-pin, the posterior tomas®-RPE eyelet should be placed in exactly the same way as in the anterior region. A stable steel wire (Ø 1.5 mm) should be adjusted and welded between the posterior tomas®-RPE eyelet and the posterior retention legs or the body of the RPE screw (Fig. 8). Wires should definitely not be welded to the body of the RPE screw as this may impair the screw function.

##### 4.3 Fitting of the RPE appliance (dentist)

First, it should be checked if the RPE appliance fits well on the model and in the mouth. Small corrections should be made if necessary. This is the case especially for the direction of insertion of the tomas®-pin. Small alterations can still be made if the tomas®-RPE eyelet should need realigning. The tomas®-RPE eyelet should not exert any pressure on the mucosa.

Once the point of insertion for the tomas®-pin has been anesthetized, the thickness of the mucosa in the direction of insertion is measured. The RPE appliance can be fitted and fixed.

Please observe the Instructions for use for the insertion of the tomas®-pin (REF 989-534-00). To fix the RPE appliance skeletally, screw the tomas®-pin through the tomas®-RPE eyelet. Push the tip of the tomas®-pin through the middle of the tomas®-RPE eyelet into the mucosa until the surface of the bone can be felt. Then check that the direction of insertion is correct and begin screwing the pin in place. Stop screwing the pin when the insertion instrument comes into contact with the tomas®-RPE eyelet (Fig. 9).

##### 4.4 Therapy with the RPE appliance

Once the tomas®-pin has been inserted, expansion can begin using the RPE screw. Intervals for opening the RPE screw depend on the individual case. Once the active opening phase of the screw is completed, the appliance should be left in the mouth for a suitable amount of time to retain the results achieved. The length of the retention phase depends also on the individual patient.

##### 4.5 Removal of the RPE appliance

The RPE appliance can be removed after a suitable period of retention. First, the tomas®-pin are removed, see Instructions for use (REF 989-534-00). The RPE appliance is then removed. .

#### 5. Composition

The composition is included in the materials list; please refer to the catalog or [www.dentaurum.com](http://www.dentaurum.com).

#### 6. Safety instructions

The product should not be used if there is a known allergic reaction to one or more of the material components.

#### 7. Information for single use products

The tomas®-RPE eyelet is designed for single use only. Reconditioning of a tomas®-RPE eyelet that has already been used (recycling) and its reuse on a patient are not permitted.

#### 8. Further information

Should the dental professional and/or the patient become aware of serious problems arising from the use of the product, it is important that the manufacturer and the competent authority in the country in which the dental professional and/or the patient is resident are informed accordingly.

#### 9. Quality information

Dentaurum ensures faultless quality of its products. These recommendations are based upon Dentaurum's own experiences. The user is solely responsible for the processing of the products. Responsibility for failures cannot be taken, as Dentaurum has no influence on the processing on site.

#### 10. Explanation of symbols used on the label

 Refer to the label. Additional information can be found at [www.dentaurum.com](http://www.dentaurum.com) (Explanation of symbols REF 989-313-00).

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