REF 094-782-00

Gebrauchsanweisung	S. 2
Instructions for use	P. 20
Mode d'emploi	P. 39
Modo de empleo	P. 58

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Instructions for use

Contents

1	Declaration of Conformity	22
2	Safety Notices	23
2.1	Intended Use	23
2.2	Warning Symbols	23
2.3	Warranty and Liability	23
2.4	Owner Responsibility	23
2.5	Operating Safety Notices	24
2.6	Disposal	24
2.7	Delivery	24
3	Application and Description	25
3.1	Operation	25
4	Installing the Appliance	25
4.1	Unboxing and Transport	25
4.2	Where to Place the Appliance	25
4.3	Connecting to the Mains	25
4.4	First-Time Operation	25
5	Using the Dublitherm compact	26
5.1	Turning on the Appliance	26
5.2	Menu Overview and Operating Method	27
5.3	Lid Locking Mechanism	28
5.4	Filling the Unit	28
5.5	Notes on Working with Duplicating Gel	29
5.6	Starting the Melting Process	30
5.7	Pouring the Melted Gel	31
5.8	Cancelling the Melting Process	31
5.9	Setting the Melting Temperature	31
5.10	Setting the Pouring Temperature	32
5.11	Setting the Temperature Display °C or °F	33
5.12	Setting the Timer Function	34
5.13	Activating the Timer	35
5.14	Stopping the Timer	35
6	Care and Maintenance	36
6.1	Care and Maintenance of the Mixing paddles	36
6.2	Cleaning the Appliance	36
6.3	Replacing the Fuse	36
7	Troubleshooting	37
8	Replacement Parts	37
9	Block Diagram	38
10	Technical Information	38

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hereby declares that the design and construction of the laboratory equipment described below, including the version marketed by us, comply with the basic regulations governing safety and health as stated in the EEC Guidelines. This declaration will become invalid if the laboratory equipment is modified or altered in any way without our prior consent.

Description of unit:	Dublitherm Compact	
	Duplicating machine	
	REF 094-782-00	
Start with Unit No.:	0041 - 2011	
EEC guidelines:	2006/42/EC	Guideline for machine
	2006/95/EC	Electrical equipment used within certain voltage limits
	2004/108/EC	Electromagnetic compatibility
Applied unified standards:	EN 61010-1:2001 EN 61010-2-051:2003 EN 61000-3-2:2000 EN 61326-1:2006	

Date and manufacturers signature: 04.04.2011
Signatory:

- i.V. Dipl.Ing. (FH) K. Merkle -
Production Manager Mechanic

Printing Date: 18.04.11

2 Safety Notices

 **Caution!** Please read carefully the instructions for use before installing and using the appliance. Do not turn on the appliance until you have read the instructions!

2.1 Intended Use

The **Dublitherm compact** is designed for the smooth melting of dental duplicating gels in the laboratory. Using this unit any other way than specified is regarded as misuse. Dentaureum GmbH & Co. KG is not liable for damage resulting from misuse. Intended use includes adherence to these instructions and providing regular inspection and maintenance.

Final deactivation and disposal of any Dentaureum product must be in accordance to specific regulations of the country of residence. Please direct any enquiries about correct disposal to Dentaureum or the regional dental depot.

2.2 Warning Symbols

The following warning symbols are used within these instructions:

 **Warning** Indicates potential danger to life and health. Ignoring this sign can critically affect a person's health and may result in life-threatening injuries.

 **Caution** Indicates a potentially dangerous situation. Ignoring this sign can result in minor injury or damage to the equipment.

 **Note** Indicates general appliance information.

2.3 Warranty and Liability

Our general conditions of sale and delivery apply. Warranty and liability claims for personal or mechanical damages are invalid when found to have resulted from any of the following:

- Improper installation, operation, assembly or maintenance
- Misuse
- Operation with faulty safety equipment, improperly installed or inoperable safety and protective precautions
- Non-compliance with the instructions for use
- Inadequate monitoring of working parts
- Unauthorised structural changes to the appliance
- Improper repairs

2.4 Owner Responsibility

The owner is responsible for ensuring that only employees who

- have been trained to use the appliance and are familiar with the safety precautions
- have read and understood the safety notices and instructions for use
- have been instructed in the relevant accident precautions

operate the appliance.

2.5 Operating Safety Notices

- Before using the appliance, all users must first read the basic health and safety regulations, operating safety notices, and the instructions for use.
- Keep the workspace clean and tidy. An untidy workspace increases the risk of accidents.
- For your own safety, use only materials and accessories recommended in the instructions for use.
- The intended safety precautions may be impaired if the appliance is not used according to the manufacturer's instructions.
- Always ensure the appliance is switched off and disconnected from the mains before maintenance.

Warning Keep hands clear of rotating mixing apparatus. Risk of injury.

Caution Hot surface! The appliance may get hot, particularly on the top surface and around the melting pot. Touching surfaces where this symbol is seen may result in burns.

2.6 Disposal

i Note

This symbol indicates the product may not be disposed of as normal domestic waste. Current law forbids commercial customers from disposing of electronic waste in communal collection areas. Further information is available from Dentaaurum or your dental depot.

2.7 Delivery

1	Duplicating Unit Dublitherm compact – 230 V	REF 094-782-00
1	Instructions for use including Declaration of Conformity	
1	Mains Cable	REF 907-027-10
1	Spare Fuse 4 A TH 5x20	REF 907-681-20

3 Application and Description

3.1 Operation

The **Dublitherm compact** duplicating unit is a universal dental melting appliance. It can be used to melt all conventional duplicating materials.

The built-in industrial motor is very powerful.

Other **Dublitherm compact** characteristics:

- Compact unit with 2.5 kg duplicating gel capacity
- Fast heating and cooling times
- Gentle processing of materials
- Will remelt duplicating gel set in the melting pot
- Special combi-paddles
- Low energy consumption
- Exact temperature control
- Easy to use
- Pre-select start time delay between 30 minutes and 72 hours

4 Installing the Appliance

4.1 Unboxing and Transport

To remove the **Dublitherm compact** from the packaging or move it around the laboratory grasp the underside. Do not use the lid to lift the unit out of the box.

4.2 Where to Place the Appliance

Position the unit in a clean and dust free area. Place the unit so the main switch on the back is accessible at all times.

Caution Do not place furnaces or other heat-emitting appliances next to the **Dublitherm compact**. Do not cover the air vents.

4.3 Connecting to the Mains

The voltage marked on the rating plate and the mains supply must be identical.

- Connect the unit to a 230 V 50/ 60 Hz socket
- The socket must be fused.

Caution All maintenance to power supply sockets and cables must be performed by qualified electricians.

4.4 First-Time Operation

After correct unboxing and installation, the **Dublitherm compact** can be turned on for the first time.

5 Using the Dublitherm compact


5.1 Turning on the Appliance

Press the main switch at the back of the appliance to turn on the power. The unit will now run a self-test. The software version is shown and all LEDs light up at once.

Main switch on the back of the unit.

The appliance is now running and ready for use. The display now shows **Stby** (Standby)

① Heat

② Cool

③ Pour

④ Timer

5.2 Menu Overview and Operating Method

Once switched on the appliance is in Operating Mode.

Press the „Menu“ button once to change features such as timer, temperature unit, melting and pouring temperature.

You are now in Parameter Mode.


In this mode, the chosen parameter can be dialled up using the arrow keys $\uparrow$ or $\downarrow$.

Once you have chosen the parameter to change, press the „Menu“ button once more to enter Editing Mode. Change the value displayed using the arrow keys $\uparrow$ or $\downarrow$. Press the „Menu“ button again and the edited value is saved to memory.

To exit the Edit or Parameter Modes and return to Operating Mode without saving, simply press and hold the „Menu“ button for longer than 3 seconds.

i Note

Section 5.6 and 5.11 cover these programming steps in more detail.

5.3 Lid Locking Mechanism

The unit's lid features a locking mechanism. This prevents accidental opening during use or duplicating gel forcing the lid open.

The locking mechanism activates automatically when the lid is closed. To open the lid, pull the locking knob outwards.

By pulling the knob outwards and turning it 90° to the right, the lock is latched and the lid can be opened without prior unlocking.

Pulling out the locking mechanism

Latching the knob in place - lid is permanently open

5.4 Filling the Unit

The unit should be filled with small pieces of duplicating gel for further processing.

The combi-paddles and powerful motor enable the unit to process even cooled and set gel still in the tank.

Note Please observe the maximum capacity of the tank to ensure the best possible results.

Melting chopped gel

To fill the tank, first add 1000 g finely chopped gel. This equates to a filling when the top mixing paddle is still just visible (note the groove). Close the lid tightly.

After starting the melting process (see Section 5.5) the gel is heated to 40 °C. Once this temperature is reached, the mixing paddles start and the level of gel will fall. Chopped gel may now be added to the tank. The tank holds a maximum of 2500 g of gel (indicated by the bottom edge of the uppermost mixing paddle).

Fill level: The groove is situated on the lower mixing paddle only and marks the maximum fill level.

Melting set gel:

The mixing paddles start when the unit reaches 40 °C. The tank will hold no more than 2500 g of gel, as stated above (see groove on the lower mixing paddle).

Caution

Avoid overfilling. Maximum first load = 1000 g chopped gel
Maximum amount of set gel = 2500 g
Observe groove.

5.5 Notes on Working with Duplicating Gel

When first filling the tank with gel such as Dubliform, Duplikat, Dublitol® or Dubliplast®, chop the gel into small pieces. This increases filling capacity and shortens melting time.

Apart from duplicating gels, duplicating concentrates may also be used.


Note

Observe mixing ratio!

Gel kept melted for longer periods must be checked regularly for moisture loss. Refill with small amounts of water.

Replace used gel completely with new material. Do not mix old and new gel. Clean the tank regularly to prevent build up of plaster or investment residue.

5.6 Starting the Melting Process



Press the Start/Timer button once to start the melting process.

Ensure the lid is closed properly.

First, the gel is heated to 40 °C, the LED labelled „Heating“ lights red, and the display shows the current temperature in the tank.

When the temperature reaches 40 °C, the mixing paddles start automatically.

Note Press the Menu button during the melting process to display both the pre-set melting and pouring temperatures for 3 seconds each. After a further 3 seconds, the display returns to show the current temperature.

Once the pre-set melting temperature is reached, the „Heating“ LED turns off and the „Cooling“ LED lights blue. The heating elements are turned off and a ventilator cools the melting tank by a cooling channel.



5.7 Pouring the Melted Gel

Pull the outlet valve to release the gel. Any gel left in the valve is pressed back into the tank on closing.

5.8 Cancelling the Melting Process

To stop or cancel the melting process, press and hold the Stop/Reset button for 3 seconds. The unit returns to Operating Mode and the display shows **Stby** (Standby).


5.9 Setting the Melting Temperature

Press Menu once to enter Parameter Mode.

The display shows the last selected parameter. Navigate to the **HEAT** parameter using the arrow buttons **↑** or **↓**.

Press Menu once more to edit the temperature value using the arrow keys **↑** and **↓**. Choose any temperature between 90...95 °C.



Confirm the selected temperature and exit Editing Mode by pressing Menu once more. Pressing Menu for longer than 3 seconds returns you to Operating Mode **Stby**.

i Note

Should several parameters need editing, press Menu just once. This confirms the entered value but keeps you in Editing Mode.


5.10 Setting the Pouring Temperature

Press Menu once to enter Parameter Mode.



The display shows the last selected parameter (for example, **HEAT** for Melting temperature). Navigate to the **Hold** parameter with the arrow buttons **↑** or **↓**.

Press Menu once more to edit the temperature value with the arrow keys **↑** or **↓**. Choose any temperature between 45...60 °C.



Confirm the selected temperature and exit editing mode by pressing Menu once more. You are now in Parameter Mode. Pressing Menu for longer than 3 seconds returns you to Operating Mode **5t64**.

5.11 Setting the Temperature Display °C oder °F

Press Menu once to enter Parameter Mode. The display shows the last selected parameter (for example, *HEAL* for Melting temperature). Navigate to the *d.5P* parameter with the arrow buttons $\uparrow$ or $\downarrow$.

Press Menu once more to edit the temperature display with the arrow keys $\uparrow$ and $\downarrow$. Choose either °C (Celsius) or °F (Fahrenheit).



Confirm the selected temperature unit and exit Editing Mode by pressing Menu once more. You are now in Parameter Mode. Pressing Menu for longer than 3 seconds returns you to Operating Mode *5664*.

5.12 Setting the Timer Function


The Dublitherm compact has an integrated timer and can provide melted gel at any predetermined time.

The delay-time can be set between 30 minutes (the display shows 0.30) and 72 hours (the display shows 72.00) in 15-minute stages.

Press Menu once to enter Parameter Mode. The display shows the last selected parameter (for example, **HEAT** for melting temperature). Navigate to the **Endn** parameter with the arrow buttons **↑** or **↓**.

Press Menu once more to edit the delay-time value with the arrow keys **↑** and **↓**. Choose a delay-time from 30 minutes to 72 hours.

For example, the display shows 15 hours.

Confirm the selected delay time and exit Editing Mode by pressing Menu once more. You are now in Parameter Mode. Pressing Menu for longer than 3 seconds returns you to Operating Mode **Stby**.


5.13 Activating the Timer

Once the timer is set (see Section 5.11), fill the tank, close the lid, and ensure the unit is connected to the mains and switched on.

Activate the timer by pressing the Start/ Timer button for at least 3 seconds.

The „Timer“ LED now flashes in time with the dot on the time display. The timer counts down.

Note Should several parameters need editing, press Menu just once. This confirms the entered value but keeps you in Editing Mode.

The melting process begins automatically when the count down stops.

5.14 Stopping the Timer

To stop the countdown, press and hold the Stop/ Reset button for longer than 3 seconds. You are now in Operating Mode and the display shows **5664**.

6 Care and Maintenance

- ⚠ Warning** Always switch off the appliance and disconnect from the mains supply when servicing or during maintenance. Never work alone!
Always ensure a second person who is informed about the effects of electricity is at hand during service and maintenance work.
- ⚠ Caution** The mixing paddles have sharp edges. Danger of injury. Danger of cuts through contact. Do not reach into moving parts. Danger of injury.

6.1 Care and Maintenance of the Mixing paddles

The mixing paddles can be removed for cleaning. To remove, turn the top knurled nut counter-clockwise and pull the mixing paddles upwards. Wash the mixing paddles in warm water.

6.2 Cleaning the Appliance

Clean the appliance regularly with a damp cloth. Do not use aggressive cleaning agents or sharp objects for cleaning.

- ⚠ Caution** Never clean with compressed air, water, or solvents.

6.3 Replacing the Fuse

To replace a fuse, first switch off the appliance and disconnect from the mains supply. Remove the fuse box from the back panel and replace the faulty fuse with a new 4.0 A TH fuse.

Main switch with fuse box on the rear panel.

7 Troubleshooting

The following errors may occur:

Error	Cause	Corrective Measures
- Machine does not mix - No LED display	(1) appliance is not plugged in at the mains (2) appliance is not switched on (3) faulty fuse	(1) plug in at the mains (2) switch the machine on (3) change the fuse (see Section 6.3!)
- Appliance not mixing - No LED lighting and no display in the control pad.	(4) Overheating, power supply cut	(4) allow the appliance to cool down
- Appliance will not start	Lid not closed properly	Close lid

Error Message	Cause	Orrective Measures
CO	Lid not closed properly	Close lid
E-OL	Motor overloaded	Turn off appliance and allow to cool down. Restart the appliance.
E-SE	Temperature sensor fault	Contact Dentaaurum Service Personell


If you are unable to correct a fault on the basis of this guide, please contact your service technician or contact Dentaaurum Service direct:

Tel: +49 72 31/803-211 · e-mail: angelika.stoehr@dentaaurum.de

8 Replacement Parts

Mains Cable	REF 907-027-10
Fuse, 4.0 A TH	REF 907-681-20
Teflon housing	REF 971-503-00
Control PCB	REF 907-681-05
Mixing paddle, Front	REF 907-696-01
Mixing paddle, Rear	REF 907-696-10
Knurled Nut for Mixing paddles	REF 907-696-60

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