


BE HAPPY!
JUST SMILE
SO FAST &
GENTLE.

discovery[®]
SMART sl



IT'S EASY #CLICK&GO

D
DENTAURUM
1886


It's great that you've taken the first steps towards your perfect smile with our help. Our new, innovative brackets will help make treatment as pleasant and as short as possible.

What is a bracket?

Brackets are small squares that are bonded directly onto your teeth. They are used to move your teeth to the right position with the help of an arch.


Standard bracket:

The arch is fixed to the bracket using a small elastic ring or a wire – the so-called ligature – to keep it in place. The ligature is regularly changed along with the arch.



Self-ligating discovery® smart sl bracket:

These brackets have an integrated Click&Go clip which holds the arch firmly in place. The teeth can move more quickly and easily into the correct position and it's easier to keep the brackets clean.



GENTLE,
STABLE &
EFFICIENT
#CLICK&GO

FAST AND
RELIABLE
#CLICK&GO





ALL FEAT
MANY BE
SMILE, BR
CLICK&

SMART TREATMENT. QUICK RESULTS.

discovery® smart sl is our latest smart bracket from the successful discovery® premium line. It was especially developed to make your treatment as fast and as pleasant as possible. Your brackets were made for you at the Dentaurem headquarters in Ispringen in the heart of the Northern Black Forest in Germany.



Find out more
online!



URES
NEFITS
RACKETS
&GO



You can find tips
on caring for your
brackets here:



YOUR ADVANTAGES AT A GLANCE

○ Shorter treatment

Your teeth will be moved efficiently so they reach the right position quickly.

○ Fewer recalls and less time at the dentist's

- There are no elastics to change so you have fewer recalls.
- Arches can be changed quickly thanks to the click mechanism – you spend less time at the dentist's. Simply **#CLICK&GO**.

○ More comfortable to wear

- Tooth movement is gentle and efficient as friction is reduced and movement forces are controlled.
- The brackets are small and flat so they are comfortable to wear and cause little bother.

○ High levels of material tolerability

The base of the bracket is laser-structured and formed to fit on the tooth structure. This makes the material more tolerable and also gives the brackets a secure hold and stability.

○ Easier oral hygiene

As there are no wires or elastics, fewer bits of food get caught so cleaning is easier for you.

PACKED HYGIENICALLY

Your brackets are cleaned thoroughly and then packed hygienically into blister packaging. They're not opened at the orthodontist's until you get them.



Your dental office:



Dentaaurum GmbH & Co. KG



Turnstr. 31
75228 Ispringen · Germany



www.dentaaurum.com

D
DENTAURUM
1886