



Orthodontics Catalog



Edition 24

Orthodontics Catalog

Digital orthodontics	
Lingual system	
Self-ligating brackets	
Metal brackets	
Aesthetic brackets	
Nickel-free products	
Bracket accessories	
Buccal tubes	
Bands	
Adhesives	
Wires Arches Accessories	
Class II appliances	
Intra- and extra-oral	
Anchorage systems	
Supplies	
Expansion screws	
Orthodontic acrylics	
Pliers Instruments	
Equipment	
Knowledge Marketing	
Courses, webinars and events	
Service	

*"What we do today
determines how the
world shall go".*

Boris Pasternak

Customer number:

Please enter your Dentaaurum customer number here. For faster service, we kindly ask you to use your customer number when ordering.

Dear Customer,

Here at Dentaaurum, the oldest dental company in the world, it has been our aim since 1886 to develop and produce products of the highest quality and reliability. We have set strict standards to ensure we offer products with excellent characteristics and a high level of biocompatibility. Our focus is always firmly placed on our customers, their demands and their satisfaction. We base our decisions and our production processes on aspects that are of central importance to us: conserving resources, reducing the burden on the environment and the sourcing of raw materials in an ethically correct manner. This is the path we have taken since as far back as 1989; our environmental management system has been certified since 1996.

We have succeeded in setting new standards worldwide for many Dentaaurum products. Our affiliates in Benelux, France, Italy, Spain, Switzerland, USA, Canada and Australia and a global network of distributors sell our high-quality products in more than 130 countries and offer a first-class service locally.

We offer events and training worldwide, in person and via webinars, to keep you up to date with our products as well as developments and findings in dental orthodontics.

Since March 1, 2022, there are three members in our management team. This means three-fold expertise and experience in a variety of business segments, ensuring that Dentaaurum will continue its market success.

We are convinced that the key to the company's future lies in the success of our digital products and have developed many new products in this field in the past years. Alongside prime4me® Aligner, the innovative 2-splint treatment concept, the prime4me® RETAIN3R has also become an important part of our digital product range. We are looking forward to continuing along the path of digitization together with you. We will, however, not be neglecting our analog products as they remain indispensable despite advances in digitization.

Best regards from Ispringen,

Mark Stephen Pace

Claudia Stöhrle

Ralph Dittes

Dentaaurum Management

Dentaurum

Germany | Benelux | España | France | Italia | Switzerland | Australia | Canada | USA
and in more than 130 countries worldwide.



Dental technologies that set standards.

Dentaurum develops, produces and sells dental products and services worldwide. Dentaurum has played an active role in the dental sector for more than 137 years as an owner-managed family company with its headquarters in Ispringen (near Pforzheim, northern Black Forest), Germany.

Experience Dentaurum by taking a virtual tour in a 360° perspective. There is a total of five 360° viewpoints which give you a fantastic overview.





Orthodontics

Implantology

Prosthetics

Ceramics



Photo © Christian Frennig®

High-tech in-house.



There is a high level of professional knowledge in our company that has been built up over the years in our own research and development facilities, both in Germany and in France. Highly qualified employees work together in interdisciplinary teams to find answers to the challenges the future poses. At the same time, long-standing co-operations with experts from universities and clinics contribute to finding new developments and innovations.

A further result of these efforts: a comprehensive product portfolio which is one of Dentaureum's strengths. No other dental company has such an extensive range of products offering a total of more than 8500 articles.





Dentaurum Online Shop.

Order easily online.




shop.dentaurum.com

Your advantages:

- Detailed product information.
- Product availability is automatically checked.
- Individualized list of articles frequently ordered.
- Your terms and prices are displayed.
- Secure internet connection with SSL encryption.
- Determine the delivery date of your order.
- Up-to-date information on our product portfolio and respective CE marking.

- | | |
|---|-----------------------|
| ■ shop.dentaurum.de | ➤ Germany and Austria |
| ■ shop.dentaurum.ch | ➤ Switzerland |
| ■ shop.dentaurum.be | ➤ Benelux |
| ■ shop.dentaurum.es | ➤ Spain |
| ■ shop.dentaurum.fr | ➤ France |
| ■ shop.dentaurum.it | ➤ Italy |
| ■ shop.dentaurum.us | ➤ USA |
| ■ shop.dentaurum.ca | ➤ Canada |
| ■ shop.dentaurum.com.au | ➤ Australia |



DENTAURUM
DIGITAL

Innovative order platform for
patient-specific products and services.



digital.dentaurum.de



Order products and services for the individual patient quickly and easily. Products are manufactured to a high level of quality thanks to professional quality management and validated manufacturing processes as well as industry standard for machines and processes.

Products available:

- **prime4me® Aligner:** Orthodontic aligner therapy from Dentaureum.
- **prime4me® model printing:** Base models and dental arch models for orthodontics.
- **prime4me® RETAIN3R:** Unique fit using a three-dimensional construction design.

Your advantages:

- Easy, intuitive ordering for all Dentaureum Digital products.
- Easy file upload.
- Display of current order status.
- Easy approval of prime4me® RETAIN3R design and prime4me® Aligner treatment proposal.
- Filing of all open and completed orders.
- Online order-related communication with an experienced support team via the order platform.

Availability may vary depending on the country. Please ask your local representative.

Your direct contact.



Chat with us.



Hello, we can answer your questions quickly and easily.



Click on the chat symbol on our website, and send us your questions. We're happy to chat with you.



Give it a try.
We look forward to hearing from you!

D
DENTAURUM

Telephone	+4972 31/803-0
Fax	+4972 31/803-295
Email	info@dentaurum.com

D
DENTAURUM
AUSTRALIA

Telephone	+61 2 97 43 48 11
Fax	+61 2 97 43 45 99
Email	info@dentaurum.com.au

D
DENTAURUM
CANADA

Telephone	+1-800-363-9973
Fax	+1-888-771-3640
Email	info@dentaurum.ca

D
DENTAURUM
USA

Telephone	+1-800-523-3946
Fax	+1-800-553-6389
Email	info@dentaurum-us.com



Contents

	Digital orthodontics	16–41
	prime4me® Aligner	26
	prime4me® model print	26
	prime4me® RETAIN3R	27
	Expansion screws without retention legs for digital workflow	28–37
	OnyxCeph ^{3TM}	38–40
	Lingual system	42–51
	discovery® delight lingual brackets/tubes	46–47
	Arches	48
	Instruments	49
	Accessories	50
	Self-ligating brackets	52–59
	discovery® sl	56
	discovery® sl 2.0	57
	Instruments and accessories	58
	Metal brackets	60–97
	discovery®	68–73
	discovery® smart	74–77
	topic	78–80
	equilibrium® mini	81–85
	equilibrium® 2	86–94
	equilibrium® ti	95–96
	Aesthetic brackets	98–109
	discovery® pearl ceramic brackets	104–108
	Accessories	109
	Nickel-free products	110–113
	Bracket accessories	114–115
	Buccal tubes	116–135
	Bonding technique	120–125
	Banding technique	126–133
	Prewelded parts	134–135
	Bands	136–171
	dentaform® Snap	142–146
	dentaform®	147–156
	Standard	157–167
	Accessories	168–170


	Adhesives	172–181
	Adhesives and accessories	178–181
	Wires Arches Accessories	182–237
	Delivery program for arches / Straight wires / Wires on coils	192–199
	Fixed appliance technique / Arches / Straight wires / Wires on coils / Retainers / Accessories	200–226
	Removable appliance technique / Straight wires / Wires on coils / Accessories	227–235
	Palatal and lingual arch systems	236–237
	Class II appliances	238–255
	SUS ³ – Sabbagh Universal Spring	244–245
	Herbst bite jumping hinges / Sets and accessories	246–250
	amda® / Starter set, appliance and accessories	251–253
	tomas® / Pins	254
	Intra- and extra-oral	256–271
	Extra-oral appliances	260–262
	Accessories / Extra-oral	262–263
	Elastics and auxiliaries / Intra-oral	264–267
	Oral screens	268
	Anchorage systems	272–293
	tomas® / Tool Sets	278
	tomas® / Pins	279–280
	tomas® / Instruments and accessories	281–289
	amda® / Starter set	290
	amda® / Appliance and accessories	291–292
	Patient consultation material	293
	Supplies	294–309
	Impression taking	298
	Model fabrication	299–300
	Appliance containers and accessories	301–302
	Orthodontic study and demonstration models	303–307
	Other	308
	Expansion screws	310–341
	Standard expansion screws	314–316
	Special expansion screws	317–329
	Variety expansion screws	330–333
	RPE screws	334–339
	Accessories	340–341

	Orthodontic acrylics	342–367
	Orthocryl® LC / Starter sets / Single packs / Accessories	348–349
	Orthocryl® / Assortments	350
	Orthocryl® / Single packs	351–352
	Disco glitter	353
	Decals	354–357
	Accessories	358–367
	Pliers Instruments	368–401
	Premium-Line pliers	372–383
	EQ-Line pliers	386–389
	Instruments	390–397
	Accessories	398–400
	Equipment	402–425
	Equipment and accessories	410–424
	Knowledge Marketing	426–439
	Books / Apps / Downloads	428–439
	Courses, webinars and events	440–445
	CDC – Centrum Dentale Communicatie / Courses	442–445
	Service	446–519
	Contact	11
	Bracket systems – information	448
	Pre-welding service – banding technique	450–461
	Bands – comparison charts	462–469
	Wires – information	470–471
	Materials list	472–479
	Conversion charts	480
	Environment and sustainability	482–483
	Alphabetical index	484–492
	Numerical index	494–514
	Disclaimer / Image credits	519

prime4me® Aligner	Page	26
prime4me® model print	Page	26
prime4me® RETAIN3R	Page	27
Expansion screws without retention legs, for digital workflow	Page	28
OnyxCeph ³ ™	Page	38




Digital orthodontics.



DENTAURUM
**ONLINE
SHOP**
shop.dentaurum.com





DENTAURUM
DIGITAL


Innovative order platform for
patient-specific products and services.





passion and precision.





DENTAURUM
DIGITAL

Innovative order platform for
patient-specific products and services.

2 splint system:

Innovative splint concept for the best
treatment results.

» **Effective | Precise | Aesthetic | Competent**







Orthodontic aligner therapy
from Dentaaurum.

Advantages:

- ” Digital case planning with OnyxCeph³™.
- ” Aligner covers part of the gingiva for better transfer of forces.
- ” Transparent aligner material for almost invisible therapy.
- ” Personal consultation from an experienced support team.







DENTAURUM
DIGITAL

Innovative order platform for
patient-specific products and services.



prime4me® model print.

Printed base models and dental arch models for orthodontics are produced with maximum precision, providing an ideal basis for the fabrication of orthodontic appliances. As always, Dentaureum stays true to its own demands on quality: top industry standard for machines and validated processes together with highly qualified personnel.






prime4me® RETAIN3R offers numerous benefits.

The retainer is designed according to the customer's wishes using a unique three-dimensional construction design, guaranteeing a perfect fit. The digital design also works for situations where space is limited. The material, titanium grade 5, is biocompatible, making the prime4me® RETAIN3R highly suitable for patients with a nickel allergy.

Once the design has been approved, the retainer is milled one-to-one from a titanium blank. This ensures the patient will find the retainer comfortable to wear as its lingual fit is ideal and the adhesive surfaces are smaller. This, in turn, also helps improve and ease oral hygiene.







Dental imaging software.

With OnyxCeph^{3TM}, Image Instruments set the benchmark for dental imaging software. The comprehensive solutions for archiving, diagnostics, treatment planning and patient education are used by customers in more than 99 countries around the world – with an impressive list of successes.


Treatment planning and areas of application:




Module Ortho Apps 3D




Module FA_Bonding 3D




Module Kylix 3D



Module Bonding Trays 3D



Module Retainer 3D



Module Aligner 3D

OnyxCeph^{3TM} – Program versions and optional modules for the ideal combination with laboratory scanners or TRIOS.

	Program version*	
	3D PRO REF 410-004-00	3D LAB REF 410-005-00
SQL database		
Network licence (concurrent clients server)	5 1	20 1
Stand-alone licence (client server)	1 1	1 1
Imaging Module 2D		
Add image (incl. interfaces to digital X-ray)	■	■
Adjust image (classify and align image)	■	■
Clone finding	■	■
Edit	■	■
Digitize (cephalometric tracing)	■	■
Combine	■	■
Mirror	■	■
CO CR Correction	■	■
Treatment Simulation (orthodontics, maxillofacial surgery)	■	■
Ricketts V.T.O.	■	■
Image Comparison	■	■
Copy Save Send Show Print	■	■
Imaging Module 3D		
Add image (incl. DICOM extraction)	■	■
Adjust image (classify and align image)	■	■
Clone finding	■	■
Digitize (cephalometric tracing)	■	■
Combine	■	■
Mirror	■	■
Treatment Simulation	□	□
Cast Adjust	■	■
Segmentation	■	■
FA_Bonding	■	■
V.T.O. 3D (virtual set-up)	□	■
Wire_Bonding 3D	□	□
Kylix 3D (indirect transfer)	■	■
Bonding Trays 3D (indirect transfer)	□	□
Bonding Jigs 3D (indirect transfer)	■	■
Bracket Erase 3D	□	□
Aligner 3D	□	□
Retainer 3D	□	□
Bite Splint 3D	□	□
Ortho Apps 3D	□	□
Waefer Creation 3D	□	□
TADmatch™	□	□
3D Export (STL,OBJ, ZPR, inter alia)		■
Presentation (Slide Show, Media)	■	■
Letter Writing	■	■
Online forms (DHTML)	■	■
Internet browser	■	■
Recycling Bin (Lost & Found)	■	■
Web Reports	■	■
Container Exchange	■	■
Image and Result Export	■	■
Findings-Lexicon	■	■

* plus annual support fees
□ optional components against extra charge

Dentaurum Digital



prime4me® Aligner

After checking your data, our team will suggest a suitable aligner package for you.

Designation	REF	Quantity
prime4me® Aligner smart	430-140-20	1 case
prime4me® Aligner advanced	430-141-21	1 case
prime4me® Aligner performance	430-142-21	1 case

NEW

prime4me® base model

Printed base models for orthodontics.

REF 430-110-10	1 piece
----------------	---------



prime4me® dental arch model

Printed dental arch models for orthodontics.

REF 430-110-30	1 piece
----------------	---------



Order now:
digital.dentaurum.de

3D separating medium

Forms a separating layer between the acrylic models and Orthocryl® / Orthocryl® LC.
Color: green

REF 162-801-00

100 ml

**prime4me® RETAIN3R maxilla**

Unique fit using a three-dimensional construction design and highest biocompatibility: milled retainer in grade 5 titanium. Can be ordered together with prime4me® RETAIN3R optionally: Transfer tray for prime4me® RETAIN3R for placement in the mouth.

REF 430-120-30

1 piece

**prime4me® RETAIN3R mandible**

Unique fit using a three-dimensional construction design and highest biocompatibility: milled retainer in grade 5 titanium. Can be ordered together with prime4me® RETAIN3R optionally: Transfer tray for prime4me® RETAIN3R for placement in the mouth.

REF 430-120-40

1 piece


**Transfer tray**

Transfer tray for prime4me® RETAIN3R for placement in the mouth.

Can only be ordered together with prime4me® RETAIN3R.

REF 430-120-90

1 piece

NEW**Orthodontic demonstration model
prime4me® RETAIN3R**

For patient demonstration.

Fitted with:

- prime4me® base model
- prime4me® RETAIN3R upper jaw
- transfer tray for prime4me® RETAIN3R

REF 892-000-29

1 piece

NEW

Download data sets for
digital orthodontics:



Expansion screws without retention legs, for digital workflow



hyrax® click Mini – 7

RPE screw with click feature to prevent unwanted reversing

For effective and reliable rapid palatal expansion by means of fixed appliances.

- optimal screw locking during the active treatment phase
 - reliable activation / deactivation of the screw with clearly perceptible click after a $\frac{1}{4}$ turn
 - no accidental reversing when pulling out safety key during activation
 - built-in stops prevent the screws from coming apart once maximum expansion has been obtained
 - small dimensions with maximum expansion
 - offset body for earliest possible angling of the retention legs
 - permanently laser marked arrow on the body of the screw indicates the opening direction and max. expansion
- Pack (1 piece) with 1 safety key, pack (10 pieces) with 10 safety keys.




Screw body: 9.8 mm

Maximum expansion: 7.0 mm

4 x $\frac{1}{4}$ turn = 0.8 mm

Material: stainless steel

Retention legs	REF	Quantity	REF	Quantity
without	602-830-12	1 piece	602-830-32	10 pieces



hyrax® click Medium – 10

RPE screw with click feature to prevent unwanted reversing

For effective and reliable rapid palatal expansion by means of fixed appliances.

- optimal screw locking during the active treatment phase
 - reliable activation / deactivation of the screw with clearly perceptible click after a $\frac{1}{4}$ turn
 - no accidental reversing when pulling out safety key during activation
 - built-in stops prevent the screws from coming apart once maximum expansion has been obtained
 - small dimensions with maximum expansion
 - offset body for earliest possible angling of the retention legs
 - permanently laser marked arrow on the body of the screw indicates the opening direction and max. expansion
- Pack (1 piece) with 1 safety key, pack (10 pieces) with 10 safety keys.



Screw body: 12.8 mm

Maximum expansion: 10.0 mm

4 x $\frac{1}{4}$ turn = 0.8 mm

Material: stainless steel

Retention legs	REF	Quantity	REF	Quantity
without	602-832-12	1 piece	602-832-32	10 pieces



Expansion screws without retention legs, for digital workflow

hyrax® click Maxi – 12**RPE screw with click feature to prevent unwanted reversing****For effective and reliable rapid palatal expansion by means of fixed appliances.**

- optimal screw locking during the active treatment phase
 - reliable activation / deactivation of the screw with clearly perceptible click after a $\frac{1}{4}$ turn
 - no accidental reversing when pulling out safety key during activation
 - built-in stops prevent the screws from coming apart once maximum expansion has been obtained
 - small dimensions with maximum expansion
 - offset body for earliest possible angling of the retention legs
 - permanently laser marked arrow on the body of the screw indicates the opening direction and max. expansion
- Pack (1 piece) with 1 safety key, pack (10 pieces) with 10 safety keys.

Screw body: 14.8 mm



Maximum expansion: 12.0 mm

4 x $\frac{1}{4}$ turn = 0.8 mm

Material: stainless steel



Retention legs	REF	Quantity	REF	Quantity
without	602-834-12	1 piece	602-834-32	10 pieces



Expansion screws without retention legs, for digital workflow



hyrax® Mini – 7

RPE screw

For rapid palatal expansion by means of fixed appliances.

- built-in stops prevent the screws from coming apart once maximum expansion has been obtained
 - small dimensions with maximum expansion
 - offset body for earliest possible angling of the retention legs
 - the integrated spindle friction prevents unwanted reversing
 - permanently laser marked arrow on the body of the screw indicates the opening direction and max. expansion
 - screw can be tied down once desired expansion has been reached
- Pack (1 piece) with 1 safety key, pack (10 pieces) with 10 safety keys.




Screw body: 9.4 mm

Maximum expansion: 7.0 mm

4 x 1/4 turn = 0.8 mm

Material: stainless steel

Retention legs	REF	Quantity	REF	Quantity
without	602-800-12	1 piece	602-800-32	10 pieces



hyrax® Medium – 10

RPE screw

For rapid palatal expansion by means of fixed appliances.

- built-in stops prevent the screws from coming apart once maximum expansion has been obtained
- small dimensions with maximum expansion
- offset body for earliest possible angling of the retention legs
- the integrated spindle friction prevents unwanted reversing
- permanently laser marked arrow on the body of the screw indicates the opening direction and max. expansion
- screw can be tied down once desired expansion has been reached

Pack (1 piece) with 1 safety key, pack (10 pieces) with 10 safety keys.



Screw body: 12.4 mm

Maximum expansion: 10.0 mm

4 x 1/4 turn = 0.8 mm

Material: stainless steel

Retention legs	REF	Quantity	REF	Quantity
without	602-808-12	1 piece	602-808-32	10 pieces



Expansion screws without retention legs, for digital workflow

hyrax® Maxi – 12

RPE screw

For rapid palatal expansion by means of fixed appliances.

- built-in stops prevent the screws from coming apart once maximum expansion has been obtained
- small dimensions with maximum expansion
- offset body for earliest possible angling of the retention legs
- the integrated spindle friction prevents unwanted reversing
- permanently laser marked arrow on the body of the screw indicates the opening direction and max. expansion
- screw can be tied down once desired expansion has been reached

Pack (1 piece) with 1 safety key, pack (10 pieces) with 10 safety keys.

Screw body: 14.4 mm



Maximum expansion: 12.0 mm

4 x 1/4 turn = 0.8 mm

Material: stainless steel



Retention legs	REF	Quantity	REF	Quantity
without	602-811-12	1 piece	602-811-32	10 pieces



hyrax® neo Mini – 7

RPE screw

For rapid palatal expansion by means of a fixed appliance.

- (no active) expansion stop
- small dimensions with maximum expansion
- TAD-optimized design
- the integrated spindle friction prevents unwanted reversing
- permanently laser marked arrow on the body of the screw indicates the opening direction and maximum expansion
- screw can be tied down once desired expansion has been reached

Pack (1 piece) with 1 safety key, pack (10 pieces) with 10 safety keys.

Screw body: 11.2 mm



Maximum expansion: 7.0 mm

4 x 1/4 turn = 0.8 mm

Material: stainless steel



Retention legs	REF	Quantity	REF	Quantity
without	602-845-12	1 piece	602-845-32	10 pieces



Expansion screws without retention legs, for digital workflow



hyrax® neo Medium – 10

RPE screw

For rapid palatal expansion by means of a fixed appliance.

- (no active) expansion stop
 - small dimensions with maximum expansion
 - TAD-optimized design
 - the integrated spindle friction prevents unwanted reversing
 - permanently laser marked arrow on the body of the screw indicates the opening direction and maximum expansion
 - screw can be tied down once desired expansion has been reached
- Pack (1 piece) with 1 safety key, pack (10 pieces) with 10 safety keys.



Screw body: 12.4 mm

Maximum expansion: 10.0 mm

4 x 1/4 turn = 0.8 mm

Material: stainless steel

Retention legs	REF	Quantity	REF	Quantity
without	602-846-12	1 piece	602-846-32	10 pieces



hyrax® neo Maxi – 12

RPE screw

For rapid palatal expansion by means of a fixed appliance.

- (no active) expansion stop
 - small dimensions with maximum expansion
 - TAD-optimized design
 - the integrated spindle friction prevents unwanted reversing
 - permanently laser marked arrow on the body of the screw indicates the opening direction and maximum expansion
 - screw can be tied down once desired expansion has been reached
- Pack (1 piece) with 1 safety key, pack (10 pieces) with 10 safety keys.



Screw body: 14.4 mm

Maximum expansion: 12.0 mm

4 x 1/4 turn = 0.8 mm

Material: stainless steel

Retention legs	REF	Quantity	REF	Quantity
without	602-847-12	1 piece	602-847-32	10 pieces



Expansion screws without retention legs, for digital workflow

Titan hyrax® Maxi – 12

RPE screw of titanium

For rapid palatal expansion by means of fixed appliances.

The screw is anchored to the teeth using cast titanium and laser-welded cap splints or alternatively using acrylic cap splints fabricated from Orthocryl® acrylic.

Pack (1 piece) with 1 safety key, pack (10 pieces) with 10 safety keys.

Screw body: 14.4 mm




Maximum expansion: 12.0 mm

4 x 1/4 turn = 0.8 mm

Material: titanium



Retention legs	REF	Quantity	REF	Quantity
without	602-810-12	1 piece	602-810-32	10 pieces



Expansion screws without retention legs, for digital workflow



Variety click Mini – 8

Monoguided special expansion screw with click feature to prevent unwanted reversing. Primarily for treatment of malocclusions in the mandible using fixed appliances, e.g. for distalization, mesialization or for space closure by moving individual teeth. Also for early palatal expansion.

- optimal screw locking during the active treatment phase
- small dimensions with maximum expansion
- reliable activation / deactivation of the screw with clearly perceptible click after a $\frac{1}{4}$ turn
- permanently laser marked arrow on the body of the screw indicates the opening direction and max. expansion

Pack (1 piece) with 1 safety key, pack (10 pieces) with 10 safety keys.

Assembly dimension: 12.9 mm

Screw body: 9.6 mm

Maximum expansion: 8.0 mm

4 x $\frac{1}{4}$ turn = 0.8 mm

Material: stainless steel

Retention legs	REF	Quantity	REF	Quantity
without	602-836-12	1 piece	602-836-32	10 pieces



Variety click SP (small palate) Mini – 8

Monoguided special expansion screw with click feature to prevent unwanted reversing. For the fixed appliance technique. Can be widely used in the upper and lower jaw, for distalization, mesialization of posterior teeth, or for space closure by moving individual teeth. Also for early palatal expansion.

- optimal screw locking during the active treatment phase
- especially suitable for very narrow jaws
- retention legs can be bent very close to the screw body
- small dimensions with maximum expansion
- reliable activation / deactivation of the screw with clearly perceptible click after a $\frac{1}{4}$ turn
- permanently laser marked arrow on the body of the screw indicates the opening direction and max. expansion

Pack (1 piece) with 1 safety key, pack (10 pieces) with 10 safety keys.

Assembly dimension: 12.9 mm

Screw body: 9.6 mm



Maximum expansion: 8.0 mm

4 x $\frac{1}{4}$ turn = 0.8 mm

Material: stainless steel

Retention legs	REF	Quantity	REF	Quantity
without	602-838-12	1 piece	602-838-32	10 pieces

Dimensions Variety click / Variety click SP Mini – 8



Expansion screws without retention legs, for digital workflow

Variety click Maxi – 12

Monoguided special expansion screw with click feature to prevent unwanted reversing. Primarily for treatment of malocclusions in the mandible using fixed appliances, e.g. for distalization, mesialization or for space closure by moving individual teeth. Also for early palatal expansion.

- optimal screw locking during the active treatment phase
 - small dimensions with maximum expansion
 - reliable activation / deactivation of the screw with clearly perceptible click after a $\frac{1}{4}$ turn
 - permanently laser marked arrow on the body of the screw indicates the opening direction and max. expansion
- Pack (1 piece) with 1 safety key, pack (10 pieces) with 10 safety keys.

Assembly dimension: 17.2 mm

Screw body: 9.6 mm

Maximum expansion: 12.0 mm

4 x $\frac{1}{4}$ turn = 0.8 mm

Material: stainless steel



Retention legs	REF	Quantity	REF	Quantity
without	602-837-12	1 piece	602-837-32	10 pieces

Variety click SP (small palate) Maxi – 12

Monoguided special expansion screw with click feature to prevent unwanted reversing. Primarily for treatment of malocclusions in the mandible using fixed appliances, e.g. for distalization, mesialization or for space closure by moving individual teeth. Also for early palatal expansion.

- optimal screw locking during the active treatment phase
- especially suitable for very narrow jaws
- retention legs can be bent very close to the screw body
- small dimensions with maximum expansion
- reliable activation / deactivation of the screw with clearly perceptible click after a $\frac{1}{4}$ turn
- permanently laser marked arrow on the body of the screw indicates the opening direction and max. expansion

Pack (1 piece) with 1 safety key, pack (10 pieces) with 10 safety keys.

Assembly dimension: 17.2 mm

Screw body: 9.6 mm

Maximum expansion: 12.0 mm

4 x $\frac{1}{4}$ turn = 0.8 mm

Material: stainless steel





Retention legs	REF	Quantity	REF	Quantity
without	602-839-12	1 piece	602-839-32	10 pieces

Variety  click



Variety  click SP

Dimensions Variety click / Variety click SP Maxi – 12



Expansion screws without retention legs, for digital workflow



Variety Mini – 8

Monoguided special expansion screw. Primarily for treatment of malocclusions in the mandible using fixed appliances, e.g. for distalization, mesialization or for space closure by moving individual teeth. Also for early palatal expansion.

- small dimensions with maximum expansion
- the integrated spindle friction prevents unwanted reversing
- permanently laser marked arrow on the body of the screw indicates the opening direction and max. expansion
- screw can be tied down once desired expansion has been reached

Pack (1 piece) with 1 safety key, pack (10 pieces) with 10 safety keys.

Assembly dimension: 12.9 mm

Screw body: 9.6 mm

Maximum expansion: 8.0 mm

4 x 1/4 turn = 0.8 mm

Material: stainless steel

Retention legs	REF	Quantity	REF	Quantity
without	602-815-12	1 piece	602-815-32	10 pieces



Variety SP (small palate) Mini – 8

Monoguided special expansion screw. Primarily for treatment of malocclusions in the mandible using fixed appliances, e.g. for distalization, mesialization or for space closure by moving individual teeth. Also for early palatal expansion.

- especially suitable for very narrow jaws
- retention legs can be bent very close to the screw body
- small dimensions with maximum expansion
- the integrated spindle friction prevents unwanted reversing
- permanently laser marked arrow on the body of the screw indicates the opening direction and max. expansion
- screw can be tied down once desired expansion has been reached

Pack (1 piece) with 1 safety key, pack (10 pieces) with 10 safety keys.

Assembly dimension: 12.9 mm

Screw body: 9.6 mm



Maximum expansion: 8.0 mm

4 x 1/4 turn = 0.8 mm

Material: stainless steel

Retention legs	REF	Quantity	REF	Quantity
without	602-819-12	1 piece	602-819-32	10 pieces

Dimensions Variety / Variety SP Mini – 8



Expansion screws without retention legs, for digital workflow

Variety Maxi – 12

Monoguided special expansion screw. Primarily for treatment of malocclusions in the mandible using fixed appliances, e.g. for distalization, mesialization or for space closure by moving individual teeth. Also for early palatal expansion.

- small dimensions with maximum expansion
- the integrated spindle friction prevents unwanted reversing
- permanently laser marked arrow on the body of the screw indicates the opening direction and max. expansion
- screw can be tied down once desired expansion has been reached

Pack (1 piece) with 1 safety key, pack (10 pieces) with 10 safety keys.

Assembly dimension: 14.9 mm

Screw body: 9.6 mm

Maximum expansion: 12.0 mm

4 x 1/4 turn = 0.8 mm

Material: stainless steel



Retention legs	REF	Quantity	REF	Quantity
without	602-806-12	1 piece	602-806-32	10 pieces

Variety SP (small palate) Maxi – 12

Monoguided special expansion screw. Primarily for treatment of malocclusions in the mandible using fixed appliances, e.g. for distalization, mesialization or for space closure by moving individual teeth. Also for early palatal expansion.

- especially suitable for very narrow jaws
- retention legs can be bent very close to the screw body
- small dimensions with maximum expansion
- the integrated spindle friction prevents unwanted reversing
- permanently laser marked arrow on the body of the screw indicates the opening direction and max. expansion
- screw can be tied down once desired expansion has been reached

Pack (1 piece) with 1 safety key, pack (10 pieces) with 10 safety keys.

Assembly dimension: 14.9 mm

Screw body: 9.6 mm

Maximum expansion: 12.0 mm

4 x 1/4 turn = 0.8 mm

Material: stainless steel



Retention legs	REF	Quantity	REF	Quantity
without	602-809-12	1 piece	602-809-32	10 pieces



Variety



Variety SP



Dimensions Variety / Variety SP Maxi – 12



OnyxCeph^{3TM}


 OnyxCeph^{3TM}	OnyxCeph^{3TM} Basic (purchase license) Software for archiving and patient consultation incl. 5 network working stations for simultaneous operation.	REF 410-001-00	1 piece
 OnyxCeph^{3TM}	OnyxCeph^{3TM} 2D Pro (purchase license) Software for archiving, diagnosis, treatment planning and patient consultation in 2D incl. 5 network working stations for simultaneous operation.	REF 410-002-00	1 piece
 OnyxCeph^{3TM}	OnyxCeph^{3TM} 3D Pro (purchase license) Software for archiving, diagnosis, treatment planning and patient consultation in 3D incl. 5 network working stations for simultaneous operation.	REF 410-004-00	1 piece
 OnyxCeph^{3TM}	OnyxCeph^{3TM} 3D Lab (purchase license) Software for archiving, diagnosis, treatment planning and patient consultation in 3D incl. Option Virtual Set-up (V.T.O.3D) and a total of 20 network working places for simultaneous operation.	REF 410-005-00	1 piece
 OnyxCeph^{3TM}	Additional server license Second network license (LAN 5/1) / same practice name.	REF 410-000-01	1 piece
 OnyxCeph^{3TM}	Additional standalone license Notebook license (1/1) for same practice name.	REF 410-000-02	1 piece
 OnyxCeph^{3TM}	Additional site license Additional site license in network.	REF 410-000-03	1 piece
 OnyxCeph^{3TM}	User administration program Permission settings depending on user login.	REF 410-000-04	1 piece
 OnyxCeph^{3TM}	Statistics program Statistical evaluation of findings.	REF 410-000-05	1 piece

Module Waefer Creation

Creation of surgical splints.

REF 410-000-06

1 piece




Module V.T.O.3D (Virtual Treatment Objective)

Creation of virtual set ups.

REF 410-000-07

1 piece



Module Bonding Trays 3D

Construction of bracket transfer trays.

REF 410-000-11

1 piece




OnyxCeph^{3TM} module Bite Splint 3D

Construction of occlusal splints, anti-snoring splints, sports guards and other splints.

REF 410-000-13

1 piece




Module Ortho Apps 3D


Construction of RPE appliances, lingual arch appliances, Herbst appliances and others.

REF 410-000-14

1 piece



OnyxCeph^{3TM}




Module TADmatch™

Module from OnyxCeph^{3TM} for the virtual positioning of tomas®-pin and as a basis for the fabrication of insertion templates.

REF 410-000-15

1 piece




Module Wire_Bonding

Positioning of wires / brackets on set-up.

REF 410-000-08

1 piece



Module Bracket Erase 3D

Erase virtual bracket.

REF 410-000-12

1 piece




Module Aligner 3D

Planning of aligners.

REF 410-000-09

1 piece



Module Retainer 3D

Planning of retainers.

REF 410-000-10

1 piece



20 horizontal lines for writing.

Lingual brackets	Page	46
Lingual tubes	Page	47
Arches	Page	48
Instruments	Page	49
Accessories	Page	50



Lingual system.

Surprisingly different.



DENTAURUM
**ONLINE
SHOP**
shop.dentaurum.com





Incredibly effective.

discovery® delight dramatically reduces the number of components required for treatment while ensuring a free selection of all lingual treatment options, from 2D to 3D treatment.

In collaboration with Prof. Dr. Rubens Demichieri, Montevideo, a lingual system that addresses patient needs and offers all the advantages of the trusted discovery® bracket family emerged:

- Highest quality.
- Virtually invisible.
- Short treatment time.
- Very comfortable to wear.
- Perfect design.






Lingual brackets / Assortments

discovery® delight

Slot: 0.51 mm x 0.51 mm / 20 x 20 - with vertical hook

Tooth	16	15	14	13	12	11	21	22	23	24	25	26
Width	2.0 mm	1.5 mm	1.5 mm	1.5 mm	1.5 mm	1.5 mm	1.5 mm	1.5 mm	1.5 mm	1.5 mm	1.5 mm	2.0 mm
REF	724-102-51	790-452-00	790-452-00	790-450-00	790-450-00	790-450-00	790-450-00	790-450-00	790-450-00	790-452-00	790-452-00	724-102-51
Quantity	1 piece	1 piece	1 piece	1 piece	1 piece	1 piece	1 piece	1 piece	1 piece	1 piece	1 piece	1 piece
REF	724-103-51	790-453-00	790-453-00	790-451-00	790-451-00	790-451-00	790-451-00	790-451-00	790-451-00	790-453-00	790-453-00	724-103-51
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces



Tooth	46	45	44	43	42	41	31	32	33	34	35	36
Width	2.0 mm	1.5 mm	1.5 mm	1.5 mm	1.5 mm	1.5 mm	1.5 mm	1.5 mm	1.5 mm	1.5 mm	1.5 mm	2.0 mm
REF	724-102-51	790-452-00	790-452-00	790-450-00	790-450-00	790-450-00	790-450-00	790-450-00	790-450-00	790-452-00	790-452-00	724-102-51
Quantity	1 piece	1 piece	1 piece	1 piece	1 piece	1 piece	1 piece	1 piece	1 piece	1 piece	1 piece	1 piece
REF	724-103-51	790-453-00	790-453-00	790-451-00	790-451-00	790-451-00	790-451-00	790-451-00	790-451-00	790-453-00	790-453-00	724-103-51
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Assortments:

discovery® delight social six pack	1 case	6 anterior brackets (REF 790-450-00)	REF 891-406-00
discovery® delight single jaw	1 case	6 anterior brackets (REF 790-450-00) + 4 posterior brackets (REF 790-452-00)	REF 891-407-00
discovery® delight single jaw + tubes	1 case	6 anterior brackets (REF 790-450-00) + 4 posterior brackets (REF 790-452-00) + 2 molar tubes (REF 724-102-51)	REF 891-408-00
discovery® delight complete case	1 case	12 anterior brackets (REF 790-450-00) + 8 posterior brackets (REF 790-452-00) + 4 molar tubes (REF 724-102-51)	REF 891-409-00

Lingual brackets and lingual tubes

discovery® delight anterior bracket

with vertical hook

Tooth	Torque	Distal Offset	In/Out	Hook	Width	REF	Quantity
13-11 / 21-23 / 43-41 / 31-33				•	1.5 mm	790-450-00	1 piece
13-11 / 21-23 / 43-41 / 31-33				•	1.5 mm	790-451-00	10 pieces

Slot: 0.51 mm x 0.51 mm / 20 x 20



discovery® delight posterior bracket, open

with vertical hook

Tooth	Torque	Distal Offset	In/Out	Hook	Width	REF	Quantity
15-14 / 24-25 / 45-44 / 34-35				•	1.5 mm	790-452-00	1 piece
15-14 / 24-25 / 45-44 / 34-35				•	1.5 mm	790-453-00	10 pieces

Slot: 0.51 mm x 0.51 mm / 20 x 20



discovery® delight posterior bracket, closed

non-convertible · with vertical hook

Tooth	Torque	Distal Offset	In/Out	Hook	Width	REF	Quantity
15-14 / 24-25 / 45-44 / 34-35				•	2.0 mm	724-100-51	1 piece
15-14 / 24-25 / 45-44 / 34-35				•	2.0 mm	724-101-51	10 pieces

Slot: 0.51 mm x 0.51 mm / 20 x 20



discovery® delight molar tube / Bonding technique

non-convertible · with vertical hook


Tooth	Torque	Distal Offset	In/Out	Hook	Width	REF	Quantity
16 / 26 / 46 / 36				•	2.0 mm	724-102-51	1 piece
16 / 26 / 46 / 36				•	2.0 mm	724-103-51	10 pieces

Slot: 0.51 mm x 0.51 mm / 20 x 20



Arches

Nickel titanium



rematitan® LITE lingual arch, round

Strength category: super elastic

ø	Length	Diagonal M.	REF	Quantity
0.25 mm / 10			766-760-00	10 pieces
0.30 mm / 12			766-761-00	10 pieces
0.35 mm / 14			766-762-00	10 pieces
0.40 mm / 16			766-763-00	10 pieces
0.45 mm / 18			766-764-00	10 pieces



rematitan® LITE lingual arch, rectangular

Strength category: super elastic

∅	Length	Diagonal M.	REF	Quantity
0.41 mm x 0.41 mm / 16 x 16		0.54 mm	766-765-00	10 pieces
0.48 mm x 0.48 mm / 19 x 19		0.62 mm	766-766-00	10 pieces

Stainless steel (Nickel chromium)



remanium® lingual arch, round

Strength category: spring hard · 1800 N/mm² - 2000 N/mm²


ø	Length	Diagonal M.	REF	Quantity
0.40 mm / 16			766-767-00	10 pieces
0.45 mm / 18			766-768-00	10 pieces



remanium® lingual arch, rectangular

Strength category: spring hard · 1800 N/mm² - 2000 N/mm²

∅	Length	Diagonal M.	REF	Quantity
0.46 mm x 0.46 mm / 18 x 18			766-769-00	10 pieces
0.48 mm x 0.48 mm / 19 x 19			766-770-00	10 pieces



dentaflex® lingual arch, round, 6-strand twisted

Strength category: super spring hard · 2500 N/mm² - 2700 N/mm²

ø	Length	Diagonal M.	REF	Quantity
0.38 mm / 15			766-771-00	10 pieces

Lingual bending pliers

sterilizable

For bending steps and compensation bends in archwires for the lingual technique.

Steps: 1.5 mm; 2.0 mm; 2.5 mm and 3.0 mm

Max. wire thickness: ø 0.5 mm, hard

Premium-Line · 10 year guarantee

REF 003-001-00

1 piece



Mosquito forceps, curved

sterilizable

For placing and ligating Dentalastics® plastic ligatures, especially for the lingual technique.

REF 000-732-00

1 piece



Crimping tweezers

sterilizable

Delicate tweezer design for easy bracket and buccal tube placement.

REF 025-280-00

1 piece



Positioning tool, lingual

sterilizable

REF 023-278-00

1 piece



Dental mirror with reference lines (assembled), lingual

sterilizable

REF 023-279-00

1 piece



Dental mirror with reference lines (without grip), lingual

sterilizable

REF 023-279-10

1 piece



Bracket removing pliers, angled

sterilizable

For easy removal of lingual brackets and brackets in the posterior region. Inserts made of hard metal.

Premium-Line · 10 year guarantee

REF 004-347-00

1 piece



Accessories



SUPER Elasto-Force plastic chains

Consistent traction force over a long period of time. High elasticity with excellent elastic recovery.

Color	Inner ø	Length	Connector	REF	Quantity
transparent	1.3 mm	4.6 m	•	774-220-01	1 roll
transparent	1.3 mm	4.6 m		774-221-01	1 roll



**Orthodontic demonstration model
discovery® delight**

Fitted with:

- discovery® delight lingual brackets
 - discovery® delight molar tubes
 - rematitan® LITE lingual arch
 - dentaflex® lingual arch
- Maxilla and mandible made of plastic, hinged.

REF 892-000-22	1 piece
----------------	---------



discovery® sl / Assortments

Roth 18/22	Page	56
------------	------	----

discovery® sl 2.0 / Assortment

McLaughlin-Bennett-Trevisi* 22	Page	57
--------------------------------	------	----

Instruments and accessories

Page	58
------	----

* The Dentaaurum version of this prescription is not claimed to be a duplication of any other, nor does Dentaaurum imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.



FDI System

18	17	16	15	14	13	12	11	21	22	23	24	25	26	27	28		
8	7	6	5	4	3	2	1	1	2	3	4	5	6	7	8		
R	8	7	6	5	4	3	2	1	1	2	3	4	5	6	7	8	L
	48	47	46	45	44	43	42	41	31	32	33	34	35	36	37	38	

Self-ligating brackets.

Dynamic and functional.



DENTAURUM
**ONLINE
SHOP**
shop.dentaurum.com



The logo for the discovery® sl self-ligating bracket system. It features the word "discovery" in a black, italicized sans-serif font, followed by a registered trademark symbol (®). To the right of "discovery" is the letters "sl" in a red, italicized sans-serif font. A thin black curved line arches over the text.The logo for the discovery® sl 2.0 self-ligating bracket system. It features the word "discovery" in a black, italicized sans-serif font, followed by a registered trademark symbol (®). Below "discovery" is the text "sl 2.0" in a blue, italicized sans-serif font. A thin black curved line arches over the text.

Orthodontic treatment at a high standard.

The self-ligating bracket systems discovery® sl from Dentaaurum have much to offer: for the patient, the brackets are easy to handle and are comfortable to wear; for the dental office, they are easy to handle and offer effective treatment with only short periods of time in the chair.

discovery® sl and discovery® sl 2.0 still rank among the smallest and flattest sl-brackets on the market. At the heart of the design lies a closing mechanism in the form of a lid. The lid has a smooth surface making it comfortable for the patient to wear and is designed to be closed during treatment, thereby offering a high level of hygiene. An instrument is available specifically for opening and closing the lid. The sl brackets from Dentaaurum impress both patients and orthodontists and fulfill all of the requirements expected of a modern bracket system.


discovery® sl

One of the smallest sl-brackets in the world.



discovery® sl 2.0

Proven design for your success.



Self-ligating brackets


discovery® sl / Assortments

Roth 18

Slot: 0.46 mm x 0.64 mm / 18 x 25 · Cuspid brackets with hook

Torque	- 7°	- 7°	- 2°	+ 8°	+ 12°	+ 12°	+ 8°	- 2°	- 7°	- 7°
Angulation	0°	0°	+ 13°	+ 9°	+ 5°	+ 5°	+ 9°	+ 13°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	2.6 mm	2.6 mm	2.7 mm	2.4 mm	2.6 mm	2.6 mm	2.4 mm	2.7 mm	2.6 mm	2.6 mm
REF	718-101-15	718-101-14	718-141-13	718-101-12	718-101-11	718-101-21	718-101-22	718-141-23	718-101-24	718-101-25
Quantity	1 piece	1 piece	1 piece	1 piece	1 piece	1 piece	1 piece	1 piece	1 piece	1 piece

Even easier bracket identification with FDI-marking on the base!



Torque	- 22°	- 17°	- 11°	- 1°	- 1°	- 1°	- 1°	- 11°	- 17°	- 22°
Angulation	0°	0°	+ 7°	0°	0°	0°	0°	+ 7°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	1.3 mm	1.3 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	2.6 mm	2.6 mm	2.7 mm	2.4 mm	2.4 mm	2.4 mm	2.4 mm	2.7 mm	2.6 mm	2.6 mm
REF	718-101-45	718-101-44	718-141-43	718-101-42	718-101-41	718-101-31	718-101-32	718-141-33	718-101-34	718-101-35
Quantity	1 piece	1 piece	1 piece	1 piece	1 piece	1 piece	1 piece	1 piece	1 piece	1 piece

Assortments:


1 case	20 brackets	REF 718-101-01
10 cases	200 brackets	REF 718-101-10

Roth 22

Slot: 0.56 mm x 0.64 mm / 22 x 25 · Cuspid brackets with hook

Torque	- 7°	- 7°	- 2°	+ 8°	+ 12°	+ 12°	+ 8°	- 2°	- 7°	- 7°
Angulation	0°	0°	+ 13°	+ 9°	+ 5°	+ 5°	+ 9°	+ 13°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	2.6 mm	2.6 mm	2.7 mm	2.4 mm	2.6 mm	2.6 mm	2.4 mm	2.7 mm	2.6 mm	2.6 mm
REF	722-101-15	722-101-14	722-141-13	722-101-12	722-101-11	722-101-21	722-101-22	722-141-23	722-101-24	722-101-25
Quantity	1 piece	1 piece	1 piece	1 piece	1 piece	1 piece	1 piece	1 piece	1 piece	1 piece

Even easier bracket identification with FDI-marking on the base!



Torque	- 22°	- 17°	- 11°	- 1°	- 1°	- 1°	- 1°	- 11°	- 17°	- 22°
Angulation	0°	0°	+ 7°	0°	0°	0°	0°	+ 7°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	1.3 mm	1.3 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	2.6 mm	2.6 mm	2.7 mm	2.4 mm	2.4 mm	2.4 mm	2.4 mm	2.7 mm	2.6 mm	2.6 mm
REF	722-101-45	722-101-44	722-141-43	722-101-42	722-101-41	722-101-31	722-101-32	722-141-33	722-101-34	722-101-35
Quantity	1 piece	1 piece	1 piece	1 piece	1 piece	1 piece	1 piece	1 piece	1 piece	1 piece

Assortments:


1 case	20 brackets	REF 722-101-01
10 cases	200 brackets	REF 722-101-10

McLaughlin-Bennett-Trevisi** 22

Slot: 0.56 mm x 0.70 mm / 22 x 28 · Cuspid brackets with hook

Torque	- 7°	- 7°	0°	+ 10°	+ 17°	+ 17°	+ 10°	0°	- 7°	- 7°
Angulation	0°	0°	+ 8°	+ 8°	+ 4°	+ 4°	+ 8°	+ 8°	0°	0°
In/Out	1.2 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	1.2 mm
Width	2.6 mm	2.6 mm	2.6 mm	2.6 mm	2.6 mm	2.6 mm	2.6 mm	2.6 mm	2.6 mm	2.6 mm
REF	722-102-15	722-102-14	722-142-13	722-102-12	722-102-11	722-102-21	722-102-22	722-142-23	722-102-24	722-102-25
Quantity	1 piece	1 piece	1 piece	1 piece	1 piece	1 piece	1 piece	1 piece	1 piece	1 piece

Even easier bracket identification with FDI-marking on the base!



Torque	- 17°	- 12°	0°	- 6°	- 6°	- 6°	- 6°	0°	- 12°	- 17°
Angulation	+ 2°	+ 2°	+ 3°	0°	0°	0°	0°	+ 3°	+ 2°	+ 2°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	1.3 mm	1.3 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	2.6 mm	2.6 mm	2.6 mm	2.6 mm	2.6 mm	2.6 mm	2.6 mm	2.6 mm	2.6 mm	2.6 mm
REF	722-102-45	722-102-44	722-142-43	722-102-42	722-102-41	722-102-31	722-102-32	722-142-33	722-102-34	722-102-35
Quantity	1 piece	1 piece	1 piece	1 piece	1 piece	1 piece	1 piece	1 piece	1 piece	1 piece

Assortments:

1 case	20 brackets	REF 722-102-01
--------	-------------	-----------------------

** The Dentaureum version of this prescription is not claimed to be a duplication of any other, nor does Dentaureum imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.



Self-ligating brackets

Instruments and accessories



sl-instrument, for discovery® sl

sterilizable

For opening and closing self-ligating brackets **discovery® sl** (Roth system), **dynamique® m** and **dynamique® c**, and for easier ligation of the archwire.

REF 023-276-00

1 piece



sl-instrument, for discovery® sl 2.0

sterilizable

For opening and closing self-ligating brackets **discovery® sl 2.0** (McLaughlin-Bennett-Trevisi* system) and easier ligation of the archwire.

REF 023-277-00

1 piece

* The Dentaaurum version of this prescription is not claimed to be a duplication of any other, nor does Dentaaurum imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.



Orthodontic demonstration model discovery® sl

Angle Class I, normal occlusion

Roth 22

Fitted with:

- discovery® sl self-ligating metal brackets
- Ortho-Cast M-Series mini buccal tubes
- Equire stainless steel arch, maxilla
- rematitan® sl arch, mandible

Maxilla and mandible made of plastic, hinged.

REF 892-000-17

1 piece



Orthodontic demonstration model discovery® sl 2.0

Angle Class I, normal occlusion

McLaughlin-Bennett-Trevisi* 22

Fitted with:

- discovery® sl 2.0 self-ligating metal brackets
- Ortho-Cast M-Series and M-Series mini buccal tubes
- TensiC® arches

Maxilla and mandible made of plastic, hinged.

REF 892-000-19

1 piece

* The Dentaaurum version of this prescription is not claimed to be a duplication of any other, nor does Dentaaurum imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.

Metal brackets

discovery® / Assortments

Roth 18/22	Page	68
Ricketts® – IV. Dimension 18	Page	70
Ricketts® Universal 18	Page	70
Hilgers 18	Page	71
Standard Edgewise 18/22	Page	71

discovery® / Single packs – bonding technique

Roth 18/22	Page	73
------------	------	----

discovery® smart / Assortments

Roth 18/22	Page	74
McLaughlin-Bennett-Trevisi* 18/22	Page	76

topic / Assortments

Roth 18/22	Page	78
------------	------	----

topic / Single packs – bonding technique

Roth 18/22	Page	80
------------	------	----

equilibrium® mini / Assortments

Roth 18/22	Page	81
McLaughlin-Bennett-Trevisi* 18/22	Page	83

equilibrium® mini / Single packs – bonding technique

McLaughlin-Bennett-Trevisi* 18/22	Page	85
-----------------------------------	------	----

equilibrium® 2 / Assortments

Roth 18/22	Page	86
McLaughlin-Bennett-Trevisi* 18/22	Page	88
Ricketts® Universal 18	Page	90
Standard Edgewise 18/22	Page	90

equilibrium® 2 / Single packs – bonding technique

Roth 18/22	Page	92
Ricketts® Universal 18	Page	92

equilibrium® 2 / Single packs – banding technique

Roth 18/22	Page	93
Roth, McLaughlin-Bennett-Trevisi* 18/22	Page	93
McLaughlin-Bennett-Trevisi* 18/22	Page	93
Standard Edgewise 18/22	Page	94
Individual systems 18/22	Page	94

equilibrium® ti / Assortments

Roth 18/22	Page	95
McLaughlin-Bennett-Trevisi* 18/22	Page	96

Bracket accessories

Page	114
------	-----

* The Dentaaurum version of this prescription is not claimed to be a duplication of any other, nor does Dentaaurum imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.



Metal brackets.

Always a good choice.



DENTAURUM
**ONLINE
SHOP**
shop.dentaurum.com



Metal brackets



High performance for precision.

Our high-quality brackets are designed and manufactured using state-of-the-art design and production procedures. The highest precision is essential for ensuring safe, reliable use in orthodontic treatment. The diversity of Dentaurem metal brackets is unique in the world and provides application options for all indications. The metal brackets are manufactured in one piece from highly corrosion-resistant, medical-grade stainless steel alloys. This makes them particularly biocompatible for the patient, as no solder is used. The wide-ranging program also includes nickel-free* brackets made from a special CoCr alloy and pure titanium brackets, which are ideal for use with patients with allergies.



* see Materials list Orthodontics

discovery®

The sum of experience.



discovery® smart

Small and brilliant.



topic

The nickel-free* alternative.



equilibrium® mini

The mini bracket.



equilibrium® 2

Harmony in function and comfort.



equilibrium® ti

Outstanding biocompatibility.



FDI System

	18	17	16	15	14	13	12	11		21	22	23	24	25	26	27	28	
	8	7	6	5	4	3	2	1		1	2	3	4	5	6	7	8	
R	8	7	6	5	4	3	2	1		1	2	3	4	5	6	7	8	L
	48	47	46	45	44	43	42	41		31	32	33	34	35	36	37	38	

The discovery® family.



The discovery® family.

The discovery® bracket family comprises the premium bracket line of the Dentaaurum metal bracket range. They are manufactured completely as one-piece brackets using the MIM technique – Metal Injection Molding – in state-of-the-art facilities. The hooks on the canine and, if desired, on the first and second premolar brackets, are directly integrated into the mold and are also produced in one piece together with the body of the bracket. The innovative mushroom shape provides retention for more elastics. Our diverse system offers the right bracket for every therapy.

Innovative production procedures guarantee maximum success.

The bases of all brackets of the discovery® family are ideally adapted to the tooth contour using 3D CAD (Computer Aided Design), allowing optimal handling for the user when bonding the bracket to the tooth. The anatomical shape facilitates positioning the bracket on the tooth and also ensures it is firmly retained on the tooth.



Optimal bond strength and easy debonding.

The patented laser-structured base with FDI notation guarantees reliable retention during treatment. The ideal distribution of micro and macro retention produces an optimal bond between tooth/adhesive and the bracket. Evenly-spaced retentions support the bonding of the brackets on the tooth surface. A debonding channel ensures very easy debonding of conventionally ligated brackets. The channel allows the bracket face to flex, ensuring safe removal from the tooth.

Individuality combined with easy handling.

All discovery® brackets have been designed so that there is sufficient space for multiple ligatures. In addition, rounded bracket tie wings facilitate ligation and provide intra-oral comfort without sharp edges for the patient. All discovery® brackets include bio-compatible laser marking for easier identification and positioning.

Diverse combination options.

As all brackets of the discovery® range have the same design, there is no limit to the individuality of the treatment. All brackets can be easily combined with one another. A low In/Out of 0.7 mm allows individual combination of the brackets within the family without using any compensation bends.



The equilibrium® family.

equilibrium® **mini**

equilibrium® 2

equilibrium® ti

The equilibrium® family.

The equilibrium® series is Dentaurem's classic bracket series. Reliable results in patient treatment over many years have proven that the attractive and functional design is successful.

The brackets from the equilibrium® family are formed and milled from rod material and undergo many fully-automated processing steps. Like many other Dentaurem brackets, they are genuine one-piece brackets.

equilibrium® 2 and equilibrium® mini brackets are produced in one piece – body and base – from high-strength stainless steel*. equilibrium® ti brackets are made from titanium** and are nickel-free. Brackets from the equilibrium® family are particularly bio-compatible as no solders are used.



* Material number acc. to DIN 1.4404

** Material number acc. to DIN 3.7035



Optimal bond strength and easy debonding.

The patented laser-structured base with FDI notation guarantees reliable retention during treatment with equilibrium® brackets. The ideal distribution of micro and macro retention produces an optimal bond between tooth/adhesive and the bracket. Uniform retention prevents premature debonding of the brackets.

Easy handling and reliable treatment.

All equilibrium® brackets include biocompatible laser marking for easier identification and positioning. FDI notation is also laser marked on the base. This system of symbols aids in assigning each bracket to a specific tooth. The laser marking is abrasion resistant and does not contain any harmful coloring agents.


Proven design.

The design of equilibrium® brackets provides excellent bracket tie wing undercuts, which also allow the use of multiple ligatures at any time. The mesial and distal ends of the bracket have been designed to make it easier to grip and place the brackets. The design of the rhomboid base with additional central marking facilitates placement on the clinical crown. The low In/Out values enable an ideal application of force.

Roth 18

Slot: 0.46 mm x 0.76 mm / 18 x 30 · Cuspid brackets with hook

Torque	- 7°	- 7°	- 2°	+ 8°	+ 12°	+ 12°	+ 8°	- 2°	- 7°	- 7°
Angulation	0°	0°	+ 13°	+ 9°	+ 5°	+ 5°	+ 9°	+ 13°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	3.4 mm	3.4 mm	3.4 mm	3.2 mm	3.4 mm	3.4 mm	3.2 mm	3.4 mm	3.4 mm	3.4 mm
REF	790-106-00	790-106-00	790-104-00	790-102-00	790-163-00	790-164-00	790-103-00	790-105-00	790-106-00	790-106-00
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces



Torque	- 22°	- 17°	- 11°	- 1°	- 1°	- 1°	- 1°	- 11°	- 17°	- 22°
Angulation	0°	0°	+ 7°	0°	0°	0°	0°	+ 7°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	1.3 mm	1.3 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	3.3 mm	3.3 mm	3.1 mm	2.6 mm	2.6 mm	2.6 mm	2.6 mm	3.1 mm	3.3 mm	3.3 mm
REF	790-111-00	790-110-00	790-108-00	790-107-00	790-107-00	790-107-00	790-107-00	790-109-00	790-110-00	790-111-00
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Assortments:

1 case	20 brackets	REF 891-180-00
10 cases	200 brackets	REF 891-181-00
50 cases	1000 brackets	REF 891-182-00


Note:

Other brackets with Roth values see single packs.

Roth 18

Slot: 0.46 mm x 0.76 mm / 18 x 30 · Cuspid, 1st bicuspid and 2nd bicuspid brackets with hook

Torque	- 7°	- 7°	- 2°	+ 8°	+ 12°	+ 12°	+ 8°	- 2°	- 7°	- 7°
Angulation	0°	0°	+ 13°	+ 9°	+ 5°	+ 5°	+ 9°	+ 13°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	3.4 mm	3.4 mm	3.4 mm	3.2 mm	3.4 mm	3.4 mm	3.2 mm	3.4 mm	3.4 mm	3.4 mm
REF	790-141-00	790-141-00	790-104-00	790-102-00	790-163-00	790-164-00	790-103-00	790-105-00	790-142-00	790-142-00
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces



Torque	- 22°	- 17°	- 11°	- 1°	- 1°	- 1°	- 1°	- 11°	- 17°	- 22°
Angulation	0°	0°	+ 7°	0°	0°	0°	0°	+ 7°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	1.3 mm	1.3 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	3.3 mm	3.3 mm	3.1 mm	2.6 mm	2.6 mm	2.6 mm	2.6 mm	3.1 mm	3.3 mm	3.3 mm
REF	790-145-00	790-143-00	790-108-00	790-107-00	790-107-00	790-107-00	790-107-00	790-109-00	790-144-00	790-146-00
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Assortments:

1 case	20 brackets	REF 891-295-00
10 cases	200 brackets	REF 891-296-00
50 cases	1000 brackets	REF 891-297-00


Note:

Other brackets with Roth values see single packs.

Roth 22

Slot: 0.56 mm x 0.76 mm / 22 x 30 · Cuspid brackets with hook

Torque	- 7°	- 7°	- 2°	+ 8°	+ 12°	+ 12°	+ 8°	- 2°	- 7°	- 7°
Angulation	0°	0°	+ 13°	+ 9°	+ 5°	+ 5°	+ 9°	+ 13°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	3.4 mm	3.4 mm	3.4 mm	3.2 mm	3.4 mm	3.4 mm	3.2 mm	3.4 mm	3.4 mm	3.4 mm
REF	790-118-00	790-118-00	790-116-00	790-114-00	790-165-00	790-166-00	790-115-00	790-117-00	790-118-00	790-118-00
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces



Torque	- 22°	- 17°	- 11°	- 1°	- 1°	- 1°	- 1°	- 11°	- 17°	- 22°
Angulation	0°	0°	+ 7°	0°	0°	0°	0°	+ 7°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	1.3 mm	1.3 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	3.3 mm	3.3 mm	3.1 mm	2.6 mm	2.6 mm	2.6 mm	2.6 mm	3.1 mm	3.3 mm	3.3 mm
REF	790-123-00	790-122-00	790-120-00	790-119-00	790-119-00	790-119-00	790-119-00	790-121-00	790-122-00	790-123-00
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Assortments:

1 case	20 brackets	REF 891-183-00
10 cases	200 brackets	REF 891-184-00
50 cases	1000 brackets	REF 891-185-00


Note:

Other brackets with Roth values see single packs.

Roth 22

Slot: 0.56 mm x 0.76 mm / 22 x 30 · Cuspid, 1st bicuspid and 2nd bicuspid brackets with hook

Torque	- 7°	- 7°	- 2°	+ 8°	+ 12°	+ 12°	+ 8°	- 2°	- 7°	- 7°
Angulation	0°	0°	+ 13°	+ 9°	+ 5°	+ 5°	+ 9°	+ 13°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	3.4 mm	3.4 mm	3.4 mm	3.4 mm	3.2 mm	3.4 mm	3.2 mm	3.4 mm	3.4 mm	3.4 mm
REF	790-147-00	790-147-00	790-116-00	790-114-00	790-165-00	790-166-00	790-115-00	790-117-00	790-148-00	790-148-00
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces



Torque	- 22°	- 17°	- 11°	- 1°	- 1°	- 1°	- 1°	- 11°	- 17°	- 22°
Angulation	0°	0°	+ 7°	0°	0°	0°	0°	+ 7°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	1.3 mm	1.3 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	3.3 mm	3.3 mm	3.1 mm	2.6 mm	2.6 mm	2.6 mm	2.6 mm	3.1 mm	3.3 mm	3.3 mm
REF	790-151-00	790-149-00	790-120-00	790-119-00	790-119-00	790-119-00	790-119-00	790-121-00	790-150-00	790-152-00
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Assortments:

1 case	20 brackets	REF 891-298-00
10 cases	200 brackets	REF 891-299-00
50 cases	1000 brackets	REF 891-300-00

Note:


Other brackets with Roth values see single packs.

discovery® / Assortments

Ricketts® - IV. dimension 18

Slot: 0.46 mm x 0.76 mm / 18 x 30 - Cuspid brackets with hook

Torque	- 7°	- 7°	+ 7°	+ 14°	+ 22°	+ 22°	+ 14°	+ 7°	- 7°	- 7°
Angulation	0°	0°	+ 5°	+ 8°	+ 5°	+ 5°	+ 8°	+ 5°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm
Width	3.4 mm	3.4 mm	3.4 mm	3.2 mm	3.8 mm	3.8 mm	3.2 mm	3.4 mm	3.4 mm	3.4 mm
REF	790-106-00	790-106-00	790-134-00	790-132-00	790-130-00	790-131-00	790-133-00	790-135-00	790-106-00	790-106-00
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces



Torque	- 17°	- 11°	+ 7°	- 1°	- 1°	- 1°	- 1°	+ 7°	- 11°	- 17°
Angulation	0°	0°	+ 5°	0°	0°	0°	0°	+ 5°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	1.3 mm	1.3 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	3.3 mm	3.3 mm	3.1 mm	2.6 mm	2.6 mm	2.6 mm	2.6 mm	3.1 mm	3.3 mm	3.3 mm
REF	790-110-00	790-128-00	790-138-00	790-107-00	790-107-00	790-107-00	790-107-00	790-139-00	790-128-00	790-110-00
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces


Assortments:

1 case	20 brackets	REF 891-375-00
10 cases	200 brackets	REF 891-376-00
50 cases	1000 brackets	REF 891-377-00

Ricketts® Universal 18

Slot: 0.46 mm x 0.76 mm / 18 x 30 - Cuspid brackets with hook

Torque	0°	0°	+ 7°	+ 14°	+ 22°	+ 22°	+ 14°	+ 7°	0°	0°
Angulation	0°	0°	+ 5°	+ 8°	+ 5°	+ 5°	+ 8°	+ 5°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm
Width	3.3 mm	3.3 mm	3.4 mm	3.2 mm	3.8 mm	3.8 mm	3.2 mm	3.4 mm	3.3 mm	3.3 mm
REF	790-136-00	790-136-00	790-134-00	790-132-00	790-130-00	790-131-00	790-133-00	790-135-00	790-136-00	790-136-00
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces



Torque	- 14°	0°	+ 7°	0°	0°	0°	0°	+ 7°	0°	- 14°
Angulation	0°	0°	+ 5°	0°	0°	0°	0°	+ 5°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm
Width	3.3 mm	3.3 mm	3.1 mm	2.6 mm	2.6 mm	2.6 mm	2.6 mm	3.1 mm	3.3 mm	3.3 mm
REF	790-140-00	790-136-00	790-138-00	790-137-00	790-137-00	790-137-00	790-137-00	790-139-00	790-136-00	790-140-00
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces


Assortments:

1 case	20 brackets	REF 891-188-00
10 cases	200 brackets	REF 891-189-00
50 cases	1000 brackets	REF 891-196-00

Hilgers 18

Slot: 0.46 mm x 0.76 mm / 18 x 30 · Cuspid brackets with hook

Torque	- 7°	- 7°	+ 7°	+ 14°	+ 22°	+ 22°	+ 14°	+ 7°	- 7°	- 7°
Angulation	0°	0°	+ 5°	+ 8°	+ 5°	+ 5°	+ 8°	+ 5°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm
Width	3.4 mm	3.4 mm	3.4 mm	3.2 mm	3.8 mm	3.8 mm	3.2 mm	3.4 mm	3.4 mm	3.4 mm
REF	790-106-00	790-106-00	790-134-00	790-132-00	790-130-00	790-131-00	790-133-00	790-135-00	790-106-00	790-106-00
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces



Torque	- 17°	- 11°	+ 7°	- 1°	- 1°	- 1°	- 1°	+ 7°	- 11°	- 17°
Angulation	0°	0°	+ 5°	0°	0°	0°	0°	+ 5°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	1.3 mm	1.3 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	3.3 mm	3.3 mm	3.1 mm	2.6 mm	2.6 mm	2.6 mm	2.6 mm	3.1 mm	3.3 mm	3.3 mm
REF	790-110-00	790-128-00	790-138-00	790-107-00	790-107-00	790-107-00	790-107-00	790-139-00	790-128-00	790-110-00
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Assortments:


1 case	20 brackets	REF 891-261-00
10 cases	200 brackets	REF 891-262-00
50 cases	1000 brackets	REF 891-263-00

Standard Edgewise 18

Slot: 0.46 mm x 0.76 mm / 18 x 30 · without hook

Torque	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
Angulation	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm
Width	2.8 mm	2.8 mm	2.8 mm	2.6 mm	2.6 mm	3.0 mm	3.0 mm	2.6 mm	2.8 mm	2.8 mm
REF	718-205-15	718-205-14	718-205-13	718-205-12	718-205-11	718-205-21	718-205-22	718-205-23	718-205-24	718-205-25
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Even easier bracket identification with FDI-marking on the base!



Torque	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
Angulation	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm
Width	2.8 mm	2.8 mm	2.8 mm	2.1 mm	2.1 mm	2.1 mm	2.1 mm	2.8 mm	2.8 mm	2.8 mm
REF	718-205-45	718-205-44	718-205-43	718-205-42	718-205-41	718-205-31	718-205-32	718-205-33	718-205-34	718-205-35
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Assortments:

1 case	20 brackets	REF 718-205-01
10 cases	200 brackets	REF 718-205-10
50 cases	1000 brackets	REF 718-205-50


discovery® / Assortments

Standard Edgewise 22

Slot: 0.56 mm x 0.76 mm / 22 x 30 - without hook

Torque	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
Angulation	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm
Width	2.8 mm	2.8 mm	2.8 mm	2.6 mm	3.0 mm	3.0 mm	2.6 mm	2.8 mm	2.8 mm	2.8 mm
REF	722-205-15	722-205-14	722-205-13	722-205-12	722-205-11	722-205-21	722-205-22	722-205-23	722-205-24	722-205-25
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Even easier bracket identification with FDI-marking on the base!



Torque	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
Angulation	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm
Width	2.8 mm	2.8 mm	2.8 mm	2.1 mm	2.1 mm	2.1 mm	2.1 mm	2.8 mm	2.8 mm	2.8 mm
REF	722-205-45	722-205-44	722-205-43	722-205-42	722-205-41	722-205-31	722-205-32	722-205-33	722-205-34	722-205-35
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Assortments:

1 case	20 brackets	REF 722-205-01
10 cases	200 brackets	REF 722-205-10
50 cases	1000 brackets	REF 722-205-50

discovery® / Single packs – bonding technique

Roth 18/22

Tooth	Torque	Angulation	In/Out	Hook	Width	REF (18)	REF (22)	Quantity
11	+ 12°	+ 5°	0.7 mm		3.8 mm	790-100-00	790-112-00	10 pieces
21	+ 12°	+ 5°	0.7 mm		3.8 mm	790-101-00	790-113-00	10 pieces
13	- 2°	+ 9°	0.7 mm	•	3.4 mm	790-155-00	790-157-00	10 pieces
23	- 2°	+ 9°	0.7 mm	•	3.4 mm	790-156-00	790-158-00	10 pieces



Slot 18: 0.46 mm x 0.76 mm / 18 x 30
Slot 22: 0.56 mm x 0.76 mm / 22 x 30



Metal brackets


discovery® smart / Assortments

Roth 18

Slot: 0.46 mm x 0.76 mm / 18 x 30 · Cuspid brackets with hook

Torque	- 7°	- 7°	- 2°	+ 8°	+ 12°	+ 12°	+ 8°	- 2°	- 7°	- 7°
Angulation	0°	0°	+ 13°	+ 9°	+ 5°	+ 5°	+ 9°	+ 13°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	2.8 mm	2.8 mm	2.8 mm	2.4 mm	2.8 mm	2.8 mm	2.4 mm	2.8 mm	2.8 mm	2.8 mm
REF	718-901-15	718-901-14	718-941-13	718-901-12	718-901-11	718-901-21	718-901-22	718-941-23	718-901-24	718-901-25
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Even easier bracket identification with FDI-marking on the base!



Torque	- 22°	- 17°	- 11°	- 1°	- 1°	- 1°	- 1°	- 11°	- 17°	- 22°
Angulation	0°	0°	+ 7°	0°	0°	0°	0°	+ 7°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	1.3 mm	1.3 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	2.8 mm	2.8 mm	2.8 mm	2.3 mm	2.3 mm	2.3 mm	2.3 mm	2.8 mm	2.8 mm	2.8 mm
REF	718-901-45	718-901-44	718-941-43	718-901-42	718-901-41	718-901-31	718-901-32	718-941-33	718-901-34	718-901-35
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Assortments:


1 case	20 brackets	REF 718-901-01
10 cases	200 brackets	REF 718-901-10
50 cases	1000 brackets	REF 718-901-50

Roth 18

Slot: 0.46 mm x 0.76 mm / 18 x 30 · Cuspid, 1st bicuspid and 2nd bicuspid brackets with hook

Torque	- 7°	- 7°	- 2°	+ 8°	+ 12°	+ 12°	+ 8°	- 2°	- 7°	- 7°
Angulation	0°	0°	+ 13°	+ 9°	+ 5°	+ 5°	+ 9°	+ 13°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	2.8 mm	2.8 mm	2.8 mm	2.4 mm	2.8 mm	2.8 mm	2.4 mm	2.8 mm	2.8 mm	2.8 mm
REF	718-941-15	718-941-14	718-941-13	718-901-12	718-901-11	718-901-21	718-901-22	718-941-23	718-941-24	718-941-25
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Even easier bracket identification with FDI-marking on the base!



Torque	- 22°	- 17°	- 11°	- 1°	- 1°	- 1°	- 1°	- 11°	- 17°	- 22°
Angulation	0°	0°	+ 7°	0°	0°	0°	0°	+ 7°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	1.3 mm	1.3 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	2.8 mm	2.8 mm	2.8 mm	2.3 mm	2.3 mm	2.3 mm	2.3 mm	2.8 mm	2.8 mm	2.8 mm
REF	718-941-45	718-941-44	718-941-43	718-901-42	718-901-41	718-901-31	718-901-32	718-941-33	718-941-34	718-941-35
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Assortments:


1 case	20 brackets	REF 718-941-01
10 cases	200 brackets	REF 718-941-10
50 cases	1000 brackets	REF 718-941-50

Roth 22

Slot: 0.56 mm x 0.76 mm / 22 x 30 · Cuspid brackets with hook

Torque	- 7°	- 7°	- 2°	+ 8°	+ 12°	+ 12°	+ 8°	- 2°	- 7°	- 7°
Angulation	0°	0°	+ 13°	+ 9°	+ 5°	+ 5°	+ 9°	+ 13°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	2.8 mm	2.8 mm	2.8 mm	2.4 mm	2.8 mm	2.8 mm	2.4 mm	2.8 mm	2.8 mm	2.8 mm
REF	722-901-15	722-901-14	722-941-13	722-901-12	722-901-11	722-901-21	722-901-22	722-941-23	722-901-24	722-901-25
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Even easier bracket identification with FDI-marking on the base!



Torque	- 22°	- 17°	- 11°	- 1°	- 1°	- 1°	- 1°	- 11°	- 17°	- 22°
Angulation	0°	0°	+ 7°	0°	0°	0°	0°	+ 7°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	1.3 mm	1.3 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	2.8 mm	2.8 mm	2.8 mm	2.3 mm	2.3 mm	2.3 mm	2.3 mm	2.8 mm	2.8 mm	2.8 mm
REF	722-901-45	722-901-44	722-941-43	722-901-42	722-901-41	722-901-31	722-901-32	722-941-33	722-901-34	722-901-35
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Assortments:


1 case	20 brackets	REF 722-901-01
10 cases	200 brackets	REF 722-901-10
50 cases	1000 brackets	REF 722-901-50

Roth 22

Slot: 0.56 mm x 0.76 mm / 22 x 30 · Cuspid, 1st bicuspid and 2nd bicuspid brackets with hook

Torque	- 7°	- 7°	- 2°	+ 8°	+ 12°	+ 12°	+ 8°	- 2°	- 7°	- 7°
Angulation	0°	0°	+ 13°	+ 9°	+ 5°	+ 5°	+ 9°	+ 13°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	2.8 mm	2.8 mm	2.8 mm	2.4 mm	2.8 mm	2.8 mm	2.4 mm	2.8 mm	2.8 mm	2.8 mm
REF	722-941-15	722-941-14	722-941-13	722-901-12	722-901-11	722-901-21	722-901-22	722-941-23	722-941-24	722-941-25
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Even easier bracket identification with FDI-marking on the base!



Torque	- 22°	- 17°	- 11°	- 1°	- 1°	- 1°	- 1°	- 11°	- 17°	- 22°
Angulation	0°	0°	+ 7°	0°	0°	0°	0°	+ 7°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	1.3 mm	1.3 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	2.8 mm	2.8 mm	2.8 mm	2.3 mm	2.3 mm	2.3 mm	2.3 mm	2.8 mm	2.8 mm	2.8 mm
REF	722-941-45	722-941-44	722-941-43	722-901-42	722-901-41	722-901-31	722-901-32	722-941-33	722-941-34	722-941-35
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Assortments:


1 case	20 brackets	REF 722-941-01
10 cases	200 brackets	REF 722-941-10
50 cases	1000 brackets	REF 722-941-50

McLaughlin-Bennett-Trevisi** 18

Slot: 0.46 mm x 0.76 mm / 18 x 30 · Cuspid brackets with hook

Torque	- 7°	- 7°	0°	+ 10°	+ 17°	+ 17°	+ 10°	0°	- 7°	- 7°
Angulation	0°	0°	+ 8°	+ 8°	+ 4°	+ 4°	+ 8°	+ 8°	0°	0°
In/Out	1.2 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	1.2 mm
Width	2.8 mm	2.8 mm	2.8 mm	2.4 mm	2.8 mm	2.8 mm	2.4 mm	2.8 mm	2.8 mm	2.8 mm
REF	718-902-15	718-902-14	718-942-13	718-902-12	718-902-11	718-902-21	718-902-22	718-942-23	718-902-24	718-902-25
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Even easier bracket identification with FDI-marking on the base!



Torque	- 17°	- 12°	0°	- 6°	- 6°	- 6°	- 6°	0°	- 12°	- 17°
Angulation	+ 2°	+ 2°	+ 3°	0°	0°	0°	0°	+ 3°	+ 2°	+ 2°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	1.3 mm	1.3 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	2.8 mm	2.8 mm	2.8 mm	2.3 mm	2.3 mm	2.3 mm	2.3 mm	2.8 mm	2.8 mm	2.8 mm
REF	718-902-45	718-902-44	718-942-43	718-902-42	718-902-41	718-902-31	718-902-32	718-942-33	718-902-34	718-902-35
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Assortments:

1 case	20 brackets	REF 718-902-01
10 cases	200 brackets	REF 718-902-10
50 cases	1000 brackets	REF 718-902-50


** The Dentaaurum version of this prescription is not claimed to be a duplication of any other, nor does Dentaaurum imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.

McLaughlin-Bennett-Trevisi** 18

Slot: 0.46 mm x 0.76 mm / 18 x 30 · Cuspid, 1st bicuspid and 2nd bicuspid brackets with hook

Torque	- 7°	- 7°	0°	+ 10°	+ 17°	+ 17°	+ 10°	0°	- 7°	- 7°
Angulation	0°	0°	+ 8°	+ 8°	+ 4°	+ 4°	+ 8°	+ 8°	0°	0°
In/Out	1.2 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	1.2 mm
Width	2.8 mm	2.8 mm	2.8 mm	2.4 mm	2.8 mm	2.8 mm	2.4 mm	2.8 mm	2.8 mm	2.8 mm
REF	718-942-15	718-942-14	718-942-13	718-902-12	718-902-11	718-902-21	718-902-22	718-942-23	718-942-24	718-942-25
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Even easier bracket identification with FDI-marking on the base!



Torque	- 17°	- 12°	0°	- 6°	- 6°	- 6°	- 6°	0°	- 12°	- 17°
Angulation	+ 2°	+ 2°	+ 3°	0°	0°	0°	0°	+ 3°	+ 2°	+ 2°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	1.3 mm	1.3 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	2.8 mm	2.8 mm	2.8 mm	2.3 mm	2.3 mm	2.3 mm	2.3 mm	2.8 mm	2.8 mm	2.8 mm
REF	718-942-45	718-942-44	718-942-43	718-902-42	718-902-41	718-902-31	718-902-32	718-942-33	718-942-34	718-942-35
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Assortments:

1 case	20 brackets	REF 718-942-01
10 cases	200 brackets	REF 718-942-10
50 cases	1000 brackets	REF 718-942-50


** The Dentaaurum version of this prescription is not claimed to be a duplication of any other, nor does Dentaaurum imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.

McLaughlin-Bennett-Trevisi** 22

Slot: 0.56 mm x 0.76 mm / 22 x 30 · Cuspid brackets with hook

Torque	- 7°	- 7°	0°	+ 10°	+ 17°	+ 17°	+ 10°	0°	- 7°	- 7°
Angulation	0°	0°	+ 8°	+ 8°	+ 4°	+ 4°	+ 8°	+ 8°	0°	0°
In/Out	1.2 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	1.2 mm
Width	2.8 mm	2.8 mm	2.8 mm	2.4 mm	2.8 mm	2.8 mm	2.4 mm	2.8 mm	2.8 mm	2.8 mm
REF	722-902-15	722-902-14	722-942-13	722-902-12	722-902-11	722-902-21	722-902-22	722-942-23	722-902-24	722-902-25
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Even easier bracket identification with FDI-marking on the base!



Torque	- 17°	- 12°	0°	- 6°	- 6°	- 6°	- 6°	0°	- 12°	- 17°
Angulation	+ 2°	+ 2°	+ 3°	0°	0°	0°	0°	+ 3°	+ 2°	+ 2°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	1.3 mm	1.3 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	2.8 mm	2.8 mm	2.8 mm	2.3 mm	2.3 mm	2.3 mm	2.3 mm	2.8 mm	2.8 mm	2.8 mm
REF	722-902-45	722-902-44	722-942-43	722-902-42	722-902-41	722-902-31	722-902-32	722-942-33	722-902-34	722-902-35
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Assortments:

1 case	20 brackets	REF 722-902-01
10 cases	200 brackets	REF 722-902-10
50 cases	1000 brackets	REF 722-902-50


** The Dentaureum version of this prescription is not claimed to be a duplication of any other, nor does Dentaureum imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.

McLaughlin-Bennett-Trevisi** 22

Slot: 0.56 mm x 0.76 mm / 22 x 30 · Cuspid, 1st bicuspid and 2nd bicuspid brackets with hook

Torque	- 7°	- 7°	0°	+ 10°	+ 17°	+ 17°	+ 10°	0°	- 7°	- 7°
Angulation	0°	0°	+ 8°	+ 8°	+ 4°	+ 4°	+ 8°	+ 8°	0°	0°
In/Out	1.2 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	1.2 mm
Width	2.8 mm	2.8 mm	2.8 mm	2.4 mm	2.8 mm	2.8 mm	2.4 mm	2.8 mm	2.8 mm	2.8 mm
REF	722-942-15	722-942-14	722-942-13	722-902-12	722-902-11	722-902-21	722-902-22	722-942-23	722-942-24	722-942-25
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Even easier bracket identification with FDI-marking on the base!



Torque	- 17°	- 12°	0°	- 6°	- 6°	- 6°	- 6°	0°	- 12°	- 17°
Angulation	+ 2°	+ 2°	+ 3°	0°	0°	0°	0°	+ 3°	+ 2°	+ 2°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	1.3 mm	1.3 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	2.8 mm	2.8 mm	2.8 mm	2.3 mm	2.3 mm	2.3 mm	2.3 mm	2.8 mm	2.8 mm	2.8 mm
REF	722-942-45	722-942-44	722-942-43	722-902-42	722-902-41	722-902-31	722-902-32	722-942-33	722-942-34	722-942-35
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Assortments:

1 case	20 brackets	REF 722-942-01
10 cases	200 brackets	REF 722-942-10
50 cases	1000 brackets	REF 722-942-50


** The Dentaureum version of this prescription is not claimed to be a duplication of any other, nor does Dentaureum imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.

topic / Assortments

Roth 18

Slot: 0.46 mm x 0.76 mm / 18 x 30 · Cuspid brackets with hook · nickel-free*

Torque	- 7°	- 7°	- 2°	+ 8°	+ 12°	+ 12°	+ 8°	- 2°	- 7°	- 7°
Angulation	0°	0°	+ 13°	+ 9°	+ 5°	+ 5°	+ 9°	+ 13°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	3.4 mm	3.4 mm	3.4 mm	3.2 mm	3.4 mm	3.4 mm	3.2 mm	3.4 mm	3.4 mm	3.4 mm
REF	790-306-00	790-306-00	790-304-00	790-302-00	790-363-00	790-364-00	790-303-00	790-305-00	790-306-00	790-306-00
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces



Torque	- 22°	- 17°	- 11°	- 1°	- 1°	- 1°	- 1°	- 11°	- 17°	- 22°
Angulation	0°	0°	+ 7°	0°	0°	0°	0°	+ 7°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	1.3 mm	1.3 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	3.4 mm	3.4 mm	3.1 mm	2.7 mm	2.7 mm	2.7 mm	2.7 mm	3.1 mm	3.4 mm	3.4 mm
REF	790-311-00	790-310-00	790-308-00	790-307-00	790-307-00	790-307-00	790-307-00	790-309-00	790-310-00	790-311-00
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Assortments:

1 case	20 brackets	REF 891-333-00
10 cases	200 brackets	REF 891-334-00
50 cases	1000 brackets	REF 891-335-00

Note:


Other brackets with Roth values see single packs.

* see Materials list Orthodontics

Roth 18

Slot: 0.46 mm x 0.76 mm / 18 x 30 · Cuspid, 1st bicuspid and 2nd bicuspid brackets with hook · nickel-free*

Torque	- 7°	- 7°	- 2°	+ 8°	+ 12°	+ 12°	+ 8°	- 2°	- 7°	- 7°
Angulation	0°	0°	+ 13°	+ 9°	+ 5°	+ 5°	+ 9°	+ 13°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	3.4 mm	3.4 mm	3.4 mm	3.2 mm	3.4 mm	3.4 mm	3.2 mm	3.4 mm	3.4 mm	3.4 mm
REF	790-341-00	790-341-00	790-304-00	790-302-00	790-363-00	790-364-00	790-303-00	790-305-00	790-342-00	790-342-00
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces



Torque	- 22°	- 17°	- 11°	- 1°	- 1°	- 1°	- 1°	- 11°	- 17°	- 22°
Angulation	0°	0°	+ 7°	0°	0°	0°	0°	+ 7°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	1.3 mm	1.3 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	3.4 mm	3.4 mm	3.1 mm	2.7 mm	2.7 mm	2.7 mm	2.7 mm	3.1 mm	3.4 mm	3.4 mm
REF	790-345-00	790-343-00	790-308-00	790-307-00	790-307-00	790-307-00	790-307-00	790-309-00	790-344-00	790-346-00
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Assortments:

1 case	20 brackets	REF 891-336-00
10 cases	200 brackets	REF 891-337-00
50 cases	1000 brackets	REF 891-338-00

Note:


Other brackets with Roth values see single packs.

* see Materials list Orthodontics

Roth 22

Slot: 0.56 mm x 0.76 mm / 22 x 30 · Cuspid brackets with hook · nickel-free*

Torque	- 7°	- 7°	- 2°	+ 8°	+ 12°	+ 12°	+ 8°	- 2°	- 7°	- 7°
Angulation	0°	0°	+ 13°	+ 9°	+ 5°	+ 5°	+ 9°	+ 13°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	3.4 mm	3.4 mm	3.4 mm	3.2 mm	3.4 mm	3.4 mm	3.2 mm	3.4 mm	3.4 mm	3.4 mm
REF	790-318-00	790-318-00	790-316-00	790-314-00	790-365-00	790-366-00	790-315-00	790-317-00	790-318-00	790-318-00
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces



Torque	- 22°	- 17°	- 11°	- 1°	- 1°	- 1°	- 1°	- 11°	- 17°	- 22°
Angulation	0°	0°	+ 7°	0°	0°	0°	0°	+ 7°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	1.3 mm	1.3 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	3.4 mm	3.4 mm	3.1 mm	2.7 mm	2.7 mm	2.7 mm	2.7 mm	3.1 mm	3.4 mm	3.4 mm
REF	790-323-00	790-322-00	790-320-00	790-319-00	790-319-00	790-319-00	790-319-00	790-321-00	790-322-00	790-323-00
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Assortments:

1 case	20 brackets	REF 891-339-00
10 cases	200 brackets	REF 891-340-00
50 cases	1000 brackets	REF 891-341-00

Note:


Other brackets with Roth values see single packs.

* see Materials list Orthodontics

Roth 22

Slot: 0.56 mm x 0.76 mm / 22 x 30 · Cuspid, 1st bicuspid and 2nd bicuspid brackets with hook · nickel-free*

Torque	- 7°	- 7°	- 2°	+ 8°	+ 12°	+ 12°	+ 8°	- 2°	- 7°	- 7°
Angulation	0°	0°	+ 13°	+ 9°	+ 5°	+ 5°	+ 9°	+ 13°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	3.4 mm	3.4 mm	3.4 mm	3.2 mm	3.4 mm	3.4 mm	3.2 mm	3.4 mm	3.4 mm	3.4 mm
REF	790-347-00	790-347-00	790-316-00	790-314-00	790-365-00	790-366-00	790-315-00	790-317-00	790-348-00	790-348-00
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces



Torque	- 22°	- 17°	- 11°	- 1°	- 1°	- 1°	- 1°	- 11°	- 17°	- 22°
Angulation	0°	0°	+ 7°	0°	0°	0°	0°	+ 7°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	1.3 mm	1.3 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	3.4 mm	3.4 mm	3.1 mm	2.7 mm	2.7 mm	2.7 mm	2.7 mm	3.1 mm	3.4 mm	3.4 mm
REF	790-351-00	790-349-00	790-320-00	790-319-00	790-319-00	790-319-00	790-319-00	790-321-00	790-350-00	790-352-00
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Assortments:

1 case	20 brackets	REF 891-342-00
10 cases	200 brackets	REF 891-343-00
50 cases	1000 brackets	REF 891-344-00

Note:

Other brackets with Roth values see single packs.

* see Materials list Orthodontics

topic / Single packs – bonding technique



Roth 18/22

nickel-free*

Tooth	Torque	Angulation	In/Out	Hook	Width	REF (18)	REF (22)	Quantity
11	+ 12°	+ 5°	0.7 mm		3.9 mm	790-300-00	790-312-00	10 pieces
21	+ 12°	+ 5°	0.7 mm		3.9 mm	790-301-00	790-313-00	10 pieces

Slot 18: 0.46 mm x 0.76 mm / 18 x 30

Slot 22: 0.56 mm x 0.76 mm / 22 x 30


* see Materials list Orthodontics

Roth 18

Slot: 0.46 mm x 0.76 mm / 18 x 30 · Cuspid brackets with hook

Torque	- 7°	- 7°	- 2°	+ 8°	+ 12°	+ 12°	+ 8°	- 2°	- 7°	- 7°
Angulation	0°	0°	+ 13°	+ 9°	+ 5°	+ 5°	+ 9°	+ 13°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	2.4 mm	2.4 mm	2.4 mm	2.2 mm	2.4 mm	2.4 mm	2.2 mm	2.4 mm	2.4 mm	2.4 mm
REF	718-401-15	718-401-14	718-441-13	718-401-12	718-401-11	718-401-21	718-401-22	718-441-23	718-401-24	718-401-25
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Even easier bracket identification with FDI-marking on the base!



Torque	- 22°	- 17°	- 11°	- 1°	- 1°	- 1°	- 1°	- 11°	- 17°	- 22°
Angulation	0°	0°	+ 7°	0°	0°	0°	0°	+ 7°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	1.3 mm	1.3 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	2.4 mm	2.4 mm	2.4 mm	2.2 mm	2.2 mm	2.2 mm	2.2 mm	2.4 mm	2.4 mm	2.4 mm
REF	718-401-45	718-401-44	718-441-43	718-401-42	718-401-41	718-401-31	718-401-32	718-441-33	718-401-34	718-401-35
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Assortments:


1 case	20 brackets	REF 718-401-01
10 cases	200 brackets	REF 718-401-10
50 cases	1000 brackets	REF 718-401-50

Roth 18

Slot: 0.46 mm x 0.76 mm / 18 x 30 · Cuspid, 1st bicuspid and 2nd bicuspid brackets with hook

Torque	- 7°	- 7°	- 2°	+ 8°	+ 12°	+ 12°	+ 8°	- 2°	- 7°	- 7°
Angulation	0°	0°	+ 13°	+ 9°	+ 5°	+ 5°	+ 9°	+ 13°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	2.4 mm	2.4 mm	2.4 mm	2.2 mm	2.4 mm	2.4 mm	2.2 mm	2.4 mm	2.4 mm	2.4 mm
REF	718-441-15	718-441-14	718-441-13	718-401-12	718-401-11	718-401-21	718-401-22	718-441-23	718-441-24	718-441-25
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Even easier bracket identification with FDI-marking on the base!



Torque	- 22°	- 17°	- 11°	- 1°	- 1°	- 1°	- 1°	- 11°	- 17°	- 22°
Angulation	0°	0°	+ 7°	0°	0°	0°	0°	+ 7°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	1.3 mm	1.3 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	2.4 mm	2.4 mm	2.4 mm	2.2 mm	2.2 mm	2.2 mm	2.2 mm	2.4 mm	2.4 mm	2.4 mm
REF	718-441-45	718-441-44	718-441-43	718-401-42	718-401-41	718-401-31	718-401-32	718-441-33	718-441-34	718-441-35
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Assortments:

1 case	20 brackets	REF 718-441-01
10 cases	200 brackets	REF 718-441-10
50 cases	1000 brackets	REF 718-441-50

Metal brackets


equilibrium® mini / Assortments

Roth 22

Slot: 0.56 mm x 0.76 mm / 22 x 30 · Cuspid brackets with hook

Torque	- 7°	- 7°	- 2°	+ 8°	+ 12°	+ 12°	+ 8°	- 2°	- 7°	- 7°
Angulation	0°	0°	+ 13°	+ 9°	+ 5°	+ 5°	+ 9°	+ 13°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	2.4 mm	2.4 mm	2.4 mm	2.2 mm	2.4 mm	2.4 mm	2.2 mm	2.4 mm	2.4 mm	2.4 mm
REF	722-401-15	722-401-14	722-441-13	722-401-12	722-401-11	722-401-21	722-401-22	722-441-23	722-401-24	722-401-25
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Even easier bracket identification with FDI-marking on the base!



Torque	- 22°	- 17°	- 11°	- 1°	- 1°	- 1°	- 1°	- 11°	- 17°	- 22°
Angulation	0°	0°	+ 7°	0°	0°	0°	0°	+ 7°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	1.3 mm	1.3 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	2.4 mm	2.4 mm	2.4 mm	2.2 mm	2.2 mm	2.2 mm	2.2 mm	2.4 mm	2.4 mm	2.4 mm
REF	722-401-45	722-401-44	722-441-43	722-401-42	722-401-41	722-401-31	722-401-32	722-441-33	722-401-34	722-401-35
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Assortments:


1 case	20 brackets	REF 722-401-01
10 cases	200 brackets	REF 722-401-10
50 cases	1000 brackets	REF 722-401-50

Roth 22

Slot: 0.56 mm x 0.76 mm / 22 x 30 · Cuspid, 1st bicuspid and 2nd bicuspid brackets with hook

Torque	- 7°	- 7°	- 2°	+ 8°	+ 12°	+ 12°	+ 8°	- 2°	- 7°	- 7°
Angulation	0°	0°	+ 13°	+ 9°	+ 5°	+ 5°	+ 9°	+ 13°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	2.4 mm	2.4 mm	2.4 mm	2.2 mm	2.4 mm	2.4 mm	2.2 mm	2.4 mm	2.4 mm	2.4 mm
REF	722-441-15	722-441-14	722-441-13	722-401-12	722-401-11	722-401-21	722-401-22	722-441-23	722-441-24	722-441-25
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Even easier bracket identification with FDI-marking on the base!



Torque	- 22°	- 17°	- 11°	- 1°	- 1°	- 1°	- 1°	- 11°	- 17°	- 22°
Angulation	0°	0°	+ 7°	0°	0°	0°	0°	+ 7°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	1.3 mm	1.3 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	2.4 mm	2.4 mm	2.4 mm	2.2 mm	2.2 mm	2.2 mm	2.2 mm	2.4 mm	2.4 mm	2.4 mm
REF	722-441-45	722-441-44	722-441-43	722-401-42	722-401-41	722-401-31	722-401-32	722-441-33	722-441-34	722-441-35
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Assortments:


1 case	20 brackets	REF 722-441-01
10 cases	200 brackets	REF 722-441-10
50 cases	1000 brackets	REF 722-441-50

McLaughlin-Bennett-Trevisi** 18

Slot: 0.46 mm x 0.76 mm / 18 x 30 · Cuspid brackets with hook

Torque	- 7°	- 7°	0°	+ 10°	+ 17°	+ 17°	+ 10°	0°	- 7°	- 7°
Angulation	0°	0°	+ 8°	+ 8°	+ 4°	+ 4°	+ 8°	+ 8°	0°	0°
In/Out	1.2 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	1.2 mm
Width	2.4 mm	2.4 mm	2.4 mm	2.2 mm	2.4 mm	2.4 mm	2.2 mm	2.4 mm	2.4 mm	2.4 mm
REF	718-402-15	718-402-14	718-442-13	718-402-12	718-402-11	718-402-21	718-402-22	718-442-23	718-402-24	718-402-25
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Even easier bracket identification with FDI-marking on the base!



Torque	- 17°	- 12°	0°	- 6°	- 6°	- 6°	- 6°	0°	- 12°	- 17°
Angulation	+ 2°	+ 2°	+ 3°	0°	0°	0°	0°	+ 3°	+ 2°	+ 2°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	1.3 mm	1.3 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	2.4 mm	2.4 mm	2.4 mm	2.2 mm	2.2 mm	2.2 mm	2.2 mm	2.4 mm	2.4 mm	2.4 mm
REF	718-402-45	718-402-44	718-442-43	718-402-42	718-402-41	718-402-31	718-402-32	718-442-33	718-402-34	718-402-35
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Assortments:

1 case	20 brackets	REF 718-402-01
10 cases	200 brackets	REF 718-402-10
50 cases	1000 brackets	REF 718-402-50

Note:

Other brackets with McLaughlin-Bennett-Trevisi** values see single packs.


** The Dentaurem version of this prescription is not claimed to be a duplication of any other, nor does Dentaurem imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.

McLaughlin-Bennett-Trevisi** 18

Slot: 0.46 mm x 0.76 mm / 18 x 30 · Cuspid, 1st bicuspid and 2nd bicuspid brackets with hook

Torque	- 7°	- 7°	0°	+ 10°	+ 17°	+ 17°	+ 10°	0°	- 7°	- 7°
Angulation	0°	0°	+ 8°	+ 8°	+ 4°	+ 4°	+ 8°	+ 8°	0°	0°
In/Out	1.2 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	1.2 mm
Width	2.4 mm	2.4 mm	2.4 mm	2.2 mm	2.4 mm	2.4 mm	2.2 mm	2.4 mm	2.4 mm	2.4 mm
REF	718-442-15	718-442-14	718-442-13	718-402-12	718-402-11	718-402-21	718-402-22	718-442-23	718-442-24	718-442-25
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Even easier bracket identification with FDI-marking on the base!



Torque	- 17°	- 12°	0°	- 6°	- 6°	- 6°	- 6°	0°	- 12°	- 17°
Angulation	+ 2°	+ 2°	+ 3°	0°	0°	0°	0°	+ 3°	+ 2°	+ 2°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	1.3 mm	1.3 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	2.4 mm	2.4 mm	2.4 mm	2.2 mm	2.2 mm	2.2 mm	2.2 mm	2.4 mm	2.4 mm	2.4 mm
REF	718-442-45	718-442-44	718-442-43	718-402-42	718-402-41	718-402-31	718-402-32	718-442-33	718-442-34	718-442-35
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Assortments:

1 case	20 brackets	REF 718-442-01
10 cases	200 brackets	REF 718-442-10
50 cases	1000 brackets	REF 718-442-50

Note:

Other brackets with McLaughlin-Bennett-Trevisi** values see single packs.

** The Dentaurem version of this prescription is not claimed to be a duplication of any other, nor does Dentaurem imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.


equilibrium® mini / Assortments

McLaughlin-Bennett-Trevisi** 22

Slot: 0.56 mm x 0.76 mm / 22 x 30 · Cuspid brackets with hook

Torque	- 7°	- 7°	0°	+ 10°	+ 17°	+ 17°	+ 10°	0°	- 7°	- 7°
Angulation	0°	0°	+ 8°	+ 8°	+ 4°	+ 4°	+ 8°	+ 8°	0°	0°
In/Out	1.2 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	1.2 mm
Width	2.4 mm	2.4 mm	2.4 mm	2.2 mm	2.4 mm	2.4 mm	2.2 mm	2.4 mm	2.4 mm	2.4 mm
REF	722-402-15	722-402-14	722-442-13	722-402-12	722-402-11	722-402-21	722-402-22	722-442-23	722-402-24	722-402-25
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Even easier bracket identification with FDI-marking on the base!



Torque	- 17°	- 12°	0°	- 6°	- 6°	- 6°	- 6°	0°	- 12°	- 17°
Angulation	+ 2°	+ 2°	+ 3°	0°	0°	0°	0°	+ 3°	+ 2°	+ 2°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	1.3 mm	1.3 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	2.4 mm	2.4 mm	2.4 mm	2.2 mm	2.2 mm	2.2 mm	2.2 mm	2.4 mm	2.4 mm	2.4 mm
REF	722-402-45	722-402-44	722-442-43	722-402-42	722-402-41	722-402-31	722-402-32	722-442-33	722-402-34	722-402-35
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Assortments:

1 case	20 brackets	REF 722-402-01
10 cases	200 brackets	REF 722-402-10
50 cases	1000 brackets	REF 722-402-50

Note:

Other brackets with McLaughlin-Bennett-Trevisi** values see single packs.


** The Dentauro version of this prescription is not claimed to be a duplication of any other, nor does Dentauro imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.

McLaughlin-Bennett-Trevisi** 22

Slot: 0.56 mm x 0.76 mm / 22 x 30 · Cuspid, 1st bicuspid and 2nd bicuspid brackets with hook

Torque	- 7°	- 7°	0°	+ 10°	+ 17°	+ 17°	+ 10°	0°	- 7°	- 7°
Angulation	0°	0°	+ 8°	+ 8°	+ 4°	+ 4°	+ 8°	+ 8°	0°	0°
In/Out	1.2 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	1.2 mm
Width	2.4 mm	2.4 mm	2.4 mm	2.2 mm	2.4 mm	2.4 mm	2.2 mm	2.4 mm	2.4 mm	2.4 mm
REF	722-442-15	722-442-14	722-442-13	722-402-12	722-402-11	722-402-21	722-402-22	722-442-23	722-442-24	722-442-25
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Even easier bracket identification with FDI-marking on the base!



Torque	- 17°	- 12°	0°	- 6°	- 6°	- 6°	- 6°	0°	- 12°	- 17°
Angulation	+ 2°	+ 2°	+ 3°	0°	0°	0°	0°	+ 3°	+ 2°	+ 2°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	1.3 mm	1.3 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	2.4 mm	2.4 mm	2.4 mm	2.2 mm	2.2 mm	2.2 mm	2.2 mm	2.4 mm	2.4 mm	2.4 mm
REF	722-442-45	722-442-44	722-442-43	722-402-42	722-402-41	722-402-31	722-402-32	722-442-33	722-442-34	722-442-35
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Assortments:

1 case	20 brackets	REF 722-442-01
10 cases	200 brackets	REF 722-442-10
50 cases	1000 brackets	REF 722-442-50

Note:

Other brackets with McLaughlin-Bennett-Trevisi** values see single packs.

** The Dentauro version of this prescription is not claimed to be a duplication of any other, nor does Dentauro imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.

equilibrium® mini / Single packs – bonding technique

McLaughlin-Bennett-Trevisi** 18/22

Tooth	Torque	Angulation	In/Out	Hook	Width	REF (18)	REF (22)	Quantity
11	+ 17°	+ 4°	0.7 mm		2.8 mm		722-462-11	10 pieces
21	+ 17°	+ 4°	0.7 mm		2.8 mm		722-462-21	10 pieces
12	+ 10°	+ 8°	1.3 mm		2.4 mm		722-462-12	10 pieces
22	+ 10°	+ 8°	1.3 mm		2.4 mm		722-462-22	10 pieces
15	- 7°	0°	0.7 mm		2.4 mm	718-452-15		10 pieces
25	- 7°	0°	0.7 mm		2.4 mm	718-452-25		10 pieces
15	- 7°	0°	0.7 mm	•	2.4 mm	718-462-15		10 pieces
25	- 7°	0°	0.7 mm	•	2.4 mm	718-462-25		10 pieces



Even easier bracket identification with FDI-marking on the base!

Slot 18: 0.46 mm x 0.76 mm / 18 x 30

Slot 22: 0.56 mm x 0.76 mm / 22 x 30

** The Dentaaurum version of this prescription is not claimed to be a duplication of any other, nor does Dentaaurum imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.




equilibrium® 2 / Assortments

Roth 18

Slot: 0.46 mm x 0.76 mm / 18 x 30 · Cuspid brackets with hook

Torque	- 7°	- 7°	- 2°	+ 8°	+ 12°	+ 12°	+ 8°	- 2°	- 7°	- 7°
Angulation	0°	0°	+ 13°	+ 9°	+ 5°	+ 5°	+ 9°	+ 13°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	3.1 mm	3.1 mm	3.1 mm	3.0 mm	3.4 mm	3.4 mm	3.0 mm	3.1 mm	3.1 mm	3.1 mm
REF	718-301-15	718-301-14	718-341-13	718-301-12	718-301-11	718-301-21	718-301-22	718-341-23	718-301-24	718-301-25
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Even easier bracket identification with FDI-marking on the base!



Torque	- 22°	- 17°	- 11°	- 1°	- 1°	- 1°	- 1°	- 11°	- 17°	- 22°
Angulation	0°	0°	+ 7°	0°	0°	0°	0°	+ 7°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	1.3 mm	1.3 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	3.1 mm	3.1 mm	3.1 mm	2.4 mm	2.4 mm	2.4 mm	2.4 mm	3.1 mm	3.1 mm	3.1 mm
REF	718-301-45	718-301-44	718-341-43	718-301-42	718-301-41	718-301-31	718-301-32	718-341-33	718-301-34	718-301-35
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Assortments:

1 case	20 brackets	REF 718-301-01
10 cases	200 brackets	REF 718-301-10
50 cases	1000 brackets	REF 718-301-50

Note:


Other brackets with Roth values see single packs.

Roth 18

Slot: 0.46 mm x 0.76 mm / 18 x 30 · Cuspid, 1st bicuspid and 2nd bicuspid brackets with hook

Torque	- 7°	- 7°	- 2°	+ 8°	+ 12°	+ 12°	+ 8°	- 2°	- 7°	- 7°
Angulation	0°	0°	+ 13°	+ 9°	+ 5°	+ 5°	+ 9°	+ 13°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	3.1 mm	3.1 mm	3.1 mm	3.0 mm	3.4 mm	3.4 mm	3.0 mm	3.1 mm	3.1 mm	3.1 mm
REF	718-341-15	718-341-14	718-341-13	718-301-12	718-301-11	718-301-21	718-301-22	718-341-23	718-341-24	718-341-25
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Even easier bracket identification with FDI-marking on the base!



Torque	- 22°	- 17°	- 11°	- 1°	- 1°	- 1°	- 1°	- 11°	- 17°	- 22°
Angulation	0°	0°	+ 7°	0°	0°	0°	0°	+ 7°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	1.3 mm	1.3 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	3.1 mm	3.1 mm	3.1 mm	2.4 mm	2.4 mm	2.4 mm	2.4 mm	3.1 mm	3.1 mm	3.1 mm
REF	718-341-45	718-341-44	718-341-43	718-301-42	718-301-41	718-301-31	718-301-32	718-341-33	718-341-34	718-341-35
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Assortments:

1 case	20 brackets	REF 718-341-01
10 cases	200 brackets	REF 718-341-10
50 cases	1000 brackets	REF 718-341-50

Note:


Other brackets with Roth values see single packs.

Roth 22

Slot: 0.56 mm x 0.76 mm / 22 x 30 · Cuspid brackets with hook

Torque	- 7°	- 7°	- 2°	+ 8°	+ 12°	+ 12°	+ 8°	- 2°	- 7°	- 7°
Angulation	0°	0°	+ 13°	+ 9°	+ 5°	+ 5°	+ 9°	+ 13°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	3.1 mm	3.1 mm	3.1 mm	3.0 mm	3.4 mm	3.4 mm	3.0 mm	3.1 mm	3.1 mm	3.1 mm
REF	722-301-15	722-301-14	722-341-13	722-301-12	722-301-11	722-301-21	722-301-22	722-341-23	722-301-24	722-301-25
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Even easier bracket identification with FDI-marking on the base!



Torque	- 22°	- 17°	- 11°	- 1°	- 1°	- 1°	- 1°	- 11°	- 17°	- 22°
Angulation	0°	0°	+ 7°	0°	0°	0°	0°	+ 7°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	1.3 mm	1.3 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	3.1 mm	3.1 mm	3.1 mm	2.4 mm	2.4 mm	2.4 mm	2.4 mm	3.1 mm	3.1 mm	3.1 mm
REF	722-301-45	722-301-44	722-341-43	722-301-42	722-301-41	722-301-31	722-301-32	722-341-33	722-301-34	722-301-35
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Assortments:

1 case	20 brackets	REF 722-301-01
10 cases	200 brackets	REF 722-301-10
50 cases	1000 brackets	REF 722-301-50

Note:


Other brackets with Roth values see single packs.

Roth 22

Slot: 0.56 mm x 0.76 mm / 22 x 30 · Cuspid, 1st bicuspid and 2nd bicuspid brackets with hook

Torque	- 7°	- 7°	- 2°	+ 8°	+ 12°	+ 12°	+ 8°	- 2°	- 7°	- 7°
Angulation	0°	0°	+ 13°	+ 9°	+ 5°	+ 5°	+ 9°	+ 13°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	3.1 mm	3.1 mm	3.1 mm	3.0 mm	3.4 mm	3.4 mm	3.0 mm	3.1 mm	3.1 mm	3.1 mm
REF	722-341-15	722-341-14	722-341-13	722-301-12	722-301-11	722-301-21	722-301-22	722-341-23	722-341-24	722-341-25
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Even easier bracket identification with FDI-marking on the base!



Torque	- 22°	- 17°	- 11°	- 1°	- 1°	- 1°	- 1°	- 11°	- 17°	- 22°
Angulation	0°	0°	+ 7°	0°	0°	0°	0°	+ 7°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	1.3 mm	1.3 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	3.1 mm	3.1 mm	3.1 mm	2.4 mm	2.4 mm	2.4 mm	2.4 mm	3.1 mm	3.1 mm	3.1 mm
REF	722-341-45	722-341-44	722-341-43	722-301-42	722-301-41	722-301-31	722-301-32	722-341-33	722-341-34	722-341-35
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Assortments:

1 case	20 brackets	REF 722-341-01
10 cases	200 brackets	REF 722-341-10
50 cases	1000 brackets	REF 722-341-50

Note:

Other brackets with Roth values see single packs.


equilibrium® 2 / Assortments

McLaughlin-Bennett-Trevisi** 18

Slot: 0.46 mm x 0.76 mm / 18 x 30 · Cuspid brackets with hook

Torque	- 7°	- 7°	0°	+ 10°	+ 17°	+ 17°	+ 10°	0°	- 7°	- 7°
Angulation	0°	0°	+ 8°	+ 8°	+ 4°	+ 4°	+ 8°	+ 8°	0°	0°
In/Out	1.2 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	1.2 mm
Width	3.1 mm	3.1 mm	3.1 mm	3.0 mm	3.4 mm	3.4 mm	3.0 mm	3.1 mm	3.1 mm	3.1 mm
REF	718-302-15	718-302-14	718-342-13	718-302-12	718-302-11	718-302-21	718-302-22	718-342-23	718-302-24	718-302-25
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Even easier bracket identification with FDI-marking on the base!



Torque	- 17°	- 12°	0°	- 6°	- 6°	- 6°	- 6°	0°	- 12°	- 17°
Angulation	+ 2°	+ 2°	+ 3°	0°	0°	0°	0°	+ 3°	+ 2°	+ 2°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	1.3 mm	1.3 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	3.1 mm	3.1 mm	3.1 mm	2.4 mm	2.4 mm	2.4 mm	2.4 mm	3.1 mm	3.1 mm	3.1 mm
REF	718-302-45	718-302-44	718-342-43	718-302-42	718-302-41	718-302-31	718-302-32	718-342-33	718-302-34	718-302-35
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Assortments:

1 case	20 brackets	REF 718-302-01
10 cases	200 brackets	REF 718-302-10
50 cases	1000 brackets	REF 718-302-50


** The Dentaaurum version of this prescription is not claimed to be a duplication of any other, nor does Dentaaurum imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.

McLaughlin-Bennett-Trevisi** 18

Slot: 0.46 mm x 0.76 mm / 18 x 30 · Cuspid, 1st bicuspid and 2nd bicuspid brackets with hook

Torque	- 7°	- 7°	0°	+ 10°	+ 17°	+ 17°	+ 10°	0°	- 7°	- 7°
Angulation	0°	0°	+ 8°	+ 8°	+ 4°	+ 4°	+ 8°	+ 8°	0°	0°
In/Out	1.2 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	1.2 mm
Width	3.1 mm	3.1 mm	3.1 mm	3.0 mm	3.4 mm	3.4 mm	3.0 mm	3.1 mm	3.1 mm	3.1 mm
REF	718-342-15	718-342-14	718-342-13	718-302-12	718-302-11	718-302-21	718-302-22	718-342-23	718-342-24	718-342-25
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Even easier bracket identification with FDI-marking on the base!



Torque	- 17°	- 12°	0°	- 6°	- 6°	- 6°	- 6°	0°	- 12°	- 17°
Angulation	+ 2°	+ 2°	+ 3°	0°	0°	0°	0°	+ 3°	+ 2°	+ 2°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	1.3 mm	1.3 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	3.1 mm	3.1 mm	3.1 mm	2.4 mm	2.4 mm	2.4 mm	2.4 mm	3.1 mm	3.1 mm	3.1 mm
REF	718-342-45	718-342-44	718-342-43	718-302-42	718-302-41	718-302-31	718-302-32	718-342-33	718-342-34	718-342-35
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Assortments:

1 case	20 brackets	REF 718-342-01
10 cases	200 brackets	REF 718-342-10
50 cases	1000 brackets	REF 718-342-50


** The Dentaaurum version of this prescription is not claimed to be a duplication of any other, nor does Dentaaurum imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.

McLaughlin-Bennett-Trevisi** 22

Slot: 0.56 mm x 0.76 mm / 22 x 30 · Cuspid brackets with hook

Torque	- 7°	- 7°	0°	+ 10°	+ 17°	+ 17°	+ 10°	0°	- 7°	- 7°
Angulation	0°	0°	+ 8°	+ 8°	+ 4°	+ 4°	+ 8°	+ 8°	0°	0°
In/Out	1.2 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	1.2 mm
Width	3.1 mm	3.1 mm	3.1 mm	3.0 mm	3.4 mm	3.4 mm	3.0 mm	3.1 mm	3.1 mm	3.1 mm
REF	722-302-15	722-302-14	722-342-13	722-302-12	722-302-11	722-302-21	722-302-22	722-342-23	722-302-24	722-302-25
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Even easier bracket identification with FDI-marking on the base!



Torque	- 17°	- 12°	0°	- 6°	- 6°	- 6°	- 6°	0°	- 12°	- 17°
Angulation	+ 2°	+ 2°	+ 3°	0°	0°	0°	0°	+ 3°	+ 2°	+ 2°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	1.3 mm	1.3 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	3.1 mm	3.1 mm	3.1 mm	2.4 mm	2.4 mm	2.4 mm	2.4 mm	3.1 mm	3.1 mm	3.1 mm
REF	722-302-45	722-302-44	722-342-43	722-302-42	722-302-41	722-302-31	722-302-32	722-342-33	722-302-34	722-302-35
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Assortments:

1 case	20 brackets	REF 722-302-01
10 cases	200 brackets	REF 722-302-10
50 cases	1000 brackets	REF 722-302-50


** The Dentaureum version of this prescription is not claimed to be a duplication of any other, nor does Dentaureum imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.

McLaughlin-Bennett-Trevisi** 22

Slot: 0.56 mm x 0.76 mm / 22 x 30 · Cuspid, 1st bicuspid and 2nd bicuspid brackets with hook

Torque	- 7°	- 7°	0°	+ 10°	+ 17°	+ 17°	+ 10°	0°	- 7°	- 7°
Angulation	0°	0°	+ 8°	+ 8°	+ 4°	+ 4°	+ 8°	+ 8°	0°	0°
In/Out	1.2 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	1.2 mm
Width	3.1 mm	3.1 mm	3.1 mm	3.0 mm	3.4 mm	3.4 mm	3.0 mm	3.1 mm	3.1 mm	3.1 mm
REF	722-342-15	722-342-14	722-342-13	722-302-12	722-302-11	722-302-21	722-302-22	722-342-23	722-342-24	722-342-25
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Even easier bracket identification with FDI-marking on the base!



Torque	- 17°	- 12°	0°	- 6°	- 6°	- 6°	- 6°	0°	- 12°	- 17°
Angulation	+ 2°	+ 2°	+ 3°	0°	0°	0°	0°	+ 3°	+ 2°	+ 2°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	1.3 mm	1.3 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	3.1 mm	3.1 mm	3.1 mm	2.4 mm	2.4 mm	2.4 mm	2.4 mm	3.1 mm	3.1 mm	3.1 mm
REF	722-342-45	722-342-44	722-342-43	722-302-42	722-302-41	722-302-31	722-302-32	722-342-33	722-342-34	722-342-35
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Assortments:

1 case	20 brackets	REF 722-342-01
10 cases	200 brackets	REF 722-342-10
50 cases	1000 brackets	REF 722-342-50

** The Dentaureum version of this prescription is not claimed to be a duplication of any other, nor does Dentaureum imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.


equilibrium® 2 / Assortments

Ricketts® Universal 18

Slot: 0.46 mm x 0.76 mm / 18 x 30 · without hook

Torque	0°	0°	+ 7°	+ 14°	+ 22°	+ 22°	+ 14°	+ 7°	0°	0°
Angulation	0°	0°	+ 5°	+ 8°	0°	0°	+ 8°	+ 5°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm
Width	3.1 mm	3.1 mm	3.1 mm	2.9 mm	3.3 mm	3.3 mm	2.9 mm	3.1 mm	3.1 mm	3.1 mm
REF	718-303-15	718-303-14	718-303-13	718-303-12	718-303-11	718-303-21	718-303-22	718-303-23	718-303-24	718-303-25
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Even easier bracket identification with FDI-marking on the base!



Torque	0°	0°	+ 7°	0°	0°	0°	0°	+ 7°	0°	0°
Angulation	0°	0°	+ 5°	0°	0°	0°	0°	+ 5°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm
Width	3.1 mm	3.1 mm	2.7 mm	2.3 mm	2.3 mm	2.3 mm	2.3 mm	2.7 mm	3.1 mm	3.1 mm
REF	718-303-45	718-303-44	718-303-43	718-303-42	718-303-41	718-303-31	718-303-32	718-303-33	718-303-34	718-303-35
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Assortments:

1 case	20 brackets	REF 718-303-01
10 cases	200 brackets	REF 718-303-10
50 cases	1000 brackets	REF 718-303-50

Note:


Other brackets with Ricketts® Universal values see single packs.

Standard Edgewise 18

Slot: 0.46 mm x 0.76 mm / 18 x 30 · without hook

Torque	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
Angulation	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm
Width	3.1 mm	3.1 mm	3.1 mm	2.9 mm	2.9 mm	3.3 mm	3.3 mm	2.9 mm	3.1 mm	3.1 mm
REF	718-305-15	718-305-14	718-305-13	718-305-12	718-305-11	718-305-21	718-305-22	718-305-23	718-305-24	718-305-25
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Even easier bracket identification with FDI-marking on the base!



Torque	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
Angulation	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm
Width	3.1 mm	3.1 mm	3.1 mm	2.3 mm	2.3 mm	2.3 mm	2.3 mm	3.1 mm	3.1 mm	3.1 mm
REF	718-305-45	718-305-44	718-305-43	718-305-42	718-305-41	718-305-31	718-305-32	718-305-33	718-305-34	718-305-35
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Assortments:


1 case	20 brackets	REF 718-305-01
10 cases	200 brackets	REF 718-305-10
50 cases	1000 brackets	REF 718-305-50

Standard Edgewise 22

Slot: 0.56 mm x 0.76 mm / 22 x 30 - without hook

Torque	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
Angulation	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm
Width	3.1 mm	3.1 mm	3.1 mm	2.9 mm	3.3 mm	3.3 mm	2.9 mm	3.1 mm	3.1 mm	3.1 mm
REF	722-305-15	722-305-14	722-305-13	722-305-12	722-305-11	722-305-21	722-305-22	722-305-23	722-305-24	722-305-25
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Even easier bracket identification with FDI-marking on the base!



Torque	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
Angulation	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm	0.7 mm
Width	3.1 mm	3.1 mm	3.1 mm	2.3 mm	2.3 mm	2.3 mm	2.3 mm	3.1 mm	3.1 mm	3.1 mm
REF	722-305-45	722-305-44	722-305-43	722-305-42	722-305-41	722-305-31	722-305-32	722-305-33	722-305-34	722-305-35
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Assortments:

1 case	20 brackets	REF 722-305-01
10 cases	200 brackets	REF 722-305-10
50 cases	1000 brackets	REF 722-305-50



Metal brackets

equilibrium® 2 / Single packs – bonding technique



Roth 18/22

Tooth	Torque	Angulation	In/Out	Hook	Width	REF (18)	REF (22)	Quantity
13	- 2°	+ 13°	0.7 mm		3.1 mm		722-301-13	10 pieces
23	- 2°	+ 13°	0.7 mm		3.1 mm		722-301-23	10 pieces



Tooth	Torque	Angulation	In/Out	Hook	Width	REF (18)	REF (22)	Quantity
43	- 11°	+ 7°	0.7 mm		3.1 mm	718-301-43	722-301-43	10 pieces
33	- 11°	+ 7°	0.7 mm		3.1 mm	718-301-33	722-301-33	10 pieces

Even easier bracket identification with FDI-marking on the base!
 Slot 18: 0.46 mm x 0.76 mm / 18 x 30
 Slot 22: 0.56 mm x 0.76 mm / 22 x 30



Ricketts® Universal 18

Tooth	Torque	Angulation	In/Out	Hook	Width	REF (18)	REF (22)	Quantity
13	+ 7°	+ 5°	0.7 mm	•	3.1 mm	718-343-13		10 pieces
23	+ 7°	+ 5°	0.7 mm	•	3.1 mm	718-343-23		10 pieces



Tooth	Torque	Angulation	In/Out	Hook	Width	REF (18)	REF (22)	Quantity
43	+ 7°	+ 5°	0.7 mm	•	2.7 mm	718-343-43		10 pieces
33	+ 7°	+ 5°	0.7 mm	•	2.7 mm	718-343-33		10 pieces

Even easier bracket identification with FDI-marking on the base!
 Slot 18: 0.46 mm x 0.76 mm / 18 x 30

equilibrium® 2 / Single packs – banding technique

Roth 18/22

Tooth	Torque	Angulation	In/Out	Hook	Width	REF (18)	REF (22)	Quantity
13	- 2°	+ 13°	0.7 mm		3.1 mm	790-259-02	790-263-02	10 pieces
13	- 2°	+ 13°	0.7 mm	•	3.1 mm	790-204-02	790-216-02	10 pieces
23	- 2°	+ 13°	0.7 mm		3.1 mm	790-260-02	790-264-02	10 pieces
23	- 2°	+ 13°	0.7 mm	•	3.1 mm	790-205-02	790-217-02	10 pieces



Tooth	Torque	Angulation	In/Out	Hook	Width	REF (18)	REF (22)	Quantity
43	- 11°	+ 7°	0.7 mm		3.1 mm	790-261-02	790-265-02	10 pieces
43	- 11°	+ 7°	0.7 mm	•	3.1 mm	790-208-02	790-220-02	10 pieces
33	- 11°	+ 7°	0.7 mm		3.1 mm	790-262-02	790-266-02	10 pieces
33	- 11°	+ 7°	0.7 mm	•	3.1 mm	790-209-02	790-221-02	10 pieces
44 / 34	- 17°	0°	0.7 mm		3.1 mm	790-210-02	790-222-02	10 pieces
44	- 17°	0°	0.7 mm	•	3.1 mm	790-251-02	790-255-02	10 pieces
34	- 17°	0°	0.7 mm	•	3.1 mm	790-252-02	790-256-02	10 pieces
45 / 35	- 22°	0°	0.7 mm		3.1 mm	790-211-02	790-223-02	10 pieces
45	- 22°	0°	0.7 mm	•	3.1 mm	790-253-02	790-257-02	10 pieces
35	- 22°	0°	0.7 mm	•	3.1 mm	790-254-02	790-258-02	10 pieces



Slot 18: 0.46 mm x 0.76 mm / 18 x 30
Slot 22: 0.56 mm x 0.76 mm / 22 x 30

Roth, McLaughlin-Bennett-Trevisi** 18/22

Tooth	Torque	Angulation	In/Out	Hook	Width	REF (18)	REF (22)	Quantity
15-14 / 24-25	- 7°	0°	0.7 mm		3.1 mm	790-206-02	790-218-02	10 pieces
15-14	- 7°	0°	0.7 mm	•	3.1 mm	790-249-02	790-243-02	10 pieces
24-25	- 7°	0°	0.7 mm	•	3.1 mm	790-250-02	790-244-02	10 pieces



Slot 18: 0.46 mm x 0.76 mm / 18 x 30
Slot 22: 0.56 mm x 0.76 mm / 22 x 30

** The Dentaaurum version of this prescription is not claimed to be a duplication of any other, nor does Dentaaurum imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.

McLaughlin-Bennett-Trevisi** 18/22

Tooth	Torque	Angulation	In/Out	Hook	Width	REF (18)	REF (22)	Quantity
13	0°	+ 8°	0.7 mm	•	3.1 mm	790-274-02	790-228-02	10 pieces
13	- 7°	+ 8°	0.7 mm		3.1 mm		790-230-02	10 pieces
23	0°	+ 8°	0.7 mm	•	3.1 mm	790-275-02	790-229-02	10 pieces
23	- 7°	+ 8°	0.7 mm		3.1 mm		790-231-02	10 pieces

Tooth	Torque	Angulation	In/Out	Hook	Width	REF (18)	REF (22)	Quantity
43	0°	+ 3°	0.7 mm	•	3.1 mm	790-278-02	790-234-02	10 pieces
43	- 6°	+ 3°	0.7 mm		3.1 mm		790-236-02	10 pieces
33	0°	+ 3°	0.7 mm	•	3.1 mm	790-279-02	790-235-02	10 pieces
33	- 6°	+ 3°	0.7 mm		3.1 mm		790-237-02	10 pieces



Slot 18: 0.46 mm x 0.76 mm / 18 x 30
Slot 22: 0.56 mm x 0.76 mm / 22 x 30

** The Dentaaurum version of this prescription is not claimed to be a duplication of any other, nor does Dentaaurum imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.

equilibrium® 2 / Single packs – banding technique



Standard Edgewise 18/22

Tooth	Torque	Angulation	In/Out	Hook	Width	REF (18)	REF (22)	Quantity
42-41 / 31-32	0°	0°	0.7 mm		2.4 mm	790-603-02	790-607-02	10 pieces



Tooth	Torque	Angulation	In/Out	Hook	Width	REF (18)	REF (22)	Quantity
15-13 / 23-25 / 45-43 / 33-35	0°	0°	0.7 mm		3.1 mm	790-602-02	790-606-02	10 pieces

Slot 18: 0.46 mm x 0.76 mm / 18 x 30
Slot 22: 0.56 mm x 0.76 mm / 22 x 30



Individual systems 18/22

Tooth	Torque	Angulation	In/Out	Hook	Width	REF (18)	REF (22)	Quantity
13 / 24-25	+ 7°	0°	0.7 mm	•	3.1 mm	790-290-02		10 pieces
15-14 / 23	+ 7°	0°	0.7 mm	•	3.1 mm	790-291-02		10 pieces



Tooth	Torque	Angulation	In/Out	Hook	Width	REF (18)	REF (22)	Quantity
45-43	- 11°	0°	0.7 mm	•	3.1 mm		790-294-02	10 pieces
33-35	- 11°	0°	0.7 mm	•	3.1 mm		790-295-02	10 pieces


Slot 18: 0.46 mm x 0.76 mm / 18 x 30
Slot 22: 0.56 mm x 0.76 mm / 22 x 30

Roth 18

Slot: 0.46 mm x 0.76 mm / 18 x 30 · Cuspid brackets with hook

Torque	- 7°	- 7°	- 2°	+ 8°	+ 12°	+ 12°	+ 8°	- 2°	- 7°	- 7°
Angulation	0°	0°	+ 13°	+ 9°	+ 5°	+ 5°	+ 9°	+ 13°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	3.1 mm	3.1 mm	3.1 mm	3.0 mm	3.4 mm	3.4 mm	3.0 mm	3.1 mm	3.1 mm	3.1 mm
REF	718-501-15	718-501-14	718-541-13	718-501-12	718-501-11	718-501-21	718-501-22	718-541-23	718-501-24	718-501-25
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Even easier bracket identification with FDI-marking on the base!



Torque	- 22°	- 17°	- 11°	- 1°	- 1°	- 1°	- 1°	- 11°	- 17°	- 22°
Angulation	0°	0°	+ 7°	0°	0°	0°	0°	+ 7°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	1.3 mm	1.3 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	3.1 mm	3.1 mm	3.1 mm	2.4 mm	2.4 mm	2.4 mm	2.4 mm	3.1 mm	3.1 mm	3.1 mm
REF	718-501-45	718-501-44	718-541-43	718-501-42	718-501-41	718-501-31	718-501-32	718-541-33	718-501-34	718-501-35
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Assortments:

1 case	20 brackets	REF 718-501-01
10 cases	200 brackets	REF 718-501-10
50 cases	1000 brackets	REF 718-501-50

Note:


1 case assortment REF 718-501-02 is equivalent to REF 718-501-01, with additional 4 buccal tubes with vertical hook (1 x 723-028-00 / 1 x 723-032-00 / 2 x 723-026-00).

Roth 22

Slot: 0.56 mm x 0.76 mm / 22 x 30 · Cuspid brackets with hook

Torque	- 7°	- 7°	- 2°	+ 8°	+ 12°	+ 12°	+ 8°	- 2°	- 7°	- 7°
Angulation	0°	0°	+ 13°	+ 9°	+ 5°	+ 5°	+ 9°	+ 13°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	3.1 mm	3.1 mm	3.1 mm	3.0 mm	3.4 mm	3.4 mm	3.0 mm	3.1 mm	3.1 mm	3.1 mm
REF	722-501-15	722-501-14	722-541-13	722-501-12	722-501-11	722-501-21	722-501-22	722-541-23	722-501-24	722-501-25
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Even easier bracket identification with FDI-marking on the base!



Torque	- 22°	- 17°	- 11°	- 1°	- 1°	- 1°	- 1°	- 11°	- 17°	- 22°
Angulation	0°	0°	+ 7°	0°	0°	0°	0°	+ 7°	0°	0°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	1.3 mm	1.3 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	3.1 mm	3.1 mm	3.1 mm	2.4 mm	2.4 mm	2.4 mm	2.4 mm	3.1 mm	3.1 mm	3.1 mm
REF	722-501-45	722-501-44	722-541-43	722-501-42	722-501-41	722-501-31	722-501-32	722-541-33	722-501-34	722-501-35
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Assortments:

1 case	20 brackets	REF 722-501-01
10 cases	200 brackets	REF 722-501-10
50 cases	1000 brackets	REF 722-501-50

Note:

1 case assortment REF 722-501-02 is equivalent to REF 722-501-01, with additional 4 buccal tubes with vertical hook (1 x 723-029-00 / 1 x 723-033-00 / 2 x 723-027-00).


equilibrium® ti / Assortments

McLaughlin-Bennett-Trevisi** 18

Slot: 0.46 mm x 0.76 mm / 18 x 30 · Cuspid brackets with hook

Torque	- 7°	- 7°	0°	+ 10°	+ 17°	+ 17°	+ 10°	0°	- 7°	- 7°
Angulation	0°	0°	+ 8°	+ 8°	+ 4°	+ 4°	+ 8°	+ 8°	0°	0°
In/Out	1.2 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	1.2 mm
Width	3.1 mm	3.1 mm	3.1 mm	3.0 mm	3.4 mm	3.4 mm	3.0 mm	3.1 mm	3.1 mm	3.1 mm
REF	718-502-15	718-502-14	718-542-13	718-502-12	718-502-11	718-502-21	718-502-22	718-542-23	718-502-24	718-502-25
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Even easier bracket identification with FDI-marking on the base!



Torque	- 17°	- 12°	0°	- 6°	- 6°	- 6°	- 6°	0°	- 12°	- 17°
Angulation	+ 2°	+ 2°	+ 3°	0°	0°	0°	0°	+ 3°	+ 2°	+ 2°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	1.3 mm	1.3 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	3.1 mm	3.1 mm	3.1 mm	2.4 mm	2.4 mm	2.4 mm	2.4 mm	3.1 mm	3.1 mm	3.1 mm
REF	718-502-45	718-502-44	718-542-43	718-502-42	718-502-41	718-502-31	718-502-32	718-542-33	718-502-34	718-502-35
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Assortments:

1 case	20 brackets	REF 718-502-01
10 cases	200 brackets	REF 718-502-10
50 cases	1000 brackets	REF 718-502-50

Note:

1 case assortment REF 718-502-02 is equivalent to REF 718-502-01, with additional 4 buccal tubes with vertical hook (1 x 723-028-00 / 1 x 723-032-00 / 2 x 723-030-00).


** The Dentaaurum version of this prescription is not claimed to be a duplication of any other, nor does Dentaaurum imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.

McLaughlin-Bennett-Trevisi** 22

Slot: 0.56 mm x 0.76 mm / 22 x 30 · Cuspid brackets with hook

Torque	- 7°	- 7°	0°	+ 10°	+ 17°	+ 17°	+ 10°	0°	- 7°	- 7°
Angulation	0°	0°	+ 8°	+ 8°	+ 4°	+ 4°	+ 8°	+ 8°	0°	0°
In/Out	1.2 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	1.3 mm	0.7 mm	0.7 mm	1.2 mm
Width	3.1 mm	3.1 mm	3.1 mm	3.0 mm	3.4 mm	3.4 mm	3.0 mm	3.1 mm	3.1 mm	3.1 mm
REF	722-502-15	722-502-14	722-542-13	722-502-12	722-502-11	722-502-21	722-502-22	722-542-23	722-502-24	722-502-25
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Even easier bracket identification with FDI-marking on the base!



Torque	- 17°	- 12°	0°	- 6°	- 6°	- 6°	- 6°	0°	- 12°	- 17°
Angulation	+ 2°	+ 2°	+ 3°	0°	0°	0°	0°	+ 3°	+ 2°	+ 2°
In/Out	0.7 mm	0.7 mm	0.7 mm	1.3 mm	1.3 mm	1.3 mm	1.3 mm	0.7 mm	0.7 mm	0.7 mm
Width	3.1 mm	3.1 mm	3.1 mm	2.4 mm	2.4 mm	2.4 mm	2.4 mm	3.1 mm	3.1 mm	3.1 mm
REF	722-502-45	722-502-44	722-542-43	722-502-42	722-502-41	722-502-31	722-502-32	722-542-33	722-502-34	722-502-35
Quantity	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces	10 pieces

Assortments:

1 case	20 brackets	REF 722-502-01
10 cases	200 brackets	REF 722-502-10
50 cases	1000 brackets	REF 722-502-50

Note:

1 case assortment REF 722-502-02 is equivalent to REF 722-502-01, with additional 4 buccal tubes with vertical hook (1 x 723-029-00 / 1 x 723-033-00 / 2 x 723-031-00).

** The Dentaaurum version of this prescription is not claimed to be a duplication of any other, nor does Dentaaurum imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.



Handwriting practice lines consisting of 20 horizontal grey lines.

discovery® pearl ceramic brackets / Assortments

Roth 18/22	Page	104
McLaughlin-Bennett-Trevisi* 18/22	Page	105

discovery® pearl ceramic brackets / Single packs

Roth 18/22	Page	106
McLaughlin-Bennett-Trevisi* 18/22	Page	107

discovery® pearl & smart combined assortments

Roth 22	Page	108
McLaughlin-Bennett-Trevisi* 22	Page	108

discovery® pearl accessories

Page	109
------	-----

Nickel-free products

Page	110
------	-----

Bracket accessories


Page	114
------	-----

* The Dentaurem version of this prescription is not claimed to be a duplication of any other, nor does Dentaurem imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.



Aesthetic brackets.

For an even more beautiful smile.



DENTAURUM
**ONLINE
SHOP**
shop.dentaurum.com





Rediscover aesthetics.

The highest quality standards.

Brackets in the wide Dentaureum range, whether metal or ceramic, are manufactured in state-of-the-art design and production procedures on site at Ispringen, Germany. All these processes have been optimized and perfectly coordinated over the years. A very high level of precision is essential to ensure safe, reliable use in orthodontic treatment. With innovative production technologies, such as the CIM process (Ceramic Injection Molding), the most diverse ceramic products can be manufactured using injection molding. Dentaureum seeks to continually develop the aesthetic materials to enable us to react optimally to the increased demand in the ceramic bracket sector.



discovery® pearl

Fascinatingly natural.



FDI System

	18	17	16	15	14	13	12	11		21	22	23	24	25	26	27	28	
	8	7	6	5	4	3	2	1		1	2	3	4	5	6	7	8	
R	8	7	6	5	4	3	2	1		1	2	3	4	5	6	7	8	L
	48	47	46	45	44	43	42	41		31	32	33	34	35	36	37	38	



fascinatingly natural.



Maximum quality and performance naturally without compromise.

discovery® pearl transcends all the previous boundaries of aesthetics, quality and performance. When it is bonded to the tooth, discovery® pearl adapts perfectly to the natural environment of the tooth due to its shade and high degree of translucency.

The bracket is made completely from ceramic without any bothersome metal inserts in the slot. It was designed in one piece and is manufactured in state-of-the-art production facilities using the CIM process (Ceramic Injection Molding).

The base of the bracket is optimally adapted to the contour of the tooth using CAD (Computer Aided Design), which allows it to be easily positioned on the tooth. A clear system consisting of colored positioning guides and FDI notation on the base facilitates tooth-specific assignment of each bracket. The positioning guides also help the operator align and place the bracket very easily.

The innovative laser-structured base ensures a strong bond, enables safe and reliable debonding and guarantees that discovery® pearl will produce optimal treatment outcomes and meet the ever increasing demands for aesthetic tooth correction.



Reliable product characteristics:

Proven bracket geometry of the discovery® family

ensures easy combination between discovery® metal and ceramic brackets.

Perfect shade adaptation

to the natural environment of the tooth – as if the bracket were not there.

Colored positioning guides

facilitate positioning on the tooth.

Optimal size

compact in design for optimal treatment & patient comfort.

Anatomical base

for a perfect fit on the tooth.

Laser-structured base with FDI notation

for optimal bonding and debonding characteristics.

Now also available as 1 case assortment without positioning guide.



Aesthetic brackets

discovery® pearl ceramic brackets / Assortments



Roth 18/22

Cuspid brackets with hook

Tooth	with positioning guide		without positioning guide		Quantity
	REF (18)	REF (22)	REF (18)	REF (22)	
13-11 / 21-23	718-701-03	722-701-03	718-701-73	722-701-73	6 pieces (1 bracket each)
15-11 / 21-25	718-701-05	722-701-05	718-701-75	722-701-75	10 pieces (1 bracket each)

Tooth	with positioning guide		without positioning guide		Quantity
	REF (18)	REF (22)	REF (18)	REF (22)	
13-11 / 21-23 / 43-41 / 31-33	718-701-63	722-701-63	718-701-83	722-701-83	12 pieces (1 bracket each)
15-11 / 21-25 / 45-41 / 31-35	718-701-65	722-701-65	718-701-85	722-701-85	20 pieces (1 bracket each)

Slot 18: 0.46 mm x 0.76 mm / 18 x 30

Slot 22: 0.56 mm x 0.76 mm / 22 x 30



Roth 18/22

Cuspid, 1st bicuspid and 2nd bicuspid brackets with hook

Tooth	with positioning guide		without positioning guide		Quantity
	REF (18)	REF (22)	REF (18)	REF (22)	
15-11 / 21-25	718-741-60	722-741-60	718-741-80	722-741-80	10 pieces (1 bracket each)

Tooth	with positioning guide		without positioning guide		Quantity
	REF (18)	REF (22)	REF (18)	REF (22)	
15-11 / 21-25 / 45-41 / 31-35	718-741-65	722-741-65	718-741-85	722-741-85	20 pieces (1 bracket each)

Slot 18: 0.46 mm x 0.76 mm / 18 x 30

Slot 22: 0.56 mm x 0.76 mm / 22 x 30

discovery® pearl ceramic brackets / Assortments

McLaughlin-Bennett-Trevisi** 18/22

Cuspid brackets with hook

Tooth	with positioning guide		without positioning guide		Quantity
	REF (18)	REF (22)	REF (18)	REF (22)	
13-11 / 21-23	718-702-03	722-702-03	718-702-73	722-702-73	6 pieces (1 bracket each)
15-11 / 21-25	718-702-05	722-702-05	718-702-75	722-702-75	10 pieces (1 bracket each)

Tooth	with positioning guide		without positioning guide		Quantity
	REF (18)	REF (22)	REF (18)	REF (22)	
13-11 / 21-23 / 43-41 / 31-33	718-702-63	722-702-63	718-702-83	722-702-83	12 pieces (1 bracket each)
15-11 / 21-25 / 45-41 / 31-35	718-702-65	722-702-65	718-702-85	722-702-85	20 pieces (1 bracket each)

Slot 18: 0.46 mm x 0.76 mm / 18 x 30

Slot 22: 0.56 mm x 0.76 mm / 22 x 30

** The Dentaaurum version of this prescription is not claimed to be a duplication of any other, nor does Dentaaurum imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.



McLaughlin-Bennett-Trevisi** 18/22

Cuspid, 1st bicuspid and 2nd bicuspid brackets with hook

Tooth	with positioning guide		without positioning guide		Quantity
	REF (18)	REF (22)	REF (18)	REF (22)	
15-11 / 21-25	718-742-60	722-742-60	718-742-80	722-742-80	10 pieces (1 bracket each)

Tooth	with positioning guide		without positioning guide		Quantity
	REF (18)	REF (22)	REF (18)	REF (22)	
15-11 / 21-25 / 45-41 / 31-35	718-742-65	722-742-65	718-742-85	722-742-85	20 pieces (1 bracket each)

Slot 18: 0.46 mm x 0.76 mm / 18 x 30

Slot 22: 0.56 mm x 0.76 mm / 22 x 30

** The Dentaaurum version of this prescription is not claimed to be a duplication of any other, nor does Dentaaurum imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.



discovery® pearl ceramic brackets / Single packs



Roth 18/22

with positioning guide

Tooth	Torque	Angulation	In/Out	Hook	Width	REF (18)	REF (22)	Quantity
11	+ 12°	+ 5°	0.7 mm		3.1 mm	718-701-11	722-701-11	1 piece
21	+ 12°	+ 5°	0.7 mm		3.1 mm	718-701-21	722-701-21	1 piece
12	+ 8°	+ 9°	1.3 mm		2.7 mm	718-701-12	722-701-12	1 piece
22	+ 8°	+ 9°	1.3 mm		2.7 mm	718-701-22	722-701-22	1 piece
13	- 2°	+ 13°	0.7 mm	•	3.1 mm	718-741-13	722-741-13	1 piece
23	- 2°	+ 13°	0.7 mm	•	3.1 mm	718-741-23	722-741-23	1 piece
14	- 7°	0°	0.7 mm		3.1 mm	718-701-14	722-701-14	1 piece
24	- 7°	0°	0.7 mm		3.1 mm	718-701-24	722-701-24	1 piece
14	- 7°	0°	0.7 mm	•	3.1 mm	718-741-14	722-741-14	1 piece
24	- 7°	0°	0.7 mm	•	3.1 mm	718-741-24	722-741-24	1 piece
15	- 7°	0°	0.7 mm		3.1 mm	718-701-15	722-701-15	1 piece
25	- 7°	0°	0.7 mm		3.1 mm	718-701-25	722-701-25	1 piece
15	- 7°	0°	0.7 mm	•	3.1 mm	718-741-15	722-741-15	1 piece
25	- 7°	0°	0.7 mm	•	3.1 mm	718-741-25	722-741-25	1 piece



Tooth	Torque	Angulation	In/Out	Hook	Width	REF (18)	REF (22)	Quantity
41	- 1°	0°	1.3 mm		2.7 mm	718-701-41	722-701-41	1 piece
31	- 1°	0°	1.3 mm		2.7 mm	718-701-31	722-701-31	1 piece
42	- 1°	0°	1.3 mm		2.7 mm	718-701-42	722-701-42	1 piece
32	- 1°	0°	1.3 mm		2.7 mm	718-701-32	722-701-32	1 piece
43	- 11°	+ 7°	0.7 mm	•	3.1 mm	718-741-43	722-741-43	1 piece
33	- 11°	+ 7°	0.7 mm	•	3.1 mm	718-741-33	722-741-33	1 piece
44	- 17°	0°	0.7 mm		3.1 mm	718-701-44	722-701-44	1 piece
34	- 17°	0°	0.7 mm		3.1 mm	718-701-34	722-701-34	1 piece
44	- 17°	0°	0.7 mm	•	3.1 mm	718-741-44	722-741-44	1 piece
34	- 17°	0°	0.7 mm	•	3.1 mm	718-741-34	722-741-34	1 piece
45	- 22°	0°	0.7 mm		3.1 mm	718-701-45	722-701-45	1 piece
35	- 22°	0°	0.7 mm		3.1 mm	718-701-35	722-701-35	1 piece
45	- 22°	0°	0.7 mm	•	3.1 mm	718-741-45	722-741-45	1 piece
35	- 22°	0°	0.7 mm	•	3.1 mm	718-741-35	722-741-35	1 piece

Slot 18: 0.46 mm x 0.76 mm / 18 x 30
Slot 22: 0.56 mm x 0.76 mm / 22 x 30

discovery® pearl ceramic brackets / Single packs

McLaughlin-Bennett-Trevisi** 18/22

with positioning guide



Tooth	Torque	Angulation	In/Out	Hook	Width	REF (18)	REF (22)	Quantity
11	+ 17°	+ 4°	0.7 mm		3.1 mm	718-702-11	722-702-11	1 piece
21	+ 17°	+ 4°	0.7 mm		3.1 mm	718-702-21	722-702-21	1 piece
12	+ 10°	+ 8°	1.3 mm		2.7 mm	718-702-12	722-702-12	1 piece
22	+ 10°	+ 8°	1.3 mm		2.7 mm	718-702-22	722-702-22	1 piece
13	0°	+ 8°	0.7 mm	•	3.1 mm	718-742-13	722-742-13	1 piece
23	0°	+ 8°	0.7 mm	•	3.1 mm	718-742-23	722-742-23	1 piece
14	- 7°	0°	0.7 mm		3.1 mm	718-702-14	722-702-14	1 piece
24	- 7°	0°	0.7 mm		3.1 mm	718-702-24	722-702-24	1 piece
14	- 7°	0°	0.7 mm	•	3.1 mm	718-742-14	722-742-14	1 piece
24	- 7°	0°	0.7 mm	•	3.1 mm	718-742-24	722-742-24	1 piece
15	- 7°	0°	1.2 mm		3.1 mm	718-702-15	722-702-15	1 piece
25	- 7°	0°	1.2 mm		3.1 mm	718-702-25	722-702-25	1 piece
15	- 7°	0°	1.2 mm	•	3.1 mm	718-742-15	722-742-15	1 piece
25	- 7°	0°	1.2 mm	•	3.1 mm	718-742-25	722-742-25	1 piece

Tooth	Torque	Angulation	In/Out	Hook	Width	REF (18)	REF (22)	Quantity
41	- 6°	0°	1.3 mm		2.7 mm	718-702-41	722-702-41	1 piece
31	- 6°	0°	1.3 mm		2.7 mm	718-702-31	722-702-31	1 piece
42	- 6°	0°	1.3 mm		2.7 mm	718-702-42	722-702-42	1 piece
32	- 6°	0°	1.3 mm		2.7 mm	718-702-32	722-702-32	1 piece
43	0°	+ 3°	0.7 mm	•	3.1 mm	718-742-43	722-742-43	1 piece
33	0°	+ 3°	0.7 mm	•	3.1 mm	718-742-33	722-742-33	1 piece
44	- 12°	+ 2°	0.7 mm		3.1 mm	718-702-44	722-702-44	1 piece
34	- 12°	+ 2°	0.7 mm		3.1 mm	718-702-34	722-702-34	1 piece
44	- 12°	+ 2°	0.7 mm	•	3.1 mm	718-742-44	722-742-44	1 piece
34	- 12°	+ 2°	0.7 mm	•	3.1 mm	718-742-34	722-742-34	1 piece
45	- 17°	+ 2°	0.7 mm		3.1 mm	718-702-45	722-702-45	1 piece
35	- 17°	+ 2°	0.7 mm		3.1 mm	718-702-35	722-702-35	1 piece
45	- 17°	+ 2°	0.7 mm	•	3.1 mm	718-742-45	722-742-45	1 piece
35	- 17°	+ 2°	0.7 mm	•	3.1 mm	718-742-35	722-742-35	1 piece



Slot 18: 0.46 mm x 0.76 mm / 18 x 30

Slot 22: 0.56 mm x 0.76 mm / 22 x 30

** The Dentaaurum version of this prescription is not claimed to be a duplication of any other, nor does Dentaaurum imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.

Aesthetic brackets

discovery® pearl & smart combined assortments



Roth 22

Cuspid brackets with hook · with positioning guide

Brackets	Tooth	REF (22)	Quantity
discovery® pearl	15-11 / 21-25	891-404-00	20 pieces (1 bracket each)
discovery® smart	45-41 / 31-35		

Slot 22: 0.56 mm x 0.76 mm / 22 x 30



McLaughlin-Bennett-Trevisi** 22

Cuspid brackets with hook · with positioning guide

Brackets	Tooth	REF (22)	Quantity
discovery® pearl	15-11 / 21-25	891-400-00	20 pieces (1 bracket each)
discovery® smart	45-41 / 31-35		

Slot 22: 0.56 mm x 0.76 mm / 22 x 30

** The Dentaaurum version of this prescription is not claimed to be a duplication of any other, nor does Dentaaurum imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.

discovery® pearl ceramic brackets / Accessories

discovery® pearl debonding instrument sterilizable

For easy and safe removal of **discovery® pearl** ceramic brackets.

REF 019-001-00

1 piece




Dentalastics® pearl

Plastic ligatures

Especially delicate plastic ligatures, instead of wire ligatures, for ligating archwires in brackets.

Contents: 50 sticks x 20 rings = 1000 pieces

Color	Inner ø	REF	Quantity
pearl	1.3 mm	774-563-01	1 pack



SUPER Elasto-Force plastic chains

Consistent traction force over a long period of time. High elasticity with excellent elastic recovery.

Color	Inner ø	Length	REF	Quantity
pearl	1.3 mm	4.6 m	774-236-01	1 roll



Orthodontic demonstration model discovery® pearl

McLaughlin-Bennett-Trevisi* 22

Fitted with:

- discovery® pearl ceramic brackets
- translucent arch pearl
- Dentalastics® pearl plastic ligatures

Maxilla made of plastic.

REF 892-000-18

1 piece

* The Dentaureum version of this prescription is not claimed to be a duplication of any other, nor does Dentaureum imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.



Orthodontic demonstration model discovery® pearl

Roth 22

Fitted with:

- discovery® pearl ceramic brackets, maxilla / mandible 5-5
- Ortho-Cast M-Series buccal tubes, maxilla / mandible 6
- translucent arch pearl, maxilla
- Tensic® arch, mandible
- Dentalastics® pearl plastic ligatures

Maxilla and mandible made of plastic, hinged.

REF 892-000-21

1 piece



Nickel-free products

You can find further information at
www.dentaurum.de/nickel-free.



Nickel-free products.

The nickel-free therapy.







DENTAURUM
**ONLINE
SHOP**
shop.dentaurum.com



Product overview.

Dentaurum offers a wide selection of nickel-free products for holistic therapy.

Products	Chapter	
topic bracket made of remanium® star*	Metal brackets	
equilibrium® ti bracket made of pure titanium	Metal brackets	
discovery® pearl bracket made of aluminum oxide	Aesthetic brackets	

Products	Chapter	
Noninium® facebow**	Intra- and extra-oral	
rematitan® buccal tubes made of pure titanium	Buccal tubes/ Bonding technique	
rematitan® SPECIAL wires and rematitan® SPECIAL arches	Wires Arches Accessories	
rematitan® & Titan hyrax® · Expansion screws made of pure titanium	Expansion screws	
prime4me® RETAIN3R · Milled from titanium grade 5	Digital orthodontics	

* see Materials list Orthodontics

**Nickel ≤ 0.2% is due to the process of manufacturing

Bracket accessories



Bite-Bumper™

Bonding technique
Material: Plastic
Palatal bite raiser for raising a deep overbite
in Class II treatment.

REF 705-000-00

10 pieces



Titanium button with chain acc. to Watted

Bonding technique
With laser-structured titanium base and titanium chain.
For aligning retracted anteriors and cuspids.
Length: 5 cm

REF 750-003-51

1 piece



Lingual buttons

Bonding technique
With laser-structured base, flat base, square. For incisors.
Dimensions (H x W): 3.5 mm x 3.5 mm

REF 750-002-51

10 pieces



Lingual buttons

Bonding technique
With laser-structured base, contoured base, round.
For cuspids and bicuspid.
ø 3.0 mm

REF 750-001-51

10 pieces



Lingual cleats

Bonding technique
With mesh base for bonding, flat.
For additional elastics or traction springs.
Dimensions (H x W): 3.8 mm x 3.8 mm

REF 750-602-50

10 pieces



Mesh base

Bonding technique
Anatomically contoured.

Indication	Base	Dimensions (H x W)	REF	Quantity
for incisors	flat	3.8 mm x 3.8 mm	799-202-00	10 pieces
for buccal tubes and 2nd molars	small	3.0 mm x 6.5 mm	799-206-00	10 pieces

Protection wax

Protects mucosa and lips in first days after bracket bonding until patient has become used to the appliance. Also for relief from irritation during treatment.

Form of delivery:

50 small boxes with 5 bars of protection wax each.

Color: white

REF 124-305-01

50 packs

**Bracket tray**

With self sticking foil to presort brackets.

REF 891-002-00

50 pieces

**Bracket container**

For storage of previously sorted brackets.

REF 891-001-00

1 piece

**Storage box**

For sorting brackets individually.

REF 766-003-10

1 piece



Bonding technique	Page	120
Bonding technique / Titanium tubes	Page	125
Banding technique	Page	126
Prewelded parts	Page	134



Note

Tube dimensions	Main tube
	Slot: 0.46 mm x 0.64 mm / 18 x 25 or Slot: 0.56 mm x 0.70 mm / 22 x 28
	Auxiliary tube
	Slot: 0.46 mm x 0.64 mm / 18 x 25 e.g. for segmental arches in the segment arch technique (max. recommended wire thickness: 0.43 mm x 0.64 mm / 17 x 25)
	Lip bumper/Facebow tubes
	ø 1.15 mm / 45
	In/Out
	0.3 mm
Hook	Hooks are positioned mesially unless stated otherwise. Easy identification:
	+ = maxilla - = mandible

Images of the buccal tubes are examples only.


Information can be found easily with Dentaaurum's color coding in the catalog

Maxillary
 Mandibular
 Maxillary + Mandibular

Buccal tubes.

Perfect in function and design.





Perfect in function and design.

Just like the Dentaurem brackets, the buccal tubes are manufactured using state-of-the-art design and production techniques. Depending on the material, the tubes are available in different variations ranging from single, non-convertible tubes to more complex triple combinations. This allows more flexibility when planning therapy. Many buccal tubes are available for the adhesive technique as well as the banding technique. The buccal tubes can be welded to the bands of your choice in our pre-welding service.

The latest buccal tube generation, the Ortho-Cast M-Series, is particularly easy to handle. Thanks to contouring of the base and a groove for the positioning tweezers, the tubes are extremely easy to place. The funnel-shaped opening ensures easy archwire insertion. The M-Series is also available as Ortho-Cast M-Series mini for the second molar.



Ortho-Cast M-Series mini

The delicate buccal tube.



Ortho-Cast M-Series

The latest generation.



Ortho-Cast

The MIM buccal tube.



Ortho

The classic buccal tube.



rematitan®

The only buccal tube made of pure titanium.



Buccal tubes

Bonding technique

single round



Ortho

non-convertible · with lip bumper/facebow tube · Quantity: 10 pieces

Tooth	Torque	Distal Offset	Hook	Width	REF	System
17-16 / 26-27 / 47-46 / 36-37		0°		6.0 mm	725-203-51	

Slot: 0.9 mm / 36

single rectangular



Ortho-Cast M-Series mini 18/22

non-convertible · Quantity: 10 pieces

Tooth	Torque	Distal Offset	Hook	Width	REF (18)	REF (22)	System
17	- 14°	0°	•	2.5 mm	818-411-17	822-411-17	Roth
27	- 14°	0°	•	2.5 mm	818-411-27	822-411-27	Roth
17	- 14°	0°	•	2.5 mm	818-412-17	822-412-17	McLaughlin-Bennett-Trevisi**
27	- 14°	0°	•	2.5 mm	818-412-27	822-412-27	McLaughlin-Bennett-Trevisi**
17	0°	0°	•	2.5 mm	818-413-17		Ricketts® Universal
27	0°	0°	•	2.5 mm	818-413-27		Ricketts® Universal
17	0°	0°	•	2.5 mm	818-415-17	822-415-17	Standard Edgewise
27	0°	0°	•	2.5 mm	818-415-27	822-415-27	Standard Edgewise



Tooth	Torque	Distal Offset	Hook	Width	REF (18)	REF (22)	System
47	- 30°	0°	•	2.5 mm	818-411-47	822-411-47	Roth
37	- 30°	0°	•	2.5 mm	818-411-37	822-411-37	Roth
47	- 10°	0°	•	2.5 mm	818-412-47	822-412-47	McLaughlin-Bennett-Trevisi**
37	- 10°	0°	•	2.5 mm	818-412-37	822-412-37	McLaughlin-Bennett-Trevisi**
46	- 20°	0°	•	2.5 mm	818-412-46	822-412-46	McLaughlin-Bennett-Trevisi**
36	- 20°	0°	•	2.5 mm	818-412-36	822-412-36	McLaughlin-Bennett-Trevisi**
47	0°	0°	•	2.5 mm	818-413-47		Ricketts® Universal
37	0°	0°	•	2.5 mm	818-413-37		Ricketts® Universal
47	0°	0°	•	2.5 mm	818-415-47	822-415-47	Standard Edgewise
37	0°	0°	•	2.5 mm	818-415-37	822-415-37	Standard Edgewise

Slot 18: 0.46 mm x 0.64 mm / 18 x 25

Slot 22: 0.56 mm x 0.70 mm / 22 x 28

** The Dentaaurum version of this prescription is not claimed to be a duplication of any other, nor does Dentaaurum imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.

single rectangular

Ortho-Cast M-Series 18/22

non-convertible · Quantity: 10 pieces

Tooth	Torque	Distal Offset	Hook	Width	REF (18)	REF (22)	System
16	- 14°	14°	•	3.3 mm	818-211-16	822-211-16	Roth
26	- 14°	14°	•	3.3 mm	818-211-26	822-211-26	Roth
16	- 14°	10°	•	3.3 mm	818-212-16	822-212-16	McLaughlin-Bennett-Trevisi**
26	- 14°	10°	•	3.3 mm	818-212-26	822-212-26	McLaughlin-Bennett-Trevisi**
16	0°	0°	•	3.3 mm	818-215-16	822-215-16	Standard Edgewise
26	0°	0°	•	3.3 mm	818-215-26	822-215-26	Standard Edgewise



Tooth	Torque	Distal Offset	Hook	Width	REF (18)	REF (22)	System
46	- 30°	4°	•	3.3 mm	818-211-46	822-211-46	Roth
36	- 30°	4°	•	3.3 mm	818-211-36	822-211-36	Roth
47	- 10°	0°	•	3.3 mm	818-212-47	822-212-47	McLaughlin-Bennett-Trevisi**
37	- 10°	0°	•	3.3 mm	818-212-37	822-212-37	McLaughlin-Bennett-Trevisi**
46	- 20°	0°	•	3.3 mm	818-212-46	822-212-46	McLaughlin-Bennett-Trevisi**
36	- 20°	0°	•	3.3 mm	818-212-36	822-212-36	McLaughlin-Bennett-Trevisi**
46	0°	0°	•	3.3 mm	818-215-46	822-215-46	Standard Edgewise
36	0°	0°	•	3.3 mm	818-215-36	822-215-36	Standard Edgewise
46	- 25°	0°	•	3.3 mm	818-219-46	822-219-46	
36	- 25°	0°	•	3.3 mm	818-219-36	822-219-36	



Slot 18: 0.46 mm x 0.64 mm / 18 x 25

Slot 22: 0.56 mm x 0.70 mm / 22 x 28

** The Dentaureum version of this prescription is not claimed to be a duplication of any other, nor does Dentaureum imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.

Ortho-Cast M-Series 18/22

convertible · Quantity: 10 pieces

Tooth	Torque	Distal Offset	Hook	Width	REF (18)	REF (22)	System
16	- 14°	14°	•	3.3 mm	818-111-16	822-111-16	Roth
26	- 14°	14°	•	3.3 mm	818-111-26	822-111-26	Roth
16	- 14°	10°	•	3.3 mm	818-112-16	822-112-16	McLaughlin-Bennett-Trevisi**
26	- 14°	10°	•	3.3 mm	818-112-26	822-112-26	McLaughlin-Bennett-Trevisi**



Tooth	Torque	Distal Offset	Hook	Width	REF (18)	REF (22)	System
46	- 30°	4°	•	3.3 mm	818-111-46	822-111-46	Roth
36	- 30°	4°	•	3.3 mm	818-111-36	822-111-36	Roth
46	- 20°	0°	•	3.3 mm	818-112-46	822-112-46	McLaughlin-Bennett-Trevisi**
36	- 20°	0°	•	3.3 mm	818-112-36	822-112-36	McLaughlin-Bennett-Trevisi**



Slot 18: 0.46 mm x 0.64 mm / 18 x 25

Slot 22: 0.56 mm x 0.70 mm / 22 x 28

** The Dentaureum version of this prescription is not claimed to be a duplication of any other, nor does Dentaureum imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.

Bonding technique

single rectangular



Ortho-Cast 18/22

non-convertible · Quantity: 10 pieces

Tooth	Torque	Distal Offset	Hook	Width	REF (18)	REF (22)	System
17-16 / 36-37	0°	6°	•	4.0 mm	724-092-51	724-094-51	
26-27 / 47-46	0°	6°	•	4.0 mm	724-093-51	724-095-51	
17-16 / 26-27 / 47-46 / 36-37	0°	6°		4.0 mm	724-090-51	724-091-51	

Slot 18: 0.46 mm x 0.64 mm / 18 x 25
Slot 22: 0.56 mm x 0.70 mm / 22 x 28



Ortho-Cast 18

non-convertible · with larger base · Quantity: 10 pieces

Tooth	Torque	Distal Offset	Hook	Width	REF (18)	REF (22)	System
17-16 / 26-27 / 47-46 / 36-37	0°	6°		4.0 mm	724-090-50		

Slot 18: 0.46 mm x 0.64 mm / 18 x 25



Ortho-Cast 18/22

non-convertible · distal notch · Quantity: 10 pieces

Tooth	Torque	Distal Offset	Hook	Width	REF (18)	REF (22)	System
17-16 / 36-37	0°	6°	•	3.3 mm	724-607-51	724-610-51	
26-27 / 47-46	0°	6°	•	3.3 mm	724-608-51	724-611-51	
17-16 / 26-27 / 47-46 / 36-37	0°	6°		3.3 mm	724-609-51	724-612-51	

Slot 18: 0.46 mm x 0.64 mm / 18 x 25
Slot 22: 0.56 mm x 0.70 mm / 22 x 28



Ortho 18/22

non-convertible · distal extension · Quantity: 10 pieces

Tooth	Torque	Distal Offset	Hook	Width	REF (18)	REF (22)	System
17-16 / 36-37	0°	0°	•	6.0 mm	725-164-50	725-166-50	Ricketts® Universal 18, Standard Edgewise
26-27 / 47-46	0°	0°	•	6.0 mm	725-165-50	725-167-50	Ricketts® Universal 18, Standard Edgewise
17-16 / 36-37	0°	6°	•	6.0 mm	725-168-50	725-170-50	
26-27 / 47-46	0°	6°	•	6.0 mm	725-169-50	725-171-50	

Slot 18: 0.46 mm x 0.64 mm / 18 x 25
Slot 22: 0.56 mm x 0.70 mm / 22 x 28

single rectangular

Ortho-Cast 18

convertible · Quantity: 10 pieces

Tooth	Torque	Distal Offset	Hook	Width	REF (18)	REF (22)	System
17-16	- 10°	8°	•	4.0 mm	724-048-51		
26-27	- 10°	8°	•	4.0 mm	724-049-51		

Slot 18: 0.46 mm x 0.64 mm / 18 x 25



double rectangular

Ortho-Cast M-Series 18/22

non-convertible · Quantity: 10 pieces

Tooth	Torque	Distal Offset	Hook	Width	REF (18)	REF (22)	System
16	- 14°	14°	•	3.3 mm	818-221-16	822-221-16	Roth
26	- 14°	14°	•	3.3 mm	818-221-26	822-221-26	Roth
16	- 14°	10°	•	3.3 mm	818-222-16	822-222-16	McLaughlin-Bennett-Trevisi**
26	- 14°	10°	•	3.3 mm	818-222-26	822-222-26	McLaughlin-Bennett-Trevisi**

Tooth	Torque	Distal Offset	Hook	Width	REF (18)	REF (22)	System
46	- 30°	4°	•	3.3 mm	818-221-46	822-221-46	Roth
36	- 30°	4°	•	3.3 mm	818-221-36	822-221-36	Roth
46	- 20°	0°	•	3.3 mm	818-222-46	822-222-46	McLaughlin-Bennett-Trevisi**
36	- 20°	0°	•	3.3 mm	818-222-36	822-222-36	McLaughlin-Bennett-Trevisi**

Slot 18: 0.46 mm x 0.64 mm / 18 x 25 · Slot auxiliary tube: 0.46 mm x 0.64 mm / 18 x 25

Slot 22: 0.56 mm x 0.70 mm / 22 x 28 · Slot auxiliary tube: 0.46 mm x 0.64 mm / 18 x 25

** The Dentaaurum version of this prescription is not claimed to be a duplication of any other, nor does Dentaaurum imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.



Ortho-Cast M-Series 18/22

convertible · Quantity: 10 pieces

Tooth	Torque	Distal Offset	Hook	Width	REF (18)	REF (22)	System
16	- 14°	14°	•	3.3 mm	818-121-16	822-121-16	Roth
26	- 14°	14°	•	3.3 mm	818-121-26	822-121-26	Roth
16	- 14°	10°	•	3.3 mm	818-122-16	822-122-16	McLaughlin-Bennett-Trevisi**
26	- 14°	10°	•	3.3 mm	818-122-26	822-122-26	McLaughlin-Bennett-Trevisi**

Tooth	Torque	Distal Offset	Hook	Width	REF (18)	REF (22)	System
46	- 30°	4°	•	3.3 mm	818-121-46	822-121-46	Roth
36	- 30°	4°	•	3.3 mm	818-121-36	822-121-36	Roth
46	- 20°	0°	•	3.3 mm	818-122-46	822-122-46	McLaughlin-Bennett-Trevisi**
36	- 20°	0°	•	3.3 mm	818-122-36	822-122-36	McLaughlin-Bennett-Trevisi**

Slot 18: 0.46 mm x 0.64 mm / 18 x 25 · Slot auxiliary tube: 0.46 mm x 0.64 mm / 18 x 25

Slot 22: 0.56 mm x 0.70 mm / 22 x 28 · Slot auxiliary tube: 0.46 mm x 0.64 mm / 18 x 25

** The Dentaaurum version of this prescription is not claimed to be a duplication of any other, nor does Dentaaurum imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.



Bonding technique

double rectangular



Ortho-Cast 18/22

non-convertible · with vertical hook · **Quantity: 10 pieces**

Tooth	Torque	Distal Offset	Hook	Width	REF (18)	REF (22)	System
17-16 / 26-27 / 47-46 / 36-37	0°	0°	•	3.5 mm	724-018-51	724-019-51	Ricketts® Universal 18, Standard Edgewise

Slot 18: 0.46 mm x 0.64 mm / 18 x 25 · Slot auxiliary tube: 0.46 mm x 0.64 mm / 18 x 25
Slot 22: 0.56 mm x 0.70 mm / 22 x 28 · Slot auxiliary tube: 0.56 mm x 0.70 mm / 22 x 28

double combination



Ortho-Cast 18/22

convertible · with lip bumper/facebow tube gingival · **Quantity: 10 pieces**

Tooth	Torque	Distal Offset	Hook	Width	REF (18)	REF (22)	System
47-46	- 20°	0°	•	4.0 mm	724-240-51	724-222-51	McLaughlin-Bennett- Trevisi** for the 1st molar
36-37	- 20°	0°	•	4.0 mm	724-241-51	724-223-51	McLaughlin-Bennett- Trevisi** for the 1st molar

Slot 18: 0.46 mm x 0.64 mm / 18 x 25 · Slot round tube: 1.15 mm / 45 (gingival)

Slot 22: 0.56 mm x 0.70 mm / 22 x 28 · Slot round tube: 1.15 mm / 45 (gingival)

** The Dentaaurum version of this prescription is not claimed to be a duplication of any other, nor does Dentaaurum imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.

double rectangular

rematitan® 18/22

non-convertible · Quantity: 10 pieces

Tooth	Torque	Distal Offset	Hook	Width	REF (18)	REF (22)	System
17-16 / 26-27	- 10°	0°		4.0 mm	723-020-00	723-021-00	

Tooth	Torque	Distal Offset	Hook	Width	REF (18)	REF (22)	System
47-46 / 36-37	- 22°	0°		4.0 mm	723-022-00	723-023-00	

Slot 18: 0.46 mm x 0.64 mm / 18 x 25 · Slot auxiliary tube: 0.46 mm x 0.64 mm / 18 x 25
 Slot 22: 0.56 mm x 0.70 mm / 22 x 28 · Slot auxiliary tube: 0.56 mm x 0.70 mm / 22 x 28



rematitan® 18/22

non-convertible · with vertical hook · Quantity: 10 pieces

Tooth	Torque	Distal Offset	Hook	Width	REF (18)	REF (22)	System
17-16	- 14°	10°	•	4.0 mm	723-028-00	723-029-00	McLaughlin-Bennett-Trevisi**
26-27	- 14°	10°	•	4.0 mm	723-032-00	723-033-00	McLaughlin-Bennett-Trevisi**

Tooth	Torque	Distal Offset	Hook	Width	REF (18)	REF (22)	System
47-46 / 36-37	- 20°	0°	•	4.0 mm	723-030-00	723-031-00	McLaughlin-Bennett-Trevisi** for the 1st molar
47-46 / 36-37	- 22°	0°	•	4.0 mm	723-026-00	723-027-00	

Tooth	Torque	Distal Offset	Hook	Width	REF (18)	REF (22)	System
17-16 / 26-27 / 47-46 / 36-37	- 10°	0°	•	4.0 mm	723-024-00	723-025-00	

Slot 18: 0.46 mm x 0.64 mm / 18 x 25 · Slot auxiliary tube: 0.46 mm x 0.64 mm / 18 x 25

Slot 22: 0.56 mm x 0.70 mm / 22 x 28 · Slot auxiliary tube: 0.56 mm x 0.70 mm / 22 x 28

** The Dentaaurum version of this prescription is not claimed to be a duplication of any other, nor does Dentaaurum imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.



Buccal tubes

Banding technique

single round



Ortho

non-convertible · with lip bumper/facebow tube · Quantity: 10 pieces

Tooth	Torque	Distal Offset	Hook	Width		REF	System
17-16 / 26-27 / 47-46 / 36-37		0°		6.0 mm		725-203-00	

Slot: 0.9 mm / 36



Ortho

non-convertible · with lip bumper/facebow tube · Quantity: 10 pieces

Tooth	Torque	Distal Offset	Hook	Width		REF	System
17-16 / 26-27 / 47-46 / 36-37		0°		6.0 mm		725-206-00	

Slot: 1.15 mm / 45

double round



Ortho

non-convertible · with lip bumper/facebow tube · Quantity: 10 pieces

Tooth	Torque	Distal Offset	Hook	Width		REF	System
17-16 / 26-27 / 47-46 / 36-37		0°		4.6 mm		726-011-00	

Slot: 0.9 mm / 36 · Slot round tube: 1.15 mm / 45 (occlusal)

single rectangular

Ortho-Cast M-Series 18/22

non-convertible · Quantity: 10 pieces

Tooth	Torque	Distal Offset	Hook	Width	REF (18)	REF (22)	System
17-16	- 14°	10°	•	3.3 mm	818-212-10	822-212-10	McLaughlin-Bennett-Trevisi**
26-27	- 14°	10°	•	3.3 mm	818-212-20	822-212-20	McLaughlin-Bennett-Trevisi**
17-16	- 10°	8°	•	3.3 mm	818-217-10	822-217-10	
26-27	- 10°	8°	•	3.3 mm	818-217-20	822-217-20	



Tooth	Torque	Distal Offset	Hook	Width	REF (18)	REF (22)	System
47	- 10°	0°	•	3.3 mm	818-512-40	822-512-40	McLaughlin-Bennett-Trevisi**
37	- 10°	0°	•	3.3 mm	818-512-30	822-512-30	McLaughlin-Bennett-Trevisi**
47-46	- 25°	8°	•	3.3 mm	818-217-40	822-217-40	
36-37	- 25°	8°	•	3.3 mm	818-217-30	822-217-30	
47-46	- 30°	12°	•	3.3 mm	818-218-40	822-218-40	
36-37	- 30°	12°	•	3.3 mm	818-218-30	822-218-30	



Slot 18: 0.46 mm x 0.64 mm / 18 x 25

Slot 22: 0.56 mm x 0.70 mm / 22 x 28

** The Dentaaurum version of this prescription is not claimed to be a duplication of any other, nor does Dentaaurum imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.

Ortho-Cast M-Series 18/22

convertible · Quantity: 10 pieces

Tooth	Torque	Distal Offset	Hook	Width	REF (18)	REF (22)	System
17-16	- 14°	14°	•	3.3 mm	818-111-10	822-111-10	Roth
26-27	- 14°	14°	•	3.3 mm	818-111-20	822-111-20	Roth
17-16	- 14°	10°	•	3.3 mm	818-112-10	822-112-10	McLaughlin-Bennett-Trevisi**
26-27	- 14°	10°	•	3.3 mm	818-112-20	822-112-20	McLaughlin-Bennett-Trevisi**



Tooth	Torque	Distal Offset	Hook	Width	REF (18)	REF (22)	System
47-46	- 30°	4°	•	3.3 mm	818-111-40	822-111-40	Roth
36-37	- 30°	4°	•	3.3 mm	818-111-30	822-111-30	Roth
47-46	- 20°	0°	•	3.3 mm	818-112-40	822-112-40	McLaughlin-Bennett-Trevisi** for the 1st molar
36-37	- 20°	0°	•	3.3 mm	818-112-30	822-112-30	McLaughlin-Bennett-Trevisi** for the 1st molar



Slot 18: 0.46 mm x 0.64 mm / 18 x 25

Slot 22: 0.56 mm x 0.70 mm / 22 x 28

** The Dentaaurum version of this prescription is not claimed to be a duplication of any other, nor does Dentaaurum imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.

Buccal tubes

Banding technique

single rectangular



Ortho-Cast 18/22

non-convertible · Quantity: 10 pieces

Tooth	Torque	Distal Offset	Hook	Width	REF (18)	REF (22)	System
17-16 / 36-37	0°	6°	•	4.0 mm	724-092-00	724-094-00	
26-27 / 47-46	0°	6°	•	4.0 mm	724-093-00	724-095-00	
17-16 / 26-27 / 47-46 / 36-37	0°	6°		4.0 mm	724-090-00	724-091-00	

Slot 18: 0.46 mm x 0.64 mm / 18 x 25
Slot 22: 0.56 mm x 0.70 mm / 22 x 28



Ortho-Cast 18

convertible · Quantity: 10 pieces

Tooth	Torque	Distal Offset	Hook	Width	REF (18)	REF (22)	System
17-16	- 10°	8°	•	4.0 mm	724-048-00		
26-27	- 10°	8°	•	4.0 mm	724-049-00		

Slot 18: 0.46 mm x 0.64 mm / 18 x 25



Ortho-Cast 18/22

non-convertible · distal notch · Quantity: 10 pieces

Tooth	Torque	Distal Offset	Hook	Width	REF (18)	REF (22)	System
17-16 / 36-37	0°	6°	•	3.3 mm	724-607-00	724-610-00	
26-27 / 47-46	0°	6°	•	3.3 mm	724-608-00	724-611-00	
17-16 / 26-27 / 47-46 / 36-37	0°	6°		3.3 mm	724-609-00	724-612-00	

Slot 18: 0.46 mm x 0.64 mm / 18 x 25
Slot 22: 0.56 mm x 0.70 mm / 22 x 28



Ortho 18/22

non-convertible · distal extension · Quantity: 10 pieces

Tooth	Torque	Distal Offset	Hook	Width	REF (18)	REF (22)	System
17-16 / 36-37	0°	0°	•	6.0 mm	725-164-00	725-166-00	Ricketts® Universal 18, Standard Edgewise
26-27 / 47-46	0°	0°	•	6.0 mm	725-165-00	725-167-00	Ricketts® Universal 18, Standard Edgewise
17-16 / 36-37	0°	6°	•	6.0 mm	725-168-00	725-170-00	
26-27 / 47-46	0°	6°	•	6.0 mm	725-169-00	725-171-00	

Slot 18: 0.46 mm x 0.64 mm / 18 x 25
Slot 22: 0.56 mm x 0.70 mm / 22 x 28

double rectangular

Ortho-Cast M-Series 18/22

non-convertible · Quantity: 10 pieces

Tooth	Torque	Distal Offset	Hook	Width	REF (18)	REF (22)	System
17-16	- 14°	14°	•	3.3 mm	818-221-10	822-221-10	Roth
26-27	- 14°	14°	•	3.3 mm	818-221-20	822-221-20	Roth
17-16	- 14°	10°	•	3.3 mm	818-222-10	822-222-10	McLaughlin-Bennett-Trevisi**
26-27	- 14°	10°	•	3.3 mm	818-222-20	822-222-20	McLaughlin-Bennett-Trevisi**
17-16	- 9°	10°	•	3.6 mm	818-227-10	822-227-10	
26-27	- 9°	10°	•	3.6 mm	818-227-20	822-227-20	



Tooth	Torque	Distal Offset	Hook	Width	REF (18)	REF (22)	System
47-46	- 30°	4°	•	3.3 mm	818-221-40	822-221-40	Roth
36-37	- 30°	4°	•	3.3 mm	818-221-30	822-221-30	Roth
47-46	- 20°	0°	•	3.3 mm	818-222-40	822-222-40	McLaughlin-Bennett-Trevisi** for the 1st molar
36-37	- 20°	0°	•	3.3 mm	818-222-30	822-222-30	McLaughlin-Bennett-Trevisi** for the 1st molar
47-46	- 22°	12°	•	3.6 mm	818-224-40	822-224-40	Ricketts® 18 - IV. dimension for the 1st molar
36-37	- 22°	12°	•	3.6 mm	818-224-30	822-224-30	Ricketts® 18 - IV. dimension for the 1st molar



Slot 18: 0.46 mm x 0.64 mm / 18 x 25 · Slot auxiliary tube: 0.46 mm x 0.64 mm / 18 x 25

Slot 22: 0.56 mm x 0.70 mm / 22 x 28 · Slot auxiliary tube: 0.46 mm x 0.64 mm / 18 x 25

** The Dentaureum version of this prescription is not claimed to be a duplication of any other, nor does Dentaureum imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.

Ortho-Cast M-Series 18/22

convertible · Quantity: 10 pieces

Tooth	Torque	Distal Offset	Hook	Width	REF (18)	REF (22)	System
17-16	- 14°	14°	•	3.3 mm	818-121-10	822-121-10	Roth
26-27	- 14°	14°	•	3.3 mm	818-121-20	822-121-20	Roth
17-16	- 14°	10°	•	3.3 mm	818-122-10	822-122-10	McLaughlin-Bennett-Trevisi**
26-27	- 14°	10°	•	3.3 mm	818-122-20	822-122-20	McLaughlin-Bennett-Trevisi**



Tooth	Torque	Distal Offset	Hook	Width	REF (18)	REF (22)	System
47-46	0°	0°	•	3.6 mm	818-125-40	822-125-40	Ricketts® Universal 18, Standard Edgewise
36-37	0°	0°	•	3.6 mm	818-125-30	822-125-30	Ricketts® Universal 18, Standard Edgewise



Slot 18: 0.46 mm x 0.64 mm / 18 x 25 · Slot auxiliary tube: 0.46 mm x 0.64 mm / 18 x 25

Slot 22: 0.56 mm x 0.70 mm / 22 x 28 · Slot auxiliary tube: 0.46 mm x 0.64 mm / 18 x 25

** The Dentaureum version of this prescription is not claimed to be a duplication of any other, nor does Dentaureum imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.

Banding technique

double rectangular



Ortho-Cast M-Series 18/22

convertible · Auxiliary tube placed in gingival/buccal direction · Quantity: 10 pieces

Tooth	Torque	Distal Offset	Hook	Width	REF (18)	REF (22)	System
47-46	- 30°	4°	•	3.6 mm	818-151-40	822-151-40	Roth
36-37	- 30°	4°	•	3.6 mm	818-151-30	822-151-30	Roth
47-46	- 20°	0°	•	3.6 mm	818-152-40	822-152-40	McLaughlin-Bennett-Trevisi** for the 1st molar
36-37	- 20°	0°	•	3.6 mm	818-152-30	822-152-30	McLaughlin-Bennett-Trevisi** for the 1st molar

Slot 18: 0.46 mm x 0.64 mm / 18 x 25 · Slot auxiliary tube: 0.46 mm x 0.64 mm / 18 x 25

Slot 22: 0.56 mm x 0.70 mm / 22 x 28 · Slot auxiliary tube: 0.46 mm x 0.64 mm / 18 x 25

** The Dentaaurum version of this prescription is not claimed to be a duplication of any other, nor does Dentaaurum imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.



Ortho-Cast 18/22

non-convertible · with vertical hook · Quantity: 10 pieces

Tooth	Torque	Distal Offset	Hook	Width	REF (18)	REF (22)	System
47-46 / 36-37	- 30°	0°	•	4.0 mm	724-616-00	724-617-00	
Tooth	Torque	Distal Offset	Hook	Width	REF (18)	REF (22)	System
17-16 / 26-27 / 47-46 / 36-37	0°	0°	•	3.5 mm	724-018-00	724-019-00	Ricketts® Universal 18, Standard Edgewise

Slot 18: 0.46 mm x 0.64 mm / 18 x 25 · Slot auxiliary tube: 0.46 mm x 0.64 mm / 18 x 25

Slot 22: 0.56 mm x 0.70 mm / 22 x 28 · Slot auxiliary tube: 0.56 mm x 0.70 mm / 22 x 28



Ortho-Cast 18/22

convertible · Quantity: 10 pieces

Tooth	Torque	Distal Offset	Hook	Width	REF (18)	REF (22)	System
47-46	- 25°	8°	•	4.0 mm	724-056-00	724-156-00	
36-37	- 25°	8°	•	4.0 mm	724-057-00	724-157-00	

Slot 18: 0.46 mm x 0.64 mm / 18 x 25 · Slot auxiliary tube: 0.46 mm x 0.64 mm / 18 x 25

Slot 22: 0.56 mm x 0.70 mm / 22 x 28 · Slot auxiliary tube: 0.46 mm x 0.64 mm / 18 x 25

double combination

Ortho-Cast M-Series 18/22

convertible · with lip bumper/facebow tube occlusal · **Quantity: 10 pieces**

Tooth	Torque	Distal Offset	Hook	Width	REF (18)	REF (22)	System
17-16	- 14°	14°	•	3.6 mm	818-141-10	822-141-10	Roth
26-27	- 14°	14°	•	3.6 mm	818-141-20	822-141-20	Roth
17-16	- 14°	10°	•	3.6 mm		822-142-10	McLaughlin-Bennett-Trevisi**
26-27	- 14°	10°	•	3.6 mm		822-142-20	McLaughlin-Bennett-Trevisi**

Slot 18: 0.46 mm x 0.64 mm / 18 x 25 · Slot round tube: 1.15 mm / 45 (occlusal)

Slot 22: 0.56 mm x 0.70 mm / 22 x 28 · Slot round tube: 1.15 mm / 45 (occlusal)

** The Dentaurem version of this prescription is not claimed to be a duplication of any other, nor does Dentaurem imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.



Ortho-Cast M-Series 18/22

convertible · with lip bumper/facebow tube gingival · **Quantity: 10 pieces**

Tooth	Torque	Distal Offset	Hook	Width	REF (18)	REF (22)	System
47-46	- 25°	8°	•	3.6 mm	818-147-40	822-147-40	
36-37	- 25°	8°	•	3.6 mm	818-147-30	822-147-30	

Slot 18: 0.46 mm x 0.64 mm / 18 x 25 · Slot round tube: 1.15 mm / 45 (gingival)

Slot 22: 0.56 mm x 0.70 mm / 22 x 28 · Slot round tube: 1.15 mm / 45 (gingival)



Ortho-Cast 18/22

convertible · with lip bumper/facebow tube gingival · **Quantity: 10 pieces**

Tooth	Torque	Distal Offset	Hook	Width	REF (18)	REF (22)	System
47-46	- 20°	0°	•	4.0 mm	724-240-00	724-222-00	McLaughlin-Bennett-Trevisi** for the 1st molar
36-37	- 20°	0°	•	4.0 mm	724-241-00	724-223-00	McLaughlin-Bennett-Trevisi** for the 1st molar

Slot 18: 0.46 mm x 0.64 mm / 18 x 25 · Slot round tube: 1.15 mm / 45 (gingival)

Slot 22: 0.56 mm x 0.70 mm / 22 x 28 · Slot round tube: 1.15 mm / 45 (gingival)

** The Dentaurem version of this prescription is not claimed to be a duplication of any other, nor does Dentaurem imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.



Ortho 18/22

non-convertible · with lip bumper/facebow tube occlusal · **Quantity: 10 pieces**

Tooth	Torque	Distal Offset	Hook	Width	REF (18)	REF (22)	System
17-16	0°	0°	•	4.6 mm	727-365-00	727-385-00	Ricketts® Universal 18, Standard Edgewise
26-27	0°	0°	•	4.6 mm	727-366-00	727-386-00	Ricketts® Universal 18, Standard Edgewise
17-16 / 26-27	0°	0°		4.6 mm	727-265-00	727-266-00	Ricketts® Universal 18, Standard Edgewise
17-16	0°	6°	•	4.6 mm	727-367-00	727-387-00	
26-27	0°	6°	•	4.6 mm	727-368-00	727-388-00	

Slot 18: 0.46 mm x 0.64 mm / 18 x 25 · Slot round tube: 1.15 mm / 45 (occlusal)

Slot 22: 0.56 mm x 0.70 mm / 22 x 28 · Slot round tube: 1.15 mm / 45 (occlusal)



Banding technique

double combination



Ortho 18/22

non-convertible · with lip bumper/facebow tube gingival · Quantity: 10 pieces

Tooth	Torque	Distal Offset	Hook	Width	REF (18)	REF (22)	System
47-46	0°	0°	•	4.6 mm	727-267-00	727-269-00	Ricketts® Universal 18, Standard Edgewise
36-37	0°	0°	•	4.6 mm	727-268-00	727-270-00	Ricketts® Universal 18, Standard Edgewise

Slot 18: 0.46 mm x 0.64 mm / 18 x 25 · Slot round tube: 1.15 mm / 45 (gingival)
Slot 22: 0.56 mm x 0.70 mm / 22 x 28 · Slot round tube: 1.15 mm / 45 (gingival)

triple combination



Ortho-Cast M-Series 18/22

convertible · with lip bumper/facebow tube occlusal · Quantity: 10 pieces

Tooth	Torque	Distal Offset	Hook	Width	REF (18)	REF (22)	System
17-16	- 14°	14°	•	3.6 mm	818-131-10	822-131-10	Roth
26-27	- 14°	14°	•	3.6 mm	818-131-20	822-131-20	Roth
17-16	- 14°	10°	•	3.6 mm	818-132-10	822-132-10	McLaughlin-Bennett-Trevisi**
26-27	- 14°	10°	•	3.6 mm	818-132-20	822-132-20	McLaughlin-Bennett-Trevisi**
17-16	- 10°	8°	•	3.6 mm		822-137-10	
26-27	- 10°	8°	•	3.6 mm		822-137-20	

Slot 18: 0.46 mm x 0.64 mm / 18 x 25 · Slot auxiliary tube: 0.46 mm x 0.64 mm / 18 x 25 · Slot round tube: 1.15 mm / 45 (occlusal)
Slot 22: 0.56 mm x 0.70 mm / 22 x 28 · Slot auxiliary tube: 0.46 mm x 0.64 mm / 18 x 25 · Slot round tube: 1.15 mm / 45 (occlusal)

** The Dentaurem version of this prescription is not claimed to be a duplication of any other, nor does Dentaurem imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.



Ortho-Cast 18/22

non-convertible · with lip bumper/facebow tube central · Quantity: 10 pieces

Tooth	Torque	Distal Offset	Hook	Width	REF (18)	REF (22)	System
17-16	0°	0°	•	3.6 mm	724-600-00		Ricketts® Universal 18, Standard Edgewise
26-27	0°	0°	•	3.6 mm	724-601-00		Ricketts® Universal 18, Standard Edgewise
17-16	0°	6°	•	3.6 mm	724-602-00	724-613-00	
26-27	0°	6°	•	3.6 mm	724-603-00	724-614-00	
17-16	- 9°	10°	•	3.6 mm	724-618-00	724-620-00	
26-27	- 9°	10°	•	3.6 mm	724-619-00	724-621-00	

Slot 18: 0.46 mm x 0.64 mm / 18 x 25 · Slot auxiliary tube: 0.46 mm x 0.64 mm / 18 x 25 · Slot round tube: 1.15 mm / 45 (central)
Slot 22: 0.56 mm x 0.70 mm / 22 x 28 · Slot auxiliary tube: 0.56 mm x 0.70 mm / 22 x 28 · Slot round tube: 1.15 mm / 45 (central)



Tooth	Torque	Distal Offset	Hook	Width	REF (18)	REF (22)	System
47-46	- 25°	8°	•	4.0 mm	724-060-00	724-160-00	
36-37	- 25°	8°	•	4.0 mm	724-061-00	724-161-00	

Slot 18: 0.46 mm x 0.64 mm / 18 x 25 · Slot auxiliary tube: 0.46 mm x 0.64 mm / 18 x 25 · Slot round tube: 1.15 mm / 45 (central)
Slot 22: 0.56 mm x 0.70 mm / 22 x 28 · Slot auxiliary tube: 0.46 mm x 0.64 mm / 18 x 25 · Slot round tube: 1.15 mm / 45 (central)



Tooth	Torque	Distal Offset	Hook	Width	REF (18)	REF (22)	System
17-16 / 26-27 / 47-46 / 36-37	0°	0°		3.6 mm	724-604-00	724-615-00	Ricketts® Universal 18, Standard Edgewise
17-16 / 26-27 / 47-46 / 36-37	0°	6°		3.6 mm	724-605-00		

Slot 18: 0.46 mm x 0.64 mm / 18 x 25 · Slot auxiliary tube: 0.46 mm x 0.64 mm / 18 x 25 · Slot round tube: 1.15 mm / 45 (central)
Slot 22: 0.56 mm x 0.70 mm / 22 x 28 · Slot auxiliary tube: 0.56 mm x 0.70 mm / 22 x 28 · Slot round tube: 1.15 mm / 45 (central)

triple combination

Ortho-Cast 18

non-convertible · with mesial/distal hook · with lip bumper/facebow tube central · **Quantity: 10 pieces**

Tooth	Torque	Distal Offset	Hook	Width	REF (18)	REF (22)	System
17-16 / 26-27 / 47-46 / 36-37	0°	0°	•	3.6 mm	724-606-00		Ricketts® Universal 18, Standard Edgewise



Slot 18: 0.46 mm x 0.64 mm / 18 x 25 · Slot auxiliary tube: 0.46 mm x 0.64 mm / 18 x 25 · Slot round tube: 1.15 mm / 45 (central)

Ortho-Cast 18/22

convertible · with lip bumper/facebow tube occlusal · **Quantity: 10 pieces**

Tooth	Torque	Distal Offset	Hook	Width	REF (18)	REF (22)	System
17-16	0°	8°	•	4.0 mm	724-006-00	724-106-00	
26-27	0°	8°	•	4.0 mm	724-007-00	724-107-00	



Slot 18: 0.46 mm x 0.64 mm / 18 x 25 · Slot auxiliary tube: 0.46 mm x 0.64 mm / 18 x 25 · Slot round tube: 1.15 mm / 45 (occlusal)
Slot 22: 0.56 mm x 0.70 mm / 22 x 28 · Slot auxiliary tube: 0.46 mm x 0.64 mm / 18 x 25 · Slot round tube: 1.15 mm / 45 (occlusal)



Buccal tubes

Prewelded parts



Lingual / Palatal sheath

for Orthorama®, Lingual/Palatal bars, Goshgarian, Quad Helix · open

Tooth	Torque	Distal Offset	In/Out	Hook	Width	REF	Quantity
17-16 / 26-27 / 47-46 / 36-37		0°	0.25 mm		3.5 mm	728-010-00	10 pieces

Slot: 2 x 0.9 mm / 36



Lingual/Palatal sheaths

for Orthorama®, Lingual/Palatal bars, Goshgarian, Quad Helix · closed

Tooth	Torque	Distal Offset	In/Out	Hook	Width	REF	Quantity
17-16 / 26-27 / 47-46 / 36-37		0°	0.25 mm		3.5 mm	728-110-00	10 pieces

Slot: 2 x 0.9 mm / 36



Lingual/Palatal sheaths

for Orthorama®, Lingual/Palatal bars, Goshgarian, Quad Helix · open · with mesial/distal hook

Tooth	Torque	Distal Offset	In/Out	Hook	Width	REF	Quantity
17-16 / 36-37		0°	0.25 mm	•	3.5 mm	728-011-00	10 pieces
26-27 / 47-46		0°	0.25 mm	•	3.5 mm	728-012-00	10 pieces

Slot: 2 x 0.9 mm / 36



Lingual/Palatal sheaths

for Orthorama®, Lingual/Palatal bars, Goshgarian, Quad Helix · closed · with mesial/distal hook

Tooth	Torque	Distal Offset	In/Out	Hook	Width	REF	Quantity
17-16 / 36-37		0°	0.25 mm	•	3.5 mm	728-111-00	10 pieces
26-27 / 47-46		0°	0.25 mm	•	3.5 mm	728-112-00	10 pieces

Slot: 2 x 0.9 mm / 36



Lingual/Palatal sheaths

for Orthorama®, Lingual/Palatal bars, Goshgarian, Quad Helix · open · with distal hook

Tooth	Torque	Distal Offset	In/Out	Hook	Width	REF	Quantity
17-16 / 36-37		0°	0.25 mm	•	3.5 mm	728-016-00	10 pieces
26-27 / 47-46		0°	0.25 mm	•	3.5 mm	728-017-00	10 pieces

Slot: 2 x 0.9 mm / 36



Lingual/Palatal sheaths easy in

for Orthorama®, Lingual/Palatal bars, Goshgarian, Quad Helix · closed · with insertion groove

An insertion groove facilitates the insertion and removal of Orthorama®, lingual and palatal arches, Goshgarian and Quad Helix into/from the lingual palatal sheaths easy-in. The anterior teeth pose no hindrance. The insertion groove should always face the occlusal plane during welding so that the arches can be inserted horizontally in the occlusal plane.

Tooth	Torque	Distal Offset	In/Out	Hook	Width	REF	Quantity
17-16 / 36-37		0°	0.25 mm		3.5 mm	728-113-00	10 pieces
26-27 / 47-46		0°	0.25 mm		3.5 mm	728-114-00	10 pieces

Slot: 2 x 0.9 mm / 36

Button / Banding technique

For elastics or Dentalastics®.

Version	REF	Quantity
round base, long neck	750-002-00	10 pieces
broad base, short neck	750-101-00	10 pieces



Band seating lugs / Banding technique

For all preformed bands.


Form of delivery	REF	Quantity
on strips	750-201-00	100 pieces



Double hooks / Banding technique

For additional elastics or tension springs.

Version	Form of delivery	REF	Quantity
narrow, for cuspid and bicuspid bands	on strips	750-602-00	10 pieces
narrow, for cuspid and bicuspid bands	single	750-602-10	10 pieces
wide, for molar bands	on strips	750-603-00	10 pieces
wide, for molar bands	single	750-603-10	10 pieces



Hooks with seating lug / Banding technique

For elastics, ligatures, etc.

Form of delivery	REF	Quantity
single	750-401-00	100 pieces



dentaform® Snap

Assortments	Page	142
Single packs	Page	143

dentaform®

Assortments	Page	147
Single packs	Page	149

Standard

Assortments	Page	157
Single packs	Page	161

Accessories	Page	168
-------------	------	-----



FDI System

18	17	16	15	14	13	12	11	21	22	23	24	25	26	27	28		
8	7	6	5	4	3	2	1	1	2	3	4	5	6	7	8		
R	8	7	6	5	4	3	2	1	1	2	3	4	5	6	7	8	L
	48	47	46	45	44	43	42	41	31	32	33	34	35	36	37	38	

Bands.

Modern band systems for satisfied orthodontists and patients.



If you wish to order individually pre-welded bands, please contact your local representative.

**DENTAURUM
SERVICE**



Best quality.

Easy and precise placement, elasticity and anatomical fit are fundamental requirements for an orthodontic band. In band production, Dentaaurum places the highest value on the best quality and fit, while using the best materials. All bands carry a permanent, biocompatible laser marking. The bands have a high-luster polishing on the outside and are roughened by a special process on the inside to improve adhesion.

The anatomically contoured dentaform® bands are made of a very flexible material that can be adapted perfectly to the tooth. dentaform® Snap bands consist of a harder material for optimal stability and permanent fit. The high edge stability prevents unwanted deformation of the band periphery. In addition, we offer standard bands for universal use.

dentaform® Snap

The bands with the super snap effect.

dentaform®

The bands which are easy to adapt.

Standard

The universal standard bands.



Pre-welding service – banding technique.

Reliable and quick.

Optimal service.

If required, all Dentaureum brackets, buccal tubes and palatal auxiliaries can be welded to Dentaureum bands by our welding service. Our competent team will assist you with your order and ensure quick processing. Computer-assisted recording of individual data enables fast and reliable delivery.



Quick and easy ordering.

- **Ask for the pre-welding order form.**

To ensure quick processing, please always use our welding service forms.

- **Personal card for welding service with ordering variant.**

Your personal card is drawn up to your own requirements and specifications. We will be pleased to assist you at any time.

Pre-welded goods from stock.

Use our pre-welded goods on stock* to further shorten delivery time.

- **Finished bands.**

You will receive your welded bands packed in mini grip bags with a complete statement of the contents.

* Maxilla/mandible first molar, right/left in band quality dentaform® Snap/dentaform® welded acc. to Roth 18/22 and maxilla/mandible first and second premolar right/left in band quality Standard welded acc. to Roth 18/22.

dentaform® Snap / Assortments

- with super snap effect
- hard, stiff quality
- anatomically shaped
- very accurate fit
- ideal for the extra-oral technique (facebow) or SUS
- secure fit and optimal adhesion due to rough interior surface
- delivery in a tray (fits in band cabinet REF 890-414-00)
- individual welding, e.g. acc. to McLaughlin-Bennett-Trevisi* 18/22, is possible on request

We recommend our dental vernier (REF 042-751-00) for determining the largest mesio-distal diameter of the tooth crown.

Upper first molar bands

Mini

Tooth	Band sizes (Contents per size)	welded		unwelded	Quantity
		REF (Roth 18) ¹⁾	REF (Roth 22) ²⁾	REF	
16	No. 5-36 (2)	878-364-18	878-364-22	878-064-00	64 pieces
26	No. 5-36 (2)	879-364-18	879-364-22	879-064-00	64 pieces

Medium

Tooth	Band sizes (Contents per size)	welded		unwelded	Quantity
		REF (Roth 18) ¹⁾	REF (Roth 22) ²⁾	REF	
16	No. 5-8, 30-36 (2), 9-18 (4), 19-20 (6), 21-25, 28-29 (8), 26-27 (10)	878-350-18	878-350-22	878-150-00	150 pieces
26	No. 5-8, 30-36 (2), 9-18 (4), 19-20 (6), 21-25, 28-29 (8), 26-27 (10)	879-350-18	879-350-22	879-150-00	150 pieces

¹⁾ (16) welded with buccal tubes 818-131-10 and double hook 750-602-10
 (26) welded with buccal tubes 818-131-20 and double hook 750-602-10

²⁾ (16) welded with buccal tubes 822-131-10 and double hook 750-602-10
 (26) welded with buccal tubes 822-131-20 and double hook 750-602-10

* The Dentaform version of this prescription is not claimed to be a duplication of any other, nor does Dentaform imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.

Lower first molar bands

Mini

Tooth	Band sizes (Contents per size)	welded		unwelded	Quantity
		REF (Roth 18) ¹⁾	REF (Roth 22) ²⁾	REF	
46	No. 5-36 (2)	880-364-18	880-364-22	880-264-00	64 pieces
36	No. 5-36 (2)	881-364-18	881-364-22	881-264-00	64 pieces

Medium

Tooth	Band sizes (Contents per size)	welded		unwelded	Quantity
		REF (Roth 18) ¹⁾	REF (Roth 22) ²⁾	REF	
46	No. 5-8, 30-36 (2), 9-18 (4), 19-20 (6), 21-25, 28-29 (8), 26-27 (10)	880-350-18	880-350-22	880-250-00	150 pieces
36	No. 5-8, 30-36 (2), 9-18 (4), 19-20 (6), 21-25, 28-29 (8), 26-27 (10)	881-350-18	881-350-22	881-250-00	150 pieces

¹⁾ (36) welded with buccal tubes 818-151-30 and double hook 750-602-10
 (46) welded with buccal tubes 818-151-40 and double hook 750-602-10

²⁾ (36) welded with buccal tubes 822-151-30 and double hook 750-602-10
 (46) welded with buccal tubes 822-151-40 and double hook 750-602-10

* The Dentaform version of this prescription is not claimed to be a duplication of any other, nor does Dentaform imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.

dentaform® Snap / Single packs

- with super snap effect
- hard, stiff quality
- anatomically shaped
- very accurate fit
- ideal for the extra-oral technique (facebow) or SUS
- secure fit and optimal adhesion due to rough interior surface
- individual welding, e.g. acc. to McLaughlin-Bennett-Trevisi* 18/22, is possible on request

We recommend our dental vernier (REF 042-751-00) for determining the largest mesio-distal diameter of the tooth crown.

Upper first molar bands, right

Tooth	Band size	Band thickness	ø mesial-distal	welded		Quantity	unwelded	Quantity
				REF (Roth 18) ¹⁾	REF (Roth 22) ²⁾		REF	
16	1	0.12 mm	8.78 mm	878-301-18	878-301-22	1 piece	878-001-00	10 pieces
16	2	0.12 mm	8.89 mm	878-302-18	878-302-22	1 piece	878-002-00	10 pieces
16	3	0.12 mm	8.99 mm	878-303-18	878-303-22	1 piece	878-003-00	10 pieces
16	4	0.12 mm	9.09 mm	878-304-18	878-304-22	1 piece	878-004-00	10 pieces
16	5	0.12 mm	9.20 mm	878-305-18	878-305-22	1 piece	878-005-00	10 pieces
16	6	0.12 mm	9.30 mm	878-306-18	878-306-22	1 piece	878-006-00	10 pieces
16	7	0.12 mm	9.40 mm	878-307-18	878-307-22	1 piece	878-007-00	10 pieces
16	8	0.12 mm	9.50 mm	878-308-18	878-308-22	1 piece	878-008-00	10 pieces
16	9	0.15 mm	9.61 mm	878-309-18	878-309-22	1 piece	878-009-00	10 pieces
16	10	0.15 mm	9.71 mm	878-310-18	878-310-22	1 piece	878-010-00	10 pieces
16	11	0.15 mm	9.81 mm	878-311-18	878-311-22	1 piece	878-011-00	10 pieces
16	12	0.15 mm	9.92 mm	878-312-18	878-312-22	1 piece	878-012-00	10 pieces
16	13	0.15 mm	10.02 mm	878-313-18	878-313-22	1 piece	878-013-00	10 pieces
16	14	0.15 mm	10.12 mm	878-314-18	878-314-22	1 piece	878-014-00	10 pieces
16	15	0.15 mm	10.23 mm	878-315-18	878-315-22	1 piece	878-015-00	10 pieces
16	16	0.15 mm	10.33 mm	878-316-18	878-316-22	1 piece	878-016-00	10 pieces
16	17	0.15 mm	10.43 mm	878-317-18	878-317-22	1 piece	878-017-00	10 pieces
16	18	0.15 mm	10.54 mm	878-318-18	878-318-22	1 piece	878-018-00	10 pieces
16	19	0.15 mm	10.64 mm	878-319-18	878-319-22	1 piece	878-019-00	10 pieces
16	20	0.15 mm	10.74 mm	878-320-18	878-320-22	1 piece	878-020-00	10 pieces
16	21	0.15 mm	10.85 mm	878-321-18	878-321-22	1 piece	878-021-00	10 pieces
16	22	0.15 mm	10.95 mm	878-322-18	878-322-22	1 piece	878-022-00	10 pieces
16	23	0.15 mm	11.05 mm	878-323-18	878-323-22	1 piece	878-023-00	10 pieces
16	24	0.15 mm	11.16 mm	878-324-18	878-324-22	1 piece	878-024-00	10 pieces
16	25	0.18 mm	11.26 mm	878-325-18	878-325-22	1 piece	878-025-00	10 pieces
16	26	0.18 mm	11.36 mm	878-326-18	878-326-22	1 piece	878-026-00	10 pieces
16	27	0.18 mm	11.47 mm	878-327-18	878-327-22	1 piece	878-027-00	10 pieces
16	28	0.18 mm	11.57 mm	878-328-18	878-328-22	1 piece	878-028-00	10 pieces
16	29	0.18 mm	11.67 mm	878-329-18	878-329-22	1 piece	878-029-00	10 pieces
16	30	0.18 mm	11.78 mm	878-330-18	878-330-22	1 piece	878-030-00	10 pieces
16	31	0.18 mm	11.88 mm	878-331-18	878-331-22	1 piece	878-031-00	10 pieces
16	32	0.18 mm	11.98 mm	878-332-18	878-332-22	1 piece	878-032-00	10 pieces
16	33	0.18 mm	12.08 mm	878-333-18	878-333-22	1 piece	878-033-00	10 pieces
16	34	0.18 mm	12.19 mm	878-334-18	878-334-22	1 piece	878-034-00	10 pieces
16	35	0.18 mm	12.29 mm	878-335-18	878-335-22	1 piece	878-035-00	10 pieces
16	36	0.18 mm	12.39 mm	878-336-18	878-336-22	1 piece	878-036-00	10 pieces



¹⁾ welded with buccal tubes 818-131-10 and double hook 750-602-10

²⁾ welded with buccal tubes 822-131-10 and double hook 750-602-10

* The Dentaform version of this prescription is not claimed to be a duplication of any other, nor does Dentaform imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.

dentaform® Snap / Single packs

- with super snap effect
- hard, stiff quality
- anatomically shaped
- very accurate fit
- ideal for the extra-oral technique (facebow) or SUS
- secure fit and optimal adhesion due to rough interior surface
- individual welding, e.g. acc. to McLaughlin-Bennett-Trevisi* 18/22, is possible on request

We recommend our dental vernier (REF 042-751-00) for determining the largest mesio-distal diameter of the tooth crown.

Upper first molar bands, left

Tooth	Band size	Band thickness	ø mesial-distal	welded		Quantity	unwelded	Quantity
				REF (Roth 18) ¹⁾	REF (Roth 22) ²⁾		REF	
26	1	0.12 mm	8.78 mm	879-301-18	879-301-22	1 piece	879-001-00	10 pieces
26	2	0.12 mm	8.89 mm	879-302-18	879-302-22	1 piece	879-002-00	10 pieces
26	3	0.12 mm	8.99 mm	879-303-18	879-303-22	1 piece	879-003-00	10 pieces
26	4	0.12 mm	9.09 mm	879-304-18	879-304-22	1 piece	879-004-00	10 pieces
26	5	0.12 mm	9.20 mm	879-305-18	879-305-22	1 piece	879-005-00	10 pieces
26	6	0.12 mm	9.30 mm	879-306-18	879-306-22	1 piece	879-006-00	10 pieces
26	7	0.12 mm	9.40 mm	879-307-18	879-307-22	1 piece	879-007-00	10 pieces
26	8	0.12 mm	9.50 mm	879-308-18	879-308-22	1 piece	879-008-00	10 pieces
26	9	0.15 mm	9.61 mm	879-309-18	879-309-22	1 piece	879-009-00	10 pieces
26	10	0.15 mm	9.71 mm	879-310-18	879-310-22	1 piece	879-010-00	10 pieces
26	11	0.15 mm	9.81 mm	879-311-18	879-311-22	1 piece	879-011-00	10 pieces
26	12	0.15 mm	9.92 mm	879-312-18	879-312-22	1 piece	879-012-00	10 pieces
26	13	0.15 mm	10.02 mm	879-313-18	879-313-22	1 piece	879-013-00	10 pieces
26	14	0.15 mm	10.12 mm	879-314-18	879-314-22	1 piece	879-014-00	10 pieces
26	15	0.15 mm	10.23 mm	879-315-18	879-315-22	1 piece	879-015-00	10 pieces
26	16	0.15 mm	10.33 mm	879-316-18	879-316-22	1 piece	879-016-00	10 pieces
26	17	0.15 mm	10.43 mm	879-317-18	879-317-22	1 piece	879-017-00	10 pieces
26	18	0.15 mm	10.54 mm	879-318-18	879-318-22	1 piece	879-018-00	10 pieces
26	19	0.15 mm	10.64 mm	879-319-18	879-319-22	1 piece	879-019-00	10 pieces
26	20	0.15 mm	10.74 mm	879-320-18	879-320-22	1 piece	879-020-00	10 pieces
26	21	0.15 mm	10.85 mm	879-321-18	879-321-22	1 piece	879-021-00	10 pieces
26	22	0.15 mm	10.95 mm	879-322-18	879-322-22	1 piece	879-022-00	10 pieces
26	23	0.15 mm	11.05 mm	879-323-18	879-323-22	1 piece	879-023-00	10 pieces
26	24	0.15 mm	11.16 mm	879-324-18	879-324-22	1 piece	879-024-00	10 pieces
26	25	0.18 mm	11.26 mm	879-325-18	879-325-22	1 piece	879-025-00	10 pieces
26	26	0.18 mm	11.36 mm	879-326-18	879-326-22	1 piece	879-026-00	10 pieces
26	27	0.18 mm	11.47 mm	879-327-18	879-327-22	1 piece	879-027-00	10 pieces
26	28	0.18 mm	11.57 mm	879-328-18	879-328-22	1 piece	879-028-00	10 pieces
26	29	0.18 mm	11.67 mm	879-329-18	879-329-22	1 piece	879-029-00	10 pieces
26	30	0.18 mm	11.78 mm	879-330-18	879-330-22	1 piece	879-030-00	10 pieces
26	31	0.18 mm	11.88 mm	879-331-18	879-331-22	1 piece	879-031-00	10 pieces
26	32	0.18 mm	11.98 mm	879-332-18	879-332-22	1 piece	879-032-00	10 pieces
26	33	0.18 mm	12.08 mm	879-333-18	879-333-22	1 piece	879-033-00	10 pieces
26	34	0.18 mm	12.19 mm	879-334-18	879-334-22	1 piece	879-034-00	10 pieces
26	35	0.18 mm	12.29 mm	879-335-18	879-335-22	1 piece	879-035-00	10 pieces
26	36	0.18 mm	12.39 mm	879-336-18	879-336-22	1 piece	879-036-00	10 pieces

¹⁾ welded with buccal tubes 818-131-20 and double hook 750-602-10

²⁾ welded with buccal tubes 822-131-20 and double hook 750-602-10



* The Dentaform version of this prescription is not claimed to be a duplication of any other, nor does Dentaform imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.

dentaform® Snap / Single packs

- with super snap effect
- hard, stiff quality
- anatomically shaped
- very accurate fit
- ideal for the extra-oral technique (facebow) or SUS
- secure fit and optimal adhesion due to rough interior surface
- individual welding, e.g. acc. to McLaughlin-Bennett-Trevisi* 18/22, is possible on request

We recommend our dental vernier (REF 042-751-00) for determining the largest mesio-distal diameter of the tooth crown.

Lower first molar bands, right

Tooth	Band size	Band thickness	ø mesial-distal	welded		Quantity	unwelded	Quantity
				REF (Roth 18) ¹⁾	REF (Roth 22) ²⁾		REF	
46	1	0.12 mm	8.91 mm	880-301-18	880-301-22	1 piece	880-201-00	10 pieces
46	2	0.12 mm	9.03 mm	880-302-18	880-302-22	1 piece	880-202-00	10 pieces
46	3	0.12 mm	9.15 mm	880-303-18	880-303-22	1 piece	880-203-00	10 pieces
46	4	0.12 mm	9.28 mm	880-304-18	880-304-22	1 piece	880-204-00	10 pieces
46	5	0.12 mm	9.40 mm	880-305-18	880-305-22	1 piece	880-205-00	10 pieces
46	6	0.12 mm	9.52 mm	880-306-18	880-306-22	1 piece	880-206-00	10 pieces
46	7	0.12 mm	9.65 mm	880-307-18	880-307-22	1 piece	880-207-00	10 pieces
46	8	0.12 mm	9.77 mm	880-308-18	880-308-22	1 piece	880-208-00	10 pieces
46	9	0.15 mm	9.89 mm	880-309-18	880-309-22	1 piece	880-209-00	10 pieces
46	10	0.15 mm	10.01 mm	880-310-18	880-310-22	1 piece	880-210-00	10 pieces
46	11	0.15 mm	10.14 mm	880-311-18	880-311-22	1 piece	880-211-00	10 pieces
46	12	0.15 mm	10.25 mm	880-312-18	880-312-22	1 piece	880-212-00	10 pieces
46	13	0.15 mm	10.38 mm	880-313-18	880-313-22	1 piece	880-213-00	10 pieces
46	14	0.15 mm	10.50 mm	880-314-18	880-314-22	1 piece	880-214-00	10 pieces
46	15	0.15 mm	10.63 mm	880-315-18	880-315-22	1 piece	880-215-00	10 pieces
46	16	0.15 mm	10.75 mm	880-316-18	880-316-22	1 piece	880-216-00	10 pieces
46	17	0.15 mm	10.87 mm	880-317-18	880-317-22	1 piece	880-217-00	10 pieces
46	18	0.15 mm	10.99 mm	880-318-18	880-318-22	1 piece	880-218-00	10 pieces
46	19	0.15 mm	11.12 mm	880-319-18	880-319-22	1 piece	880-219-00	10 pieces
46	20	0.15 mm	11.24 mm	880-320-18	880-320-22	1 piece	880-220-00	10 pieces
46	21	0.15 mm	11.36 mm	880-321-18	880-321-22	1 piece	880-221-00	10 pieces
46	22	0.15 mm	11.48 mm	880-322-18	880-322-22	1 piece	880-222-00	10 pieces
46	23	0.15 mm	11.61 mm	880-323-18	880-323-22	1 piece	880-223-00	10 pieces
46	24	0.15 mm	11.73 mm	880-324-18	880-324-22	1 piece	880-224-00	10 pieces
46	25	0.18 mm	11.85 mm	880-325-18	880-325-22	1 piece	880-225-00	10 pieces
46	26	0.18 mm	11.97 mm	880-326-18	880-326-22	1 piece	880-226-00	10 pieces
46	27	0.18 mm	12.10 mm	880-327-18	880-327-22	1 piece	880-227-00	10 pieces
46	28	0.18 mm	12.22 mm	880-328-18	880-328-22	1 piece	880-228-00	10 pieces
46	29	0.18 mm	12.34 mm	880-329-18	880-329-22	1 piece	880-229-00	10 pieces
46	30	0.18 mm	12.46 mm	880-330-18	880-330-22	1 piece	880-230-00	10 pieces
46	31	0.18 mm	12.59 mm	880-331-18	880-331-22	1 piece	880-231-00	10 pieces
46	32	0.18 mm	12.71 mm	880-332-18	880-332-22	1 piece	880-232-00	10 pieces
46	33	0.18 mm	12.83 mm	880-333-18	880-333-22	1 piece	880-233-00	10 pieces
46	34	0.18 mm	12.96 mm	880-334-18	880-334-22	1 piece	880-234-00	10 pieces
46	35	0.18 mm	13.08 mm	880-335-18	880-335-22	1 piece	880-235-00	10 pieces
46	36	0.18 mm	13.20 mm	880-336-18	880-336-22	1 piece	880-236-00	10 pieces

¹⁾ welded with buccal tubes 818-151-40 and double hook 750-602-10

²⁾ welded with buccal tubes 822-151-40 and double hook 750-602-10

* The Dentaform version of this prescription is not claimed to be a duplication of any other, nor does Dentaform imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.



dentaform® Snap / Single packs

- with super snap effect
- hard, stiff quality
- anatomically shaped
- very accurate fit
- ideal for the extra-oral technique (facebow) or SUS
- secure fit and optimal adhesion due to rough interior surface
- individual welding, e.g. acc. to McLaughlin-Bennett-Trevisi* 18/22, is possible on request

We recommend our dental vernier (REF 042-751-00) for determining the largest mesio-distal diameter of the tooth crown.

Lower first molar bands, left

Tooth	Band size	Band thickness	ø mesial-distal	welded		Quantity	unwelded	Quantity
				REF (Roth 18) ¹⁾	REF (Roth 22) ²⁾		REF	
36	1	0.12 mm	8.91 mm	881-301-18	881-301-22	1 piece	881-201-00	10 pieces
36	2	0.12 mm	9.03 mm	881-302-18	881-302-22	1 piece	881-202-00	10 pieces
36	3	0.12 mm	9.15 mm	881-303-18	881-303-22	1 piece	881-203-00	10 pieces
36	4	0.12 mm	9.28 mm	881-304-18	881-304-22	1 piece	881-204-00	10 pieces
36	5	0.12 mm	9.40 mm	881-305-18	881-305-22	1 piece	881-205-00	10 pieces
36	6	0.12 mm	9.52 mm	881-306-18	881-306-22	1 piece	881-206-00	10 pieces
36	7	0.12 mm	9.65 mm	881-307-18	881-307-22	1 piece	881-207-00	10 pieces
36	8	0.12 mm	9.77 mm	881-308-18	881-308-22	1 piece	881-208-00	10 pieces
36	9	0.15 mm	9.89 mm	881-309-18	881-309-22	1 piece	881-209-00	10 pieces
36	10	0.15 mm	10.01 mm	881-310-18	881-310-22	1 piece	881-210-00	10 pieces
36	11	0.15 mm	10.14 mm	881-311-18	881-311-22	1 piece	881-211-00	10 pieces
36	12	0.15 mm	10.25 mm	881-312-18	881-312-22	1 piece	881-212-00	10 pieces
36	13	0.15 mm	10.38 mm	881-313-18	881-313-22	1 piece	881-213-00	10 pieces
36	14	0.15 mm	10.50 mm	881-314-18	881-314-22	1 piece	881-214-00	10 pieces
36	15	0.15 mm	10.63 mm	881-315-18	881-315-22	1 piece	881-215-00	10 pieces
36	16	0.15 mm	10.75 mm	881-316-18	881-316-22	1 piece	881-216-00	10 pieces
36	17	0.15 mm	10.87 mm	881-317-18	881-317-22	1 piece	881-217-00	10 pieces
36	18	0.15 mm	10.99 mm	881-318-18	881-318-22	1 piece	881-218-00	10 pieces
36	19	0.15 mm	11.12 mm	881-319-18	881-319-22	1 piece	881-219-00	10 pieces
36	20	0.15 mm	11.24 mm	881-320-18	881-320-22	1 piece	881-220-00	10 pieces
36	21	0.15 mm	11.36 mm	881-321-18	881-321-22	1 piece	881-221-00	10 pieces
36	22	0.15 mm	11.48 mm	881-322-18	881-322-22	1 piece	881-222-00	10 pieces
36	23	0.15 mm	11.61 mm	881-323-18	881-323-22	1 piece	881-223-00	10 pieces
36	24	0.15 mm	11.73 mm	881-324-18	881-324-22	1 piece	881-224-00	10 pieces
36	25	0.18 mm	11.85 mm	881-325-18	881-325-22	1 piece	881-225-00	10 pieces
36	26	0.18 mm	11.97 mm	881-326-18	881-326-22	1 piece	881-226-00	10 pieces
36	27	0.18 mm	12.10 mm	881-327-18	881-327-22	1 piece	881-227-00	10 pieces
36	28	0.18 mm	12.22 mm	881-328-18	881-328-22	1 piece	881-228-00	10 pieces
36	29	0.18 mm	12.34 mm	881-329-18	881-329-22	1 piece	881-229-00	10 pieces
36	30	0.18 mm	12.46 mm	881-330-18	881-330-22	1 piece	881-230-00	10 pieces
36	31	0.18 mm	12.59 mm	881-331-18	881-331-22	1 piece	881-231-00	10 pieces
36	32	0.18 mm	12.71 mm	881-332-18	881-332-22	1 piece	881-232-00	10 pieces
36	33	0.18 mm	12.83 mm	881-333-18	881-333-22	1 piece	881-233-00	10 pieces
36	34	0.18 mm	12.96 mm	881-334-18	881-334-22	1 piece	881-234-00	10 pieces
36	35	0.18 mm	13.08 mm	881-335-18	881-335-22	1 piece	881-235-00	10 pieces
36	36	0.18 mm	13.20 mm	881-336-18	881-336-22	1 piece	881-236-00	10 pieces



¹⁾ welded with buccal tubes 818-151-30 and double hook 750-602-10

²⁾ welded with buccal tubes 822-151-30 and double hook 750-602-10

* The Dentaform version of this prescription is not claimed to be a duplication of any other, nor does Dentaform imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.

- easy adaptation
- slightly softer band quality (compared with dentaform® Snap)
- anatomically shaped
- very accurate fit
- unwelded bands (with or without double hook) ideal for the Herbst technique (welding service for double bands available free of charge on request)
- secure fit and optimal adhesion due to rough interior surface
- delivery in a tray (fits in band cabinet REF 890-414-00)
- individual welding, e.g. acc. to McLaughlin-Bennett-Trevisi* 18/22, is possible on request

We recommend our dental vernier (REF 042-751-00) for determining the largest mesio-distal diameter of the tooth crown.

Upper first molar bands

Mini

Tooth	Band sizes (Contents per size)	welded		unwelded	Quantity
		REF (Roth 18) ¹⁾	REF (Roth 22) ²⁾	REF	
16	No. 5-36 (2)	882-364-18	882-364-22	882-064-00	64 pieces
26	No. 5-36 (2)	883-364-18	883-364-22	883-064-00	64 pieces

Medium

Tooth	Band sizes (Contents per size)	welded		unwelded	Quantity
		REF (Roth 18) ¹⁾	REF (Roth 22) ²⁾	REF	
16	No. 5-8, 30-36 (2), 9-18 (4), 19-20 (6), 21-25, 28-29 (8), 26-27 (10)	882-350-18	882-350-22	882-150-00	150 pieces
26	No. 5-8, 30-36 (2), 9-18 (4), 19-20 (6), 21-25, 28-29 (8), 26-27 (10)	883-350-18	883-350-22	883-150-00	150 pieces

Maxi

Tooth	Band sizes (Contents per size)	welded		unwelded	Quantity
		REF (Roth 18)	REF (Roth 22)	REF	
16	No. 5-8, 30-36 (4), 9-18 (8), 19-20 (12), 21-25, 28-29 (16), 26-27 (20)			882-300-00	300 pieces
26	No. 5-8, 30-36 (4), 9-18 (8), 19-20 (12), 21-25, 28-29 (16), 26-27 (20)			883-300-00	300 pieces

¹⁾(16) welded with buccal tubes 818-131-10 and double hook 750-602-10

(26) welded with buccal tubes 818-131-20 and double hook 750-602-10

²⁾(16) welded with buccal tubes 822-131-10 and double hook 750-602-10

(26) welded with buccal tubes 822-131-20 and double hook 750-602-10

* The Dentaform version of this prescription is not claimed to be a duplication of any other, nor does Dentaform imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.

Upper second molar bands

Mini

Tooth	Band sizes (Contents per size)	welded		unwelded	Quantity
		REF (Roth 18)	REF (Roth 22)	REF	
17	No. 1-32 (2)			886-064-00	64 pieces
27	No. 1-32 (2)			887-064-00	64 pieces

Medium

Tooth	Band sizes (Contents per size)	welded		unwelded	Quantity
		REF (Roth 18)	REF (Roth 22)	REF	
17	No. 1-4, 28-32 (2), 5-10, 21-27 (4), 11-12, 19-20 (6), 13, 18 (8), 14-17 (10)			886-150-00	150 pieces
27	No. 1-4, 28-32 (2), 5-10, 21-27 (4), 11-12, 19-20 (6), 13, 18 (8), 14-17 (10)			887-150-00	150 pieces

* The Dentaform version of this prescription is not claimed to be a duplication of any other, nor does Dentaform imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.

dentaform® / Assortments

- easy adaptation
- slightly softer band quality (compared with dentaform® Snap)
- anatomically shaped
- very accurate fit
- unwelded bands (with or without double hook) ideal for the Herbst technique (welding service for double bands available free of charge on request)
- secure fit and optimal adhesion due to rough interior surface
- delivery in a tray (fits in band cabinet REF 890-414-00)
- individual welding, e.g. acc. to McLaughlin-Bennett-Trevisi* 18/22, is possible on request

We recommend our dental vernier (REF 042-751-00) for determining the largest mesio-distal diameter of the tooth crown.

Lower first molar bands

Mini

Tooth	Band sizes (Contents per size)	welded		unwelded	Quantity
		REF (Roth 18) ¹⁾	REF (Roth 22) ²⁾	REF	
46	No. 5-36 (2)	884-364-18	884-364-22	884-064-00	64 pieces
36	No. 5-36 (2)	885-364-18	885-364-22	885-064-00	64 pieces

Medium

Tooth	Band sizes (Contents per size)	welded		unwelded	Quantity
		REF (Roth 18) ¹⁾	REF (Roth 22) ²⁾	REF	
46	No. 5-8, 30-36 (2), 9-18 (4), 19-20 (6), 21-25, 28-29 (8), 26-27 (10)	884-350-18	884-350-22	884-150-00	150 pieces
36	No. 5-8, 30-36 (2), 9-18 (4), 19-20 (6), 21-25, 28-29 (8), 26-27 (10)	885-350-18	885-350-22	885-150-00	150 pieces

Maxi

Tooth	Band sizes (Contents per size)	welded		unwelded	Quantity
		REF (Roth 18)	REF (Roth 22)	REF	
46	No. 5-8, 30-36 (4), 9-18 (8), 19-20 (12), 21-25, 28-29 (16), 26-27 (20)			884-300-00	300 pieces
36	No. 5-8, 30-36 (4), 9-18 (8), 19-20 (12), 21-25, 28-29 (16), 26-27 (20)			885-300-00	300 pieces

¹⁾(36) welded with buccal tubes 818-151-30 and double hook 750-602-10

(46) welded with buccal tubes 818-151-40 and double hook 750-602-10

²⁾(36) welded with buccal tubes 822-151-30 and double hook 750-602-10

(46) welded with buccal tubes 822-151-40 and double hook 750-602-10

* The Dentaform version of this prescription is not claimed to be a duplication of any other, nor does Dentaform imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.

Lower second molar bands

Mini

Tooth	Band sizes (Contents per size)	welded		unwelded	Quantity
		REF (Roth 18)	REF (Roth 22)	REF	
47	No. 1-32 (2)			888-064-00	64 pieces
37	No. 1-32 (2)			889-064-00	64 pieces

Medium

Tooth	Band sizes (Contents per size)	welded		unwelded	Quantity
		REF (Roth 18)	REF (Roth 22)	REF	
47	No. 1-4, 28-32 (2), 5-10, 21-27 (4), 11-12, 19-20 (6), 13, 18 (8), 14-17 (10)			888-150-00	150 pieces
37	No. 1-4, 28-32 (2), 5-10, 21-27 (4), 11-12, 19-20 (6), 13, 18 (8), 14-17 (10)			889-150-00	150 pieces

* The Dentaform version of this prescription is not claimed to be a duplication of any other, nor does Dentaform imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.

- easy adaptation
- slightly softer band quality (compared with dentaform® Snap)
- anatomically shaped
- very accurate fit
- unwelded bands (with or without double hook) ideal for the Herbst technique (welding service for double bands available free of charge on request)
- secure fit and optimal adhesion due to rough interior surface
- individual welding, e.g. acc. to McLaughlin-Bennett-Trevisi* 18/22, is possible on request

We recommend our dental vernier (REF 042-751-00) for determining the largest mesio-distal diameter of the tooth crown.

Upper first molar bands, right

Tooth	Band size	Band thickness	ø mesial-distal	welded		Quantity	unwelded	Quantity
				REF (Roth 18) ¹⁾	REF (Roth 22) ²⁾		REF	
16	1	0.12 mm	8.78 mm	882-301-18	882-301-22	1 piece	882-001-00	10 pieces
16	2	0.12 mm	8.89 mm	882-302-18	882-302-22	1 piece	882-002-00	10 pieces
16	3	0.12 mm	8.99 mm	882-303-18	882-303-22	1 piece	882-003-00	10 pieces
16	4	0.12 mm	9.09 mm	882-304-18	882-304-22	1 piece	882-004-00	10 pieces
16	5	0.12 mm	9.20 mm	882-305-18	882-305-22	1 piece	882-005-00	10 pieces
16	6	0.12 mm	9.30 mm	882-306-18	882-306-22	1 piece	882-006-00	10 pieces
16	7	0.12 mm	9.40 mm	882-307-18	882-307-22	1 piece	882-007-00	10 pieces
16	8	0.12 mm	9.50 mm	882-308-18	882-308-22	1 piece	882-008-00	10 pieces
16	9	0.15 mm	9.61 mm	882-309-18	882-309-22	1 piece	882-009-00	10 pieces
16	10	0.15 mm	9.71 mm	882-310-18	882-310-22	1 piece	882-010-00	10 pieces
16	11	0.15 mm	9.81 mm	882-311-18	882-311-22	1 piece	882-011-00	10 pieces
16	12	0.15 mm	9.92 mm	882-312-18	882-312-22	1 piece	882-012-00	10 pieces
16	13	0.15 mm	10.02 mm	882-313-18	882-313-22	1 piece	882-013-00	10 pieces
16	14	0.15 mm	10.12 mm	882-314-18	882-314-22	1 piece	882-014-00	10 pieces
16	15	0.15 mm	10.23 mm	882-315-18	882-315-22	1 piece	882-015-00	10 pieces
16	16	0.15 mm	10.33 mm	882-316-18	882-316-22	1 piece	882-016-00	10 pieces
16	17	0.15 mm	10.43 mm	882-317-18	882-317-22	1 piece	882-017-00	10 pieces
16	18	0.15 mm	10.54 mm	882-318-18	882-318-22	1 piece	882-018-00	10 pieces
16	19	0.15 mm	10.64 mm	882-319-18	882-319-22	1 piece	882-019-00	10 pieces
16	20	0.15 mm	10.74 mm	882-320-18	882-320-22	1 piece	882-020-00	10 pieces
16	21	0.15 mm	10.85 mm	882-321-18	882-321-22	1 piece	882-021-00	10 pieces
16	22	0.15 mm	10.95 mm	882-322-18	882-322-22	1 piece	882-022-00	10 pieces
16	23	0.15 mm	11.05 mm	882-323-18	882-323-22	1 piece	882-023-00	10 pieces
16	24	0.15 mm	11.16 mm	882-324-18	882-324-22	1 piece	882-024-00	10 pieces
16	25	0.18 mm	11.26 mm	882-325-18	882-325-22	1 piece	882-025-00	10 pieces
16	26	0.18 mm	11.36 mm	882-326-18	882-326-22	1 piece	882-026-00	10 pieces
16	27	0.18 mm	11.47 mm	882-327-18	882-327-22	1 piece	882-027-00	10 pieces
16	28	0.18 mm	11.57 mm	882-328-18	882-328-22	1 piece	882-028-00	10 pieces
16	29	0.18 mm	11.67 mm	882-329-18	882-329-22	1 piece	882-029-00	10 pieces
16	30	0.18 mm	11.78 mm	882-330-18	882-330-22	1 piece	882-030-00	10 pieces
16	31	0.18 mm	11.88 mm	882-331-18	882-331-22	1 piece	882-031-00	10 pieces
16	32	0.18 mm	11.98 mm	882-332-18	882-332-22	1 piece	882-032-00	10 pieces
16	33	0.18 mm	12.08 mm	882-333-18	882-333-22	1 piece	882-033-00	10 pieces
16	34	0.18 mm	12.19 mm	882-334-18	882-334-22	1 piece	882-034-00	10 pieces
16	35	0.18 mm	12.29 mm	882-335-18	882-335-22	1 piece	882-035-00	10 pieces
16	36	0.18 mm	12.39 mm	882-336-18	882-336-22	1 piece	882-036-00	10 pieces

¹⁾ welded with buccal tubes 818-131-10 and double hook 750-602-10

²⁾ welded with buccal tubes 822-131-10 and double hook 750-602-10



* The Dentaform version of this prescription is not claimed to be a duplication of any other, nor does Dentaform imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.

dentaform® / Single packs

- easy adaptation
- slightly softer band quality (compared with dentaform® Snap)
- anatomically shaped
- very accurate fit
- unwelded bands (with or without double hook) ideal for the Herbst technique (welding service for double bands available free of charge on request)
- secure fit and optimal adhesion due to rough interior surface
- individual welding, e.g. acc. to McLaughlin-Bennett-Trevisi* 18/22, is possible on request

We recommend our dental vernier (REF 042-751-00) for determining the largest mesio-distal diameter of the tooth crown.

Upper first molar bands, left

Tooth	Band size	Band thickness	ø mesial-distal	welded		Quantity	unwelded	Quantity
				REF (Roth 18) ¹⁾	REF (Roth 22) ²⁾		REF	
26	1	0.12 mm	8.78 mm	883-301-18	883-301-22	1 piece	883-001-00	10 pieces
26	2	0.12 mm	8.89 mm	883-302-18	883-302-22	1 piece	883-002-00	10 pieces
26	3	0.12 mm	8.99 mm	883-303-18	883-303-22	1 piece	883-003-00	10 pieces
26	4	0.12 mm	9.09 mm	883-304-18	883-304-22	1 piece	883-004-00	10 pieces
26	5	0.12 mm	9.20 mm	883-305-18	883-305-22	1 piece	883-005-00	10 pieces
26	6	0.12 mm	9.30 mm	883-306-18	883-306-22	1 piece	883-006-00	10 pieces
26	7	0.12 mm	9.40 mm	883-307-18	883-307-22	1 piece	883-007-00	10 pieces
26	8	0.12 mm	9.50 mm	883-308-18	883-308-22	1 piece	883-008-00	10 pieces
26	9	0.15 mm	9.61 mm	883-309-18	883-309-22	1 piece	883-009-00	10 pieces
26	10	0.15 mm	9.71 mm	883-310-18	883-310-22	1 piece	883-010-00	10 pieces
26	11	0.15 mm	9.81 mm	883-311-18	883-311-22	1 piece	883-011-00	10 pieces
26	12	0.15 mm	9.92 mm	883-312-18	883-312-22	1 piece	883-012-00	10 pieces
26	13	0.15 mm	10.02 mm	883-313-18	883-313-22	1 piece	883-013-00	10 pieces
26	14	0.15 mm	10.12 mm	883-314-18	883-314-22	1 piece	883-014-00	10 pieces
26	15	0.15 mm	10.23 mm	883-315-18	883-315-22	1 piece	883-015-00	10 pieces
26	16	0.15 mm	10.33 mm	883-316-18	883-316-22	1 piece	883-016-00	10 pieces
26	17	0.15 mm	10.43 mm	883-317-18	883-317-22	1 piece	883-017-00	10 pieces
26	18	0.15 mm	10.54 mm	883-318-18	883-318-22	1 piece	883-018-00	10 pieces
26	19	0.15 mm	10.64 mm	883-319-18	883-319-22	1 piece	883-019-00	10 pieces
26	20	0.15 mm	10.74 mm	883-320-18	883-320-22	1 piece	883-020-00	10 pieces
26	21	0.15 mm	10.85 mm	883-321-18	883-321-22	1 piece	883-021-00	10 pieces
26	22	0.15 mm	10.95 mm	883-322-18	883-322-22	1 piece	883-022-00	10 pieces
26	23	0.15 mm	11.05 mm	883-323-18	883-323-22	1 piece	883-023-00	10 pieces
26	24	0.15 mm	11.16 mm	883-324-18	883-324-22	1 piece	883-024-00	10 pieces
26	25	0.18 mm	11.26 mm	883-325-18	883-325-22	1 piece	883-025-00	10 pieces
26	26	0.18 mm	11.36 mm	883-326-18	883-326-22	1 piece	883-026-00	10 pieces
26	27	0.18 mm	11.47 mm	883-327-18	883-327-22	1 piece	883-027-00	10 pieces
26	28	0.18 mm	11.57 mm	883-328-18	883-328-22	1 piece	883-028-00	10 pieces
26	29	0.18 mm	11.67 mm	883-329-18	883-329-22	1 piece	883-029-00	10 pieces
26	30	0.18 mm	11.78 mm	883-330-18	883-330-22	1 piece	883-030-00	10 pieces
26	31	0.18 mm	11.88 mm	883-331-18	883-331-22	1 piece	883-031-00	10 pieces
26	32	0.18 mm	11.98 mm	883-332-18	883-332-22	1 piece	883-032-00	10 pieces
26	33	0.18 mm	12.08 mm	883-333-18	883-333-22	1 piece	883-033-00	10 pieces
26	34	0.18 mm	12.19 mm	883-334-18	883-334-22	1 piece	883-034-00	10 pieces
26	35	0.18 mm	12.29 mm	883-335-18	883-335-22	1 piece	883-035-00	10 pieces
26	36	0.18 mm	12.39 mm	883-336-18	883-336-22	1 piece	883-036-00	10 pieces

¹⁾ welded with buccal tubes 818-131-20 and double hook 750-602-10

²⁾ welded with buccal tubes 822-131-20 and double hook 750-602-10



* The Dentaform version of this prescription is not claimed to be a duplication of any other, nor does Dentaform imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.

- easy adaptation
- slightly softer band quality (compared with dentaform® Snap)
- anatomically shaped
- very accurate fit
- unwelded bands (with or without double hook) ideal for the Herbst technique (welding service for double bands available free of charge on request)
- secure fit and optimal adhesion due to rough interior surface
- individual welding, e.g. acc. to McLaughlin-Bennett-Trevisi* 18/22, is possible on request

We recommend our dental vernier (REF 042-751-00) for determining the largest mesio-distal diameter of the tooth crown.

Lower first molar bands, right

Tooth	Band size	Band thickness	ø mesial-distal	welded		Quantity	unwelded	Quantity
				REF (Roth 18) ¹⁾	REF (Roth 22) ²⁾		REF	
46	1	0.12 mm	8.91 mm	884-301-18	884-301-22	1 piece	884-001-00	10 pieces
46	2	0.12 mm	9.03 mm	884-302-18	884-302-22	1 piece	884-002-00	10 pieces
46	3	0.12 mm	9.15 mm	884-303-18	884-303-22	1 piece	884-003-00	10 pieces
46	4	0.12 mm	9.28 mm	884-304-18	884-304-22	1 piece	884-004-00	10 pieces
46	5	0.12 mm	9.40 mm	884-305-18	884-305-22	1 piece	884-005-00	10 pieces
46	6	0.12 mm	9.52 mm	884-306-18	884-306-22	1 piece	884-006-00	10 pieces
46	7	0.12 mm	9.65 mm	884-307-18	884-307-22	1 piece	884-007-00	10 pieces
46	8	0.12 mm	9.77 mm	884-308-18	884-308-22	1 piece	884-008-00	10 pieces
46	9	0.15 mm	9.89 mm	884-309-18	884-309-22	1 piece	884-009-00	10 pieces
46	10	0.15 mm	10.01 mm	884-310-18	884-310-22	1 piece	884-010-00	10 pieces
46	11	0.15 mm	10.14 mm	884-311-18	884-311-22	1 piece	884-011-00	10 pieces
46	12	0.15 mm	10.25 mm	884-312-18	884-312-22	1 piece	884-012-00	10 pieces
46	13	0.15 mm	10.38 mm	884-313-18	884-313-22	1 piece	884-013-00	10 pieces
46	14	0.15 mm	10.50 mm	884-314-18	884-314-22	1 piece	884-014-00	10 pieces
46	15	0.15 mm	10.63 mm	884-315-18	884-315-22	1 piece	884-015-00	10 pieces
46	16	0.15 mm	10.75 mm	884-316-18	884-316-22	1 piece	884-016-00	10 pieces
46	17	0.15 mm	10.87 mm	884-317-18	884-317-22	1 piece	884-017-00	10 pieces
46	18	0.15 mm	10.99 mm	884-318-18	884-318-22	1 piece	884-018-00	10 pieces
46	19	0.15 mm	11.12 mm	884-319-18	884-319-22	1 piece	884-019-00	10 pieces
46	20	0.15 mm	11.24 mm	884-320-18	884-320-22	1 piece	884-020-00	10 pieces
46	21	0.15 mm	11.36 mm	884-321-18	884-321-22	1 piece	884-021-00	10 pieces
46	22	0.15 mm	11.48 mm	884-322-18	884-322-22	1 piece	884-022-00	10 pieces
46	23	0.15 mm	11.61 mm	884-323-18	884-323-22	1 piece	884-023-00	10 pieces
46	24	0.15 mm	11.73 mm	884-324-18	884-324-22	1 piece	884-024-00	10 pieces
46	25	0.18 mm	11.85 mm	884-325-18	884-325-22	1 piece	884-025-00	10 pieces
46	26	0.18 mm	11.97 mm	884-326-18	884-326-22	1 piece	884-026-00	10 pieces
46	27	0.18 mm	12.10 mm	884-327-18	884-327-22	1 piece	884-027-00	10 pieces
46	28	0.18 mm	12.22 mm	884-328-18	884-328-22	1 piece	884-028-00	10 pieces
46	29	0.18 mm	12.34 mm	884-329-18	884-329-22	1 piece	884-029-00	10 pieces
46	30	0.18 mm	12.46 mm	884-330-18	884-330-22	1 piece	884-030-00	10 pieces
46	31	0.18 mm	12.59 mm	884-331-18	884-331-22	1 piece	884-031-00	10 pieces
46	32	0.18 mm	12.71 mm	884-332-18	884-332-22	1 piece	884-032-00	10 pieces
46	33	0.18 mm	12.83 mm	884-333-18	884-333-22	1 piece	884-033-00	10 pieces
46	34	0.18 mm	12.96 mm	884-334-18	884-334-22	1 piece	884-034-00	10 pieces
46	35	0.18 mm	13.08 mm	884-335-18	884-335-22	1 piece	884-035-00	10 pieces
46	36	0.18 mm	13.20 mm	884-336-18	884-336-22	1 piece	884-036-00	10 pieces

¹⁾ welded with buccal tubes 818-151-40 and double hook 750-602-10

²⁾ welded with buccal tubes 822-151-40 and double hook 750-602-10

* The Dentaform version of this prescription is not claimed to be a duplication of any other, nor does Dentaform imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.



dentaform® / Single packs

- easy adaptation
- slightly softer band quality (compared with dentaform® Snap)
- anatomically shaped
- very accurate fit
- unwelded bands (with or without double hook) ideal for the Herbst technique (welding service for double bands available free of charge on request)
- secure fit and optimal adhesion due to rough interior surface
- individual welding, e.g. acc. to McLaughlin-Bennett-Trevisi* 18/22, is possible on request

We recommend our dental vernier (REF 042-751-00) for determining the largest mesio-distal diameter of the tooth crown.

Lower first molar bands, left

Tooth	Band size	Band thickness	ø mesial-distal	welded		Quantity	unwelded	Quantity
				REF (Roth 18) ¹⁾	REF (Roth 22) ²⁾		REF	
36	1	0.12 mm	8.91 mm	885-301-18	885-301-22	1 piece	885-001-00	10 pieces
36	2	0.12 mm	9.03 mm	885-302-18	885-302-22	1 piece	885-002-00	10 pieces
36	3	0.12 mm	9.15 mm	885-303-18	885-303-22	1 piece	885-003-00	10 pieces
36	4	0.12 mm	9.28 mm	885-304-18	885-304-22	1 piece	885-004-00	10 pieces
36	5	0.12 mm	9.40 mm	885-305-18	885-305-22	1 piece	885-005-00	10 pieces
36	6	0.12 mm	9.52 mm	885-306-18	885-306-22	1 piece	885-006-00	10 pieces
36	7	0.12 mm	9.65 mm	885-307-18	885-307-22	1 piece	885-007-00	10 pieces
36	8	0.12 mm	9.77 mm	885-308-18	885-308-22	1 piece	885-008-00	10 pieces
36	9	0.15 mm	9.89 mm	885-309-18	885-309-22	1 piece	885-009-00	10 pieces
36	10	0.15 mm	10.01 mm	885-310-18	885-310-22	1 piece	885-010-00	10 pieces
36	11	0.15 mm	10.14 mm	885-311-18	885-311-22	1 piece	885-011-00	10 pieces
36	12	0.15 mm	10.25 mm	885-312-18	885-312-22	1 piece	885-012-00	10 pieces
36	13	0.15 mm	10.38 mm	885-313-18	885-313-22	1 piece	885-013-00	10 pieces
36	14	0.15 mm	10.50 mm	885-314-18	885-314-22	1 piece	885-014-00	10 pieces
36	15	0.15 mm	10.63 mm	885-315-18	885-315-22	1 piece	885-015-00	10 pieces
36	16	0.15 mm	10.75 mm	885-316-18	885-316-22	1 piece	885-016-00	10 pieces
36	17	0.15 mm	10.87 mm	885-317-18	885-317-22	1 piece	885-017-00	10 pieces
36	18	0.15 mm	10.99 mm	885-318-18	885-318-22	1 piece	885-018-00	10 pieces
36	19	0.15 mm	11.12 mm	885-319-18	885-319-22	1 piece	885-019-00	10 pieces
36	20	0.15 mm	11.24 mm	885-320-18	885-320-22	1 piece	885-020-00	10 pieces
36	21	0.15 mm	11.36 mm	885-321-18	885-321-22	1 piece	885-021-00	10 pieces
36	22	0.15 mm	11.48 mm	885-322-18	885-322-22	1 piece	885-022-00	10 pieces
36	23	0.15 mm	11.61 mm	885-323-18	885-323-22	1 piece	885-023-00	10 pieces
36	24	0.15 mm	11.73 mm	885-324-18	885-324-22	1 piece	885-024-00	10 pieces
36	25	0.18 mm	11.85 mm	885-325-18	885-325-22	1 piece	885-025-00	10 pieces
36	26	0.18 mm	11.97 mm	885-326-18	885-326-22	1 piece	885-026-00	10 pieces
36	27	0.18 mm	12.10 mm	885-327-18	885-327-22	1 piece	885-027-00	10 pieces
36	28	0.18 mm	12.22 mm	885-328-18	885-328-22	1 piece	885-028-00	10 pieces
36	29	0.18 mm	12.34 mm	885-329-18	885-329-22	1 piece	885-029-00	10 pieces
36	30	0.18 mm	12.46 mm	885-330-18	885-330-22	1 piece	885-030-00	10 pieces
36	31	0.18 mm	12.59 mm	885-331-18	885-331-22	1 piece	885-031-00	10 pieces
36	32	0.18 mm	12.71 mm	885-332-18	885-332-22	1 piece	885-032-00	10 pieces
36	33	0.18 mm	12.83 mm	885-333-18	885-333-22	1 piece	885-033-00	10 pieces
36	34	0.18 mm	12.96 mm	885-334-18	885-334-22	1 piece	885-034-00	10 pieces
36	35	0.18 mm	13.08 mm	885-335-18	885-335-22	1 piece	885-035-00	10 pieces
36	36	0.18 mm	13.20 mm	885-336-18	885-336-22	1 piece	885-036-00	10 pieces

¹⁾ welded with buccal tubes 818-151-30 and double hook 750-602-10

²⁾ welded with buccal tubes 822-151-30 and double hook 750-602-10

* The Dentaform version of this prescription is not claimed to be a duplication of any other, nor does Dentaform imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.



- easy adaptation
- slightly softer band quality (compared with dentaform® Snap)
- anatomically shaped
- very accurate fit
- secure fit and optimal adhesion due to rough interior surface
- individual welding, e.g. acc. to McLaughlin-Bennett-Trevisi* 18/22, is possible on request

We recommend our dental vernier (REF 042-751-00) for determining the largest mesio-distal diameter of the tooth crown.

Upper second molar bands, right

Tooth	Band size	Band thickness	Ø mesial-distal	welded		Quantity	unwelded	Quantity
				REF (Roth 18)	REF (Roth 22)		REF	
17	1	0.12 mm	8.10 mm				886-001-00	10 pieces
17	2	0.12 mm	8.22 mm				886-002-00	10 pieces
17	3	0.12 mm	8.33 mm				886-003-00	10 pieces
17	4	0.12 mm	8.45 mm				886-004-00	10 pieces
17	5	0.12 mm	8.56 mm				886-005-00	10 pieces
17	6	0.12 mm	8.68 mm				886-006-00	10 pieces
17	7	0.15 mm	8.80 mm				886-007-00	10 pieces
17	8	0.15 mm	8.91 mm				886-008-00	10 pieces
17	9	0.15 mm	9.03 mm				886-009-00	10 pieces
17	10	0.15 mm	9.14 mm				886-010-00	10 pieces
17	11	0.15 mm	9.26 mm				886-011-00	10 pieces
17	12	0.15 mm	9.38 mm				886-012-00	10 pieces
17	13	0.15 mm	9.49 mm				886-013-00	10 pieces
17	14	0.15 mm	9.61 mm				886-014-00	10 pieces
17	15	0.15 mm	9.72 mm				886-015-00	10 pieces
17	16	0.15 mm	9.84 mm				886-016-00	10 pieces
17	17	0.15 mm	9.96 mm				886-017-00	10 pieces
17	18	0.15 mm	10.07 mm				886-018-00	10 pieces
17	19	0.15 mm	10.19 mm				886-019-00	10 pieces
17	20	0.15 mm	10.30 mm				886-020-00	10 pieces
17	21	0.15 mm	10.42 mm				886-021-00	10 pieces
17	22	0.15 mm	10.54 mm				886-022-00	10 pieces
17	23	0.18 mm	10.65 mm				886-023-00	10 pieces
17	24	0.18 mm	10.77 mm				886-024-00	10 pieces
17	25	0.18 mm	10.88 mm				886-025-00	10 pieces
17	26	0.18 mm	11.00 mm				886-026-00	10 pieces
17	27	0.18 mm	11.12 mm				886-027-00	10 pieces
17	28	0.18 mm	11.23 mm				886-028-00	10 pieces
17	29	0.18 mm	11.35 mm				886-029-00	10 pieces
17	30	0.18 mm	11.46 mm				886-030-00	10 pieces
17	31	0.18 mm	11.58 mm				886-031-00	10 pieces
17	32	0.18 mm	11.70 mm				886-032-00	10 pieces

* The Dentaform version of this prescription is not claimed to be a duplication of any other, nor does Dentaform imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.



dentaform® / Single packs

- easy adaptation
- slightly softer band quality (compared with dentaform® Snap)
- anatomically shaped
- very accurate fit
- secure fit and optimal adhesion due to rough interior surface
- individual welding, e.g. acc. to McLaughlin-Bennett-Trevisi* 18/22, is possible on request

We recommend our dental vernier (REF 042-751-00) for determining the largest mesio-distal diameter of the tooth crown.

Upper second molar bands, left

Tooth	Band size	Band thickness	Ø mesial-distal	welded		Quantity	unwelded	Quantity
				REF (Roth 18)	REF (Roth 22)		REF	
27	1	0.12 mm	8.10 mm				887-001-00	10 pieces
27	2	0.12 mm	8.22 mm				887-002-00	10 pieces
27	3	0.12 mm	8.33 mm				887-003-00	10 pieces
27	4	0.12 mm	8.45 mm				887-004-00	10 pieces
27	5	0.12 mm	8.56 mm				887-005-00	10 pieces
27	6	0.12 mm	8.68 mm				887-006-00	10 pieces
27	7	0.15 mm	8.80 mm				887-007-00	10 pieces
27	8	0.15 mm	8.91 mm				887-008-00	10 pieces
27	9	0.15 mm	9.03 mm				887-009-00	10 pieces
27	10	0.15 mm	9.14 mm				887-010-00	10 pieces
27	11	0.15 mm	9.26 mm				887-011-00	10 pieces
27	12	0.15 mm	9.38 mm				887-012-00	10 pieces
27	13	0.15 mm	9.49 mm				887-013-00	10 pieces
27	14	0.15 mm	9.61 mm				887-014-00	10 pieces
27	15	0.15 mm	9.72 mm				887-015-00	10 pieces
27	16	0.15 mm	9.84 mm				887-016-00	10 pieces
27	17	0.15 mm	9.96 mm				887-017-00	10 pieces
27	18	0.15 mm	10.07 mm				887-018-00	10 pieces
27	19	0.15 mm	10.19 mm				887-019-00	10 pieces
27	20	0.15 mm	10.30 mm				887-020-00	10 pieces
27	21	0.15 mm	10.42 mm				887-021-00	10 pieces
27	22	0.15 mm	10.54 mm				887-022-00	10 pieces
27	23	0.18 mm	10.65 mm				887-023-00	10 pieces
27	24	0.18 mm	10.77 mm				887-024-00	10 pieces
27	25	0.18 mm	10.88 mm				887-025-00	10 pieces
27	26	0.18 mm	11.00 mm				887-026-00	10 pieces
27	27	0.18 mm	11.12 mm				887-027-00	10 pieces
27	28	0.18 mm	11.23 mm				887-028-00	10 pieces
27	29	0.18 mm	11.35 mm				887-029-00	10 pieces
27	30	0.18 mm	11.46 mm				887-030-00	10 pieces
27	31	0.18 mm	11.58 mm				887-031-00	10 pieces
27	32	0.18 mm	11.70 mm				887-032-00	10 pieces

* The Dentaform version of this prescription is not claimed to be a duplication of any other, nor does Dentaform imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.



- easy adaptation
- slightly softer band quality (compared with dentaform® Snap)
- anatomically shaped
- very accurate fit
- secure fit and optimal adhesion due to rough interior surface
- individual welding, e.g. acc. to McLaughlin-Bennett-Trevisi* 18/22, is possible on request

We recommend our dental vernier (REF 042-751-00) for determining the largest mesio-distal diameter of the tooth crown.

Lower second molar bands, right

Tooth	Band size	Band thickness	Ø mesial-distal	welded		Quantity	unwelded	Quantity
				REF (Roth 18)	REF (Roth 22)		REF	
47	1	0.12 mm	9.88 mm				888-001-00	10 pieces
47	2	0.12 mm	9.92 mm				888-002-00	10 pieces
47	3	0.12 mm	10.04 mm				888-003-00	10 pieces
47	4	0.12 mm	10.17 mm				888-004-00	10 pieces
47	5	0.12 mm	10.29 mm				888-005-00	10 pieces
47	6	0.12 mm	10.41 mm				888-006-00	10 pieces
47	7	0.15 mm	10.53 mm				888-007-00	10 pieces
47	8	0.15 mm	10.65 mm				888-008-00	10 pieces
47	9	0.15 mm	10.78 mm				888-009-00	10 pieces
47	10	0.15 mm	10.90 mm				888-010-00	10 pieces
47	11	0.15 mm	11.02 mm				888-011-00	10 pieces
47	12	0.15 mm	11.14 mm				888-012-00	10 pieces
47	13	0.15 mm	11.26 mm				888-013-00	10 pieces
47	14	0.15 mm	11.39 mm				888-014-00	10 pieces
47	15	0.15 mm	11.51 mm				888-015-00	10 pieces
47	16	0.15 mm	11.63 mm				888-016-00	10 pieces
47	17	0.15 mm	11.75 mm				888-017-00	10 pieces
47	18	0.15 mm	11.87 mm				888-018-00	10 pieces
47	19	0.15 mm	11.99 mm				888-019-00	10 pieces
47	20	0.15 mm	12.12 mm				888-020-00	10 pieces
47	21	0.15 mm	12.24 mm				888-021-00	10 pieces
47	22	0.15 mm	12.36 mm				888-022-00	10 pieces
47	23	0.18 mm	12.48 mm				888-023-00	10 pieces
47	24	0.18 mm	12.61 mm				888-024-00	10 pieces
47	25	0.18 mm	12.73 mm				888-025-00	10 pieces
47	26	0.18 mm	12.85 mm				888-026-00	10 pieces
47	27	0.18 mm	12.97 mm				888-027-00	10 pieces
47	28	0.18 mm	13.09 mm				888-028-00	10 pieces
47	29	0.18 mm	13.22 mm				888-029-00	10 pieces
47	30	0.18 mm	13.34 mm				888-030-00	10 pieces
47	31	0.18 mm	13.46 mm				888-031-00	10 pieces
47	32	0.18 mm	13.59 mm				888-032-00	10 pieces

* The Dentaform version of this prescription is not claimed to be a duplication of any other, nor does Dentaform imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.



dentaform® / Single packs

- easy adaptation
- slightly softer band quality (compared with dentaform® Snap)
- anatomically shaped
- very accurate fit
- secure fit and optimal adhesion due to rough interior surface
- individual welding, e.g. acc. to McLaughlin-Bennett-Trevisi* 18/22, is possible on request

We recommend our dental vernier (REF 042-751-00) for determining the largest mesio-distal diameter of the tooth crown.

Lower second molar bands, left

Tooth	Band size	Band thickness	Ø mesial-distal	welded		Quantity	unwelded	Quantity
				REF (Roth 18)	REF (Roth 22)		REF	
37	1	0.12 mm	9.88 mm				889-001-00	10 pieces
37	2	0.12 mm	9.92 mm				889-002-00	10 pieces
37	3	0.12 mm	10.04 mm				889-003-00	10 pieces
37	4	0.12 mm	10.17 mm				889-004-00	10 pieces
37	5	0.12 mm	10.29 mm				889-005-00	10 pieces
37	6	0.12 mm	10.41 mm				889-006-00	10 pieces
37	7	0.15 mm	10.53 mm				889-007-00	10 pieces
37	8	0.15 mm	10.65 mm				889-008-00	10 pieces
37	9	0.15 mm	10.78 mm				889-009-00	10 pieces
37	10	0.15 mm	10.90 mm				889-010-00	10 pieces
37	11	0.15 mm	11.02 mm				889-011-00	10 pieces
37	12	0.15 mm	11.14 mm				889-012-00	10 pieces
37	13	0.15 mm	11.26 mm				889-013-00	10 pieces
37	14	0.15 mm	11.39 mm				889-014-00	10 pieces
37	15	0.15 mm	11.51 mm				889-015-00	10 pieces
37	16	0.15 mm	11.63 mm				889-016-00	10 pieces
37	17	0.15 mm	11.75 mm				889-017-00	10 pieces
37	18	0.15 mm	11.87 mm				889-018-00	10 pieces
37	19	0.15 mm	11.99 mm				889-019-00	10 pieces
37	20	0.15 mm	12.12 mm				889-020-00	10 pieces
37	21	0.15 mm	12.24 mm				889-021-00	10 pieces
37	22	0.15 mm	12.36 mm				889-022-00	10 pieces
37	23	0.18 mm	12.48 mm				889-023-00	10 pieces
37	24	0.18 mm	12.61 mm				889-024-00	10 pieces
37	25	0.18 mm	12.73 mm				889-025-00	10 pieces
37	26	0.18 mm	12.85 mm				889-026-00	10 pieces
37	27	0.18 mm	12.97 mm				889-027-00	10 pieces
37	28	0.18 mm	13.09 mm				889-028-00	10 pieces
37	29	0.18 mm	13.22 mm				889-029-00	10 pieces
37	30	0.18 mm	13.34 mm				889-030-00	10 pieces
37	31	0.18 mm	13.46 mm				889-031-00	10 pieces
37	32	0.18 mm	13.59 mm				889-032-00	10 pieces

* The Dentaform version of this prescription is not claimed to be a duplication of any other, nor does Dentaform imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.



Standard bands / Assortments

- universally applicable for right and left halves of the jaw
- unwelded bands (with or without double hook) ideal for the Herbst technique (welding service for double bands available free of charge on request)
- secure fit and optimal adhesion due to rough interior surface
- delivery in a tray (fits in band cabinet REF 890-414-00)
- individual welding, e.g. acc. to McLaughlin-Bennett-Trevisi* 18/22, is possible on request

We recommend our dental vernier (REF 042-751-00) for determining the largest mesio-distal diameter of the tooth crown.

Upper first and second bicuspid bands

Mini

Tooth	Band sizes (Contents per size)	welded		unwelded	Quantity
		REF (Roth 18) ¹⁾	REF (Roth 22) ²⁾	REF	
15-14 / 24-25	No. 1-32 (2)	860-064-18	860-064-22	864-064-00	64 pieces

Medium

Tooth	Band sizes (Contents per size)	welded		unwelded	Quantity
		REF (Roth 18) ¹⁾	REF (Roth 22) ²⁾	REF	
15-14 / 24-25	No. 1-4, 28-32 (2), 5-10, 21-27 (4), 11-12, 19-20 (6), 13, 18 (8), 14-17 (10)	860-150-18	860-150-22	864-150-00	150 pieces

Maxi

Tooth	Band sizes (Contents per size)	welded		unwelded	Quantity
		REF (Roth 18) ¹⁾	REF (Roth 22) ²⁾	REF	
15-14 / 24-25	No. 1-4, 28-32 (4), 5-10, 21-27 (8), 11-12, 19-20 (12), 13, 18 (16), 14-17 (20)			864-300-00	300 pieces

¹⁾ welded with equilibrium® 2 bracket 790-206-02 and double hook 750-602-10

²⁾ welded with equilibrium® 2 bracket 790-218-02 and double hook 750-602-10

* The Dentaureum version of this prescription is not claimed to be a duplication of any other, nor does Dentaureum imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.

Lower first bicuspid bands

Mini

Tooth	Band sizes (Contents per size)	welded		unwelded	Quantity
		REF (Roth 18) ¹⁾	REF (Roth 22) ²⁾	REF	
44 / 34	No. 1-32 (2)	861-064-18	861-064-22		64 pieces

Medium

Tooth	Band sizes (Contents per size)	welded		unwelded	Quantity
		REF (Roth 18) ¹⁾	REF (Roth 22) ²⁾	REF	
44 / 34	No. 1-4, 28-32 (2), 5-10, 21-27 (4), 11-12, 19-20 (6), 13, 18 (8), 14-17 (10)	861-150-18	861-150-22		150 pieces

¹⁾ welded with equilibrium® 2 bracket 790-210-02 and double hook 750-602-10

²⁾ welded with equilibrium® 2 bracket 790-222-02 and double hook 750-602-10

* The Dentaureum version of this prescription is not claimed to be a duplication of any other, nor does Dentaureum imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.

Standard bands / Assortments

- universally applicable for right and left halves of the jaw
- unwelded bands (with or without double hook) ideal for the Herbst technique (welding service for double bands available free of charge on request)
- secure fit and optimal adhesion due to rough interior surface
- delivery in a tray (fits in band cabinet REF 890-414-00)
- individual welding, e.g. acc. to McLaughlin-Bennett-Trevisi* 18/22, is possible on request

We recommend our dental vernier (REF 042-751-00) for determining the largest mesio-distal diameter of the tooth crown.

Lower second bicuspid bands

Mini

Tooth	Band sizes (Contents per size)	welded		unwelded	Quantity
		REF (Roth 18) ¹⁾	REF (Roth 22) ²⁾	REF	
45 / 35	No. 1-32 (2)	862-064-18	862-064-22		64 pieces

Medium

Tooth	Band sizes (Contents per size)	welded		unwelded	Quantity
		REF (Roth 18) ¹⁾	REF (Roth 22) ²⁾	REF	
45 / 35	No. 1-4, 28-32 (2), 5-10, 21-27 (4), 11-12, 19-20 (6), 13, 18 (8), 14-17 (10)	862-150-18	862-150-22		150 pieces

¹⁾ welded with equilibrium® 2 bracket 790-211-02 and double hook 750-602-10

²⁾ welded with equilibrium® 2 bracket 790-223-02 and double hook 750-602-10

* The Dentaurem version of this prescription is not claimed to be a duplication of any other, nor does Dentaurem imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.

Lower first and second bicuspid bands

Mini

Tooth	Band sizes (Contents per size)	welded		unwelded	Quantity
		REF (Roth 18)	REF (Roth 22)	REF	
45-44 / 34-35	No. 1-32 (2)			865-064-00	64 pieces

Medium

Tooth	Band sizes (Contents per size)	welded		unwelded	Quantity
		REF (Roth 18)	REF (Roth 22)	REF	
45-44 / 34-35	No. 1-4, 28-32 (2), 5-10, 21-27 (4), 11-12, 19-20 (6), 13, 18 (8), 14-17 (10)			865-150-00	150 pieces

Maxi

Tooth	Band sizes (Contents per size)	welded		unwelded	Quantity
		REF (Roth 18)	REF (Roth 22)	REF	
45-44 / 34-35	No. 1-4, 28-32 (4), 5-10, 21-27 (8), 11-12, 19-20 (12), 13, 18 (16), 14-17 (20)			865-300-00	300 pieces

* The Dentaurem version of this prescription is not claimed to be a duplication of any other, nor does Dentaurem imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.

Standard bands / Assortments

- universally applicable for right and left halves of the jaw
- secure fit and optimal adhesion due to rough interior surface
- delivery in a tray (fits in band cabinet REF 890-414-00)
- individual welding, e.g. acc. to McLaughlin-Bennett-Trevisi* 18/22, is possible on request

We recommend our dental vernier (REF 042-751-00) for determining the largest mesio-distal diameter of the tooth crown.

Upper first and second molar bands

Mini

Tooth	Band sizes (Contents per size)	welded		unwelded	Quantity
		REF (Roth 18)	REF (Roth 22)	REF	
17-16 / 26-27	No. 1-32 (2)			872-064-00	64 pieces

Medium

Tooth	Band sizes (Contents per size)	welded		unwelded	Quantity
		REF (Roth 18)	REF (Roth 22)	REF	
17-16 / 26-27	No. 1-8, 30-32 (2), 9-18 (4), 19-20 (6), 21-25, 28-29 (8), 26-27 (10)			872-150-00	150 pieces

* The Dentaureum version of this prescription is not claimed to be a duplication of any other, nor does Dentaureum imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.

Lower first and second molar bands

Mini

Tooth	Band sizes (Contents per size)	welded		unwelded	Quantity
		REF (Roth 18)	REF (Roth 22)	REF	
47-46 / 36-37	No. 1-32 (2)			873-064-00	64 pieces

Medium

Tooth	Band sizes (Contents per size)	welded		unwelded	Quantity
		REF (Roth 18)	REF (Roth 22)	REF	
47-46 / 36-37	No. 1-8, 30-32 (2), 9-18 (4), 19-20 (6), 21-25, 28-29 (8), 26-27 (10)			873-150-00	150 pieces

* The Dentaureum version of this prescription is not claimed to be a duplication of any other, nor does Dentaureum imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.



Standard bands / Assortments

- universally applicable for right and left halves of the jaw
- unwelded bands (with or without double hook) ideal for the Herbst technique (welding service for double bands available free of charge on request)
- secure fit and optimal adhesion due to rough interior surface
- delivery in a tray (fits in band cabinet REF 890-414-00)
- individual welding, e.g. acc. to McLaughlin-Bennett-Trevisi* 18/22, is possible on request

We recommend our dental vernier (REF 042-751-00) for determining the largest mesio-distal diameter of the tooth crown.

Upper cuspid bands

Mini

Tooth	Band sizes (Contents per size)	welded		unwelded	Quantity
		REF (Roth 18)	REF (Roth 22)	REF	
13 / 23	No. 1-20 (2)			855-050-00	50 pieces

Medium

Tooth	Band sizes (Contents per size)	welded		unwelded	Quantity
		REF (Roth 18)	REF (Roth 22)	REF	
13 / 23	No. 1-5; 19-20 (4), 6-8; 8.5; 9.5; 10.5; 11.5; 12.5; 16-18 (6), 9; 10; 11; 12; 13-15 (8)			855-150-00	150 pieces

* The Dentaaurum version of this prescription is not claimed to be a duplication of any other, nor does Dentaaurum imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.

Lower cuspid bands

Mini

Tooth	Band sizes (Contents per size)	welded		unwelded	Quantity
		REF (Roth 18)	REF (Roth 22)	REF	
43 / 33	No. 1-16 (2)			857-042-00	42 pieces

Medium

Tooth	Band sizes (Contents per size)	welded		unwelded	Quantity
		REF (Roth 18)	REF (Roth 22)	REF	
43 / 33	No. 1-2; 5.5; 6.5; 15-16 (4), 7.5; 8.5; 9.5 (6), 3-5; 12-14 (8), 6; 7; 8; 9; 10-11 (10)			857-150-00	150 pieces

* The Dentaaurum version of this prescription is not claimed to be a duplication of any other, nor does Dentaaurum imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.

Standard bands / Single packs

- universally applicable for right and left halves of the jaw
- unwelded bands (with or without double hook) ideal for the Herbst technique (welding service for double bands available free of charge on request)
- secure fit and optimal adhesion due to rough interior surface
- individual welding, e.g. acc. to McLaughlin-Bennett-Trevisi* 18/22, is possible on request

We recommend our dental vernier (REF 042-751-00) for determining the largest mesio-distal diameter of the tooth crown.

Upper first and second bicuspid bands, right and left

Tooth	Band size	Band thickness	Ø mesial-distal	welded		Quantity	unwelded	Quantity
				REF (Roth 18) ¹⁾	REF (Roth 22) ²⁾		REF	
15-14 / 24-25	1	0.12 mm	5.52 mm	860-101-18	860-101-22	1 piece	860-101-00	10 pieces
15-14 / 24-25	2	0.12 mm	5.57 mm	860-102-18	860-102-22	1 piece	860-102-00	10 pieces
15-14 / 24-25	3	0.12 mm	5.62 mm	860-103-18	860-103-22	1 piece	860-103-00	10 pieces
15-14 / 24-25	4	0.12 mm	5.67 mm	860-104-18	860-104-22	1 piece	860-104-00	10 pieces
15-14 / 24-25	5	0.12 mm	5.75 mm	860-105-18	860-105-22	1 piece	860-105-00	10 pieces
15-14 / 24-25	6	0.12 mm	5.82 mm	860-106-18	860-106-22	1 piece	860-106-00	10 pieces
15-14 / 24-25	7	0.12 mm	5.89 mm	860-107-18	860-107-22	1 piece	860-107-00	10 pieces
15-14 / 24-25	8	0.12 mm	5.94 mm	860-108-18	860-108-22	1 piece	860-108-00	10 pieces
15-14 / 24-25	9	0.12 mm	6.01 mm	860-109-18	860-109-22	1 piece	860-109-00	10 pieces
15-14 / 24-25	10	0.12 mm	6.08 mm	860-110-18	860-110-22	1 piece	860-110-00	10 pieces
15-14 / 24-25	11	0.12 mm	6.13 mm	860-111-18	860-111-22	1 piece	860-111-00	10 pieces
15-14 / 24-25	12	0.12 mm	6.20 mm	860-112-18	860-112-22	1 piece	860-112-00	10 pieces
15-14 / 24-25	13	0.12 mm	6.24 mm	860-113-18	860-113-22	1 piece	860-113-00	10 pieces
15-14 / 24-25	14	0.12 mm	6.30 mm	860-114-18	860-114-22	1 piece	860-114-00	10 pieces
15-14 / 24-25	15	0.12 mm	6.34 mm	860-115-18	860-115-22	1 piece	860-115-00	10 pieces
15-14 / 24-25	16	0.12 mm	6.42 mm	860-116-18	860-116-22	1 piece	860-116-00	10 pieces
15-14 / 24-25	17	0.12 mm	6.48 mm	860-117-18	860-117-22	1 piece	860-117-00	10 pieces
15-14 / 24-25	18	0.12 mm	6.56 mm	860-118-18	860-118-22	1 piece	860-118-00	10 pieces
15-14 / 24-25	19	0.12 mm	6.64 mm	860-119-18	860-119-22	1 piece	860-119-00	10 pieces
15-14 / 24-25	20	0.12 mm	6.71 mm	860-120-18	860-120-22	1 piece	860-120-00	10 pieces
15-14 / 24-25	21	0.12 mm	6.80 mm	860-121-18	860-121-22	1 piece	860-121-00	10 pieces
15-14 / 24-25	22	0.12 mm	6.90 mm	860-122-18	860-122-22	1 piece	860-122-00	10 pieces
15-14 / 24-25	23	0.12 mm	7.00 mm	860-123-18	860-123-22	1 piece	860-123-00	10 pieces
15-14 / 24-25	24	0.12 mm	7.11 mm	860-124-18	860-124-22	1 piece	860-124-00	10 pieces
15-14 / 24-25	25	0.12 mm	7.18 mm	860-125-18	860-125-22	1 piece	860-125-00	10 pieces
15-14 / 24-25	26	0.12 mm	7.26 mm	860-126-18	860-126-22	1 piece	860-126-00	10 pieces
15-14 / 24-25	27	0.12 mm	7.35 mm	860-127-18	860-127-22	1 piece	860-127-00	10 pieces
15-14 / 24-25	28	0.12 mm	7.43 mm	860-128-18	860-128-22	1 piece	860-128-00	10 pieces
15-14 / 24-25	29	0.12 mm	7.51 mm	860-129-18	860-129-22	1 piece	860-129-00	10 pieces
15-14 / 24-25	30	0.12 mm	7.61 mm	860-130-18	860-130-22	1 piece	860-130-00	10 pieces
15-14 / 24-25	31	0.12 mm	7.71 mm	860-131-18	860-131-22	1 piece	860-131-00	10 pieces
15-14 / 24-25	32	0.12 mm	7.80 mm	860-132-18	860-132-22	1 piece	860-132-00	10 pieces

¹⁾ welded with equilibrium® 2 bracket 790-206-02 and double hook 750-602-10

²⁾ welded with equilibrium® 2 bracket 790-218-02 and double hook 750-602-10



* The Dentaurem version of this prescription is not claimed to be a duplication of any other, nor does Dentaurem imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.

Bands

Standard bands / Single packs

- universally applicable for right and left halves of the jaw
- unwelded bands (with or without double hook) ideal for the Herbst technique (welding service for double bands available free of charge on request)
- secure fit and optimal adhesion due to rough interior surface
- individual welding, e.g. acc. to McLaughlin-Bennett-Trevisi* 18/22, is possible on request

We recommend our dental vernier (REF 042-751-00) for determining the largest mesio-distal diameter of the tooth crown.

Lower first bicuspid bands, right and left

Tooth	Band size	Band thickness	Ø mesial-distal	welded		Quantity	unwelded	Quantity
				REF (Roth 18) ¹⁾	REF (Roth 22) ²⁾		REF	
44 / 34	1	0.12 mm	6.34 mm	861-101-18	861-101-22	1 piece	861-101-00	10 pieces
44 / 34	2	0.12 mm	6.45 mm	861-102-18	861-102-22	1 piece	861-102-00	10 pieces
44 / 34	3	0.12 mm	6.56 mm	861-103-18	861-103-22	1 piece	861-103-00	10 pieces
44 / 34	4	0.12 mm	6.68 mm	861-104-18	861-104-22	1 piece	861-104-00	10 pieces
44 / 34	5	0.12 mm	6.80 mm	861-105-18	861-105-22	1 piece	861-105-00	10 pieces
44 / 34	6	0.12 mm	6.90 mm	861-106-18	861-106-22	1 piece	861-106-00	10 pieces
44 / 34	7	0.12 mm	7.02 mm	861-107-18	861-107-22	1 piece	861-107-00	10 pieces
44 / 34	8	0.12 mm	7.13 mm	861-108-18	861-108-22	1 piece	861-108-00	10 pieces
44 / 34	9	0.12 mm	7.25 mm	861-109-18	861-109-22	1 piece	861-109-00	10 pieces
44 / 34	10	0.12 mm	7.37 mm	861-110-18	861-110-22	1 piece	861-110-00	10 pieces
44 / 34	11	0.12 mm	7.50 mm	861-111-18	861-111-22	1 piece	861-111-00	10 pieces
44 / 34	12	0.12 mm	7.63 mm	861-112-18	861-112-22	1 piece	861-112-00	10 pieces
44 / 34	13	0.12 mm	7.75 mm	861-113-18	861-113-22	1 piece	861-113-00	10 pieces
44 / 34	14	0.12 mm	7.88 mm	861-114-18	861-114-22	1 piece	861-114-00	10 pieces
44 / 34	15	0.12 mm	8.05 mm	861-115-18	861-115-22	1 piece	861-115-00	10 pieces
44 / 34	16	0.12 mm	8.13 mm	861-116-18	861-116-22	1 piece	861-116-00	10 pieces
44 / 34	17	0.12 mm	8.26 mm	861-117-18	861-117-22	1 piece	861-117-00	10 pieces
44 / 34	18	0.12 mm	8.40 mm	861-118-18	861-118-22	1 piece	861-118-00	10 pieces
44 / 34	19	0.12 mm	8.53 mm	861-119-18	861-119-22	1 piece	861-119-00	10 pieces
44 / 34	20	0.12 mm	8.70 mm	861-120-18	861-120-22	1 piece	861-120-00	10 pieces
44 / 34	21	0.12 mm	8.84 mm	861-121-18	861-121-22	1 piece	861-121-00	10 pieces
44 / 34	22	0.12 mm	9.02 mm	861-122-18	861-122-22	1 piece	861-122-00	10 pieces
44 / 34	23	0.12 mm	9.17 mm	861-123-18	861-123-22	1 piece	861-123-00	10 pieces
44 / 34	24	0.12 mm	9.28 mm	861-124-18	861-124-22	1 piece	861-124-00	10 pieces
44 / 34	25	0.12 mm	9.38 mm	861-125-18	861-125-22	1 piece	861-125-00	10 pieces
44 / 34	26	0.12 mm	9.48 mm	861-126-18	861-126-22	1 piece	861-126-00	10 pieces
44 / 34	27	0.12 mm	9.59 mm	861-127-18	861-127-22	1 piece	861-127-00	10 pieces
44 / 34	28	0.12 mm	9.70 mm	861-128-18	861-128-22	1 piece	861-128-00	10 pieces
44 / 34	29	0.12 mm	9.82 mm	861-129-18	861-129-22	1 piece	861-129-00	10 pieces
44 / 34	30	0.12 mm	9.96 mm	861-130-18	861-130-22	1 piece	861-130-00	10 pieces
44 / 34	31	0.12 mm	10.11 mm	861-131-18	861-131-22	1 piece	861-131-00	10 pieces
44 / 34	32	0.12 mm	10.28 mm	861-132-18	861-132-22	1 piece	861-132-00	10 pieces

¹⁾ welded with equilibrium* 2 bracket 790-210-02 and double hook 750-602-10

²⁾ welded with equilibrium* 2 bracket 790-222-02 and double hook 750-602-10



* The Dentaurem version of this prescription is not claimed to be a duplication of any other, nor does Dentaurem imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevi.

Standard bands / Single packs

- universally applicable for right and left halves of the jaw
- unwelded bands (with or without double hook) ideal for the Herbst technique (welding service for double bands available free of charge on request)
- secure fit and optimal adhesion due to rough interior surface
- individual welding, e.g. acc. to McLaughlin-Bennett-Trevisi* 18/22, is possible on request

We recommend our dental vernier (REF 042-751-00) for determining the largest mesio-distal diameter of the tooth crown.

Lower second bicuspid bands, right and left

Tooth	Band size	Band thickness	Ø mesial-distal	welded		Quantity	unwelded	Quantity
				REF (Roth 18) ¹⁾	REF (Roth 22) ²⁾		REF	
45 / 35	1	0.12 mm	6.34 mm	862-101-18	862-101-22	1 piece	861-101-00	10 pieces
45 / 35	2	0.12 mm	6.45 mm	862-102-18	862-102-22	1 piece	861-102-00	10 pieces
45 / 35	3	0.12 mm	6.56 mm	862-103-18	862-103-22	1 piece	861-103-00	10 pieces
45 / 35	4	0.12 mm	6.68 mm	862-104-18	862-104-22	1 piece	861-104-00	10 pieces
45 / 35	5	0.12 mm	6.80 mm	862-105-18	862-105-22	1 piece	861-105-00	10 pieces
45 / 35	6	0.12 mm	6.90 mm	862-106-18	862-106-22	1 piece	861-106-00	10 pieces
45 / 35	7	0.12 mm	7.02 mm	862-107-18	862-107-22	1 piece	861-107-00	10 pieces
45 / 35	8	0.12 mm	7.13 mm	862-108-18	862-108-22	1 piece	861-108-00	10 pieces
45 / 35	9	0.12 mm	7.25 mm	862-109-18	862-109-22	1 piece	861-109-00	10 pieces
45 / 35	10	0.12 mm	7.37 mm	862-110-18	862-110-22	1 piece	861-110-00	10 pieces
45 / 35	11	0.12 mm	7.50 mm	862-111-18	862-111-22	1 piece	861-111-00	10 pieces
45 / 35	12	0.12 mm	7.63 mm	862-112-18	862-112-22	1 piece	861-112-00	10 pieces
45 / 35	13	0.12 mm	7.75 mm	862-113-18	862-113-22	1 piece	861-113-00	10 pieces
45 / 35	14	0.12 mm	7.88 mm	862-114-18	862-114-22	1 piece	861-114-00	10 pieces
45 / 35	15	0.12 mm	8.05 mm	862-115-18	862-115-22	1 piece	861-115-00	10 pieces
45 / 35	16	0.12 mm	8.13 mm	862-116-18	862-116-22	1 piece	861-116-00	10 pieces
45 / 35	17	0.12 mm	8.26 mm	862-117-18	862-117-22	1 piece	861-117-00	10 pieces
45 / 35	18	0.12 mm	8.40 mm	862-118-18	862-118-22	1 piece	861-118-00	10 pieces
45 / 35	19	0.12 mm	8.53 mm	862-119-18	862-119-22	1 piece	861-119-00	10 pieces
45 / 35	20	0.12 mm	8.70 mm	862-120-18	862-120-22	1 piece	861-120-00	10 pieces
45 / 35	21	0.12 mm	8.84 mm	862-121-18	862-121-22	1 piece	861-121-00	10 pieces
45 / 35	22	0.12 mm	9.02 mm	862-122-18	862-122-22	1 piece	861-122-00	10 pieces
45 / 35	23	0.12 mm	9.17 mm	862-123-18	862-123-22	1 piece	861-123-00	10 pieces
45 / 35	24	0.12 mm	9.28 mm	862-124-18	862-124-22	1 piece	861-124-00	10 pieces
45 / 35	25	0.12 mm	9.38 mm	862-125-18	862-125-22	1 piece	861-125-00	10 pieces
45 / 35	26	0.12 mm	9.48 mm	862-126-18	862-126-22	1 piece	861-126-00	10 pieces
45 / 35	27	0.12 mm	9.59 mm	862-127-18	862-127-22	1 piece	861-127-00	10 pieces
45 / 35	28	0.12 mm	9.70 mm	862-128-18	862-128-22	1 piece	861-128-00	10 pieces
45 / 35	29	0.12 mm	9.82 mm	862-129-18	862-129-22	1 piece	861-129-00	10 pieces
45 / 35	30	0.12 mm	9.96 mm	862-130-18	862-130-22	1 piece	861-130-00	10 pieces
45 / 35	31	0.12 mm	10.11 mm	862-131-18	862-131-22	1 piece	861-131-00	10 pieces
45 / 35	32	0.12 mm	10.28 mm	862-132-18	862-132-22	1 piece	861-132-00	10 pieces

¹⁾ welded with equilibrium® 2 bracket 790-211-02 and double hook 750-602-10

²⁾ welded with equilibrium® 2 bracket 790-223-02 and double hook 750-602-10



* The Dentaurem version of this prescription is not claimed to be a duplication of any other, nor does Dentaurem imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.

Bands

Standard bands / Single packs

- universally applicable for right and left halves of the jaw
- secure fit and optimal adhesion due to rough interior surface
- individual welding, e.g. acc. to McLaughlin-Bennett-Trevisi* 18/22, is possible on request

We recommend our dental vernier (REF 042-751-00) for determining the largest mesio-distal diameter of the tooth crown.

Upper first and second molar bands, right and left

Tooth	Band size	Band thickness	Ø mesial-distal	welded		Quantity	unwelded	Quantity
				REF (Roth 18)	REF (Roth 22)		REF	
17-16 / 26-27	1	0.15 mm	7.94 mm				870-001-00	10 pieces
17-16 / 26-27	2	0.15 mm	8.06 mm				870-002-00	10 pieces
17-16 / 26-27	3	0.15 mm	8.18 mm				870-003-00	10 pieces
17-16 / 26-27	4	0.15 mm	8.30 mm				870-004-00	10 pieces
17-16 / 26-27	5	0.15 mm	8.42 mm				870-005-00	10 pieces
17-16 / 26-27	6	0.15 mm	8.54 mm				870-006-00	10 pieces
17-16 / 26-27	7	0.15 mm	8.66 mm				870-007-00	10 pieces
17-16 / 26-27	8	0.15 mm	8.78 mm				870-008-00	10 pieces
17-16 / 26-27	9	0.15 mm	8.90 mm				870-009-00	10 pieces
17-16 / 26-27	10	0.15 mm	9.02 mm				870-010-00	10 pieces
17-16 / 26-27	11	0.15 mm	9.14 mm				870-011-00	10 pieces
17-16 / 26-27	12	0.15 mm	9.26 mm				870-012-00	10 pieces
17-16 / 26-27	13	0.15 mm	9.38 mm				870-013-00	10 pieces
17-16 / 26-27	14	0.15 mm	9.50 mm				870-014-00	10 pieces
17-16 / 26-27	15	0.15 mm	9.62 mm				870-015-00	10 pieces
17-16 / 26-27	16	0.15 mm	9.74 mm				870-016-00	10 pieces
17-16 / 26-27	17	0.15 mm	9.86 mm				870-017-00	10 pieces
17-16 / 26-27	18	0.15 mm	9.98 mm				870-018-00	10 pieces
17-16 / 26-27	19	0.15 mm	10.10 mm				870-019-00	10 pieces
17-16 / 26-27	20	0.15 mm	10.22 mm				870-020-00	10 pieces
17-16 / 26-27	21	0.15 mm	10.34 mm				870-021-00	10 pieces
17-16 / 26-27	22	0.15 mm	10.46 mm				870-022-00	10 pieces
17-16 / 26-27	23	0.15 mm	10.58 mm				870-023-00	10 pieces
17-16 / 26-27	24	0.15 mm	10.70 mm				870-024-00	10 pieces
17-16 / 26-27	25	0.15 mm	10.82 mm				870-025-00	10 pieces
17-16 / 26-27	26	0.15 mm	10.94 mm				870-026-00	10 pieces
17-16 / 26-27	27	0.15 mm	11.06 mm				870-027-00	10 pieces
17-16 / 26-27	28	0.15 mm	11.18 mm				870-028-00	10 pieces
17-16 / 26-27	29	0.15 mm	11.30 mm				870-029-00	10 pieces
17-16 / 26-27	30	0.15 mm	11.42 mm				870-030-00	10 pieces
17-16 / 26-27	31	0.15 mm	11.54 mm				870-031-00	10 pieces
17-16 / 26-27	32	0.15 mm	11.66 mm				870-032-00	10 pieces

* The Dentaaurum version of this prescription is not claimed to be a duplication of any other, nor does Dentaaurum imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.



Standard bands / Single packs

- universally applicable for right and left halves of the jaw
- secure fit and optimal adhesion due to rough interior surface
- individual welding, e.g. acc. to McLaughlin-Bennett-Trevisi* 18/22, is possible on request

We recommend our dental vernier (REF 042-751-00) for determining the largest mesio-distal diameter of the tooth crown.

Lower first and second molar bands, right and left

Tooth	Band size	Band thickness	Ø mesial-distal	welded		Quantity	unwelded	Quantity
				REF (Roth 18)	REF (Roth 22)		REF	
47-46 / 36-37	1	0.15 mm	8.94 mm				871-001-00	10 pieces
47-46 / 36-37	2	0.15 mm	9.06 mm				871-002-00	10 pieces
47-46 / 36-37	3	0.15 mm	9.18 mm				871-003-00	10 pieces
47-46 / 36-37	4	0.15 mm	9.30 mm				871-004-00	10 pieces
47-46 / 36-37	5	0.15 mm	9.42 mm				871-005-00	10 pieces
47-46 / 36-37	6	0.15 mm	9.54 mm				871-006-00	10 pieces
47-46 / 36-37	7	0.15 mm	9.56 mm				871-007-00	10 pieces
47-46 / 36-37	8	0.15 mm	9.78 mm				871-008-00	10 pieces
47-46 / 36-37	9	0.15 mm	9.90 mm				871-009-00	10 pieces
47-46 / 36-37	10	0.15 mm	10.02 mm				871-010-00	10 pieces
47-46 / 36-37	11	0.15 mm	10.14 mm				871-011-00	10 pieces
47-46 / 36-37	12	0.15 mm	10.25 mm				871-012-00	10 pieces
47-46 / 36-37	13	0.15 mm	10.38 mm				871-013-00	10 pieces
47-46 / 36-37	14	0.15 mm	10.50 mm				871-014-00	10 pieces
47-46 / 36-37	15	0.15 mm	10.62 mm				871-015-00	10 pieces
47-46 / 36-37	16	0.15 mm	10.74 mm				871-016-00	10 pieces
47-46 / 36-37	17	0.15 mm	10.86 mm				871-017-00	10 pieces
47-46 / 36-37	18	0.15 mm	10.98 mm				871-018-00	10 pieces
47-46 / 36-37	19	0.15 mm	11.10 mm				871-019-00	10 pieces
47-46 / 36-37	20	0.15 mm	11.22 mm				871-020-00	10 pieces
47-46 / 36-37	21	0.15 mm	11.34 mm				871-021-00	10 pieces
47-46 / 36-37	22	0.15 mm	11.46 mm				871-022-00	10 pieces
47-46 / 36-37	23	0.15 mm	11.58 mm				871-023-00	10 pieces
47-46 / 36-37	24	0.15 mm	11.70 mm				871-024-00	10 pieces
47-46 / 36-37	25	0.15 mm	11.82 mm				871-025-00	10 pieces
47-46 / 36-37	26	0.15 mm	11.94 mm				871-026-00	10 pieces
47-46 / 36-37	27	0.15 mm	12.06 mm				871-027-00	10 pieces
47-46 / 36-37	28	0.15 mm	12.18 mm				871-028-00	10 pieces
47-46 / 36-37	29	0.15 mm	12.30 mm				871-029-00	10 pieces
47-46 / 36-37	30	0.15 mm	12.42 mm				871-030-00	10 pieces
47-46 / 36-37	31	0.15 mm	12.54 mm				871-031-00	10 pieces
47-46 / 36-37	32	0.15 mm	12.66 mm				871-032-00	10 pieces

* The Dentaureum version of this prescription is not claimed to be a duplication of any other, nor does Dentaureum imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.



Bands

Standard bands / Single packs

- universally applicable for right and left halves of the jaw
- unwelded bands (with or without double hook) ideal for the Herbst technique (welding service for double bands available free of charge on request)
- secure fit and optimal adhesion due to rough interior surface
- individual welding, e.g. acc. to McLaughlin-Bennett-Trevisi* 18/22, is possible on request

We recommend our dental vernier (REF 042-751-00) for determining the largest mesio-distal diameter of the tooth crown.

Upper cuspid bands, right and left

Tooth	Band size	Band thickness	ø mesial-distal	welded		Quantity	unwelded	Quantity
				REF (Roth 18)	REF (Roth 22)		REF	
13 / 23	1	0.12 mm	7.10 mm				850-101-00	10 pieces
13 / 23	2	0.12 mm	7.20 mm				850-102-00	10 pieces
13 / 23	3	0.12 mm	7.30 mm				850-103-00	10 pieces
13 / 23	4	0.12 mm	7.50 mm				850-104-00	10 pieces
13 / 23	5	0.12 mm	7.60 mm				850-105-00	10 pieces
13 / 23	6	0.12 mm	7.80 mm				850-106-00	10 pieces
13 / 23	7	0.12 mm	7.90 mm				850-107-00	10 pieces
13 / 23	8	0.12 mm	8.00 mm				850-108-00	10 pieces
13 / 23	8.5	0.12 mm	8.10 mm				850-121-00	10 pieces
13 / 23	9	0.12 mm	8.20 mm				850-109-00	10 pieces
13 / 23	9.5	0.12 mm	8.30 mm				850-122-00	10 pieces
13 / 23	10	0.12 mm	8.40 mm				850-110-00	10 pieces
13 / 23	10.5	0.12 mm	8.50 mm				850-123-00	10 pieces
13 / 23	11	0.12 mm	8.60 mm				850-111-00	10 pieces
13 / 23	11.5	0.12 mm	8.70 mm				850-124-00	10 pieces
13 / 23	12	0.12 mm	8.80 mm				850-112-00	10 pieces
13 / 23	12.5	0.12 mm	8.90 mm				850-125-00	10 pieces
13 / 23	13	0.12 mm	9.10 mm				850-113-00	10 pieces
13 / 23	14	0.12 mm	9.20 mm				850-114-00	10 pieces
13 / 23	15	0.12 mm	9.40 mm				850-115-00	10 pieces
13 / 23	16	0.12 mm	9.50 mm				850-116-00	10 pieces
13 / 23	17	0.12 mm	9.80 mm				850-117-00	10 pieces
13 / 23	18	0.12 mm	9.90 mm				850-118-00	10 pieces
13 / 23	19	0.12 mm	10.00 mm				850-119-00	10 pieces
13 / 23	20	0.12 mm	10.10 mm				850-120-00	10 pieces

* The Dentaureum version of this prescription is not claimed to be a duplication of any other, nor does Dentaureum imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.



Standard bands / Single packs

- universally applicable for right and left halves of the jaw
- unwelded bands (with or without double hook) ideal for the Herbst technique (welding service for double bands available free of charge on request)
- secure fit and optimal adhesion due to rough interior surface
- individual welding, e.g. acc. to McLaughlin-Bennett-Trevisi* 18/22, is possible on request

We recommend our dental vernier (REF 042-751-00) for determining the largest mesio-distal diameter of the tooth crown.

Lower cuspid bands, right and left

Tooth	Band size	Band thickness	Ø mesial-distal	welded		Quantity	unwelded	Quantity
				REF (Roth 18)	REF (Roth 22)		REF	
43 / 33	1	0.12 mm	6.20 mm				851-101-00	10 pieces
43 / 33	2	0.12 mm	6.30 mm				851-102-00	10 pieces
43 / 33	3	0.12 mm	6.50 mm				851-103-00	10 pieces
43 / 33	4	0.12 mm	6.60 mm				851-104-00	10 pieces
43 / 33	5	0.12 mm	6.80 mm				851-105-00	10 pieces
43 / 33	5.5	0.12 mm	6.90 mm				851-117-00	10 pieces
43 / 33	6	0.12 mm	7.00 mm				851-106-00	10 pieces
43 / 33	6.5	0.12 mm	7.10 mm				851-118-00	10 pieces
43 / 33	7	0.12 mm	7.20 mm				851-107-00	10 pieces
43 / 33	7.5	0.12 mm	7.30 mm				851-119-00	10 pieces
43 / 33	8	0.12 mm	7.40 mm				851-108-00	10 pieces
43 / 33	8.5	0.12 mm	7.50 mm				851-120-00	10 pieces
43 / 33	9	0.12 mm	7.60 mm				851-109-00	10 pieces
43 / 33	9.5	0.12 mm	7.80 mm				851-121-00	10 pieces
43 / 33	10	0.12 mm	7.90 mm				851-110-00	10 pieces
43 / 33	11	0.12 mm	8.00 mm				851-111-00	10 pieces
43 / 33	12	0.12 mm	8.20 mm				851-112-00	10 pieces
43 / 33	13	0.12 mm	8.30 mm				851-113-00	10 pieces
43 / 33	14	0.12 mm	8.40 mm				851-114-00	10 pieces
43 / 33	15	0.12 mm	8.60 mm				851-115-00	10 pieces
43 / 33	16	0.12 mm	8.70 mm				851-116-00	10 pieces

* The Dentaaurum version of this prescription is not claimed to be a duplication of any other, nor does Dentaaurum imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.



Accessories



Band material

on roll

Soft, upper side highly polished, underside mat.

Thickness x width	Length	REF	Quantity
0.10 mm x 3.80 mm / 4 x 150	10 m	812-038-00	1 piece
0.12 mm x 4.55 mm / 5 x 180	10 m	813-045-00	1 piece
0.15 mm x 4.55 mm / 6 x 180	10 m	814-045-00	1 piece



Matrix band

on roll

Material No.: 1.4310

spring-hard

Thickness x width	Length	REF	Quantity
0.05 mm x 5.00 mm	1 m	831-601-00	1 piece
0.05 mm x 5.00 mm	10 m	831-701-00	1 piece
0.05 mm x 6.00 mm	1 m	831-606-00	1 piece
0.05 mm x 7.00 mm	1 m	831-603-00	1 piece
0.05 mm x 7.00 mm	10 m	831-703-00	1 piece

Band cabinet**Standard model plastic**

For storing preformed bands, 10 trays (without contents or vacuum-formed inserts).

Dimensions (H x W x D)	REF	Quantity
31 cm x 27.5 cm x 25 cm	890-414-00	1 piece

**Trays with lids and inserts**

For band cabinet (REF 890-414-00).

Marking	REF	Quantity
Standard upper cuspid bands	890-420-00	1 piece
Standard lower cuspid bands	890-421-00	1 piece
Standard upper bicuspid bands	890-422-00	1 piece
Standard lower bicuspid bands	890-423-00	1 piece
Standard upper molar bands	890-424-00	1 piece
Standard lower molar bands	890-425-00	1 piece
dentaform® upper first molar bands, right	890-426-00	1 piece
dentaform® upper first molar bands, left	890-427-00	1 piece
dentaform® upper second molar bands, right	890-428-00	1 piece
dentaform® upper second molar bands, left	890-429-00	1 piece
dentaform® lower first molar bands, right	890-430-00	1 piece
dentaform® lower first molar bands, left	890-431-00	1 piece
dentaform® lower second molar bands, right	890-432-00	1 piece
dentaform® lower second molar bands, left	890-433-00	1 piece
Direct bonding technique	890-412-00	1 piece
Drawer with plain insert	890-406-00	1 piece



Accessories



Stopper

For band cabinet (REF 890-414-00).

REF 890-414-10

1 piece



Head-end framework

For customized assembly of the band cabinet (REF 890-414-00).

Without stopper (REF 890-414-10).

REF 890-402-00

1 piece



Extension framework

For expanding the capacity of the band cabinet individually (REF 890-414-00).

Incl. stopper (REF 890-414-10).

REF 890-407-00

1 piece



Base plate

For customized assembly of the band cabinet (REF 890-414-00).

Without stopper (REF 890-414-10).

REF 890-403-00

1 piece



Adhesives.

For perfect bonding results.



DENTAURUM
**ONLINE
SHOP**
shop.dentaurum.com





Our adhesives for your dental practice.

The range of adhesives from Dentaureum gives users a wide selection for all bonding tasks in the orthodontic practice.

The light-curing adhesives CONTEC Ic and CONTEC lcr have been especially designed for bonding fixed orthodontic appliances – CONTEC Ic is for bonding metal, ceramic and plastic brackets and buccal tubes, CONTEC lcr is for bonding orthodontic retainers using adhesive. Both adhesives have a very high bond strength. An alternative to the light-curing adhesive CONTEC Ic is the chemically-curing bracket adhesive CONTEC go. CONTEC Ic contains no BisGMA (bisphenol A glycerolate dimethacrylate), TEGDMA (triethylene glycol dimethacrylate) or HEMA (hydroxyethyl methacrylate). CONTEC lcr and CONTEC go are free of TEGDMA and HEMA. The adhesives are therefore also suitable for patients with an allergy to methacrylates.

CONTEC

Best compatibility, no TEGDMA or HEMA.



Adhesives – areas of application.

	CONTEC lc	CONTEC lcr	CONTEC go
discovery® sl	+++	+	+++
discovery® sl 2.0	+++	+	+++
discovery®	+++	+	+++
discovery® smart	+++	+	+++
discovery® pearl	+++	+	+++
discovery® delight	+++	++	+
topic	+++	+	+++
equilibrium® 2 / ti	+++	+	+++
equilibrium® mini	+++	+	+++
Retainer	++	+++	—
Buccal tubes	++	+	++

+++ highly suitable
 ++ very suitable
 + suitable with restrictions
 — not suitable

Adhesives and accessories



CONTEC Ic set

Light-curing.

For bonding metal brackets and buccal tubes, ceramic brackets and plastic brackets.

For approx. 800 brackets (40 treatment cases).

A plastic primer must be used to bond plastic brackets.

Content:

- 3 x CONTEC Ic bracket adhesive, each 3 g (REF 163-111-00)
- 1 x CONTEC Ic / lcr primer, 5 ml (REF 163-113-10)
- 2 x CONTEC etching gel, each 2.5 ml (REF 163-115-00)
- 6 x CONTEC disposable tube for etching gel
- 50 x CONTEC disposable brush
- 1 x CONTEC combination instrument
- 1 x CONTEC mixing pad

When used with metal and ceramic brackets, CONTEC Ic is free of BisGMA, TEGDMA and HEMA.

CONTEC Ic, used with plastic brackets, is free of TEGDMA and HEMA.

REF 163-110-00

1 set



CONTEC Ic bracket adhesive

Light-curing.

For bonding metal brackets and buccal tubes, ceramic brackets and plastic brackets.

A plastic primer must be used to bond plastic brackets.

When used with metal and ceramic brackets, CONTEC Ic is free of BisGMA, TEGDMA and HEMA.

CONTEC Ic, used with plastic brackets, is free of TEGDMA and HEMA.

REF 163-111-00

2 pieces (3 g each)



CONTEC lcr set

Light-curing.

For bonding of retainers.

For approx. 80 retainers, braided or approx. 250 retainer 3-3.

Content:

- 2 x CONTEC lcr retainer adhesive, each 3.4 g (REF 163-131-00)
- 1 x CONTEC Ic / lcr primer, 5 ml (REF 163-113-10)
- 2 x CONTEC etching gel, each 2.5 ml (REF 163-115-00)
- 6 x CONTEC lcr disposable tube
- 6 x CONTEC disposable tube for etching gel
- 50 x CONTEC disposable brush
- 1 x CONTEC combination instrument
- 1 x CONTEC mixing pad

CONTEC lcr is free of TEGDMA and HEMA.

REF 163-130-00

1 set

CONTEC lcr retainer adhesive

Light-curing.

For bonding retainers.

Without disposable tube.

CONTEC lcr is free of TEGDMA and HEMA.

REF 163-131-00

3.4 g



CONTEC lc / lcr primer

Light-curing.

For CONTEC lc bracket adhesive (REF 163-111-00)
and CONTEC lcr retainer adhesive (REF 163-131-00).

REF 163-113-10

5 ml



CONTEC etching gel

35% phosphoric acid

For etching dental enamel.

Suitable for all CONTEC adhesives.

Without disposable tubes.

REF 163-115-00

2 pieces (2.5 ml each)



Adhesives and accessories



CONTEC go set

Chemically-curing

For bonding metal brackets and buccal tubes, ceramic brackets and plastic brackets.

For approx. 600 brackets (30 treatment cases).

Content:

- 3 x CONTEC go bracket adhesive, each 2.5 g (REF 163-121-00)
- 1 x CONTEC go primer, 5 ml (REF 163-114-00)
- 2 x CONTEC etching gel, je 2,5 ml (REF 163-115-00)
- 6 x CONTEC disposable tube for etching gel
- 50 x CONTEC disposable brush
- 1 x CONTEC combination instrument
- 1 x CONTEC mixing pad

CONTEC go is free of TEGDMA and HEMA.

Only compatible with REF 163-114-00 CONTEC go primer and REF 163-121-00 CONTEC go bracket adhesive. No other Dentaureum adhesive is compatible and may, therefore, not be used in combination.

REF 163-107-00

1 set

NEW



CONTEC go bracket adhesive

Chemically-curing

For bonding metal brackets and buccal tubes, ceramic brackets and plastic brackets.

CONTEC go is free of TEGDMA and HEMA.

Only compatible with REF 163-107-00 CONTEC go set and REF 163-114-00 CONTEC go primer. No other Dentaureum adhesive is compatible and may, therefore, not be used in combination.

REF 163-121-00

2 pieces (each 2.5 g)

NEW



CONTEC go primer

Chemically-curing

For bonding metal brackets and buccal tubes, ceramic brackets and plastic brackets

Only compatible with REF 163-107-00 CONTEC go set and REF 163-121-00 CONTEC go bracket adhesive. No other Dentaureum adhesive is compatible and may, therefore, not be used in combination.

REF 163-114-00

5 ml

NEW



CONTEC etching gel

35% phosphoric acid

For etching dental enamel.

Suitable for all CONTEC adhesives.

Without disposable tubes.

REF 163-115-00

2 pieces (2.5 ml each)

Adhesive residue remover

For gentle removal of residual adhesive and cement without damaging the enamel.

Use water cooling!

Recommended speed: 5000 min⁻¹ - max. 10000 min⁻¹

Form	REF	Quantity
conical	123-700-00	10 pieces
radial	123-701-00	10 pieces
pear, large	123-702-00	10 pieces
pear, small	123-703-00	10 pieces



Carbide bur sterilizable

For the removal of adhesive residue.

Shaft ø: 2.35 mm

Recommended speed: 40000 min⁻¹

Form	REF	Quantity	REF	Quantity
fine	123-604-00	1 piece	123-604-30	10 pieces
coarse	123-603-00	1 piece	123-603-30	10 pieces



Product range

Arches / Straight wires / Wires on coils	Page	192
--	------	-----

Fixed appliance technique

Arches

Nickel titanium	Page	200
Thermo-active nickel titanium	Page	207
Stainless steel (NiCr)	Page	211
Cobalt-chrome	Page	213
Nickel-free* titanium molybdenum	Page	214
Straight wires	Page	217
Wires on coils	Page	219
Retainer	Page	220
Accessories	Page	222

Removable appliance technique

Straight wires	Page	227
Wires on coils	Page	229
Accessories	Page	231


Palatal and lingual arch systems	Page	236
----------------------------------	------	-----

* see Materials list Orthodontics




Wires | Arches Accessories.

Highest level of precision
for reliable results.



DENTAURUM
**ONLINE
SHOP**
shop.dentaurum.com





Arches for every kind of treatment.

The choice of arch during therapy has a decisive influence on the optimal progress of an orthodontic treatment. That is why Dentaureum offers a wide range of high-quality arches made of different materials for the treatment of all kinds of situations. The depth of the program with respect to shapes, dimensions and strengths is unique as well.

Dentaureum's assortment not only includes wires for the fixed technique, but also for the orthodontic laboratory – from clasp wires to pre-contoured wire elements for the immediate use on removable appliances.

In addition, the wire and arch program is complemented by highly biocompatible materials like rematitan® SPECIAL. These materials are the perfect complement to the nickel-free brackets and in combination enable a complete nickel-free treatment of sensitive allergy patients.

**rematitan® LITE /
rematitan® LITE White**
Nickel titanium.

rematitan® LITE Spee®
Nickel titanium.

Tensic® / Tensic® White
Thermo-active nickel titanium.

Equire Thermo-Active
Thermo-active nickel titanium.

dentaflex®
Stainless steel (twisted and braided).

remanium® / Equire
Stainless steel (NiCr).

remaloy®
Cobalt-chrome.

Noninium®
Nickel-free* stainless steel wire.

rematitan® SPECIAL
Titanium molybdenum.

* Nickel trace elements of $\leq 0.2\%$ are due to the process of manufacturing

Suggested applications for Dentaureum wires and arches.

Wire material	Initial phase			Intermediate phase				Final phase	
	Levelling	Torque	Rotation	Dentoalveolar compensation	Torque	Rotation	Retraction with loops	Dentoalveolar compensation	Torque
Nickel titanium									
rematitan® LITE round	+		+						
rematitan® LITE rectangular		+		+	+	+		+	+
rematitan® sl round	+		+						
Nickel titanium (white coated)									
rematitan® LITE White round	+		+			+			
rematitan® LITE White rectangular		+			+	+			
Thermo-active nickel titanium									
Tensic® round	+		+			+			
Tensic® rectangular	+	+			+	+			
Equire Thermo-Active round	+		+			+			
Equire Thermo-Active rectangular	+	+			+	+			
Thermo-active nickel titanium (white coated)									
Tensic® White round	+		+			+			
Tensic® White rectangular	+	+			+	+			
Stainless steel (twisted and braided)									
dentaflex® 3 strands round	+								
dentaflex® 6 strands round	+								
dentaflex® 3 strands rectangular		+							
dentaflex® 8 strands rectangular	+	+							
Stainless steel									
remanium® round				+					
remanium® rectangular				+	+		+	+	+
Equire round				+					
Equire rectangular				+	+		+	+	+
Cobalt-chrome									
remaloy® rectangular					+		+		
Nickel-free titanium molybdenum									
rematitan® SPECIAL round				+					
rematitan® SPECIAL rectangular		+			+		+		+


Strength categories.

Strength category	N/mm ²	rematitan® LITE	Tensic®	Equire	dentaflex®	remanium®	remaloy®	Noninium®	rematitan® SPECIAL
middle plus 0	1100 – < 1400								+
hard 1	1400 – < 1600					+	+	+	
hard plus 2	1600 – < 1800				+		+		
spring hard 3	1800 – < 2000					+	+	+	
spring hard plus 4	2000 – < 2300			+		+			
extra spring hard 5	2000 – < 2500					+			
super spring hard 6	2500 – < 2700				+				
super elastic 7	–	+	+	+					





Arch forms in comparison.

Illustrations are for reference purposes only. The actual product may deviate from the illustration.



- Available in rematitan® LITE only.
- Same shape for maxilla and mandible.



¹ Damon® is a registered trademark of Ormco BV.

² Available in rematitan® LITE only.

³ The Dentaaurum version of this prescription is not claimed to be a duplication of any other, nor does Dentaaurum imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.

A brief overview of the program.

rematitan® LITE – nickel titanium.

rematitan® LITE arches are made from a super elastic, nickel titanium alloy. The super elasticity allows for reversibility and an even and gentle transfer of force. Due to its elastic recovery, rematitan® LITE arches are ideal in the levelling phase.

Available in the strength categories 0 1 2 3 4 5 6 7

Tensic® – thermo-active nickel titanium.

Tensic® arches are super elastic and very easy to bend in martensitic state at room temperature. The transformation temperature, which is close to body temperature, is responsible for the transfer of constant, gentle forces, which prevents the risk of root resorption during the levelling phase.

Available in the strength categories 0 1 2 3 4 5 6 7

Equire – thermo-active nickel-titanium (American Style $\triangle$ Damon® Style¹).

Equire nickel titanium arches are super elastic. Apart from their arch form they are identical to Tensic® arches and may be used throughout the entire treatment, even when the wire is frequently replaced. Equire nickel titanium arches are easy to form and weld in a martensitic state at room temperature and have smooth, low friction surfaces. The arch form complies with the American orthodontic ideal size.

Available in the strength categories 0 1 2 3 4 5 6 7

dentaflex® – stainless steel (twisted and braided).

dentaflex® wires and arches are available as triple-stranded and twisted wires. They are slightly stronger than the 6-strand, coaxial version and offer a well proven alternative to NiTi wires. dentaflex® 8-strand braided arches are available in addition. These particular arches are especially useful for torque during the initial stages.

Available in the strength categories 0 1 2 3 4 5 6 7

remanium® – stainless steel.

remanium® wires and arches are the best solution when stronger forces are required or when precise torque control is required towards the end of treatment. The wide spectrum of strength categories available in this series allows the use of remanium® wires and arches throughout the entire treatment, even when the wire is frequently replaced. They are easy to form and weld and have smooth surfaces with low friction.

Available in the strength categories 0 1 2 3 4 5 6 7

Equire – stainless steel (American Style $\triangle$ Damon® Style¹).

Equire stainless steel arches are best suitable when forces are required such as those produced by remanium® arches. They may be used throughout the entire treatment, even when the wire is frequently replaced. Equire arches are easy to form, can be welded and have smooth, low friction surfaces. The arch shape complies with the American orthodontic ideal size.

Available in the strength categories 0 1 2 3 4 5 6 7

remaloy® – cobalt-chrome.

remaloy® wires and arches are made of a cobalt chrome alloy, which can be tempered and are extremely break-resistant. This makes bending loops with fingers or pliers much easier. The surfaces are smooth, and low in friction.

Available in the strength categories 0 1 2 3 4 5 6 7

Noninium® – nickel-free² stainless steel wire.

Noninium® wires are especially designed for patients who are sensitized to nickel as they do not contain nickel². They are very break-resistant and have outstanding mechanical properties.

Available in the strength categories 0 1 2 3 4 5 6 7

rematitan® SPECIAL – titanium molybdenum.

The rematitan® SPECIAL wires and arches are made of a titanium molybdenum nickel-free² alloy, and are therefore suitable for patients who are sensitized to nickel. This alloy's modulus of elasticity, which is much lower compared to that of stainless steel wires, creates significantly lower forces but at the same time enables high elasticity and easy bending properties.

Available in the strength categories 0 1 2 3 4 5 6 7

¹ Damon® is a registered trademark of Ormco BV

² Nickel trace elements of $\leq 0.2\%$ are due to the process of manufacturing

Delivery program arches

Nickel titanium, super elastic

ø	REF (Maxillary)	REF (Mandibular)	Quantity
---	--------------------	---------------------	----------

rematitan® LITE ideal arch, round

12	766-058-00	766-059-00	10 pieces
14	766-080-00	766-081-00	10 pieces
14	766-080-01	766-081-01	50 pieces
16	766-082-00	766-083-00	10 pieces
16	766-082-01	766-083-01	50 pieces
18	766-084-00	766-085-00	10 pieces
18	766-084-01	766-085-01	50 pieces
20	768-100-00	768-101-00	10 pieces
20	768-100-01	768-101-01	50 pieces

rematitan® LITE ideal arch, rectangular

16 x 16	766-088-00	766-089-00	10 pieces
16 x 16	766-088-01	766-089-01	50 pieces
16 x 22	766-090-00	766-091-00	10 pieces
16 x 22	766-090-01	766-091-01	50 pieces
17 x 25	766-092-00	766-093-00	10 pieces
18 x 18	768-102-00	768-103-00	10 pieces
18 x 18	768-102-01	768-103-01	50 pieces
18 x 25	768-104-00	768-105-00	10 pieces
18 x 25	768-104-01	768-105-01	50 pieces
19 x 25	768-106-00	768-107-00	10 pieces
19 x 25	768-106-01	768-107-01	50 pieces
20 x 20	768-108-00	768-109-00	10 pieces
20 x 20	768-108-01	768-109-01	50 pieces
21 x 25	768-110-00	768-111-00	10 pieces
21 x 25	768-110-01	768-111-01	50 pieces

rematitan® LITE ideal arch, round, with dimple

12	766-060-00	766-061-00	10 pieces
14	766-062-00	766-063-00	10 pieces
16	766-064-00	766-065-00	10 pieces
18	766-066-00	766-067-00	10 pieces
20	768-112-00	768-113-00	10 pieces

rematitan® LITE ideal arch, rectangular, with dimple

16 x 16	766-068-00	766-069-00	10 pieces
16 x 22	766-070-00	766-071-00	10 pieces
17 x 25	766-072-00	766-073-00	10 pieces
18 x 18	766-074-00	766-075-00	10 pieces
18 x 25	766-076-00	766-077-00	10 pieces
19 x 25	766-078-00	766-079-00	10 pieces
20 x 20	768-114-00	768-115-00	10 pieces
21 x 25	766-086-00	766-087-00	10 pieces

rematitan® LITE ideal arch, round, narrow form, with dimple			
12	766-160-00	766-161-00	10 pieces
14	766-162-00	766-163-00	10 pieces
16	766-164-00	766-165-00	10 pieces
18	766-166-00	766-167-00	10 pieces
20	768-116-00	768-117-00	10 pieces

rematitan® LITE ideal arch, rectangular, narrow form, with dimple

16 x 16	766-168-00	766-169-00	10 pieces
16 x 22	766-170-00	766-171-00	10 pieces
17 x 25	766-172-00	766-173-00	10 pieces
18 x 18	766-174-00	766-175-00	10 pieces
18 x 25	766-176-00	766-177-00	10 pieces
19 x 25	766-178-00	766-179-00	10 pieces
20 x 20	768-118-00	768-119-00	10 pieces
21 x 25	768-120-00	768-121-00	10 pieces

rematitan® LITE Spee® retraction arch, round, with dimple

12	766-014-00	766-015-00	10 pieces
14	766-018-00	766-019-00	10 pieces
16	766-050-00	766-051-00	10 pieces
18	766-052-00	766-053-00	10 pieces

rematitan® LITE Spee® retraction arch, rectangular, with dimple

16 x 22	766-054-00	766-055-00	10 pieces
17 x 25	766-056-00	766-057-00	10 pieces
18 x 18	766-032-00	766-033-00	10 pieces
18 x 25	766-034-00	766-037-00	10 pieces
19 x 25	766-036-00	766-039-00	10 pieces

Nickel titanium, super elastic

☑	REF	Quantity
rematitan® LITE arch, arch form I, pointed, rectangular		
16 x 16	766-180-00	10 pieces
16 x 22	766-181-00	10 pieces
17 x 25	766-182-00	10 pieces

rematitan® LITE arch, arch form II, square, rectangular		
16 x 16	766-183-00	10 pieces
16 x 22	766-184-00	10 pieces
17 x 25	766-185-00	10 pieces

rematitan® LITE arch, arch form III, oval, rectangular		
16 x 16	766-186-00	10 pieces
16 x 22	766-187-00	10 pieces
17 x 25	766-188-00	10 pieces

rematitan® LITE Ricketts® PENTA arch, normal, rectangular		
16 x 16	766-250-00	10 pieces
16 x 22	766-251-00	10 pieces
17 x 25	766-252-00	10 pieces

rematitan® LITE Ricketts® PENTA arch, pointed, rectangular		
16 x 16	766-253-00	10 pieces
16 x 22	766-254-00	10 pieces
17 x 25	766-255-00	10 pieces

rematitan® LITE Ricketts® PENTA arch, oval, rectangular		
16 x 16	766-256-00	10 pieces
16 x 22	766-257-00	10 pieces
17 x 25	766-258-00	10 pieces

rematitan® LITE Ricketts® PENTA arch, extra pointed, rectangular		
16 x 16	766-259-00	10 pieces
16 x 22	766-260-00	10 pieces
17 x 25	766-261-00	10 pieces

rematitan® LITE Ricketts® PENTA arch, oval narrow, rectangular		
16 x 16	766-262-00	10 pieces
16 x 22	766-263-00	10 pieces
17 x 25	766-264-00	10 pieces

Thermo-active nickel titanium, super elastic

ø	REF (Maxillary)	REF (Mandibular)	Quantity
Tensic® ideal arch, round			
12	766-700-00	766-701-00	10 pieces
14	766-702-00	766-703-00	10 pieces
16	766-704-00	766-705-00	10 pieces
18	766-706-00	766-707-00	10 pieces

Tensic® ideal arch, rectangular			
16 x 16	766-710-00	766-711-00	10 pieces
16 x 22	766-712-00	766-713-00	10 pieces
17 x 25	766-714-00	766-715-00	10 pieces
18 x 25	766-716-00	766-717-00	10 pieces
19 x 25	766-718-00	766-719-00	10 pieces
21 x 25	766-720-00	766-721-00	10 pieces

Equire thermo-active, preformed ideal arch, round, arch form: American Style ± Damon Style*			
12	766-800-00	766-801-00	10 pieces
14	766-802-00	766-803-00	10 pieces
16	766-804-00	766-805-00	10 pieces
18	766-806-00	766-807-00	10 pieces

Equire thermo-active, preformed ideal arch, rectangular, arch form: American Style ± Damon Style*			
16 x 16	766-808-00	766-809-00	10 pieces
16 x 22	766-810-00	766-811-00	10 pieces
17 x 25	766-812-00	766-813-00	10 pieces
18 x 25	766-814-00	766-815-00	10 pieces
19 x 25	766-816-00	766-817-00	10 pieces
21 x 25	766-818-00	766-819-00	10 pieces

*Damon® is a registered trademark of Ormco BV

Delivery program arches

Stainless steel (twisted and braided), super spring hard

ø	REF (Maxillary)	REF (Mandibular)	Quantity
---	--------------------	---------------------	----------

dentaflex® ideal arch, round, 3-strand twisted

15	766-120-00	766-121-00	10 pieces
15	766-120-01	766-121-01	50 pieces
18	766-122-00	766-123-00	10 pieces
18	766-122-01	766-123-01	50 pieces
20	766-124-00	766-125-00	10 pieces

dentaflex® ideal arch, rectangular, 3-strand twisted

16 x 16	766-130-00	766-131-00	10 pieces
16 x 22	766-132-00	766-133-00	10 pieces
16 x 22	766-132-01	766-133-01	50 pieces

dentaflex® ideal arch, rectangular, 8-strand braided

16 x 16	766-100-00	766-200-00	10 pieces
16 x 16	766-100-01	766-200-01	50 pieces
16 x 22	766-101-00	766-201-00	10 pieces
16 x 22	766-101-01	766-201-01	50 pieces
17 x 25	766-103-00	766-203-00	10 pieces
18 x 25	766-104-00	766-204-00	10 pieces
19 x 25	766-105-00	766-205-00	10 pieces

Stainless steel (twisted and braided), hard plus

ø	REF (Maxillary)	REF (Mandibular)	Quantity
---	--------------------	---------------------	----------

dentaflex® ideal arch, round, 6-strand co-axial

15	766-300-00	766-400-00	10 pieces
15	766-300-01	766-400-01	50 pieces
18	766-301-00	766-401-00	10 pieces
18	766-301-01	766-401-01	50 pieces

Stainless steel (nickel chromium), spring hard plus

ø	REF (Maxillary)	REF (Mandibular)	Quantity
---	--------------------	---------------------	----------

remanium® ideal arch, round

14	765-200-00	765-300-00	10 pieces
16	765-201-00	765-301-00	10 pieces
16	765-201-01	765-301-01	50 pieces
18	765-202-00	765-302-00	10 pieces
18	765-202-01	765-302-01	50 pieces
20	765-203-00	765-303-00	10 pieces

remanium® ideal arch, rectangular

16 x 16	767-100-00	767-200-00	10 pieces
16 x 16	767-100-01	767-200-01	50 pieces
16 x 22	767-101-00	767-201-00	10 pieces
16 x 22	767-101-01	767-201-01	50 pieces
17 x 22	767-102-00	767-202-00	10 pieces
17 x 25	767-103-00	767-203-00	10 pieces
17 x 25	767-103-01	767-203-01	50 pieces
18 x 25	767-104-00	767-204-00	10 pieces
19 x 25	767-105-00	767-205-00	10 pieces
19 x 25	767-105-01	767-205-01	50 pieces
21.5 x 25	767-106-00	767-206-00	10 pieces
21.5 x 28	767-107-00	767-207-00	10 pieces

Equire preformed ideal arch, round, arch form:

American Style ± Damon Style*

12	769-000-00	769-001-00	10 pieces
14	769-002-00	769-003-00	10 pieces
16	769-004-00	769-005-00	10 pieces
18	769-006-00	769-007-00	10 pieces

Equire preformed ideal arch, rectangular, arch form:

American Style ± Damon Style*

16 x 16	769-008-00	769-009-00	10 pieces
16 x 22	769-010-00	769-011-00	10 pieces
17 x 25	769-012-00	769-013-00	10 pieces
18 x 25	769-014-00	769-015-00	10 pieces
19 x 25	769-016-00	769-017-00	10 pieces
21 x 25	769-018-00	769-019-00	10 pieces

*Damon® is a registered trademark of Ormco BV

Cobalt chromium, hard plus

□	REF (Maxillary)	REF (Mandibular)	Quantity
---	--------------------	---------------------	----------

remaloy® ideal arch, rectangular

16 x 16	769-100-00	769-200-00	10 pieces
16 x 22	769-101-00	769-201-00	10 pieces

Nickel-free titanium molybdenum, middle plus

ø	REF (Maxillary)	REF (Mandibular)	Quantity
---	--------------------	---------------------	----------

rematitan® SPECIAL ideal arch, round

16	766-605-00	766-606-00	10 pieces
18	766-607-00	766-608-00	10 pieces

rematitan® SPECIAL ideal arch, rectangular

16 x 16	766-631-00	766-632-00	10 pieces
16 x 22	766-609-00	766-610-00	10 pieces
17 x 25	766-611-00	766-612-00	10 pieces
18 x 18	766-613-00	766-614-00	10 pieces
18 x 25	766-615-00	766-616-00	10 pieces
19 x 25	766-617-00	766-618-00	10 pieces
21 x 25	766-619-00	766-620-00	10 pieces

Nickel titanium (white coated), super elastic

ø	REF (Maxillary)	REF (Mandibular)	Quantity
---	--------------------	---------------------	----------

rematitan® LITE White ideal arch, round

16	768-014-00	768-015-00	3 pieces
18	768-016-00	768-017-00	3 pieces

rematitan® LITE White ideal arch, rectangular

16 x 16	768-018-00	768-019-00	3 pieces
16 x 22	768-020-00	768-021-00	3 pieces
17 x 25	768-022-00	768-023-00	3 pieces

Thermo-active nickel titanium (white coated),
super elastic

ø	REF (Maxillary)	REF (Mandibular)	Quantity
---	--------------------	---------------------	----------

Tensic® White ideal arch, round

16	768-026-00	768-027-00	3 pieces
18	768-028-00	768-029-00	3 pieces

Tensic® White ideal arch, rectangular

16 x 16	768-030-00	768-031-00	3 pieces
16 x 22	768-032-00	768-033-00	3 pieces
17 x 25	768-034-00	768-035-00	3 pieces
18 x 25	768-038-00	768-039-00	3 pieces
19 x 25	768-040-00	768-041-00	3 pieces

Delivery program straight wires

Stainless steel (twisted and braided)

ø	REF	Quantity
---	-----	----------

dentaflex® straight wire, round, 3-strand twisted, super spring hard

0.38 mm / 15	545-738-00	10 pieces
0.45 mm / 18	545-745-00	10 pieces
0.50 mm / 20	545-750-00	10 pieces

dentaflex® straight wire, round, 6-strand co-axial, hard plus

0.38 mm / 15	545-638-00	10 pieces
0.45 mm / 18	545-645-00	10 pieces

rematain® flat retainer wire, 6-strand braided

10 x 29	530-000-00	10 pieces
---------	------------	-----------

Stainless steel (nickel chromium)

ø	REF	Quantity
---	-----	----------

remanium® straight wire, round, extra spring hard

0.35 mm / 14	535-035-00	25 pieces
0.40 mm / 16	535-040-00	25 pieces
0.45 mm / 18	535-045-00	25 pieces

remanium® straight wire, rectangular, spring hard

16 x 16	537-340-00	10 pieces
16 x 22	537-310-00	10 pieces
17 x 22	537-311-00	10 pieces
17 x 25	537-313-00	10 pieces
18 x 25	537-314-00	10 pieces
19 x 19	535-048-00	25 pieces
20 x 20	535-051-00	25 pieces
21.5 x 28	537-324-00	10 pieces

remanium® straight wire, rectangular, spring hard plus

16 x 16	537-440-00	10 pieces
16 x 22	537-410-00	10 pieces
17 x 22	537-411-00	10 pieces
17 x 25	537-413-00	10 pieces
18 x 25	537-414-00	10 pieces
21.5 x 25	537-423-00	10 pieces
21.5 x 28	537-424-00	10 pieces

remanium® straight wire, round, spring hard

0.40 mm / 16	525-040-00	25 pieces
0.50 mm / 20	535-002-00	25 pieces
0.60 mm / 23	535-003-00	25 pieces
0.70 mm / 28	527-070-00	25 pieces
0.70 mm / 28	527-005-00	900 pieces
0.70 mm / 28	527-004-00	1100 pieces
0.75 mm / 30	535-005-00	25 pieces
0.80 mm / 31	527-080-00	25 pieces
0.80 mm / 31	527-006-00	875 pieces
0.90 mm / 36	527-090-00	25 pieces
1.00 mm / 39	527-100-00	25 pieces
1.15 mm / 45	535-009-00	25 pieces

Cobalt chromium

☑	REF	Quantity
---	-----	----------

remaloy® straight wire, rectangular, hard plus

16 x 16	537-540-00	10 pieces
16 x 22	537-510-00	10 pieces

remaloy® straight wire, round, hard

0.70 mm / 28	528-070-00	10 pieces
0.80 mm / 31	528-080-00	10 pieces
0.90 mm / 36	528-090-00	10 pieces
1.00 mm / 39	528-100-00	10 pieces
1.10 mm / 43	528-110-00	10 pieces
1.20 mm / 47	528-120-00	10 pieces
1.30 mm / 51	528-130-00	10 pieces
1.50 mm / 59	528-150-00	10 pieces

remaloy® straight wire, half-round, spring hard

1.50 mm x 0.75 mm / 59 x 30	528-155-00	10 pieces
1.75 mm x 0.90 mm / 69 x 36	528-158-00	10 pieces

remaloy® straight wire, rectangular rounded, hard plus

1.92 mm x 0.90 mm / 76 x 36	528-159-00	10 pieces
-----------------------------	------------	-----------

Delivery program straight wires

Nickel-free titanium molybdenum

ø	REF	Quantity
---	-----	----------

rematitan® SPECIAL straight wire, round, middle plus

0.40 mm / 16	766-600-00	10 pieces
0.45 mm / 18	766-601-00	10 pieces
0.80 mm / 32	766-629-00	10 pieces

rematitan® SPECIAL straight wire, rectangular, middle plus

16 x 16	766-633-00	10 pieces
16 x 22	766-602-00	10 pieces
17 x 25	766-603-00	10 pieces
19 x 25	766-604-00	10 pieces

Titanium

ø	REF	Quantity
---	-----	----------

Titanium retainer wire, 3-strand twisted, Ti1, round

0.44 mm / 17	528-001-01	10 pieces
0.50 mm / 20	528-000-01	10 pieces

Titanium retainer wire, 3-strand twisted, Ti5, round

0.44 mm / 17	528-001-00	10 pieces
0.50 mm / 20	528-000-00	10 pieces

Gold

ø	REF	Quantity
---	-----	----------

Gold retainer, 3-strand twisted, round

0.50 mm / 20	529-000-00	1 piece
--------------	-------------------	---------

Dimensions

Inch		mm
10	round	0.25 mm
12	round	0.30 mm
14	round	0.35 mm
15	round	0.38 mm
16	round	0.40 mm
18	round	0.45 mm
20	round	0.50 mm
23	round	0.60 mm
28	round	0.70 mm
32	round	0.80 mm
36	round	0.90 mm
39	round	1.00 mm
43	round	1.10 mm
47	round	1.20 mm
51	round	1.30 mm
59	round	1.50 mm

Inch		mm
16 x 16	rectangular	0.41 mm x 0.41 mm
16 x 22	rectangular	0.41 mm x 0.56 mm
17 x 22	rectangular	0.43 mm x 0.56 mm
17 x 25	rectangular	0.43 mm x 0.64 mm
18 x 25	rectangular	0.46 mm x 0.64 mm
19 x 25	rectangular	0.48 mm x 0.64 mm
21 x 25	rectangular	0.53 mm x 0.64 mm
21.5 x 25	rectangular	0.55 mm x 0.64 mm
21.5 x 28	rectangular	0.55 mm x 0.70 mm

Delivery program wires on coils

Nickel titanium

ø	REF	Quantity
rematitan® LITE wire on coils, round, super elastic		
0.30 mm / 12	766-094-00	4.5 m
0.35 mm / 14	766-095-00	4.5 m
0.40 mm / 16	766-196-00	4.5 m

Stainless steel (twisted and braided)

ø	REF	Quantity
dentaflex® wire on coils, round, 3-strand twisted, super spring hard		
0.38 mm / 15	545-738-01	4 m
0.45 mm / 18	545-745-01	4 m

dentaflex® wire on coils, round, 6-strand co-axial, super spring hard		
0.38 mm / 15	545-638-01	4 m
0.45 mm / 18	545-645-01	4 m

Stainless steel (nickel chromium)

ø	REF	Quantity
remanium® wire on coils, round, extra spring hard		
0.35 mm / 14	530-035-00	7.5 m
0.40 mm / 16	530-040-00	7.5 m
0.45 mm / 18	530-045-00	7.5 m

remanium® ligature wire, round, soft		
0.25 mm / 10	500-025-00	80 m
0.25 mm / 10	501-025-00	1160 m
0.30 mm / 12	500-030-00	60 m
0.40 mm / 16	500-040-00	30 m
0.50 mm / 20	503-050-00	50 m
0.50 mm / 20	501-050-00	320 m

remanium® laboratory coil, round, hard		
0.70 mm / 28	513-070-00	30 m
0.80 mm / 31	513-080-00	20 m
0.90 mm / 36	513-090-00	10 m
1.00 mm / 39	513-100-00	10 m
1.20 mm / 47	513-120-00	10 m
1.50 mm / 59	513-150-00	10 m

remanium® clinical coil, round, hard		
0.70 mm / 28	514-070-00	165 m
0.80 mm / 31	514-080-00	125 m
0.90 mm / 36	514-090-00	100 m
1.00 mm / 39	514-100-00	80 m

remanium® laboratory coil, round, spring hard		
0.40 mm / 16	521-040-00	30 m
0.50 mm / 20	523-050-00	50 m
0.60 mm / 23	523-060-00	40 m
0.70 mm / 28	523-070-00	30 m
0.80 mm / 31	523-080-00	20 m
0.90 mm / 36	523-090-00	10 m
1.00 mm / 39	523-100-00	10 m
1.10 mm / 43	523-110-00	10 m
1.20 mm / 47	523-120-00	10 m
1.50 mm / 59	523-150-00	10 m

remanium® clinical coil, round, spring hard		
0.60 mm / 23	524-060-00	225 m
0.70 mm / 28	524-070-00	165 m
0.80 mm / 31	524-080-00	125 m
0.90 mm / 36	524-090-00	100 m
1.00 mm / 39	524-100-00	80 m
1.10 mm / 43	524-110-00	70 m

remanium® wire on coils, half-round, hard		
1.30 mm x 0.65 mm / 51 x 26	308-513-00	10 m
1.50 mm x 0.75 mm / 59 x 30	308-515-00	10 m
1.75 mm x 0.90 mm / 69 x 36	308-518-00	10 m
2.00 mm x 1.00 mm / 79 x 40	308-520-00	10 m

remanium® clinical coil, half-round, hard		
1.30 mm x 0.65 mm / 51 x 26	308-713-00	90 m
1.50 mm x 0.75 mm / 59 x 30	308-715-00	65 m
1.75 mm x 0.90 mm / 69 x 36	308-718-00	50 m

remanium® wire on coils, oval, hard		
2.40 mm x 1.40 mm / 95 x 55	385-624-00	10 m
3.00 mm x 1.50 mm / 118 x 59	385-630-00	10 m

Nickel-free* stainless steel

ø	REF	Quantity
Noninium® wire on coils, round, hard		
0.70 mm / 28	520-070-10	30 m
0.80 mm / 31	520-080-10	20 m
0.90 mm / 36	520-090-10	10 m
1.20 mm / 47	520-120-10	10 m

Noninium® wire on coils, round, spring hard		
0.60 mm / 23	520-062-10	40 m
0.70 mm / 28	520-072-10	30 m
0.80 mm / 31	520-082-10	20 m
0.90 mm / 36	520-092-10	10 m

* Nickel trace elements of ≤ 0.05 % are due to the process of manufacturing


Titanium

ø	REF	Quantity
Titanium retainer wire, 3-strand twisted, Ti1, round		
0.50 mm / 20	528-000-02	2 m



Fixed appliance technique / Arches

Nickel titanium



rematitan® LITE ideal arch, round


Strength category: super elastic

Color-coded packaging: orange

Mid-line marking: Maxillary 3 dots

ø	Length	Diagonal M.	REF	Quantity
0.30 mm / 12			766-058-00	10 pieces
0.35 mm / 14			766-080-00	10 pieces
0.35 mm / 14			766-080-01	50 pieces
0.40 mm / 16			766-082-00	10 pieces
0.40 mm / 16			766-082-01	50 pieces
0.45 mm / 18			766-084-00	10 pieces
0.45 mm / 18			766-084-01	50 pieces
0.50 mm / 20			768-100-00	10 pieces
0.50 mm / 20			768-100-01	50 pieces

Mid-line marking: Mandibular 1 dot



ø	Length	Diagonal M.	REF	Quantity
0.30 mm / 12			766-059-00	10 pieces
0.35 mm / 14			766-081-00	10 pieces
0.35 mm / 14			766-081-01	50 pieces
0.40 mm / 16			766-083-00	10 pieces
0.40 mm / 16			766-083-01	50 pieces
0.45 mm / 18			766-085-00	10 pieces
0.45 mm / 18			766-085-01	50 pieces
0.50 mm / 20			768-101-00	10 pieces
0.50 mm / 20			768-101-01	50 pieces

**rematitan® LITE ideal arch, rectangular**

Strength category: super elastic

Color-coded packaging: orange

Mid-line marking: Maxillary 3 dots

	Length	Diagonal M.	REF	Quantity
0.41 mm x 0.41 mm / 16 x 16		0.54 mm	766-088-00	10 pieces
0.41 mm x 0.41 mm / 16 x 16		0.54 mm	766-088-01	50 pieces
0.41 mm x 0.56 mm / 16 x 22		0.62 mm	766-090-00	10 pieces
0.41 mm x 0.56 mm / 16 x 22		0.62 mm	766-090-01	50 pieces
0.43 mm x 0.64 mm / 17 x 25		0.70 mm	766-092-00	10 pieces
0.46 mm x 0.46 mm / 18 x 18		0.59 mm	768-102-00	10 pieces
0.46 mm x 0.46 mm / 18 x 18		0.59 mm	768-102-01	50 pieces
0.46 mm x 0.64 mm / 18 x 25		0.71 mm	768-104-00	10 pieces
0.46 mm x 0.64 mm / 18 x 25		0.71 mm	768-104-01	50 pieces
0.48 mm x 0.64 mm / 19 x 25		0.73 mm	768-106-00	10 pieces
0.48 mm x 0.64 mm / 19 x 25		0.73 mm	768-106-01	50 pieces
0.51 mm x 0.51 mm / 20 x 20		0.64 mm	768-108-00	10 pieces
0.51 mm x 0.51 mm / 20 x 20		0.64 mm	768-108-01	50 pieces
0.53 mm x 0.64 mm / 21 x 25		0.76 mm	768-110-00	10 pieces
0.53 mm x 0.64 mm / 21 x 25		0.76 mm	768-110-01	50 pieces


Mid-line marking: Mandibular 1 dot



	Length	Diagonal M.	REF	Quantity
0.41 mm x 0.41 mm / 16 x 16		0.54 mm	766-089-00	10 pieces
0.41 mm x 0.41 mm / 16 x 16		0.54 mm	766-089-01	50 pieces
0.41 mm x 0.56 mm / 16 x 22		0.62 mm	766-091-00	10 pieces
0.41 mm x 0.56 mm / 16 x 22		0.62 mm	766-091-01	50 pieces
0.43 mm x 0.64 mm / 17 x 25		0.70 mm	766-093-00	10 pieces
0.46 mm x 0.46 mm / 18 x 18		0.59 mm	768-103-00	10 pieces
0.46 mm x 0.46 mm / 18 x 18		0.59 mm	768-103-01	50 pieces
0.46 mm x 0.64 mm / 18 x 25		0.71 mm	768-105-00	10 pieces
0.46 mm x 0.64 mm / 18 x 25		0.71 mm	768-105-01	50 pieces
0.48 mm x 0.64 mm / 19 x 25		0.73 mm	768-107-00	10 pieces
0.48 mm x 0.64 mm / 19 x 25		0.73 mm	768-107-01	50 pieces
0.51 mm x 0.51 mm / 20 x 20		0.64 mm	768-109-00	10 pieces
0.51 mm x 0.51 mm / 20 x 20		0.64 mm	768-109-01	50 pieces
0.53 mm x 0.64 mm / 21 x 25		0.76 mm	768-111-00	10 pieces
0.53 mm x 0.64 mm / 21 x 25		0.76 mm	768-111-01	50 pieces

Fixed appliance technique / Arches

Nickel titanium




rematitan® LITE ideal arch, round

Strength category: super elastic

Color-coded packaging: orange

Mid-line marking: with dimple

ø	Length	Diagonal M.	REF	Quantity
0.30 mm / 12			766-060-00	10 pieces
0.35 mm / 14			766-062-00	10 pieces
0.40 mm / 16			766-064-00	10 pieces
0.45 mm / 18			766-066-00	10 pieces
0.50 mm / 20			768-112-00	10 pieces



ø	Length	Diagonal M.	REF	Quantity
0.30 mm / 12			766-061-00	10 pieces
0.35 mm / 14			766-063-00	10 pieces
0.40 mm / 16			766-065-00	10 pieces
0.45 mm / 18			766-067-00	10 pieces
0.50 mm / 20			768-113-00	10 pieces




rematitan® LITE ideal arch, rectangular

Strength category: super elastic


Color-coded packaging: orange

Mid-line marking: with dimple

∅	Length	Diagonal M.	REF	Quantity
0.41 mm x 0.41 mm / 16 x 16		0.54 mm	766-068-00	10 pieces
0.41 mm x 0.56 mm / 16 x 22		0.62 mm	766-070-00	10 pieces
0.43 mm x 0.64 mm / 17 x 25		0.70 mm	766-072-00	10 pieces
0.46 mm x 0.46 mm / 18 x 18		0.59 mm	766-074-00	10 pieces
0.46 mm x 0.64 mm / 18 x 25		0.71 mm	766-076-00	10 pieces
0.48 mm x 0.64 mm / 19 x 25		0.73 mm	766-078-00	10 pieces
0.51 mm x 0.51 mm / 20 x 20		0.64 mm	768-114-00	10 pieces
0.53 mm x 0.64 mm / 21 x 25		0.76 mm	766-086-00	10 pieces




∅	Length	Diagonal M.	REF	Quantity
0.41 mm x 0.41 mm / 16 x 16		0.54 mm	766-069-00	10 pieces
0.41 mm x 0.56 mm / 16 x 22		0.62 mm	766-071-00	10 pieces
0.43 mm x 0.64 mm / 17 x 25		0.70 mm	766-073-00	10 pieces
0.46 mm x 0.46 mm / 18 x 18		0.59 mm	766-075-00	10 pieces
0.46 mm x 0.64 mm / 18 x 25		0.71 mm	766-077-00	10 pieces
0.48 mm x 0.64 mm / 19 x 25		0.73 mm	766-079-00	10 pieces
0.51 mm x 0.51 mm / 20 x 20		0.64 mm	768-115-00	10 pieces
0.53 mm x 0.64 mm / 21 x 25		0.76 mm	766-087-00	10 pieces

**rematitan® LITE ideal arch, round****For narrow jaws. Narrow arch form.**

Strength category: super elastic

Color-coded packaging: orange

Mid-line marking: with dimple



ø	Length	Diagonal M.	REF	Quantity
0.30 mm / 12			766-160-00	10 pieces
0.35 mm / 14			766-162-00	10 pieces
0.40 mm / 16			766-164-00	10 pieces
0.45 mm / 18			766-166-00	10 pieces
0.50 mm / 20			768-116-00	10 pieces

ø	Length	Diagonal M.	REF	Quantity
0.30 mm / 12			766-161-00	10 pieces
0.35 mm / 14			766-163-00	10 pieces
0.40 mm / 16			766-165-00	10 pieces
0.45 mm / 18			766-167-00	10 pieces
0.50 mm / 20			768-117-00	10 pieces

**rematitan® LITE ideal arch, rectangular****For narrow jaws. Narrow arch form.**

Strength category: super elastic

Color-coded packaging: orange

Mid-line marking: with dimple



▣	Length	Diagonal M.	REF	Quantity
0.41 mm x 0.41 mm / 16 x 16		0.54 mm	766-168-00	10 pieces
0.41 mm x 0.56 mm / 16 x 22		0.62 mm	766-170-00	10 pieces
0.43 mm x 0.64 mm / 17 x 25		0.70 mm	766-172-00	10 pieces
0.46 mm x 0.46 mm / 18 x 18		0.59 mm	766-174-00	10 pieces
0.46 mm x 0.64 mm / 18 x 25		0.71 mm	766-176-00	10 pieces
0.48 mm x 0.64 mm / 19 x 25		0.73 mm	766-178-00	10 pieces
0.51 mm x 0.51 mm / 20 x 20		0.64 mm	768-118-00	10 pieces
0.53 mm x 0.64 mm / 21 x 25		0.76 mm	768-120-00	10 pieces

▣	Length	Diagonal M.	REF	Quantity
0.41 mm x 0.41 mm / 16 x 16		0.54 mm	766-169-00	10 pieces
0.41 mm x 0.56 mm / 16 x 22		0.62 mm	766-171-00	10 pieces
0.43 mm x 0.64 mm / 17 x 25		0.70 mm	766-173-00	10 pieces
0.46 mm x 0.46 mm / 18 x 18		0.59 mm	766-175-00	10 pieces
0.46 mm x 0.64 mm / 18 x 25		0.71 mm	766-177-00	10 pieces
0.48 mm x 0.64 mm / 19 x 25		0.73 mm	766-179-00	10 pieces
0.51 mm x 0.51 mm / 20 x 20		0.64 mm	768-119-00	10 pieces
0.53 mm x 0.64 mm / 21 x 25		0.76 mm	768-121-00	10 pieces

Fixed appliance technique / Arches

Nickel titanium




rematitan® LITE Spee® retraction arch, round

Strength category: super elastic

Mid-line marking: with dimple

ø	Length	Diagonal M.	REF	Quantity
0.30 mm / 12			766-014-00	10 pieces
0.35 mm / 14			766-018-00	10 pieces
0.40 mm / 16			766-050-00	10 pieces
0.45 mm / 18			766-052-00	10 pieces



ø	Length	Diagonal M.	REF	Quantity
0.30 mm / 12			766-015-00	10 pieces
0.35 mm / 14			766-019-00	10 pieces
0.40 mm / 16			766-051-00	10 pieces
0.45 mm / 18			766-053-00	10 pieces



rematitan® LITE Spee® retraction arch, rectangular

Strength category: super elastic

Mid-line marking: with dimple

∅	Length	Diagonal M.	REF	Quantity
0.41 mm x 0.56 mm / 16 x 22		0.62 mm	766-054-00	10 pieces
0.43 mm x 0.64 mm / 17 x 25		0.70 mm	766-056-00	10 pieces
0.46 mm x 0.46 mm / 18 x 18		0.59 mm	766-032-00	10 pieces
0.46 mm x 0.64 mm / 18 x 25		0.71 mm	766-034-00	10 pieces
0.48 mm x 0.64 mm / 19 x 25		0.73 mm	766-036-00	10 pieces



∅	Length	Diagonal M.	REF	Quantity
0.41 mm x 0.56 mm / 16 x 22		0.62 mm	766-055-00	10 pieces
0.43 mm x 0.64 mm / 17 x 25		0.70 mm	766-057-00	10 pieces
0.46 mm x 0.46 mm / 18 x 18		0.59 mm	766-033-00	10 pieces
0.46 mm x 0.64 mm / 18 x 25		0.71 mm	766-037-00	10 pieces
0.48 mm x 0.64 mm / 19 x 25		0.73 mm	766-039-00	10 pieces


rematitan® LITE arch, arch form I, pointed
For McLaughlin-Bennett-Trevisi technique**

Strength category: super elastic

Color-coded packaging: orange

Mid-line marking: Maxillary / Mandibular 1 dot, etched

▣	Length	Diagonal M.	REF	Quantity
0.41 mm x 0.41 mm / 16 x 16		0.54 mm	766-180-00	10 pieces
0.41 mm x 0.56 mm / 16 x 22		0.62 mm	766-181-00	10 pieces
0.43 mm x 0.64 mm / 17 x 25		0.70 mm	766-182-00	10 pieces

** The Dentaaurum version of this prescription is not claimed to be a duplication of any other, nor does Dentaaurum imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.


rematitan® LITE arch, arch form II, square
For McLaughlin-Bennett-Trevisi technique**

Strength category: super elastic

Color-coded packaging: orange

Mid-line marking: Maxillary / Mandibular 1 dot, etched

▣	Length	Diagonal M.	REF	Quantity
0.41 mm x 0.41 mm / 16 x 16		0.54 mm	766-183-00	10 pieces
0.41 mm x 0.56 mm / 16 x 22		0.62 mm	766-184-00	10 pieces
0.43 mm x 0.64 mm / 17 x 25		0.70 mm	766-185-00	10 pieces

** The Dentaaurum version of this prescription is not claimed to be a duplication of any other, nor does Dentaaurum imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.


rematitan® LITE arch, arch form III, oval
For McLaughlin-Bennett-Trevisi technique**

Strength category: super elastic

Color-coded packaging: orange

Mid-line marking: Maxillary / Mandibular 1 dot, etched

▣	Length	Diagonal M.	REF	Quantity
0.41 mm x 0.41 mm / 16 x 16		0.54 mm	766-186-00	10 pieces
0.41 mm x 0.56 mm / 16 x 22		0.62 mm	766-187-00	10 pieces
0.43 mm x 0.64 mm / 17 x 25		0.70 mm	766-188-00	10 pieces

** The Dentaaurum version of this prescription is not claimed to be a duplication of any other, nor does Dentaaurum imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.


rematitan® LITE Ricketts® PENTA arch, normal

Strength category: super elastic

Color-coded packaging: orange

Mid-line marking: Maxillary / Mandibular 1 dot, etched

▣	Length	Diagonal M.	REF	Quantity
0.41 mm x 0.41 mm / 16 x 16		0.54 mm	766-250-00	10 pieces
0.41 mm x 0.56 mm / 16 x 22		0.62 mm	766-251-00	10 pieces
0.43 mm x 0.64 mm / 17 x 25		0.70 mm	766-252-00	10 pieces

Fixed appliance technique / Arches

Nickel titanium



rematitan® LITE Ricketts® PENTA arch, pointed

Strength category: super elastic

Color-coded packaging: orange

Mid-line marking: Maxillary / Mandibular 1 dot, etched

▣	Length	Diagonal M.	REF	Quantity
0.41 mm x 0.41 mm / 16 x 16		0.54 mm	766-253-00	10 pieces
0.41 mm x 0.56 mm / 16 x 22		0.62 mm	766-254-00	10 pieces
0.43 mm x 0.64 mm / 17 x 25		0.70 mm	766-255-00	10 pieces



rematitan® LITE Ricketts® PENTA arch, oval

Strength category: super elastic

Color-coded packaging: orange

Mid-line marking: Maxillary / Mandibular 1 dot, etched

▣	Length	Diagonal M.	REF	Quantity
0.41 mm x 0.41 mm / 16 x 16		0.54 mm	766-256-00	10 pieces
0.41 mm x 0.56 mm / 16 x 22		0.62 mm	766-257-00	10 pieces
0.43 mm x 0.64 mm / 17 x 25		0.70 mm	766-258-00	10 pieces



rematitan® LITE Ricketts® PENTA arch, extra pointed

Strength category: super elastic

Color-coded packaging: orange

Mid-line marking: Maxillary / Mandibular 1 dot, etched

▣	Length	Diagonal M.	REF	Quantity
0.41 mm x 0.41 mm / 16 x 16		0.54 mm	766-259-00	10 pieces
0.41 mm x 0.56 mm / 16 x 22		0.62 mm	766-260-00	10 pieces
0.43 mm x 0.64 mm / 17 x 25		0.70 mm	766-261-00	10 pieces




rematitan® LITE Ricketts® PENTA arch, oval narrow

Strength category: super elastic

Color-coded packaging: orange

Mid-line marking: Maxillary / Mandibular 1 dot, etched

▣	Length	Diagonal M.	REF	Quantity
0.41 mm x 0.41 mm / 16 x 16		0.54 mm	766-262-00	10 pieces
0.41 mm x 0.56 mm / 16 x 22		0.62 mm	766-263-00	10 pieces
0.43 mm x 0.64 mm / 17 x 25		0.70 mm	766-264-00	10 pieces

**Tensic® ideal arch, round**

Strength category: super elastic


Color-coded packaging: orange

Mid-line marking: Maxillary 3 dots

ø	Length	Diagonal M.	REF	Quantity
0.30 mm / 12			766-700-00	10 pieces
0.35 mm / 14			766-702-00	10 pieces
0.40 mm / 16			766-704-00	10 pieces
0.45 mm / 18			766-706-00	10 pieces

Mid-line marking: Mandibular 1 dot

ø	Length	Diagonal M.	REF	Quantity
0.30 mm / 12			766-701-00	10 pieces
0.35 mm / 14			766-703-00	10 pieces
0.40 mm / 16			766-705-00	10 pieces
0.45 mm / 18			766-707-00	10 pieces

**Tensic® ideal arch, rectangular**

Strength category: super elastic

Color-coded packaging: orange

Mid-line marking: Maxillary 3 dots

▣	Length	Diagonal M.	REF	Quantity
0.41 mm x 0.41 mm / 16 x 16		0.54 mm	766-710-00	10 pieces
0.41 mm x 0.56 mm / 16 x 22		0.62 mm	766-712-00	10 pieces
0.43 mm x 0.64 mm / 17 x 25		0.70 mm	766-714-00	10 pieces
0.46 mm x 0.64 mm / 18 x 25		0.71 mm	766-716-00	10 pieces
0.48 mm x 0.64 mm / 19 x 25		0.73 mm	766-718-00	10 pieces
0.53 mm x 0.64 mm / 21 x 25		0.76 mm	766-720-00	10 pieces


Mid-line marking: Mandibular 1 dot

▣	Length	Diagonal M.	REF	Quantity
0.41 mm x 0.41 mm / 16 x 16		0.54 mm	766-711-00	10 pieces
0.41 mm x 0.56 mm / 16 x 22		0.62 mm	766-713-00	10 pieces
0.43 mm x 0.64 mm / 17 x 25		0.70 mm	766-715-00	10 pieces
0.46 mm x 0.64 mm / 18 x 25		0.71 mm	766-717-00	10 pieces
0.48 mm x 0.64 mm / 19 x 25		0.73 mm	766-719-00	10 pieces
0.53 mm x 0.64 mm / 21 x 25		0.76 mm	766-721-00	10 pieces



Fixed appliance technique / Arches

Thermo-active nickel titanium



Equire thermo-active, preformed ideal arch, round

Arch form: American Style $\hat{=}$ Damon Style*

Strength category: super elastic


Mid-line marking: Maxillary 3 dots, etched

Ø	Length	Diagonal M.	REF	Quantity
0.30 mm / 12			766-800-00	10 pieces
0.35 mm / 14			766-802-00	10 pieces
0.40 mm / 16			766-804-00	10 pieces
0.45 mm / 18			766-806-00	10 pieces

Mid-line marking: Mandibular 1 dot, etched

Ø	Length	Diagonal M.	REF	Quantity
0.30 mm / 12			766-801-00	10 pieces
0.35 mm / 14			766-803-00	10 pieces
0.40 mm / 16			766-805-00	10 pieces
0.45 mm / 18			766-807-00	10 pieces

*Damon® is a registered trademark of Ormco BV



Equire thermo-active, preformed ideal arch, rectangular

Arch form: American Style $\hat{=}$ Damon Style*

Strength category: super elastic

Mid-line marking: Maxillary 3 dots, etched

∅	Length	Diagonal M.	REF	Quantity
0.41 mm x 0.41 mm / 16 x 16		0.54 mm	766-808-00	10 pieces
0.41 mm x 0.56 mm / 16 x 22		0.62 mm	766-810-00	10 pieces
0.43 mm x 0.64 mm / 17 x 25		0.70 mm	766-812-00	10 pieces
0.46 mm x 0.64 mm / 18 x 25		0.71 mm	766-814-00	10 pieces
0.48 mm x 0.64 mm / 19 x 25		0.73 mm	766-816-00	10 pieces
0.53 mm x 0.64 mm / 21 x 25		0.76 mm	766-818-00	10 pieces

Mid-line marking: Mandibular 1 dot, etched

∅	Length	Diagonal M.	REF	Quantity
0.41 mm x 0.41 mm / 16 x 16		0.54 mm	766-809-00	10 pieces
0.41 mm x 0.56 mm / 16 x 22		0.62 mm	766-811-00	10 pieces
0.43 mm x 0.64 mm / 17 x 25		0.70 mm	766-813-00	10 pieces
0.46 mm x 0.64 mm / 18 x 25		0.71 mm	766-815-00	10 pieces
0.48 mm x 0.64 mm / 19 x 25		0.73 mm	766-817-00	10 pieces
0.53 mm x 0.64 mm / 21 x 25		0.76 mm	766-819-00	10 pieces

*Damon® is a registered trademark of Ormco BV



Stainless steel (twisted and braided)


dentaflex® ideal arch, round, 3-strand twisted

Strength category: super spring hard · 2500 N/mm² - 2700 N/mm²

Color-coded packaging: red

Mid-line marking: Maxillary 3 dots

ø	Length	Diagonal M.	REF	Quantity
0.38 mm / 15			766-120-00	10 pieces
0.38 mm / 15			766-120-01	50 pieces
0.45 mm / 18			766-122-00	10 pieces
0.45 mm / 18			766-122-01	50 pieces
0.50 mm / 20			766-124-00	10 pieces

Mid-line marking: Mandibular 1 dot

ø	Length	Diagonal M.	REF	Quantity
0.38 mm / 15			766-121-00	10 pieces
0.38 mm / 15			766-121-01	50 pieces
0.45 mm / 18			766-123-00	10 pieces
0.45 mm / 18			766-123-01	50 pieces
0.50 mm / 20			766-125-00	10 pieces


dentaflex® ideal arch, rectangular, 3-strand twisted

Strength category: super spring hard · 2500 N/mm² - 2700 N/mm²

Color-coded packaging: red

Mid-line marking: Maxillary 3 dots

∅	Length	Diagonal M.	REF	Quantity
0.41 mm x 0.41 mm / 16 x 16		0.54 mm	766-130-00	10 pieces
0.41 mm x 0.56 mm / 16 x 22		0.62 mm	766-132-00	10 pieces
0.41 mm x 0.56 mm / 16 x 22		0.62 mm	766-132-01	50 pieces

Mid-line marking: Mandibular 1 dot

∅	Length	Diagonal M.	REF	Quantity
0.41 mm x 0.41 mm / 16 x 16		0.54 mm	766-131-00	10 pieces
0.41 mm x 0.56 mm / 16 x 22		0.62 mm	766-133-00	10 pieces
0.41 mm x 0.56 mm / 16 x 22		0.62 mm	766-133-01	50 pieces



Fixed appliance technique / Arches

Stainless steel (twisted and braided)



dentaflex® ideal arch, rectangular, 8-strand braided

Strength category: super spring hard · 2500 N/mm² - 2700 N/mm²


Color-coded packaging: red

Mid-line marking: Maxillary 3 dots

Ø	Length	Diagonal M.	REF	Quantity
0.41 mm x 0.41 mm / 16 x 16		0.54 mm	766-100-00	10 pieces
0.41 mm x 0.41 mm / 16 x 16		0.54 mm	766-100-01	50 pieces
0.41 mm x 0.56 mm / 16 x 22		0.62 mm	766-101-00	10 pieces
0.41 mm x 0.56 mm / 16 x 22		0.62 mm	766-101-01	50 pieces
0.43 mm x 0.64 mm / 17 x 25		0.70 mm	766-103-00	10 pieces
0.46 mm x 0.64 mm / 18 x 25		0.71 mm	766-104-00	10 pieces
0.48 mm x 0.64 mm / 19 x 25		0.73 mm	766-105-00	10 pieces

Mid-line marking: Mandibular 1 dot

Ø	Length	Diagonal M.	REF	Quantity
0.41 mm x 0.41 mm / 16 x 16		0.54 mm	766-200-00	10 pieces
0.41 mm x 0.41 mm / 16 x 16		0.54 mm	766-200-01	50 pieces
0.41 mm x 0.56 mm / 16 x 22		0.62 mm	766-201-00	10 pieces
0.41 mm x 0.56 mm / 16 x 22		0.62 mm	766-201-01	50 pieces
0.43 mm x 0.64 mm / 17 x 25		0.70 mm	766-203-00	10 pieces
0.46 mm x 0.64 mm / 18 x 25		0.71 mm	766-204-00	10 pieces
0.48 mm x 0.64 mm / 19 x 25		0.73 mm	766-205-00	10 pieces



dentaflex® ideal arch, round, 6-strand co-axial

Strength category: hard plus · 1600 N/mm² - 1800 N/mm²


Color-coded packaging: red

Mid-line marking: Maxillary 3 dots


Ø	Length	Diagonal M.	REF	Quantity
0.38 mm / 15			766-300-00	10 pieces
0.38 mm / 15			766-300-01	50 pieces
0.45 mm / 18			766-301-00	10 pieces
0.45 mm / 18			766-301-01	50 pieces

Mid-line marking: Mandibular 1 dot

Ø	Length	Diagonal M.	REF	Quantity
0.38 mm / 15			766-400-00	10 pieces
0.38 mm / 15			766-400-01	50 pieces
0.45 mm / 18			766-401-00	10 pieces
0.45 mm / 18			766-401-01	50 pieces



Stainless steel (nickel chromium)


remanium® ideal arch, round


Strength category: spring hard plus · 2000 N/mm² - 2300 N/mm²

Color-coded packaging: red

Mid-line marking: Maxillary 3 dots

ø	Length	Diagonal M.	REF	Quantity
0.35 mm / 14			765-200-00	10 pieces
0.40 mm / 16			765-201-00	10 pieces
0.40 mm / 16			765-201-01	50 pieces
0.45 mm / 18			765-202-00	10 pieces
0.45 mm / 18			765-202-01	50 pieces
0.50 mm / 20			765-203-00	10 pieces

Mid-line marking: Mandibular 1 dot



ø	Length	Diagonal M.	REF	Quantity
0.35 mm / 14			765-300-00	10 pieces
0.40 mm / 16			765-301-00	10 pieces
0.40 mm / 16			765-301-01	50 pieces
0.45 mm / 18			765-302-00	10 pieces
0.45 mm / 18			765-302-01	50 pieces
0.50 mm / 20			765-303-00	10 pieces



Fixed appliance technique / Arches

Stainless steel (nickel chromium)



remanium® ideal arch, rectangular

Strength category: spring hard plus · 2000 N/mm² - 2300 N/mm²

Color-coded packaging: red


Mid-line marking: Maxillary 3 dots

∅	Length	Diagonal M.	REF	Quantity
0.41 mm x 0.41 mm / 16 x 16		0.54 mm	767-100-00	10 pieces
0.41 mm x 0.41 mm / 16 x 16		0.54 mm	767-100-01	50 pieces
0.41 mm x 0.56 mm / 16 x 22		0.62 mm	767-101-00	10 pieces
0.41 mm x 0.56 mm / 16 x 22		0.62 mm	767-101-01	50 pieces
0.43 mm x 0.56 mm / 17 x 22		0.64 mm	767-102-00	10 pieces
0.43 mm x 0.64 mm / 17 x 25		0.70 mm	767-103-00	10 pieces
0.43 mm x 0.64 mm / 17 x 25		0.70 mm	767-103-01	50 pieces
0.46 mm x 0.64 mm / 18 x 25		0.71 mm	767-104-00	10 pieces
0.48 mm x 0.64 mm / 19 x 25		0.73 mm	767-105-00	10 pieces
0.48 mm x 0.64 mm / 19 x 25		0.73 mm	767-105-01	50 pieces
0.55 mm x 0.64 mm / 21.5 x 25		0.76 mm	767-106-00	10 pieces
0.55 mm x 0.70 mm / 21.5 x 28		0.81 mm	767-107-00	10 pieces

Mid-line marking: Mandibular 1 dot



∅	Length	Diagonal M.	REF	Quantity
0.41 mm x 0.41 mm / 16 x 16		0.54 mm	767-200-00	10 pieces
0.41 mm x 0.41 mm / 16 x 16		0.54 mm	767-200-01	50 pieces
0.41 mm x 0.56 mm / 16 x 22		0.62 mm	767-201-00	10 pieces
0.41 mm x 0.56 mm / 16 x 22		0.62 mm	767-201-01	50 pieces
0.43 mm x 0.56 mm / 17 x 22		0.64 mm	767-202-00	10 pieces
0.43 mm x 0.64 mm / 17 x 25		0.70 mm	767-203-00	10 pieces
0.43 mm x 0.64 mm / 17 x 25		0.70 mm	767-203-01	50 pieces
0.46 mm x 0.64 mm / 18 x 25		0.71 mm	767-204-00	10 pieces
0.48 mm x 0.64 mm / 19 x 25		0.73 mm	767-205-00	10 pieces
0.48 mm x 0.64 mm / 19 x 25		0.73 mm	767-205-01	50 pieces
0.55 mm x 0.64 mm / 21.5 x 25		0.76 mm	767-206-00	10 pieces
0.55 mm x 0.70 mm / 21.5 x 28		0.81 mm	767-207-00	10 pieces



Equire preformed ideal arch, round


Arch form: American Style ≙ Damon Style*

Strength category: spring hard plus · 2000 N/mm² - 2300 N/mm²

Mid-line marking: Maxillary 3 dots, etched

∅	Length	Diagonal M.	REF	Quantity
0.30 mm / 12			769-000-00	10 pieces
0.35 mm / 14			769-002-00	10 pieces
0.40 mm / 16			769-004-00	10 pieces
0.45 mm / 18			769-006-00	10 pieces

Mid-line marking: Mandibular 1 dot, etched



∅	Length	Diagonal M.	REF	Quantity
0.30 mm / 12			769-001-00	10 pieces
0.35 mm / 14			769-003-00	10 pieces
0.40 mm / 16			769-005-00	10 pieces
0.45 mm / 18			769-007-00	10 pieces

*Damon® is a registered trademark of Ormco BV

Stainless steel (nickel chromium)


Equire preformed ideal arch, rectangular
Arch form: American Style $\hat{=}$ Damon Style*

Strength category: spring hard plus · 2000 N/mm² - 2300 N/mm²

Mid-line marking: Maxillary 3 dots, etched

Length	Diagonal M.	REF	Quantity
0.41 mm x 0.41 mm / 16 x 16	0.54 mm	769-008-00	10 pieces
0.41 mm x 0.56 mm / 16 x 22	0.62 mm	769-010-00	10 pieces
0.43 mm x 0.64 mm / 17 x 25	0.70 mm	769-012-00	10 pieces
0.46 mm x 0.64 mm / 18 x 25	0.71 mm	769-014-00	10 pieces
0.48 mm x 0.64 mm / 19 x 25	0.73 mm	769-016-00	10 pieces
0.53 mm x 0.64 mm / 21 x 25	0.76 mm	769-018-00	10 pieces

Mid-line marking: Mandibular 1 dot, etched



Length	Diagonal M.	REF	Quantity
0.41 mm x 0.41 mm / 16 x 16	0.54 mm	769-009-00	10 pieces
0.41 mm x 0.56 mm / 16 x 22	0.62 mm	769-011-00	10 pieces
0.43 mm x 0.64 mm / 17 x 25	0.70 mm	769-013-00	10 pieces
0.46 mm x 0.64 mm / 18 x 25	0.71 mm	769-015-00	10 pieces
0.48 mm x 0.64 mm / 19 x 25	0.73 mm	769-017-00	10 pieces
0.53 mm x 0.64 mm / 21 x 25	0.76 mm	769-019-00	10 pieces

*Damon® is a registered trademark of Ormco BV

CoCr


remaloy® ideal arch, rectangular

Strength category: hard plus · 1600 N/mm² - 1800 N/mm²

Color-coded packaging: blue

Mid-line marking: Maxillary 3 dots

Length	Diagonal M.	REF	Quantity
0.41 mm x 0.41 mm / 16 x 16	0.54 mm	769-100-00	10 pieces
0.41 mm x 0.56 mm / 16 x 22	0.62 mm	769-101-00	10 pieces


Mid-line marking: Mandibular 1 dot



Length	Diagonal M.	REF	Quantity
0.41 mm x 0.41 mm / 16 x 16	0.54 mm	769-200-00	10 pieces
0.41 mm x 0.56 mm / 16 x 22	0.62 mm	769-201-00	10 pieces

Fixed appliance technique / Arches

Nickel-free titanium molybdenum



rematitan® SPECIAL ideal arch, round

Strength category: middle plus - 1100 N/mm² - 1400 N/mm²


Color-coded packaging: yellow

Mid-line marking: Maxillary 3 dots

ø	Length	Diagonal M.	REF	Quantity
0.40 mm / 16			766-605-00	10 pieces
0.45 mm / 18			766-607-00	10 pieces

Mid-line marking: Mandibular 1 dot

ø	Length	Diagonal M.	REF	Quantity
0.40 mm / 16			766-606-00	10 pieces
0.45 mm / 18			766-608-00	10 pieces



rematitan® SPECIAL ideal arch, rectangular

Strength category: middle plus - 1100 N/mm² - 1400 N/mm²

Color-coded packaging: yellow



Mid-line marking: Maxillary 3 dots

∅	Length	Diagonal M.	REF	Quantity
0.41 mm x 0.41 mm / 16 x 16		0.54 mm	766-631-00	10 pieces
0.41 mm x 0.56 mm / 16 x 22		0.62 mm	766-609-00	10 pieces
0.43 mm x 0.64 mm / 17 x 25		0.70 mm	766-611-00	10 pieces
0.46 mm x 0.46 mm / 18 x 18		0.59 mm	766-613-00	10 pieces
0.46 mm x 0.64 mm / 18 x 25		0.71 mm	766-615-00	10 pieces
0.48 mm x 0.64 mm / 19 x 25		0.73 mm	766-617-00	10 pieces
0.53 mm x 0.64 mm / 21 x 25		0.76 mm	766-619-00	10 pieces

Mid-line marking: Mandibular 1 dot

∅	Length	Diagonal M.	REF	Quantity
0.41 mm x 0.41 mm / 16 x 16		0.54 mm	766-632-00	10 pieces
0.41 mm x 0.56 mm / 16 x 22		0.62 mm	766-610-00	10 pieces
0.43 mm x 0.64 mm / 17 x 25		0.70 mm	766-612-00	10 pieces
0.46 mm x 0.46 mm / 18 x 18		0.59 mm	766-614-00	10 pieces
0.46 mm x 0.64 mm / 18 x 25		0.71 mm	766-616-00	10 pieces
0.48 mm x 0.64 mm / 19 x 25		0.73 mm	766-618-00	10 pieces
0.53 mm x 0.64 mm / 21 x 25		0.76 mm	766-620-00	10 pieces




rematitan® LITE White ideal arch, round


Strength category: super elastic

Color-coded packaging: orange

Mid-line marking: none

ø	Length	Diagonal M.	REF	Quantity
0.40 mm / 16			768-014-00	3 pieces
0.45 mm / 18			768-016-00	3 pieces

ø	Length	Diagonal M.	REF	Quantity
0.40 mm / 16			768-015-00	3 pieces
0.45 mm / 18			768-017-00	3 pieces


rematitan® LITE White ideal arch, rectangular

Strength category: super elastic

Color-coded packaging: orange

Mid-line marking: none


▣	Length	Diagonal M.	REF	Quantity
0.41 mm x 0.41 mm / 16 x 16		0.54 mm	768-018-00	3 pieces
0.41 mm x 0.56 mm / 16 x 22		0.62 mm	768-020-00	3 pieces
0.43 mm x 0.64 mm / 17 x 25		0.70 mm	768-022-00	3 pieces

▣	Length	Diagonal M.	REF	Quantity
0.41 mm x 0.41 mm / 16 x 16		0.54 mm	768-019-00	3 pieces
0.41 mm x 0.56 mm / 16 x 22		0.62 mm	768-021-00	3 pieces
0.43 mm x 0.64 mm / 17 x 25		0.70 mm	768-023-00	3 pieces



Fixed appliance technique / Arches

Thermo-active nickel titanium (white coated)




Tensic® White ideal arch, round

Strength category: super elastic


Color-coded packaging: orange

Mid-line marking: none

Ø	Length	Diagonal M.	REF	Quantity
0.40 mm / 16			768-026-00	3 pieces
0.45 mm / 18			768-028-00	3 pieces



Ø	Length	Diagonal M.	REF	Quantity
0.40 mm / 16			768-027-00	3 pieces
0.45 mm / 18			768-029-00	3 pieces



Tensic® White ideal arch, rectangular

Strength category: super elastic

Color-coded packaging: orange

Mid-line marking: none

▣	Length	Diagonal M.	REF	Quantity
0.41 mm x 0.41 mm / 16 x 16		0.54 mm	768-030-00	3 pieces
0.41 mm x 0.56 mm / 16 x 22		0.62 mm	768-032-00	3 pieces
0.43 mm x 0.64 mm / 17 x 25		0.70 mm	768-034-00	3 pieces
0.46 mm x 0.64 mm / 18 x 25		0.71 mm	768-038-00	3 pieces
0.48 mm x 0.64 mm / 19 x 25		0.73 mm	768-040-00	3 pieces



▣	Length	Diagonal M.	REF	Quantity
0.41 mm x 0.41 mm / 16 x 16		0.54 mm	768-031-00	3 pieces
0.41 mm x 0.56 mm / 16 x 22		0.62 mm	768-033-00	3 pieces
0.43 mm x 0.64 mm / 17 x 25		0.70 mm	768-035-00	3 pieces
0.46 mm x 0.64 mm / 18 x 25		0.71 mm	768-039-00	3 pieces
0.48 mm x 0.64 mm / 19 x 25		0.73 mm	768-041-00	3 pieces

Fixed appliance technique / Straight wires

Stainless steel (twisted and braided)

dentaflex® straight wire, round, 3-strand twisted

Strength category: super spring hard · 2500 N/mm² - 2700 N/mm²

Color-coded packaging: red

ø	Length	Diagonal M.	REF	Quantity
0.38 mm / 15	each 380 mm		545-738-00	10 pieces
0.45 mm / 18	each 380 mm		545-745-00	10 pieces
0.50 mm / 20	each 380 mm		545-750-00	10 pieces

dentaflex® straight wire, round, 6-strand co-axial

Strength category: hard plus · 1600 N/mm² - 1800 N/mm²

Color-coded packaging: red

ø	Length	Diagonal M.	REF	Quantity
0.38 mm / 15	each 380 mm		545-638-00	10 pieces
0.45 mm / 18	each 380 mm		545-645-00	10 pieces

Stainless steel (nickel chromium)

remanium® straight wire, round

Strength category: extra spring hard · 2300 N/mm² - 2500 N/mm²

Color-coded packaging: red

ø	Length	Diagonal M.	REF	Quantity
0.35 mm / 14	each 380 mm		535-035-00	25 pieces
0.40 mm / 16	each 380 mm		535-040-00	25 pieces
0.45 mm / 18	each 380 mm		535-045-00	25 pieces

remanium® straight wire, rectangular

Strength category: spring hard · 1800 N/mm² - 2000 N/mm²

Color-coded packaging: red

▣	Length	Diagonal M.	REF	Quantity
0.41 mm x 0.41 mm / 16 x 16	each 380 mm	0.54 mm	537-340-00	10 pieces
0.41 mm x 0.56 mm / 16 x 22	each 380 mm	0.62 mm	537-310-00	10 pieces
0.43 mm x 0.56 mm / 17 x 22	each 380 mm	0.64 mm	537-311-00	10 pieces
0.43 mm x 0.64 mm / 17 x 25	each 380 mm	0.70 mm	537-313-00	10 pieces
0.46 mm x 0.64 mm / 18 x 25	each 380 mm	0.71 mm	537-314-00	10 pieces
0.48 mm x 0.48 mm / 19 x 19	each 380 mm	0.62 mm	535-048-00	25 pieces
0.51 mm x 0.51 mm / 20 x 20	each 380 mm	0.64 mm	535-051-00	25 pieces
0.55 mm x 0.70 mm / 21.5 x 28	each 380 mm	0.81 mm	537-324-00	10 pieces

Fixed appliance technique / Straight wires

Stainless steel (nickel chromium)

remanium® straight wire, rectangular

Strength category: spring hard plus · 2000 N/mm² - 2300 N/mm²

Color-coded packaging: red

▣	Length	Diagonal M.	REF	Quantity
0.41 mm x 0.41 mm / 16 x 16	each 380 mm	0.54 mm	537-440-00	10 pieces
0.41 mm x 0.56 mm / 16 x 22	each 380 mm	0.62 mm	537-410-00	10 pieces
0.43 mm x 0.56 mm / 17 x 22	each 380 mm	0.64 mm	537-411-00	10 pieces
0.43 mm x 0.64 mm / 17 x 25	each 380 mm	0.70 mm	537-413-00	10 pieces
0.46 mm x 0.64 mm / 18 x 25	each 380 mm	0.71 mm	537-414-00	10 pieces
0.55 mm x 0.64 mm / 21.5 x 25	each 380 mm	0.76 mm	537-423-00	10 pieces
0.55 mm x 0.70 mm / 21.5 x 28	each 380 mm	0.81 mm	537-424-00	10 pieces



CoCr

remaloy® straight wire, rectangular

Strength category: hard plus · 1600 N/mm² - 1800 N/mm²

Color-coded packaging: blue

▣	Length	Diagonal M.	REF	Quantity
0.41 mm x 0.41 mm / 16 x 16	each 380 mm	0.54 mm	537-540-00	10 pieces
0.41 mm x 0.56 mm / 16 x 22	each 380 mm	0.62 mm	537-510-00	10 pieces



Nickel-free titanium molybdenum

rematitan® SPECIAL straight wire, round

Strength category: middle plus · 1100 N/mm² - 1400 N/mm²

Color-coded packaging: yellow

ø	Length	Diagonal M.	REF	Quantity
0.40 mm / 16	each 380 mm		766-600-00	10 pieces
0.45 mm / 18	each 380 mm		766-601-00	10 pieces
0.80 mm / 32	each 360 mm		766-629-00	10 pieces



rematitan® SPECIAL straight wire, rectangular

Strength category: middle plus · 1100 N/mm² - 1400 N/mm²

Color-coded packaging: yellow

▣	Length	Diagonal M.	REF	Quantity
0.41 mm x 0.41 mm / 16 x 16	each 360 mm	0.54 mm	766-633-00	10 pieces
0.41 mm x 0.56 mm / 16 x 22	each 380 mm	0.62 mm	766-602-00	10 pieces
0.43 mm x 0.64 mm / 17 x 25	each 380 mm	0.70 mm	766-603-00	10 pieces
0.48 mm x 0.64 mm / 19 x 25	each 380 mm	0.73 mm	766-604-00	10 pieces



Nickel titanium


rematitan® LITE wire on coils, round

For gentle levelling in the initial phase immediately after bonding.

Strength category: super elastic

ø	Length	Diagonal M.	REF	Quantity
0.30 mm / 12	4.5 m		766-094-00	1 piece
0.35 mm / 14	4.5 m		766-095-00	1 piece
0.40 mm / 16	4.5 m		766-196-00	1 piece

Stainless steel (twisted and braided)


dentaflex® wire on coils, round, 3-strand twisted

Strength category: super spring hard · 2500 N/mm² - 2700 N/mm²

ø	Length	Diagonal M.	REF	Quantity
0.38 mm / 15	4 m		545-738-01	1 piece
0.45 mm / 18	4 m		545-745-01	1 piece


dentaflex® wire on coils, round, 6-strand co-axial

Strength category: super spring hard · 2500 N/mm² - 2700 N/mm²

ø	Length	Diagonal M.	REF	Quantity
0.38 mm / 15	4 m		545-638-01	1 piece
0.45 mm / 18	4 m		545-645-01	1 piece

Stainless steel (nickel chromium)


remanium® wire on coils, round

Strength category: extra spring hard · 2300 N/mm² - 2500 N/mm²

ø	Length	Diagonal M.	REF	Quantity
0.35 mm / 14	7.5 m		530-035-00	1 piece
0.40 mm / 16	7.5 m		530-040-00	1 piece
0.45 mm / 18	7.5 m		530-045-00	1 piece



Fixed appliance technique / Retainer wires

Stainless steel (nickel chromium)

rematain® flat retainer wire

For maxilla and mandible, 6-strand braided.

Easy to adapt, very flat for high level of patient comfort

	Length	Diagonal M.	REF	Quantity
0.25 mm x 0.73 mm / 10 x 29	150 mm	0.77 mm	530-000-00	10 pieces

Titanium

Titanium retainer wire, Ti1

For maxilla and mandible, 3-strand twisted.

Grade 1 indicates a very low oxygen content, relative degree of softness and easy adaptation when cold.

ø	Length	Diagonal M.	REF	Quantity
0.44 mm / 17	50 mm		528-001-01	10 pieces
0.50 mm / 20	50 mm		528-000-01	10 pieces

Titanium retainer wire, Ti5

For maxilla and mandible, 3-strand twisted.

Grade 5 indicates a very high mechanical strength.

ø	Length	Diagonal M.	REF	Quantity
0.44 mm / 17	50 mm		528-001-00	10 pieces
0.50 mm / 20	50 mm		528-000-00	10 pieces

Titanium retainer wire, Ti1

For maxilla and mandible, 3-strand twisted.

Grade 1 indicates a very low oxygen content, relative degree of softness and easy adaptation when cold.

ø	Length	Diagonal M.	REF	Quantity
0.50 mm / 20	2 m		528-000-02	1 piece

Gold

Gold retainer wire

For maxilla and mandible, 3-strand twisted.

ø	Length	Diagonal M.	REF	Quantity
0.50 mm / 20	160 mm		529-000-00	1 piece

Fixed appliance technique / Retainer wires

remanium® lingual retainer, assortment

For lower cuspids 3-3, anatomically preformed.

1 assortment = 20 pieces

Material: stainless steel

ø 0.80 mm / 32

REF 728-500-00	1 assortment
Contents:	
4 x lingual retainer, 22 mm	REF 728-501-00
4 x lingual retainer, 24 mm	REF 728-502-00
4 x lingual retainer, 26 mm	REF 728-503-00
4 x lingual retainer, 28 mm	REF 728-504-00
4 x lingual retainer, 30 mm	REF 728-505-00



remanium® lingual retainer

For lower cuspids 3-3, anatomically preformed.

The stated length conforms with the transverse distance of both lower cuspids.

Material: stainless steel



ø	Length	REF	Quantity
0.80 mm / 32	22 mm	728-501-00	5 pieces
0.80 mm / 32	24 mm	728-502-00	5 pieces
0.80 mm / 32	26 mm	728-503-00	5 pieces
0.80 mm / 32	28 mm	728-504-00	5 pieces
0.80 mm / 32	30 mm	728-505-00	5 pieces

prime4me® RETAIN3R maxilla*

Unique fit using a three-dimensional construction design and highest biocompatibility: milled retainer in grade 5 titanium.

Can be ordered together with prime4me® RETAIN3R optionally: Transfer tray for prime4me® RETAIN3R for placement in the mouth.

REF 430-120-30	1 piece
----------------	---------



prime4me® RETAIN3R mandible*

Unique fit using a three-dimensional construction design and highest biocompatibility: milled retainer in grade 5 titanium.

Can be ordered together with prime4me® RETAIN3R optionally: Transfer tray for prime4me® RETAIN3R for placement in the mouth.

REF 430-120-40	1 piece
----------------	---------



Fixed appliance technique / Accessories



remanium® ligature wire

Material: stainless steel

Strength category: soft · $\leq 800 \text{ N/mm}^2$

ø	Length	REF	Quantity
0.25 mm / 10	80 m	500-025-00	1 piece
0.25 mm / 10	1160 m	501-025-00	1 piece
0.30 mm / 12	60 m	500-030-00	1 piece
0.40 mm / 16	30 m	500-040-00	1 piece
0.50 mm / 20	50 m	503-050-00	1 piece
0.50 mm / 20	320 m	501-050-00	1 piece




remanium® preformed ligature

Material: stainless steel

Strength category: soft · $\leq 800 \text{ N/mm}^2$

Version	ø	REF	Quantity
short	0.25 mm / 10	751-000-00	1000 pieces
long	0.25 mm / 10	751-001-00	1000 pieces



remanium® Kobayashi ligature

With hook for elastics.

Material: stainless steel

Strength category: soft · $\leq 800 \text{ N/mm}^2$

Version	ø	REF	Quantity
short	0.30 mm / 12	751-005-00	100 pieces
short	0.35 mm / 14	751-004-00	100 pieces
long	0.30 mm / 12	751-003-00	100 pieces



Ligature White, short

Material: stainless steel

Strength category: soft · $\leq 800 \text{ N/mm}^2$

ø	REF	Quantity
0.30 mm / 12	768-036-00	100 pieces

Kobayashi White ligature, short

With hook for elastics.

Material: stainless steel

Strength category: soft · $\leq 800 \text{ N/mm}^2$



ø	REF	Quantity
0.40 mm / 16	768-037-00	20 pieces

remanium® tension spring

Material: stainless steel

Outer ø	Inner ø	Length	REF	Quantity
1.15 mm / 45	0.70 mm / 28	1 m	758-165-00	1 piece
1.21 mm / 48	0.75 mm / 30	1 m	758-175-00	1 piece
1.66 mm / 66	1.20 mm / 47	10 mm	758-180-00	6 pieces



rematitan® LITE tension spring

Advantages:

- high elasticity with virtually constant restoring force
- fatigue-free, low, uniform transfer of forces

With eyelets for easy engagement.

Material: nickel titanium

Strength category: super elastic



Outer ø	Inner ø	Length with eyelets	Wire strength ø	REF	Quantity
1.26 mm / 50	0.76 mm / 30	9 mm	0.25 mm / 10	758-160-00	10 pieces
1.26 mm / 50	0.76 mm / 30	12 mm	0.25 mm / 10	758-161-00	10 pieces



Fixed appliance technique / Accessories



Cuspid retraction spring

acc. to Ladanyi

For distalizing cuspids.

Material: stainless steel


Quadrant	REF	Quantity
upper right / lower left	758-861-00	10 pieces
upper left / lower right	758-862-00	10 pieces



remanium® compression spring

Material: stainless steel

Outer ø	Inner ø	Length	REF	Quantity
1.16 mm / 45	0.70 mm / 28	1 m	758-365-00	1 piece
1.21 mm / 48	0.75 mm / 30	1 m	758-375-00	1 piece
1.66 mm / 66	1.20 mm / 47	15 mm	758-380-00	6 pieces



rematitan® LITE compression spring

Advantages:

- high elasticity with virtually constant restoring force
- fatigue-free, low, uniform transfer of forces

May be individually shortened.

Material: nickel titanium

Strength category: super elastic

Inner ø	Length	Wire strength ø	Force	REF	Quantity
0.76 mm / 30	18 cm	0.25 mm / 10	light	758-360-00	3 pieces
1.14 mm / 45	18 cm	0.25 mm / 10	light	758-367-00	3 pieces
0.76 mm / 30	18 cm	0.30 mm / 12	medium	758-361-00	3 pieces
1.00 mm / 40	18 cm	0.30 mm / 12	medium	758-368-00	3 pieces
0.90 mm / 36	18 cm	0.35 mm / 14	heavy	758-362-00	3 pieces



Jaw fracture splint

With hooks.

Material: stainless steel

REF 402-253-00

1 m

**Sliding tubes, round**

Closed with hook. For lip bumper etc.

Material: stainless steel

Inner ø	Length	REF	Quantity
0.60 mm / 23	2 mm	480-803-00	10 pieces
0.90 mm / 36	2 mm	480-809-00	10 pieces

**Sliding tubes, rectangular**

Slotted with hook. Can be placed directly on the archwire and crimped in place with removing ligature. Universal application for the 18/22 technique.

Material: stainless steel

Quadrant	REF	Quantity
upper right / lower left	480-619-00	10 pieces
upper left / lower right	480-620-00	10 pieces

**Sliding tubes, rectangular**

Slotted with wire ø 0.7 mm, spring hard.

Can be placed directly on the arch wire and crimped in place without removing ligature.

Universal application for the 18/22 technique.

Material: stainless steel

Quadrant	REF	Quantity
upper / lower right and left	480-621-00	10 pieces
upper / lower right and left	480-621-02	2 pieces

**Crimping pliers**

sterilizable

Tip with V-shaped notch for holding and crimping sliding/clamping and stop tubes.

Premium-Line · 10 year guarantee

REF 003-711-00

1 piece



Fixed appliance technique / Accessories



Cross tube

Transverse connection between two small rectangular tubes for joining two wires without welding. Ideal for stainless steel or titanium connections.

Horizontal tube:

Overall length: 4.0 mm

Inner $\square$ 0.46 mm x 0.64 mm / 18 x 25

Vertical tube:

Overall length: 2.0 mm

Inner $\square$ 0.56 mm x 0.70 mm / 22 x 28

Material: stainless steel

REF 480-000-00

10 pieces



Cross tube

Transverse connection between two small rectangular tubes for joining two wires without welding. Ideal for stainless steel or titanium connections.

Horizontal tube:

Overall length: 4.0 mm

Inner $\square$: 0.56 mm x 0.70 mm / 22 x 28

Vertical tube:

Overall length: 2.0 mm

Inner $\square$: 0.56 mm x 0.70 mm / 22 x 28

Material: stainless steel

REF 480-100-00

10 pieces



Rectangular tube

For use with suprastructures in orthodontic implantology for holding wires and springs.

Material: stainless steel

Inner $\square$	Length	REF	Quantity
0.56 mm x 0.70 mm / 22 x 28	5 mm	480-022-00	10 pieces



Stop tubes, slotted

Prevents slipping within the bracket buccal tube assembly with nickel titanium arches or as stop for compression springs.

Material: stainless steel

$\square$	REF	Quantity
0.46 mm x 0.64 mm / 18 x 25	480-617-01	100 pieces
0.56 mm x 0.70 mm / 22 x 28	480-618-01	100 pieces



Archwire holder for preformed arches

For 11 different arches.

Especially hygienic due to protective cover.

Material: acrylic

Dimensions: 25.3 cm x 11.0 cm x 10.3 cm (W x H x D)

REF 055-104-00

1 piece

Stainless steel (nickel chromium)



remanium® straight wire, round

Strength category: spring hard · 1800 N/mm² - 2000 N/mm²

Color-coded packaging: red

ø	Length	Diagonal M.	REF	Quantity
0.40 mm / 16	each 380 mm		525-040-00	25 pieces
0.50 mm / 20	each 380 mm		535-002-00	25 pieces
0.60 mm / 23	each 380 mm		535-003-00	25 pieces
0.70 mm / 28	each 300 mm		527-070-00	25 pieces
0.70 mm / 28	each 190 mm		527-005-00	900 pieces
0.70 mm / 28	each 150 mm		527-004-00	1100 pieces
0.75 mm / 30	each 380 mm		535-005-00	25 pieces
0.80 mm / 31	each 300 mm		527-080-00	25 pieces
0.80 mm / 31	each 150 mm		527-006-00	875 pieces
0.90 mm / 36	each 300 mm		527-090-00	25 pieces
1.00 mm / 39	each 300 mm		527-100-00	25 pieces
1.15 mm / 45	each 380 mm		535-009-00	25 pieces

CoCr



remaloy® straight wire, round

Recommended for Quad Helix, Crozat appliances, palatal and lingual bars etc.

Can be heat treated · 450 °C / 842 °F for about 5 min

Strength category: hard · 1400 N/mm² - 1600 N/mm²

Color-coded packaging: blue

ø	Length	Diagonal M.	REF	Quantity
0.70 mm / 28	each 380 mm		528-070-00	10 pieces
0.80 mm / 31	each 380 mm		528-080-00	10 pieces
0.90 mm / 36	each 380 mm		528-090-00	10 pieces
1.00 mm / 39	each 380 mm		528-100-00	10 pieces
1.10 mm / 43	each 380 mm		528-110-00	10 pieces
1.20 mm / 47	each 380 mm		528-120-00	10 pieces
1.30 mm / 51	each 380 mm		528-130-00	10 pieces
1.50 mm / 59	each 380 mm		528-150-00	10 pieces

Removable appliance technique / Straight wires

CoCr



remaloy® straight wire, half-round

Can be heat treated · 450 °C / 842 °F for about 5 min

Strength category: spring hard · 1800 N/mm² - 2000 N/mm²

Color-coded packaging: blue

ø	Length	Diagonal M.	REF	Quantity
1.50 mm x 0.75 mm / 59 x 30	each 380 mm		528-155-00	10 pieces
1.75 mm x 0.90 mm / 69 x 36	each 380 mm		528-158-00	10 pieces



remaloy® straight wire, rectangular rounded

With rounded edges, for use as TPA or reinforcing bow.

Can be heat treated · 450 °C / 842 °F for about 5 min

Strength category: hard plus · 1600 N/mm² - 1800 N/mm²

Color-coded packaging: blue

ø	Length	Diagonal M.	REF	Quantity
1.92 mm x 0.90 mm / 76 x 36	each 380 mm		528-159-00	10 pieces

Stainless steel (nickel chromium)


remanium® laboratory coil, round

Strength category: hard - 1400 N/mm² - 1600 N/mm²

ø	Length	Diagonal M.	REF	Quantity
0.70 mm / 28	30 m		513-070-00	1 piece
0.80 mm / 31	20 m		513-080-00	1 piece
0.90 mm / 36	10 m		513-090-00	1 piece
1.00 mm / 39	10 m		513-100-00	1 piece
1.20 mm / 47	10 m		513-120-00	1 piece
1.50 mm / 59	10 m		513-150-00	1 piece


remanium® clinical coil, round

Strength category: hard - 1400 N/mm² - 1600 N/mm²

ø	Length	Diagonal M.	REF	Quantity
0.70 mm / 28	165 m		514-070-00	1 piece
0.80 mm / 31	125 m		514-080-00	1 piece
0.90 mm / 36	100 m		514-090-00	1 piece
1.00 mm / 39	80 m		514-100-00	1 piece


remanium® laboratory coil, round

Strength category: spring hard - 1800 N/mm² - 2000 N/mm²

ø	Length	Diagonal M.	REF	Quantity
0.40 mm / 16	30 m		521-040-00	1 piece
0.50 mm / 20	50 m		523-050-00	1 piece
0.60 mm / 23	40 m		523-060-00	1 piece
0.70 mm / 28	30 m		523-070-00	1 piece
0.80 mm / 31	20 m		523-080-00	1 piece
0.90 mm / 36	10 m		523-090-00	1 piece
1.00 mm / 39	10 m		523-100-00	1 piece
1.10 mm / 43	10 m		523-110-00	1 piece
1.20 mm / 47	10 m		523-120-00	1 piece
1.50 mm / 59	10 m		523-150-00	1 piece


remanium® clinical coil, round

Strength category: spring hard - 1800 N/mm² - 2000 N/mm²

ø	Length	Diagonal M.	REF	Quantity
0.60 mm / 23	225 m		524-060-00	1 piece
0.70 mm / 28	165 m		524-070-00	1 piece
0.80 mm / 31	125 m		524-080-00	1 piece
0.90 mm / 36	100 m		524-090-00	1 piece
1.00 mm / 39	80 m		524-100-00	1 piece
1.10 mm / 43	70 m		524-110-00	1 piece

Removable appliance technique / Wires on coils

Stainless steel (nickel chromium)



remanium® wire on coils, half-round

Strength category: hard · 1400 N/mm² - 1600 N/mm²

ø	Length	Diagonal M.	REF	Quantity
1.30 mm x 0.65 mm / 51 x 26	10 m		308-513-00	1 piece
1.50 mm x 0.75 mm / 59 x 30	10 m		308-515-00	1 piece
1.75 mm x 0.90 mm / 69 x 36	10 m		308-518-00	1 piece
2.00 mm x 1.00 mm / 79 x 40	10 m		308-520-00	1 piece



remanium® clinical coil, half-round

Strength category: hard · 1400 N/mm² - 1600 N/mm²

ø	Length	Diagonal M.	REF	Quantity
1.30 mm x 0.65 mm / 51 x 26	90 m		308-713-00	1 piece
1.50 mm x 0.75 mm / 59 x 30	65 m		308-715-00	1 piece
1.75 mm x 0.90 mm / 69 x 36	50 m		308-718-00	1 piece



remanium® wire on coils, oval

Strength category: hard · 1400 N/mm² - 1600 N/mm²

ø	Length	Diagonal M.	REF	Quantity
2.40 mm x 1.40 mm / 95 x 55	10 m		385-624-00	1 piece
3.00 mm x 1.50 mm / 118 x 59	10 m		385-630-00	1 piece

Nickel-free* stainless steel



Noninium® wire on coils, round

Strength category: hard · 1400 N/mm² - 1600 N/mm²

ø	Length	Diagonal M.	REF	Quantity
0.70 mm / 28	30 m		520-070-10	1 piece
0.80 mm / 31	20 m		520-080-10	1 piece
0.90 mm / 36	10 m		520-090-10	1 piece
1.20 mm / 47	10 m		520-120-10	1 piece

* Nickel trace elements of ≤ 0.05 % are due to the process of manufacturing



Noninium® wire on coils, round

Strength category: spring hard · 1800 N/mm² - 2000 N/mm²

ø	Length	Diagonal M.	REF	Quantity
0.60 mm / 23	40 m		520-062-10	1 piece
0.70 mm / 28	30 m		520-072-10	1 piece
0.80 mm / 31	20 m		520-082-10	1 piece
0.90 mm / 36	10 m		520-092-10	1 piece

* Nickel trace elements of ≤ 0.05 % are due to the process of manufacturing


Removable appliance technique / Accessories

Noninium® standard labial arch

Material: nickel-free* stainless steel

Strength category: spring hard · 1800 N/mm² - 2000 N/mm²

ø	Length mesial-mesial	REF	Quantity
0.80 mm / 31	28 mm	621-028-01	10 pieces
0.80 mm / 31	31 mm	621-031-01	10 pieces
0.80 mm / 31	34 mm	621-034-01	10 pieces
0.80 mm / 31	37 mm	621-037-01	10 pieces
0.80 mm / 31	40 mm	621-040-01	10 pieces
0.80 mm / 31	43 mm	621-043-01	10 pieces
0.80 mm / 31	46 mm	621-046-01	10 pieces



* Nickel trace elements of ≤ 0.05 % are due to the process of manufacturing

Noninium® Adams clasp, assortment

Strength category: spring hard · 1800 N/mm² - 2000 N/mm²

1 assortment = 250 pieces

Material: nickel-free* stainless steel

ø 0.70 mm / 28

REF 621-500-01	1 assortment
Contents:	
30 x Adams clasp, 7 mm	REF 621-404-01
30 x Adams clasp, 8 mm	REF 621-405-01
40 x Adams clasp, 9 mm	REF 621-406-01
60 x Adams clasp, 10 mm	REF 621-407-01
60 x Adams clasp, 11 mm	REF 621-408-01
30 x Adams clasp, 12 mm	REF 621-409-01
1 x ruler	REF 766-000-00



* Nickel trace elements of ≤ 0.05 % are due to the process of manufacturing

Noninium® Adams clasp

Material: nickel-free* stainless steel

Strength category: spring hard · 1800 N/mm² - 2000 N/mm²

ø	Length mesial-distal	REF	Quantity
0.70 mm / 28	7 mm	621-404-01	10 pieces
0.70 mm / 28	8 mm	621-405-01	10 pieces
0.70 mm / 28	9 mm	621-406-01	10 pieces
0.70 mm / 28	10 mm	621-407-01	10 pieces
0.70 mm / 28	11 mm	621-408-01	10 pieces
0.70 mm / 28	12 mm	621-409-01	10 pieces



* Nickel trace elements of ≤ 0.05 % are due to the process of manufacturing

Removable appliance technique / Accessories



Noninium® triangular clasp

Material: nickel-free* stainless steel

Strength category: hard · 1400 N/mm² - 1600 N/mm²

ø	Length	REF	Quantity
0.70 mm / 28	each 36 mm	621-300-01	100 pieces

* Nickel trace elements of ≤ 0.05 % are due to the process of manufacturing



Noninium® arrow clasp

Material: nickel-free* stainless steel

Strength category: spring hard · 1800 N/mm² - 2000 N/mm²

ø	Length	REF	Quantity
0.70 mm / 28	each 41 mm	621-200-01	100 pieces

* Nickel trace elements of ≤ 0.05 % are due to the process of manufacturing



Noninium® ball retainer clasp

Material: nickel-free* stainless steel

Strength category: hard · 1400 N/mm² - 1600 N/mm²

ø	Length	REF	Quantity
0.70 mm / 28	each 38 mm	620-207-02	100 pieces
0.80 mm / 32	each 38 mm	620-208-02	100 pieces
0.90 mm / 36	each 38 mm	620-209-02	100 pieces

* Nickel trace elements of ≤ 0.05 % are due to the process of manufacturing



remanium® ball retainer clasp

Material: stainless steel

Strength category: hard · 1400 N/mm² - 1600 N/mm²

ø	Length	REF	Quantity
0.70 mm / 28	each 38 mm	620-107-00	100 pieces
0.80 mm / 31	each 38 mm	620-108-00	100 pieces
0.90 mm / 36	each 38 mm	620-109-00	100 pieces
1.00 mm / 39	each 38 mm	620-110-00	100 pieces

Removable appliance technique / Accessories

Guide pin

acc. to G. Müller®

Guide pin between upper and lower plates intended for bite correction of Angle Class II + Class III.

Material: stainless steel

REF 607-204-00	2 pieces
REF 607-204-30	20 pieces


Guide pin

acc. to Hinz

Material: stainless steel

ø	REF	Quantity	REF	Quantity
1.50 mm / 59	607-202-00	1 piece	607-202-30	10 pieces


Lingual bar

Anatomically shaped e.g. for occlusion plates.

Material: stainless steel

Length	Dimensions	REF	Quantity
6 cm	3.0 mm x 1.5 mm	311-306-30	10 pieces
7 cm	3.0 mm x 1.5 mm	311-307-30	10 pieces
8 cm	3.0 mm x 1.5 mm	311-308-30	10 pieces



Removable appliance technique / Accessories



Strengthened

Twisted and rolled flat.

Material: stainless steel

Length	☑	REF	Quantity
10 m	1.8 mm x 0.8 mm	312-106-00	1 piece
10 m	2.4 mm x 0.8 mm	312-107-00	1 piece



Stainless steel tube

Round

Material: stainless steel

Inner ø	Length	REF	Quantity
0.50 mm / 20	30 cm	480-001-00	1 piece
0.60 mm / 23	30 cm	480-003-00	1 piece
0.70 mm / 28	30 cm	480-005-00	1 piece
0.80 mm / 31	30 cm	480-007-00	1 piece
0.90 mm / 36	30 cm	480-009-00	1 piece
1.00 mm / 39	30 cm	480-011-00	1 piece
1.10 mm / 43	30 cm	480-013-00	1 piece
1.20 mm / 47	30 cm	480-015-00	1 piece
1.30 mm / 51	30 cm	480-017-00	1 piece



remanium® Goshgarian palatal bars, assortment

Strength category: hard · 1400 N/mm² - 1600 N/mm²

Reverse Loop – mesial

The inverted loop is better for suppressing the gag reflex.

Enables controllable movements of upper molars.

1 assortment = 30 pieces

Material: stainless steel

ø 0.90 mm / 36

REF 728-020-00	1 assortment
Contents:	
5 x palatal bars, 37 mm	REF 728-021-00
5 x palatal bars, 40 mm	REF 728-022-00
5 x palatal bars, 43 mm	REF 728-023-00
5 x palatal bars, 46 mm	REF 728-024-00
5 x palatal bars, 49 mm	REF 728-025-00
5 x palatal bars, 52 mm	REF 728-026-00


remanium® Goshgarian palatal bars

Reverse Loop – mesial

The inverted loop is better for suppressing the gag reflex.
Enables controllable movements of upper molars.

Material: stainless steel

Strength category: hard · 1400 N/mm² - 1600 N/mm²



ø	Size	REF	Quantity
0.90 mm / 36	1 = 37 mm	728-021-00	10 pieces
0.90 mm / 36	2 = 40 mm	728-022-00	10 pieces
0.90 mm / 36	3 = 43 mm	728-023-00	10 pieces
0.90 mm / 36	4 = 46 mm	728-024-00	10 pieces
0.90 mm / 36	5 = 49 mm	728-025-00	10 pieces
0.90 mm / 36	6 = 52 mm	728-026-00	10 pieces


remanium® Goshgarian palatal bars

Standard Loop – distal

Enables controllable movement of upper molars.

Material: stainless steel

Strength category: hard · 1400 N/mm² - 1600 N/mm²



ø	Size	REF	Quantity
0.90 mm / 36	1 = 34 mm	728-005-00	10 pieces
0.90 mm / 36	2 = 37 mm	728-006-00	10 pieces
0.90 mm / 36	3 = 41 mm	728-007-00	10 pieces
0.90 mm / 36	4 = 45 mm	728-008-00	10 pieces
0.90 mm / 36	5 = 47 mm	728-013-00	10 pieces
0.90 mm / 36	6 = 49 mm	728-014-00	10 pieces
0.90 mm / 36	7 = 52 mm	728-015-00	10 pieces


remanium® Quad Helix

Removable; for insertion into palatal sheaths on the upper molars.

Application: for expansion of posterior segments

Material: stainless steel

Strength category: spring hard · 1800 N/mm² - 2000 N/mm²



ø	Length	Size	REF	Quantity
0.90 mm / 36	17.3 mm	1 = 36 mm	728-100-01	10 pieces
0.90 mm / 36	19.0 mm	2 = 38 mm	728-101-01	10 pieces
0.90 mm / 36	20.7 mm	3 = 40 mm	728-102-01	10 pieces

Palatal and lingual arch systems




Orthorama® palatal arches, assortment

1 assortment = 10 pieces

Material: CoCr

REF 728-205-00	1 assortment
Contents:	
2 x Orthorama® palatal arch, 40 mm	REF 728-200-00
2 x Orthorama® palatal arch, 45 mm	REF 728-201-00
2 x Orthorama® palatal arch, 50 mm	REF 728-202-00
2 x Orthorama® palatal arch, 55 mm	REF 728-203-00
2 x Orthorama® palatal arch, 60 mm	REF 728-204-00
1 x ruler	REF 766-000-00



Orthorama® palatal arch

Material: CoCr

ø	Size	REF	Quantity
0.90 mm / 36	1 = 40 mm	728-200-00	10 pieces
0.90 mm / 36	2 = 45 mm	728-201-00	10 pieces
0.90 mm / 36	3 = 50 mm	728-202-00	10 pieces
0.90 mm / 36	4 = 55 mm	728-203-00	10 pieces
0.90 mm / 36	5 = 60 mm	728-204-00	10 pieces




Orthorama® lingual arches, assortment

1 assortment = 18 pieces

Material: CoCr

REF 728-215-00	1 assortment
Contents:	
2 x Orthorama® lingual arch, 48 mm	REF 728-206-00
2 x Orthorama® lingual arch, 51 mm	REF 728-207-00
2 x Orthorama® lingual arch, 54 mm	REF 728-208-00
2 x Orthorama® lingual arch, 57 mm	REF 728-209-00
2 x Orthorama® lingual arch, 60 mm	REF 728-210-00
2 x Orthorama® lingual arch, 63 mm	REF 728-211-00
2 x Orthorama® lingual arch, 66 mm	REF 728-212-00
2 x Orthorama® lingual arch, 69 mm	REF 728-213-00
2 x Orthorama® lingual arch, 72 mm	REF 728-214-00
1 x ruler	REF 766-000-00



Orthorama® lingual arch

Material: CoCr

ø	Size	REF	Quantity
0.90 mm / 36	1 = 48 mm	728-206-00	10 pieces
0.90 mm / 36	2 = 51 mm	728-207-00	10 pieces
0.90 mm / 36	3 = 54 mm	728-208-00	10 pieces
0.90 mm / 36	4 = 57 mm	728-209-00	10 pieces
0.90 mm / 36	5 = 60 mm	728-210-00	10 pieces
0.90 mm / 36	6 = 63 mm	728-211-00	10 pieces
0.90 mm / 36	7 = 66 mm	728-212-00	10 pieces
0.90 mm / 36	8 = 69 mm	728-213-00	10 pieces
0.90 mm / 36	9 = 72 mm	728-214-00	10 pieces

Palatal and lingual arch systems

Orthorama® Multi-Action lingual arches assortment

1 assortment = 8 pieces

Material: CoCr


REF 728-222-00	1 assortment
Contents:	
2 x Orthorama® Multi-Action lingual arch, 53 mm	REF 728-218-00
2 x Orthorama® Multi-Action lingual arch, 57 mm	REF 728-219-00
2 x Orthorama® Multi-Action lingual arch, 61 mm	REF 728-220-00
2 x Orthorama® Multi-Action lingual arch, 65 mm	REF 728-221-00
1 x ruler	REF 766-000-00



Orthorama® Multi-Action lingual arch

Material: CoCr

ø	Size	REF	Quantity
0.90 mm / 36	1 = 53 mm	728-218-00	10 pieces
0.90 mm / 36	2 = 57 mm	728-219-00	10 pieces
0.90 mm / 36	3 = 61 mm	728-220-00	10 pieces
0.90 mm / 36	4 = 65 mm	728-221-00	10 pieces



Orthorama® sectional arch

Material: CoCr

Version	REF	Quantity
with Omega loop	728-216-00	10 pieces
without Omega loop	728-217-00	10 pieces



Orthorama® transfer system

The indirect transfer system permits the precise fitting of an Orthorama® palatal or lingual arch onto a plaster model, as well as the identical transfer of the palatal or lingual sheath's position in the mouth onto the model.

Material: stainless steel

REF 728-400-00	10 pairs
Contents:	
10 x transfer jigs	REF 728-401-00
10 x retention band	REF 728-402-00



Orthorama® transfer jig

Material: stainless steel

REF 728-401-00	10 pieces
----------------	-----------



Orthorama® retention band

Material: stainless steel

REF 728-402-00	10 pieces
----------------	-----------



SUS ³ – Sabbagh Universal Spring	Page	244
Herbst bite jumping hinge / Sets	Page	246
Herbst bite jumping hinge / Accessories	Page	248
amda® / Starter Set	Page	251
amda® / Appliance and accessories	Page	252
tomas® / Pins	Page	254



Class II appliances.

Achieve Class I successfully.




DENTAURUM
**ONLINE
SHOP**
shop.dentaurum.com






The ideal path
to achieving Class I.



SUS³ – first-class treatment for your Class II patients.

The SUS³ Sabbagh Universal Spring is a telescope appliance with spring for universal intermaxillary use when treating Class II. It provides a perfect alternative for patients with compliance issues. The SUS³ acts with constant, gentle forces, mainly horizontal, on the temporomandibular joint, and can help to avoid the need for some extractions and dysgnathic surgery.





Indications:

- Distoclusion/Class II (unilateral/bilateral).
- Molar distalization (headgear substitute).
- Dentoalveolar compensation of the occlusion (elastics substitute).
- Midline displacement.
- Gap closure (aplasia)/anchorage (alternative to mini screws).
- Alternative to dysgnathia operation.
- Temporomandibular joint disorder and craniomandibular dysfunction (repositioning effect).
- Therapy for snoring and sleep apnea.

The ideal solution for patients with:

- poor compliance.
- minimal residual growth.
- hindrance to using removable appliances.
- illnesses to the primary respiratory tract or asthma.
- allergy to components in acrylic.
- disability.

Herbst bite jumping hinges.

Safe, proven, effective.



The classic correction of distocclusions.

Dentaurum's program consists of four versions of the proven Herbst bite jumping hinge. Each of the versions has a high mechanical stability and an excellent finish which guarantee adequate safety reserves for an effective therapy, even when under great stress.





Herbst TS.

This Herbst version is fitted with an inner telescope. This prevents the hinge from sliding apart.

The new laser markings on the telescopic element indicating size and side guarantee clear allocation to the jaw.

Herbst I.

Classic retention hinge for welding to bands.

The Dentaureum Herbst I bite jumping hinge is the classic version that is welded onto bands or cast splints, depending on the aim of the treatment and the experience of the orthodontist. This means maximum flexibility and enables individual treatment.

Herbst II.

Modification of classic appliance for fixing onto the archwire of multi-bracket appliances.

Herbst IV.

Modified version with ball and socket joint and C-clips for fixing. Greater lateral freedom of movement.

This clever solution opens up a new dimension of movement in the mandible during treatment with Herbst bite jumping hinges since the ball and socket joint on Herbst IV greatly improves patient comfort.

The C-clips mean screws need not be screwed in place on the Dentaureum Herbst IV retention hinge. This makes it easier to handle and offers in particular more security.

Class II appliances

SUS³ – Sabbagh Universal Spring



SUS³ starter set



REF 607-160-00	1 set
-----------------------	-------

Contents:

2 x SUS ³ telescope element	REF 607-157-00
2 x SUS ³ telescope rod, assembled	REF 607-158-00
2 x SUS ³ ball retainer clasp, ø 1.0 mm	REF 607-154-00
2 x SUS ³ turbo spring (stainless steel), length 27 mm	REF 607-156-00
2 x spacer ring, crimpable, length 1 mm	REF 607-103-11
2 x spacer ring, crimpable, length 2 mm	REF 607-104-11
1 x hexagon socket key 1.5	REF 607-144-00



SUS³ set

An Allen key is not included in the delivery.



REF 607-165-00	1 set
-----------------------	-------

Contents:

2 x SUS ³ telescope element	REF 607-157-00
2 x SUS ³ telescope rod, assembled	REF 607-158-00
2 x SUS ³ ball retainer clasp, ø 1.0 mm	REF 607-154-00
2 x SUS ³ turbo spring (stainless steel), length 27 mm	REF 607-156-00
2 x spacer ring, crimpable, length 1 mm	REF 607-103-11
2 x spacer ring, crimpable, length 2 mm	REF 607-104-11



SUS³ telescope rod, assembled

REF 607-158-00	2 pieces
-----------------------	----------

Contents:

2 x SUS ³ arch adapter / fixing screw (1 pair)	REF 607-153-00
2 x SUS ³ telescope rod	REF 607-159-00



SUS³ telescope rod

REF 607-159-00	2 pieces
-----------------------	----------



SUS³ arch adapter / fixing screw

REF 607-153-00	1 pair
-----------------------	--------



SUS³ telescope element

REF 607-157-00	1 piece
-----------------------	---------

SUS³ – Sabbagh Universal Spring

SUS³ turbo spring (stainless steel)

Length: 27 mm

REF 607-156-00

6 pieces



SUS³ compression spring (stainless steel)

Length: 15 mm

REF 607-151-00

6 pieces



SUS³ distance spring

Length: 10 mm

REF 607-152-00

6 pieces



Crimping pliers, for Herbst TS / SUS sterilizable

For clamping slotted spacers for Herbst TS and SUS.

Premium-Line · 10 year guarantee

REF 003-710-00

1 piece



SUS³ ball retainer clasp

Shaft ø: 1 mm

Head ø: 2.1 mm

Length: 18 mm

REF 607-154-00

10 pieces



Spacer ring, crimpable

Length	for	REF	Quantity
1 mm	Herbst TS, SUS	607-103-11	10 pieces
2 mm	Herbst TS, SUS	607-104-11	10 pieces
3 mm	Herbst TS, SUS	607-105-11	10 pieces
4 mm	Herbst TS, SUS	607-107-11	10 pieces



Hexagon socket key 1.5

sterilizable

Key width	for	REF	Quantity
1.5 mm	Herbst TS / I / II, SUS	607-144-00	1 piece



Class II appliances

Herbst bite jumping hinges / Sets



Herbst TS set

Telescope system. For welding onto bands or cast splints.

Contents:

- 1 x telescopic tube, right
- 1 x telescopic tube, left
- 4 x hexagon socket screw
- 4 x F-base

A hex key is not included in the set!

Length	REF	Quantity
20 mm	607-126-20	1 set
22 mm	607-126-22	1 set
24 mm	607-126-24	1 set

Please note that we offer a welding service for double bands for the Herbst technique free of charge. Please contact your local medical product advisor.



Herbst I set

With hexagon socket screws.

For welding onto bands or cast splints.

A hex key is not included in the set!



REF 607-120-00	1 set
-----------------------	-------

Contents:

1 x telescopic tube, right	REF 607-121-00
1 x telescopic tube, left	REF 607-122-00
2 x telescopic rod	REF 607-124-00
4 x hexagon socket screw	REF 607-122-30
4 x F-base	REF 607-127-30

Please note that we offer a welding service for double bands for the Herbst technique free of charge. Please contact your local medical product advisor.



Herbst II set

With hexagon socket screws.

To fix onto arch wires with the multi-bracket technique.

A hex key is not included in the set!



REF 607-108-00	1 set
-----------------------	-------

Contents:

1 x telescopic tube, right	REF 607-121-00
1 x telescopic tube, left	REF 607-122-00
2 x telescopic rod	REF 607-124-00
2 x hexagon socket screw, short	REF 607-108-10
2 x hexagon socket screw, long	REF 607-108-20
4 x pivot	REF 607-111-00

Herbst bite jumping hinges / Sets

Herbst IV set

With ball head.

For welding onto bands or cast splints. Big lateral movement. It is retained on the base using lock washers.

No screw retention and consequently no lost screws.

Great patient comfort.

Contents:

- 2 x telescopic tube with ball and socket joint
- 2 x telescopic rod with ball head
- 2 x coupling for ball head
- 4 x base
- 10 x C-clip, lock washer



REF 607-115-00

1 set

Please note that we offer a welding service for double bands for the Herbst technique free of charge. Please contact your local medical product advisor.



Class II appliances

Herbst bite jumping hinges / Accessories

Hexagon socket screw



Version	for	REF	Quantity
with base	Herbst TS / I	607-121-30	10 pairs
without base	Herbst TS / I	607-122-30	10 pieces
with short shank and tip	Herbst II	607-108-10	10 pieces
with long shank and tip	Herbst II	607-108-20	10 pieces

Slotted screw



Version	for	REF	Quantity
with base	Herbst TS / I	607-101-00	10 pairs
without base	Herbst TS / I	607-102-00	10 pieces
with short shank and tip	Herbst II	607-109-00	10 pieces
with long shank and tip	Herbst II	607-110-00	10 pieces

Lock washer



Designation	for	REF	Quantity
C-clip	Herbst IV	607-116-00	10 pieces

Base / pivot

Designation	for	REF	Quantity
F-base, form concave	Herbst TS / I	607-127-30	10 pieces
Base, form cylindrical	Herbst TS / I	607-117-00	10 pieces
Base, curved inner shape	Herbst IV	607-118-00	10 pieces
Pivot	Herbst II	607-111-00	10 pieces



Telescopic tube

Designation	for	REF	Quantity
Telescopic tube, right	Herbst I / II	607-121-00	1 piece
Telescopic tube, left	Herbst I / II	607-122-00	1 piece
Telescopic tube with ball joint	Herbst IV	607-125-00	1 piece



Telescopic rod

Designation	for	REF	Quantity
Telescopic rod	Herbst I / II	607-124-00	1 piece
Telescopic rod with ball head	Herbst IV	607-128-00	1 piece



Class II appliances

Herbst bite jumping hinges / Accessories



Spacer ring

For Herbst I to IV.

Length	for	REF	Quantity
1 mm	Herbst I / II / IV	607-103-00	10 pieces
2 mm	Herbst I / II / IV	607-104-00	10 pieces
3 mm	Herbst I / II / IV	607-105-00	10 pieces
4 mm	Herbst I / II / IV	607-107-00	10 pieces
5 mm	Herbst I / II / IV	607-112-00	10 pieces



Spacer ring, crimpable

Length	for	REF	Quantity
1 mm	Herbst TS, SUS	607-103-11	10 pieces
2 mm	Herbst TS, SUS	607-104-11	10 pieces
3 mm	Herbst TS, SUS	607-105-11	10 pieces
4 mm	Herbst TS, SUS	607-107-11	10 pieces



Crimping pliers, for Herbst TS / SUS

sterilizable

For clamping slotted spacers for Herbst TS and SUS.

Premium-Line - 10 year guarantee

REF 003-710-00

1 piece



Hexagon socket key 1.5

sterilizable

Key width	for	REF	Quantity
1.5 mm	Herbst TS / I / II, SUS	607-144-00	1 piece



Hexagon socket key

Only intended for use in the laboratory!

Key width	for	REF	Quantity
1.5 mm	Herbst TS / I	607-123-10	1 piece

amda® starter set (assembled)

Skeletally anchored appliance for the distalization of molars with the tomas®-pin EP.



REF 607-180-00	1 set
----------------	-------

Contents:

1 x amda® telescope (2 pieces)	REF 607-181-02
1 x amda® connector, pair 7 mm	REF 607-184-07
1 x amda® connector, pair 10 mm	REF 607-184-10
1 x amda® connector, pair 13 mm	REF 607-184-13
1 x amda® palatal arch	REF 607-185-01
1 x hexagon socket key 0.9	REF 607-129-00

Class II appliances

amda® / Appliance and accessories



amda® telescope

Encapsulated stainless steel spring with two stop screws
(REF 302-013-01).

Spring force: max. 4.7 N

Spring travel: max. 4.0 mm

REF 607-181-02	2 pieces
REF 607-181-10	10 pieces



amda® connector, left

Elements to connect the amda® appliance and the palatal sheath on the molars.


Height	REF	Quantity
7 mm	607-182-07	10 pieces
10 mm	607-182-10	10 pieces
13 mm	607-182-13	10 pieces



amda® connector, right

Elements to connect the amda® appliance and the palatal sheath on the molars.

Height	REF	Quantity
7 mm	607-183-07	10 pieces
10 mm	607-183-10	10 pieces
13 mm	607-183-13	10 pieces



amda® connector, pair

Elements to connect the amda® appliance and the palatal sheath on the molars.

Height	REF	Quantity
7 mm	607-184-07	1 pair
10 mm	607-184-10	1 pair
13 mm	607-184-13	1 pair

amda® palatal arch

(incl. 2 x tomas®-abutments EP, REF 302-027-00).

Palatal arch for amda® for fixing the amda® telescopes.

REF 607-185-01	1 piece
REF 607-185-10	10 pieces



Hexagon socket key 0.9

sterilizable

Key width: 0.9 mm

REF 607-129-00	1 piece
----------------	---------



Orthodontic demonstration model tomas® / amda®

Fitted with:

- tomas®-pin EP
- amda®
- lingual / palatal sheaths easy in,
welded to standard molar bands

Maxilla made of plastic.

REF 892-000-26	1 piece
----------------	---------



tomas® / Pins



tomas®-pin EP

self-drilling

sterile

Titanium grade 5

Thread ø: 1.6 mm

Core ø: 1.2 mm

Head type: mushroom-shaped

Length	Color-coded packaging	REF	Quantity
6 mm	white	302-206-00	1 piece
8 mm	grey	302-208-00	1 piece
10 mm	black	302-210-00	1 piece



tomas®-pin EP N

self-drilling

non-sterile

sterilizable

Titanium grade 5

Thread ø: 1.6 mm

Core ø: 1.2 mm

Head type: mushroom-shaped

Length	REF	Quantity	REF	Quantity
6 mm	302-206-10	1 piece	302-206-30	10 pieces
8 mm	302-208-10	1 piece	302-208-30	10 pieces
10 mm	302-210-10	1 piece	302-210-30	10 pieces

Extra-oral appliances	Page	260
Accessories / Extra-oral	Page	262
Elastics and auxiliaries / Intra-oral	Page	264
Oral screens	Page	268




Intra- and extra-oral.

Appliances for every treatment technique.



DENTAURUM
**ONLINE
SHOP**
shop.dentaurum.com





Complete delivery program for all areas of treatment.

The range of intra-oral and extra-oral products from Dentaureum extends from simple elastics or separators for intra-oral use to innovative appliances such as the FaceFormer.

During the development, production and selection of all intra- and extra-oral products, great care is taken to comply with the Dentaureum quality standards. Not only does this guarantee a high degree of security for the practitioner and patient, but it also enables time-saving and effective treatment.



Elastic auxiliaries


Available in many colors, diameters and with various traction forces.

Dentalastics® separators

For quick and easy separation of teeth.



Extra-oral appliances



Noninium® facebow, medium

nickel-free

(Nickel trace elements of $\leq 0.2\%$ are due to the process of manufacturing)

Break resistant laser welded connection between outer and inner bow.

Safety Note:

Facebows must always be used with headgear and neck bands with safety modules!

Inner bow ø	Safety caps	REF	Quantity
1.15 mm / 45	white	740-000-00	1 piece




Standard facebow, mini

Break resistant laser welded connection between outer and inner bow.

Safety Note:

Facebows must always be used with headgear and neck bands with safety modules!

Inner bow ø	Stop loops	Safety caps	REF	Quantity
1.15 mm / 45		white	740-019-00	1 piece
1.15 mm / 45	97 mm	red	741-049-00	1 piece



Standard facebow, medium

Break resistant laser welded connection between outer and inner bow.

Safety Note:

Facebows must always be used with headgear and neck bands with safety modules!

Inner bow ø	Safety caps	REF	Quantity
1.00 mm / 40	white	740-010-00	1 piece
1.15 mm / 45	white	740-020-00	1 piece

Standard facebow, maxi

Break resistant laser welded connection between outer and inner bow.

Safety Note:

Facebows must always be used with headgear and neck bands with safety modules!

Inner bow ø	Safety caps	REF	Quantity
1.15 mm / 45	white	740-021-00	1 piece



Facebow, medium, with stop loops

Facebow: without coating


Inner bow ø: 1.15 mm / 45

Break resistant laser welded connection between outer and inner bow.

Safety Note:

Facebows must always be used with headgear and neck bands with safety modules!

Outer bow size	Inner bow size	Stop loops	Safety caps / Size	Cuspid hooks	REF	Quantity
49 mm	29 mm	83 mm	yellow / Size 1		741-040-00	1 piece
49 mm	32 mm	90 mm	grey / Size 2		741-045-00	1 piece
49 mm	35 mm	97 mm	red / Size 3		741-050-00	1 piece
49 mm	38 mm	104 mm	green / Size 4		741-055-00	1 piece
49 mm	41 mm	111 mm	blue / Size 5		741-060-00	1 piece
49 mm	35 mm	97 mm	red / Size 3	•	741-150-00	1 piece



Accessories / Extra-oral



Dentalastics® extra-oral plastic chain

Color	Length	REF	Quantity
transparent	118 mm	774-300-00	10 pieces



Extra-oral latex elastics

Highest elasticity and fatigue resistance.

The elastic packs are marked according to size and strength with symbols.

Every type can therefore be quickly identified.

Strength:

Strength ratings are based on the elastic being stretched 3 times the listed inner diameter.

Content:

30 x 100 rings = 3000 pieces or 30 x 50 rings = 1500 pieces



Traction force	Inner ø	REF	Quantity
light 2.2 N (225 g)	6.4 mm / 1/4 inch	772-203-02	3000 pieces
light 2.2 N (225 g)	8.0 mm / 5/16 inch	772-204-02	3000 pieces
light 2.2 N (225 g)	9.5 mm / 3/8 inch	772-205-02	3000 pieces
strong 4.0 N (400 g)	8.0 mm / 5/16 inch	772-208-02	1500 pieces
strong 4.0 N (400 g)	12.7 mm / 1/2 inch	772-210-02	1500 pieces



Lip bumper

With plastic shield and hook for intermaxillary elastics.

Indications:

- for distalization of molars
- for stabilization of molars in the late mixed dentition
- with labial tipping of the anterior teeth caused by tongue pressure

for tubes	REF	Quantity
ø 0.90 mm / 36	748-009-00	1 piece
ø 1.15 mm / 45	748-114-00	1 piece

Stop tube

For facebow and lip bumper.

Inner ø	Length	REF	Quantity
1.15 mm / 45	2 mm	480-914-00	100 pieces



Stop screw

acc. to Ladanyi

For facebow and lip bumper.

Inner ø	REF	Quantity
1.15 mm / 45	480-918-02	10 pieces



Key for hexagon socket screws

For stop screws acc. to Ladanyi.

Only intended for use in the laboratory!

Key width	REF	Quantity
0.9 mm	480-924-02	1 piece



Activator tubes

E.g. for Teuscher activator. These tubes can be built directly into an activator to hold a facebow.

Used for attaching a facebow.

Inner ø	REF	Quantity
1.15 mm / 45	607-300-00	10 pieces
1.30 mm / 51	607-301-00	10 pieces



Elastics and auxiliaries / Intra-oral



Elastic placer

For seating and removing intra-oral elastics.
Color: assorted colors

REF 772-400-00

50 pieces



Intra-oral latex elastics

Highest elasticity and fatigue resistance.

The elastic packs are marked according to size and strength with symbols.


Every type can therefore be quickly identified.

Strength:

Strength ratings are based on the elastic being stretched 3 times the listed inner diameter.

Content:

10 x 100 pieces = 1000 pieces



Color	Traction force	Inner ø	REF	Quantity
natural	light 0.7 N (71 g)	3.2 mm / 1/8 inch	772-103-00	1000 pieces
natural	light 0.7 N (71 g)	4.8 mm / 3/16 inch	772-104-00	1000 pieces
natural	light 0.7 N (71 g)	6.4 mm / 1/4 inch	772-106-00	1000 pieces
natural	light 0.7 N (71 g)	8.0 mm / 5/16 inch	772-107-00	1000 pieces
natural	light 0.7 N (71 g)	9.5 mm / 3/8 inch	772-109-00	1000 pieces
natural	medium 1.3 N (128 g)	3.2 mm / 1/8 inch	772-113-00	1000 pieces
natural	medium 1.3 N (128 g)	4.8 mm / 3/16 inch	772-114-00	1000 pieces
natural	medium 1.3 N (128 g)	6.4 mm / 1/4 inch	772-116-00	1000 pieces
natural	medium 1.3 N (128 g)	8.0 mm / 5/16 inch	772-117-00	1000 pieces
natural	medium 1.3 N (128 g)	9.5 mm / 3/8 inch	772-119-00	1000 pieces
natural	strong 1.8 N (184 g)	3.2 mm / 1/8 inch	772-123-00	1000 pieces
natural	strong 1.8 N (184 g)	4.8 mm / 3/16 inch	772-124-00	1000 pieces
natural	strong 1.8 N (184 g)	6.4 mm / 1/4 inch	772-125-00	1000 pieces
natural	strong 1.8 N (184 g)	8.0 mm / 5/16 inch	772-126-00	1000 pieces
natural	strong 1.8 N (184 g)	9.5 mm / 3/8 inch	772-127-00	1000 pieces

Intra-oral non-latex elastics

Highest elasticity and fatigue resistance.

The elastic packs are marked according to size and strength with symbols.

Every type can therefore be quickly identified.

Strength:

Strength ratings are based on the elastic being stretched 3 times the listed inner diameter.

Content:

30 x 100 pieces = 3000 pieces



Color	Traction force	Inner ø	REF	Quantity
transparent	medium 1.3 N (128 g)	3.2 mm / 1/8 inch	772-313-01	3000 pieces
transparent	medium 1.3 N (128 g)	4.8 mm / 3/16 inch	772-314-01	3000 pieces
transparent	medium 1.3 N (128 g)	6.4 mm / 1/4 inch	772-316-01	3000 pieces
transparent	medium 1.3 N (128 g)	8.0 mm / 5/16 inch	772-317-01	3000 pieces
transparent	strong 1.8 N (184 g)	3.2 mm / 1/8 inch	772-323-01	3000 pieces
transparent	strong 1.8 N (184 g)	4.8 mm / 3/16 inch	772-324-01	3000 pieces
transparent	strong 1.8 N (184 g)	6.4 mm / 1/4 inch	772-325-01	3000 pieces
transparent	strong 1.8 N (184 g)	8.0 mm / 5/16 inch	772-326-01	3000 pieces



Container for intra-oral elastics

For mini patient packs. With wall mount.

Dimensions (W x H x D): 508 mm x 169 mm x 107 mm

Color: white

REF 772-000-00

1 piece



Intra-oral latex elastics

Special sizes

Cuspid to cuspid, for facebows with cuspid hooks. Highest elasticity and fatigue resistance.

Strength ratings are based on the elastic being stretched 3 times the listed inner diameter.

Traction force	Inner ø	REF	Quantity
medium 1.3 N (128 g)	16 mm / 5/8 inch	772-100-00	1000 pieces
medium 1.3 N (128 g)	19 mm / 3/4 inch	772-200-00	1000 pieces




Elastics and auxiliaries / Intra-oral



Dentalastics® plastic ligatures

For ligating arch wires to brackets, instead of wire ligatures.

Color	Inner ø	Height	REF	Quantity
transparent	1.3 mm	0.9 mm	774-002-01	1000 pieces



Dentalastics® separators

For quick and easy separation of teeth. Radiopaque (ideal setting 60 kV and 3 mA).

Color	Inner ø	Indication	REF	Quantity
green	1.5 mm	for anteriors	774-100-01	1000 pieces
blue	2.1 mm	for posteriors	774-200-01	1000 pieces



SUPER Elasto-Force plastic chains

Consistent traction force over a long period of time. High elasticity with excellent elastic recovery.

Color	Inner ø	Length	Connector	REF	Quantity
transparent	1.3 mm	4.6 m	•	774-220-01	1 roll
pink	1.3 mm	4.6 m	•	774-231-01	1 roll
blue	1.3 mm	4.6 m	•	774-232-01	1 roll
red	1.3 mm	4.6 m	•	774-233-01	1 roll
silver	1.3 mm	4.6 m	•	774-235-01	1 roll
transparent	1.3 mm	4.6 m		774-221-01	1 roll
blue	1.3 mm	4.6 m		774-226-01	1 roll
silver	1.3 mm	4.6 m		774-229-01	1 roll
pearl	1.3 mm	4.6 m		774-236-01	1 roll



Dentalastics® Personal / assortment


Colored ligatures

Especially delicate plastic ligatures, instead of wire ligatures, for ligating archwires in brackets.

Content:

10 colors (without yellow or orange), 5 sticks each, total 50 sticks x 20 rings = 1000 pieces

Color	Inner ø	Height	REF	Quantity
10 colors	1.3 mm	0.7 mm	774-561-04	1 assortment



Dentalastics® Personal / Single packs

Colored ligatures

Especially delicate plastic ligatures, instead of wire ligatures, for ligating archwires in brackets.

Contents: 50 sticks x 20 rings = 1000 pieces

Color	Inner ø	Height	REF	Quantity
silver	1.3 mm	0.7 mm	774-549-01	1 pack
transparent	1.3 mm	0.7 mm	774-551-02	1 pack
black	1.3 mm	0.7 mm	774-552-01	1 pack
white	1.3 mm	0.7 mm	774-553-01	1 pack
pink	1.3 mm	0.7 mm	774-554-01	1 pack
blue	1.3 mm	0.7 mm	774-555-01	1 pack
red	1.3 mm	0.7 mm	774-556-01	1 pack
purple	1.3 mm	0.7 mm	774-557-01	1 pack
yellow	1.3 mm	0.7 mm	774-558-01	1 pack
green	1.3 mm	0.7 mm	774-559-01	1 pack
orange	1.3 mm	0.7 mm	774-560-01	1 pack
pearl	1.3 mm	0.7 mm	774-563-01	1 pack



Oral screens



Oral screen Ulmer model

Indications:

- incorrect functions of the tongue (pressing of the tongue, incorrect position of the tongue)
- biting the lips
- sucking the lips
- thumb sucking
- weak muscle structure of the lips

Advantages:

- no pressure points due to the flexible material
- more comfortable to wear
- better acceptance by the patient
- easy to monitor the performance of the myofunctional exercises

Special appliance containers see Supplies!

Size	REF	Quantity
Mini / Size 1 = for preschool age patients	074-000-00	1 piece
Maxi / Size 2 = for patients from school age to adults	074-001-00	1 piece



FaceFormer ONE blue with hygiene box

acc. to Berndsen

FaceFormer therapy targets the cause and is holistic. Regular FaceFormer training actively tones tissues, muscles and skin in the mouth and face. It strengthens and harmonizes muscle functions that are particularly important for chewing, swallowing, speaking, breathing, tooth position as well as facial expressions and tightening of the skin.

Indications:

- difficulty in swallowing (dysphagia)
- incorrect functions of the tongue
- dysfunctional mastication
- speech disorders with functional causes
- curing thumb/pacifier sucking, lip biting etc.
- paralysis of the oral and facial areas
- following an operation on the face, jaw or tongue
- malalignment of the teeth, anomalies in the jaw
- grinding the teeth (bruxism), temporomandibular problems
- sleep disorders, snoring, sleep apnea
- hearing loss, tinnitus
- dizziness, Menière's disease
- problems compensating pressure in the middle ear
- activation of nasal breathing, curing mouth breathing
- etc.

Color	Form	REF	Quantity
transparent blue	anatomical	074-101-01	1 piece



A series of 20 horizontal lines for writing, spanning the width of the page.




A series of 20 horizontal lines for writing, spanning the width of the page.

tomas® / Tool sets	Page	278
tomas® / Pins	Page	279
tomas® / Instruments and accessories	Page	281
amda® / Starter set	Page	290
amda® / Appliance and accessories	Page	291
Patient consultation material	Page	293





Anchorage systems.
innovative – comprehensive – efficient.



DENTAURUM
**ONLINE
SHOP**
shop.dentaurum.com





The complete system for all indications.

The concept.

The tomas® system has been one of the leading skeletal anchorage systems for many years. Users value its easy handling and high reliability. By integrating the new tomas®-pin EP and a selection of abutments for direct and indirect coupling, the system can be used in the buccal and the palatal region. Be it distalization, mesialization, intrusion, palatal expansion or indirect anchorage: tomas® has something to offer both the beginner and the professional in skeletal anchorage systems.



Indications.

The tomas®-pin SD and the tomas®-pin EP are used to provide temporary orthodontic anchorage for the following treatment examples:

- Distalization and mesialization of teeth.
- Uprighting of molars.
- Intrusion of teeth.
- Space closure in class I occlusion.
- Sliding mechanics in class II.
- Preventing the protrusion of incisors.
- Palatal expansion.
- Oligodontia.
- Changing the tooth position as part of preprosthetic treatment.
- Anchoring temporary crowns.


The advantages:

- Quick and easy to fit.
- Coordinated accessories.
- Effective around the clock – day and night.
- No patient compliance required.
- Perfect biomechanics.





advanced
molar
distalization
appliance.



advanced molar distalization appliance.

Together with the tomas® system, amda® offers a safe system for the distalization of upper molars. As a skeletal anchorage device it can be used for the bi- or unilateral distalization or the mesialization of upper molars and for the retraction of the upper anterior teeth depending on how it is fitted.

The appliance serves as a skeletally anchored transpalatal arch during the retention phase. amda® is more than merely a distalization appliance!

The amda® system was designed by Prof. Dr. med. dent. Moschos A. Papadopoulos. Dentaureum developed the system in close cooperation with him. Dr. Papadopoulos is a Professor and Chairman of the Department of Orthodontics at the School of Dentistry, Aristotle University in Thessaloniki, Greece.



Indications:

- Distalization of molars.
- Mesialization of molars.
- Retraction of anterior teeth.
- Combination of distalization and mesialization.
- Retention of molars.

The advantages:

- No patient compliance required.
- In situ – 24 hours.
- Hygienic and comfortable for the patient.
- No contamination of the spring.
- Easy to fit/Easy to remove.



Tool sets



tomas®-tool set XL

The tomas®-tool set XL contains all instruments for manual and mechanical insertion of a tomas®-pin.

REF 302-151-00	1 set
Contents:	
1 x tomas®-drill SD 1.1, short pilot drill	REF 302-103-00
1 x tomas®-guide driver	REF 302-080-00
1 x tomas®-driver	REF 302-004-50
1 x tomas®-applicator	REF 302-004-20
1 x tomas®-wheel, ø 10 mm	REF 302-004-30
1 x tomas®-applicator, long version	REF 302-004-70
1 x tomas®-torque ratchet	REF 302-004-41
1 x tomas®-screwdriver	REF 302-004-10
1 x tomas®-tool tray	REF 302-156-00



tomas®-tool set S

The tomas®-tool set S contains the instruments required to insert a tomas®-pin.

REF 302-251-00	1 set
Contents:	
1 x tomas®-drill SD 1.1, short pilot drill	REF 302-103-00
1 x tomas®-guide driver	REF 302-080-00
1 x tomas®-screwdriver	REF 302-004-10
1 x tomas®-tool tray	REF 302-156-00

tomas®-pin SD
self-drilling
sterile

Titanium grade 5
Thread $\varnothing$: 1.6 mm
Core $\varnothing$: 1.2 mm
Head type: cross-slot

Length	Color-coded packaging	REF	Quantity
6 mm	yellow	302-106-00	1 piece
8 mm	blue	302-108-00	1 piece
10 mm	green	302-110-00	1 piece



tomas®-pin EP
self-drilling
sterile


Titanium grade 5
Thread $\varnothing$: 1.6 mm
Core $\varnothing$: 1.2 mm
Head type: mushroom-shaped

Length	Color-coded packaging	REF	Quantity
6 mm	white	302-206-00	1 piece
8 mm	grey	302-208-00	1 piece
10 mm	black	302-210-00	1 piece



ISO 19023

	SD / SD N	EP / EP N
d₁	1.6 mm	1.6 mm
d₂	1.2 mm	1.2 mm
d₃	2.8 mm	2.8 mm
k₁	2.2 mm	2.9 mm
k₂	2.0 mm	2.0 mm
l₁	10.2 mm / 12.2 mm / 14.2 mm	10.9 mm / 12.9 mm / 14.9 mm
l₂	6.0 mm / 8.0 mm / 10.0 mm	6.0 mm / 8.0 mm / 10.0 mm



Pins, non-sterile



tomas®-pin SD N
self-drilling
non-sterile
sterilizable

Titanium grade 5
Thread ø: 1.6 mm
Core ø: 1.2 mm
Head type: cross-slot

Length	REF	Quantity	REF	Quantity
6 mm	302-106-10	1 piece	302-106-30	10 pieces
8 mm	302-108-10	1 piece	302-108-30	10 pieces
10 mm	302-110-10	1 piece	302-110-30	10 pieces



tomas®-pin EP N
self-drilling
non-sterile
sterilizable

Titanium grade 5
Thread ø: 1.6 mm
Core ø: 1.2 mm
Head type: mushroom-shaped

Length	REF	Quantity	REF	Quantity
6 mm	302-206-10	1 piece	302-206-30	10 pieces
8 mm	302-208-10	1 piece	302-208-30	10 pieces
10 mm	302-210-10	1 piece	302-210-30	10 pieces

tomas®-drill SD 1.1, short pilot drill

sterilizable

The drilling must be undertaken at 800 min^{-1} (max. 1500 min^{-1}) with a sterile, physiological saline solution.
For perforating the compact bone ($\geq 2.0 \text{ mm}$) when using tomas®-pin SD and EP.



Length	ø	Color-coded packaging	REF	Quantity
4 mm	1.1 mm	yellow	302-103-00	1 piece

tomas®-applicator

sterilizable

For manual insertion and explantation of tomas®-pin.
Long version especially suitable for palatal insertion.

Length	REF	Quantity
18 mm	302-004-20	1 piece
28 mm	302-004-70	1 piece



tomas®-wheel

sterilizable

(nut)

For manual insertion and explantation of the tomas®-pin in combination with the tomas®-applicator (REF 302-004-20 / 302-004-70).

ø: 10.0 mm

REF 302-004-30	1 piece
----------------	---------



Anchorage systems

Instruments and accessories



tomas®-torque ratchet

For manual insertion and explantation of the tomas®-pin in combination with the tomas®-applicator (REF 302-004-20 / 302-004-70); torque should be limited to 20 Ncm.

REF 302-004-41

1 piece



tomas®-driver sterilizable

(mechanical screwing aid)

For mechanical insertion of the tomas®-pin with the aid of a lockable torque instrument. The setting should be limited to 20 Ncm and 25 min⁻¹.

Length: 22 mm

REF 302-004-50

1 piece



tomas®-guide driver sterilizable

(power-driven insertion aid)

Long version.

Handpiece-held insertion instrument coordinated with tomas®-guide tubes. For the power-driven insertion of the tomas®-pin with the help of a contra-angle handpiece with restricted torque. Torque restriction should be set at 20 Ncm and 25 min⁻¹. Highly suitable for palatal insertion of the tomas®-pin.

Length: 27 mm

REF 302-080-00

1 piece



Manually turned unit mod. acc. to Pauls, with adjustable torque from 0 - 40 Ncm sterilizable

(Manually turned unit)

Manually turned unit with adjustable torque from 0 Ncm - 40 Ncm. Torque should not exceed 20 Ncm.

Can be fitted to a contra-angle handpiece as cost-effective alternative to a surgical tool or electrically-powered equipment.

REF 302-004-81

1 piece



tomas®-screwdriver sterilizable

For manual insertion and explantation of tomas®-pin.

Contains tomas®-screwdriver key 1.3 (hex key REF 302-004-16) and tomas®-screwdriver key 3.0 (REF 302-004-14).

REF 302-004-10

1 piece



tomas®-tool tray sterilizable

For storage of tomas®-components.

REF 302-156-00

1 piece

tomas®-abutment kit

(set of coupling elements)

The tomas®-abutment kit contains abutments for the tomas®-pin SD and the tomas®-pin EP.

For skeletal anchorage of appliances for mesialization, distalization and intrusion of teeth and for maxillary expansion. Appliances can be made in the dental practice or in a laboratory.

REF 302-050-00	1 set
Contents:	
6 x tomas®-abutment universal	REF 302-025-00
6 x tomas®-abutment tube 1.1	REF 302-025-11
2 x tomas®-abutment tube 1.5	REF 302-025-15
2 x tomas®-abutment tube square 22	REF 302-025-22
2 x tomas®-abutment wire 6	REF 302-025-06
4 x tomas®-abutment wire 12	REF 302-025-12
2 x tomas®-abutment plain	REF 302-026-00
2 x tomas®-abutment EP	REF 302-027-00
50 x tomas®-transfer cap	REF 302-028-01
10 x tomas®-laboratory pin EP	REF 302-029-01
10 x tomas®-laboratory pin SD	REF 302-030-01
1 x tomas®-box	REF 302-023-10



tomas®-abutment universal

(coupling element)

Indirect coupling between tomas®-pin and conventional multi-bracket systems over one or more teeth. Direct coupling to skeletally-anchored appliances. Temporary fixation of wires with 1.1 mm ϕ by means of snap effect (only in connection with tomas®-pin SD!).

REF 302-025-00	2 pieces
----------------	----------



tomas®-abutment tube 1.1

(coupling element with tube)

Abutment with rotatable tube for direct coupling between tomas®-pin and various appliances, e.g. for bi- or unilateral mesialization, distalization, intrusion etc.

The tomas®-abutment wire 6 (REF 302-025-06) and tomas®-abutment wire 12 (REF 302-025-12) can be used in addition.

Inner ϕ tube	Width tube	REF	Quantity
1.15 mm / 45	3.0 mm	302-025-11	1 piece



tomas®-abutment tube 1.5

(coupling element with tube)

Abutment with rotatable tube for direct coupling between the tomas®-pin and all hyrax® RPE screws.

Can be used for hybrid RPE appliance or skeletal RPE appliance using Pin-First® approach.

Inner ϕ tube	Width tube	REF	Quantity
1.55 mm / 61	3.0 mm	302-025-15	1 piece



tomas®-RPE eyelet

(coupling element for RPE screws)

Eyelet for direct coupling between tomas®-pin with all hyrax® RPE screws with two supports for easier placement on working model.

Can be used for hybrid RPE appliance or skeletal RPE appliance using Pin-Last® approach.

Outer ϕ eyelet	Inner ϕ eyelet	Thickness eyelet	REF	Quantity
5.0 mm	top 3.0 mm bottom 2.3 mm	1.5 mm	302-602-00	12 pieces + 2 holders



Instruments and accessories



tomas®-abutment tube square 22

(coupling element with tube)

Abutment with rotatable square tube for direct coupling between the tomas®-pin and various appliances, e.g. for intrusion etc.

Inner $\square$ tube	Width tube	REF	Quantity
0.56 mm x 0.70 mm / 22 x 28	4.0 mm	302-025-22	1 piece



tomas®-abutment double tube 1.5

Concept by Dr. Sohani

(coupling element with 2 tubes)

Abutment with rotatable double tube for direct coupling between tomas®-pin and various appliances, e.g. for bi- or unilateral mesialization, distalization, intrusion etc. The tomas®-abutment U-wire (REF 302-025-21) can be used in addition.

Inner ϕ tube	Width tube	REF	Quantity
1.55 mm / 61	2 x 3.0 mm	302-025-30	1 piece



tomas®-abutment U-wire

(coupling element with wire)

Abutment with rotatable bilateral wire for direct coupling between the tomas®-pin and various appliances, e.g. for bi- or unilateral mesialization, distalization etc. The tomas®-abutment double tube 1.5 (REF 302-025-30) can be used in addition.

Wire ϕ	Wire length	REF	Quantity
1.10 mm / 43	2 x 80 mm	302-025-21	1 piece



tomas®-abutments median

Combination of tomas®-abutment U-wire (REF 302-025-21) and tomas®-abutment double tube 1.5 (REF 302-025-30). Combination is intended for two tomas®-pins which are placed sagittally to the mid-palatal suture. The distance between the tomas®-pins must be at least 6 mm.

Inner ϕ tube	Width tube	Wire ϕ	Wire length	REF	Quantity
1.55 mm / 61	2 x 3.0 mm	1.10 mm / 43	2 x 80 mm	302-025-51	1 piece



tomas®-abutment wire 6

(coupling element with wire)

Abutment with rotatable unilateral wire for direct coupling between the tomas®-pin and various appliances, e.g. for bi- or unilateral mesialization, distalization etc. The tomas®-abutment tube 1.1 (REF 302-025-11) can be used in addition.

Wire ϕ	Wire length	REF	Quantity
1.10 mm / 43	60 mm	302-025-06	1 piece

tomas®-abutment wire 12

(coupling element with wire)

Abutment with rotatable unilateral wire for direct coupling between the tomas®-pin and various appliances, e.g. for bi- or unilateral mesialization, distalization etc. The tomas®-abutment tube 1.1 (REF 302-025-11) can be used in addition.

Wire ø	Wire length	REF	Quantity
1.10 mm / 43	120 mm	302-025-12	1 piece



tomas®-abutment plain

(coupling element)

For welding individual elements and for anchoring temporary crowns.

REF 302-026-00	1 piece
-----------------------	---------



tomas®-abutment EP

(coupling element)

Direct coupling between the tomas®-pin EP and distalization appliances.

Inner ø tube	Width tube	Inner ø eyelet	Thickness eyelet	REF	Quantity
1.05 mm / 41	2.0 mm	2.4 mm / 94	0.8 mm	302-027-00	2 pieces



tomas®-transfer cap

Transfer of oral status to a dental cast.

REF 302-028-01	2 pieces
-----------------------	----------



tomas®-laboratory pin EP

Laboratory analog for the tomas®-pin EP for manufacture of orthodontic appliances with skeletal anchorage in a laboratory.

REF 302-029-01	1 piece
-----------------------	---------



tomas®-laboratory pin SD

Laboratory analog for the tomas®-pin SD for manufacture of orthodontic appliances with skeletal anchorage in a laboratory.

REF 302-030-01	1 piece
-----------------------	---------



Instruments and accessories



tomas®-auxiliary kit

(set of accessories)

The tomas®-auxiliary kit contains accessories for mesializing, distalizing and uprighting of teeth in connection with the tomas®-pin.

REF 302-001-20

1 set

Contents:

10 x tomas®-coil spring, light (8 mm)	REF 302-012-00
10 x tomas®-coil spring, medium (8 mm)	REF 302-012-10
10 x tomas®-coil spring, heavy (8 mm)	REF 302-012-20
10 x tomas®-uprighting spring	REF 302-009-00
3 x rematitan® LITE compression spring	REF 758-361-00
10 x tomas®-stop screw	REF 302-013-01
1 x hexagon socket key 0.9	REF 607-129-00
10 x tomas®-T-wire	REF 302-024-00
10 x tomas®-power arm, square	REF 302-015-00
10 x tomas®-power arm, round	REF 302-019-00
10 x tomas®-crimp hook, Upper left / Lower right	REF 302-020-00
10 x tomas®-crimp hook, Upper right / Lower left	REF 302-020-10
10 x tomas®-aligner hook, Upper left / Lower right	REF 302-021-00
10 x tomas®-aligner hook, Upper right / Lower left	REF 302-021-10
1 x tomas®-box	REF 302-023-10



tomas®-coil spring

Nickel titanium tension springs with eyelets for linking to the tomas®-pin and hooks.

Length: 8 mm (measured from the center of the two eyelets)

Inner ø: 0.90 mm / 36

Strength category: super elastic

Traction force	Wire strength ø	Color coded packaging	REF	Quantity
light 1.38 N (140 g)	0.21 mm / 8	blue	302-012-00	10 pieces
medium 2.01 N (200 g)	0.23 mm / 9	yellow	302-012-10	10 pieces
heavy 2.74 N (270 g)	0.23 mm / 9	red	302-012-20	10 pieces



SUPER Elasto-Force plastic chains

Consistent traction force over a long period of time. High elasticity with excellent elastic recovery.

Color	Inner ø	Length	REF	Quantity
transparent	1.3 mm	4.6 m	774-221-01	1 roll



tomas®-uprighting spring

For uprighting and simultaneous intrusion or extrusion of molars.

REF 302-009-00

10 pieces

rematitan® LITE compression spring

Advantages:


- high elasticity with virtually constant restoring force
- fatigue-free, low, uniform transfer of forces

May be individually shortened.

Material: nickel titanium

Strength category: super elastic

Inner ø	Length	Wire strength ø	Force	REF	Quantity
0.76 mm / 30	18 cm	0.30 mm / 12	medium	758-361-00	3 pieces



tomas®-stop screw

Without hexagon socket key.

REF 302-013-01	10 pieces
-----------------------	-----------



Hexagon socket key 0.9 sterilizable

Key width: 0.9 mm

REF 607-129-00	1 piece
-----------------------	---------



tomas®-slotted stop

Suitable as stop for tension springs.

☑	REF	Quantity
0.46 mm x 0.64 mm / 18 x 25	302-022-18	100 pieces
0.56 mm x 0.70 mm / 22 x 28	302-022-22	100 pieces



tomas®-hook

For fixing elastic elements (plastic chains, springs)
to the tomas®-pin SD head.

REF 302-009-10	10 pieces
-----------------------	-----------



tomas®-T-wire

Made of remanium® for direct and indirect connection with the tomas®-pin SD.

Material: stainless steel

Strength category: spring hard · 1800 N/mm² - 2000 N/mm²

☑	REF	Quantity
0.55 mm x 0.70 mm / 21.5 x 28	302-024-00	10 pieces



Instruments and accessories



tomas®-cross tube 18

Transverse connection between two small rectangular tubes for joining two wires without welding. Ideal for the connection between tomas®-pin SD and arch.

Horizontal tube:

Overall length: 4.0 mm

Inner $\square$: 0.46 mm x 0.64 mm / 18 x 25

Vertical tube:

Overall length: 2.0 mm

Inner $\square$: 0.56 mm x 0.70 mm / 22 x 28

REF 302-014-18

10 pieces



tomas®-cross tube 22

Transverse connection between two small rectangular tubes for joining two wires without welding. Ideal for the connection between tomas®-pin SD and arch.

Horizontal tube:

Overall length: 4.0 mm

Inner $\square$: 0.56 mm x 0.70 mm / 22 x 28

Vertical tube:

Overall length: 2.0 mm

Inner $\square$: 0.56 mm x 0.70 mm / 22 x 28

REF 302-014-22

10 pieces



tomas®-power arm, square

Strength category: spring hard · 1800 N/mm² - 2000 N/mm²

Crimping tube slotted, with spring hard wire, square

($\square$ 0.55 mm x 0.70 mm / 21.5 x 28). Can be placed directly on the arch wire, without loss of ligation of the archwire.

Universal application for the slot 18/22 technique. Ideal for coupling with the tomas®-pin SD and for bending individual hooks.

Dimensions:

Tube length: 2.0 mm

Wire length: 25.0 mm

REF 302-015-00

10 pieces



tomas®-power arm, round

Strength category: spring hard · 1800 N/mm² - 2000 N/mm²

Crimping tube slotted, with spring hard wire, round

($\varnothing$ 0.7 mm / 28). Can be placed directly on the arch wire, without loss of ligation of the archwire. Universal application for the slot 18/22 technique. For bending individual hooks.

Dimensions:

Tube length: 2.0 mm

Wire length: 25.0 mm

REF 302-019-00

10 pieces

Crimping pliers

sterilizable

Tip with V-shaped notch for holding and crimping sliding/clamping and stop tubes.

Premium-Line · 10 year guarantee

REF 003-711-00

1 piece



tomas®-crimp hook

Crimping hook slotted, with 4 hooks to engage springs and plastic chains.

Universal application for the slot 18/22 technique.

Strength category: spring hard · 1800 N/mm² - 2000 N/mm²



Quadrant	REF	Quantity
upper left / lower right	302-020-00	10 pieces
upper right / lower left	302-020-10	10 pieces

tomas®-aligner hook

(bonding hook)

Hook for bonding on teeth or aligner for direct coupling with tomas®-pin.

Special adhesive (Bond Aligner™ from Reliance) is required to bond onto plastic of aligner.

Strength category: spring hard · 1800 N/mm² - 2000 N/mm²



Quadrant	REF	Quantity
upper left / lower right	302-021-00	10 pieces
upper right / lower left	302-021-10	10 pieces

tomas®-box

Plastic box with 2 inserts.

REF 302-023-10

1 piece



amda® / Starter set



amda® starter set (assembled)

Skeletally anchored appliance for the distalization of molars with the tomas®-pin EP.



REF 607-180-00	1 set
Contents:	
1 x amda® telescope (2 pieces)	REF 607-181-02
1 x amda® connector, pair 7 mm	REF 607-184-07
1 x amda® connector, pair 10 mm	REF 607-184-10
1 x amda® connector, pair 13 mm	REF 607-184-13
1 x amda® palatal arch	REF 607-185-01
1 x hexagon socket key 0.9	REF 607-129-00

amda® telescope

Encapsulated stainless steel spring with two stop screws
(REF 302-013-01).

Spring force: max. 4.7 N

Spring travel: max. 4.0 mm

REF 607-181-02	2 pieces
REF 607-181-10	10 pieces



amda® connector, left

Elements to connect the amda® appliance and the palatal sheath on the molars.

Height	REF	Quantity
7 mm	607-182-07	10 pieces
10 mm	607-182-10	10 pieces
13 mm	607-182-13	10 pieces



amda® connector, right

Elements to connect the amda® appliance and the palatal sheath on the molars.

Height	REF	Quantity
7 mm	607-183-07	10 pieces
10 mm	607-183-10	10 pieces
13 mm	607-183-13	10 pieces



amda® connector, pair

Elements to connect the amda® appliance and the palatal sheath on the molars.

Height	REF	Quantity
7 mm	607-184-07	1 pair
10 mm	607-184-10	1 pair
13 mm	607-184-13	1 pair



amda® / Appliance and accessories



amda® palatal arch

(incl. 2 x tomas®-abutments EP, REF 302-027-00).

Palatal arch for amda® for fixing the amda® telescopes.

REF 607-185-01	1 piece
REF 607-185-10	10 pieces



Hexagon socket key 0.9

sterilizable

Key width: 0.9 mm

REF 607-129-00	1 piece
----------------	---------

tomas®-demonstration pin SD 08 self-drilling

For patient demonstration.

Scale 5 : 1

REF 302-007-20

1 piece



tomas®-demonstration pin EP 08 self-drilling

For patient demonstration.

Scale 5 : 1

REF 302-007-30

1 piece



Orthodontic demonstration model

tomas® / amda®

Fitted with:

- tomas®-pin EP
- amda®
- lingual / palatal sheaths easy in,
welded to standard molar bands

Maxilla made of plastic.

REF 892-000-26

1 piece



Impression-taking	Page	298
Model fabrication	Page	299
Appliance containers and accessories	Page	301
Orthodontic study and demonstration models	Page	303
Other	Page	308



Supplies.

Reliability at great value.



DENTAURUM
**ONLINE
SHOP**
shop.dentaurum.com





Professional materials for the orthodontic practice.

Our wide product range for supplies ranges from impression materials to orthodontic tooth brushes. Every question from your patient can be answered with the help of our demonstration models and orthodontic study models.

Whether mini, maxi or magnum, your patients will be delighted with our appliance containers.

With the Ortho Brush, it is easy to clean multi-bracket appliances reliably and quickly.



Orthalgenat® Duo

The universal orthodontic impression material.

Flexaponal®

Special wax according to Uhlig.

Orthalgenat® Ex

The perfect cleaning material for impression trays.

Rapidur®

White, high-strength synthetic dental stone.

Orthodontic study and demonstration models

Wide range of models.

Appliance container

In all sizes and colors.

Impression-taking



Orthalgenat® Duo

Impression material incl. measuring set.

Flavor	Color	REF	Quantity
mint	green	150-501-20	1 kg
mint	green	150-509-20	9 kg
mint	green	150-520-20	20 kg

Packaging: 1 kg (box), 9 x 1 kg and 20 x 1 kg (bag)
DIN EN ISO 21563



Dosing bottle

For liquids.

Capacity	REF	Quantity
500 ml	150-000-00	1 piece



Orthalgenat® Ex

Impression tray cleaner. For removing alginate and plaster residue. Concentrate for manual or ultrasonic cleaning.

REF 150-010-00	1.5 kg
----------------	--------

Flexaponal®, special wax

acc. to Uhlig

- for bite registration (as bite wafer or bite rim)
- for adhesive borders on impression trays
- soft, flexible and adhesive

Thickness: 1.0 mm



Dimensions (L x W)	Color	REF	Quantity
120 mm x 15 mm	red	124-003-00	50 pieces
200 mm x 10 mm	red	124-002-00	50 pieces
200 mm x 10 mm	white	124-202-00	50 pieces

Rapidur®

Synthetic dental stone

- extra hard, low expansion
- quick-setting

Application:

suitable for all areas of orthodontics and dental prosthetics

Norm: DIN EN ISO 6873

Class 4

Setting expansion	max. 0.15 %
Working time	min. 3 min
Compression resistance	min. 35 MPa
Setting time	max. 30 min
Color	white

REF 164-450-00

6 kg

**Mixing bowl**

For plaster, impression and investment materials, easy to clean.

Capacity	REF	Quantity
400 ml	164-651-00	1 piece
600 ml	164-650-00	1 piece

**Model gloss**

Environmentally friendly impregnation liquid for plaster demonstration models (e.g. Rapidur®, REF 164-450-00) for shiny, dust-repellent surfaces.

REF 166-450-00

4.5 kg

**Ultrasonic cleaning solution**

To remove	REF	Quantity
Polishing paste, wax, grease	094-610-00	1000 ml
Plaster and plaster-bonded investment material	094-611-00	1000 ml
Calculus, cement, etching residues	094-612-00	1000 ml



Model fabrication



Orthodontic base unit

acc. to Hinz

For fixing three-dimensionally oriented models. The uniform height (6 cm) of the orthodontic models is determined by the base unit, which is easy to use.

REF 050-500-00

1 piece



Base form with hinge

Made of impact-resistant plastic. Eliminates time consuming plaster model preparation. Suitable for orthodontic base unit acc. to Hinz.

Material	Color	REF	Quantity
Plastic	white	050-501-00	100 pairs



Bar with hinge

For base form (REF 050-501-00).

Suitable for orthodontic base unit acc. to Hinz.

Material	Color	REF	Quantity
Plastic	white	050-503-01	200 pieces

NEW



Orthodontic model tray

With colored labels.

For 5 sets of models, stackable.

Dimensions (W x H x D): 94 mm x 83 mm x 330 mm

Material: Plastic

Color: white

REF 050-502-00

10 pieces

Appliance containers and accessories

Appliance container, mini

The mini version is especially suitable for single appliances and oral screens (e.g. the Ulmer model).

Material: Plastic (dimensionally stable and impact resistant)

Dimensions (W x H x D): 7.00 cm x 3.30 cm x 4.80 cm

Size	Color	REF	Quantity
mini	pearl blue	163-602-01	10 pieces
mini	turquoise	163-610-01	10 pieces
mini	neon coral	163-611-01	10 pieces
mini	neon green	163-612-01	10 pieces

**Appliance container, maxi**

The maxi version is especially suitable for bimaxillary appliances (e.g. activator, bionator).

Material: Plastic (dimensionally stable and impact resistant)

Dimensions (W x H x D): 7.80 cm x 4.50 cm x 5.60 cm

Size	Color	REF	Quantity
maxi	ruby-red	163-701-01	10 pieces
maxi	pearl blue	163-702-01	10 pieces
maxi	pink	163-703-01	10 pieces
maxi	black	163-707-01	10 pieces
maxi	turquoise	163-710-01	10 pieces
maxi	neon coral	163-711-01	10 pieces
maxi	neon green	163-712-01	10 pieces
maxi	neon yellow	163-713-01	10 pieces

**Appliance containers, magnum**

The magnum version is especially suitable for extremely large appliances (e.g. Fraenkel appliances, orthodontic appliances for adult treatment, or dentures).

Material: Plastic (dimensionally stable and impact resistant)

Dimensions (W x H x D): 6.80 cm x 4.70 cm x 7.80 cm

Size	Color	REF	Quantity
magnum	ruby-red	163-801-01	10 pieces
magnum	pearl blue	163-802-01	10 pieces



Appliance containers and accessories



Appliance container, football

Universal size for removable upper appliances, bimaxillary appliances and function regulators incl. snap hook.
Material: Plastic (dimensionally stable and impact resistant)

Color	REF	Quantity
black-white	163-805-01	10 pieces



Snap hook

For attaching the appliance container to a belt loop.

REF 163-801-00	10 pieces
----------------	-----------



Cord with safety lock

For appliance containers, printed with Dentaaurum.

Color	REF	Quantity
a variety of colors	163-807-00	25 pieces

Orthodontic study and demonstration models

Typodont / Dentoform

Study and test device for fixed appliances. The unit mechanically corresponds to a simple articulator. Depending upon the task, the metal/plastic teeth are set in the wax sockets, banded or bonded with brackets and the arch wire fitted.

Availability:

Unit without teeth and accessories.

REF 050-801-00

1 piece

**Metal/Plastic teeth**

For Typodont / Dentoform unit (REF 050-801-00).

Tooth crowns made of plastic, tooth roots made of metal.

REF 050-803-01

28 pieces

**Wax bases**

For Typodont / Dentoform unit (REF 050-801-00).

Ideal occlusion (maxillary and mandibular incl. 2 wax sheets).

REF 050-802-00

1 pair

**Orthodontic study model****Angle Class I, normal occlusion**

Maxilla and mandible made of plastic.

REF 892-001-00

1 piece

**Orthodontic study model****Angle Class II/1, distal occlusion with deep bite**

Maxilla and mandible made of plastic.

REF 892-002-01

1 piece



Orthodontic study and demonstration models



Orthodontic study model

Angle Class II/2, distal occlusion with deep overbite
Maxilla and mandible made of plastic.

REF 892-002-02

1 piece



Orthodontic study model

Angle Class III, mesial occlusion, prognathic dentition
Maxilla and mandible made of plastic.

REF 892-003-00

1 piece



Orthodontic demonstration model discovery®

Angle Class I, normal occlusion

Roth 22

Fitted with:

- discovery® metal brackets
- Ortho-Cast M-Series buccal tubes and M-Series mini buccal tubes
- dentaform® bands
- rematitan® LITE arches

Maxilla and mandible made of plastic, hinged.

REF 892-000-07

1 piece



Orthodontic demonstration model equilibrium® mini / equilibrium® 2

Angle Class I, normal occlusion

Roth 22

Fitted with:

- equilibrium® mini and equilibrium® 2 metal brackets
- Ortho-Cast M-Series mini buccal tubes and M-Series mini buccal tubes
- dentaform® bands
- rematitan® LITE arches

Maxilla and mandible made of plastic, hinged.

REF 892-000-14

1 piece

Orthodontic study and demonstration models

**Orthodontic demonstration model
discovery® delight****Fitted with:**

- discovery® delight lingual brackets
- discovery® delight molar tubes
- rematitan® LITE lingual arch
- dentaflex® lingual arch

Maxilla and mandible made of plastic, hinged.



REF 892-000-22

1 piece

**Orthodontic demonstration model
discovery® smart****Angle Class I, normal occlusion**

McLaughlin-Bennett-Trevisi* 22

Fitted with:

- discovery® smart metal brackets
- Ortho-Cast M-Series and M-Series mini buccal tubes
- rematitan® LITE arches

Maxilla and mandible made of plastic, hinged.

REF 892-000-20

1 piece

* The Dentaurum version of this prescription is not claimed to be a duplication of any other, nor does Dentaurum imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.

**Orthodontic demonstration model
discovery® pearl**

McLaughlin-Bennett-Trevisi* 22

Fitted with:

- discovery® pearl ceramic brackets
- translucent arch pearl
- Dentalastics® pearl plastic ligatures

Maxilla made of plastic.

REF 892-000-18

1 piece

* The Dentaurum version of this prescription is not claimed to be a duplication of any other, nor does Dentaurum imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.

**Orthodontic demonstration model
discovery® pearl**

Roth 22

Fitted with:

- discovery® pearl ceramic brackets, maxilla / mandible 5-5
- Ortho-Cast M-Series buccal tubes, maxilla / mandible 6
- translucent arch pearl, maxilla
- Tensic® arch, mandible
- Dentalastics® pearl plastic ligatures

Maxilla and mandible made of plastic, hinged.

REF 892-000-21

1 piece



Orthodontic study and demonstration models



Orthodontic demonstration model discovery® sl

Angle Class I, normal occlusion

Roth 22

Fitted with:

- discovery® sl self-ligating metal brackets
- Ortho-Cast M-Series mini buccal tubes
- Equire stainless steel arch, maxilla
- rematitan® sl arch, mandible

Maxilla and mandible made of plastic, hinged.

REF 892-000-17

1 piece



Orthodontic demonstration model discovery® sl 2.0

Angle Class I, normal occlusion

McLaughlin-Bennett-Trevisi* 22

Fitted with:

- discovery® sl 2.0 self-ligating metal brackets
- Ortho-Cast M-Series and M-Series mini buccal tubes
- Tensic® arches

Maxilla and mandible made of plastic, hinged.

REF 892-000-19

1 piece

* The Dentaaurum version of this prescription is not claimed to be a duplication of any other, nor does Dentaaurum imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.



Orthodontic demonstration model prime4me® RETAIN3R

For patient demonstration.

Fitted with:

- prime4me® base model
- prime4me® RETAIN3R upper jaw
- transfer tray for prime4me® RETAIN3R

REF 892-000-29

1 piece

NEW



Orthodontic demonstration model tomas® / amda®

Fitted with:

- tomas®-pin EP
- amda®
- lingual / palatal sheaths easy in, welded to standard molar bands

Maxilla made of plastic.

REF 892-000-26

1 piece



Orthodontic study and demonstration models

tomas®-demonstration pin SD 08 self-drilling

For patient demonstration.

Scale 5 : 1

REF 302-007-20

1 piece

**tomas®-demonstration pin EP 08 self-drilling**

For patient demonstration.

Scale 5 : 1

REF 302-007-30

1 piece



Other



Ortho Brush
Orthodontic tooth brush.
Ideal for cleaning interproximal areas,
fixed brackets and band appliances.


REF 163-300-00	10 pieces
----------------	-----------

Standard expansion screws	Page	314
Special expansion screws	Page	317
Variety expansion screws	Page	330
RPE screws	Page	334
Accessories	Page	340



Expansion screws.

Perfect in form, function and design.



DENTAURUM
**ONLINE
SHOP**
shop.dentaurum.com



Perfect in form, function and design.



A wide range of products.

Dentaurum's selection of stainless steel expansion screws ranges from standard expansion screws to more complex special expansion screws. We have extended our range of special expansion screws with the Neo-bite jumping appliance acc. to F.M. Sander* – a modern, patented three-dimensional Bertoni expansion screw with guide ridges. Titan hyrax® and rematitan® expansion screws, the biocompatible alternatives in pure titanium, round off this wide selection.

Perfect in every detail.

Dentaurum's expansion screws achieve wide expansion despite their small size. The screws are guaranteed to open evenly and not reverse inadvertently thanks to the friction on the thread which has been carefully calculated, and an integrated stop mechanism. High precision manufacture and excellent quality ensure that the screw construction is strong and stable. Additional holes in the body ensure optimal retention in the acrylic.



Dentaurum expansion screws

The right solution for every type of treatment.

Neo-bite jumping appliance acc. to F.M. Sander*

Latest generation of the Sander appliance.

hyrax®

Special expansion screw for palatal expansion.

hyrax® click

The RPE screw with special feature
to prevent unwanted reversing.

hyrax® neo

The alternative to the classic hyrax® screw.

Titan hyrax®

The world's first RPE screw in titanium.

Variety

The versatile, multi-functional expansion screw.

Variety click

The special expansion screw with special feature
to prevent unwanted reversing.

Expansion screws

Standard expansion screws



Expansion screw Mini

Material: stainless steel


- with **cone-shaped or wide** expansion screw holder
 - skeleton design
 - holes in body of screws for better retention in acrylic
 - built-in stops prevent screws from coming apart once maximum expansion has been achieved
- With expansion screw key, standard (straight, small).

Indication: for transverse expansion and distalization in the mandible or for especially narrow maxillae

Directional arrow: Plastic, white

Maximum expansion	4 x 1/4 turn	REF (cone-shaped)	REF (wide)	Quantity
6.5 mm	= 0.7 mm	600-300-10	600-320-10	1 piece
6.5 mm	= 0.7 mm	600-300-30	600-320-30	10 pieces
6.5 mm	= 0.7 mm	600-300-60	600-320-60	100 pieces

→ 10 mm →



Expansion screw Medium

Material: stainless steel


- with **cone-shaped or wide** expansion screw holder
 - skeleton design
 - holes in body of screws for better retention in acrylic
 - built-in stops prevent screws from coming apart once maximum expansion has been achieved
- With expansion screw key straight, standard (straight, small).

Indication: for transverse expansion and distalization in maxilla and mandible

Directional arrow: Plastic, neon yellow

Maximum expansion	4 x 1/4 turn	REF (cone-shaped)	REF (wide)	Quantity
7.5 mm	= 0.8 mm	600-301-10	600-321-10	1 piece
7.5 mm	= 0.8 mm	600-301-30	600-321-30	10 pieces
7.5 mm	= 0.8 mm	600-301-60	600-321-60	100 pieces

→ 11 mm →



Expansion screw Maxi

Material: stainless steel

- with **cone-shaped** or **wide** expansion screw holder
- skeleton design
- holes in body of screws for better retention in acrylic
- built-in stops prevent screws from coming apart once maximum expansion has been achieved

With expansion screw key straight, standard (straight, small).


Indication: for transverse expansion in the maxilla

Directional arrow: Plastic, neon yellow

Maximum expansion	4 x 1/4 turn	REF (cone-shaped)	REF (wide)	Quantity
7.5 mm	= 0.8 mm	600-302-10	600-322-10	1 piece
7.5 mm	= 0.8 mm	600-302-30	600-322-30	10 pieces
7.5 mm	= 0.8 mm	600-302-60	600-322-60	100 pieces



11 mm



Expansion screw Magnum


Material: stainless steel

- with **cone-shaped** expansion screw holder
 - skeleton design
 - version: extra strong, with larger screw body
 - holes in body of screws for better retention in acrylic
 - built-in stops prevent screws from coming apart
- With expansion screw key straight, standard (straight, small).



Indication: for transverse expansion in the maxilla and bimaxillary functional gnatho-orthopedic appliances

Directional arrow: Plastic, neon yellow

Maximum expansion	4 x 1/4 turn	REF	Quantity
11.0 mm	= 0.8 mm	600-303-10	1 piece
11.0 mm	= 0.8 mm	600-303-30	10 pieces
11.0 mm	= 0.8 mm	600-303-60	100 pieces



14,5 mm



Expansion screws

Standard expansion screws



rematitan® expansion screw Medium

Material: titanium

- with **cone-shaped** expansion screw holder
- skeleton design
- holes in body of screws for better retention in acrylic
- built-in stops prevent screws from coming apart


With expansion screw key, standard (straight, small).

Indication: for transverse expansion and distalization in maxilla and mandible

Directional arrow: Plastic, pearl-blue

Maximum expansion	4 x 1/4 turn	REF	Quantity
7.5 mm	= 0.8 mm	600-400-10	1 piece
7.5 mm	= 0.8 mm	600-400-30	10 pieces

→ 11 mm →



Trapezoidal expansion screw

Material: stainless steel

- with **cone-shaped** expansion screw holder
- skeleton design
- holes in body of screws for better retention in acrylic
- built-in stops prevent screws from coming apart



With expansion screw key, standard (straight, small).

Indication: for transverse expansion of narrow maxillae and as activation screw for bimaxillary functional gnatho-orthopedic appliances

Directional arrow: Plastic, neon yellow

Maximum expansion	4 x 1/4 turn	REF	Quantity
7.5 mm	= 0.8 mm	600-310-30	10 pieces

→ 14,5 mm →
→ 11 mm →



Sectional screw Mini

Material: stainless steel

- with **wide** expansion screw holder
- skeleton design
- guide pins: U-shaped, **straight**

With expansion screw key, standard (straight, small).


Indication: for movement of individual teeth

Directional arrow: printed, blue-yellow

Maximum expansion	4 x 1/4 turn	REF	Quantity
4.0 mm	= 0.35 mm	600-502-10	1 piece
4.0 mm	= 0.35 mm	600-502-30	10 pieces
4.0 mm	= 0.35 mm	600-502-60	100 pieces



← 9,5 mm →



Sectional screw Mini

Material: stainless steel

- with **wide** expansion screw holder
- skeleton design
- guide pins: U-shaped, **angulated facing down**

With expansion screw key, standard (straight, small).


Indication: for movement of individual teeth

Directional arrow: printed, blue-yellow

Maximum expansion	4 x 1/4 turn	REF	Quantity
4.0 mm	= 0.35 mm	600-503-10	1 piece
4.0 mm	= 0.35 mm	600-503-30	10 pieces
4.0 mm	= 0.35 mm	600-503-60	100 pieces



← 9,5 mm →



Special expansion screws



Sectional screw Medium

Material: stainless steel



- with **cone-shaped** expansion screw holder
 - skeleton design
 - guide pins: **straight**, with mushroom head retention
 - holes in body of screws for better retention in acrylic
 - built-in stops prevent screws from coming apart
- With expansion screw key, standard (straight, small).

Indication: for distalization and segmentation of appliances, Y plates

Directional arrow: Plastic, neon yellow

Maximum expansion	4 x 1/4 turn	REF	Quantity
6.0 mm	= 0.4 mm	600-500-10	1 piece
6.0 mm	= 0.4 mm	600-500-30	10 pieces
6.0 mm	= 0.4 mm	600-500-60	100 pieces

→ 12,5 mm ←



Sectional screw Medium

Material: stainless steel


- with **cone-shaped** expansion screw holder
 - skeleton design
 - guide pins: U-shaped, **angulated facing down**
 - holes in body of screws for better retention in acrylic
- With expansion screw key, standard (straight, small).

Indication: for distalization and segmentation of appliances, protrusion of an anterior segment

Directional arrow: Plastic, neon yellow

Maximum expansion	4 x 1/4 turn	REF	Quantity
6.0 mm	= 0.4 mm	600-501-10	1 piece
6.0 mm	= 0.4 mm	600-501-30	10 pieces
6.0 mm	= 0.4 mm	600-501-60	100 pieces

→ 13 mm ←



Retraction screw Extra Mini

Material: stainless steel

- with **cone-shaped** expansion screw holder
- screw body: closed
- guide pins: U-shaped, straight



With expansion screw key, standard (straight, small).

Indication: for space closing and movement of individual teeth

Directional arrow: printed, blue-yellow

Maximum retraction	4 x 1/4 turn	REF	Quantity
3.0 mm	= 0.4 mm	600-602-30	10 pieces

← 13,5 mm →



Retraction screw Mini

Material: stainless steel

- with **cone-shaped** expansion screw holder
- screw body: closed
- guide pins: U-shaped, straight



With expansion screw key, standard (straight, small).

Indication: for space closing and movement of individual teeth

Directional arrow: printed, blue-yellow

Maximum retraction	4 x 1/4 turn	REF	Quantity
2.5 mm	= 0.4 mm	600-601-30	10 pieces

← 13 mm →



Expansion screws

Special expansion screws



Retraction screw Medium

Material: stainless steel


- with **cone-shaped** expansion screw holder
- screw body: closed
- guide pins: U-shaped, straight

With expansion screw key, standard (straight, small).

Indication: for space closing and movement of individual teeth

Directional arrow: printed, blue-yellow

Maximum retraction	4 x 1/4 turn	REF	Quantity
4.5 mm	= 0.4 mm	600-600-30	10 pieces



Neo-Bertoni three dimensional screw – three movements (type A)

Material: stainless steel


- anatomically angled multi-dimensional screw
- three independently controlled movements
- same total expansion per segment

With expansion screw key, standard (straight, small).

Indication: for combined transverse and protrusion movements with Y-plates

Directional arrow: Plastic, neon yellow

Maximum expansion	4 x 1/4 turn	REF	Quantity
3.5 mm / segment	= 0.40 mm / segment	600-131-10	1 piece



NEW

Neo-Bertoni three dimensional screw – two movements (type B)

Material: stainless steel



- anatomically angled multi-dimensional screw
- two independently controlled movements (symmetrical transverse expansion)
- same total expansion per segment

With expansion screw key, standard (straight, small).

Indication: for combined transverse and protrusion movements with Y plates

Directional arrow: Plastic, neon yellow

Maximum expansion	4 x 1/4 turn	REF	Quantity
3.5 mm / segment	= 0.40 mm / segment	600-132-10	1 piece



NEW



Neo-Bertoni 3D

Material: stainless steel


- anatomical, individually adjustable three dimensional screw
- front segment can be adapted individually to the palate
- two independently controlled movements (symmetrical transverse expansion)
- same total expansion per segment

With expansion screw key, standard (straight, small).

Indication: for combined transverse and protrusion movements with Y plates

Directional arrow: Plastic, neon yellow


Maximum expansion	4 x 1/4 turn	REF	Quantity
3.5 mm / segment	= 0.40 mm / segment	600-133-10	1 piece



NEW



Special expansion screws



Three dimensional screw acc. to Bertoni – three movements (type A)

Material: stainless steel



- strengthened guides
- three independently controlled movements
- same total expansion per segment

With expansion screw key straight, standard (straight, small).

Indication: for combined transverse and protrusion movements with Y plates

Directional arrow: printed, blue-yellow

Maximum expansion	4 x 1/4 turn	REF	Quantity
4.0 mm / segment	= 0.35 mm / segment	602-605-10	1 piece



Three dimensional screw acc. to Bertoni – two movements (type B)

Material: stainless steel


- strengthened guides
- two independently controlled movements (symmetrical transverse expansion)
- same total expansion per segment

With expansion screw key straight, standard (straight, small).

Indication: for combined transverse and protrusion movements with Y plates

Directional arrow: printed, blue-yellow

Maximum expansion	4 x 1/4 turn	REF	Quantity
4.0 mm / segment	= 0.35 mm / segment	602-606-10	1 piece



Three dimensional screw acc. to Steiner

Material: stainless steel

- strengthened guides
- individually adjustable middle bar for adapting to different palatal vault curvatures
- two independently controlled movements (symmetrical transverse expansion)
- same total expansion per segment


With expansion screw key, standard (straight, small).

Indication: for combined transverse and protrusion movements with Y plates

Directional arrow: printed, blue-yellow



Maximum expansion	4 x 1/4 turn	REF	Quantity
4.0 mm / segment	= 0.35 mm / segment	602-604-10	1 piece



Expansion screws

Special expansion screws



Fan-type expansion screw

Material: stainless steel

- with **wide** expansion screw holder
- plate sections are opened fanwise

With expansion screw key straight, standard (straight, small).

Indication: for maxillary anterior or posterior expansion

Directional arrow: printed, blue-yellow

Maximum expansion	4 x 1/4 turn	REF	Quantity
7.0 mm	= 0.8 mm	606-600-00	1 piece
7.0 mm	= 0.8 mm	606-600-30	10 pieces

→ 12 mm →



Mandibular bow screw acc. to Müller

Material: stainless steel

- with **cone-shaped** expansion screw holder
- minimal fan-type opening of the arch screw, according to the anterior arch


With expansion screw key straight, standard (straight, small).

Indication: for transverse expansion of the mandibular anterior arch

Directional arrow: Plastic, neon yellow

Maximum expansion	4 x 1/4 turn	REF	Quantity
5.0 mm	= 0.8 mm	606-701-00	1 piece
5.0 mm	= 0.8 mm	606-701-30	10 pieces

← 20 mm →




Neo-bite jumping appliance acc. to F.M. Sander*, set

Material: stainless steel

Patented, three-dimensional Bertoni expansion screw with guide ridges for forward thrust acc. to F.M. Sander*. The latest generation of the Sander appliance allows, alongside class II/1 treatment, corrections in class II/2 thanks to the sagittal expansion element.

REF 600-120-00	1 set
Contents:	
1 x Neo-bite jumping appliance acc. to F.M. Sander	REF 600-121-10
1 x Neo-bite jumping appliance mounting template	REF 600-122-00
1 x Neo-bite jumping appliance guide template (for inclined plane)	REF 600-123-00




NEW

Neo-bite jumping appliance acc. to F.M. Sander*

Material: stainless steel

Directional arrow: Plastic, neon yellow

Maximum expansion	4 x 1/4 turn	REF	Quantity
4.0 mm / segment	= 0.35 mm / segment	600-121-10	1 piece



NEW



Neo-bite jumping appliance mounting template

acc. to Sander & Höbel

The angle to the occlusal plane can be adapted to each individual patient in combination with the slots at 55°, 60° and 65° in the expansion screw holder of the neo-bite jumping appliance.

REF 600-122-00	1 piece
----------------	---------

NEW



Neo-bite jumping appliance guide template

acc. to Sander & Höbel

To produce an inclined plane in the mandible.

REF 600-123-00	10 pieces
----------------	-----------

NEW



Expansion screws

Special expansion screws



Flexible forward thrust double plate system acc. to Schaneng, set

Material: stainless steel

- flexible intermaxillary element
- maxillary screws with intermaxillary guide ridge for forward thrust double plate. Double loop, wax coated

With expansion screw key, standard (straight, small).

Indication: for treatment of Angle Class II/1 (distoclusion), for advancement of the mandible

REF 600-000-00	1 set
----------------	-------

Contents:

3 x screw for maxilla with guide ridge for forward thrust double plate	REF 600-001-00
1 x adapter, narrow	REF 600-003-00
1 x adapter, wide	REF 600-004-00
1 x locating template	REF 600-005-00



Guide ridge for flexible forward thrust double plate system acc. to Schaneng

With double loop, wax coated.

Material: stainless steel

Length of guide ridge: 14.0 mm

Version	REF	Quantity	REF	Quantity
with palatal expansion screw	600-001-00	1 piece	600-001-30	10 pieces
without palatal expansion screw	600-002-00	1 piece	600-002-30	10 pieces



Guide ridge adapter for flexible forward thrust double plate system acc. to Schaneng

To produce the diagonal level in the mandibular.

Adapter	Width	REF	Quantity
narrow	14.0 mm	600-003-00	10 pieces
wide	18.0 mm	600-004-00	10 pieces



Locating template acc. to Schaneng

Locating template to produce the 60° angle of inclination.

Material: stainless steel

REF 600-005-00	1 piece
----------------	---------

Expansion screw for split activator

Material: stainless steel

- with **cone-shaped** expansion screw holder
- strengthened guides

With expansion screw key, standard (straight, small).



Indication: for separate expansion of maxilla and mandible

Directional arrow: printed, blue-yellow

Maximum expansion	4 x 1/4 turn	REF	Quantity
7.0 mm	= 0.8 mm	600-610-00	1 piece
7.0 mm	= 0.8 mm	600-610-30	10 pieces



→ 14,5 mm →



Expansion screw for progenia activator acc. to Weise

Material: stainless steel

- without expansion screw holder
- for fabricating horizontally split activators, which allow connection and shifting of both split planes
- stable connection of the maxillary plate with the mandibular plate using a strong U-shaped bar

With expansion screw key, standard (straight, small).



Indication: for the treatment of Angle Class III (progenia), for distal shifting of the mandible

Directional arrow: printed, blue-yellow

Maximum expansion	4 x 1/4 turn	REF	Quantity
5.0 mm	= 0.9 mm	603-000-00	1 piece



→ 14 mm →



Expansion screws

Special expansion screws




Piston spring screws, set

Material: stainless steel

- delicate design
- constant, elastic and finely controllable pressure on the tooth to be moved
- a protrusion of up to 7 mm can be attained according to the progress of treatment using three different lengths of screws (4 / 6 / 8 mm)

Indication: for movement of individual teeth

4 x 1/4 turn = 0.35 mm



REF 609-500-00

1 set

Contents:

10 x piston spring screw, 4 mm	REF 609-004-00
10 x piston spring screw, 6 mm	REF 609-006-00
10 x piston spring screw, 8 mm	REF 609-008-00
10 x positioning screw	REF 609-100-00
10 x serrated housing nut	REF 609-200-30
1 x reamer	REF 609-300-01
1 x screwdriver	REF 609-400-00



Piston spring screw

Material: stainless steel

- delicate design
- constant, elastic and finely controllable pressure on the tooth to be moved

Indication: for movement of individual teeth

Length	4 x 1/4 turn	REF	Quantity
4.0 mm	= 0.35 mm	609-004-00	10 pieces
6.0 mm	= 0.35 mm	609-006-00	10 pieces
8.0 mm	= 0.35 mm	609-008-00	10 pieces



Positioning screw

For piston spring screw.

To position the piston spring screw into the acrylic.

Material: stainless steel

Length	REF	Quantity
12.0 mm	609-100-00	10 pieces



Serrated housing nut

For piston spring screw.

To position into acrylic or to weld onto wire elements.

Material: stainless steel

Length	REF	Quantity
3.0 mm	609-200-30	10 pieces

Reamer

For piston spring screw.

Material: stainless steel

REF 609-300-01

1 piece

NEW



Screwdriver

For piston spring screw.

To activate the piston spring screw.

Material: aluminum

REF 609-400-00

1 piece



Variety expansion screws



Variety click Mini – 8

Monoguided special expansion screw with click feature to prevent unwanted reversing. Primarily for treatment of malocclusions in the mandible using fixed appliances, e.g. for distalization, mesialization or for space closure by moving individual teeth. Also for early palatal expansion.

- optimal screw locking during the active treatment phase
- small dimensions with maximum expansion
- reliable activation / deactivation of the screw with clearly perceptible click after a $\frac{1}{4}$ turn
- break resistant retention legs and stable guidance throughout the entire expansion
- permanently laser marked arrow on the body of the screw indicates the opening direction and max. expansion

Pack (1 piece) with 1 safety key, pack (10 pieces) with 10 safety keys.

Assembly dimension: 12.9 mm

Screw body: 9.6 mm

4 x $\frac{1}{4}$ turn = 0.8 mm

Material: stainless steel

Maximum expansion	REF (without retention legs)	REF (straight retention legs)	Quantity
8.0 mm	602-836-12	602-836-10	1 piece
8.0 mm	602-836-32	602-836-30	10 pieces



Variety click SP (small palate) Mini – 8

Monoguided special expansion screw with click feature to prevent unwanted reversing. For the fixed appliance technique. Can be widely used in the upper and lower jaw, for distalization, mesialization of posterior teeth, or for space closure by moving individual teeth. Also for early palatal expansion.

- optimal screw locking during the active treatment phase
- especially suitable for very narrow jaws
- retention legs can be bent very close to the screw body
- small dimensions with maximum expansion
- reliable activation / deactivation of the screw with clearly perceptible click after a $\frac{1}{4}$ turn
- break resistant retention legs and stable guidance throughout the entire expansion
- permanently laser marked arrow on the body of the screw indicates the opening direction and max. expansion

Pack (1 piece) with 1 safety key, pack (10 pieces) with 10 safety keys.

Assembly dimension: 12.9 mm


Screw body: 9.6 mm

4 x $\frac{1}{4}$ turn = 0.8 mm

Material: stainless steel

Maximum expansion	REF (without retention legs)	REF (straight retention legs)	Quantity
8.0 mm	602-838-12	602-838-10	1 piece
8.0 mm	602-838-32	602-838-30	10 pieces

Dimensions Variety click / Variety click SP Mini – 8



Variety click Maxi – 12

Monoguided special expansion screw with click feature to prevent unwanted reversing. Primarily for treatment of malocclusions in the mandible using fixed appliances, e.g. for distalization, mesialization or for space closure by moving individual teeth. Also for early palatal expansion.

- optimal screw locking during the active treatment phase
 - small dimensions with maximum expansion
 - reliable activation / deactivation of the screw with clearly perceptible click after a $\frac{1}{4}$ turn
 - break resistant retention legs and stable guidance throughout the entire expansion
 - permanently laser marked arrow on the body of the screw indicates the opening direction and max. expansion
- Pack (1 piece) with 1 safety key, pack (10 pieces) with 10 safety keys.

Assembly dimension: 17.2 mm

Screw body: 9.6 mm

4 x $\frac{1}{4}$ turn = 0.8 mm

Material: stainless steel

Maximum expansion	REF (without retention legs)	REF (straight retention legs)	Quantity
12.0 mm	602-837-12	602-837-10	1 piece
12.0 mm	602-837-32	602-837-30	10 pieces



Variety click SP (small palate) Maxi – 12

Monoguided special expansion screw with click feature to prevent unwanted reversing. Primarily for treatment of malocclusions in the mandible using fixed appliances, e.g. for distalization, mesialization or for space closure by moving individual teeth. Also for early palatal expansion.

- optimal screw locking during the active treatment phase
 - especially suitable for very narrow jaws
 - retention legs can be bent very close to the screw body
 - small dimensions with maximum expansion
 - reliable activation / deactivation of the screw with clearly perceptible click after a $\frac{1}{4}$ turn
 - break resistant retention legs and stable guidance throughout the entire expansion
 - permanently laser marked arrow on the body of the screw indicates the opening direction and max. expansion
- Pack (1 piece) with 1 safety key, pack (10 pieces) with 10 safety keys.

Assembly dimension: 17.2 mm

Screw body: 9.6 mm

4 x $\frac{1}{4}$ turn = 0.8 mm

Material: stainless steel

Maximum expansion	REF (without retention legs)	REF (straight retention legs)	Quantity
12.0 mm	602-839-12	602-839-10	1 piece
12.0 mm	602-839-32	602-839-30	10 pieces



Variety click

Variety click SP

Dimensions Variety click / Variety click SP Maxi – 12

Variety expansion screws



Variety Mini – 8

Monoguided special expansion screw. Primarily for treatment of malocclusions in the mandible using fixed appliances, e.g. for distalization, mesialization or for space closure by moving individual teeth.

Also for early palatal expansion.

- small dimensions with maximum expansion
- break resistant retention legs and stable guidance throughout the entire expansion
- the integrated spindle friction prevents unwanted reversing
- permanently laser marked arrow on the body of the screw indicates the opening direction and max. expansion
- screw can be tied down once desired expansion has been reached

Pack (1 piece) with 1 safety key, pack (10 pieces) with 10 safety keys.

Assembly dimension: 12.9 mm

Screw body: 9.6 mm

4 x 1/4 turn = 0.8 mm

Material: stainless steel

Maximum expansion	REF (without retention legs)	REF (straight retention legs)	Quantity
8.0 mm	602-815-12	602-815-10	1 piece
8.0 mm	602-815-32	602-815-30	10 pieces



Variety SP (small palate) Mini – 8

Monoguided special expansion screw. Primarily for treatment of malocclusions in the mandible using fixed appliances, e.g. for distalization, mesialization or for space closure by moving individual teeth.

Also for early palatal expansion.

- especially suitable for very narrow jaws
- retention legs can be bent very close to the screw body
- small dimensions with maximum expansion
- break resistant retention legs and stable guidance throughout the entire expansion
- the integrated spindle friction prevents unwanted reversing
- permanently laser marked arrow on the body of the screw indicates the opening direction and max. expansion
- screw can be tied down once desired expansion has been reached

Pack (1 piece) with 1 safety key, pack (10 pieces) with 10 safety keys.

Assembly dimension: 12.9 mm


Screw body: 9.6 mm

4 x 1/4 turn = 0.8 mm

Material: stainless steel

Maximum expansion	REF (without retention legs)	REF (straight retention legs)	Quantity
8.0 mm	602-819-12	602-819-10	1 piece
8.0 mm	602-819-32	602-819-30	10 pieces

Dimensions Variety / Variety SP Mini – 8



Variety Maxi – 12

Monoguided special expansion screw. Primarily for treatment of malocclusions in the mandible using fixed appliances, e.g. for distalization, mesialization or for space closure by moving individual teeth. Also for early palatal expansion.

- small dimensions with maximum expansion
- break resistant retention legs and stable guidance throughout the entire expansion
- the integrated spindle friction prevents unwanted reversing
- permanently laser marked arrow on the body of the screw indicates the opening direction and max. expansion
- screw can be tied down once desired expansion has been reached

Pack (1 piece) with 1 safety key, pack (10 pieces) with 10 safety keys.

Assembly dimension: 14.9 mm

Screw body: 9.6 mm

4 x 1/4 turn = 0.8 mm

Material: stainless steel



Maximum expansion	REF (without retention legs)	REF (straight retention legs)	Quantity
12.0 mm	602-806-12	602-806-10	1 piece
12.0 mm	602-806-32	602-806-30	10 pieces

Variety SP (small palate) Maxi – 12

Monoguided special expansion screw. Primarily for treatment of malocclusions in the mandible using fixed appliances, e.g. for distalization, mesialization or for space closure by moving individual teeth. Also for early palatal expansion.

- especially suitable for very narrow jaws
- retention legs can be bent very close to the screw body
- small dimensions with maximum expansion
- break resistant retention legs and stable guidance throughout the entire expansion
- the integrated spindle friction prevents unwanted reversing
- permanently laser marked arrow on the body of the screw indicates the opening direction and max. expansion
- screw can be tied down once desired expansion has been reached

Pack (1 piece) with 1 safety key, pack (10 pieces) with 10 safety keys.

Assembly dimension: 14.9 mm

Screw body: 9.6 mm

4 x 1/4 turn = 0.8 mm

Material: stainless steel



Maximum expansion	REF (without retention legs)	REF (straight retention legs)	Quantity
12.0 mm	602-809-12	602-809-10	1 piece
12.0 mm	602-809-32	602-809-30	10 pieces



Variety



Variety SP



Dimensions Variety / Variety SP Maxi – 12



Expansion screws

RPE screws



Titan hyrax® Maxi – 12

RPE screw of titanium

For rapid palatal expansion by means of fixed appliances.

The screw is anchored to the teeth using cast titanium and laser-welded cap splints or alternatively using acrylic cap splints fabricated from Orthocryl® acrylic.




Pack (1 piece) with 1 safety key, pack (10 pieces) with 10 safety keys.

Screw body: 14.4 mm

4 x 1/4 turn = 0.8 mm

Material: titanium

Maximum expansion	REF (without retention legs)	REF (straight retention legs)	Quantity
12.0 mm	602-810-12	602-810-10	1 piece
12.0 mm	602-810-32	602-810-30	10 pieces



rematitan® straight wire

Titanium grade 1

For laser welding.

Norm: DIN EN ISO 28319, DIN EN ISO 5832-2

ø	Length	REF	Quantity
1.00 mm / 39	each 100 mm	528-041-00	10 pieces
1.20 mm / 47	each 100 mm	528-042-00	10 pieces



hyrax® click Mini – 7

RPE screw with click feature to prevent unwanted reversing

For effective and reliable rapid palatal expansion by means of fixed appliances.

- optimal screw locking during the active treatment phase
- reliable activation / deactivation of the screw with clearly perceptible click after a $\frac{1}{4}$ turn
- no accidental reversing when pulling out safety key during activation
- built-in stops prevent the screws from coming apart once maximum expansion has been obtained
- small dimensions with maximum expansion
- break resistant retention legs, stable guidance throughout complete expansion
- offset body for earliest possible angling of the retention legs
- permanently laser marked arrow on the body of the screw indicates the opening direction and max. expansion



Pack (1 piece) with 1 safety key, pack (10 pieces) with 10 safety keys.

Screw body: 9.8 mm

4 x $\frac{1}{4}$ turn = 0.8 mm

Material: stainless steel

Maximum expansion	REF (without retention legs)	REF (straight retention legs)	REF (angulated retention legs)	Quantity
7.0 mm	602-830-12	602-830-10	602-831-10	1 piece
7.0 mm	602-830-32	602-830-30	602-831-30	10 pieces



hyrax® click Medium – 10

RPE screw with click feature to prevent unwanted reversing

For effective and reliable rapid palatal expansion by means of fixed appliances.

- optimal screw locking during the active treatment phase
- reliable activation / deactivation of the screw with clearly perceptible click after a $\frac{1}{4}$ turn
- no accidental reversing when pulling out safety key during activation
- built-in stops prevent the screws from coming apart once maximum expansion has been obtained
- small dimensions with maximum expansion
- break resistant retention legs, stable guidance throughout complete expansion
- offset body for earliest possible angling of the retention legs
- permanently laser marked arrow on the body of the screw indicates the opening direction and max. expansion



Pack (1 piece) with 1 safety key, pack (10 pieces) with 10 safety keys.

Screw body: 12.8 mm

4 x $\frac{1}{4}$ turn = 0.8 mm

Material: stainless steel

Maximum expansion	REF (without retention legs)	REF (straight retention legs)	REF (angulated retention legs)	Quantity
10.0 mm	602-832-12	602-832-10	602-833-10	1 piece
10.0 mm	602-832-32	602-832-30	602-833-30	10 pieces



Expansion screws

RPE screws



hyrax® click Maxi – 12

RPE screw with click feature to prevent unwanted reversing

For effective and reliable rapid palatal expansion by means of fixed appliances.

- optimal screw locking during the active treatment phase
- reliable activation / deactivation of the screw with clearly perceptible click after a $\frac{1}{4}$ turn
- no accidental reversing when pulling out safety key during activation
- built-in stops prevent the screws from coming apart once maximum expansion has been obtained
- small dimensions with maximum expansion
- break resistant retention legs, stable guidance throughout complete expansion
- offset body for earliest possible angling of the retention legs
- permanently laser marked arrow on the body of the screw indicates the opening direction and max. expansion


Pack (1 piece) with 1 safety key, pack (10 pieces) with 10 safety keys.

Screw body: 14.8 mm

4 x $\frac{1}{4}$ turn = 0.8 mm

Material: stainless steel

Maximum expansion	REF (without retention legs)	REF (straight retention legs)	REF (angulated retention legs)	Quantity
12.0 mm	602-834-12	602-834-10	602-835-10	1 piece
12.0 mm	602-834-32	602-834-30	602-835-30	10 pieces



tomas®-RPE eyelet

(coupling element for RPE screws)

Eyelet for direct coupling between tomas®-pin with all hyrax® RPE screws with two supports for easier placement on working model.

Can be used for hybrid RPE appliance or skeletal RPE appliance using Pin-Last® approach.

Outer ø eyelet	Inner ø eyelet	Thickness eyelet	REF	Quantity
5.0 mm	top 3.0 mm bottom 2.3 mm	1.5 mm	302-602-00	12 pieces + 2 holders

hyrax® Mini – 7

RPE screw

For rapid palatal expansion by means of fixed appliances.


- built-in stops prevent the screws from coming apart once maximum expansion has been obtained
 - small dimensions with maximum expansion
 - break resistant retention legs, stable guidance throughout complete expansion
 - offset body for earliest possible angling of the retention legs
 - the integrated spindle friction prevents unwanted reversing
 - permanently laser marked arrow on the body of the screw indicates the opening direction and max. expansion
 - screw can be tied down once desired expansion has been reached
- Pack (1 piece) with 1 safety key, pack (10 pieces) with 10 safety keys.

Screw body: 9.4 mm

4 x 1/4 turn = 0.8 mm

Material: stainless steel

Maximum expansion	REF (without retention legs)	REF (straight retention legs)	REF (angulated retention legs)	Quantity
7.0 mm	602-800-12	602-800-10	602-801-10	1 piece
7.0 mm	602-800-32	602-800-30	602-801-30	10 pieces



hyrax® Medium – 10

RPE screw

For rapid palatal expansion by means of fixed appliances.



- built-in stops prevent the screws from coming apart once maximum expansion has been obtained
 - small dimensions with maximum expansion
 - break resistant retention legs, stable guidance throughout complete expansion
 - offset body for earliest possible angling of the retention legs
 - the integrated spindle friction prevents unwanted reversing
 - permanently laser marked arrow on the body of the screw indicates the opening direction and max. expansion
 - screw can be tied down once desired expansion has been reached
- Pack (1 piece) with 1 safety key, pack (10 pieces) with 10 safety keys.

Screw body: 12.4 mm

4 x 1/4 turn = 0.8 mm

Material: stainless steel

Maximum expansion	REF (without retention legs)	REF (straight retention legs)	REF (angulated retention legs)	Quantity
10.0 mm	602-808-12	602-808-10	602-840-10	1 piece
10.0 mm	602-808-32	602-808-30	602-840-30	10 pieces



RPE screws



hyrax® Maxi – 12

RPE screw

For rapid palatal expansion by means of fixed appliances.

- built-in stops prevent the screws from coming apart once maximum expansion has been obtained
- small dimensions with maximum expansion
- break resistant retention legs, stable guidance throughout complete expansion
- offset body for earliest possible angling of the retention legs
- the integrated spindle friction prevents unwanted reversing
- permanently laser marked arrow on the body of the screw indicates the opening direction and max. expansion
- screw can be tied down once desired expansion has been reached


Pack (1 piece) with 1 safety key, pack (10 pieces) with 10 safety keys.

Screw body: 14.4 mm

4 x 1/4 turn = 0.8 mm

Material: stainless steel

Maximum expansion	REF (without retention legs)	REF (straight retention legs)	REF (angulated retention legs)	Quantity
12.0 mm	602-811-12	602-811-10	602-812-10	1 piece
12.0 mm	602-811-32	602-811-30	602-812-30	10 pieces



hyrax® neo Mini – 7

RPE screw

For rapid palatal expansion by means of a fixed appliance.

- (no active) expansion stop
- small dimensions with maximum expansion
- TAD-optimized design
- break resistant retention legs, stable guidance throughout complete expansion
- the integrated spindle friction prevents unwanted reversing
- permanently laser marked arrow on the body of the screw indicates the opening direction and maximum expansion
- screw can be tied down once desired expansion has been reached



Pack (1 piece) with 1 safety key, pack (10 pieces) with 10 safety keys.

Screw body: 11.2 mm

4 x 1/4 turn = 0.8 mm

Material: stainless steel

Maximum expansion	REF (without retention legs)	REF (straight retention legs)	Quantity
7.0 mm	602-845-12	602-845-10	1 piece
7.0 mm	602-845-32	602-845-30	10 pieces



hyrax® neo Medium – 10

RPE screw

For rapid palatal expansion by means of a fixed appliance.



- (no active) expansion stop
 - small dimensions with maximum expansion
 - TAD-optimized design
 - break resistant retention legs, stable guidance throughout complete expansion
 - the integrated spindle friction prevents unwanted reversing
 - permanently laser marked arrow on the body of the screw indicates the opening direction and maximum expansion
 - screw can be tied down once desired expansion has been reached
- Pack (1 piece) with 1 safety key, pack (10 pieces) with 10 safety keys.

Screw body: 12.4 mm

4 x 1/4 turn = 0.8 mm

Material: stainless steel

Maximum expansion	REF (without retention legs)	REF (straight retention legs)	Quantity
10.0 mm	602-846-12	602-846-10	1 piece
10.0 mm	602-846-32	602-846-30	10 pieces



hyrax® neo Maxi – 12

RPE screw

For rapid palatal expansion by means of a fixed appliance.


- (no active) expansion stop
 - small dimensions with maximum expansion
 - TAD-optimized design
 - break resistant retention legs, stable guidance throughout complete expansion
 - the integrated spindle friction prevents unwanted reversing
 - permanently laser marked arrow on the body of the screw indicates the opening direction and maximum expansion
 - screw can be tied down once desired expansion has been reached
- Pack (1 piece) with 1 safety key, pack (10 pieces) with 10 safety keys.

Screw body: 14.4 mm

4 x 1/4 turn = 0.8 mm

Material: stainless steel


Maximum expansion	REF (without retention legs)	REF (straight retention legs)	Quantity
12.0 mm	602-847-12	602-847-10	1 piece
12.0 mm	602-847-32	602-847-30	10 pieces



Accessories

Arrow to indicate the opening direction

For expansion screws.



Color	REF	Quantity
yellow	611-130-00	50 pieces
white	611-131-00	50 pieces
pearl-blue	611-132-00	50 pieces



Universal bending key

Easy bending of the retention legs of all hyrax® and Variety screws without using any additional pliers (leverage principle).
Material: stainless steel

REF 611-122-00	1 piece
----------------	---------



Safety key

Can be used with all Dentaaurum hyrax® and Variety screws.
For safe intra-oral activation of the expansion screws in fixed appliances.

REF 611-120-01	10 pieces
----------------	-----------



Expansion screw key, standard straight, small

Can be used with all Dentaaurum expansion screws.
To activate all expansion screws in removable appliances (outside of the mouth).
Material: stainless steel

REF 611-115-00	1 piece
REF 611-115-30	10 pieces
REF 611-115-60	100 pieces



Expansion screw key, large with elongated handle

Can be used with all Dentaaurum expansion screws.
To activate all expansion screws in removable appliances (outside of the mouth).
Material: stainless steel

REF 611-116-00	1 piece
----------------	---------

Expansion screw sample kit

The kit includes a selection from Dentaureum's large range of expansion screws and is for promotional and demonstration purposes only. The content of this kit is not certified for use on patients!

Filled with 27 expansion screws that are on stock, 2 expansion screw keys and 1 universal bending key. **1 set = 30 pieces.**

REF 612-000-00

1 set



Wedge gauge

For checking the gap width on palatal expanders.

REF 043-720-00

1 piece




Orthocryl® LC	Page	348
Orthocryl® / Assortments	Page	350
Orthocryl® / Single packs	Page	351
Disco glitter	Page	353
Decals	Page	354
Accessories	Page	358




Orthodontic acrylics.

Vibrant colors – fascinating varieties.



DENTAURUM
**ONLINE
SHOP**
shop.dentaurum.com






The versatile acrylic system.

Orthocryl® has been a fixture in orthodontic laboratories for 60 years. For generations, it has combined excellent product characteristics and easy handling.

There are no limits to your creativity, regardless of whether you choose black, white or multi-color – Orthocryl® offers maximum freedom in design with its variety of lively colors, glitter and decals. Apart from the conventional cold-curing acrylic, Orthocryl® is also available as light-curing material. Orthocryl® LC is especially designed to meet the requirements for fabricating orthodontic appliances. It is light-curing, time-saving and easy to use.



Orthocryl® LC

The light-curing acrylic.

Orthocryl®

Suited for the salt-and-pepper technique.

Orthocryl® Neon

Vivid colors.

Orthocryl® EQ

Suited for the doughing technique.

Orthocryl® black & white

High-contrast, versatile.

A wide range of field-proven accessories


All from one source.



Orthocryl® LC



The light-curing acrylic for orthodontists.



Anything but ordinary.

Orthocryl® LC is an acrylic with many advantages – it is light-curing, time-saving and easy to use. Orthocryl® LC is especially designed to meet the requirements for fabricating orthodontic appliances. But it can also be used to make drill templates for implantology, bite plates and anti-snoring devices within seconds and with little effort.

Orthocryl® LC is anything but ordinary and that makes it delightfully different. A true premium quality product with undeniable advantages.

Advantages:

- Extremely easy to use.
- Saves time and material.
- Optimal flowability.
- Wide range of applications.
- Completely odorless.
- Biocompatible.


Easy to use and flexible.

Orthocryl® LC is extremely easy to use. The product is delivered ready to use in black, opaque cartridges. Time-consuming preparation steps have become superfluous: the acrylic does not have to be doughed or salt-and-peppered anymore.

The material can be directly applied from the cartridge to the prepared model. The application injector helps to apply the right dosage in the right spot. The easy handling saves both a great amount of time when processing and finishing the appliance and a lot of material.

Due to its optimal flowability, Orthocryl® LC can be easily adapted to the prepared model.

Orthocryl® LC is available in three colors. Adding Disco glitter or decals from the extensive Orthocryl® program makes each appliance unique.



Orthocryl® LC



Orthocryl® LC starter set

Light-curing acrylic.

REF 160-400-00	1 set
----------------	-------

Contents:

2 x Orthocryl® LC, clear	REF 160-401-00
2 x Orthocryl® LC, transparent pink	REF 160-402-00
2 x Orthocryl® LC, blue	REF 160-404-00
1 x injector for Orthocryl® LC	REF 160-420-00
50 x cartridge tip for Orthocryl® LC	REF 160-421-00



Orthocryl® LC starter set, clear

Light-curing acrylic.

REF 160-400-20	1 set
----------------	-------

Contents:

6 x Orthocryl® LC, clear	REF 160-401-00
1 x injector for Orthocryl® LC	REF 160-420-00
50 x cartridge tip for Orthocryl® LC	REF 160-421-00

Orthocryl® LC

Light-curing acrylic in opaque cartridge.

Color	REF	Quantity
clear	160-401-00	30 g
transparent pink	160-402-00	30 g
blue	160-404-00	30 g



Bite paste

For Orthocryl® LC.
Color: transparent

REF 160-410-00

50 g



Injector

For Orthocryl® LC.

REF 160-420-00

1 piece



Cartridge tip

For Orthocryl® LC.

REF 160-421-00

50 pieces



Cleaning liquid

For Orthocryl® LC.

REF 160-430-00

100 ml



Stand for cartridges

For Orthocryl® LC. Fits 6 cartridges.

Material: Plastic

Do not use disinfectants containing ethanol,
phenol or alcohol.

Dimensions: approx. 170 mm x 110 mm x 110 mm (W x H x D)

Color: transparent, clear

REF 160-422-00

1 piece



Orthodentic acrylics

Orthocryl® / Assortments



Orthocryl® starter set, clear

REF 161-912-00	1 set
----------------	-------

Contents:

1 x Orthocryl® powder (polymer), clear	REF 160-112-00
1 x Orthocryl® liquid (monomer), clear	REF 161-100-00
1 x separating medium	REF 162-800-00
2 x spray bottle	REF 162-100-00
1 x measuring beaker	REF 162-200-00
1 x funnel	REF 162-600-00
1 x spray nozzle for powder (polymer)	REF 162-751-00
1 x spray nozzle for liquid (monomer)	REF 162-752-00



Orthocryl® set, small

Small set 200 g powder (polymer) and 100 ml liquid (monomer)

Color	REF	Quantity
transparent pink	160-501-00	1 set
clear	160-502-00	1 set



Orthocryl® multi-colored palette, assortment

REF 160-600-00	1 assortment
----------------	--------------

Contents:

1 x Orthocryl® powder (polymer), clear	REF 160-112-00
1 x Orthocryl® liquid (monomer), red	REF 161-127-00
1 x Orthocryl® liquid (monomer), green	REF 161-128-00
1 x Orthocryl® liquid (monomer), blue	REF 161-129-00
1 x Orthocryl® liquid (monomer), yellow	REF 161-130-00
1 x funnel	REF 162-600-00



Orthocryl® Neon assortment

Brightly colored acrylic

Liquid (monomer)

Brightly colored Orthocryl® liquid (monomer) for producing brightly colored orthodontic appliances. Cold curing acrylic resin with wide range of processing options permitting the efficient fabrication of orthodontic appliances. Easy processing, using both the doughing and the salt-and-pepper technique in conjunction with the clear Orthocryl® powder (polymer, REF 160-112-00 / 160-212-00). Disco glitter can also be used to enhance the effect of the colors.

REF 160-700-00	1 assortment
----------------	--------------

Contents:

1 x Orthocryl® powder (polymer), clear	REF 160-112-00
1 x Orthocryl® liquid (monomer), neon yellow	REF 161-135-00
1 x Orthocryl® liquid (monomer), neon orange	REF 161-136-00
1 x Orthocryl® liquid (monomer), neon pink	REF 161-137-00
1 x Orthocryl® liquid (monomer), neon green	REF 161-138-00
1 x Orthocryl® liquid (monomer), neon blue	REF 161-139-00
1 x funnel	REF 162-600-00

Orthocryl®

Powder (polymer)

Cold curing acrylic resin with wide range of processing options permitting the efficient fabrication of orthodontic appliances. **Optimal stability using the salt-and-pepper technique.**

Also very suitable for the doughing technique.

Color	REF	Quantity	REF	Quantity
clear	160-112-00	1 kg	160-212-00	10 kg



Orthocryl® EQ

Powder (polymer)

Cold curing acrylic resin with wide range of processing options permitting the efficient fabrication of orthodontic appliances. **Optimal flowability using the doughing technique.**

Also very suitable for the salt-and-pepper technique.

Color	REF	Quantity	REF	Quantity
clear	160-300-00	1 kg	160-301-00	10 kg



Orthocryl® black & white

Powder (polymer), opaque

Cold curing acrylic resin with wide range of processing options permitting the efficient fabrication of orthodontic appliances. **Optimal stability using the salt-and-pepper technique.**

Also very suitable for the doughing technique.

Color	REF	Quantity	REF	Quantity
black	160-011-00	200 g	160-013-00	1 kg
white	160-012-00	200 g	160-014-00	1 kg



Orthocryl® Neon

Brightly colored acrylic

Liquid (monomer)

Brightly colored Orthocryl® liquid (monomer) for producing brightly colored orthodontic appliances. Cold curing acrylic resin with wide range of processing options permitting the efficient fabrication of orthodontic appliances. Easy processing, using both the doughing and the salt-and-pepper technique in conjunction with the clear Orthocryl® powder (polymer, REF 160-112-00 / 160-212-00). Disco glitter can also be used to enhance the effect of the colors.

Color	REF	Quantity
neon yellow	161-135-00	250 ml
neon orange	161-136-00	250 ml
neon pink	161-137-00	250 ml
neon green	161-138-00	250 ml
neon blue	161-139-00	250 ml



Orthodontic acrylics

Orthocryl® / Single packs



Orthocryl®

Liquid (monomer)

Cold curing acrylic resin with wide range of processing options permitting the efficient fabrication of orthodontic appliances.

Color	REF	Quantity	REF	Quantity	REF	Quantity
red	161-127-00	250 ml				
green	161-128-00	250 ml				
blue	161-129-00	250 ml				
yellow	161-130-00	250 ml				
turquoise	161-131-00	250 ml				
pink	161-132-00	250 ml				
emerald green	161-133-00	250 ml				
violet	161-134-00	250 ml				
clear			161-100-00	500 ml	161-150-00	1000 ml
transparent pink			161-300-00	500 ml	161-350-00	1000 ml



Orthocryl®

Color concentrate

For tinting clear Orthocryl® liquid (monomer, REF 161-100-00 / 161-150-00). With one dispensing pipette.

Color	Mixing ratio	REF	Quantity
red	1 : 10	161-620-00	100 ml
blue	1 : 50	161-622-00	100 ml

Mixing ratio:

	Liquid, clear	Color concentrate	Liquid, clear	Color concentrate
1 : 10 (red)	100 ml	10 ml	50 ml	5 ml
1 : 50 (blue)	100 ml	2 ml	50 ml	1 ml

Orthocryl® Disco glitter

For coloring the removable appliances. The glitter is simply added as desired to the Orthocryl® powder (polymer) to give the appliance an attractive luster. No change in the processing of Orthocryl® is necessary. The appliances can be finished and polished without difficulty.

Mixing ratio:

0.5 : 100 (by volume glitter : powder)

Color	REF	Quantity
iris	161-605-00	50 g
red	161-606-00	50 g
blue	161-607-00	50 g
silver	161-608-00	50 g
gold	161-609-00	50 g
mother-of-pearl	161-611-00	50 g



Decals



Decals, assortment

For placing in orthodontic appliances. For polymerizing into Orthocryl® acrylic, allowing removable orthodontic appliances to be fabricated to suit individual patient preferences.

Contents:

Assortment consists of 20 individual designs
20 pieces each = 400 pieces, plus tweezers.

REF 160-100-00	1 assortment
Contents:	
20 x Crocodile	REF 160-100-01
20 x Frog	REF 160-100-02
20 x Parrot	REF 160-100-04
20 x Ladybug	REF 160-100-05
20 x Dragon	REF 160-100-06
20 x Fish	REF 160-100-07
20 x Formula 1 (Indy racecar)	REF 160-100-09
20 x Elephant	REF 160-100-10
20 x Heart	REF 160-100-11
20 x Duck	REF 160-100-12
20 x Dinosaur	REF 160-100-13
20 x Motorcycle	REF 160-100-14
20 x Snail	REF 160-100-15
20 x Teddy bear	REF 160-100-16
20 x Mushrooms	REF 160-100-18
20 x Car	REF 160-100-19
20 x Raven	REF 160-100-27
20 x Devil	REF 160-100-29
20 x Soccer player	REF 160-100-31
20 x Mouse	REF 160-100-45



Decal Football

Sheet with two different sizes of design

For placing in orthodontic appliances. For polymerizing into Orthocryl® acrylic, allowing removable orthodontic appliances to be fabricated to suit individual patient preferences.

Motif	REF	Quantity
Football	160-100-39	1 sheet = 20 pieces



Decal Beaver

For placing in orthodontic appliances. For polymerizing into Orthocryl® acrylic, allowing removable orthodontic appliances to be fabricated to suit individual patient preferences.

Motif	REF	Quantity
Beaver	160-100-40	1 sheet = 20 pieces



Decals Black Edition

For placing in orthodontic appliances. For polymerizing into Orthocryl® acrylic, allowing removable orthodontic appliances to be fabricated to suit individual patient preferences.

Motif	REF	Quantity
Black Edition	160-100-41	1 sheet = 12 motifs

Decal

For placing in orthodontic appliances. For polymerizing into Orthocryl® acrylic, allowing removable orthodontic appliances to be fabricated to suit individual patient preferences.

Motif	REF	Quantity
Hello Beaver	160-100-42	2 sheets = 20 pieces
Horse	160-100-43	2 sheets = 20 pieces
Pink Pony	160-100-44	2 sheets = 20 pieces
Unicorn	160-100-46	2 sheets = 20 pieces

**Decal**

For placing in orthodontic appliances. For polymerizing into Orthocryl® acrylic, allowing removable orthodontic appliances to be fabricated to suit individual patient preferences.

Motif	REF	Quantity
Crocodile	160-100-01	20 pieces
Frog	160-100-02	20 pieces
Parrot	160-100-04	20 pieces



Decals

Decal

For placing in orthodontic appliances. For polymerizing into Orthocryl® acrylic, allowing removable orthodontic appliances to be fabricated to suit individual patient preferences.



Motif	REF	Quantity
Ladybug	160-100-05	20 pieces
Dragon	160-100-06	20 pieces
Fish	160-100-07	20 pieces
Formula 1 (Indy Racecar)	160-100-09	20 pieces
Elephant	160-100-10	20 pieces
Heart	160-100-11	20 pieces
Duck	160-100-12	20 pieces
Dinosaur	160-100-13	20 pieces
Motorcycle	160-100-14	20 pieces
Snail	160-100-15	20 pieces

Decal

For placing in orthodontic appliances. For polymerizing into Orthocryl® acrylic, allowing removable orthodontic appliances to be fabricated to suit individual patient preferences.

Motif	REF	Quantity
Teddybear	160-100-16	20 pieces
Mushrooms	160-100-18	20 pieces
Car	160-100-19	20 pieces
Raven	160-100-27	20 pieces
Spider	160-100-28	20 pieces
Devil	160-100-29	20 pieces
Graffiti	160-100-30	20 pieces
Soccer Player	160-100-31	20 pieces
Mouse	160-100-45	20 pieces
Dolphin	160-100-34	20 pieces



Accessories



Spray bottle

For powder (polymer) and liquid (monomer).

Capacity	REF	Quantity
100 ml	162-100-00	1 piece



Spray nozzle for Orthocryl® liquid (monomer)

Inner ø	Length	REF	Quantity
0.40 mm	10 mm	162-752-00	12 pieces



Spray nozzle for Orthocryl® powder (polymer)

REF 162-751-00	1 piece
----------------	---------



Mixing bowl

Flexible, with suction base. For all mixing jobs in the practice and laboratory.

Do not use with silicone materials!

Material: silicone

Capacity	REF	Quantity
5 ml	162-501-00	3 pieces
25 ml	162-502-00	1 piece
75 ml	162-503-00	1 piece



Funnel

REF 162-600-00	1 piece
----------------	---------



Measuring beaker

REF 162-200-00	1 piece
----------------	---------



Dispensing pipette

For Orthocryl® color concentrates.

REF 162-101-00	10 pieces
----------------	-----------

Separating medium

Special alginate-based separating medium that forms a thin, abrasive-resistant separating film.

Separates:

- plaster from acrylic
- plaster from plaster

Color: pink

REF 162-800-00

1000 ml



3D separating medium

Forms a separating layer between the acrylic models and Orthocryl® / Orthocryl® LC.

Color: green

REF 162-801-00

100 ml



Modeling wax, pink

For acrylic work, construction bites etc.

- pliable
- very easy to model
- firm up to approx. 37 °C / 98.6 °F
- easy fixation of teeth



Melting point	Dimensions (L x W)	Thickness	REF	Quantity
60 °C / 140 °F	175 mm x 80 mm	1.50 mm	119-150-00	500 g

Thermo wax

- heat-resistant wax for blocking-out
- for covering springs used in the Orthocryl® technique
- prevents acrylic from becoming milky in appearance

Melting point	Color	REF	Quantity
87 °C / 188.6 °F	pink	120-170-00	50 g



Sticky wax

- good adhesive effect on all materials
- ideal for the fixation of retentive or support elements, repair and soldering work

Melting point	Color	REF	Quantity
71 °C / 159.8 °F	red	122-300-00	60 g



Accessories

Bur for acrylics

with mill toothing cut, straight

- cutting properties

Application: acrylics

Maximum speed: 15000 min⁻¹



ISO number coding	Size	Form	REF	Quantity
310.104.266.171.040	small	bud, rounded, long	180-201-00	1 piece
310.104.266.171.050	medium	bud, rounded, long	180-202-00	1 piece
310.104.260.171.070	medium	bud, rounded	180-304-00	1 piece
310.104.260.171.080	large	bud, rounded	180-305-00	1 piece



Super hard metal bur

with cross toothing cut, coarse

- high precision milled
- best cutting properties

Application: plaster, acrylics

Recommended speed: 15000 min⁻¹ - 20000 min⁻¹

Maximum speed: 30000 min⁻¹

ISO number coding	Form	REF	Quantity
500.104.263.190.070	bud, rounded, slender	180-308-00	1 piece

**Super hard metal bur and cutter
with cross toothing cut, coarse**

- high precision milled
- best cutting properties

Application: acrylics, cast partial denture

Recommended speed: 20000 min⁻¹ - 25000 min⁻¹

Maximum speed: 50000 min⁻¹

ISO number coding	Form	REF	Quantity
500.104.116.190.023	cylindrical, only side cutting	123-580-00	1 piece
500.104.137.190.023	cylindrical, domed end	123-581-00	1 piece
500.104.292.190.023	torpedo, cylindrical	123-582-00	1 piece
500.104.001.190.023	spherical (round), ø 2.1 mm	123-583-00	1 piece
500.104.274.190.060	bullet	123-584-00	1 piece
500.104.194.190.023	conical, domed end	123-585-00	1 piece



**Super hard metal bur
with cross toothing cut, coarse**

- high precision milled
- best cutting properties

Application: plaster, acrylics

Recommended speed: 20000 min⁻¹ - 25000 min⁻¹

Maximum speed: 30000 min⁻¹

ISO number coding	Form	REF	Quantity
500.104.237.220.060	pear	123-594-00	1 piece
500.104.194.220.060	conical, domed end	123-596-00	1 piece



Accessories

Super hard metal bur

with mill toothing cut, wound right

- tungsten carbide
- precision cut
- best cutting properties


Application: CoCr alloys, acrylics

Recommended speed: 20000 min⁻¹ - 25000 min⁻¹

Maximum speed: 30000 min⁻¹



ISO number coding	Form	REF	Quantity
500.104.274.175.060	bullet	123-597-00	1 piece
500.104.194.175.045	conical, domed end	123-598-00	1 piece



Saw

For orthodontics / prosthetics.

Small saw with replaceable blades.

Contents:

Saw incl. 3 saw blades

REF 180-702-00	1 piece
-----------------------	---------



Saw blade

For saw.

Length	REF	Quantity
70 mm	180-703-00	5 pieces



Polishing brush

For polishing machines.

Application: orthodontics, crowns and bridges, plastic prosthetics, cast partial denture

Maximum speed: 2800 min⁻¹

ø	REF	Quantity	REF	Quantity
50 mm	141-700-00	1 piece	141-700-10	10 pieces



Polishing brush

Combination Chungking bristles / cotton, for polishing machines.

Application: orthodontics, crowns and bridges, plastic prosthetics

Maximum speed: 1500 min⁻¹

ø	REF	Quantity
80 mm	141-710-00	1 piece

NEW

Polishing brush

Short Chungking bristles, for polishing machines.

Application: orthodontics, crowns and bridges, cast partial denture

Maximum speed: 6000 min⁻¹

Number of brush ranges	ø	REF	Quantity
2	60 mm	141-711-00	10 pieces
3	60 mm	141-712-00	10 pieces

**Palatal cup brush**

Chungking cup bristles, for polishing machines.

Application: orthodontics, plastic prosthetics

Maximum speed: 3000 min⁻¹

ø	REF	Quantity	REF	Quantity
30 mm	141-713-10	1 piece	141-713-00	10 pieces

**Cloth polishing brush, mounted**

For handpiece.

Application: orthodontics, crowns and bridges, plastic prosthetics, cast partial denture

Recommended speed: 1500 min⁻¹ - 4000 min⁻¹

Maximum speed: 5000 min⁻¹

ø	REF	Quantity
20 mm	181-100-01	6 pieces

**Cloth polishing brush**

For polishing machines.

Application: orthodontics, crowns and bridges, plastic prosthetics, cast partial denture

Maximum speed: 5000 min⁻¹

ø	REF	Quantity
100 mm	181-200-00	1 piece

**Chamois leather polishing brush**

For high shine polish, for polishing machines.

Application: orthodontics, crowns and bridges, plastic prosthetics, cast partial denture

Maximum speed: 1500 min⁻¹

ø	REF	Quantity
90 mm	181-201-00	1 piece



Accessories



Edelweiss polishing stick

Ideally suitable for high-gloss finish on acrylic.

REF 190-100-00	80 g
----------------	------



Edelweiss acrylic polishing liquid

Ideally suitable for high-gloss finish on acrylic.

REF 190-200-00	1000 g
----------------	--------

Silicon polisher

- non-damaging to metal parts, clasps

Application: orthodontics appliances, pre-polishing of prosthetic acrylics

Recommended speed: 5000 min⁻¹ - 7000 min⁻¹

Maximum speed: 15000 min⁻¹



ISO number coding	Form	REF	Quantity
658.104.275.523.100	torpedo	138-623-00	10 pieces



Mandrel

Version	Maximum speed	REF	Quantity
for emery cloth	16000 min ⁻¹	139-500-00	10 pieces



Emery cloth

Grain size	Width	REF	Quantity
120 = fine	15 mm	140-015-00	25 m

FKO-split-post fixator

acc. to Gordon

This split-post fixator is designed for precise model positioning and holding. Special conical arms allow for the removal and replacement of models.

REF 072-004-00

1 piece



Polyclav®

Pressure vessel for polymerization of cold-curing polymers, induction-compatible.

With safety valve.

Without pressure hose.

Capacity: 4.5 l

REF 070-010-00

1 piece



Pressure hose

For Polyclav® (REF 070-010-00).

Connection plug nipple to fit inlet valve and standard coupling for compressed air sockets.

Length: 1.5 m

Pmax: 8.0 bar

REF 070-015-00

1 piece



Heating plate

For Polyclav® (REF 070-000-00 and 070-010-00).

Power voltage	REF	Quantity
230 V~	070-900-00	1 piece

Technical data

Output	125 W
Thermostat for water temperature	43 °C / 109.4 °F
Power frequency	50 Hz / 60 Hz
Weight	2.6 kg
Dimensions (W x H x D)	240 mm x 240 mm x 75 mm



Accessories



Suction apparatus with charcoal filter

For protection against vapor during processing of acrylics in the laboratory.

- controlled by a photoelectric sensor (automatic switch on/off)
- includes hour meter
- easy-change filter (replace the filter cartridge, REF 072-001-40, after approx. 50 operating hours)

Power voltage	REF	Quantity
230 V~	072-001-00	1 piece
115 V~	072-002-00	1 piece

Technical data

Output	100 W
Power frequency	50 Hz
Weight	6 kg
Dimensions (W x H x D)	220 mm x 260 mm x 260 mm



activeblue cleaning unit

With activefluid cleaning fluid.

For the gentle removal of plaque from orthodontic removable appliances.

Included:

- 1 x cleaning fluid, 1000 ml (REF 150-020-00)
- 1 x cleaning cup, white
- 1 x cleaning pins, 75 g (REF 098-001-00)
- 1 x tweezers
- 1 x holding magnet
- 1 x power cable (230 V)

Technical data

Power voltage	100 V~ - 240 V~ at 50 Hz / 60 Hz
Turning rotations	1200 min ⁻¹ / 1600 min ⁻¹ / 2000 min ⁻¹
Electric fuse	T 1.6 A L 250 V
Noise level	< 70 dB (A)
Weight	3 kg
Dimensions (W x H x D)	150 mm x 240 mm x 280 mm

REF 098-000-00	1 piece
-----------------------	---------



Cleaning pins for activeblue

For activeblue cleaning unit.

REF 098-001-00	1 container 75 g each
-----------------------	-----------------------

activefluid cleaning fluid

For activeblue cleaning unit.

REF 150-020-00	1000 ml
----------------	---------



Dermatril®, disposable rubber gloves

For protection during work with polymethyl methacrylates such as Orthocryl®.

Made from special nitrile, free from powder and natural latex.

Layer thickness	Size	REF	Quantity
0.10 mm	S	162-907-00	100 pieces
0.10 mm	M	162-908-00	100 pieces
0.10 mm	L	162-909-00	100 pieces



Premium-Line pliers / Fixed appliance technique	Page 372
--	----------

Premium-Line pliers / Removable appliance technique	Page 380
--	----------



Premium-Line pliers.

Premium quality – premium function.



DENTAURUM
**ONLINE
SHOP**
shop.dentaurum.com






One program to cover all areas.

Whether for orthodontics or for prosthetics, the Dentaureum pliers program includes instruments to serve all standard requirements. The user can always trust in consistent quality and reliability. This ensures quicker, more effective use.

Made in Germany.

Highly experienced craftsmanship and state-of-the-art machinery have been used to produce these high-quality pliers which remain true to the tag line "Made in Germany". For this reason, Dentaureum's Premium-Line pliers are extra long-lasting*.






PREMIUM
LINE

10 year guarantee!**



*For instructions on how to clean and take care of the Dentaaurum pliers and instruments, please refer to our separate information sheet REF 989-590-20. Repair service available on manufacturer's site (e.g. wear).

** The 10 year guarantee is on material and workmanship, loosening of cutting inserts, breakage and corrosion of pliers' body, providing the instruments have been used correctly and properly maintained. The guarantee does not apply if the equipment has been manipulated in any way.

Premium-Line pliers / Fixed appliance technique



Ortho-Cast M-Series cap removal pliers sterilizable

With replaceable blade.

Only suitable for use with Ortho-Cast M-Series buccal tubes.

Components supplied:

Pliers, fitted with a blade plus two replacement blades and one hex key.

Premium-Line · 10 year guarantee

REF 003-805-00

1 piece



Blade sterilizable

For Ortho-Cast M-Series cap removal pliers, REF 003-805-00.

Components supplied:

10 blades + 2 screws

REF 038-805-10

10 pieces



Special bracket debonding pliers, 45° angled sterilizable

For the removal of brackets and adhesive residues.

45° angled cutting surface, specially ground to fit purpose.

Hard metal inserts.

Premium-Line · 10 year guarantee

REF 004-350-00

1 piece



Bracket removing pliers, straight sterilizable

For removing brackets between bracket base and adhesive.

With hard metal inserts.

Length: 13.0 cm

Premium-Line · 10 year guarantee

REF 004-346-00

1 piece



Bracket removing pliers sterilizable

For removing brackets between bracket base and adhesive.

Premium-Line · 10 year guarantee

REF 003-349-00

1 piece



Plastic insert sterilizable

For bracket removing pliers, REF 003-349-00 and 045-005-00.

REF 038-356-00

4 pieces

Premium-Line pliers / Fixed appliance technique

Bracket removing pliers, angled sterilizable

For easy removal of lingual brackets and brackets in the posterior region. Inserts made of hard metal.

Premium-Line · 10 year guarantee

REF 004-347-00

1 piece



Oliver-Jones band removing pliers sterilizable

With replaceable plastic tip. With hard metal insert.

Premium-Line · 10 year guarantee

REF 026-347-00

1 piece



Plastic tip sterilizable

For Oliver-Jones band removing pliers, REF 026-347-00 and 045-004-00.

REF 039-357-00

10 pieces



Johnson contouring pliers sterilizable

Curvature of the beaks only on one side. For contouring and adapting bands and crowns.

Premium-Line · 10 year guarantee

REF 003-114-00

1 piece



Separating pliers sterilizable

For expanding Dentalastics® separating rings before positioning.

Premium-Line · 10 year guarantee

REF 015-160-00

1 piece



How pliers, straight sterilizable

Universal pliers for the multi-bracket technique. Due to the slender working tips these pliers are perfect for positioning bands with prewelded brackets. Also in the lower anterior area. **Not suitable for debonding brackets.**

Premium-Line · 10 year guarantee

REF 003-112-00

1 piece



Premium-Line pliers / Fixed appliance technique



How pliers, curved sterilizable

These curved pliers are especially suited for the posterior region. **Not suitable for debonding brackets.**

Premium-Line · 10 year guarantee

REF 003-111-00

1 piece



Side cutter Maxi sterilizable

Inserts made of hard metal.

We do not recommend using the front third of the cutters for cutting wires with the maximum indicated thickness!

Max. wire thickness: $\varnothing$ 1.0 mm, hard

Premium-Line · 10 year guarantee

REF 004-266-00

1 piece



Distal end cutter Mini sterilizable

For shortening the archwire ends in the multi-bracket technique. Delicate design. Wire holder.

Inserts made of hard metal.

Length: 12.0 cm

Max. wire thickness: $\boxtimes$ 0.56 mm x 0.70 mm / 22 x 28

Premium-Line · 10 year guarantee

REF 003-701-00

1 piece



Distal end cutter Mini sterilizable

With longer handle for easier handling when shortening archwire ends. Delicate design. Wire holder.

Inserts made of hard metal.

Length: 14.0 cm

Max. wire thickness: $\boxtimes$ 0.56 mm x 0.70 mm / 22 x 28

Premium-Line · 10 year guarantee

REF 003-702-00

1 piece



Distal end cutter Maxi sterilizable

For shortening the archwire ends in the multi-bracket technique. Robust design. Wire holder.

Inserts made of hard metal.

Length: 12.0 cm

Max. wire thickness: $\boxtimes$ 0.56 mm x 0.70 mm / 22 x 28

Premium-Line · 10 year guarantee

REF 003-700-00

1 piece



Flush Cutter Mini, with silicone O-ring sterilizable

For shortening round wires in the multi-bracket technique.

Elastic wire holder due to silicone O-ring, REF 037-356-00.

Inserts made of hard metal.

Length: 11.7 cm

Max. wire thickness: $\varnothing$ 0.50 mm / 20

Premium-Line · 10 year guarantee

REF 003-705-00

1 piece

Premium-Line pliers / Fixed appliance technique

Flush Cutter Maxi, with silicone O-ring sterilizable

With longer handle for easier handling when shortening round wires in the multi-bracket technique.
Elastic wire holder due to silicone O-ring, REF 037-356-00.
Inserts made of hard metal.

Length: 14.2 cm

Max. wire thickness: $\varnothing$ 0.50 mm / 20

Premium-Line · 10 year guarantee

REF 003-706-00

1 piece



Silicone O-ring sterilizable

For Flush Cutter Mini and Maxi, REF 003-705-00 and 003-706-00.

REF 037-356-00

10 pieces



Ligature cutter Mini sterilizable

Accuracy right to the tip for smooth cutting of soft ligature wires. With hard metal inserts.
We do not recommend using the front third of the cutters for cutting wires with the maximum indicated thickness!

Max. wire thickness: $\varnothing$ 0.4 mm, soft

Premium-Line · 10 year guarantee

REF 014-151-00

1 piece



Ligature cutter 45° sterilizable

45° angled cutting tips for easy cutting of ligatures, especially in the posterior areas. With hard metal inserts.
We do not recommend using the front third of the cutters for cutting wires with the maximum indicated thickness!

Max. wire thickness: $\varnothing$ 0.4 mm, soft

Premium-Line · 10 year guarantee

REF 014-152-00

1 piece



Side cutter Medium sterilizable

Black DLC coating. Smooth cuts guaranteed by accurate cutting edges. Inserts made of hard metal.
We do not recommend using the front third of the cutters for cutting wires with the maximum indicated thickness!

Max. wire thickness: $\varnothing$ 0.7 mm, hard

Premium-Line · 10 year guarantee

REF 016-153-00

1 piece



Light Wire pliers sterilizable

For especially fine wire bending in the light wire technique.
Max. wire thickness: $\varnothing$ 0.5 mm, spring hard

Premium-Line · 10 year guarantee

REF 003-130-00

1 piece



Premium-Line pliers / Fixed appliance technique



Begg Light Wire pliers

sterilizable

For especially fine wire bending in the light wire technique. Additional cutting device. Inserts made of hard metal.

Max. wire thickness: $\varnothing$ 0.5 mm, spring hard

Premium-Line · 10 year guarantee

REF 004-130-00

1 piece



Angle / Tweed ribbon arch pliers

sterilizable

For applying torque to square wires. Parallel grip prevents crushing of wire. Inserts made of hard metal.

Max. wire thickness: $\square$ 0.56 mm x 0.70 mm / 22 x 28

Premium-Line · 10 year guarantee

REF 003-442-00

1 piece



Weingart universal pliers

sterilizable

For positioning the arch and bending the arch wire ends. Accurately closing serrated tips. With hard metal inserts. Also suitable for debonding brackets.

Max. wire thickness: $\varnothing$ 0.5 mm, hard

Premium-Line · 10 year guarantee

REF 003-120-00

1 piece



Weingart universal pliers Mini

sterilizable

Very slender design. For holding and positioning the arch. Accurately closing serrated tips.

Not suitable for debonding brackets.

Max. wire thickness: $\varnothing$ 0.5 mm, hard

Premium-Line · 10 year guarantee

REF 003-121-00

1 piece



Lingual bending pliers

sterilizable

For bending steps and compensation bends in archwires for the lingual technique.

Steps: 1.5 mm; 2.0 mm; 2.5 mm and 3.0 mm

Max. wire thickness: $\varnothing$ 0.5 mm, hard

Premium-Line · 10 year guarantee

REF 003-001-00

1 piece



Nance loop bending pliers

sterilizable

With four-step beaks for bending loops of different sizes in round and rectangular wires.

Max. wire thickness: $\varnothing$ 0.7 mm, hard

Premium-Line · 10 year guarantee

REF 000-001-00

1 piece

Premium-Line pliers / Fixed appliance technique

Jarabak Light Wire pliers

sterilizable

With guiding grooves.

Max. wire thickness: $\varnothing$ 0.5 mm, spring hard

Premium-Line · 10 year guarantee

REF 003-125-00

1 piece



Tweed loop forming pliers

sterilizable

For especially fine bendings.

Max. wire thickness: $\varnothing$ 0.5 mm, hard

Premium-Line · 10 year guarantee

REF 012-350-00

1 piece



Tweed loop forming pliers

sterilizable

For especially fine bendings.

With individually replaceable hard metal tip.

Max. wire thickness: $\varnothing$ 0.5 mm, hard

Premium-Line · 10 year guarantee

REF 012-351-00

1 piece



Tip

sterilizable

For Tweed loop forming plier, REF 012-351-00.

REF 038-357-00

3 pieces



De La Rosa contouring pliers

sterilizable

With guiding grooves for forming round and square arches 16, 18, 22.

Premium-Line · 10 year guarantee

REF 004-109-00

1 piece



Lingual arch bending pliers

sterilizable

For individual construction and activation of lingual/palatal bars. Using these pliers, the bent tips of the wire will fit into all Dentaurem lingual/palatal tubes $\varnothing$ 2 x 0.9 mm.

Max. wire thickness: $\varnothing$ 0.9 mm, hard

Premium-Line · 10 year guarantee

REF 003-185-00

1 piece



Premium-Line pliers / Fixed appliance technique



Crimping pliers, for Herbst TS / SUS sterilizable

For clamping slotted spacers for Herbst TS and SUS.
Premium-Line · 10 year guarantee

REF 003-710-00

1 piece



Crimping pliers sterilizable

Tip with V-shaped notch for holding and crimping
sliding/clamping and stop tubes.

Premium-Line · 10 year guarantee

REF 003-711-00

1 piece



Angle wire bending pliers standard sterilizable

Sturdy constructed pliers with one cone and one
pyramid shaped beak.

Max. wire thickness: $\varnothing$ 0.7 mm, hard

Premium-Line · 10 year guarantee

REF 003-139-00

1 piece



Angle wire bending pliers sterilizable

Sturdy constructed pliers with one cone and one
pyramid shaped beak. Inserts made of hard metal.

Max. wire thickness: $\varnothing$ 0.7 mm, hard

Premium-Line · 10 year guarantee

REF 004-139-00

1 piece



Facebow bending pliers Magnum sterilizable

Ideal for bending the outer bows of a facebow.

Max. wire thickness: $\varnothing$ 1.8 mm

Premium-Line · 10 year guarantee

REF 003-180-01

1 piece

Premium-Line pliers / Fixed appliance technique

Hollow chop pliers Extra-Mini

sterilizable

Ideally suited for forming different loops and adjusting preshaped arches for the fixed appliance technique.

Max. wire thickness: $\varnothing$ 0.5 mm, hard

Premium-Line · 10 year guarantee

REF 003-350-00

1 piece



Hollow chop pliers Mini

sterilizable

For loop forming and bending of round and rectangular wires.

Max. wire thickness: $\varnothing$ 0.5 mm, hard

Premium-Line · 10 year guarantee

REF 003-351-00

1 piece



Hammer head pliers

sterilizable

For bending nickel-titanium arches/wires in the molar region.

The two ends fit together like a tongue and groove around the nickel titanium wire, so that the ends of the wire attain the required bend with one squeeze of the pliers.

Max. wire thickness:

$\varnothing$ 0.45 mm / 18, ∇ 0.41 mm x 0.56 mm / 16 x 22

Premium-Line · 10 year guarantee

REF 003-355-00

1 piece



Aderer 3-prong pliers Extra-Mini

sterilizable

For bending wires/arches for the fixed appliance technique.

Specially suited for bending NiTi wires/arches such as rematitan® LITE and Tensic®.

Max. wire thickness:

$\varnothing$ 0.5 mm, hard, ∇ 0.48 mm x 0.64 mm / 19 x 25

Premium-Line · 10 year guarantee

REF 003-205-00

1 piece



Premium-Line pliers / Removable appliance technique



Hollow chop pliers Medium

For bending small springs, eyelets and loops.

Max. wire thickness: $\varnothing$ 0.7 mm, hard

Premium-Line · 10 year guarantee

REF 003-053-00

1 piece



Hollow chop pliers Maxi

For bending springs, loops and space maintainers.

Max. wire thickness: $\varnothing$ 0.8 mm, hard

Premium-Line · 10 year guarantee

REF 003-233-00

1 piece



Ruhland loop forming pliers Magnum

For bending coffin springs and connecting springs, labial arches and U-loops.

Max. wire thickness: $\varnothing$ 1.1 mm, hard

Premium-Line · 10 year guarantee

REF 001-650-00

1 piece



Adams pliers Medium

With two smooth, rectangular beaks. **With guiding grooves.**

Length: 12.0 cm

Max. wire thickness: $\varnothing$ 0.7 mm, hard

Premium-Line · 10 year guarantee

REF 003-514-00

1 piece



Adams pliers Medium

With two smooth, rectangular beaks. Hard metal coated working faces.

Length: 12.0 cm

Max. wire thickness: $\varnothing$ 0.7 mm, hard

Premium-Line · 10 year guarantee

REF 003-515-00

1 piece



Adams pliers Medium

With two smooth, rectangular beaks.

Inserts made of hard metal. With guiding groove.

Length: 12.0 cm

Max. wire thickness: $\varnothing$ 0.7 mm, spring hard

Premium-Line · 10 year guarantee

REF 004-515-00

1 piece

Premium-Line pliers / Removable appliance technique

Adams pliers Maxi

With two smooth, rectangular beaks.

Length: 14.5 cm

Max. wire thickness: $\varnothing$ 0.8 mm, hard

Premium-Line · 10 year guarantee

REF 003-094-00

1 piece



Adams pliers Medium

With one smooth, rectangular and one smooth, rounded beak.

Length: 14.0 cm

Max. wire thickness: $\varnothing$ 0.7 mm, hard

Premium-Line · 10 year guarantee

REF 003-516-00

1 piece



Adams pliers Maxi

With one smooth, rectangular and one smooth, rounded beak. Inserts made of hard metal.

Length: 16.0 cm

Max. wire thickness: $\varnothing$ 0.8 mm, hard

Premium-Line · 10 year guarantee

REF 004-095-00

1 piece



Wire bending pliers, wide

Serrated tips. With guiding groove.

Max. wire thickness: $\varnothing$ 0.8 mm, hard

Premium-Line · 10 year guarantee

REF 033-452-00

1 piece



Wire bending pliers Maxi

Serrated tips. With guiding groove. Universal bending and cutting pliers. Inserts made of hard metal.

Max. wire thickness: $\varnothing$ 0.9 mm, hard

Cuts up to: $\varnothing$ 0.7 mm, hard

Premium-Line · 10 year guarantee

REF 013-522-00

1 piece



Universal wire bending pliers

With serrated tips and guiding grooves.

Hard metal coated working faces.

Max. wire thickness: $\varnothing$ 0.8 mm, hard

Premium-Line · 10 year guarantee

REF 003-290-00

1 piece



Premium-Line pliers / Removable appliance technique



Pointed beak pliers

Smooth beaks with guiding groove.

Inserts made of hard metal.

Max. wire thickness: $\varnothing$ 0.7 mm, spring hard

Premium-Line · 10 year guarantee

REF 004-455-00

1 piece



Pointed beak pliers

Serrated beaks with guiding groove. Fine pointed.

Max. wire thickness: $\varnothing$ 0.7 mm, hard

Premium-Line · 10 year guarantee

REF 035-455-00

1 piece



Young loop pliers

Universal pliers for different size loops. Ideal for labial arches.

Max. wire thickness: $\varnothing$ 0.7 mm, hard

Premium-Line · 10 year guarantee

REF 012-074-00

1 piece



Aderer 3-Prong pliers Mini

Length: 12.5 cm

Max. wire thickness: $\varnothing$ 0.7 mm, hard

Premium-Line · 10 year guarantee

REF 003-201-00

1 piece



Aderer 3-Prong pliers Medium

Length: 12.5 cm

Max. wire thickness: $\varnothing$ 0.9 mm, hard

Premium-Line · 10 year guarantee

REF 003-200-00

1 piece



Aderer 3-Prong pliers Maxi

Sturdy design.

Length: 15.0 cm

Max. wire thickness: $\varnothing$ 1.5 mm, hard

Premium-Line · 10 year guarantee

REF 003-202-00

1 piece

Premium-Line pliers / Removable appliance technique

Arrow clasp former

Max. wire thickness: $\varnothing$ 0.7 mm, hard

Premium-Line · 10 year guarantee

REF 012-520-00

1 piece




Arrow clasp bending pliers

Max. wire thickness: $\varnothing$ 0.7 mm, hard

Premium-Line · 10 year guarantee

REF 012-521-00

1 piece



EQ-Line pliers / Fixed appliance technique	Page 386
EQ-Line pliers / Removable appliance technique	Page 388
Instruments	Page 390
Accessories	Page 398



EQ-Line pliers. Dental instruments.

Good quality pliers at reasonable prices.
A versatile range for the orthodontic practice
and laboratory.



DENTAURUM
**ONLINE
SHOP**
shop.dentaurum.com



EQ-Line pliers / Fixed appliance technique



Bracket removing pliers

sterilizable

For removing brackets.

REF 045-005-00

1 piece



Plastic insert

sterilizable

For bracket removing pliers, REF 003-349-00 and 045-005-00.

REF 038-356-00

4 pieces



Distal end cutter

sterilizable

For shortening the archwire ends in the multi-bracket technique. Wire holder. Inserts made of hard metal.

Max. wire thickness: $\square$ 0.56 mm x 0.70 mm / 22 x 28

REF 045-021-00

1 piece



Side cutter Medium

sterilizable

For cutting wires and arches. Inserts made of hard metal.

We do not recommend using the front third of the cutters for cutting wires with the maximum indicated thickness!

Max. wire thickness: $\varnothing$ 0.7 mm, hard

REF 045-022-00

1 piece



Oliver-Jones band removing pliers

sterilizable

With replaceable plastic tip.

REF 045-004-00

1 piece



Plastic tip

sterilizable

For Oliver-Jones band removing pliers, REF 026-347-00 and 045-004-00.

REF 039-357-00

10 pieces



How pliers, straight

sterilizable

Universal pliers for the multi-bracket technique.

The special tip design makes these pliers suited for seating individual anterior bands.

REF 045-006-00

1 piece

EQ-Line pliers / Fixed appliance technique

Angle / Tweed ribbon arch pliers

sterilizable

For applying torque to square wires. Parallel retention prevents compression of the wire.

Max. wire thickness: ∇ 0.56 mm x 0.70 mm / 22 x 28

REF 045-000-00

1 piece



Weingart universal pliers

sterilizable

For positioning the arch and bending the arch wire ends.

Max. wire thickness: $\emptyset$ 0.5 mm, hard

REF 045-003-00

1 piece



Nance loop bending pliers

sterilizable

With four-step beaks for producing loops of different sizes in round and rectangular wires.

Max. wire thickness: $\emptyset$ 0.7 mm, spring hard

REF 045-010-00

1 piece



Jarabak Light Wire pliers

sterilizable

Max. wire thickness: $\emptyset$ 0.5 mm, spring hard

REF 045-009-00

1 piece



Tweed loop forming pliers

sterilizable

For especially fine bendings.

Max. wire thickness: $\emptyset$ 0.5 mm, hard

REF 045-001-00

1 piece



Angle wire bending pliers

sterilizable

Sturdy constructed pliers with one cone and one pyramid shaped beak.

Max. wire thickness: $\emptyset$ 0.7 mm, hard

REF 045-002-00

1 piece



EQ-Line pliers / Removable appliance technique



Loop forming pliers Medium

Max. wire thickness: $\varnothing$ 0.7 mm, hard

REF 045-012-00

1 piece



Loop forming pliers Maxi

E.g. for bending space maintainers.

Max. wire thickness: $\varnothing$ 0.8 mm, hard

REF 045-013-00

1 piece



Ruhland loop forming pliers Magnum

For bending coffin springs and connecting springs, labial arches and loops.

Max. wire thickness: $\varnothing$ 1.1 mm, hard

REF 045-014-00

1 piece



Adams pliers Medium

With two smooth rectangular beaks.

With guiding groove.

Max. wire thickness: $\varnothing$ 0.7 mm, hard

REF 045-015-00

1 piece



Adams pliers Maxi

With two smooth rectangular beaks.

With groove.

Max. wire thickness: $\varnothing$ 0.8 mm, hard

REF 045-016-00

1 piece



Adams pliers Maxi

With one smooth, rectangular and

one smooth, rounded beak.

Max. wire thickness: $\varnothing$ 0.8 mm, hard

REF 045-017-00

1 piece



Wire bending pliers Maxi

Max. wire thickness: $\varnothing$ 0.9 mm, hard

Cuts up to: $\varnothing$ 0.7 mm, hard

REF 045-019-00

1 piece

EQ-Line pliers / Removable appliance technique

Pointed beak pliers

Smooth beaks with guiding groove. Fine pointed.
Max. wire thickness: ø 0.7 mm, hard

REF 045-018-00

1 piece



Young loop pliers

Universal pliers for different size loops. Ideal for labial arches.
Max. wire thickness: ø 0.7 mm, hard

REF 045-011-00

1 piece



Aderer 3-Prong pliers Extra-Mini

Max. wire thickness: ø 0.8 mm, hard

REF 045-020-00

1 piece



Aderer 3-Prong pliers Mini

Max. wire thickness: ø 0.7 mm, hard

REF 045-007-00

1 piece



Aderer 3-Prong pliers Medium

Max. wire thickness: ø 0.9 mm, hard

REF 045-008-00

1 piece



Heavy gauge side cutter Magnum

Little force required for cutting through large size wires
due to special design.

We do not recommend using the front third of the cutters
for cutting wires with the maximum indicated thickness!

Max. wire thickness: ø 1.3 mm, hard

REF 044-177-00

1 piece



Instruments



Crimping tweezers sterilizable

For bracket placement. The blade at the end of the tweezers can be used for angulating brackets. Delicate design.

REF 025-277-00

1 piece



Crimping tweezers sterilizable

With long placement and angulation aid. For bracket placement. The blade at the end of the tweezers can be used for angulating brackets.

REF 025-276-00

1 piece



Crimping tweezers sterilizable

For posterior teeth.

Specially suited for placing brackets and buccal tubes in the posterior area.

REF 025-279-00

1 piece



Crimping tweezers sterilizable

Delicate tweezer design for easy bracket and buccal tube placement.

REF 025-280-00

1 piece



discovery® pearl debonding instrument sterilizable

For easy and safe removal of **discovery® pearl ceramic brackets**.

REF 019-001-00

1 piece



sl-instrument, for discovery® sl sterilizable

For opening and closing self-ligating brackets **discovery® sl** (Roth system), **dinamique® m** and **dinamique® c**, and for easier ligation of the archwire.

REF 023-276-00

1 piece



sl-instrument, for discovery® sl 2.0 sterilizable

For opening and closing self-ligating brackets **discovery® sl 2.0** (McLaughlin-Bennett-Trevisi* system) and easier ligation of the archwire.

REF 023-277-00

1 piece

* The Dentaaurum version of this prescription is not claimed to be a duplication of any other, nor does Dentaaurum imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.

Scaler, double-ended

sterilizable

For removing adhesive and cement residues.
With serrated band seater.

REF 027-349-00

1 piece



Separating strips*

sterilizable

Width: 4 mm

Coating	REF	Quantity
double-sided	018-001-01	12 pieces
one-sided	018-002-01	12 pieces



NEW

Scissors

For cutting thin metal bands.

Version	REF	Quantity
straight	007-024-00	1 piece
curved	008-025-00	1 piece



Band seating instrument, plastic

sterilizable to a maximum of 180 °C / 356 °F

With hard metal, triangular serrated tip.

REF 026-355-00

1 piece



Band pusher

acc. to Mershon

sterilizable

With serrated rectangular tip and anatomically formed grip
for good transfer of force when positioning bands.

REF 024-300-01

1 piece



DENTALURUM
ONLINE SHOP
shop.dentalurum.com



Instruments



Bracket positioning gauge sterilizable

Simplifies seating the brackets on the teeth in the correct distance between bracket slot and incisal edge.

With metal tips.

REF 030-390-00

1 piece



Needle holder Medium

acc. to Mathieu

sterilizable

Ideal for ligating. Convenient and practical design with serrated tips.

REF 000-030-00

1 piece



Needle holder Medium, with finger ring

acc. to Bachmann

sterilizable

Eases ligating with the middle finger under constant pull. Can also be used in ligating jaw fracture splints. With serrated tips.

REF 000-035-00

1 piece



Needle holder Mini

acc. to Mathieu

sterilizable

With especially delicate, serrated tips. For Dentalastics® plastic ligatures. **Not for wire ligatures.**

REF 000-031-00

1 piece



Mosquito forceps, without hook

sterilizable

Crown forceps or for placing Dentalastics® plastic ligatures.

REF 000-730-00

1 piece



Mosquito forceps, with hook

sterilizable

For placing Dentalastics® plastic ligatures.

REF 000-731-00

1 piece

Mosquito forceps, curved

sterilizable

For placing and ligating Dentalastics® plastic ligatures, especially for the lingual technique.

REF 000-732-00

1 piece



Dental mirror with reference lines (assembled), lingual

sterilizable

REF 023-279-00

1 piece



Dental mirror with reference lines (without grip), lingual

sterilizable

REF 023-279-10

1 piece



Positioning tool, lingual

sterilizable

REF 023-278-00

1 piece



Ligator

sterilizable

For fast ligating of the complete arch and figure-eight ligatures. Y-shaped tip simplifies the twisting of the ligatures.

Max. wire thickness: ø 0.3 mm, soft

REF 023-292-00

1 piece



Accessories (for ligator REF 023-292-00)

sterilizable

Spare part	REF	Quantity
Hard metal tip	023-293-00	1 piece
Spool	023-294-00	1 piece
Screw	023-295-00	1 piece



Instruments



Ligature tightener and adapter sterilizable

Ligature adjuster and tucker for long preformed ligatures.

REF 022-274-00

1 piece



Ligature adapter, double ended sterilizable

For bending and adapting wire ligatures.

REF 021-273-00

1 piece



Distal end bender, double-ended

acc. to Langer

sterilizable

Sturdy instrument for bending arch wire ends.

Handle ø: 8 mm

REF 019-000-00

1 piece



Arch turret without torque

For forming flat, plain arches made of rectangular wires.

For rectangular wires:

☑ 0.41 mm / 16

☑ 0.43 mm / 17

☑ 0.46 mm / 18

☑ 0.51 mm / 20

☑ 0.56 mm / 22

Color: blue

REF 004-180-00

1 piece



Arch turret with torque / 16

For forming flat or angulated (torque) rectangular arches.

☑ 0.41 mm x 0.41 mm / 16 x 16

in 0°, 7°, 10°, 13° and 16° Torque

Color: red

REF 004-186-00

1 piece

Arch turret with torque / 16 x 16

For forming flat or angulated (torque) rectangular arches.

☐ 0.41 mm x 0.41 mm / 16 x 16

in 0°, 10°, 16° and 23° Torque

Specially designed for the Ricketts® technique.

Color: green

REF 004-187-00

1 piece



Arch turret with torque / 18

For forming flat or angulated (torque) rectangular arches.

☐ 0.46 mm x 0.64 mm / 18 x 25

in 0°, 7°, 10°, 13° and 16° Torque

Color: gold

REF 004-182-00

1 piece



Arch turret with torque / 22

For forming flat or angulated (torque) rectangular arches.

☐ 0.56 mm x 0.70 mm / 22 x 28

in 0°, 7°, 10°, 13° and 16° Torque

Color: silver

REF 004-184-00

1 piece



LinguForm arch turret without torque

For forming flat, plain arches made of round or rectangular wires for the lingual technique.

ø 0.35 mm / 14, ☐ 0.35 mm / 14

ø 0.41 mm / 16, ☐ 0.41 mm / 16

ø 0.46 mm / 18, ☐ 0.46 mm / 18

ø 0.48 mm / 19, ☐ 0.48 mm / 19

ø 0.51 mm / 20, ☐ 0.51 mm / 20

Color: silver

REF 004-190-00

1 piece



Measuring gauge

With drag indicator for measuring tensile and compressive forces. Suitable for measuring on the model.

Measuring range	REF	Quantity
25 g - 250 g	040-711-00	1 piece
50 g - 500 g	040-712-00	1 piece
100 g - 1000 g	040-713-00	1 piece



Instruments



Dental vernier

acc. to Beerendonk

With fine, long tips. Suitable for measurements on the model.

REF 042-750-00

1 piece



Dental vernier (Münchner design)

With extra delicate tips. Suitable for width and distance measurements on the model.

REF 042-751-00

1 piece



Ruler

acc. to Bernklau

sterilizable

Material: stainless steel

Perfect for all diagnostic measurements, both intra-oral and on the model.

REF 044-730-01

1 piece



Wedge gauge

For checking the gap width on palatal expanders.

REF 043-720-00

1 piece



Ruler

Material: plastic

For measuring labial arches, Adams clasps, palatal bars, etc.

Color: white

REF 766-000-00

1 piece



Ruler (Münchner design)

Material: stainless steel

With $\frac{1}{2}$ mm and mm increments.

Length: 15.0 cm

REF 044-731-00

1 piece



Measuring grid

acc. to Schmuth

For parallax-free, optical measurement of the dental arch.

When looking through the grid the red sighting grid should

be covered by the corresponding black measuring lines

to ensure parallax compensation for precise measurement of distances.

REF 003-366-00

1 piece

Instrument set acc. to Korkhaus

Set with orthometer, mm disk, cm disk and bow divider with sleeve.

REF 029-360-02

1 set



Orthometer

For Instrument set acc. to Korkhaus (REF 029-360-01 and 029-360-02).

REF 029-364-00

1 piece



Bow divider, with sleeve

For Instrument set acc. to Korkhaus (REF 029-360-01 and 029-360-02).

REF 029-361-02

1 piece




Centimeter disc

For Instrument set acc. to Korkhaus (REF 029-360-01 and 029-360-02).

REF 029-362-00

1 piece




Millimeter disc

For Instrument set acc. to Korkhaus (REF 029-360-01 and 029-360-02).

REF 029-363-00

1 piece



Three-dimensional bow divider

acc. to Korkhaus

Allows the direction and extent of tooth movement to be checked by a single measurement on the model.
Not approved for intra-oral use!

REF 028-353-00

1 piece



Bow divider

Standard design. For surveying the model and measuring tooth widths by positioning on the model.
Not approved for intra-oral use!

REF 030-395-03

1 piece



Tip for bow divider


For bow divider (REF 030-395-02 and 030-395-03).

REF 030-396-02

1 piece



Accessories



Tracing-protractor template, standard

Material: plastic
With protractor and tooth symbols.
For evaluation of radiographs.
Dimensions: 207 mm x 132 mm
Color: transparent


REF 004-368-00	1 piece
----------------	---------



Tracing-protractor template

acc. to Ricketts®
Material: plastic
For evaluation of radiographs. With mm scale, protractor and cut out tooth symbols for tracing.
Dimensions: 197 mm x 74 mm
Color: transparent, clear

REF 004-369-00	1 piece
----------------	---------



Penta arch template

acc. to Ricketts®
Material: plastic
For determining the Penta archwire shape.
Dimensions: 30.3 cm x 21.6 cm
Color: transparent, clear

REF 004-373-00	1 piece
----------------	---------



Tracing foil

For radiograph tracings. One side mat, other side shiny. Tear-off pad.

Dimensions	REF	Quantity
17.5 cm x 23.5 cm	075-100-00	100 pieces
24 cm x 30 cm	075-101-00	100 pieces

Maintenance and protection oil

For the maintenance of pliers and instruments.
Application in particular before sterilisation.

REF 055-201-00

50 ml



Stand for wires and pliers

Material: plastic

With 6 tubes for straight wires.

Do not use disinfectants containing ethanol, phenol or alcohol.

Color: transparent, clear

Height: approx. 13 cm

Width: approx. 18 cm

REF 055-101-01

1 piece



Stand for pliers

Material: plastic

Do not use disinfectants containing ethanol, phenol or alcohol.

Color: transparent, clear

suitable for	Length	REF	Quantity
approx. 7 pliers	13.0 cm	055-000-00	1 piece
approx. 14 pliers	25.5 cm	055-100-00	1 piece



Stand for pliers

Material: plastic

For 15 pliers.

Do not use disinfectants containing ethanol, phenol or alcohol.

Dimensions: 30 cm x 20 cm x 6.5 cm (W x D x H)

Color: white

REF 055-105-00

1 piece



Accessories



Bur stand magnetic holder

Partly plastic.

- permanent magnets ensure firm retention of the instruments
- magnet holder without instruments
- do not use disinfectants containing ethanol, phenol or alcohol

REF 089-430-00

1 piece



Bur stand

Material: plastic

- for 12 burs or abrasive stones.
- Bur stand without instruments
- shaft ø: 2.35 mm
- do not use disinfectants containing ethanol, phenol or alcohol

Color: transparent, clear

REF 089-431-00

1 piece



Marking pencil

Wax pencil

Color: white

REF 043-700-00

6 pieces

Marking pencil

indelible



Color	REF	Quantity
black	043-800-00	1 piece
red	043-900-00	1 piece



Equipment.

For professional workmanship.



DENTAURUM
**ONLINE
SHOP**
shop.dentaurum.com





Certified quality made in Germany.

Dentaurum equipment has been developed, manufactured and tested in close collaboration with internal and external users. The focus is on easy handling and high quality, ensuring long durability.

Metal connections in orthodontics.

We offer a wide range of equipment for all types of metal connections in orthodontics, from mini gas burners to high-tech laser welding technology.



desktop Compact

Universal laser welding unit.

Laser Welder SL10

The entry-level model, ideal for orthodontics.

Spot welding unit

Assistent 3000.

Dual wheel model trimmer

Beautiful models – trimmed quickly and neatly.

Suction unit

For protection from vapor during processing.

activeblue

Gentle removal of plaque.

Spot welding unit.

For professional welding and soldering in orthodontics.

A suitable unit for everyone.

The Assistant 3000 is the right unit for welding and soldering tasks in orthodontics, working professionally and efficiently both in the orthodontic practice and in the laboratory. The function and selection of accessories available, along with the unit's ergonomic design, are all based upon years of experience in development and close collaboration with orthodontists.

Modern technology – excellent results.

The Assistant 3000 achieves excellent welding and soldering results due to the modern electronic components used together with high quality mechanical elements. For example, the electronic welding adjustments or the extra wide hand rests make work easy, comfortable and effective.

Assistant 3000 – the comfortable unit for average use.

With the Assistant 3000, excellent results can be efficiently achieved with average use.

This unit's special features, amongst others, are

- 6 welding settings and 2 soldering settings.



Dentaurum laser technology.

For professional welding in orthodontics.

Laser welding – the modern joining technique.

Laser welding allows more efficient processing. In comparison with traditional soldering, you can save up to 80 % in time. New fabrications and repairs can be carried out better and faster by using a dental laser.

Your advantages:

- gentle on materials.
- biocompatible.
- corrosion resistant.



Maximum stability.

Laser-welded connections present a significantly higher mechanical resistance and are more durable than soldered connections. This means greater security and higher quality for the orthodontist and the patients.

Broad spectrum of application.

It is impossible to imagine the world of joining techniques without laser welding. Compared to soldering, laser welding offers a much wider range of applications.

Whether you want to process precious metal alloys with a high gold content, non-precious alloys or titanium: you can cover the whole range of dental materials using the laser welding units from Dentaureum.

Ergonomics – extremely user-friendly.

A special characteristic feature of the Dentaureum dental lasers is their special ergonomic design. Working becomes a relaxed and easy affair for you.

- The dental lasers demonstrate their excellent ergonomics in the optimal inclination of the laser axis. Thereby not only do you sit relaxed in front of the microscope, but you can also hold the workpiece at an optimal angle.
- The workpiece is always flushed with argon through the argon nozzle. This ensures an oxidation-free weld with high stability.
- The desktop Compact fulfills the highest of demands in all areas of dental prosthetics and orthodontics. The Laser Welder SL10 is the ideal entry model for laser technology.

Equipment and accessories



Power voltage	REF	Quantity
200 V - 240 V	090-578-50	1 piece

Technical data

max. power	60 W
Peak pulse output	5000 W
Pulse energy	50 J
Pulse duration	0.5 ms - 20 ms
Pulse frequency	single pulse - 25 Hz
Welding spot diameter	0.2 mm - 2.0 mm, infinitely variable motorized adjustment
Memory settings	39
Pulse shaping	4 pre-set pulse shapes
Argon nozzle	1 x flexible, 1 x fixed retractable
Air nozzle	integrated in welding chamber
Extraction	integrated with suspended particle filter (H14)
Microscope	stereo microscope with 16x magnification
Lighting	2 x 3,2 W LED
max. working height	88 mm
Electric fuse	T 16 A
Power frequency	50 Hz - 60 Hz
Weight	42 kg
Dimensions (W x H x D)	510 mm x 500 mm x 630 mm

desktop Compact

Laser welding – powerful and easy to use.

More than 30 years of dental laser welding technique accumulated in a compact table unit. From delicate, detailed welding to secure joints, the desktop Compact laser can be used for all areas of application in dental prosthetics and orthodontics. The finely adjustable dosage makes welding so much easier. Until now, such compact machine measurements were not available in this laser performance category. Simple machine operation prevents fatigue during work.

Advantages:

- extremely compact table top unit
- clear and user-friendly touchscreen
- high output reserve
- ergonomically designed – large working chamber with inclining laser axis
- high-quality microscope with a wide field of vision
- integrated pulse forming – 4 pre-set pulse shapes
- two inert gas nozzles, 1x flexible and 1x fixed retractable
- large pulse adjustment range – easy operation using joysticks
- integrated extraction and cool air nozzle

Further product information is available on request.

Laser Welder SL10

The Laser Welder SL10 from Siro Lasertec complements the desktop Compact as an economical entry-level model for the laser welding technique. It is ideal for smaller dental laboratories and orthodontic practices.

Product advantages:

- excellent price-performance ratio
- large color touch display
- particularly suitable for orthodontic appliances
- welding power suitable for all fields of dentistry
- easy to operate
- with microscope
- solid housing
- adjustment of welding parameters on touch panel
- color camera shows welding process on the display
- desktop design



Power voltage	REF	Quantity
230 V~	090-594-00	1 piece

Technical data

max. power	60 W
Peak pulse output	3800 W
Pulse energy	33 J
Pulse duration	0.1 ms - 10 ms
Pulse frequency	0,5 Hz - 20 Hz
Welding spot diameter	0.2 mm - 2.0 mm, infinitely variable motorized adjustment
Protective gas	Argon 99.996 % (Argon 4.6)
Pulse shaping	no
Argon nozzle	1 x flexible
Extraction	integrated with suspended particle filter
Microscope	stereo microscope with 16x magnification
Lighting	2 x 3 W LED
max. working height	200 mm
Electric fuse	T 16 A
Power frequency	50 Hz / 60 Hz
Weight	40.5 kg
Dimensions (W x H x D)	440 mm x 480 mm x 600 mm



Equipment and accessories



Electrically adjustable stand

The answer to unsurpassed ergonomics with the Dentaaurum desktop laser-welding units.

Advantages:

- this stand offers the user the best possible ergonomics
- the electrically adjustable stand enables the technician to change to their comfort level, via a foot switch, within seconds. This offers the most comfortable viewing level, and assures full concentration on the laser work
- the stand is designed to guarantee fatigue-free laser welding, even after extended periods of time
- can be used with all Dentaaurum desktop laser-welding units

Power voltage	REF	Quantity
230 V~	090-574-00	1 piece

Technical data

Power frequency	50 Hz
min. working height	740 mm
max. working height	1090 mm



Argon fitting

For laser.

For connecting the laser to the argon bottle.
Flow rate shown in liters.

REF 090-404-00

1 assortment



Armrest cushion

For laser.

REF 090-513-10

2 pieces

CoCr welding wire

Various dimensions to be used as a filling soldering wire for all CoCr alloys.

Norm: DIN EN ISO 28319

ø	Length	REF	Quantity
0.25 mm / 10	2 m	528-215-10	1 piece
0.35 mm / 14	2 m	528-210-10	1 piece
0.50 mm / 20	2 m	528-200-10	1 piece

Composition (percentage by mass)

Co	Cr	Mo
65.0	28.0	5.5

Additional elements < 1 %: Mn, N, Si

**NiCr welding wire**

For laser welding NiCr based alloys. Specifically for high grade alloys such as remanium® CSe and G-soft.

Norm: DIN EN ISO 28319

ø	Length	REF	Quantity
0.50 mm / 20	2 m	528-220-00	1 piece

Composition (percentage by mass)

Ni	Cr	Mo	Nb
63.0	22.0	9.0	3,5

Additional elements < 1 %: Al, Co, Fe, Ti

**rematitan® wire on coil**

Titanium grade 1

For laser welding of titanium frameworks.

Recommended as filler material for laser welding, e.g. with titanium RPE appliances (Titan hyrax® Maxi – 12).

Norm: DIN EN ISO 28319, DIN EN ISO 5832-2

ø	Length	REF	Quantity
0.40 mm / 16	2 m	528-039-50	1 piece
0.70 mm / 28	2 m	528-040-50	1 piece



Equipment and accessories



rematitan® straight wire

Titanium grade 1

For laser welding.

Norm: DIN EN ISO 28319, DIN EN ISO 5832-2

ø	Length	REF	Quantity
1.00 mm / 39	each 100 mm	528-041-00	10 pieces
1.20 mm / 47	each 100 mm	528-042-00	10 pieces



rematitan® straight wire

Titanium grade 4

For laser welding.

Recommended for axes of swing locks.

Norm: DIN EN ISO 28319, DIN EN ISO 5832-2

ø	Length	REF	Quantity
1.50 mm / 59	50 mm	528-050-00	1 piece



rematitan® straight wire

Titanium grade 1

For laser welding.

Norm: DIN EN ISO 28319, DIN EN ISO 5832-2

ø rolled	Length	REF	Quantity
0.50 mm x 1.5 mm / 20 x 59	each 100 mm	528-043-00	10 pieces



Spot welding unit Assistent 3000

- button for continuous soldering
- LED display
- automatic voltage detection
- current limitation if electrode is incorrect
- power saving device (stand-by mode if not used for 10 minutes)
- error indications for:
 - overheating
 - mix-up of electrodes
- connection for hand electrode
- connection for foot switch

Hand electrodes, holder for hand electrodes and foot switch are not included!



Power voltage	REF	Quantity
110 V~ / 240 V~	078-000-00	1 piece

Technical data

Output	140 W
Welding power	up to 2000 A
Welding stages	6 stages
Soldering stages	2 stages, continuous soldering
Electric fuse	T 3.15 A
Power frequency	50 Hz / 60 Hz
Weight	8 kg
Dimensions (W x H x D)	335 mm x 130 mm x 310 mm



Equipment

Equipment and accessories



Foot switch

To be used with the spot welding units Junior 3000, Assistant 3000 and Master 3000.

REF 080-116-00

1 piece



Hand electrode

Suitable for all Dentaureum spot welding units.

REF 081-101-00

1 piece



Holder for hand electrode

For the spot welding units Junior 3000, Assistant 3000 and Master 3000.

Version	REF	Quantity
for right and left side	081-204-00	2 pieces
for right side	081-205-00	1 piece
for left side	081-206-00	1 piece



Soldering carbon for hand electrode

Suitable for all Dentaureum spot welding units.
Copper casing.

REF 081-601-00

1 piece



Electrode, notched, for hand electrode

Suitable for all Dentaureum spot welding units.
For soft-annealing, tempering and welding.

REF 081-501-00

2 pieces



Clamping electrode for hand electrode

Suitable for all Dentaureum spot welding units.
For soft-annealing, tempering and welding.

REF 081-701-00

1 piece

Welding electrode

Suitable for all Dentaaurum spot welding units.
Material: copper

Welding electrode	Version	REF	Quantity
upper	round, fine	085-000-00	1 piece
upper	wide, flat	085-100-00	1 piece
upper	round, thick	085-200-00	1 piece
lower	round, fine	086-000-00	1 piece
lower	wide, notched	086-100-00	1 piece
lower	round, thick	086-200-00	1 piece
lower	wide, flat	086-300-00	1 piece



Soldering electrode

Suitable for all Dentaaurum spot welding units.
Carbon with holder.

REF 085-300-00	1 piece
----------------	---------



Carbon insert

For soldering electrode (REF 085-300-00).

REF 085-400-00	10 pieces
----------------	-----------



Equipment

Equipment and accessories



Screw

For tightening electrodes on Dentaaurum spot welding units.

REF 084-100-00

10 pieces



File

For cleaning copper and carbon electrodes on spot welding units.

REF 083-300-00

10 pieces



Heat protection

Protects metal, plastic and ceramic parts against heat damage during soldering.

REF 380-070-00

100 g



Dentaflux® flux

- universal flux
- for all alloys
- especially for connections between precious and non-precious alloys

REF 681-100-00

50 g

Universal silver solder

- in rods
- ready fluxed
- cadmium free

Not permitted for intra-oral use.

Ø 1.0 mm - 1.1 mm

Working temperature 700 °C / 1292 °F

REF 380-604-50	1.2 g
REF 380-604-60	10 pieces 1.2 g each

Composition (percentage by mass)

Ag	Zn	Cu
59.0	24.0	16.0

**Universal silver solder**

- in coils
- without flux
- cadmium free

Not permitted for intra-oral use.

Ø 0.6 mm

Working temperature 700 °C / 1292 °F

REF 380-704-50	10 g
----------------	------

Composition (percentage by mass)

Ag	Zn	Cu
59.0	25.0	16.0

**White gold solder**

- in rods
- ready fluxed
- purity: ct 18.6
- cadmium free

Not permitted for intra-oral use.

Ø 1.0 mm - 1.1 mm

Working temperature 950 °C / 1742 °F

REF 380-600-50	1.6 g
----------------	-------

Composition (percentage by mass)

Au	Ni	Zn
77.0	16.0	6.0



Equipment and accessories



activeblue cleaning unit

With activefluid cleaning fluid.

For the gentle removal of plaque from orthodontic removable appliances.

Included:

- 1 x cleaning fluid, 1000 ml (REF 150-020-00)
- 1 x cleaning cup, white
- 1 x cleaning pins, 75 g (REF 098-001-00)
- 1 x tweezers
- 1 x holding magnet
- 1 x power cable (230 V)

Technical data

Power voltage	100 V~ - 240 V~ at 50 Hz / 60 Hz
Turning rotations	1200 min ⁻¹ / 1600 min ⁻¹ / 2000 min ⁻¹
Electric fuse	T 1.6 A L 250 V
Noise level	< 70 dB (A)
Weight	3 kg
Dimensions (W x H x D)	150 mm x 240 mm x 280 mm

REF 098-000-00

1 piece



Cleaning pins for activeblue

For activeblue cleaning unit.

REF 098-001-00

1 container 75 g each



activefluid cleaning fluid

For activeblue cleaning unit.

REF 150-020-00

1000 ml

Polyclav®

Pressure vessel for polymerization of cold-curing polymers, induction-compatible.
With safety valve.

Without pressure hose.

Capacity: 4.5 l

REF 070-010-00

1 piece

**Pressure hose**

For Polyclav® (REF 070-010-00).

Connection plug nipple to fit inlet valve and standard coupling for compressed air sockets.

Length: 1.5 m

Pmax: 8.0 bar

REF 070-015-00

1 piece

**Heating plate**

For Polyclav® (REF 070-000-00 and 070-010-00).

Power voltage	REF	Quantity
230 V~	070-900-00	1 piece

Technical data

Output	125 W
Thermostat for water temperature	43 °C / 109.4 °F
Power frequency	50 Hz / 60 Hz
Weight	2.6 kg
Dimensions (W x H x D)	240 mm x 240 mm x 75 mm

**Suction apparatus with charcoal filter**

For protection against vapor during processing of acrylics in the laboratory.

- controlled by a photoelectric sensor (automatic switch on/off)
- includes hour meter
- easy-change filter (replace the filter cartridge, REF 072-001-40, after approx. 50 operating hours)

Power voltage	REF	Quantity
230 V~	072-001-00	1 piece
115 V~	072-002-00	1 piece

Technical data

Output	100 W
Power frequency	50 Hz
Weight	6 kg
Dimensions (W x H x D)	220 mm x 260 mm x 260 mm



Equipment and accessories



Dual wheel model trimmer

For plaster models.

Coarse and fine corundum grinding wheels. Can be used on either side. Separately controlled motor switches with 450 W capacity. Individually controlled water pressure by a magnetic water valve. Work area with squaring and angulation devices.

Power voltage	REF	Quantity
230 V~	091-100-00	1 piece

Technical data

Output	2 x 450 W
Turning rotations	1400 min ⁻¹
Power frequency	50 Hz / 60 Hz
Weight	39 kg
Dimensions (W x H x D)	530 mm x 400 mm x 385 mm

Corundum disc

For Dentaureum dual wheel model trimmer (REF 091-100-00).

Grain size	REF	Quantity
46 = coarse	091-101-00	1 piece
120 = fine	091-102-00	1 piece

**Diamond disc, coarse**

For Dentaureum dual wheel model trimmer (REF 091-100-00).

- long operating life
- extremely high abrasive capacity
- minimal force required

REF 091-105-00	1 piece
-----------------------	---------

**Mini-Torch**

Ignition and flame adjustment is easy and uncomplicated (automatic, with a flick of your finger). Refill as easy as filling a lighter. Adjustable nozzle.

Technical data

Max. temperature	1300 °C / 2372 °F
Filling quantity	65 ml
Weight	280 g

REF 097-400-01	1 piece
-----------------------	---------



Equipment and accessories

	Modeling tip for the Thermomat		
	Designation	REF	Quantity
	Wax knife, small	089-602-00	1 piece
	Wax knife, large	089-605-00	1 piece
	Modeling tip, fine	089-603-00	1 piece
	Modeling tip, medium	089-604-00	1 piece
	Modeling spoon	089-606-00	1 piece



Handpiece with cable
For Thermomat (REF 089-620-00).

REF 089-612-00	1 piece
----------------	---------

o-atlas II	Page	428
Media library	Page	430
Download center	Page	432
Social media	Page	434
Patient consultation	Page	436
Build-a-Brace	Page	438



Knowledge I Marketing.

The key to success.



DENTAURUM
**ONLINE
SHOP**
shop.dentaurum.com





The reference guide for
removable orthodontic appliances.

Atlas of removable orthodontic appliances.

In a fully revised, extended edition with new images, Ursula Wirtz presents all of the important removable orthodontic appliances.

It is a standard reference work designed as a work and study aid for removable orthodontic appliances and is widely illustrated. It serves as a comprehensive reference guide for both orthodontists and dental technicians.

The o-atlas II gives clear, precise information on how to make an orthodontic appliance. The new edition is characterized by its pleasing, modern design.

- Approx. 390 pages with more than 1000 illustrations as hardback.
- 235 different orthodontic appliances are shown and described. The o-atlas II includes 44 new orthodontic appliances.
- More than 800 images have been replaced with new images of high quality.
- A 360 degree product view of the 50 most interesting orthodontic appliances is available online at www.o-atlas.com.
- Intelligent connection between print and online media in the o-atlas II.

Order now – ordering information:

Book – German	REF 903-920-10
Book – English	REF 903-920-20
Book – French	REF 903-920-30
Book – Spanish	REF 903-920-40
Book – Polish	on request



Online media library.


Videos and tutorials.

www.dentaurum.com/eng/media-library

See how it works.

You can find video clips and animations in best quality explaining step-by-step how to use our products at www.dentaurum.com.





ALL

ORTHODONTICS

PROSTHETICS

CERAMICS

IMPLANTOLOGY

COURSES & EVENTS

COMPANY

ALL ORTHODONTICS



Digital orthodontics



Dentaurum expansion screws



Dentaurum Brackets



Download center.

For easy and convenient access to Instructions for use and further documents.

www.dentaurum.com/eng/downloads

Downloads free of charge.

We want to make access to important documents on our products as easy as possible for you.

At www.dentaurum.com you will find a large selection of free downloads in our download center. Its structure is easy to follow.

- Instructions for use.
- Safety data sheets.
- Product information.
- Catalogs.
- Specialist publications.
- Patient information.
- Courses.
- Certificates.





Social media.

Always up-to-date.



www.facebook.com/DentaurumGermany



www.youtube.com/Dentaurum



www.instagram.com/dentaurum_germany

You will find current information, offers and insights into our company on our social networks.
Connect with us!



www.pinterest.com/dentaurumgmbhco



www.xing.com/companies/dentaurumgmbh&co.kg



www.kununu.com/dentaurum-co

Patient consultation.

Successful treatment with the right consultation.

Patient consultation the easy way.

For a course of treatment to be successful, patients need the right consultation and professional demonstration material.

We offer a wide range of information material to support you in your daily consultations with patients. Posters, patient information brochures, demonstration models, smartphone apps – we have them all.

These all serve to support you during daily patient consultations.





Build-a-Brace.

Modern, digital patient consultation.



Patient consultation done differently.

Dentaurum's app: Build-a-Brace is an excellent and innovative tool for the modern orthodontic practice. Let your patients design their own brace using their smartphone or tablet and then you can pass the design directly on to the dental laboratory.

Free-of-charge for smartphones and tablets – for iOS and Android.

The app is available in German and English and can be downloaded free of charge in the App Store or from Google Play.





Courses	Page	442
Webinars	Page	444



Courses, webinars and events.



Knowledge connects.

Digital. Analog. International.



Dentaurum provides the setting for courses and events in different languages, both in person and digitally. Broaden your horizon with us and learn about the latest worldwide developments and findings in orthodontics.

Overview.



More information on our courses is available at:



Dentaurum
www.dentaurum.de/kurse



Dentaurum Benelux
www.dentaurum.be



Dentaurum France
www.dentaurum.fr



Dentaurum España
www.dentaurum.es



Dentaurum Italia
<https://corsi.dentaurum.it>



Dentaurum USA
www.dentaurum.us



Dentaurum Canada
www.dentaurum.ca



Dentaurum Australia
www.dentaurum.com.au

Webinars at Dentaureum.



Informative, innovative and international – the webinar program from Dentaureum. National and international experts in the fields of orthodontics, implantology, dental prosthetics and ceramics give talks on various topics. The webinars are available live or on demand and are the ideal complement to Dentaureum's extensive course program.

The webinars can be viewed at any time worldwide.



Register now at:
www.dentaureum.com/webinare



Contact	Page	11
Bracket systems – information	Page	448
Pre-welding service – banding technique	Page	450
Bands – comparison charts	Page	462
Wires – information	Page	470
Materials list	Page	472
Conversion charts	Page	480
Environment and sustainability	Page	482
Alphabetical index	Page	484
Numerical index	Page	494



Dentaurum Service.

Added value – advice and support.



Bracket systems – information

Chart of the most frequently used “fixed appliance techniques”

Maxillary	1 1		2 2		3 3		4 4		5 5		6 6		7 7			
System	Torque	Angulation	Torque	Angulation	Torque	Angulation	Torque	Angulation	Torque	Angulation	Torque	Angulation	Distal offset = anti rotation	Torque	Angulation	Distal offset = anti rotation
Alexander	+ 14°	+ 5°	+ 7°	+ 8°	− 3°	+ 10°	− 7°	0°	− 7°	0°	− 10°	0°	15°	− 10°	0°	12°
Andrews	+ 7°	+ 5°	+ 3°	+ 9°	− 7°	+ 11°	− 7°	+ 2°	− 7°	+ 2°	− 15°	+ 1°	16°	− 15°	+ 1°	16°
McLaughlin-Bennett-Trevisi*	+ 17°	+ 4°	+ 10°	+ 8°	− 7°/0°/+ 7°	+ 8°	− 7°	0°	− 7°	0°	− 14°	0°	10°	− 14°	0°	10°
Burstone	+ 7°	+ 5°	+ 3°	+ 8°	− 7°	+ 10°	− 7°	0°	− 7°	0°	− 10°	− 5°	10°	− 10°	0°	6°
Hasund	+ 20°	+ 3°	+ 14°	+ 9°	− 2°	+ 6°	− 10°	+ 2°	− 10°	+ 2°	− 20°	+ 3°	10°	− 25°	+ 5°	6°
Hilgers	+ 22°	+ 5°	+ 14°	+ 8°	+ 7°	+ 10°	− 7°	0°	− 7°	0°	− 10°	0°	14°	− 10°	0°	10°
Ricketts® – IV. Dimension	+ 22°	0°	+ 14°	+ 8°	+ 7°	+ 5°	0°	0°	0°	0°	0°	0°	15°	− 10°	0°	12°
Ricketts® Universal	+ 22°	0°	+ 14°	+ 8°	+ 7°	+ 5°	0°	0°	0°	0°	0°	0°	0°			
Roth	+ 12°	+ 5°	+ 8°	+ 9°	− 2°	+ 9°/+ 13°	− 7°	0°	− 7°	0°	− 14°	0°	14°	− 14°	0°	14°
Standard Edgewise	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
Tweed	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°/6°	0°	0°	0°/6°
Begg	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°	6°	0°	0°	6°

Mandibular	1		2		3		4		5		6		7			
System	Torque	Angulation	Torque	Angulation	Torque	Angulation	Torque	Angulation	Torque	Angulation	Torque	Angulation	Distal offset = anti rotation	Torque	Angulation	Distal offset = anti rotation
Alexander	− 5°	0°	− 5°	0°	− 7°	+ 6°	− 11°	0°	− 17°	0°	− 22°	− 6°	5°	− 27°	0°	5°
Andrews	− 6°	+ 2°	− 6°	+ 2°	− 11°	+ 7°	− 17°	− 2°	− 22°	− 2°	− 30°	− 2°	6°	− 31°	− 2°	6°
McLaughlin-Bennett-Trevisi*	− 6°	0°	− 6°	0°	− 6°/0°/+ 6°	+ 3°	− 12°	+ 2°	− 17°	+ 2°	− 20°	0°	0°	− 10°	0°	0°
Burstone	− 1°	0°	− 1°	0°	− 11°	+ 6°	− 17°	0°	− 22°	0°	− 27°	0°	5°	− 27°	+ 2°	6°
Hasund	0°	0°	0°	+ 5°	0°	+ 5°	− 10°	+ 2°	− 15°	+ 2°	− 22°	+ 4°	0°	− 25°	+ 2°	6°
Hilgers	− 1°	0°	− 1°	0°	+ 7°	+ 6°	− 11°	0°	− 17°	0°	− 25°	− 5°	7°	− 25°	0°	6°
Ricketts® – IV. Dimension	0°	0°	0°	0°	+ 7°	+ 5°	− 7° ex 0° n-ex	0°	− 7° ex − 14° n-ex	0°	− 22°	− 5°	12°	− 27°	0°	16°
Ricketts® Universal	0°	0°	0°	0°	+ 7°	+ 5°	0°	0°	− 14°	0°	0°	− 5°	0°	0°	0°	0°
Roth	− 1°	0°	− 1°	0°	− 11°	+ 7°	− 17°	0°	− 22°	0°	− 30°	− 1°	4°	− 30°	− 1°	4°
Standard Edgewise	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
Tweed	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°/6°	0°	0°	0°/6°
Begg	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°	6°	0°	0°	6°

* The Dentaurem version of this prescription is not claimed to be a duplication of any other, nor does Dentaurem imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.

Note:

No responsibility is taken for accuracy of this information, systems are subject to change!

Pre-welding service.

For your treatment system – quick and reliable.

If required, all Dentaurem brackets, buccal tubes and palatal auxiliaries can be welded to Dentaurem bands by our welding service. Customer friendly and rapid processing of all welding orders by our highly efficient staff. Computer systems allow recording and storage of individual data and rapid delivery of the correct products.

How to process your welding orders.

- Contact us directly or via your local sales representative when ordering.
- For first orders, we draw up an individual order profile to your specifications – your system card.
- With the help of the order form Welding Service you order the items required and the quantities needed for each size. Any changes are recorded on your personal system card.


Standard and custom pre-welding service.

We provide the choice of a standard or custom pre-welding service. Brackets and buccal tubes are welded onto the band according to conventional systems.


Standard pre-welding service.

The Edgewise technique, e.g. Standard Edgewise, Roth, McLaughlin-Bennett-Trevisi* etc. are available. The system Roth 18 and 22 can be delivered from stock.

Edgewise.




Lingual/palatal attachments.



Custom pre-welding service.

Brackets and buccal tubes are welded onto the bands according to your requirements, e.g. with angulation. Please provide the precise angulation and position of the lingual/palatal attachment with your order.



Angulation / Tip

(+) = angulation of root to distal

(-) = angulation of root to mesial

Pre-welding service – banding technique

Roth 18

Slot: 0.46 mm x 0.76 mm / 18 x 30 · cuspid, 1st bicuspid and 2nd bicuspid bands are welded conventionally with equilibrium® 2 brackets

Tooth	Position on band / buccal view Position auf dem Band / bukkal		Band type / size	Brackets Buccal tubes	Torque	Angulation	In / Out mm	Distal Offset	Palatal attachment Position on band	Bondable tubes / Bukkalkröhrchen für die Adhäsivtechnik
Zahn	Bandtyp / Größe		Brackets Bukkalkröhrchen					Anti Rotation	Palatinales Teil Position auf d. Band	
<u>3 3</u> 13 - 23	R distal	gingival mesial incisal L distal	Standard Eckzahn Cuspid No. 1–25	Rechts / right REF 790-259-02 (S, H) REF 790-254-02 (S, H) Links / left REF 790-260-02 (S, H) REF 790-255-02 (S, H)	– 2°	+ 13°	0.7		REF 750-602-10 central	
<u>4 4</u> 14 - 24	R distal	gingival mesial occlusal L distal	Standard Prämolaren Bicuspid No. 1–32	Rechts / right REF 790-206-02 (S, H) REF 790-249-02 (S, H) Links / left REF 790-206-02 (S, H) REF 790-250-02 (S, H)	– 7°	0°	0.7		REF 750-602-10 central	
<u>5 5</u> 15 - 25	R distal	gingival mesial occlusal L distal	Standard Prämolaren Bicuspid No. 1–32	Rechts / right REF 790-206-02 (S, H) REF 790-249-02 (S, H) Links / left REF 790-206-02 (S, H) REF 790-250-02 (S, H)	– 7°	0°	0.7		REF 750-602-10 central	
<u>6 6</u> 16 - 26	R distal	gingival mesial occlusal L distal	dentaform® dentaform® Snap 1° Molar Bands No. 1–36 Standard Molar Bands No. 1–32	Rechts / right REF 818-131-10 Links / left REF 818-131-20	occlusal – 14°	0°	0.3	14°	REF 750-602-10 central	Rechts / right REF 818-211-16 Links / left REF 818-211-26
<u>7 7</u> 17 - 27	R distal	gingival mesial occlusal L distal	dentaform® 2° Molar Bands No. 1–32 Standard Molar Bands No. 1–32	Rechts / right REF 818-212-10 Links / left REF 818-212-20	– 14°	0°	0.3	10°	REF 750-602-10 central	Rechts / right REF 818-411-17 Links / left REF 818-411-27



Roth 18

Slot: 0.46 mm x 0.76 mm / 18 x 30 · cuspid, 1st bicuspid and 2nd bicuspid bands are welded conventionally with equilibrium® 2 brackets

Tooth	Position on band / buccal view Position auf dem Band / bukkal		Band type / size	Brackets Buccal tubes	Torque	Angulation	In / Out mm	Distal Offset	Lingual attachment Position on band	Bondable tubes / Bukkalkröhrchen für die Adhäsivtechnik
Zahn	Bandtyp / Größe		Brackets Bukkalkröhrchen					Anti Rotation	Linguales Teil Position auf d. Band	
<u>3 3</u> 43 - 33	R distal	incisal mesial L distal	Standard Eckzahn Cuspid No. 1–21	Rechts / right REF 790-261-02 (S, H) REF 790-208-02 (S, H) Links / left REF 790-262-02 (S, H) REF 790-209-02 (S, H)	– 11°	+ 7°	0.7		REF 750-602-10 central	
<u>4 4</u> 44 - 34	R distal	gingival occlusal mesial L distal	Standard Prämolaren Bicuspid No. 1–32	Rechts / right REF 790-210-02 (S, H) REF 790-251-02 (S, H) Links / left REF 790-210-02 (S, H) REF 790-252-02 (S, H)	– 17°	0°	0.7		REF 750-602-10 central	
<u>5 5</u> 45 - 35	R distal	gingival occlusal mesial L distal	Standard Prämolaren Bicuspid No. 1–32	Rechts / right REF 790-211-02 (S, H) REF 790-253-02 (S, H) Links / left REF 790-211-02 (S, H) REF 790-254-02 (S, H)	– 22°	0°	0.7		REF 750-602-10 central	
<u>6 6</u> 46 - 36	R distal	gingival occlusal mesial L distal	dentaform® dentaform® Snap 1° Molar Bands No. 1–36 Standard Molar Bands No. 1–32	Rechts / right REF 818-151-40 Links / left REF 818-151-30	gingival – 30°	0°	0.3	4°	REF 750-602-10 central	Rechts / right REF 818-211-46 Links / left REF 818-211-36
<u>7 7</u> 47 - 37	R distal	gingival occlusal mesial L distal	dentaform® 2° Molar Bands No. 1–32 Standard Molar Bands No. 1–32	Rechts / right REF 818-218-40 Links / left REF 818-218-30	– 30°	0°	0.3	12°	REF 750-602-10 central	Rechts / right REF 818-411-47 Links / left REF 818-411-37



All brackets, buccal tubes, lingual and palatal auxiliaries can be welded onto bands in our welding department.

The illustrations shown are examples of bands with pre-welded attachments and are not identical with the described system.

Pre-welding service – banding technique

Roth 22

Slot: 0.56 mm x 0.76 mm / 22 x 30 · cuspid, 1st bicuspid and 2nd bicuspid bands are welded conventionally with equilibrium® 2 brackets

Tooth	Position on band/buccal view Position auf dem Band/buccal	Band type/size	Brackets Buccal tubes	Torque	Angulation	In/Out mm	Distal Offset	Palatal attachment Position on band	Bondable tubes/ Bukkalaröhren für die Adhäsivtechnik
Zahn		Bandtyp/Größe	Brackets Bukkalaröhren				Anti Rotation	Palatinales Teil Position auf d. Band	
<u>3</u> <u>3</u> 13 · 23	R distal gingival mesial L distal incisal	Standard Eckzahn Cuspid No. 1–25	Rechts / right REF 790-263-02 (S, R) REF 790-216-02 (S, R) Links / left REF 790-264-02 (S, R) REF 790-217-02 (S, R)	– 2°	+ 13°	0.7		REF 750-602-10 central	
<u>4</u> <u>4</u> 14 · 24	R distal gingival mesial L distal occlusal	Standard Prämolaren Bicuspid No. 1–32	Rechts / right REF 790-218-02 (S, R) REF 790-243-02 (S, R) Links / left REF 790-218-02 (S, R) REF 790-244-02 (S, R)	– 7°	0°	0.7		REF 750-602-10 central	
<u>5</u> <u>5</u> 15 · 25	R distal gingival mesial L distal occlusal	Standard Prämolaren Bicuspid No. 1–32	Rechts / right REF 790-218-02 (S, R) REF 790-243-02 (S, R) Links / left REF 790-218-02 (S, R) REF 790-244-02 (S, R)	– 7°	0°	0.7		REF 750-602-10 central	
<u>6</u> <u>6</u> 16 · 26	R distal gingival mesial L distal occlusal	dentaform® dentaform® Snap 1° Molar Bands No. 1–36 Standard Molar Bands No. 1–32	Rechts / right REF 822-131-10 Links / left REF 822-131-20	occlusal – 14°	0°	0.3	14°	REF 750-602-10 central	Rechts / right REF 822-211-16 Links / left REF 822-211-26
<u>7</u> <u>7</u> 17 · 27	R distal gingival mesial L distal occlusal	dentaform® 2° Molar Bands No. 1–32 Standard Molar Bands No. 1–32	Rechts / right REF 822-212-10 Links / left REF 822-212-20	– 14°	0°	0.3	10°	REF 750-602-10 central	Rechts / right REF 822-411-17 Links / left REF 822-411-27



Roth 22

Slot: 0.56 mm x 0.76 mm / 22 x 30 · cuspid, 1st bicuspid and 2nd bicuspid bands are welded conventionally with equilibrium® 2 brackets

Tooth	Position on band/buccal view Position auf dem Band/buccal	Band type/size	Brackets Buccal tubes	Torque	Angulation	In/Out mm	Distal Offset	Lingual attachment Position on band	Bondable tubes/ Bukkalaröhren für die Adhäsivtechnik
Zahn		Bandtyp/Größe	Brackets Bukkalaröhren				Anti Rotation	Linguales Teil Position auf d. Band	
<u>3</u> <u>3</u> 43 · 33	R distal incisal mesial L distal	Standard Eckzahn Cuspid No. 1–21	Rechts / right REF 790-265-02 (S, R) REF 790-220-02 (S, R) Links / left REF 790-266-02 (S, R) REF 790-221-02 (S, R)	– 11°	+ 7°	0.7		REF 750-602-10 central	
<u>4</u> <u>4</u> 44 · 34	R distal gingival occlusal L distal	Standard Prämolaren Bicuspid No. 1–32	Rechts / right REF 790-222-02 (S, R) REF 790-255-02 (S, R) Links / left REF 790-222-02 (S, R) REF 790-256-02 (S, R)	– 17°	0°	0.7		REF 750-602-10 central	
<u>5</u> <u>5</u> 45 · 35	R distal gingival occlusal L distal	Standard Prämolaren Bicuspid No. 1–32	Rechts / right REF 790-223-02 (S, R) REF 790-257-02 (S, R) Links / left REF 790-223-02 (S, R) REF 790-258-02 (S, R)	– 22°	0°	0.7		REF 750-602-10 central	
<u>6</u> <u>6</u> 46 · 36	R distal gingival occlusal L distal	dentaform® dentaform® Snap 1° Molar Bands No. 1–36 Standard Molar Bands No. 1–32	Rechts / right REF 822-151-40 Links / left REF 822-151-30	gingival – 30°	0°	0.3	4°	REF 750-602-10 central	Rechts / right REF 822-211-46 Links / left REF 822-211-36
<u>7</u> <u>7</u> 47 · 37	R distal gingival occlusal L distal	dentaform® 2° Molar Bands No. 1–32 Standard Molar Bands No. 1–32	Rechts / right REF 822-218-40 Links / left REF 822-218-30	– 30°	0°	0.3	12°	REF 750-602-10 central	Rechts / right REF 822-411-47 Links / left REF 822-411-37



All brackets, buccal tubes, lingual and palatal auxiliaries can be welded onto bands in our welding department.

The illustrations shown are examples of bands with pre-welded attachments and are not identical with the described system.

Pre-welding service – banding technique

McLaughlin-Bennett-Trevisi* 18

Slot: 0.46 mm x 0.76 mm / 18 x 30 · cuspid, 1st bicuspid and 2nd bicuspid bands are welded conventionally with equilibrium® 2 brackets

Tooth	Position on band / buccal view Position auf dem Band / bukkal		Band type / size	Brackets Buccal tubes	Torque	Angulation	In / Out mm	Distal Offset	Palatal attachment Position on band	Bondable tubes / Bukkalt�hren f�r die Adh�sivtechnik
Zahn			Bandtyp / Gr��e	Brackets Bukkalt�hren				Anti Rotation	Palatinales Teil Position auf d. Band	
<u>3 3</u> 13 - 23	R	gingival mesial incisal	L distal	Standard Eckzahn Cuspid No. 1–25	Rechts / right REF 790-274-02 (25,34) Links / left REF 790-275-02 (25,34)	0°	+ 8°	0.7	REF 750-602-10 central	
<u>4 4</u> 14 - 24	R	gingival mesial occlusal	L distal	Standard Pr�molaren Bicuspid No. 1–32	Rechts / right REF 790-206-02 (32,34) REF 790-249-02 (32,34) Links / left REF 790-206-02 (32,34) REF 790-250-02 (32,34)	– 7°	0°	0.7	REF 750-602-10 central	
<u>5 5</u> 15 - 25	R	gingival mesial occlusal	L distal	Standard Pr�molaren Bicuspid No. 1–32	Rechts / right REF 790-206-02 (32,34) REF 790-249-02 (32,34) Links / left REF 790-206-02 (32,34) REF 790-250-02 (32,34)	– 7°	0°	0.7	REF 750-602-10 central	
<u>6 6</u> 16 - 26	R	gingival mesial occlusal	L distal	dentaform® dentaform® Snap 1° Molar Bands No. 1–36 Standard Molar Bands No. 1–32	Rechts / right REF 818-132-10 Links / left REF 818-132-20	occlusal –14°	0°	0.3	REF 750-602-10 central	Rechts / right REF 818-212-16 Links / left REF 818-212-26
<u>7 7</u> 17 - 27	R	gingival mesial occlusal	L distal	dentaform® 2° Molar Bands No. 1–32 Standard Molar Bands No. 1–32	Rechts / right REF 818-212-10 Links / left REF 818-212-20	–14°	0°	0.3	REF 750-602-10 central	Rechts / right REF 818-412-17 Links / left REF 818-412-27



McLaughlin-Bennett-Trevisi* 18

Slot: 0.46 mm x 0.76 mm / 18 x 30 · cuspid, 1st bicuspid and 2nd bicuspid bands are welded conventionally with equilibrium® 2 brackets

Tooth	Position on band/buccal view Position auf dem Band/bukkal		Band type/size	Brackets Buccal tubes	Torque	Angulation	In/Out mm	Distal Offset	Lingual attachment Position on band	Bondable tubes/ Buckalr�hren f�r die Adh�sivtechnik
Zahn			Bandtyp/Gr��e	Brackets Buckalr�hren				Anti Rotation	Linguales Teil Position auf d. Band	
3 3 43 - 33	R	incisal mesial	L distal	Standard Eckzahn Cuspid No. 1–21	Rechts / right REF 790-278-02 (21,34) Links / left REF 790-279-02 (21,34)	0°	+ 3°	0.7	REF 750-602-10 central	
4 4 44 - 34	R	occlusal mesial	L distal	Standard Pr�molaren Bicuspid No. 1–32	Rechts / right REF 790-210-02 (32,34) REF 790-251-02 (32,34) Links / left REF 790-210-02 (32,34) REF 790-252-02 (32,34)	– 17°	0°	0.7	REF 750-602-10 central	
5 5 45 - 35	R	occlusal mesial	L distal	Standard Pr�molaren Bicuspid No. 1–32	Rechts / right REF 790-210-02 (32,34) REF 790-251-02 (32,34) Links / left REF 790-210-02 (32,34) REF 790-252-02 (32,34)	– 17°	0°	0.7	REF 750-602-10 central	
6 6 46 - 36	R	occlusal mesial	L distal	dentaform® dentaform® Snap 1° Molar Bands No. 1–36 Standard Molar Bands No. 1–32	Rechts / right REF 818-152-40 Links / left REF 818-152-30	gingival –20°	0°	0.3	0° REF 750-602-10 central	Rechts / right REF 818-212-46 Links / left REF 818-212-36
7 7 47 - 37	R	occlusal mesial	L distal	dentaform® 2° Molar Bands No. 1–32 Standard Molar Bands No. 1–32	Rechts / right REF 818-512-40 Links / left REF 818-512-30	gingival –10°	0°	0.3	0° REF 750-602-10 central	Rechts / right REF 818-412-47 Links / left REF 818-412-37



All brackets, buccal tubes, lingual and palatal auxiliaries can be welded onto bands in our welding department.

The illustrations shown are examples of bands with pre-welded attachments and are not identical with the described system.

Pre-welding service – banding technique

McLaughlin-Bennett-Trevisi* 22

Slot: 0.56 mm x 0.76 mm / 22 x 30 · cuspid, 1st bicuspid and 2nd bicuspid bands are welded conventionally with equilibrium® 2 brackets

Tooth	Position on band/buccal view Position auf dem Band/buccal	Band type/size	Brackets Buccal tubes	Torque	Angulation	In/Out mm	Distal Offset	Palatal attachment Position on band	Bondable tubes/ Bukkalaröhrchen für die Adhäsivtechnik
Zahn		Bandtyp/Größe	Brackets Bukkalaröhrchen				Anti Rotation	Palatinales Teil Position auf d. Band	
<u>3</u> <u>3</u> 13 · 23	R distal gingival mesial L distal incisal	Standard Eckzahn Cuspid No. 1–25	Rechts / right REF 790-230-02 (2xR) REF 790-228-02 (2xL) Links / left REF 790-231-02 (2xR) REF 790-229-02 (2xL)	- 7°	0°	+ 8°	0.7	REF 750-602-10 central	 20 32 5445 32 5445
<u>4</u> <u>4</u> 14 · 24	R distal gingival mesial L distal occlusal	Standard Prämolaren Bicuspid No. 1–32	Rechts / right REF 790-218-02 (2xR) REF 790-243-02 (2xL) Links / left REF 790-218-02 (2xR) REF 790-244-02 (2xL)	- 7°	0°	0°	0.7	REF 750-602-10 central	 32 5445
<u>5</u> <u>5</u> 15 · 25	R distal gingival mesial L distal occlusal	Standard Prämolaren Bicuspid No. 1–32	Rechts / right REF 790-218-02 (2xR) REF 790-243-02 (2xL) Links / left REF 790-218-02 (2xR) REF 790-244-02 (2xL)	- 7°	0°	0°	0.7	REF 750-602-10 central	 32 5445
<u>6</u> <u>6</u> 16 · 26	R distal gingival mesial L distal occlusal	dentaform® dentaform® Snap 1° Molar Bands No. 1–36 Standard Molar Bands No. 1–32	Rechts / right REF 822-132-10 Links / left REF 822-132-20	- 14°	0°	0.3	10°	REF 750-602-10 central	Rechts / right REF 822-212-16 Links / left REF 822-212-26
<u>7</u> <u>7</u> 17 · 27	R distal gingival mesial L distal occlusal	dentaform® 2° Molar Bands No. 1–32 Standard Molar Bands No. 1–32	Rechts / right REF 822-212-10 Links / left REF 822-212-20	- 14°	0°	0.3	10°	REF 750-602-10 central	Rechts / right REF 822-412-17 Links / left REF 822-412-27

McLaughlin-Bennett-Trevisi* 22

Slot: 0.56 mm x 0.76 mm / 22 x 30 · cuspid, 1st bicuspid and 2nd bicuspid bands are welded conventionally with equilibrium® 2 brackets

Tooth	Position on band/buccal view Position auf dem Band/buccal	Band type/size	Brackets Buccal tubes	Torque	Angulation	In/Out mm	Distal Offset	Lingual attachment Position on band	Bondable tubes/ Bukkalaröhrchen für die Adhäsivtechnik
Zahn		Bandtyp/Größe	Brackets Bukkalaröhrchen				Anti Rotation	Linguales Teil Position auf d. Band	
<u>3</u> <u>3</u> 43 · 33	R distal incisal mesial L distal	Standard Eckzahn Cuspid No. 1–21	Rechts / right REF 790-236-02 (2xR) REF 790-240-02 (2xL) Links / left REF 790-237-02 (2xR) REF 790-235-02 (2xL)	- 6°	0°	+ 3°	0.7	REF 750-602-10 central	 16 32 5445
<u>4</u> <u>4</u> 44 · 34	R distal gingival occlusal L distal	Standard Prämolaren Bicuspid No. 1–32	Rechts / right REF 790-222-02 (2xR) REF 790-294-02 (2xL) Links / left REF 790-222-02 (2xR) REF 790-295-02 (2xL)	- 17°	- 11°	0°	0.7	REF 750-602-10 central	 32 5445
<u>5</u> <u>5</u> 45 · 35	R distal gingival occlusal L distal	Standard Prämolaren Bicuspid No. 1–32	Rechts / right REF 790-222-02 (2xR) REF 790-294-02 (2xL) Links / left REF 790-222-02 (2xR) REF 790-295-02 (2xL)	- 17°	- 11°	0°	0.7	REF 750-602-10 central	 32 5445
<u>6</u> <u>6</u> 46 · 36	R distal gingival occlusal L distal	dentaform® dentaform® Snap II 1° Molar Bands No. 1–36 Standard Molar Bands No. 1–32	Rechts / right REF 822-152-40 Links / left REF 822-152-30	- 20°	0°	0.3	0°	REF 750-602-10 central	Rechts / right REF 822-212-46 Links / left REF 822-212-36
<u>7</u> <u>7</u> 47 · 37	R distal gingival occlusal L distal	dentaform® 2° Molar Bands No. 1–32 Standard Molar Bands No. 1–32	Rechts / right REF 822-512-40 Links / left REF 822-512-30	- 10°	0°	0.3	0°	REF 750-602-10 central	Rechts / right REF 822-412-47 Links / left REF 822-412-37

All brackets, buccal tubes, lingual and palatal auxiliaries can be welded onto bands in our welding department.

The illustrations shown are examples of bands with pre-welded attachments and are not identical with the described system.

* The Dentaform version of this prescription is not claimed to be a duplication of any other, nor does Dentaform imply that it is endorsed in any way by Drs. McLaughlin, Bennett and Trevisi.

Pre-welding service – banding technique

Ricketts® 18 – IV. Dimension

Slot: 0.46 mm x 0.76 mm / 18 x 30

Tooth	Position on band/buccal view Position auf dem Band/bukkal	Band type/size Bandtyp/Größe	Brackets Bukkalkörhchen	Torque	Angulation	In/Out mm	Distal Offset	Palatal attachment Position on band Palatinales Teil Position auf d. Band	Bondable tubes / Bukkalkörhchen für die Adhäsivtechnik
Zahn							Anti Rotation		
$\frac{3}{13} \mid \frac{3}{23}$	R distal gingival mesial incisal L distal	Standard Eckzahn Cuspids No. 1–25	Rechts / right REF 790-290-02 Links / left REF 790-291-02	+ 7°	0°	0.7		REF 750-602-10 central	
$\frac{4}{14} \mid \frac{4}{24}$	R distal gingival mesial occlusal L distal	Standard Prämolaren Bicuspid No. 1–32	Rechts / right REF 790-602-02 Links / left REF 790-602-02	0°	0°	0.7		REF 750-602-10 central	
$\frac{5}{15} \mid \frac{5}{25}$	R distal gingival mesial occlusal L distal	Standard Prämolaren Bicuspid No. 1–32	Rechts / right REF 790-602-02 Links / left REF 790-602-02	0°	0°	0.7		REF 750-602-10 central	
$\frac{6}{16} \mid \frac{6}{26}$	R distal gingival mesial occlusal L distal	dentaform® dentaform® Snap 1° Molar Bands No. 1–36 Standard Molar Bands No. 1–32	Rechts / right REF 724-006-00 Links / left REF 724-007-00	occlusal 0°	0°	0.3	15°	REF 750-602-10 central	
$\frac{7}{17} \mid \frac{7}{27}$	R distal gingival mesial occlusal L distal	dentaform® 2° Molar Bands No. 1–32 Standard Molar Bands No. 1–32	Rechts / right REF 818-217-10 Links / left REF 818-217-20	–10°	0°	0.3	8°	REF 750-602-10 central	Rechts / right REF 724-048-51 Links / left REF 724-049-51



Ricketts® 18 – IV. Dimension

Slot: 0.46 mm x 0.76 mm / 18 x 30

Tooth	Position on band/buccal view Position auf dem Band/bukkal	Band type/size Bandtyp/Größe	Brackets Bukkalkörhchen	Torque	Angulation	In/Out mm	Distal Offset	Lingual attachment Position on band Linguales Teil Position auf d. Band	Bondable tubes / Bukkalkörhchen für die Adhäsivtechnik
Zahn							Anti Rotation		
$\frac{3}{43} \mid \frac{3}{33}$	R distal incisal mesial L distal	Standard Eckzahn Cuspids No. 1–21	Rechts / right REF 790-278-02 (20,34) Links / left REF 790-279-02 (20,34)	0°	+ 3°	0.7		REF 750-602-10 central	
$\frac{4}{44} \mid \frac{4}{34}$	R distal gingival occlusal mesial L distal	Standard Prämolaren Bicuspid No. 1–32	Rechts / right REF 790-602-02 Links / left REF 790-602-02	0°	0°	0.7		REF 750-602-10 central	
$\frac{5}{45} \mid \frac{5}{35}$	R distal gingival occlusal mesial L distal	Standard Prämolaren Bicuspid No. 1–32	Rechts / right REF 790-210-02 Links / left REF 790-210-02	–17°	0°	0.7		REF 750-602-10 central	
$\frac{6}{46} \mid \frac{6}{36}$	R distal gingival occlusal mesial L distal	dentaform® dentaform® Snap 1° Molar Bands No. 1–36 Standard Molar Bands No. 1–32	Rechts / right REF 724-064-00 Links / left REF 724-065-00	gingival –22°	–5°	0.3	12°	REF 750-602-10 central	Rechts / right REF 818-122-46 Links / left REF 818-122-36
$\frac{7}{47} \mid \frac{7}{37}$	R distal gingival occlusal mesial L distal	dentaform® 2° Molar Bands No. 1–32 Standard Molar Bands No. 1–32	Rechts / right REF 818-218-40 Links / left REF 818-218-30	–30°	–5°	0.3	12°	REF 750-602-10 central	Rechts / right REF 818-211-46 Links / left REF 818-211-36



All brackets, buccal tubes, lingual and palatal auxiliaries can be welded onto bands in our welding department.

The illustrations shown are examples of bands with pre-welded attachments and are not identical with the described system.

Pre-welding service – banding technique

Ricketts® Universal 18

Slot: 0.46 mm x 0.76 mm / 18 x 30

Tooth	Position on band/buccal view Position auf dem Band/bukkal	Band type / size Bandtyp / Größe	Brackets Buccal tubes	Torque	Angulation	In / Out mm	Distal Offset	Palatal attachment Position on band	
Zahn									
<u>3</u> <u>3</u> 13 - 23	R distal gingival mesial L distal incisal	Standard Eckzahn Cuspids No. 1–25	Rechts / right REF 790-290-02 Links / left REF 790-290-02	+ 7°	0°	0.7		REF 750-602-10 central	    
<u>4</u> <u>4</u> 14 - 24	R distal gingival mesial L distal occlusal	Standard Prämolaren Bicuspid No. 1–32	Rechts / right REF 790-602-02 Links / left REF 790-602-02	0°	0°	0.7		REF 750-602-10 central	
<u>5</u> <u>5</u> 15 - 25	R distal gingival mesial L distal occlusal	Standard Prämolaren Bicuspid No. 1–32	Rechts / right REF 790-602-02 Links / left REF 790-602-02	0°	0°	0.7		REF 750-602-10 central	
<u>6</u> <u>6</u> 16 - 26	R distal gingival mesial L distal occlusal	dentaform® Snap 1° Molar Bands No. 1–36 Standard Molar Bands No. 1–32	Rechts / right REF 724-600-00 Links / left REF 724-601-00	central 0°	0°	0.3	0°	REF 750-602-10 central	

Ricketts® Universal 18

Slot: 0.46 mm x 0.76 mm / 18 x 30

Tooth	Position on band/buccal view Position auf dem Band/bukkal	Band type / size Bandtyp / Größe	Brackets Buccal tubes	Torque	Angulation	In / Out mm	Distal Offset	Lingual attachment Position on band	Bondable tubes/ Bukkalehrchen für die Adhäsivtechnik
Zahn									
<u>3</u> <u>3</u> 43 - 33	R distal incisal mesial L distal	Standard Eckzahn Cuspids No. 1–21	Rechts / right REF 790-278-02 (09,10) Links / left REF 790-279-02 (09,10)	0°	+ 3°	0.7		REF 750-602-10 central	    
<u>4</u> <u>4</u> 44 - 34	R distal gingival mesial L distal occlusal	Standard Prämolaren Bicuspid No. 1–32	Rechts / right REF 790-602-02 Links / left REF 790-602-02	0°	0°	0.7		REF 750-602-10 central	
<u>5</u> <u>5</u> 45 - 35	R distal gingival mesial L distal occlusal	Standard Prämolaren Bicuspid No. 1–32	Rechts / right REF 790-210-02 Links / left REF 790-210-02	- 17°	0°	0.7		REF 750-602-10 central	
<u>6</u> <u>6</u> 46 - 36	R distal gingival mesial L distal occlusal	dentaform® Snap 1° Molar Bands No. 1–36 Standard Molar Bands No. 1–32	Rechts / right REF 724-018-00 Links / left REF 724-018-00	gingival 0°	- 5°	0.3	0°	REF 750-602-10 central	
								Rechts / right REF 724-018-51 Links / left REF 724-018-51	

All brackets, buccal tubes, lingual and palatal auxiliaries can be welded onto bands in our welding department.

The illustrations shown are examples of bands with pre-welded attachments and are not identical with the described system.

Pre-welding service – banding technique

Hilgers 18

Slot: 0.46 mm x 0.76 mm / 18 x 30

Tooth	Position on band/buccal view Position auf dem Band/bukkal	Band type/size	Brackets Buccal tubes	Torque	Angulation	In / Out mm	Distal Offset	Palatal attachment Position on band	Bondable tubes / Bukkalkörhörnchen für die Adhäsivtechnik
Zahn		Bandtyp / Größe	Brackets Bukkalkörhörnchen				Anti Rotation	Palatinales Teil Position auf d. Band	
<u>3 3</u> 13 - 23	R distal gingival mesial L distal incisal	Standard Eckzahn Cuspid No. 1–25	Rechts / right REF 790-290-02 Links / left REF 790-291-02	+ 7°	0°	0.7		REF 750-602-10 central	
<u>4 4</u> 14 - 24	R distal gingival mesial L distal occlusal	Standard Prämolaren Bicuspid No. 1–32	Rechts / right REF 790-206-02 Links / left REF 790-206-02	- 7°	0°	0.7		REF 750-602-10 central	
<u>5 5</u> 15 - 25	R distal gingival mesial L distal occlusal	Standard Prämolaren Bicuspid No. 1–32	Rechts / right REF 790-206-02 Links / left REF 790-206-02	- 7°	0°	0.7		REF 750-602-10 central	
<u>6 6</u> 16 - 26	R distal gingival mesial L distal occlusal	dentaform® dentaform® Snap 1° Molar Bands No. 1–36 Standard Molar Bands No. 1–32	Rechts / right REF 818-131-10 Links / left REF 818-131-20	occlusal -14°	0°	0.3	14°	REF 750-602-10 central	Rechts / right REF 724-266-51 Links / left REF 724-267-51
<u>7 7</u> 17 - 27	R distal gingival mesial L distal occlusal	dentaform® 2° Molar Bands No. 1–32 Standard Molar Bands No. 1–32	Rechts / right REF 818-217-10 Links / left REF 818-217-20	-10°	0°	0.3	8°	REF 750-602-10 central	Rechts / right REF 724-048-51 Links / left REF 724-049-51



Hilgers 18

Slot: 0.46 mm x 0.76 mm / 18 x 30

Tooth	Position on band/buccal view Position auf dem Band/bukkal	Band type/size	Brackets Buccal tubes	Torque	Angulation	In / Out mm	Distal Offset	Lingual attachment Position on band	Bondable tubes / Bukkalkörhörnchen für die Adhäsivtechnik
Zahn		Bandtyp / Größe	Brackets Bukkalkörhörnchen				Anti Rotation	Linguales Teil Position auf d. Band	
<u>3 3</u> 43 - 33	R distal incisal mesial L distal	Standard Eckzahn Cuspid No. 1–21	Rechts / right REF 790-278-02 Links / left REF 790-279-02	0°	+ 3°	0.7		REF 750-602-10 central	
<u>4 4</u> 44 - 34	R distal gingival occlusal L distal gingival	Standard Prämolaren Bicuspid No. 1–32	Rechts / right REF 790-210-02 Links / left REF 790-210-02	- 17°	0°	0.7		REF 750-602-10 central	
<u>5 5</u> 45 - 35	R distal gingival occlusal L distal gingival	Standard Prämolaren Bicuspid No. 1–32	Rechts / right REF 790-210-02 Links / left REF 790-210-02	- 17°	0°	0.7		REF 750-602-10 central	
<u>6 6</u> 46 - 36	R distal gingival occlusal L distal gingival	dentaform® dentaform® Snap 1° Molar Bands No. 1–36 Standard Molar Bands No. 1–32	Rechts / right REF 724-056-00 Links / left REF 724-057-00	gingival -25°	- 5°	0.3	8°	REF 750-602-10 central	Rechts / right REF 724-056-51 Links / left REF 724-057-51
<u>7 7</u> 47 - 37	R distal gingival occlusal L distal gingival	dentaform® 2° Molar Bands No. 1–32 Standard Molar Bands No. 1–32	Rechts / right REF 818-218-40 Links / left REF 818-218-30	-30°	- 5°	0.3	12°	REF 750-602-10 central	Rechts / right REF 818-211-46 Links / left REF 818-211-36



All brackets, buccal tubes, lingual and palatal auxiliaries can be welded onto bands in our welding department.

The illustrations shown are examples of bands with pre-welded attachments and are not identical with the described system.

Pre-welding service – banding technique

Standard Edgewise 18

Slot: 0.46 mm x 0.76 mm / 18 x 30 · cuspid, 1st bicuspid and 2nd bicuspid bands are welded conventionally with equilibrium® 2 brackets

Tooth	Position on band/buccal view Position auf dem Band/buccal	Band type / size	Brackets Buccal tubes	Torque	Angulation	In / Out mm	Distal Offset	Palatal attachment Position on band
Zahn		Bandtyp / Größe	Brackets Bukkaloehrrchen				Anti Rotation	Palatinales Teil Position auf d. Band
<u>3</u> <u>3</u> 13 - 23	R distal gingival mesial L distal incisal	Standard Eckzahn Cuspid No. 1–25	Rechts / right REF 790-602-02 Links / left REF 790-602-02	0°	0°	0.7		REF 750-602-10 central
<u>4</u> <u>4</u> 14 - 24	R distal gingival mesial L distal occlusal	Standard Prämolaren Bicuspid No. 1–32	Rechts / right REF 790-602-02 Links / left REF 790-602-02	0°	0°	0.7		REF 750-602-10 central
<u>5</u> <u>5</u> 15 - 25	R distal gingival mesial L distal occlusal	Standard Prämolaren Bicuspid No. 1–32	Rechts / right REF 790-602-02 Links / left REF 790-602-02	0°	0°	0.7		REF 750-602-10 central
<u>6</u> <u>6</u> 16 - 26	R distal gingival mesial L distal occlusal	dentaform® Snap 1° Molar Bands No. 1–36 Standard Molar Bands No. 1–32	Rechts / right REF 727-365-00 Links / left REF 727-366-00	0°	0°	0.3	0°	REF 750-602-10 central



Standard Edgewise 18

Slot: 0.46 mm x 0.76 mm / 18 x 30 · cuspid, 1st bicuspid and 2nd bicuspid bands are welded conventionally with equilibrium® 2 brackets

Tooth	Position on band/buccal view Position auf dem Band/buccal	Band type / size	Brackets Buccal tubes	Torque	Angulation	In / Out mm	Distal Offset	Lingual attachment Position on band	Bondable tubes/ Bukkaloehrrchen für die Adhäsivtechnik
Zahn		Bandtyp / Größe	Brackets Bukkaloehrrchen				Anti Rotation	Linguales Teil Position auf d. Band	
<u>3</u> <u>3</u> 43 - 33	R distal incisal mesial L distal	Standard Eckzahn Cuspid No. 1–21	Rechts / right REF 790-602-02 Links / left REF 790-602-02	0°	0°	0.7		REF 750-602-10 central	
<u>4</u> <u>4</u> 44 - 34	R distal gingival mesial L distal occlusal	Standard Prämolaren Bicuspid No. 1–32	Rechts / right REF 790-602-02 Links / left REF 790-602-02	0°	0°	0.7		REF 750-602-10 central	
<u>5</u> <u>5</u> 45 - 35	R distal gingival mesial L distal occlusal	Standard Prämolaren Bicuspid No. 1–32	Rechts / right REF 790-602-02 Links / left REF 790-602-02	0°	0°	0.7		REF 750-602-10 central	
<u>6</u> <u>6</u> 46 - 36	R distal gingival mesial L distal	dentaform® Snap 1° Molar Bands No. 1–36 Standard Molar Bands No. 1–32	Rechts / right REF 724-018-00 Links / left REF 724-018-00	0°	0°	0.3	0°	REF 750-602-10 central	Rechts / right REF 724-018-51 Links / left REF 724-018-51



All brackets, buccal tubes, lingual and palatal auxiliaries can be welded onto bands in our welding department.

The illustrations shown are examples of bands with pre-welded attachments and are not identical with the described system.

Pre-welding service – banding technique

Standard Edgewise 22

Slot: 0.56 mm x 0.76 mm / 22 x 30 · cuspid, 1st bicuspid and 2nd bicuspid bands are welded conventionally with equilibrium® 2 brackets

Tooth	Position on band / buccal view Position auf dem Band / bukkal	Band type / size	Brackets Buccal tubes	Torque	Angulation	In / Out mm	Distal Offset	Palatal attachment Position on band	
Zahn		Bandtyp / Größe	Brackets Bukkalkörhrrchen				Anti Rotation	Palatinales Teil Position auf d. Band	
3 3 13 - 23	R distal gingival mesial L distal incisal	Standard Eckzahn Cuspid No. 1–25	Rechts / right REF 790-606-02 Links / left REF 790-606-02	0°	0°	0.7		REF 750-602-10 central	
4 4 14 - 24	R distal gingival mesial L distal occlusal	Standard Prämolaren Bicuspid No. 1–32	Rechts / right REF 790-606-02 Links / left REF 790-606-02	0°	0°	0.7		REF 750-602-10 central	
5 5 15 - 25	R distal gingival mesial L distal occlusal	Standard Prämolaren Bicuspid No. 1–32	Rechts / right REF 790-606-02 Links / left REF 790-606-02	0°	0°	0.7		REF 750-602-10 central	
6 6 16 - 26	R distal gingival mesial L distal occlusal	dentaform® dentaform® Snap 1° Molar Bands No. 1–36 Standard Molar Bands No. 1–32	Rechts / right REF 727-385-00 Links / left REF 727-386-00	 0°	0°	0.3	0°	REF 750-602-10 central	



Standard Edgewise 22

Slot: 0.56 mm x 0.76 mm / 22 x 30 · cuspid, 1st bicuspid and 2nd bicuspid bands are welded conventionally with equilibrium® 2 brackets

Tooth	Position on band / buccal view Position auf dem Band / bukkal	Band type / size	Brackets Buccal tubes	Torque	Angulation	In / Out mm	Distal Offset	Lingual attachment Position on band	Bondable tubes / Bukkalkörhrrchen für die Adhäsivtechnik
Zahn		Bandtyp / Größe	Brackets Bukkalkörhrrchen				Anti Rotation	Linguales Teil Position auf d. Band	
3 3 43 - 33	R distal gingival mesial L distal incisal	Standard Eckzahn Cuspid No. 1–21	Rechts / right REF 790-606-02 Links / left REF 790-606-02	0°	0°	0.7		REF 750-602-10 central	
4 4 44 - 34	R distal gingival mesial L distal occlusal	Standard Prämolaren Bicuspid No. 1–32	Rechts / right REF 790-606-02 Links / left REF 790-606-02	0°	0°	0.7		REF 750-602-10 central	
5 5 45 - 35	R distal gingival mesial L distal occlusal	Standard Prämolaren Bicuspid No. 1–32	Rechts / right REF 790-606-02 Links / left REF 790-606-02	0°	0°	0.7		REF 750-602-10 central	
6 6 46 - 36	R distal gingival mesial L distal occlusal	dentaform® dentaform® Snap 1° Molar Bands No. 1–36 Standard Molar Bands No. 1–32	Rechts / right REF 724-019-00 Links / left REF 724-019-00	 0°	0°	0.3	0°	REF 750-602-10 central	Rechts / right REF 724-019-51 Links / left REF 724-019-51



All brackets, buccal tubes, lingual and palatal auxiliaries can be welded onto bands in our welding department.

The illustrations shown are examples of bands with pre-welded attachments and are not identical with the described system.

Bands – comparison charts

The following charts featuring commonly used bands according to their sizes will help you to convert to Dentaurem bands when ordering.

Upper cuspid bands

Dentaurem		American Orthodontics	Dentsply Sirona Snap-Fit	Lancer Orthodontics	Ormco	Rocky Mountain Orthodontics	3M		
1		1		00			1+		
2		2	000	0			2+		
3		3.5			1	d	4		
4		4.5		1	2		5		
5		5	1	2		d+	5+		
6		6	2		3	e	6+		
7		7	3	3	4	e+	7+		
8		8	4	4	5		8+		
8.5		8.5		5	6	f	9		
9		9	5		7		9+		
9.5		9.5		6	8	f+	10		
10		10	6	7	9		10+		
10.5		10.5	7			g	11		
11		11		8	10		11+		
11.5		11.5	8	9	11		12		
12		12	9		12	g+	12+		
12.5		12.5		10			13		
13		13	10	11	14	h	13+		
14		13.5	11	12			14		
15		14.5	12	13	16	h+	15		
16		15.5		14			16		
17		16.5	14	15	18	i+	17		
18		17	15	16	19		17+		
19		18		17		j	18+		
20		19	17	19	22	j+	19+		

Bands – comparison charts

Lower cuspid bands

Dentaurum		American Orthodontics	Dentsply Sirona Snap-Fit	Lancer Orthodontics	Ormco	Rocky Mountain Orthodontics	3M		
1		1		0		B	1		
2		2		1	1	B+	2		
3		3				C	3		
4		4	00	2	2		3+		
5		5	0		3	C+	4+		
5.5		5.5		3			5		
6		6	1	4		D	5+		
6.5		6.5		5	4				
7		7	2			D+	6		
7.5		7.5		6	5		6+		
8		8	3	7		E	7		
8.5		8.5			6		7+		
9		9	4	8	7	E+	8		
9.5		9.5		9			8+		
10		10	5		8	F	9+		
11		11	7	11	10	F+	10		
12		12	9	13	12	G	11+		
13		13	10				12+		
14		14		15	14	H	13		
15		15	12		16	H+	14+		
16		16	13	16			15+		

Bands – comparison charts

Upper 1st and 2nd bicuspid bands

Dentaurum		American Orthodontics	Dentsply Sirona Snap-Fit	Lancer Orthodontics	Ormco	Rocky Mountain Orthodontics	3M CONTOURED	3M GENERAL PURPOSE	3M PROPORTIONED
1		1	1		1			00	
2		3	2	1		G		00+	
3		4	3		2			0	
4		5	4	2	3	G+		0+	1
5		6	5					1	
6			6	3	4	H		1+	
7		8	7		6	H+		2+	
8		9	8	5			2		5
9		10	9	6	7			3	
10		12	10	7	9			4	
11			11		10	I+		4+	
12		13	12	8	11		4		9
13		14	13	9	12	J		5	
14		15	14	10	13	J+		6	12
15		16	15	11	14		6	6+	
16		17	16			K		7	13
17		18	17	12	15		7		
18		19	18	13	16		8	7+	15
19		20	19	14	17	K+		8	
20		23	20	15	19	L	10	9+	
21		24	21	16	20			10	
22		25	22	17	21	L+			20
23			23		22	M	11	10+	
24		26	24	18	23			11	22
25		27	25	19	24		12	11+	
26		29	26	20	25	MI+		12+	24
27		30	27	21	26		13	13	
28		32	28	22	28	N	14	13+	27
29			29	23	29	N+		14	
30			30	24	31			15	29
31			31	25	32				
32			32						

Bands – comparison charts

Lower 1st and 2nd bicuspid bands

Dentaurum		American Orthodontics	Dentsply Sirona Snap-Fit	Lancer Orthodontics	Ormco	Rocky Mountain Orthodontics	3M CONTOURED	3M GENERAL PURPOSE	3M PROPORTIONED
1		1	1						
2			2	1					
3		2	3		1				1
4		3	4	2		F		1+	4
5		4	5		3		1		5
6		5	6	3	4	F+		2+	6
7		6	7	4	6				8
8		9	8	6	7	G		5	10
9			9	7	8	G+			11
10			10		9	H	5	5+	12
11		10	11	8	10				
12		11	12	9	11	H+	7	7	14
13		12	13	10					16
14		13	14	11	12	I	9	8+	17
15		15	15	12		I+			
16		16	16	13	14		11	10	19
17		17	17	14	15	J	12		21
18		18	18	15	17			12	
19		19	19	16		K	14		23
20		21	20		19		15	13+	25
21		22	21	17		K+			27
22		23	22	18	21		17	15	
23		24	23	19		L	18		28
24		25	24	20	22	L+		16	29
25		27	25	21		M			31
26		28	26	22	25			18	33
27		30	27	23	26	M+			34
28		31	28		27		22	19+	35
29		32	29	24	28				36
30			30	25	29				
31			31	26	30				38
32			32	27					

Bands – comparison charts

Upper 1st molar bands

Dentaurum Standard	Dentaurum dentafirm® Snap dentafirm®		American Orthodontics CONTOURED MR2, Mini Master S	American Orthodontics CONTOURED MR1, Master Series	Dentsply Sirona Snap-Fit	Lancer Orthodontics SUPRA™	Ormco Standard	Ormco ULTIMA® Trimline	Ormco WASHBON®	Rocky Mountain Orthodontics	TP Orthodontics LABS	3M CONTOURED	3M GENERAL PURPOSE
1			1							3			
2	1				1								
	2									5			29
	3				2	1	1						29+
3			2						1				
4	4				3	2	2			6			
5			3			3	3		2	7			30
6	5												30+
	6		4		4	4	4		3				
7			5			5	5	5					
	7						6	6	4	8	1	31	
8			6			6		7	5	9	2		31
	8			1.5	5	7	7	8		10		31+	31+
9	9		7		6	8	8		6		3	32	32
10	10		8			9	9			11	5		
	11			2.5	7		10	10	8		6	33	
11	12		9	3	8	10	11	11	9		7		33
12			10			11			10	12	8	33+	33+
13	13		11	3.5	9		12	12			9	34	34
14			12			12	13	13		13	10		
	14		13	4	10	13	14	14	11	14	11	34+	34+
15	15				11		15		12	15	12		
16	16		14	4.5	12	14		15	13	16	13	35	35
17	17		15	5	13	15		16	14	17	14	35+	35+
18					16	16	17	17	15	18	15		
	18		16		14	17			16		16	36	36
19				5.5			18	18			17		36+
	19		17	6	15	18		19	17	19	18		
20	20		18		16	19	19	20		20	19	36+	
	21		19	6.5	17				18	21	20	37	37
21				7		20	20	21		22	21		
	22		20		18	21		22	19	23	22	37+	37+
22				7.5			21			24			
23	23		21		19	22	22		20	25	23		38
24	24		22	8	20	23	23	23	21	26	24	38	
	25		23		21				22	27	25		
25	26				22	24	24	24		28		38+	39
26	27				23	25	25		24	29	27		
	28				24	27		26	26	31	29	39+	
27	29		26	9.5	25		27	27	27	32	30		39+
	30		27		26	28	28	28		33	31		40
28						29		29	28	34		40	40+
29	31		28	10	27		29			35	32	40+	
30	32			10.5	28	31	30	31	29				
31	33		30	11	29		31	32	30	36		41	
32	34		32	11.5	30	32							41+

Bands – comparison charts

Lower 1st molar bands

Dentaurum Standard	Dentaurum dentiform® Snap dentiform®		American Orthodontics CONTOURED MR2, Mini Master S	American Orthodontics CONTOURED MR1, Master Series, MSH	Dentsply Sirona Snap-Fit	Lancer Orthodontics SUPRA™	Ormco Standard	Ormco ULTIMA® Trimline	Ormco WASHBON®	Rocky Mountain Orthodontics	TP Orthodontics LABS	3M CONTOURED	3M GENERAL PURPOSE
	1									2			28+
	2									3			
	3				1								29+
	4									4			
1			1				1						
2	5				2								30
3			2			1	2		1	5			
4	6		3		3	2							
	7						3		2	6			30+
5			4			3							
	8				4	4	4		3	7			
6	9			1.5		5	5					31	31
7	10		5		5	5	6	5	4	8	1		
8	11		7	2	6	7	7	7		9	2	31+	
	12		8			7	8		6	10	4		32
9				2.5		8		8	7		5	32	
10	13		9		7		9	9		11	6		32+
11	14		10		8	9	10		8		7	32+	
12				3		10	11	10	9	12	8	33	33
	15		11		9				10				
	16		12		10	11	12	11		13	9		33+
13				3.5		12	13	12		14	10	33+	
14	17		13		11		14	13	11	15	11	34	34
15	18			4	12	14		14	13	16	12	34+	34+
	19		14	4.5	13		16	15		17	13	35	35
16			15			15			14		14		
	20		16	5	14	16	17	16	15	18	15	35+	35+
17	21		17		15			17	16		16		
18			18			17	18	18		19	17		
	22		19		16	18			17	20	18	36	36
19	23		20	6	17		19	19		21	19		
20				6.5		19			18	22	20	36+	36+
21	24			7	18	20	20	20		23	21		
	25		21		19			21	19	24		37	37
22	26		22	7.5	20	21	21	22			22		
	27					22	22		20			37+	37+
23	28		23		21		23	23	21	25	23		
	29			8	22	23		24	22	26	24	38	38
24	30				23	24	24			27	25		
25			25			25	25	26	24	29	27		38+
26			26	8.5		26	26		25	30	28	39	39
27			27		24			27	26	31		39+	
	31		28	9	25	27	27		27	32			39+
28	32			9.5		28	28	28		33		40	
29			29			29	29	30	28				40
30				10.5	26	30		31		34		40+	40+
31	33		30		27		30		29	35			
32	34		32		29	31		33		36		41	41
	36												

Bands – comparison charts

Upper 2nd molar bands

Dentaurum dentaform®		American Orthodontics	Dentsply Sirona Snap-Fit	Lancer Orthodontics	Ormco	Rocky Mountain Orthodontics	TP Orthodontics LABS	3M
								28
1			1					28+
2			2					29
3			3					29+
4			4	1				
5			5		1			30
6				2	2			30+
7			6	3				
8			7	4	3			31
9			8	5				31+
10			9		4		1	
11			10	6	5		2	
		1.5		7		7		32
12			11		6		3	
		2		8	7		4	32+
13			12	9			5	
		2.5			8	8	6	33
14			13	10	9		7	
		3		11	10	9	8	33+
15			14				9	
		3.5		12	11	10	10	34
16			15	13	12		11	
					13		12	34+
		4.5		14		11	13	
17			16	15	14	12		
18		5	17		15	13	14	35
19			18	16	16	14	15	
20			19	17		15	16	35+
		5.5			17		17	
21			20	18		16	18	36
22			21	19	18	17	19	
		6.5				18	20	36+
23			22	20	19	19	21	
24			23	21		20	22	37
25		7.5	24		20	21		
				22	21	22	23	37+
26		8	25	23	22		24	
						23	25	38
27			26	24	23	24		
28		8.5	27	25	24	25	26	38+
					25		27	
29		9	28	26	26	26	28	39
30			29	27	27	27	29	
		9.5				28	30	39+
31			30	28	28	29	31	
				29		30		40
32		10				31	32	
				30	29			40+
		10.5		31	30	32		
					31			41
				32			33	
		11.5			32			
							34	
							35	

Bands – comparison charts

Lower 2nd molar bands

Dentaurum dentaurum®		American Orthodontics	Dentsply Sirona Snap-Fit	Lancer Orthodontics	Omco	Rocky Mountain Orthodontics	TP Orthodontics LABS	3M
								28
				29.5				28+
								29
1			3	30	1			29+
2			4		2			30
3				30.5				30+
4			5		3		1	
5		1.5	6	31			2	31
				31.5	4		3	
				32	5			31+
6		2				7		
7			7		6		4	32
8		2.5	8		7		5	
9			9	32.5			6	32+
10		3	10		8	8	7	
			11	33	9		8	33
		3.5			10	9		
11			12	33.5			9	33+
						10	10	
12			13	34	11		11	34
13		4	14	34.5	12			
		4.5			13	11	12	34+
14						12	13	
15			15	35	14	13	14	35
		5			15	14	15	
16		5.5	16	35.5	16	15	16	35+
							17	
17		6	17	36	17	16	18	36
18		6.5	18	36.5		17	19	
19		7		37	18	18	20	36+
20			19			19	21	
		7.5			19	20		37
21			20	37.5		21	22	37+
					20	22		
22		8	21	38	21		23	38
23			22		22	23	24	
24			23	38.5	23	24	25	38+
					24	25	26	
25			24	39	25		27	39
26		8.5	25	39.5	26	26	28	
27			26		27	27		
		9				28		39+
28		9.5	27		28	29		
29		10				30		40
30			28			31		
					29			40+
31		10.5	29			32		
					30			41
32		11	30		31			
					32			

Wires – information

Comparison between Dentaurem wires/ Brand names of other competitive products

Dentaurem	American Orthodontics	Dentsply Sirona	FORESTADENT	Lancer Orthodontics	Leone	MASEL
dentaflex® twisted wires	TWIST	Wildcat® Tricat	Forestaflex®	Swiss Strand™	Twist	Triple Twist R™
dentaflex® co-axial wires	CO-AX	PentaCat® Hexacat	Forestaflex®	Memorex	Flex	Penta-One™
dentaflex® braided wires	Streight™ Woven	Quad-Cat Multibraid Arches	Forestaflex®	Simplicity Wire Braided NiTi	Twist	Octava™
remanium® stainless steel wires	SS STAINLESS STEEL Wire	Nubryte Gold Super Gold Stainless Steel	Stahl Forestanit® Forestanit®-Super	High Tensile Stainless Steel X-Resilient Stainless Steel	Leowire® Stainless Steel	Resilient Bright Golden
Equire stainless steel wires	SS STAINLESS STEEL Wire	Nubryte Gold Super Gold Stainless Steel	Stahl Forestanit® Forestanit®-Super	High Tensile Stainless Steel	Leowire® Stainless Steel	Resilient Bright Golden
remaloy®/Redur® cobalt-based alloy	CC Chromium-Cobalt (Multiphase blue)	Coballoy	Forestalloy	Bendalloy™ Blue Ortholoy Blue	Leoloy® chromium-cobalt wire Croizat	Masiloy™
rematitan® LITE nickel titanium alloy super elastic		Sentalloy® superelastic Lowland NiTi	Titanol® - Superelastic (SE) - Low Force	Titalan™ Super Elastic Orthinol II	Memoria® superelastic	Elastinol™ superelastic Authentic Nitinol Value Brand of NiTi Econo Nickel Titanium
rematitan® SPECIAL beta titanium alloy	BETA TITANIUM	Resolve® Beta Titanium	Betaflex Titan	Titalan™ Beta-Titanium Molybdenum	Beta Memoria®	CNA BETA III™
Tensic® nickel titanium, thermo-active	Therma-Ti® Lite Titanium Memory Wire	Neo Sentalloy® Heat-Activation	Titanol® Martensitic BioStarter BioTorque BioFinisher Martensitic	Titalan LT™ Thermodynamic	Thermomemoria®	Therm-A-Form Ultra™ CV SuperForce™
Equire nickel titanium, thermo-active	Therma-Ti® Lite Titanium Memory Wire	Neo Sentalloy® Heat-Activation	Titanol® Budget Titanol® Martensitic BioStarter BioTorque BioFinisher Martensitic	Titalan LT™ Thermodynamic	Thermomemoria®	Therm-A-Form Ultra™ CV SuperForce™
rematitan® LITE White nickel titanium alloy super elastic	Iconix® Aesthetic SE NiTi	Spectra™ Imagination Sentalloy® High Aesthetic	Titanol® Cosmetic	Titalan™ TS Coated	Memoria® Natura®	Orthocosmetic Elastinol™
Tensic® White nickel titanium, thermo-active			BioCosmetic®		Memoria® Mimetic rhodium coated	

Note: Some of the items listed are registered brand names e.g. TMA® from Ormco.

Comparison between Dentaureum wires / Brand names of other competitive products

Dentaureum	Ormco	Rocky Mountain Orthodontics	SCHEU-DENTAL	TP Orthodontics	World Class Orthodontics Ortho Organizers GmbH	3M
dentaflex® twisted wires	Triple Flex™	Tri-Flex™		Twisted 3-Strand Wire	Twist Stainless Steel	Hi-T™ II Twist Flex Twist
dentaflex® co-axial Wires	Respond®	Supra-Flex™		Co-Ax®	Coaxial Stainless Steel	Coaxial
dentaflex® braided wires	D-Rect® Turbo Wire	Flex-VIII™		Woven 8-Strand Wire	Braided Stainless Steel	Braided
remanium® stainless steel wires	Stainless Steel	Tru-Chrome®	Chromium	Shiny Bright Wire TP Original Stainless Steel	Stainless Steel	Permachrome Hi-T™ II Stainless Steel
Equire stainless steel wires	Stainless Steel	Tru-Chrome®	Chromium	Vela™ Shiny Bright Wire TP Original Stainless Steel	Stainless Steel	Permachrome Hi-T™ II
remaloy® / Redur® cobalt-based alloy	Azurloy™	Elgiloy® blue			Ortholoy-Blue	Blue Flexiloy®
rematitan® LITE nickel titanium alloy super elastic	Ni-Ti® Align™ SE 200 Align™ XF	Orthonol® Superelastic	Lazium®	Reflex® Superelastic	Nitanium S.E. Superelastic	Nitinol® SE Superelastic Nitinol® Classic
rematitan® SPECIAL beta titanium alloy	TMA® Titanium Molybdenum	Bendaloy™ Titanium- Molybdenum		TiMolium®	CNA™ Beta III	Beta III Titanium
rematitan® sl nickel titanium alloy super elastic					Nitanium Black Ti™	
Tensic® nickel titanium, thermo-active		Thermaloy® Thermally activated		Reflex® Heat activated	Accel™ Bio-Kinetix™ Bio-Kinetix™ Plus	Nitinol® Heat activated
Equire nickel titanium, thermo-active		Thermaloy® Thermally activated		Vela™ Reflex® Heat activated	Accel™ Bio-Kinetix™ Bio-Kinetix™ Plus	Nitinol® Heat activated
rematitan® LITE White nickel titanium alloy super elastic		Orthonol® FLI® Wire coated		Reflex® Aesthetic Wire	Nitanium Plastic Coated WOOWA Bio- Aesthetic Coated	
Tensic® White nickel titanium, thermo-active					Bio-Kinetix™ Micro Coated	

Note: Some of the items listed are registered brand names e.g. TMA® from Ormco.

Materials list

Product groups	Product	Identification no. *
Aligner	prime4me® Aligner	4.102
	prime4me® Retainer	4.102
	prime4me® Aligner Attachment Guide	4.004
Lingual system	discovery® delight – lingual brackets / molar tubes	1.005
	rematitan® LITE lingual arch	3.200
	remanium® lingual arch	1.012 or 1.013
	dentaflex® lingual arch	1.012 or 1.013
Self-ligating brackets	dinamique® c	5.000 (body), 3.300 (clip), rhodium coating
	dinamique® m	1.009 (body), 3.300 (clip)
	discovery® sl / discovery® sl 2.0	1.003, 1.009, 3.005
Metal brackets	discovery® / discovery® smart	1.005
	topic	3.004
	equilibrium® 2 / equilibrium® mini	1.005
	equilibrium® 2 / equilibrium® mini with hook	1.001, 1.005
	equilibrium® ti	3.104
	equilibrium® ti with hook	3.102, 3.104
Aesthetic brackets	discovery® pearl	5.000, 4.003 (positioning guide)
	Elegance®	1.000, 4.002 glass-fiber reinforced
	Fascination® / Fascination® 2	5.000
	Jewels	5.000
Bracket accessories	Titanium button with chain acc. to Watted	3.100, 3.101
	Lingual button for bonding technique	1.002
	Lingual double hook for bonding technique	1.005, 1.008
	Mesh bases for bonding technique	1.005
	Mesh bases for buccal tubes for bonding technique	1.000, 1.001
	Gauze	1.000
Buccal tubes w.h. = with hook w.t. = with tube w.h. & w.t. = with hook and tube	Ortho-Cast / Ortho-Cast w.h. / Ortho-Cast w.h. & w.t. / Ortho-Cast M-Series	1.005
	Ortho-Cast w.h. and central tube	1.000, 1.005
	Ortho-Cast with large base / Ortho-Cast w.h. and large base	1.000, 1.001, 1.005
	Ortho-Cast w.h. & w.t. occlusal / gingival	1.005, 1.008
	Ortho-Cast mini / Ortho-Cast sl	1.005
	Ortho-Cast NF	3.004
	Ortho / Ortho w.t.	1.000, 1.002
	Ortho w.h. / Ortho w.h. & w.t.	1.000, 1.001, 1.002
	Ortho for bonding technique / Ortho w.t. occlusal for bonding technique	1.000, 1.001, 1.002
	Ortho with mesial / distal hook and central tube	1.000, 1.002, 1.005
	rematitan®	3.100, 3.101
	rematitan® with hook	3.101, 3.102
	Lingual / palatal sheaths	1.008
	Lingual / palatal sheaths with mesial / distal hooks	1.000, 1.008
	Lingual / palatal sheaths with distal hooks	1.002, 1.008
	Lingual / Palatal sheaths easy in	1.008
	Button	1.002
	Band seating lug	1.005
	Double hook / Hook with seating lug	1.008
	Button hook	1.001
	Insert-Adapter	1.003
Bands	dentaform® Snap	1.008
	dentaform®	1.001
	Standard bicuspid and cuspid bands	1.001
	Standard 1st and 2nd molar bands	1.008
	Band material on rolls	1.001
	Matrix band on rolls	1.003
Arches	remanium® ideal arch	1.013
	Equire preformed ideal arch	1.013
	dentaflex® ideal arch	1.012 or 1.013
	Noninium® ideal arch	REF XXX-XXX-00 / 1.006
		REF XXX-XXX-10 / 1.014
	Noninium® White ideal arch	1.006, 4.008
	Tensi® ideal arch	3.203

Product groups	Product	Identification no.*
Arches	Equire thermo-active preformed ideal arch	3.200
	Tensic® White ideal arch	3.200, 4.008
	rematitan® sl ideal arch	3.200
	rematitan® LITE ideal arch / Ricketts® PENTA ideal arch	3.202
	rematitan® LITE Spee® retraction arch	3.202
	rematitan® LITE White ideal arch	3.202, 4.008
	remaloy® ideal arch	3.006
	rematitan® SPECIAL ideal arch	3.103
	Translucent ideal pearl arches	4.005 + fiber glass
Wire	remanium® straight wire / wire on coils	1.003
	dentaflex® straight wire, round, 3-strand twisted	1.003
	dentaflex® straight wire, round, 6-strand co-axial	1.012 or 1.013
	Noninium® straight wire / wire on coils	REF XXX-XXX-00 / 1.006
		REF XXX-XXX-10 / 1.014
	remaloy® straight wire	3.006
	rematitan® SPECIAL straight wire	3.103
	dentaflex® wire on coils, round, 3-strand twisted / 6-strand co-axial	1.011
	rematitan® LITE wire on coils	3.200
Retainer	rematitan® flat retainer wire	1.013
	Titanium retainer wire, Ti1	3.100
	Titanium retainer wire, Ti5	3.102
	Gold retainer wire	2.301
	remanium® lingual retainer	1.011, 1.013
	Titanium lingual retainer grade 4	3.101
	prime4me® RETAIN3R	3.102
Wires / Arches / Accessories Accessories	remanium® ligature wire	1.005
	remanium® preformed ligatures, short	1.005
	remanium® preformed ligatures, long	1.012 or 1.013
	remanium® Kobayashi ligatures	1.012 or 1.013
	Ligatures White, Kobayashi White ligatures	1.013, 4.008
	remanium® tension springs / compression springs	1.003
	rematitan® LITE tension springs / compression springs	3.200
	Cuspid retraction springs	1.000, 1.003
	Jaw fracture splint	1.000
	Sliding tubes, round with hook	1.000, 1.001
	Clamping tubes, rectangular with hook	1.005
	Clamping tubes, rectangular with wire	1.003 + 1.005
	Cross tube, rectangular tube	1.000
	Stop tubes	1.013
	Noninium® standard labial arches	1.014
	Noninium® Adams clasps / triangular clasps / arrow clasps	1.014
	Noninium® ball retainer clasps	REF XXX-XXX-01 / 1.014
		REF XXX-XXX-02 / 1.014
	remanium® ball retainer clasps / remanium® arrow clasps	1.001
	Guide pins acc. to G. Müller®	1.003
	Guide pins acc. to Hinz	1.012 or 1.013
	Lingual bar	1.003
	Strengtheners	1.005
	Stainless steel tubes	1.000
	Plastic sleeving	4.009
	remanium® Goshgarian palatal bars	1.013
	remanium® Quad Helix	1.013
	Orthorama® palatal arches / Lingual arches / Sectional arches / Multi-Action	3.006
	Orthorama® transfer jigs	1.000
	Orthorama® retention rings	1.001 + 1.008
	Grid strengtheners stainless steel	1.000
	Grid strengtheners gold-plated	1.000 gold-plated (24 carat)
Class II appliances	SUS2 / SUS3 telescope rod, assembled	1.002, 1.003
	SUS2 / SUS3 telescope rod	1.003
	SUS2 / SUS3 arch adapter / fixing screw	1.002

Materials list

Product groups	Product	Identification no.*
Class II appliances	SUS ² telescope element	1.000, 1.003
	SUS ³ telescope element	1.000, 1.003, 1.015
	SUS ² / SUS ³ turbo spring / pressure spring (stainless steel)	1.003
	SUS ² compression spring (nickel titanium)	3.200
	SUS ² / SUS ³ distance spring	1.003
	Spacer rings, crimpable	1.002
	SUS ² / SUS ³ ball retainer clasp	1.000
	Herbst I hexagon socket screw / slotted screws	1.002
	Herbst base / pivot	1.002
	Herbst telescope tubes	1.009
	Herbst telescope tubes with ball joint	1.000, 1.002, 1.003
	Herbst TS telescopic tube	1.000, 1.002, 1.003, 1.009
	Herbst telescope rods	1.002
	Herbst spacer rings	1.002
Intra- and extra-oral	Metal clip for safety modules	1.003
	Reverse-pull headgear, Tübingen model	1.000, 1.001, 1.002, 4.100, suede
	Delaire face mask	1.003, 4.100, suede
	Noninium® facebow, medium	1.006, 4.001
	Standard facebow, with stop loops	1.000, 1.003, 4.001
	Standard facebow, with stop loops, with canine hook	1.000, 1.001, 1.003, 4.001
	Chin cap, elastic	4.009, polyester, 1.003
	Chin cap, rigid, with short hooks	1.002, 1.003, 4.001, 4.100
	Chin cap, rigid, with long hooks	1.003, 4.100, suede
	Chin cap for Tübingen model	4.100, suede
	Chin cap liner	4.009, polyester
	Chin cap liner, Delaire	Suede
	J-hooks for Hickham headgear	1.011, 1.013
	Dentalastics® extra-oral plastic chain	4.009 **
	Elastics, extra-oral, latex	4.201 **
	Lip bumper	1.000, 1.001, 3.006, 4.003
	Stop tubes	1.000
	Stop screws	1.002
	Activator tubes	1.000, 1.001
	Elastic Placers	4.006
	Elastics, intra-oral, latex	4.201 **
	Elastics, intra-oral, non-latex	Thermoplastic acrylic**
	Elastic ligature	4.009
	Elasto-Force ligature thread	4.009 **
	Dentalastics® rotation wedge / plastic ligatures / separators / Personal ligatures	4.009
	Elasto-Force, SUPER Elasto-Force plastic chains	4.009
	Oral screen Ulmer Modell	4.009
	FaceFormer ONE blue with hygiene box	4.203
	OrthoSafe mouthguard	4.101
tomas® anchorage system	tomas®-pin	3.102
	tomas®-abutment universal	1.005
	tomas®-abutment tube	1.000, 1.005
	tomas®-RPE eyelet	1.000
	tomas®-abutment wire	1.000, 1.005, 3.000
	tomas®-abutment median	1.000, 1.005, 3.000
	tomas®-abutment plain	1.000, 1.005
	tomas®-abutment EP	1.000
	tomas®-compression spring	3.200
	tomas®-coil spring	3.200 + 1.013
	tomas®-uprighting spring	1.000, 1.003, 3.200
	tomas®-stop screw	1.000, 1.002
	tomas®-slotted stops	1.013
	tomas®-hook / tomas®-T-wire	1.003
	tomas®-cross tube	1.000
	tomas®-power arm	1.003, 1.005

Product groups	Product	Identification no.*
tomas® anchorage system	tomas®-crimp hook	1.005
	tomas®-aligner hook	1.005
	tomas® PI 4.0	3.102
	tomas® PI-abutment multiple	1.000
	tomas® PI-abutment screw	3.102
	tomas® PI-abutment plain	1.000
	tomas® PI-abutment plain set with wire	1.000, 3.000, 3.102
	tomas® PI-straight wire with notch	3.000
Distalization appliance – amda®	amda® telescope	1.000, 1.002, 1.003
	amda® connector	1.003
	amda® palatal arch	1.000
Expansion screws	Standard expansion screws	1.002, 1.003, 4.006
	rematitan® expansion screw Medium	3.101, 3.102, 4.006
	Trapezoidal expansion screw, sectional screw Medium	1.002, 1.003, 4.006
	Sectional screw Mini, retraction screw	1.002, 1.003
	Neo-Bertoni	1.002, 1.003, 1.005, 4.006
	Three dimensional screw acc. to Bertoni / acc. to Steiner	1.002, 1.003
	Fan-type expansion screw	1.000, 1.002
	Mandibular bow screw acc. to Müller	1.002, 1.003, 4.006
	Neo-bite jumping appliance acc. to F.M. Sander	1.002, 1.003, 1.019, 4.006
	Guide ridge for flexible forward thrust double plate system acc. to Schaneng	1.002, 1.003, 1.015
	Expansion screw for split activator	1.002, 1.003
	Expansion screw for progenia activator acc. to Weise	1.002
	Piston spring screw	1.002, 1.003
	Serrated housing nut	1.002
	Variety / Variety SP / Variety click / Variety click SP expansion screw	1.000, 1.002, 1.003
	Titan hyrax® Maxi expansion screw	3.102, 3.104
	rematitan® straight wire, round	3.100
	hyrax® / hyrax® click / hyrax® neo expansion screw	1.000, 1.002, 1.003
	Arrow to indicate the opening direction	4.006
	Expansion screw key, standard	1.003
	Expansion screw key, large	1.002 + 1.003
	Safety key for hyrax® / Variety screw	1.003, 4.100

* **Note:** If several materials are listed, this means that the product consists of all these materials.

** Also contains traces of other elements. Please ask if you have any queries.

Materials list

1. Stainless steel

Composition in weight %

Identification No.	Material		C	Si	Mn	Cr	Mo	Ni	P	S	Other	Fe
	Designation	Number (DIN, AISI)*										
1.000	Stainless steel	1.4301	≤ 0.07	≤ 1.0	≤ 2.0	17.0–19.5	–	8.0–10.5	≤ 0.045	≤ 0.03	N ≤ 0.11	Residue
1.001	Stainless steel	1.4303	≤ 0.06	≤ 1.0	≤ 2.0	17.0–19.0	–	11.0–13.0	≤ 0.045	≤ 0.03	N ≤ 0.11	Residue
1.002	Stainless steel	1.4305	≤ 0.10	≤ 1.0	≤ 2.0	17.0–19.0	–	8.0–10.0	≤ 0.045	0.15–0.35	N ≤ 0.11/ Cu ≤ 1.0	Residue
1.003	Stainless steel (dentaflex®, remanium®)	1.4310	0.05–0.15	≤ 2.0	≤ 2.0	16.0–19.0	≤ 0.8	6.0–9.5	≤ 0.045	≤ 0.015	N ≤ 0.11	Residue
1.004	Stainless steel	1.4401	≤ 0.07	≤ 1.0	≤ 2.0	16.5–18.5	2.0–2.5	10.0–13.0	≤ 0.045	≤ 0.03	N ≤ 0.11	Residue
1.005	Stainless steel (remanium®)	1.4404	≤ 0.03	≤ 1.0	≤ 2.0	16.5–18.5	2.0–2.5	10.0–13.0	≤ 0.045	≤ 0.03	N ≤ 0.11	Residue
1.006	Stainless steel (Noninium®)	1.4456	≤ 0.1	≤ 1.0	16.0–20.0	16.0–20.0	1.8–2.5	< 0.2	≤ 0.05	≤ 0.05	V ≤ 0.2 N 0.7–1.0	Residue
1.007	Stainless steel	1.4460	≤ 0.05	≤ 1.0	≤ 2.0	25.0–28.0	1.3–2.0	4.5–6.5	≤ 0.035	≤ 0.3	N 0.05–0.20	Residue
1.008	Stainless steel	1.4541	≤ 0.08	≤ 1.0	≤ 2.0	17.0–19.0	–	9.0–12.0	≤ 0.045	≤ 0.3	Ti 5 x C to 0.7	Residue
1.009	Stainless steel	1.4542	≤ 0.07	≤ 0.7	≤ 1.5	15.0–17.0	≤ 0.6	3.0–5.0	≤ 0.040	≤ 0.03	Cu 3.0–5.0 Nb 5 x C to 0.45	Residue
1.010	Stainless steel	1.4435	≤ 0.03	≤ 1.0	≤ 2.0	17.0–19.0	2.5–3.0	12.5–15.0	≤ 0.045	≤ 0.025	N ≤ 0.11	Residue
1.011	Stainless steel (dentaflex®)	AISI 302	≤ 0.15	≤ 1.0	≤ 2.0	17.0–19.0	–	8.0–10.0	≤ 0.045	≤ 0.03	–	Residue
1.012	Stainless steel (dentaflex®)	AISI 302 B	≤ 0.15	2.0–3.0	≤ 2.0	17.0–19.0	–	8.0–10.0	≤ 0.045	≤ 0.03	–	Residue
1.013	Stainless steel (dentaflex®, remanium®)	AISI 304	≤ 0.08	≤ 1.0	≤ 2.0	18.0–20.0	–	8.0–10.5	≤ 0.045	≤ 0.03	–	Residue
1.014	Stainless steel (Noninium®)	1.3808	0.15–0.25	0.2–0.6	9.5–12.5	16.5–18.0	2.7–3.7	< 0.05	0.02	0.01	Cu 0.25 N 0.45–0.55	Residue
1.015	Stainless steel	1.4568	≤ 0.03	≤ 0.70	≤ 1.0	16.0–18.0	–	6.5–7.8	≤ 0.04	≤ 0.015	Al 0.70–1.50	Residue
1.016	Stainless steel	1.4197	0.20–0.26	≤ 1.0	≤ 2.0	12.5–14.0	1.1–1.5	0.75–1.50	≤ 0.04	0.15–0.35	–	Residue
1.017	Stainless steel	1.4034	0.43–0.50	≤ 1.0	≤ 1.0	12.5–14.5	–	–	≤ 0.04	≤ 0.03	–	Residue
1.018	Stainless steel	1.4571	≤ 0.08	≤ 1.0	≤ 2.0	16.5–18.5	2.0–2.5	10.5–13.5	≤ 0.045	≤ 0.03	Ti 5 x C to 0.7	Residue
1.019	Stainless steel	1.4501	≤ 0.030	≤ 1.0	≤ 1.0	24.0–26.0	3.0–4.0	6.0–8.0	≤ 0.035	≤ 0.015	N 0.20–0.30 Cu 0.50–1.00 W 0.50–1.00	Residue

*Note: Stainless steel with foreign standards which correspond to the DIN material numbers are also acceptable.

2. Copper and precious metal alloys

Composition in weight %

Identification No.	Material		Ag	Cu	Ni	Pb	Zn	Fe	Mn	Sn	Sb	Al	Others
	Designation	Number (DIN)*											
2.000	Argentan	2.0770	—	45–48	9–11	0.5–2.0	38–45	≤ 0.05	≤ 0.5	≤ 0.3	—	—	≤ 0.1
2.001	Argentan	2.0780	—	56–58	11–13	0.3–1.5	26–33	≤ 0.05	≤ 0.5	≤ 0.3	—	—	≤ 0.1
2.002	Argentan	2.0790	—	59–63	17–19	0.3–1.5	Residue	≤ 0.03	≤ 0.7	—	—	—	≤ 0.4
2.100	Brass	2.0321	—	62.0–65.5	< 0.3	≤ 0.1	Residue	≤ 0.01	≤ 0.1	≤ 0.1	≤ 0.1	—	—
2.101	Brass	2.0360	—	59.0–62.0	< 0.3	≤ 0.1	Residue	≤ 0.01	≤ 0.1	≤ 0.1	≤ 0.1	—	—
2.200	Silver solder	2.5153 L-Ag75	74–76	Residue	—	≤ 0.02	2.0–4.0	—	—	—	—	≤ 0.005	≤ 0.1
2.201	Silver solder	2.5147 L-Ag44	43–45	29–31	—	≤ 0.02	Residue	—	—	—	—	≤ 0.005	≤ 0.1
2.202	Silver solder	2.5159 L-Ag55Sn	54–57	21–23	—	≤ 0.02	Residue	—	—	2.0–5.0	—	≤ 0.005	≤ 0.1
2.203	Silver solder	2.5151 L-Ag72	71–73	Residue	—	≤ 0.02	—	—	—	—	—	≤ 0.005	≤ 0.1
2.204	Silver solder	—	43–46	18–22	—	—	6–10	—	—	2.0–6.0	—	—	≤ 0.3
2.205	Silver solder	—	57–61	15–18	—	—	Residue	—	—	—	—	—	≤ 0.3
2.300	Gold solder	—	38.5–39.5	Residue	—	—	—	—	—	—	—	—	Au 33.0
2.301	Au-Pt alloy	—	16–17	8–10	—	—	—	—	—	—	—	—	Au 60–62 Pt 13–14

Note: The materials listed in tables 1–4 enable easy identification of chemical compositions.

Not all of the listed materials are used by Dentaaurum. You can find the most recent materials list at www.dentaaurum.com

Materials list

3. Non-precious metal alloys

Composition in weight %

Identification No.	Material		Ni	Fe	Mn	Cr	Mo	W	Ti	Si	C	Co	Al	S	O	H	N	Others
	Designation	Number (DIN)*																
3.000	CoCr alloy	—	19–23	4–6	≤ 1.0	18–22	3–5	3–5	0.1–2.0	≤ 0.5	≤ 0.03	Residue	—	≤ 0.1	—	—	—	—
3.002	CoCr alloy	—	—	27–31	≤ 1.0	28–32	4–6	—	—	≤ 0.1	≤ 0.35	31–35	—	—	—	—	—	—
3.003	CoCr alloy	—	< 1.0	—	≤ 0.3	28–32	4–6	—	—	≤ 1.0	≤ 1.0	Residue	—	—	—	—	≤ 0.3	—
3.004	CoCr alloy (topic, Ortho-Cast NF)	—	< 0.1	—	≤ 1.0	26–30	—	8–10	—	0.5–2.5	—	Residue	—	—	—	—	≤ 1.0	Nb < 1.0
3.005	CoCr alloy	—	15–16	Residue	1.5–2.0	19–21	6.5–7.5	—	—	≤ 1.2	≤ 0.15	39–41	—	≤ 0.015	—	—	—	P ≤ 0.015 Be ≤ 0.001
3.006	CoCr alloy (remaloy®)	—	Residue	5.29	0.75	18.55	3.85	3.86	0.94	0.34	0.003	45.15	—	—	—	—	—	—
3.100	Titanium grade 1 DIN EN ISO 5832-2	3.7025	—	< 0.2	—	—	—	—	> 99.5	—	< 0.08	—	—	—	< 0.18	< 0.0125	< 0.03	—
3.101	Titanium grade 4 DIN EN ISO 5832-2	3.7065	—	< 0.5	—	—	—	—	> 99.0	—	< 0.08	—	—	—	< 0.40	< 0.0125	< 0.05	—
3.102	TiAl6V4 Titanium grade 5 DIN EN ISO 5832-3	3.7165	—	< 0.3	—	—	—	—	90.0	—	< 0.08	—	6.0	—	< 0.20	< 0.0150	< 0.05	V 4.0
3.103	TiMo alloy (rematitan® SPECIAL)	—	—	—	—	—	11.5	—	78	—	—	—	—	—	—	—	—	Zr ≤ 6 Sn ≤ 4.5
3.104	Titanium grade 2 DIN EN ISO 5832-2 (equilibrium® ti, Titan hyrax®)	3.7035	—	< 0.3	—	—	—	—	> 99.3	—	< 0.08	—	—	—	< 0.25	< 0.0125	< 0.03	—
3.200	NiTi alloy (Equire)	—	50–60	≤ 0.5	—	—	—	—	Residue	—	≤ 0.1	—	≤ 0.1	—	≤ 0.1	≤ 0.01	≤ 0.01	—
3.201	NiTiCu alloy	—	50–60	≤ 0.5	—	—	—	—	Residue	—	≤ 0.1	—	≤ 0.1	—	≤ 0.1	≤ 0.017	≤ 0.01	Cu < 1.0
3.202	NiTi alloy (rematitan® LITE)	—	54.5–57	< 0.005	—	< 0.005	—	—	Residue	—	< 0.0029	< 0.005	—	—	< 0.05	< 0.005	< 0.002	Cu < 0.005 Nb < 0.005
3.203	NiTi alloy (Tensic®)	—	35–60	0–6	—	—	—	—	20–50	—	—	—	—	—	—	—	—	Cu 0–15 HF 0–40 V 0–10 Nb 0–20 B 0–1
3.300	NiCoCrMo alloy	MP35N®	33–37	≤ 1	≤ 0.15	19–21	9–10.5	—	≤ 1	≤ 0.15	≤ 0.025	Residue	—	≤ 0.01	—	—	—	P ≤ 0.015 B ≤ 0.015

4. Plastics

Identification No.	Material	
	Designation	Symbol (DIN)*
4.000	Epoxide	EP
4.001	Polyamide	PA
4.002	Polycarbonate	PC
4.003	High-density polyethylene	PE-HD
4.004	Low-density polyethylene	PE-LD
4.005	Polymethylmethacrylate	PMMA
4.006	Polypropylene	PP
4.007	Polystyrene	PS
4.008	Polytetrafluorethylene	PTFE
4.009	Polyurethane	PUR
4.010	Polyvinylchloride	PVC
4.011	Polyvinylsiloxane	
4.100	Acrylonitrile butadiene styrene	ABS
4.101	Ethylene / vinyl acetate	EVA
4.102	Polyethylene terephthalate G Copolyester	PETG
4.200	Synthetic isoprene rubber	IR
4.201	Natural rubber	NR
4.202	Silicone rubber	Q
4.203	Silicone	SI

5. Ceramic materials

Composition in weight %

Identification No.	Material		Al ₂ O ₃	ZnO	MgO	Other
	Designation	Number (DIN)*				
5.000	Aluminum oxide	—	99.99	—	—	≤ 0.01
5.001	Band cement	—	—	80–86	6–10	≤ 10

Alloy elements

Ag	silver	O	oxygen
Al	aluminum	P	phosphorus
Au	gold	Pb	lead
Be	beryllium	Pt	platinum
C	carbon	Rh	rhodium
Cr	chromium	S	sulphur
Cu	copper	Sb	antimony
Co	cobalt	Si	silicon
Fe	iron	Sn	tin
H	hydrogen	Ti	titanium
Mn	manganese	V	vanadium
Mo	molybdenum	W	tungsten
N	nitrogen	Zn	zinc
Nb	niobium	Zr	zirconia
Ni	nickel		

Note: The materials listed in tables 1– 4 enable easy identification of chemical compositions.
Not all of the listed materials are used by Dentaaurum. You can find the most recent materials list at www.dentaaurum.com

Conversion charts

Note:

- ∅ round
- ▧ rectangular

In addition to the metric measurements provided, the wires, brackets and buccal tubes are also available with an additional description in inches (DIN EN ISO 15841) e.g. 17 x 25 = .017" x .025"

1 mm	= 0.0394 inch
1 cm	= 10 mm = 0.3937 inch
1 m	= 1000 mm = 1.0936 yards
1 km	= 0.6214 miles
1 inch	= 25.40 mm
1 foot	= 12 inches = 30.48 cm
1 yard	= 3 feet = 91.44 cm
1 mile	= 1760 yards = 1.6093 km

1 mg	= 0.0154 grain
1 g	= 1000 mg = 0.0353 oz
1 kg	= 1000 g = 2.2046 lb
1 ounce	= 28.35 g
1 pound	= 16 ounces = 0.4536 kg

1 MPa	= 1 N/mm ²
1 N/mm ²	= 145.14 p. s. i.
1 p. s. i.	= 0.00689 N/mm ²

1 l	= 1000 ml = 2.11 pints
	= 1.06 quarts = 0.26 gallons

FDI System																											
	18	17	16	15	14	13	12	11		21	22	23	24	25	26	27	28										
R	8	7	6	5	4	3	2	1		1	2	3	4	5	6	7	8										
	8	7	6	5	4	3	2	1		1	2	3	4	5	6	7	8										
L	48	47	46	45	44	43	42	41		31	32	33	34	35	36	37	38										

mm ∅	x 10 ⁻³ inch	steel wire 1 kg	mm ∅	x 10 ⁻³ inch	steel wire 1 kg
0.10	4	16120 m*	0.60	24	448 m*
0.15	6	7150 m*	0.70	28	329 m*
0.20	8	4030 m*	0.80	32	251 m*
0.22	9	3335 m*	0.90	36	198 m*
0.25	10	2590 m*	1.00	40	161 m*
0.28	11	2055 m*	1.10	44	133 m*
0.30	12	1793 m*	1.20	48	112 m*
0.35	14	1315 m*	1.30	51	95 m*
0.40	16	1010 m*	1.40	55	82 m*
0.45	18	795 m*	1.50	59	71 m*
0.50	20	645 m*			

*We accept no responsibility for the accuracy of this information

inch	mm
1/64	0.0156
1/32	0.0312
1/16	0.0625
1/8	0.1250
3/16	0.1875
1/4	0.2500
5/16	0.3125
3/8	0.3750
1/2	0.5000
5/8	0.6250
3/4	0.7500
1	1.0000

°C	°F	°C	°F
-40	-40	110	230
-30	-22	120	248
-20	-4	130	266
-10	+14	140	284
0	+32	150	302
10	50	200	392
20	68	300	572
30	86	400	752
40	104	500	932
50	122	600	1112
60	140	700	1292
70	158	800	1472
80	176	900	1652
90	194	1000	1832
100	212	2000	3632

System International (SI)		
Basic size	Basic unit	Symbol
Length	meter	m
Mass	kilogram	kg
Time	second	s
Electric current	ampere	A
Thermodynamic temperature	kelvin	K
Luminous intensity	candela	cd
Amount of substance	mole	mol

1 at	= 0.981 bar = 0.981 x 10 ⁵ Pascal (Pa)
1 kp	= 9.81 Newton (N)
1 PS	= 0.746 kW
0 °C	= 273 Kelvin (K)

mm	inch	mm	inch	mm	inch	mm	inch	mm	inch
0.0254	0.001	0.4318	0.017	0.8382	0.033	1.2192	0.048	1.6256	0.064
0.0508	0.002	0.4572	0.018	0.8636	0.034	1.2446	0.049	1.6510	0.065
0.0762	0.003	0.4826	0.019	0.8890	0.035	1.2700	0.050	1.6764	0.066
0.1016	0.004	0.5080	0.020	0.9144	0.036	1.2954	0.051	1.7018	0.067
0.1270	0.005	0.5334	0.021	0.9398	0.037	1.3208	0.052	1.7272	0.068
0.1524	0.006	0.5588	0.022	0.9652	0.038	1.3462	0.053	1.7526	0.069
0.1778	0.007	0.5842	0.023	0.9906	0.039	1.3716	0.054	1.7780	0.070
0.2032	0.008	0.6069	0.024	1.0000	0.0394	1.3970	0.055	1.8034	0.071
0.2286	0.009	0.6350	0.025	1.0160	0.040	1.4224	0.056	1.8288	0.072
0.2540	0.010	0.6604	0.026	1.0414	0.041	1.4478	0.057	1.8542	0.073
0.2794	0.011	0.6858	0.027	1.0668	0.042	1.4732	0.058	1.8796	0.074
0.3048	0.012	0.7112	0.028	1.0922	0.043	1.4986	0.059	1.9050	0.075
0.3302	0.013	0.7366	0.029	1.1176	0.044	1.5240	0.060	1.9304	0.076
0.3556	0.014	0.7620	0.030	1.1430	0.045	1.5494	0.061	1.9558	0.077
0.3810	0.015	0.7874	0.031	1.1634	0.046	1.5748	0.062	1.9812	0.078
0.4064	0.016	0.8128	0.032	1.1938	0.047	1.6002	0.063	2.0000	0.0787

Environment and sustainability.



Dentaurum has implemented an environmental management system for over 25 years, making it one of the first in the dental world.

Resources are finite. In a world in which production output and consumer demands are ever increasing, companies bear responsibility for a conscious and sustainable management of our environment. To leave a habitable planet behind for future generations, habits and production techniques need to change to create an economic society that is beneficial to natural resources instead of merely consuming them. Dentaurum has been certified acc. to DIN EN ISO 14001 and EMAS for many years and is very aware of this responsibility. In all production and business processes, focus is placed on conserving resources, reducing the burden on the environment as well as sourcing raw materials in an ethically correct manner, and this has been followed since 1989.



You can find the measures we have taken to protect the environment clearly listed on our home page under [Environmental protection at Dentaurum](#).

Or you can scan the QR code.

The timeline gives you an overview of our measures and the opportunity to learn more about the areas that interest you most.



Alphabetical index

1, 2, 3	Page
3-Prong pliers	379, 382, 389
3D printing	26–27
3D retainer	27, 221
3D separating medium	27, 359
A	Page
Abutments, tomas®	283–285
Accessories, tomas®	286
Activator / bonding technique, CONTEC go	180
Activator screw	316, 327
Activator tubes	263
activeblue cleaning unit	366, 420
activefluid cleaning liquid	367, 420
Adams clasp	231
Adams pliers	380–381, 388
Aderer pliers	379, 382, 389
Adhesive (bracket adhesive)	178–180
Adhesive residue remover	181
Adhesive, CONTEC go, CONTEC lc/lcr	178, 180
Adhesives	
• CONTEC go – chemically-curing	180
• CONTEC lc/lcr – light-curing	178–179
Aesthetic brackets	
• discovery® pearl – ceramic brackets	104–108
Alginate	298
Aligner	26
amda® – distalization appliance	
• amda® connector	252, 291
• amda® palatal arch	253, 292
• amda® starter set	251, 290
• amda® telescope	252, 291
Anchorage system, tomas®	
• Instruments and accessories	281–289
• Patient consultation material	293, 306–307
• tomas®-abutment	283–285
• tomas®-abutment kit	283
• tomas®-auxiliary kit	286
• tomas®-laboratory pin	285
• tomas®-pin	254, 279–280
• tomas®-RPE eyelet	283, 336
• tomas®-tool set	278
Angle pliers	
• Ribbon arch pliers	376, 387
• Wire bending pliers	378, 387
Appliance containers	301–302
Arch template	398
Arch turret	394–395
Archwire holder	226

Archwires, pre-formed	48, 192–195, 200–216
Argon fitting for laser	412
Arrow clasp former, arrow clasp bending pliers	383
Arrow to indicate the opening direction	340
Assistant 3000 spot welding unit	415
Auxiliary kit, tomas®	286
B	Page
Ball retainer clasp	232
Band container, cabinets	169–170
Band material	168
Band pliers	373
Band pusher	391
Band removing pliers	373, 386
Band seater	391
Band seating lug / Banding technique	135
Bands	
• Accessories	168–170
• dentaform® (softer band quality)	147–156
• dentaform® Snap (hard, stiff band quality)	142–146
• Standard (higher band shape)	157–167
Base model	26
Base unit	300
Beerendonk dental vernier	396
Begg Light Wire pliers, Premium-Line	376
Bending key	340
Bernklau ruler	396
Bertoni expansion screw	320–322
Bite jumping hinge acc. to Herbst	246–247
Bite-Bumper™	114
black & white, Orthocryl®	351
Bow divider	397
Bow screw	324
Bracket accessories	114–115
Bracket adhesives	178–180
Bracket container	115
Bracket mesh bases	114
Bracket placement instruments	390
Bracket positioning gauge	392
Bracket protection wax	115
Bracket removing pliers	49, 372–373, 386
Bracket storage box	115
Bracket tray	115
Brackets / Banding technique	
• Metal brackets, equilibrium® 2 18/22	93–94
Brackets / Bonding technique	
• Ceramic brackets, discovery® pearl 18/22	104–108

Alphabetical index

• Lingual brackets, discovery® delight	46–47
• Metal brackets, discovery® 18/22	68–73
• Metal brackets, discovery® smart 18/22	74–77, 108
• Metal brackets, equilibrium® 2 18/22	86–92
• Metal brackets, equilibrium® mini 18/22	81–85
• Nickel-free* brackets, topic 18/22	78–80
• Pure titanium brackets, equilibrium® ti 18/22	95–96
• Self-ligating brackets, discovery® sl 18/22	56
• Self-ligating brackets, discovery® sl 2.0 22	57
Buccal tubes / Banding technique	
• Ortho	126, 128, 131–132
• Ortho-Cast	128, 130, 131, 132–133
• Ortho-Cast M-Series	127, 129–130, 131, 132
Buccal tubes / Banding technique	
• double combination	131–132
• double rectangular	129–130
• double round	126
• single rectangular	127–128
• single round	126
• triple combination	132–133
Buccal tubes / Bonding technique	
• Ortho	120, 122
• Ortho-Cast	122, 123, 124
• Ortho-Cast M-Series	121, 123
• Ortho-Cast M-Series mini	120
Buccal tubes / Bonding technique	
• double combination	124
• double rectangular	123–124
• double rectangular / titanium	125
• single rectangular	120–123
• single round	120
Bur stand	400
Burs	360–362
Button	
• Banding technique	135
• Bonding technique	114
• with chain / Bonding technique	114
• with hook / Bonding technique	114
Buttons / Bonding technique	114
C	Page
CAD/CAM retainer	27, 221
Cap removal pliers	372
Cap removal pliers, replacement blades	372
Carbide burs	181
Carbon insert	417
Cephalometric tracing accessories	398

Ceramic brackets	
• discovery® pearl	104–108
Chamois leather	363
Chemically-curing adhesive, CONTEC go	180
Clamping electrode	416
Cleaning liquid, activefluid	367, 420
Cleaning pins, activeblue	366, 420
Cleaning unit, activeblue	366, 420
Cloth polishing brush	363
Coil springs	223–224
Color concentrate	352
Color pencils	400
Compression gauge	395
Compression spring	224
Container for brackets	115
Container for intra-oral elastics	265
CONTEC go – adhesive	
• CONTEC go activator	180
• CONTEC go adhesive	180
• CONTEC go etching gel	180
• CONTEC go set	180
CONTEC lc – Adhesive	
• CONTEC lc bracket adhesive	178
• CONTEC lc etching gel	179
• CONTEC lc/lcr primer	179
• CONTEC lc set	178
CONTEC lcr – Adhesive	
• CONTEC lcr etching gel	179
• CONTEC lc/lcr primer	179
• CONTEC lcr retainer adhesive	179
• CONTEC lcr set	178
Cord for appliance containers	302
Corundum disc / Model trimmer	422, 423
Coupling elements, tomas®	283–285
Crimping pliers	225, 245, 250, 378
Crimping pliers, Herbst TS/SUS	245, 250, 378
Crimping tweezers	49, 390
Cross tubes	226, 288
Crozat wire	196, 227
Cuspid bands	160, 166–167
Cuspid retraction springs	224
D	Page
De La Rosa contouring pliers	377
Debonding instrument / discovery® pearl	109, 390
Decal, orthodontics	354–357
Decals	354–357
Demonstration model	
• discovery®	304
• discovery® delight	50, 305

Alphabetical index

• discovery® pearl	109, 305
• discovery® sl	58, 306
• discovery® sl 2.0	58, 306
• discovery® smart	305
• equilibrium® 2	304
• equilibrium® mini	304
• tomas® / amda®	253, 293, 306
Demonstration models	27, 50, 58, 109, 253, 293, 303–307
Dentaflex® wires	194, 196, 198, 209–210, 217, 219
Dentaflux® universal flux	418
dentaform® – bands	
• dentaform® bands / Assortments	147–148
• dentaform® bands / Single packs	149–156
dentaform® Snap – bands	
• dentaform® Snap bands / Assortments	142
• dentaform® Snap bands / Single packs	143–146
Dental arch model	26
Dental mirror with reference lines, lingual	49, 393
Dental verniers	341, 396–398
Dentalastics®	
• Chains, extra-oral	262
• Ligatures	109, 267
• Personal	109, 267
• Plastic ligatures	266
• Separators	266
Dentaurum Digital	9, 26–27, 221
Dentoform unit	303
Dermatril® disposable rubber gloves	367
desktop Compact Laser	410
Diagnostic instruments	396–398
Diamond disc / Model trimmer	422
Digital orthodontics	26–40
Disco glitter	353
discovery® – premium brackets	
• discovery® premium brackets / Assortments	68–72
• discovery® premium brackets / Single packs	73
discovery® delight – lingual brackets	
• discovery® delight lingual brackets / Assortments	46
• discovery® delight demonstration model	50, 305
• discovery® delight lingual brackets	47
• discovery® delight lingual tubes	47
discovery® pearl – ceramic brackets	
• discovery® pearl ceramic brackets / Accessories	109

• discovery® pearl ceramic brackets / Assortments	104–105
• discovery® pearl ceramic brackets / Demonstration model	109, 305
• discovery® pearl ceramic brackets / Single packs	106–107
• discovery® pearl & smart combined Assortments	108
discovery® pearl debonding instrument	109, 390
discovery® sl – self-ligating brackets	
• discovery® sl brackets 18/22 (metal)	56
• discovery® sl demonstration model	58, 306
• discovery® sl 2.0 brackets 22 (metal)	57
• discovery® sl 2.0 demonstration model	58, 306
• sl-instruments	58, 390
discovery® smart – metal brackets	
• discovery® smart metal brackets / Assortments	74–77, 108
• discovery® smart metal brackets / Demonstration model	305
Dispensing pipette for color concentrates	358
Disposable rubber gloves, Dermatril®	367
Distal end bender	394
Distal end cutter	374, 386
Distalization appliance	251–253, 290–292
Distalization appliance, amda®	
• Appliance and accessories	252–253, 291–292
• Starter set	251, 290
Dosing bottle	298
Double hooks / Banding technique	135
Double hooks / Bonding technique	114
Download STL data	27
Drill, tomas®	281
Dual wheel model trimmer, discs	422, 423
E	Page
Edelweiß	
• Acrylic dry polishing paste	364
• Acrylic polishing liquid	364
• Acrylic polishing stick	364
Elastic placers	264
Elastics	
• extra-oral, latex	262
• intra-oral, latex	264, 265
• intra-oral, non-latex	265
Electrically adjustable stand	412
Electrodes	416–417
Emery mandrel, emery cloth	364
EQ-Line pliers	386–389

Alphabetical index

equilibrium® 2 – metal brackets	
• equilibrium® 2 metal brackets / Assortments	86–91
• equilibrium® 2 metal brackets / Single packs	92–94
equilibrium® mini – mini brackets	
• equilibrium® mini mini brackets / Assortments	81–84
• equilibrium® mini mini brackets / Single packs	85
equilibrium® ti – pure titanium brackets	
• equilibrium® ti pure titanium brackets / Assortments	95–96
Equire wires	193, 194, 208, 212–213
Etching gel / bonding technique, CONTEC go, CONTEC lc/lcr	179, 180
Expansion screw key	340
Expansion screw sample kit	341
Expansion screws	
• hyrax®	30–31, 337–338
• hyrax® click	28–29, 335–336
• hyrax® neo	31–32, 338–339
• Neo-Bertoni 3D	321
• Neo-Bertoni three-dimensional screws	320–321
• Neo-bite jumping appliance acc. to F.M. Sander	325
• Retraction screws	319–320
• RPE screws	28–33, 334–339
• Special expansion screws	317–329
• Standard expansion screws	314–316
• Three-dimensional screws	322–323
• Titan hyrax®	33, 334
• Variety click expansion screws	34–35, 330–331
• Variety expansion screws	36–37, 332–333
Extra-oral plastic chain	262
Extra-oral, elastics, latex	262
Eyelet clasp	232
F	Page
Facebow	260–261
Facebow accessories	
• Activator tubes	263
• Stop screws	263
• Stop tubes	263
Facebow bending pliers	378
Facebows	
• nickel-free	260
• standard	260–261
• with stop loops	261

FaceFormer	268
Fan-type expansion screw	324
File	418
Fixator	365
FKO-Split-Post Fixator	365
Flexaponal® special wax	298
Flush Cutter	374–375
Forward thrust double plate system	325, 326
Framework for cabinet	169–170
Funnel	358
G	Page
Gauges	395
Gauges, compression and tension	395
Gold retainer wire	197, 220
Goshgarian	234–235
Guide pins	233
H	Page
Hammerhead pliers, Premium-Line	379
Hand electrodes, holder	416
Heat protection	418
Heating plate	365, 421
Herbst accessories	248–250
Herbst bite jumping hinges	246–247
Hexagon socket key	245, 250, 253, 263, 287, 292
Hexagon socket screws / Herbst	248
Hilgers 18 – bracket system	
• Metal brackets, discovery® 18	71
Hooks	
• Double hooks / Banding technique	135
• Double hooks / Bonding technique	114
• With seating lug / Banding technique	135
How pliers	373–374, 386
hyrax® screws	30–31, 334, 337–338
hyrax® click screws	28–29, 335–336
I	Page
Impression material	
• Orthalgene® Duo	298
Impression tray cleaner	
• Orthalgene® Ex	298
Isolation for printed models	27, 359
J	Page
Jarabak Light Wire pliers	377, 387
Jaw fracture splint	225
Johnson contouring pliers	373
K	Page
Keys	
• Expansion screw key	340

Alphabetical index

• Hexagon socket key	245, 250, 253, 263, 287, 292
• Safety key / hyrax® / Variety	340
• Universal bending key	340
Kobayashi ligatures	222
Korkhaus	
• Bow divider	397
• Orthometer kit	397
L	Page
Labial arch	231
Labial shield	262
Laboratory analog, tomas®	285
Laboratory pins, tomas®	285
Laser desktop Compact	410
Laser Welder SL10	411
Laser welding wires	413–414
Ligator	393
Ligature adapter and tightener	393
Ligature cutter	375
Ligature rings	109, 267
Ligature wire	222
Light Wire pliers	375–376, 377
Light-curing adhesive, CONTEC lc/lcr	178–179
Lingual arch bending pliers	378
Lingual arches	48, 236–237
Lingual bar	233
Lingual bending pliers	49, 376
Lingual brackets, discovery® delight	46–47
Lingual buttons	135
Lingual hooks	135
Lingual sheaths	134
Lingual system	
• Accessories	50
• discovery® delight – lingual brackets	47
• discovery® delight – lingual brackets / Assortments	46
• Instruments	49, 373, 376, 386, 390, 393
• Lingual arches	48
• Lingual tubes	47
• Patient consultation material	50, 305
Lingual system / Instruments	
• Bracket removing pliers	49, 373
• Crimping tweezers	49, 390
• Dental mirror with reference lines	49, 393
• Lingual bending pliers	49, 376
• Mosquito forceps	49, 393
• Positioning tool	49, 393
Lingual system Orthorama®	236–237

Lingual tubes, discovery® delight	47
LinguForm arch turret	395
Lip bumper	262
Loop closing pliers	
• Nance	376, 387
• Tweed	376, 377, 387
Loop closing pliers, Nance	376, 387
Loop forming pliers	379, 380, 388
M	Page
Maintenance and protection oil for pliers and instruments	399
Mandibular bow screw	324
Mandrel	364
Marking pencil	400
Mathieu needle holder	392
Matrix band	168
McLaughlin-Bennett-Trevisi** 18/22 – bracket system	
• Ceramic brackets, discovery® pearl 18/22	105, 107, 108
• Combined assortment, discovery® pearl & smart 22	108
• Metal brackets, discovery® smart 18/22	76–77, 108
• Metal brackets, equilibrium® 2 18/22	88–89, 93
• Mini brackets, equilibrium® mini 18/22	83–85
• Pure titanium brackets, equilibrium® ti 18/22	96
• Self-ligating brackets, discovery® sl 2.0 22	57
Measuring beaker	358
Measuring gauge	395
Measuring grid	396
Measuring instruments	341, 396–398
Mershon band pusher	391
Mesh bases (anatomically contoured) / Bonding technique	114
Metal brackets / Bonding technique	
• equilibrium® 2 18/22 – metal brackets	93–94
Metal brackets / Bonding technique	
• discovery® 18/22 – premium brackets	68–73
• discovery® delight – lingual brackets	46–47
• discovery® sl 18/22 – self-ligating brackets	56
• discovery® sl 2.0 22 – self-ligating brackets	57
• discovery® smart 18/22	74–77, 108
• equilibrium® 2 18/22 – metal brackets	86–92
• equilibrium® mini 18/22 – mini brackets	81–85
• equilibrium® ti 18/22 – pure titanium brackets	95–96
• topic 18/22 – nickel-free* brackets	78–80
Metal/Plastic teeth	303

* see Materials list Orthodontics

Alphabetical index

Mini brackets, equilibrium® mini	81–85
Mini screw	254, 279–280
Mixing bowls	299, 358
Model bases, hinges	300
Model gloss	299
Model isolation for printed models	27, 359
Model printing	26
Model separating medium for printed models	27, 359
Model tray	300
Model trimmer, discs	422, 423
Modeling tips for Thermomat	424
Modeling wax	359
Molar band removing pliers	373, 386
Molar band seater	391
Mosquito forceps	49, 392–393
Müller guide pins	233
Multi-action lingual arches	237
Münchner model dental vernier	396
N	Page
Nance pliers	376, 387
Needle holder	392
Neo-Bertoni 3D	321
Neo-Bertoni three-dimensional screws	320–321
Neo-BJA acc. to F.M. Sander	325
Nickel-free* brackets	78–80, 95–96
Nickel-free therapy	110–113
Noninium® wires	199, 230
O	Page
o-atlas II	428–429
Oliver-Jones band removing pliers	373, 386
OnyxCeph ^{3TM}	38–40
Oral screen	268
Order platform	9
Orthalgemat®	
• Duo impression material	298
• Ex impression tray cleaner	298
Ortho Brush	308
Ortho buccal tubes / Banding technique	126, 128, 131–132
Ortho buccal tubes / Bonding technique	120, 122
Ortho-Cast buccal tubes / Banding technique	128, 130, 131, 132–133
Ortho-Cast buccal tubes / Bonding technique	122, 123, 124
Ortho-Cast M-Series buccal tubes / Banding technique	127, 129–130, 131, 132
Ortho-Cast M-Series buccal tubes / Bonding technique	121, 123

Ortho-Cast M-Series cap removal pliers	372
Ortho-Cast M-Series cap removal pliers, replacement blades	372
Ortho-Cast M-Series mini buccal tubes	120
Orthocryl® – acrylics	
• Accessories	358–367
• Assortments	350
• black & white	351
• Color concentrates	352
• Decals	354–357
• Disco glitter	353
• Neon	350, 351
• Single packs	351–352
Orthocryl® LC – light-curing acrylic	
• Accessories	349
• Single packs	349
• Starter sets	348
Orthodontic acrylics	348–353
Orthodontic acrylics	
• Accessories	358–367
• Disco glitter	353
• Decals	354–357
• Orthocryl® – acrylic	350–353
• Orthocryl® LC – light-curing acrylic	348–349
Orthodontic appliance cleaner	366, 420
Orthodontic model tray	300
Orthodontic models	26
Orthodontic planning software	38–40
Orthodontic treatment planning	38–40
Orthorama®	
• Lingual arches	236–237
• Palatal arches	236
• Sectional arches	237
• Transfer system	237
P	Page
Palatal arches, palatal bars	234–235, 236
Palatal bite raiser / Bonding technique	114
Palatal expansion screw	330–339
Palatal tube, palatal sheath	134
Pilot drill, tomas®	281
Pins, tomas®	254, 279–280
Piston spring screw	328
Planning software	38–40
Plaster	299
Plastic chains	50, 109, 266
Plastic inserts / Pliers	372, 386
Pliers	49, 225, 245, 250, 289, 372–383, 386–389

Alphabetical index

Pointed beak pliers	382, 389
Polishing brushes	362–363
Polyclav® polymerizing pressure vessel	365, 421
Positioning tool, lingual	49, 393
Pre-drill, tomas®	281
Pre-formed arches	48, 192–195, 200–216
Pre-formed bands	142–156
Pre-welded parts / Banding technique	
• Band seating lug	135
• Buttons	135
• Double hooks	135
• Hooks with seating lug	135
• Lingual / Palatal sheaths	134
Premium-Line pliers	372–383
Pressure vessel, hose	365, 421
prime4me® Aligner	26
prime4me® model print	26
prime4me® RETAIN3R	27
Primer for CONTEC Ic/Icr	179
Printing service	9, 26–27
Progenia screw	327
Protection oil / Pliers, instruments	399
Protection wax	115
Protective gloves, Dermatril®	367
Q	Page
Quad Helix	235
R	Page
Rapidur® dental stone	299
Reamer	329
Rectangular tubes	226
remaloy® wires	194, 196, 213, 218, 227–228
	194, 196, 198, 211–212, 217–218, 219, 227, 229–230
rematitan® flat retainer wire	196, 220
rematitan®	
• Expansion screws	316, 334
• LITE compression springs	224
• LITE tension springs	223
• LITE wires	48, 192–193, 198, 200–206, 219
• LITE White wires	195, 215
• SPECIAL wires	195, 197, 214, 218
• Welding wires	334, 413–414
rematitan® buccal tubes	125

rematitan® LITE Spee®	192, 204
retraction arch	
	192–193, 195, 197, 198, 200–206, 214–215, 218–219
rematitan® wires	
Retainer 3D	27, 221
Retainer, milled	31, 221
Retainer, pre-made	221
Retainer, titanium, gold, stainless steel	196, 197, 199, 220, 221
Retention bands / Transfer jigs	237
Retraction arches	192, 204
Retraction screw	319–320
Ribbon arch pliers, Angle / Tweed	376, 387
Ricketts® 18 – bracket system	
• Metal brackets, discovery® 18	70
• Metal brackets, equilibrium® 2 18	90, 92
Roth 18/22 – bracket system	
• Ceramic brackets, discovery® pearl 18/22	104, 106, 108
• Metal brackets, discovery® 18/22	68–69, 73
• Metal brackets, discovery® smart 18/22	74–75, 108
• Metal brackets, equilibrium® 2 18/22	86–87, 92, 93
• Mini brackets, equilibrium® mini	81–82
• Nickel-free* brackets, topic 18/22	78–80
• Pure titanium brackets, equilibrium® ti 18/22	95
• Self-ligating brackets, discovery® sl 18/22	56
RPE eyelet, tomas®	283, 336
RPE screws	28–33, 334–339
Ruhland loop forming pliers	380, 388
Rulers	341, 396
S	Page
Sabbagh Universal Spring – SUS3	244–245
Safety key / hyrax®	340
Sander screw	325
Saw	362
Saw blades	362
Scaler	391
Schaneng forward thrust double plate system	326
Schmuth measuring grid	396
Scissors	391
Screwdriver	282
Sectional arches	237
Sectional screws	317–318
Self-ligating bracket system	
• Brackets 18/22 (metal), discovery® sl	56

Alphabetical index

• Brackets 22 (metal), discovery® sl 2.0	57
• discovery® sl demonstration model	58, 306
• discovery® sl 2.0 demonstration model	58, 306
• sl-instruments	58, 390
Separating medium	359
Separating medium for printed models	27, 359
Separating pliers, Premium-Line	373
Separating strips	391
Separators	266
Side cutter	374, 375, 386
Side cutter, heavy gauge	389
Silicone O-rings / Flush Cutter	375
Silicone polisher	364
Silver solders	419
sl-instruments	58, 390
Sliding tubes	225
Slotted screws / Herbst	248
Snap hooks	302
Soldering carbon	416
Soldering electrode	417
Solders	419
Spacer rings	245, 250
Special bracket debonding pliers	372
Special expansion screws	317–329
Spot welding unit	415
Spray nozzles, spray bottles	358
Stainless steel expansion screw	314–315, 316–325, 326, 327, 330–333, 335–339
Stainless steel tubes	234
Stand for pliers	399
Stand, electrically adjustable	412
Standard – bands	
• Standard bands / Assortments	157–160
• Standard bands / Single packs	161–167
Standard Edgewise 18/22 – bracket system	
• Metal brackets, discovery® 18/22	71–72
• Metal brackets, equilibrium® 2 18/22	90–91, 94
Standard expansion screws	314–316
Sticky wax	359
STL download	27
Stop screws / Facebow, tomas®	263, 287
Stop tubes / Facebow	263
Stop tubes / Stop arches	226
Stopper for cabinet / Bands	170
Storage box for brackets	115
Strengtheners, flat, rolled	234
Stripper	391

Study models	27, 50, 58, 109, 253, 293, 303–307
Suction apparatus with charcoal filter	366, 421
SUPER Elasto-Force plastic chains	50, 109, 266
SUS ³ – Sabbagh Universal Spring	244–245
T	Page
TAD, temporary anchorage device	254, 279–280
Teeth	303
Templates for diagnostic purposes	398
Tensic® White wires	195, 216
Tensic® wires	193, 207
Thermo arches	193, 195, 207–208, 216
Thermo wax, blocking-out wax	359
Three-dimensional screw	320–323
Titanium	
• Brackets	95–96
• Buccal tubes	125
• Buttons with chain	114
• Expansion screws	33, 316, 334
• hyrax®	33, 334
• Retainer	221
• Springs	223, 224
• Wires	48, 192–193, 195, 197, 198, 199, 200–208, 214–216, 218–219, 220, 221, 413–414
Titanium buccal tubes	125
tomas® anchorage system	
• Instruments and accessories	281–289
• Patient consultation material	293, 306–307
• tomas®-abutment	283–285
• tomas®-abutment kit	283
• tomas®-auxiliary kit	286
• tomas®-laboratory pin	285
• tomas®-pin	254, 279–280
• tomas®-RPE eyelet	283, 336
• tomas®-tool set	278
Tooth stripper	391
Tooth wax	115
Toothbrushes	308
topic – nickel-free* brackets	
• topic nickel-free* brackets / Assortments	78–79
• topic nickel-free* brackets / Single packs	80
Torch, mini	423
Torque ratchet, tomas®	282

Alphabetical index

Tracing foils	398
Tracing-protractor templates	398
Transfer jigs / Retention bands	237
Trapezoidal expansion screw	316
Trays, cabinets / Bands	169–170
Treatment planning	38–40
Triangular clasp	232
Trimmer, discs	422, 423
Tubes / Activator	263
Tubes / Banding technique	126–133
Tubes / Bonding technique	120–125
Tubes, stainless steel	234
Tweed loop forming pliers	376, 377, 387
Tweezers	390
Twisted wires	194, 196, 197, 198, 199, 209–210, 217, 219, 220
Typodont, device	303
U	Page
Ulmer Modell, oral screen	268
Ultrasonic cleaning solution	299
Universal bending key	340
Universal flux	418
Universal wire bending pliers	381
Uprighting spring	286
V	Page
Variety click expansion screws	34–35, 330–331
Variety expansion screws	36–37, 332–333
W	Page
Wax bases	303
Wax pencil	400
Waxes	115, 298, 303, 359
Wedge gauge	341, 396
Weingart pliers	376, 387
Welding electrode	417
Welding units	410–411, 415
White arches	195, 215–216
White gold solder	419
White oil / Protection oil for pliers and instruments	399
Wire bending pliers	381, 388
Wire ligatures	222
Wire, rolled flat	234
Wires	48, 192–220, 222, 227–230, 413–414
Wires and pliers, stand	399
Y	Page
Young loop bending pliers	382, 389

Numerical index

REF	Page	REF	Page	REF	Page	REF	Page
000-001-00	376	004-180-00	394	026-355-00	391	045-016-00	388
000-030-00	392	004-182-00	395	027-349-00	391	045-017-00	388
000-031-00	392	004-184-00	395	028-353-00	397	045-018-00	389
000-035-00	392	004-186-00	394	029-360-01	397	045-019-00	388
000-730-00	392	004-187-00	395	029-360-02	397	045-020-00	389
000-731-00	392	004-190-00	395	029-361-02	397	045-021-00	386
000-732-00	49, 393	004-266-00	374	029-362-00	397	045-022-00	386
001-650-00	380	004-346-00	372	029-363-00	397	050-500-00	300
003-001-00	49, 376	004-347-00	49, 373	029-364-00	397	050-501-00	300
003-053-00	380	004-350-00	372	030-390-00	392	050-502-00	300
003-094-00	381	004-368-00	398	030-395-02	397	050-503-01	300
003-111-00	374	004-369-00	398	030-395-03	397	050-801-00	303
003-112-00	373	004-373-00	398	030-396-02	397	050-802-00	303
003-114-00	373	004-455-00	382	033-452-00	381	050-803-01	303
003-120-00	376	004-515-00	380	035-455-00	382	055-000-00	399
003-121-00	376	007-024-00	391	037-356-00	374, 375	055-100-00	399
003-125-00	377	008-025-00	391	038-356-00	372, 386	055-101-01	399
003-130-00	375	012-074-00	382	038-357-00	377	055-104-00	226
003-139-00	378	012-350-00	377	038-805-10	372	055-105-00	399
003-180-01	378	012-351-00	377	039-357-00	373, 386	055-201-00	399
003-185-00	377	012-520-00	383	040-711-00	395	070-000-00	365, 421
003-200-00	382	012-521-00	383	040-712-00	395	070-010-00	365, 421
003-201-00	382	013-522-00	381	040-713-00	395	070-015-00	365, 421
003-202-00	382	014-151-00	375	042-750-00	396	070-900-00	365, 421
003-205-00	379	014-152-00	375	042-751-00	142 - 167, 396	072-001-00	366, 421
003-233-00	380	015-160-00	373	043-700-00	400	072-001-40	366, 421
003-290-00	381	016-153-00	375	043-720-00	341, 396	072-002-00	366, 421
003-349-00	372, 386	018-001-01	391	043-800-00	400	072-004-00	365
003-350-00	379	018-002-01	391	043-900-00	400	074-000-00	268
003-351-00	379	019-000-00	394	044-177-00	389	074-001-00	268
003-355-00	379	019-001-00	109, 390	044-730-01	396	074-101-01	268
003-366-00	396	021-273-00	394	044-731-00	396	075-100-00	398
003-442-00	376	022-274-00	394	045-000-00	387	075-101-00	398
003-514-00	380	023-276-00	58, 390	045-001-00	387	078-000-00	415
003-515-00	380	023-277-00	58, 390	045-002-00	387	080-116-00	416
003-516-00	381	023-278-00	49, 393	045-003-00	387	081-101-00	416
003-700-00	374	023-279-00	49, 393	045-004-00	373, 386	081-204-00	416
003-701-00	374	023-279-10	49, 393	045-005-00	372, 386	081-205-00	416
003-702-00	374	023-292-00	393	045-006-00	386	081-206-00	416
003-705-00	374, 375	023-293-00	393	045-007-00	389	081-501-00	416
003-706-00	375	023-294-00	393	045-008-00	389	081-601-00	416
003-710-00	245, 250, 378	023-295-00	393	045-009-00	387	081-701-00	416
003-711-00	225, 289, 378	024-300-01	391	045-010-00	387	083-300-00	418
003-805-00	372	025-276-00	390	045-011-00	389	084-100-00	418
004-095-00	381	025-277-00	390	045-012-00	388	085-000-00	417
004-109-00	377	025-279-00	390	045-013-00	388	085-100-00	417
004-130-00	376	025-280-00	49, 390	045-014-00	388	085-200-00	417
004-139-00	378	026-347-00	373, 386	045-015-00	388	085-300-00	417

Numerical index

REF	Page	REF	Page	REF	Page	REF	Page
085-400-00	417	123-702-00	181	160-100-34	357	161-609-00	353
086-000-00	417	123-703-00	181	160-100-39	354	161-611-00	353
086-100-00	417	124-002-00	298	160-100-40	354	161-620-00	352
086-200-00	417	124-003-00	298	160-100-41	354	161-622-00	352
086-300-00	417	124-202-00	298	160-100-42	355	161-912-00	350
089-430-00	400	124-305-01	115	160-100-43	355	162-100-00	350, 358
089-431-00	400	138-623-00	364	160-100-44	355	162-101-00	358
089-602-00	424	139-500-00	364	160-100-45	354, 357	162-200-00	350, 358
089-603-00	424	140-015-00	364	160-100-46	355	162-501-00	358
089-604-00	424	141-700-00	362	160-112-00	350, 351	162-502-00	358
089-605-00	424	141-700-10	362	160-212-00	350, 351	162-503-00	358
089-606-00	424	141-710-00	362	160-300-00	351	162-600-00	350, 358
089-612-00	424	141-711-00	363	160-301-00	351	162-751-00	350, 358
089-620-00	424	141-712-00	363	160-400-00	348	162-752-00	350, 358
090-404-00	412	141-713-00	363	160-400-20	348	162-800-00	350, 359
090-513-10	412	141-713-10	363	160-401-00	348, 349	162-801-00	27, 359
090-574-00	412	150-000-00	298	160-402-00	348, 349	162-907-00	367
090-578-50	410	150-010-00	298	160-404-00	348, 349	162-908-00	367
090-594-00	411	150-020-00	366, 367, 420	160-410-00	349	162-909-00	367
091-100-00	422, 423	150-501-20	298	160-420-00	348, 349	163-107-00	180
091-101-00	423	150-509-20	298	160-421-00	348, 349	163-110-00	178
091-102-00	423	150-520-20	298	160-422-00	349	163-111-00	178, 179
091-105-00	423	160-011-00	351	160-430-00	349	163-113-10	178, 179
094-610-00	299	160-012-00	351	160-501-00	350	163-114-00	180
094-611-00	299	160-013-00	351	160-502-00	350	163-115-00	178, 179, 180
094-612-00	299	160-014-00	351	160-600-00	350	163-121-00	180
097-400-01	423	160-100-00	354	160-700-00	350	163-130-00	178
098-000-00	366, 420	160-100-01	354, 355	161-100-00	350, 352	163-131-00	178, 179
098-001-00	366, 420	160-100-02	354, 355	161-127-00	350, 352	163-300-00	308
119-150-00	359	160-100-04	354, 355	161-128-00	350, 352	163-602-01	301
120-170-00	359	160-100-05	354, 356	161-129-00	350, 352	163-610-01	301
122-300-00	359	160-100-06	354, 356	161-130-00	350, 352	163-611-01	301
123-580-00	361	160-100-07	354, 356	161-131-00	352	163-612-01	301
123-581-00	361	160-100-09	354, 356	161-132-00	352	163-701-01	301
123-582-00	361	160-100-10	354, 356	161-133-00	352	163-702-01	301
123-583-00	361	160-100-11	354, 356	161-134-00	352	163-703-01	301
123-584-00	361	160-100-12	354, 356	161-135-00	350, 351	163-707-01	301
123-585-00	361	160-100-13	354, 356	161-136-00	350, 351	163-710-01	301
123-594-00	361	160-100-14	354, 356	161-137-00	350, 351	163-711-01	301
123-596-00	361	160-100-15	354, 356	161-138-00	350, 351	163-712-01	301
123-597-00	362	160-100-16	354, 357	161-139-00	350, 351	163-713-01	301
123-598-00	362	160-100-18	354, 357	161-150-00	352	163-801-00	302
123-603-00	181	160-100-19	354, 357	161-300-00	352	163-801-01	301
123-603-30	181	160-100-27	354, 357	161-350-00	352	163-802-01	301
123-604-00	181	160-100-28	357	161-605-00	353	163-805-01	302
123-604-30	181	160-100-29	354, 357	161-606-00	353	163-807-00	302
123-700-00	181	160-100-30	357	161-607-00	353	164-450-00	299
123-701-00	181	160-100-31	354, 357	161-608-00	353	164-650-00	299

Numerical index

REF	Page	REF	Page	REF	Page	REF	Page
164-651-00	299	302-025-15	283	380-600-50	419	480-619-00	225
166-450-00	299	302-025-21	284	380-604-50	419	480-620-00	225
180-201-00	360	302-025-22	283, 284	380-604-60	419	480-621-00	225
180-202-00	360	302-025-30	284	380-704-50	419	480-621-02	225
180-304-00	360	302-025-51	284	385-624-00	198, 230	480-803-00	225
180-305-00	360	302-026-00	283, 285	385-630-00	198, 230	480-809-00	225
180-308-00	360	302-027-00	253, 283, 285, 292	402-253-00	225	480-914-00	263
180-702-00	362	302-028-01	283, 285	410-000-01	38	480-918-02	263
180-703-00	362	302-029-01	283, 285	410-000-02	38	480-924-02	263
181-100-01	363	302-030-01	283, 285	410-000-03	38	500-025-00	198, 222
181-200-00	363	302-050-00	283	410-000-04	38	500-030-00	198, 222
181-201-00	363	302-080-00	278, 282	410-000-05	38	500-040-00	198, 222
190-100-00	364	302-103-00	278, 281	410-000-06	39	501-025-00	198, 222
190-200-00	364	302-106-00	279	410-000-07	39	501-050-00	198, 222
302-001-20	286	302-106-10	280	410-000-08	40	503-050-00	198, 222
302-004-10	278, 282	302-106-30	280	410-000-09	40	513-070-00	198, 229
302-004-14	282	302-108-00	279	410-000-10	40	513-080-00	198, 229
302-004-16	282	302-108-10	280	410-000-11	39	513-090-00	198, 229
302-004-20	278, 281, 282	302-108-30	280	410-000-12	40	513-100-00	198, 229
302-004-30	278, 281	302-110-00	279	410-000-13	39	513-120-00	198, 229
302-004-41	278, 282	302-110-10	280	410-000-14	39	513-150-00	198, 229
302-004-50	278, 282	302-110-30	280	410-000-15	40	514-070-00	198, 229
302-004-70	278, 281, 282	302-151-00	278	410-001-00	38	514-080-00	198, 229
302-004-81	282	302-156-00	278, 282	410-002-00	38	514-090-00	198, 229
302-007-20	293, 307	302-206-00	254, 279	410-004-00	38	514-100-00	198, 229
302-007-30	293, 307	302-206-10	254, 280	410-005-00	38	520-062-10	199, 230
302-009-00	286	302-206-30	254, 280	430-110-10	26	520-070-10	199, 230
302-009-10	287	302-208-00	254, 279	430-110-30	26	520-072-10	199, 230
302-012-00	286	302-208-10	254, 280	430-120-30	27, 221	520-080-10	199, 230
302-012-10	286	302-208-30	254, 280	430-120-40	27, 221	520-082-10	199, 230
302-012-20	286	302-210-00	254, 279	430-120-90	27	520-090-10	199, 230
302-013-01	252, 286, 287, 291	302-210-10	254, 280	430-140-20	26	520-092-10	199, 230
302-014-18	288	302-210-30	254, 280	430-141-21	26	520-120-10	199, 230
302-014-22	288	302-251-00	278	430-142-21	26	521-040-00	198, 229
302-015-00	286, 288	302-602-00	283, 336	480-000-00	226	523-050-00	198, 229
302-019-00	286, 288	308-513-00	198, 230	480-001-00	234	523-060-00	198, 229
302-020-00	286, 289	308-515-00	198, 230	480-003-00	234	523-070-00	198, 229
302-020-10	286, 289	308-518-00	198, 230	480-005-00	234	523-080-00	198, 229
302-021-00	286, 289	308-520-00	198, 230	480-007-00	234	523-090-00	198, 229
302-021-10	286, 289	308-713-00	198, 230	480-009-00	234	523-100-00	198, 229
302-022-18	287	308-715-00	198, 230	480-011-00	234	523-110-00	198, 229
302-022-22	287	308-718-00	198, 230	480-013-00	234	523-120-00	198, 229
302-023-10	283, 286, 289	311-306-30	233	480-015-00	234	523-150-00	198, 229
302-024-00	286, 287	311-307-30	233	480-017-00	234	524-060-00	198, 229
302-025-00	283	311-308-30	233	480-022-00	226	524-070-00	198, 229
302-025-06	283, 284	312-106-00	234	480-100-00	226	524-080-00	198, 229
302-025-11	283, 285	312-107-00	234	480-617-01	226	524-090-00	198, 229
302-025-12	283, 285	380-070-00	418	480-618-01	226	524-100-00	198, 229

Numerical index

REF	Page	REF	Page	REF	Page	REF	Page
524-110-00	198, 229	535-051-00	196, 217	600-302-60	315	602-808-30	337
525-040-00	196, 227	537-310-00	196, 217	600-303-10	315	602-808-32	30, 337
527-004-00	196, 227	537-311-00	196, 217	600-303-30	315	602-809-10	333
527-005-00	196, 227	537-313-00	196, 217	600-303-60	315	602-809-12	37, 333
527-006-00	196, 227	537-314-00	196, 217	600-310-30	316	602-809-30	333
527-070-00	196, 227	537-324-00	196, 217	600-320-10	314	602-809-32	37, 333
527-080-00	196, 227	537-340-00	196, 217	600-320-30	314	602-810-10	334
527-090-00	196, 227	537-410-00	196, 218	600-320-60	314	602-810-12	33, 334
527-100-00	196, 227	537-411-00	196, 218	600-321-10	314	602-810-30	334
528-000-00	197, 220	537-413-00	196, 218	600-321-30	314	602-810-32	33, 334
528-000-01	197, 220	537-414-00	196, 218	600-321-60	314	602-811-10	338
528-000-02	199, 220	537-423-00	196, 218	600-322-10	315	602-811-12	31, 338
528-001-00	197, 220	537-424-00	196, 218	600-322-30	315	602-811-30	338
528-001-01	197, 220	537-440-00	196, 218	600-322-60	315	602-811-32	31, 338
528-039-50	413	537-510-00	196, 218	600-400-10	316	602-812-10	338
528-040-50	413	537-540-00	196, 218	600-400-30	316	602-812-30	338
528-041-00	334, 414	545-638-00	196, 217	600-500-10	318	602-815-10	332
528-042-00	334, 414	545-638-01	198, 219	600-500-30	318	602-815-12	36, 332
528-043-00	414	545-645-00	196, 217	600-500-60	318	602-815-30	332
528-050-00	414	545-645-01	198, 219	600-501-10	318	602-815-32	36, 332
528-070-00	196, 227	545-738-00	196, 217	600-501-30	318	602-819-10	332
528-080-00	196, 227	545-738-01	198, 219	600-501-60	318	602-819-12	36, 332
528-090-00	196, 227	545-745-00	196, 217	600-502-10	317	602-819-30	332
528-100-00	196, 227	545-745-01	198, 219	600-502-30	317	602-819-32	36, 332
528-110-00	196, 227	545-750-00	196, 217	600-502-60	317	602-830-10	335
528-120-00	196, 227	600-000-00	326	600-503-10	317	602-830-12	28, 335
528-130-00	196, 227	600-001-00	326	600-503-30	317	602-830-30	335
528-150-00	196, 227	600-001-30	326	600-503-60	317	602-830-32	28, 335
528-155-00	196, 228	600-002-00	326	600-600-30	320	602-831-10	335
528-158-00	196, 228	600-002-30	326	600-601-30	319	602-831-30	335
528-159-00	196, 228	600-003-00	326	600-602-30	319	602-832-10	335
528-200-10	413	600-004-00	326	600-610-00	327	602-832-12	28, 335
528-210-10	413	600-005-00	326	600-610-30	327	602-832-30	335
528-215-10	413	600-120-00	325	602-604-10	323	602-832-32	28, 335
528-220-00	413	600-121-10	325	602-605-10	322	602-833-10	335
529-000-00	197, 220	600-122-00	325	602-606-10	322	602-833-30	335
530-000-00	196, 220	600-123-00	325	602-800-10	337	602-834-10	336
530-035-00	198, 219	600-131-10	320	602-800-12	30, 337	602-834-12	29, 336
530-040-00	198, 219	600-132-10	321	602-800-30	337	602-834-30	336
530-045-00	198, 219	600-133-10	321	602-800-32	30, 337	602-834-32	29, 336
535-002-00	196, 227	600-300-10	314	602-801-10	337	602-835-10	336
535-003-00	196, 227	600-300-30	314	602-801-30	337	602-835-30	336
535-005-00	196, 227	600-300-60	314	602-806-10	333	602-836-10	330
535-009-00	196, 227	600-301-10	314	602-806-12	37, 333	602-836-12	34, 330
535-035-00	196, 217	600-301-30	314	602-806-30	333	602-836-30	330
535-040-00	196, 217	600-301-60	314	602-806-32	37, 333	602-836-32	34, 330
535-045-00	196, 217	600-302-10	315	602-808-10	337	602-837-10	331
535-048-00	196, 217	600-302-30	315	602-808-12	30, 337	602-837-12	35, 331

Numerical index

REF	Page	REF	Page	REF	Page	REF	Page
602-837-30	331	607-117-00	249	609-004-00	328	718-101-15	56
602-837-32	35, 331	607-118-00	249	609-006-00	328	718-101-21	56
602-838-10	330	607-120-00	246	609-008-00	328	718-101-22	56
602-838-12	34, 330	607-121-00	246, 249	609-100-00	328	718-101-24	56
602-838-30	330	607-121-30	248	609-200-30	328	718-101-25	56
602-838-32	34, 330	607-122-00	246, 249	609-300-01	328, 329	718-101-31	56
602-839-10	331	607-122-30	246, 248	609-400-00	328, 329	718-101-32	56
602-839-12	35, 331	607-123-10	250	609-500-00	328	718-101-34	56
602-839-30	331	607-124-00	246, 249	611-115-00	340	718-101-35	56
602-839-32	35, 331	607-125-00	249	611-115-30	340	718-101-41	56
602-840-10	337	607-126-20	246	611-115-60	340	718-101-42	56
602-840-30	337	607-126-22	246	611-116-00	340	718-101-44	56
602-845-10	338	607-126-24	246	611-120-01	340	718-101-45	56
602-845-12	31, 338	607-127-30	246, 249	611-122-00	340	718-141-13	56
602-845-30	338	607-128-00	249	611-130-00	340	718-141-23	56
602-845-32	31, 338	607-129-00	251, 253, 286, 287, 290, 292	611-131-00	340	718-141-33	56
602-846-10	339			611-132-00	340	718-141-43	56
602-846-12	32, 339	607-144-00	244, 245, 250	612-000-00	341	718-205-01	71
602-846-30	339	607-151-00	245	620-107-00	232	718-205-10	71
602-846-32	32, 339	607-152-00	245	620-108-00	232	718-205-11	71
602-847-10	339	607-153-00	244	620-109-00	232	718-205-12	71
602-847-12	32, 339	607-154-00	244, 245	620-110-00	232	718-205-13	71
602-847-30	339	607-156-00	244, 245	620-207-02	232	718-205-14	71
602-847-32	32, 339	607-157-00	244	620-208-02	232	718-205-15	71
603-000-00	327	607-158-00	244	620-209-02	232	718-205-21	71
606-600-00	324	607-159-00	244	621-028-01	231	718-205-22	71
606-600-30	324	607-160-00	244	621-031-01	231	718-205-23	71
606-701-00	324	607-165-00	244	621-034-01	231	718-205-24	71
606-701-30	324	607-180-00	251, 290	621-037-01	231	718-205-25	71
607-101-00	248	607-181-02	251, 252, 290, 291	621-040-01	231	718-205-31	71
607-102-00	248	607-181-10	252, 291	621-043-01	231	718-205-32	71
607-103-00	250	607-182-07	252, 291	621-046-01	231	718-205-33	71
607-103-11	244, 245, 250	607-182-10	252, 291	621-200-01	232	718-205-34	71
607-104-00	250	607-182-13	252, 291	621-300-01	232	718-205-35	71
607-104-11	244, 245, 250	607-183-07	252, 291	621-404-01	231	718-205-41	71
607-105-00	250	607-183-10	252, 291	621-405-01	231	718-205-42	71
607-105-11	245, 250	607-183-13	252, 291	621-406-01	231	718-205-43	71
607-107-00	250	607-184-07	251, 252, 290, 291	621-407-01	231	718-205-44	71
607-107-11	245, 250	607-184-10	251, 252, 290, 291	621-408-01	231	718-205-45	71
607-108-00	246	607-184-13	251, 252, 290, 291	621-409-01	231	718-301-01	86
607-108-10	246, 248	607-185-01	251, 253, 290, 292	621-500-01	231	718-301-10	86
607-108-20	246, 248	607-185-10	253, 292	681-100-00	418	718-301-11	86
607-109-00	248	607-202-00	233	705-000-00	114	718-301-12	86
607-110-00	248	607-202-30	233	718-101-01	56	718-301-14	86
607-111-00	246, 249	607-204-00	233	718-101-10	56	718-301-15	86
607-112-00	250	607-204-30	233	718-101-11	56	718-301-21	86
607-115-00	247	607-300-00	263	718-101-12	56	718-301-22	86
607-116-00	248	607-301-00	263	718-101-14	56	718-301-24	86

Numerical index

REF	Page	REF	Page	REF	Page	REF	Page
718-301-25	86	718-303-45	90	718-342-43	88	718-441-24	81
718-301-31	86	718-305-01	90	718-342-44	88	718-441-25	81
718-301-32	86	718-305-10	90	718-342-45	88	718-441-33	81
718-301-34	86	718-305-11	90	718-343-13	92	718-441-34	81
718-301-35	86	718-305-12	90	718-343-23	92	718-441-35	81
718-301-41	86	718-305-13	90	718-343-43	92	718-441-43	81
718-301-42	86	718-305-14	90	718-401-01	81	718-441-44	81
718-301-44	86	718-305-15	90	718-401-10	81	718-441-45	81
718-301-45	86	718-305-21	90	718-401-11	81	718-442-01	83
718-302-01	88	718-305-22	90	718-401-12	81	718-442-10	83
718-302-10	88	718-305-23	90	718-401-14	81	718-442-13	83
718-302-11	88	718-305-24	90	718-401-15	81	718-442-14	83
718-302-12	88	718-305-25	90	718-401-21	81	718-442-15	83
718-302-14	88	718-305-31	90	718-401-22	81	718-442-23	83
718-302-15	88	718-305-32	90	718-401-24	81	718-442-24	83
718-302-21	88	718-305-33	90	718-401-25	81	718-442-25	83
718-302-22	88	718-305-34	90	718-401-31	81	718-442-33	83
718-302-24	88	718-305-35	90	718-401-32	81	718-442-34	83
718-302-25	88	718-305-41	90	718-401-34	81	718-442-35	83
718-302-31	88	718-305-42	90	718-401-35	81	718-442-43	83
718-302-32	88	718-305-43	90	718-401-41	81	718-442-44	83
718-302-34	88	718-305-44	90	718-401-42	81	718-442-45	83
718-302-35	88	718-305-45	90	718-401-44	81	718-501-01	95
718-302-41	88	718-341-01	86	718-401-45	81	718-501-02	95
718-302-42	88	718-341-10	86	718-402-01	83	718-501-10	95
718-302-44	88	718-341-13	86	718-402-10	83	718-501-11	95
718-302-45	88	718-341-14	86	718-402-11	83	718-501-12	95
718-303-01	90	718-341-15	86	718-402-12	83	718-501-14	95
718-303-10	90	718-341-23	86	718-402-14	83	718-501-15	95
718-303-11	90	718-341-24	86	718-402-15	83	718-501-21	95
718-303-12	90	718-341-25	86	718-402-21	83	718-501-22	95
718-303-13	90	718-341-33	86	718-402-22	83	718-501-24	95
718-303-14	90	718-341-34	86	718-402-24	83	718-501-25	95
718-303-15	90	718-341-35	86	718-402-25	83	718-501-31	95
718-303-21	90	718-341-43	86	718-402-31	83	718-501-32	95
718-303-22	90	718-341-44	86	718-402-32	83	718-501-34	95
718-303-23	90	718-341-45	86	718-402-34	83	718-501-35	95
718-303-24	90	718-342-01	88	718-402-35	83	718-501-41	95
718-303-25	90	718-342-10	88	718-402-41	83	718-501-42	95
718-303-31	90	718-342-13	88	718-402-42	83	718-501-44	95
718-303-32	90	718-342-14	88	718-402-44	83	718-501-45	95
718-303-33	90	718-342-15	88	718-402-45	83	718-502-01	96
718-303-34	90	718-342-23	88	718-441-01	81	718-502-02	96
718-303-35	90	718-342-24	88	718-441-10	81	718-502-10	96
718-303-41	90	718-342-25	88	718-441-13	81	718-502-11	96
718-303-42	90	718-342-33	88	718-441-14	81	718-502-12	96
718-303-43	90	718-342-34	88	718-441-15	81	718-502-14	96
718-303-44	90	718-342-35	88	718-441-23	81	718-502-15	96

Numerical index

REF	Page	REF	Page	REF	Page	REF	Page
718-502-21	96	718-702-14	107	718-742-60	105	718-941-33	74
718-502-22	96	718-702-15	107	718-742-65	105	718-941-34	74
718-502-24	96	718-702-21	107	718-742-80	105	718-941-35	74
718-502-25	96	718-702-22	107	718-742-85	105	718-941-43	74
718-502-31	96	718-702-24	107	718-901-01	74	718-941-44	74
718-502-32	96	718-702-25	107	718-901-10	74	718-941-45	74
718-502-34	96	718-702-31	107	718-901-11	74	718-942-01	76
718-502-35	96	718-702-32	107	718-901-12	74	718-942-10	76
718-502-41	96	718-702-34	107	718-901-14	74	718-942-13	76
718-502-42	96	718-702-35	107	718-901-15	74	718-942-14	76
718-502-44	96	718-702-41	107	718-901-21	74	718-942-15	76
718-502-45	96	718-702-42	107	718-901-22	74	718-942-23	76
718-541-13	95	718-702-44	107	718-901-24	74	718-942-24	76
718-541-23	95	718-702-45	107	718-901-25	74	718-942-25	76
718-541-33	95	718-702-63	105	718-901-31	74	718-942-33	76
718-541-43	95	718-702-65	105	718-901-32	74	718-942-34	76
718-542-13	96	718-702-73	105	718-901-34	74	718-942-35	76
718-542-23	96	718-702-75	105	718-901-35	74	718-942-43	76
718-542-33	96	718-702-83	105	718-901-41	74	718-942-44	76
718-542-43	96	718-702-85	105	718-901-42	74	718-942-45	76
718-701-03	104	718-741-13	106	718-901-44	74	722-101-01	56
718-701-05	104	718-741-14	106	718-901-45	74	722-101-10	56
718-701-11	106	718-741-15	106	718-902-01	76	722-101-11	56
718-701-12	106	718-741-23	106	718-902-10	76	722-101-12	56
718-701-14	106	718-741-24	106	718-902-11	76	722-101-14	56
718-701-15	106	718-741-25	106	718-902-12	76	722-101-15	56
718-701-21	106	718-741-33	106	718-902-14	76	722-101-21	56
718-701-22	106	718-741-34	106	718-902-15	76	722-101-22	56
718-701-24	106	718-741-35	106	718-902-21	76	722-101-24	56
718-701-25	106	718-741-43	106	718-902-22	76	722-101-25	56
718-701-31	106	718-741-44	106	718-902-24	76	722-101-31	56
718-701-32	106	718-741-45	106	718-902-25	76	722-101-32	56
718-701-34	106	718-741-60	104	718-902-31	76	722-101-34	56
718-701-35	106	718-741-65	104	718-902-32	76	722-101-35	56
718-701-41	106	718-741-80	104	718-902-34	76	722-101-41	56
718-701-42	106	718-741-85	104	718-902-35	76	722-101-42	56
718-701-44	106	718-742-13	107	718-902-41	76	722-101-44	56
718-701-45	106	718-742-14	107	718-902-42	76	722-101-45	56
718-701-63	104	718-742-15	107	718-902-44	76	722-102-01	57
718-701-65	104	718-742-23	107	718-902-45	76	722-102-11	57
718-701-73	104	718-742-24	107	718-941-01	74	722-102-12	57
718-701-75	104	718-742-25	107	718-941-10	74	722-102-14	57
718-701-83	104	718-742-33	107	718-941-13	74	722-102-15	57
718-701-85	104	718-742-34	107	718-941-14	74	722-102-21	57
718-702-03	105	718-742-35	107	718-941-15	74	722-102-22	57
718-702-05	105	718-742-43	107	718-941-23	74	722-102-24	57
718-702-11	107	718-742-44	107	718-941-24	74	722-102-25	57
718-702-12	107	718-742-45	107	718-941-25	74	722-102-31	57

Numerical index

REF	Page	REF	Page	REF	Page	REF	Page
722-102-32	57	722-301-25	87	722-305-44	91	722-402-01	84
722-102-34	57	722-301-31	87	722-305-45	91	722-402-10	84
722-102-35	57	722-301-32	87	722-341-01	87	722-402-11	84
722-102-41	57	722-301-34	87	722-341-10	87	722-402-12	84
722-102-42	57	722-301-35	87	722-341-13	87	722-402-14	84
722-102-44	57	722-301-41	87	722-341-14	87	722-402-15	84
722-102-45	57	722-301-42	87	722-341-15	87	722-402-21	84
722-141-13	56	722-301-43	92	722-341-23	87	722-402-22	84
722-141-23	56	722-301-44	87	722-341-24	87	722-402-24	84
722-141-33	56	722-301-45	87	722-341-25	87	722-402-25	84
722-141-43	56	722-302-01	89	722-341-33	87	722-402-31	84
722-142-13	57	722-302-10	89	722-341-34	87	722-402-32	84
722-142-23	57	722-302-11	89	722-341-35	87	722-402-34	84
722-142-33	57	722-302-12	89	722-341-43	87	722-402-35	84
722-142-43	57	722-302-14	89	722-341-44	87	722-402-41	84
722-205-01	72	722-302-15	89	722-341-45	87	722-402-42	84
722-205-10	72	722-302-21	89	722-342-01	89	722-402-44	84
722-205-11	72	722-302-22	89	722-342-10	89	722-402-45	84
722-205-12	72	722-302-24	89	722-342-13	89	722-441-01	82
722-205-13	72	722-302-25	89	722-342-14	89	722-441-10	82
722-205-14	72	722-302-31	89	722-342-15	89	722-441-13	82
722-205-15	72	722-302-32	89	722-342-23	89	722-441-14	82
722-205-21	72	722-302-34	89	722-342-24	89	722-441-15	82
722-205-22	72	722-302-35	89	722-342-25	89	722-441-23	82
722-205-23	72	722-302-41	89	722-342-33	89	722-441-24	82
722-205-24	72	722-302-42	89	722-342-34	89	722-441-25	82
722-205-25	72	722-302-44	89	722-342-35	89	722-441-33	82
722-205-31	72	722-302-45	89	722-342-43	89	722-441-34	82
722-205-32	72	722-305-01	91	722-342-44	89	722-441-35	82
722-205-33	72	722-305-10	91	722-342-45	89	722-441-43	82
722-205-34	72	722-305-11	91	722-401-01	82	722-441-44	82
722-205-35	72	722-305-12	91	722-401-10	82	722-441-45	82
722-205-41	72	722-305-13	91	722-401-11	82	722-442-01	84
722-205-42	72	722-305-14	91	722-401-12	82	722-442-10	84
722-205-43	72	722-305-15	91	722-401-14	82	722-442-13	84
722-205-44	72	722-305-21	91	722-401-15	82	722-442-14	84
722-205-45	72	722-305-22	91	722-401-21	82	722-442-15	84
722-301-01	87	722-305-23	91	722-401-22	82	722-442-23	84
722-301-10	87	722-305-24	91	722-401-24	82	722-442-24	84
722-301-11	87	722-305-25	91	722-401-25	82	722-442-25	84
722-301-12	87	722-305-31	91	722-401-31	82	722-442-33	84
722-301-13	92	722-305-32	91	722-401-32	82	722-442-34	84
722-301-14	87	722-305-33	91	722-401-34	82	722-442-35	84
722-301-15	87	722-305-34	91	722-401-35	82	722-442-43	84
722-301-21	87	722-305-35	91	722-401-41	82	722-442-44	84
722-301-22	87	722-305-41	91	722-401-42	82	722-442-45	84
722-301-23	92	722-305-42	91	722-401-44	82	722-462-11	85
722-301-24	87	722-305-43	91	722-401-45	82	722-462-12	85

Numerical index

REF	Page	REF	Page	REF	Page	REF	Page
722-462-21	85	722-701-03	104	722-741-13	106	722-901-44	75
722-462-22	85	722-701-05	104	722-741-14	106	722-901-45	75
722-501-01	95	722-701-11	106	722-741-15	106	722-902-01	77
722-501-02	95	722-701-12	106	722-741-23	106	722-902-10	77
722-501-10	95	722-701-14	106	722-741-24	106	722-902-11	77
722-501-11	95	722-701-15	106	722-741-25	106	722-902-12	77
722-501-12	95	722-701-21	106	722-741-33	106	722-902-14	77
722-501-14	95	722-701-22	106	722-741-34	106	722-902-15	77
722-501-15	95	722-701-24	106	722-741-35	106	722-902-21	77
722-501-21	95	722-701-25	106	722-741-43	106	722-902-22	77
722-501-22	95	722-701-31	106	722-741-44	106	722-902-24	77
722-501-24	95	722-701-32	106	722-741-45	106	722-902-25	77
722-501-25	95	722-701-34	106	722-741-60	104	722-902-31	77
722-501-31	95	722-701-35	106	722-741-65	104	722-902-32	77
722-501-32	95	722-701-41	106	722-741-80	104	722-902-34	77
722-501-34	95	722-701-42	106	722-741-85	104	722-902-35	77
722-501-35	95	722-701-44	106	722-742-13	107	722-902-41	77
722-501-41	95	722-701-45	106	722-742-14	107	722-902-42	77
722-501-42	95	722-701-63	104	722-742-15	107	722-902-44	77
722-501-44	95	722-701-65	104	722-742-23	107	722-902-45	77
722-501-45	95	722-701-73	104	722-742-24	107	722-941-01	75
722-502-01	96	722-701-75	104	722-742-25	107	722-941-10	75
722-502-02	96	722-701-83	104	722-742-33	107	722-941-13	75
722-502-10	96	722-701-85	104	722-742-34	107	722-941-14	75
722-502-11	96	722-702-03	105	722-742-35	107	722-941-15	75
722-502-12	96	722-702-05	105	722-742-43	107	722-941-23	75
722-502-14	96	722-702-11	107	722-742-44	107	722-941-24	75
722-502-15	96	722-702-12	107	722-742-45	107	722-941-25	75
722-502-21	96	722-702-14	107	722-742-60	105	722-941-33	75
722-502-22	96	722-702-15	107	722-742-65	105	722-941-34	75
722-502-24	96	722-702-21	107	722-742-80	105	722-941-35	75
722-502-25	96	722-702-22	107	722-742-85	105	722-941-43	75
722-502-31	96	722-702-24	107	722-901-01	75	722-941-44	75
722-502-32	96	722-702-25	107	722-901-10	75	722-941-45	75
722-502-34	96	722-702-31	107	722-901-11	75	722-942-01	77
722-502-35	96	722-702-32	107	722-901-12	75	722-942-10	77
722-502-41	96	722-702-34	107	722-901-14	75	722-942-13	77
722-502-42	96	722-702-35	107	722-901-15	75	722-942-14	77
722-502-44	96	722-702-41	107	722-901-21	75	722-942-15	77
722-502-45	96	722-702-42	107	722-901-22	75	722-942-23	77
722-541-13	95	722-702-44	107	722-901-24	75	722-942-24	77
722-541-23	95	722-702-45	107	722-901-25	75	722-942-25	77
722-541-33	95	722-702-63	105	722-901-31	75	722-942-33	77
722-541-43	95	722-702-65	105	722-901-32	75	722-942-34	77
722-542-13	96	722-702-73	105	722-901-34	75	722-942-35	77
722-542-23	96	722-702-75	105	722-901-35	75	722-942-43	77
722-542-33	96	722-702-83	105	722-901-41	75	722-942-44	77
722-542-43	96	722-702-85	105	722-901-42	75	722-942-45	77

Numerical index

REF	Page	REF	Page	REF	Page	REF	Page
723-020-00	125	724-223-51	124	727-266-00	131	728-213-00	236
723-021-00	125	724-240-00	131	727-267-00	132	728-214-00	236
723-022-00	125	724-240-51	124	727-268-00	132	728-215-00	236
723-023-00	125	724-241-00	131	727-269-00	132	728-216-00	237
723-024-00	125	724-241-51	124	727-270-00	132	728-217-00	237
723-025-00	125	724-600-00	132	727-365-00	131	728-218-00	237
723-026-00	95, 125	724-601-00	132	727-366-00	131	728-219-00	237
723-027-00	95, 125	724-602-00	132	727-367-00	131	728-220-00	237
723-028-00	95, 96, 125	724-603-00	132	727-385-00	131	728-221-00	237
723-029-00	95, 96, 125	724-606-00	133	727-386-00	131	728-222-00	237
723-030-00	96, 125	724-607-00	128	727-387-00	131	728-400-00	237
723-031-00	96, 125	724-607-51	122	727-388-00	131	728-401-00	237
723-032-00	95, 96, 125	724-608-00	128	728-005-00	235	728-402-00	237
723-033-00	95, 96, 125	724-608-51	122	728-006-00	235	728-500-00	221
724-006-00	133	724-610-00	128	728-007-00	235	728-501-00	221
724-018-00	130	724-610-51	122	728-008-00	235	728-502-00	221
724-018-51	124	724-611-00	128	728-010-00	134	728-503-00	221
724-019-00	130	724-611-51	122	728-011-00	134	728-504-00	221
724-019-51	124	724-612-00	128	728-013-00	235	728-505-00	221
724-048-00	128	724-612-51	122	728-014-00	235	740-000-00	260
724-048-51	123	724-613-00	132	728-015-00	235	740-010-00	260
724-049-00	128	724-614-00	132	728-016-00	134	740-019-00	260
724-049-51	123	724-615-00	132	728-020-00	234	740-020-00	260
724-056-00	130	724-616-00	130	728-021-00	234, 235	740-021-00	261
724-060-00	132	724-617-00	130	728-022-00	234, 235	741-040-00	261
724-061-00	132	724-618-00	132	728-023-00	234, 235	741-045-00	261
724-090-50	122	724-619-00	132	728-024-00	234, 235	741-049-00	260
724-092-00	128	724-620-00	132	728-025-00	234, 235	741-050-00	261
724-092-51	122	724-621-00	132	728-026-00	234, 235	741-055-00	261
724-093-00	128	725-164-00	128	728-100-01	235	741-060-00	261
724-093-51	122	725-164-50	122	728-101-01	235	741-150-00	261
724-094-00	128	725-165-00	128	728-102-01	235	748-009-00	262
724-094-51	122	725-165-50	122	728-110-00	134	748-114-00	262
724-095-00	128	725-166-00	128	728-111-00	134	750-001-51	114
724-095-51	122	725-166-50	122	728-113-00	134	750-002-00	135
724-100-51	47	725-167-00	128	728-200-00	236	750-002-51	114
724-101-51	47	725-167-50	122	728-201-00	236	750-003-51	114
724-102-51	46, 47	725-168-00	128	728-202-00	236	750-602-00	135
724-103-51	46, 47	725-168-50	122	728-203-00	236	750-602-10	135, 142, 152, 157, 158, 161, 163
724-106-00	133	725-170-00	128	728-204-00	236		
724-107-00	133	725-170-50	122	728-205-00	236	750-602-50	114
724-156-00	130	725-171-00	128	728-206-00	236	750-603-00	135
724-157-00	130	725-171-50	122	728-207-00	236	751-000-00	222
724-160-00	132	725-203-00	126	728-208-00	236	751-001-00	222
724-161-00	132	725-203-51	120	728-209-00	236	751-003-00	222
724-222-00	131	725-206-00	126	728-210-00	236	751-004-00	222
724-222-51	124	726-011-00	126	728-211-00	236	751-005-00	222
724-223-00	131	727-265-00	131	728-212-00	236	758-160-00	223

Numerical index

REF	Page	REF	Page	REF	Page	REF	Page
758-161-00	223	766-061-00	192, 202	766-101-01	194, 210	766-188-00	193, 205
758-165-00	223	766-062-00	192, 202	766-103-00	194, 210	766-196-00	198, 219
758-175-00	223	766-063-00	192, 202	766-104-00	194, 210	766-200-00	194, 210
758-180-00	223	766-064-00	192, 202	766-105-00	194, 210	766-200-01	194, 210
758-360-00	224	766-065-00	192, 202	766-120-00	194, 209	766-201-00	194, 210
758-361-00	224, 286, 287	766-066-00	192, 202	766-120-01	194, 209	766-201-01	194, 210
758-362-00	224	766-067-00	192, 202	766-121-00	194, 209	766-203-00	194, 210
758-365-00	224	766-068-00	192, 202	766-121-01	194, 209	766-204-00	194, 210
758-367-00	224	766-069-00	192, 202	766-122-00	194, 209	766-205-00	194, 210
758-368-00	224	766-070-00	192, 202	766-122-01	194, 209	766-250-00	193, 205
758-375-00	224	766-071-00	192, 202	766-123-00	194, 209	766-251-00	193, 205
758-380-00	224	766-072-00	192, 202	766-123-01	194, 209	766-252-00	193, 205
758-861-00	224	766-073-00	192, 202	766-124-00	194, 209	766-253-00	193, 206
758-862-00	224	766-074-00	192, 202	766-125-00	194, 209	766-254-00	193, 206
765-200-00	194, 211	766-075-00	192, 202	766-130-00	194, 209	766-255-00	193, 206
765-201-00	194, 211	766-076-00	192, 202	766-131-00	194, 209	766-256-00	193, 206
765-201-01	194, 211	766-077-00	192, 202	766-132-00	194, 209	766-257-00	193, 206
765-202-00	194, 211	766-078-00	192, 202	766-132-01	194, 209	766-258-00	193, 206
765-202-01	194, 211	766-079-00	192, 202	766-133-00	194, 209	766-259-00	193, 206
765-203-00	194, 211	766-080-00	192, 200	766-133-01	194, 209	766-260-00	193, 206
765-300-00	194, 211	766-080-01	192, 200	766-160-00	192, 203	766-261-00	193, 206
765-301-00	194, 211	766-081-00	192, 200	766-161-00	192, 203	766-262-00	193, 206
765-301-01	194, 211	766-081-01	192, 200	766-162-00	192, 203	766-263-00	193, 206
765-302-00	194, 211	766-082-00	192, 200	766-163-00	192, 203	766-264-00	193, 206
765-302-01	194, 211	766-082-01	192, 200	766-164-00	192, 203	766-300-00	194, 210
765-303-00	194, 211	766-083-00	192, 200	766-165-00	192, 203	766-300-01	194, 210
766-000-00	231, 236, 237, 396	766-083-01	192, 200	766-166-00	192, 203	766-301-00	194, 210
766-014-00	192, 204	766-084-00	192, 200	766-167-00	192, 203	766-301-01	194, 210
766-015-00	192, 204	766-084-01	192, 200	766-168-00	192, 203	766-400-00	194, 210
766-018-00	192, 204	766-085-00	192, 200	766-169-00	192, 203	766-400-01	194, 210
766-019-00	192, 204	766-085-01	192, 200	766-170-00	192, 203	766-401-00	194, 210
766-032-00	192, 204	766-086-00	192, 202	766-171-00	192, 203	766-401-01	194, 210
766-033-00	192, 204	766-087-00	192, 202	766-172-00	192, 203	766-600-00	197, 218
766-034-00	192, 204	766-088-00	192, 201	766-173-00	192, 203	766-601-00	197, 218
766-036-00	192, 204	766-088-01	192, 201	766-174-00	192, 203	766-602-00	197, 218
766-037-00	192, 204	766-089-00	192, 201	766-175-00	192, 203	766-603-00	197, 218
766-039-00	192, 204	766-089-01	192, 201	766-176-00	192, 203	766-604-00	197, 218
766-050-00	192, 204	766-090-00	192, 201	766-177-00	192, 203	766-605-00	195, 214
766-051-00	192, 204	766-090-01	192, 201	766-178-00	192, 203	766-606-00	195, 214
766-052-00	192, 204	766-091-00	192, 201	766-179-00	192, 203	766-607-00	195, 214
766-053-00	192, 204	766-091-01	192, 201	766-180-00	193, 205	766-608-00	195, 214
766-054-00	192, 204	766-092-00	192, 201	766-181-00	193, 205	766-609-00	195, 214
766-055-00	192, 204	766-093-00	192, 201	766-182-00	193, 205	766-610-00	195, 214
766-056-00	192, 204	766-094-00	198, 219	766-183-00	193, 205	766-611-00	195, 214
766-057-00	192, 204	766-095-00	198, 219	766-184-00	193, 205	766-612-00	195, 214
766-058-00	192, 200	766-100-00	194, 210	766-185-00	193, 205	766-613-00	195, 214
766-059-00	192, 200	766-100-01	194, 210	766-186-00	193, 205	766-614-00	195, 214
766-060-00	192, 202	766-101-00	194, 210	766-187-00	193, 205	766-615-00	195, 214

Numerical index

REF	Page	REF	Page	REF	Page	REF	Page
766-616-00	195, 214	766-807-00	193, 208	768-027-00	195, 216	768-121-00	192, 203
766-617-00	195, 214	766-808-00	193, 208	768-028-00	195, 216	769-000-00	194, 212
766-618-00	195, 214	766-809-00	193, 208	768-029-00	195, 216	769-001-00	194, 212
766-619-00	195, 214	766-810-00	193, 208	768-030-00	195, 216	769-002-00	194, 212
766-620-00	195, 214	766-811-00	193, 208	768-031-00	195, 216	769-003-00	194, 212
766-629-00	197, 218	766-812-00	193, 208	768-032-00	195, 216	769-004-00	194, 212
766-631-00	195, 214	766-813-00	193, 208	768-033-00	195, 216	769-005-00	194, 212
766-632-00	195, 214	766-814-00	193, 208	768-034-00	195, 216	769-006-00	194, 212
766-633-00	197, 218	766-815-00	193, 208	768-035-00	195, 216	769-007-00	194, 212
766-700-00	193, 207	766-816-00	193, 208	768-036-00	222	769-008-00	194, 213
766-701-00	193, 207	766-817-00	193, 208	768-037-00	223	769-009-00	194, 213
766-702-00	193, 207	766-818-00	193, 208	768-038-00	195, 216	769-010-00	194, 213
766-703-00	193, 207	766-819-00	193, 208	768-039-00	195, 216	769-011-00	194, 213
766-704-00	193, 207	767-100-00	194, 212	768-040-00	195, 216	769-012-00	194, 213
766-705-00	193, 207	767-100-01	194, 212	768-041-00	195, 216	769-013-00	194, 213
766-706-00	193, 207	767-101-00	194, 212	768-100-00	192, 200	769-014-00	194, 213
766-707-00	193, 207	767-101-01	194, 212	768-100-01	192, 200	769-015-00	194, 213
766-710-00	193, 207	767-102-00	194, 212	768-101-00	192, 200	769-016-00	194, 213
766-711-00	193, 207	767-103-00	194, 212	768-101-01	192, 200	769-017-00	194, 213
766-712-00	193, 207	767-103-01	194, 212	768-102-00	192, 201	769-018-00	194, 213
766-713-00	193, 207	767-104-00	194, 212	768-102-01	192, 201	769-019-00	194, 213
766-714-00	193, 207	767-105-00	194, 212	768-103-00	192, 201	769-100-00	194, 213
766-715-00	193, 207	767-105-01	194, 212	768-103-01	192, 201	769-101-00	194, 213
766-716-00	193, 207	767-106-00	194, 212	768-104-00	192, 201	769-200-00	194, 213
766-717-00	193, 207	767-107-00	194, 212	768-104-01	192, 201	769-201-00	194, 213
766-718-00	193, 207	767-200-00	194, 212	768-105-00	192, 201	772-000-00	265
766-719-00	193, 207	767-200-01	194, 212	768-105-01	192, 201	772-100-00	265
766-720-00	193, 207	767-201-00	194, 212	768-106-00	192, 201	772-103-00	264
766-721-00	193, 207	767-201-01	194, 212	768-106-01	192, 201	772-104-00	264
766-760-00	48	767-202-00	194, 212	768-107-00	192, 201	772-106-00	264
766-761-00	48	767-203-00	194, 212	768-107-01	192, 201	772-107-00	264
766-762-00	48	767-203-01	194, 212	768-108-00	192, 201	772-109-00	264
766-763-00	48	767-204-00	194, 212	768-108-01	192, 201	772-113-00	264
766-764-00	48	767-205-00	194, 212	768-109-00	192, 201	772-114-00	264
766-765-00	48	767-205-01	194, 212	768-109-01	192, 201	772-116-00	264
766-766-00	48	767-206-00	194, 212	768-110-00	192, 201	772-117-00	264
766-767-00	48	767-207-00	194, 212	768-110-01	192, 201	772-119-00	264
766-768-00	48	768-014-00	195, 215	768-111-00	192, 201	772-123-00	264
766-769-00	48	768-015-00	195, 215	768-111-01	192, 201	772-124-00	264
766-770-00	48	768-016-00	195, 215	768-112-00	192, 202	772-125-00	264
766-771-00	48	768-017-00	195, 215	768-113-00	192, 202	772-126-00	264
766-800-00	193, 208	768-018-00	195, 215	768-114-00	192, 202	772-127-00	264
766-801-00	193, 208	768-019-00	195, 215	768-115-00	192, 202	772-200-00	265
766-802-00	193, 208	768-020-00	195, 215	768-116-00	192, 203	772-203-02	262
766-803-00	193, 208	768-021-00	195, 215	768-117-00	192, 203	772-204-02	262
766-804-00	193, 208	768-022-00	195, 215	768-118-00	192, 203	772-205-02	262
766-805-00	193, 208	768-023-00	195, 215	768-119-00	192, 203	772-208-02	262
766-806-00	193, 208	768-026-00	195, 216	768-120-00	192, 203	772-210-02	262

Numerical index

REF	Page	REF	Page	REF	Page	REF	Page
772-313-01	265	790-113-00	73	790-217-02	93	790-322-00	79
772-314-01	265	790-114-00	69	790-218-02	93, 157, 161	790-323-00	79
772-316-01	265	790-115-00	69	790-220-02	93	790-341-00	78
772-317-01	265	790-116-00	69	790-221-02	93	790-342-00	78
772-323-01	265	790-117-00	69	790-222-02	93, 157, 162	790-343-00	78
772-324-01	265	790-118-00	69	790-223-02	93, 158, 163	790-344-00	78
772-325-01	265	790-119-00	69	790-228-02	93	790-345-00	78
772-326-01	265	790-120-00	69	790-229-02	93	790-346-00	78
772-400-00	264	790-121-00	69	790-230-02	93	790-347-00	79
774-002-01	266	790-122-00	69	790-231-02	93	790-348-00	79
774-100-01	266	790-123-00	69	790-234-02	93	790-349-00	79
774-200-01	266	790-128-00	70, 71	790-243-02	93	790-350-00	79
774-220-01	50, 266	790-130-00	70, 71	790-249-02	93	790-351-00	79
774-221-01	50, 266, 286	790-131-00	70, 71	790-259-02	93	790-352-00	79
774-226-01	266	790-132-00	70, 71	790-260-02	93	790-363-00	78
774-229-01	266	790-133-00	70, 71	790-261-02	93	790-364-00	78
774-231-01	266	790-134-00	70, 71	790-262-02	93	790-365-00	79
774-232-01	266	790-135-00	70, 71	790-263-02	93	790-366-00	79
774-233-01	266	790-136-00	70	790-264-02	93	790-450-00	46, 47
774-235-01	266	790-137-00	70	790-265-02	93	790-451-00	46, 47
774-236-01	109, 266	790-138-00	70, 71	790-266-02	93	790-452-00	46, 47
774-300-00	262	790-139-00	70, 71	790-274-02	93	790-453-00	46, 47
774-549-01	267	790-140-00	70	790-275-02	93	790-603-02	94
774-551-02	267	790-141-00	68	790-278-02	93	790-607-02	94
774-552-01	267	790-142-00	68	790-279-02	93	799-202-00	114
774-553-01	267	790-143-00	68	790-290-02	94	799-206-00	114
774-554-01	267	790-144-00	68	790-291-02	94	812-038-00	168
774-555-01	267	790-145-00	68	790-294-02	94	813-045-00	168
774-556-01	267	790-146-00	68	790-300-00	80	818-111-10	127
774-557-01	267	790-147-00	69	790-302-00	78	818-111-16	121
774-558-01	267	790-148-00	69	790-303-00	78	818-111-20	127
774-559-01	267	790-149-00	69	790-304-00	78	818-111-26	121
774-560-01	267	790-150-00	69	790-305-00	78	818-111-30	127
774-561-04	267	790-151-00	69	790-306-00	78	818-111-36	121
774-563-01	109, 267	790-152-00	69	790-307-00	78	818-111-40	127
790-100-00	73	790-155-00	73	790-308-00	78	818-111-46	121
790-101-00	73	790-157-00	73	790-309-00	78	818-112-10	127
790-102-00	68	790-163-00	68	790-310-00	78	818-112-16	121
790-103-00	68	790-164-00	68	790-311-00	78	818-112-20	127
790-104-00	68	790-165-00	69	790-312-00	80	818-112-26	121
790-105-00	68	790-166-00	69	790-314-00	79	818-112-40	127
790-106-00	68, 70, 71	790-204-02	93	790-315-00	79	818-112-46	121
790-107-00	68, 70, 71	790-205-02	93	790-316-00	79	818-121-10	129
790-108-00	68	790-206-02	93, 157, 161	790-317-00	79	818-121-16	123
790-109-00	68	790-208-02	93	790-318-00	79	818-121-20	129
790-110-00	68, 70, 71	790-210-02	93, 157, 162	790-319-00	79	818-121-26	123
790-111-00	68	790-211-02	93, 158, 163	790-320-00	79	818-121-36	123
790-112-00	73	790-216-02	93	790-321-00	79	818-121-46	123

Numerical index

REF	Page	REF	Page	REF	Page	REF	Page
818-122-10	129	818-222-10	129	822-122-16	123	822-219-36	121
818-122-16	123	818-222-16	123	822-122-20	129	822-219-46	121
818-122-20	129	818-222-20	129	822-122-26	123	822-221-10	129
818-122-26	123	818-222-26	123	822-122-36	123	822-221-16	123
818-122-46	123	818-222-30	129	822-122-46	123	822-221-20	129
818-125-30	129	818-222-40	129	822-125-30	129	822-221-26	123
818-125-40	129	818-222-46	123	822-125-40	129	822-221-30	129
818-131-10	132, 142, 143, 147, 149	818-224-40	129	822-131-10	132, 142, 143, 147, 149	822-221-36	123
		818-227-10	129			822-221-40	129
818-131-20	132, 142, 144, 147, 150	818-227-20	129	822-131-20	132, 142, 144, 147, 150	822-221-46	123
		818-411-17	120			822-222-10	129
818-132-10	132	818-411-27	120	822-132-10	132	822-222-16	123
818-132-20	132	818-411-37	120	822-132-20	132	822-222-20	129
818-141-10	131	818-411-47	120	822-137-10	132	822-222-26	123
818-141-20	131	818-412-17	120	822-137-20	132	822-222-30	129
818-147-40	131	818-412-27	120	822-141-10	131	822-222-36	123
818-151-30	130, 142, 146, 148, 152	818-412-37	120	822-141-20	131	822-222-40	129
		818-412-46	120	822-142-10	131	822-222-46	123
818-151-40	130, 142, 145, 148, 151	818-412-47	120	822-142-20	131	822-224-30	129
		818-413-17	120	822-147-40	131	822-224-40	129
818-152-30	130	818-413-27	120	822-151-30	130, 142, 146, 148, 152	822-227-10	129
818-152-40	130	818-415-17	120			822-227-20	129
818-211-16	121	818-415-27	120	822-151-40	130, 142, 145, 148, 151	822-411-17	120
818-211-26	121	818-512-30	127			822-411-27	120
818-211-36	121	818-512-40	127	822-152-30	130	822-411-37	120
818-211-46	121	822-111-10	127	822-152-40	130	822-411-47	120
818-212-10	127	822-111-16	121	822-211-16	121	822-412-17	120
818-212-16	121	822-111-20	127	822-211-26	121	822-412-27	120
818-212-20	127	822-111-26	121	822-211-36	121	822-412-36	120
818-212-26	121	822-111-30	127	822-211-46	121	822-412-37	120
818-212-36	121	822-111-36	121	822-212-10	127	822-412-46	120
818-212-37	121	822-111-40	127	822-212-16	121	822-412-47	120
818-212-46	121	822-111-46	121	822-212-20	127	822-415-17	120
818-212-47	121	822-112-10	127	822-212-26	121	822-415-27	120
818-215-16	121	822-112-16	121	822-212-36	121	822-415-37	120
818-215-26	121	822-112-20	127	822-212-37	121	822-415-47	120
818-217-10	127	822-112-26	121	822-212-46	121	822-512-30	127
818-217-20	127	822-112-30	127	822-212-47	121	822-512-40	127
818-217-30	127	822-112-36	121	822-215-16	121	831-601-00	168
818-217-40	127	822-112-40	127	822-215-26	121	831-603-00	168
818-221-10	129	822-112-46	121	822-215-36	121	831-606-00	168
818-221-16	123	822-121-10	129	822-215-46	121	831-701-00	168
818-221-20	129	822-121-16	123	822-217-10	127	850-101-00	166
818-221-26	123	822-121-20	129	822-217-20	127	850-102-00	166
818-221-30	129	822-121-26	123	822-217-30	127	850-103-00	166
818-221-36	123	822-121-36	123	822-217-40	127	850-104-00	166
818-221-40	129	822-121-46	123	822-218-30	127	850-105-00	166
818-221-46	123	822-122-10	129	822-218-40	127	850-106-00	166

Numerical index

REF	Page	REF	Page	REF	Page	REF	Page
850-107-00	166	860-105-22	161	860-121-22	161	861-105-18	162
850-108-00	166	860-106-00	161	860-122-00	161	861-105-22	162
850-109-00	166	860-106-18	161	860-122-18	161	861-106-00	162, 163
850-110-00	166	860-106-22	161	860-122-22	161	861-106-18	162
850-111-00	166	860-107-00	161	860-123-00	161	861-106-22	162
850-112-00	166	860-107-18	161	860-123-18	161	861-107-00	162, 163
850-113-00	166	860-107-22	161	860-123-22	161	861-107-18	162
850-114-00	166	860-108-00	161	860-124-00	161	861-107-22	162
850-121-00	166	860-108-18	161	860-124-18	161	861-108-00	162, 163
850-122-00	166	860-108-22	161	860-124-22	161	861-108-18	162
850-123-00	166	860-109-00	161	860-125-00	161	861-108-22	162
850-124-00	166	860-109-18	161	860-125-18	161	861-109-00	162, 163
850-125-00	166	860-109-22	161	860-125-22	161	861-109-18	162
851-101-00	167	860-110-00	161	860-126-00	161	861-109-22	162
851-102-00	167	860-110-18	161	860-126-18	161	861-110-00	162, 163
851-103-00	167	860-110-22	161	860-126-22	161	861-110-18	162
851-104-00	167	860-111-00	161	860-127-00	161	861-110-22	162
851-105-00	167	860-111-18	161	860-127-18	161	861-111-00	162, 163
851-106-00	167	860-111-22	161	860-127-22	161	861-111-18	162
851-107-00	167	860-112-00	161	860-128-00	161	861-111-22	162
851-108-00	167	860-112-18	161	860-128-18	161	861-112-00	162, 163
851-109-00	167	860-112-22	161	860-128-22	161	861-112-18	162
851-110-00	167	860-113-00	161	860-129-00	161	861-112-22	162
851-117-00	167	860-113-18	161	860-129-18	161	861-113-00	162, 163
851-118-00	167	860-113-22	161	860-129-22	161	861-113-18	162
851-119-00	167	860-114-00	161	860-130-00	161	861-113-22	162
851-120-00	167	860-114-18	161	860-130-18	161	861-114-00	162, 163
851-121-00	167	860-114-22	161	860-130-22	161	861-114-18	162
855-050-00	160	860-115-00	161	860-131-00	161	861-114-22	162
855-150-00	160	860-115-18	161	860-131-18	161	861-115-00	162, 163
857-042-00	160	860-115-22	161	860-131-22	161	861-115-18	162
857-150-00	160	860-116-00	161	860-150-18	157	861-115-22	162
860-064-18	157	860-116-18	161	860-150-22	157	861-116-00	162, 163
860-064-22	157	860-116-22	161	861-064-18	157	861-116-18	162
860-101-00	161	860-117-00	161	861-064-22	157	861-116-22	162
860-101-18	161	860-117-18	161	861-101-00	162, 163	861-117-00	162
860-101-22	161	860-117-22	161	861-101-18	162	861-117-18	162
860-102-00	161	860-118-00	161	861-101-22	162	861-117-22	162
860-102-18	161	860-118-18	161	861-102-00	162, 163	861-118-00	162
860-102-22	161	860-118-22	161	861-102-18	162	861-118-18	162
860-103-00	161	860-119-00	161	861-102-22	162	861-118-22	162
860-103-18	161	860-119-18	161	861-103-00	162, 163	861-119-00	162
860-103-22	161	860-119-22	161	861-103-18	162	861-119-18	162
860-104-00	161	860-120-00	161	861-103-22	162	861-119-22	162
860-104-18	161	860-120-18	161	861-104-00	162, 163	861-120-00	162
860-104-22	161	860-120-22	161	861-104-18	162	861-120-18	162
860-105-00	161	860-121-00	161	861-104-22	162	861-120-22	162
860-105-18	161	860-121-18	161	861-105-00	162, 163	861-121-00	162

Numerical index

REF	Page	REF	Page	REF	Page	REF	Page
861-121-18	162	862-107-18	163	862-131-18	163	871-008-00	165
861-121-22	162	862-107-22	163	862-131-22	163	871-009-00	165
861-122-00	162	862-108-18	163	862-150-18	158	871-010-00	165
861-122-18	162	862-108-22	163	862-150-22	158	871-011-00	165
861-122-22	162	862-109-18	163	864-064-00	157	871-012-00	165
861-123-00	162	862-109-22	163	864-150-00	157	871-013-00	165
861-123-18	162	862-110-18	163	864-300-00	157	871-014-00	165
861-123-22	162	862-110-22	163	865-064-00	158	871-015-00	165
861-124-00	162	862-111-18	163	865-150-00	158	871-016-00	165
861-124-18	162	862-111-22	163	865-300-00	158	871-017-00	165
861-124-22	162	862-112-18	163	870-001-00	164	871-018-00	165
861-125-00	162	862-112-22	163	870-002-00	164	871-019-00	165
861-125-18	162	862-113-18	163	870-003-00	164	871-020-00	165
861-125-22	162	862-113-22	163	870-004-00	164	871-021-00	165
861-126-00	162	862-114-18	163	870-005-00	164	871-022-00	165
861-126-18	162	862-114-22	163	870-006-00	164	871-023-00	165
861-126-22	162	862-115-18	163	870-007-00	164	871-024-00	165
861-127-00	162	862-115-22	163	870-008-00	164	871-025-00	165
861-127-18	162	862-116-18	163	870-009-00	164	871-026-00	165
861-127-22	162	862-116-22	163	870-010-00	164	871-027-00	165
861-128-00	162	862-117-18	163	870-011-00	164	871-028-00	165
861-128-18	162	862-117-22	163	870-012-00	164	871-029-00	165
861-128-22	162	862-118-18	163	870-013-00	164	871-030-00	165
861-129-00	162	862-118-22	163	870-014-00	164	871-031-00	165
861-129-18	162	862-119-18	163	870-015-00	164	872-064-00	159
861-129-22	162	862-119-22	163	870-016-00	164	873-064-00	159
861-130-00	162	862-120-18	163	870-017-00	164	878-001-00	143
861-130-18	162	862-120-22	163	870-018-00	164	878-002-00	143
861-130-22	162	862-121-18	163	870-019-00	164	878-003-00	143
861-131-00	162	862-121-22	163	870-020-00	164	878-004-00	143
861-131-18	162	862-122-18	163	870-021-00	164	878-005-00	143
861-131-22	162	862-122-22	163	870-022-00	164	878-006-00	143
861-150-18	157	862-123-18	163	870-023-00	164	878-007-00	143
861-150-22	157	862-123-22	163	870-024-00	164	878-008-00	143
862-064-18	158	862-124-18	163	870-025-00	164	878-009-00	143
862-064-22	158	862-124-22	163	870-026-00	164	878-010-00	143
862-101-18	163	862-125-18	163	870-027-00	164	878-011-00	143
862-101-22	163	862-125-22	163	870-028-00	164	878-012-00	143
862-102-18	163	862-126-18	163	870-029-00	164	878-013-00	143
862-102-22	163	862-126-22	163	870-030-00	164	878-014-00	143
862-103-18	163	862-127-18	163	870-031-00	164	878-015-00	143
862-103-22	163	862-127-22	163	871-001-00	165	878-016-00	143
862-104-18	163	862-128-18	163	871-002-00	165	878-017-00	143
862-104-22	163	862-128-22	163	871-003-00	165	878-018-00	143
862-105-18	163	862-129-18	163	871-004-00	165	878-064-00	142
862-105-22	163	862-129-22	163	871-005-00	165	878-301-18	143
862-106-18	163	862-130-18	163	871-006-00	165	878-301-22	143
862-106-22	163	862-130-22	163	871-007-00	165	878-302-18	143

Numerical index

REF	Page	REF	Page	REF	Page	REF	Page
878-302-22	143	878-326-22	143	879-305-18	144	879-329-18	144
878-303-18	143	878-327-18	143	879-305-22	144	879-329-22	144
878-303-22	143	878-327-22	143	879-306-18	144	879-330-18	144
878-304-18	143	878-328-18	143	879-306-22	144	879-330-22	144
878-304-22	143	878-328-22	143	879-307-18	144	879-331-18	144
878-305-18	143	878-329-18	143	879-307-22	144	879-331-22	144
878-305-22	143	878-329-22	143	879-308-18	144	879-332-18	144
878-306-18	143	878-330-18	143	879-308-22	144	879-332-22	144
878-306-22	143	878-330-22	143	879-309-18	144	879-333-18	144
878-307-18	143	878-331-18	143	879-309-22	144	879-333-22	144
878-307-22	143	878-331-22	143	879-310-18	144	879-334-18	144
878-308-18	143	878-332-18	143	879-310-22	144	879-334-22	144
878-308-22	143	878-332-22	143	879-311-18	144	879-335-18	144
878-309-18	143	878-333-18	143	879-311-22	144	879-335-22	144
878-309-22	143	878-333-22	143	879-312-18	144	879-364-18	142
878-310-18	143	878-334-18	143	879-312-22	144	879-364-22	142
878-310-22	143	878-334-22	143	879-313-18	144	880-201-00	145
878-311-18	143	878-335-18	143	879-313-22	144	880-202-00	145
878-311-22	143	878-335-22	143	879-314-18	144	880-203-00	145
878-312-18	143	878-364-18	142	879-314-22	144	880-204-00	145
878-312-22	143	878-364-22	142	879-315-18	144	880-205-00	145
878-313-18	143	879-001-00	144	879-315-22	144	880-206-00	145
878-313-22	143	879-002-00	144	879-316-18	144	880-207-00	145
878-314-18	143	879-003-00	144	879-316-22	144	880-208-00	145
878-314-22	143	879-004-00	144	879-317-18	144	880-209-00	145
878-315-18	143	879-005-00	144	879-317-22	144	880-210-00	145
878-315-22	143	879-006-00	144	879-318-18	144	880-211-00	145
878-316-18	143	879-007-00	144	879-318-22	144	880-212-00	145
878-316-22	143	879-008-00	144	879-319-18	144	880-213-00	145
878-317-18	143	879-009-00	144	879-319-22	144	880-214-00	145
878-317-22	143	879-010-00	144	879-320-18	144	880-215-00	145
878-318-18	143	879-011-00	144	879-320-22	144	880-216-00	145
878-318-22	143	879-012-00	144	879-321-18	144	880-217-00	145
878-319-18	143	879-013-00	144	879-321-22	144	880-218-00	145
878-319-22	143	879-014-00	144	879-322-18	144	880-264-00	142
878-320-18	143	879-015-00	144	879-322-22	144	880-301-18	145
878-320-22	143	879-016-00	144	879-323-18	144	880-301-22	145
878-321-18	143	879-017-00	144	879-323-22	144	880-302-18	145
878-321-22	143	879-018-00	144	879-324-18	144	880-302-22	145
878-322-18	143	879-064-00	142	879-324-22	144	880-303-18	145
878-322-22	143	879-301-18	144	879-325-18	144	880-303-22	145
878-323-18	143	879-301-22	144	879-325-22	144	880-304-18	145
878-323-22	143	879-302-18	144	879-326-18	144	880-304-22	145
878-324-18	143	879-302-22	144	879-326-22	144	880-305-18	145
878-324-22	143	879-303-18	144	879-327-18	144	880-305-22	145
878-325-18	143	879-303-22	144	879-327-22	144	880-306-18	145
878-325-22	143	879-304-18	144	879-328-18	144	880-306-22	145
878-326-18	143	879-304-22	144	879-328-22	144	880-307-18	145

Numerical index

REF	Page	REF	Page	REF	Page	REF	Page
880-307-22	145	880-331-22	145	881-310-18	146	881-334-18	146
880-308-18	145	880-332-18	145	881-310-22	146	881-334-22	146
880-308-22	145	880-332-22	145	881-311-18	146	881-335-18	146
880-309-18	145	880-333-18	145	881-311-22	146	881-335-22	146
880-309-22	145	880-333-22	145	881-312-18	146	881-364-18	142
880-310-18	145	880-334-18	145	881-312-22	146	881-364-22	142
880-310-22	145	880-334-22	145	881-313-18	146	882-001-00	149
880-311-18	145	880-335-18	145	881-313-22	146	882-002-00	149
880-311-22	145	880-335-22	145	881-314-18	146	882-003-00	149
880-312-18	145	880-364-18	142	881-314-22	146	882-004-00	149
880-312-22	145	880-364-22	142	881-315-18	146	882-005-00	149
880-313-18	145	881-201-00	146	881-315-22	146	882-006-00	149
880-313-22	145	881-202-00	146	881-316-18	146	882-007-00	149
880-314-18	145	881-203-00	146	881-316-22	146	882-008-00	149
880-314-22	145	881-204-00	146	881-317-18	146	882-009-00	149
880-315-18	145	881-205-00	146	881-317-22	146	882-010-00	149
880-315-22	145	881-206-00	146	881-318-18	146	882-011-00	149
880-316-18	145	881-207-00	146	881-318-22	146	882-012-00	149
880-316-22	145	881-208-00	146	881-319-18	146	882-013-00	149
880-317-18	145	881-209-00	146	881-319-22	146	882-014-00	149
880-317-22	145	881-210-00	146	881-320-18	146	882-015-00	149
880-318-18	145	881-211-00	146	881-320-22	146	882-016-00	149
880-318-22	145	881-212-00	146	881-321-18	146	882-017-00	149
880-319-18	145	881-213-00	146	881-321-22	146	882-018-00	149
880-319-22	145	881-214-00	146	881-322-18	146	882-064-00	147
880-320-18	145	881-215-00	146	881-322-22	146	882-150-00	147
880-320-22	145	881-216-00	146	881-323-18	146	882-301-18	149
880-321-18	145	881-217-00	146	881-323-22	146	882-301-22	149
880-321-22	145	881-218-00	146	881-324-18	146	882-302-18	149
880-322-18	145	881-264-00	142	881-324-22	146	882-302-22	149
880-322-22	145	881-301-18	146	881-325-18	146	882-303-18	149
880-323-18	145	881-301-22	146	881-325-22	146	882-303-22	149
880-323-22	145	881-302-18	146	881-326-18	146	882-304-18	149
880-324-18	145	881-302-22	146	881-326-22	146	882-304-22	149
880-324-22	145	881-303-18	146	881-327-18	146	882-305-18	149
880-325-18	145	881-303-22	146	881-327-22	146	882-305-22	149
880-325-22	145	881-304-18	146	881-328-18	146	882-306-18	149
880-326-18	145	881-304-22	146	881-328-22	146	882-306-22	149
880-326-22	145	881-305-18	146	881-329-18	146	882-307-18	149
880-327-18	145	881-305-22	146	881-329-22	146	882-307-22	149
880-327-22	145	881-306-18	146	881-330-18	146	882-308-18	149
880-328-18	145	881-306-22	146	881-330-22	146	882-308-22	149
880-328-22	145	881-307-18	146	881-331-18	146	882-309-18	149
880-329-18	145	881-307-22	146	881-331-22	146	882-309-22	149
880-329-22	145	881-308-18	146	881-332-18	146	882-310-18	149
880-330-18	145	881-308-22	146	881-332-22	146	882-310-22	149
880-330-22	145	881-309-18	146	881-333-18	146	882-311-18	149
880-331-18	145	881-309-22	146	881-333-22	146	882-311-22	149

Numerical index

REF	Page	REF	Page	REF	Page	REF	Page
882-312-18	149	882-350-18	147	883-313-18	150	883-364-18	147
882-312-22	149	882-350-22	147	883-313-22	150	883-364-22	147
882-313-18	149	882-364-18	147	883-314-18	150	884-001-00	151
882-313-22	149	882-364-22	147	883-314-22	150	884-002-00	151
882-314-18	149	883-001-00	150	883-315-18	150	884-003-00	151
882-314-22	149	883-002-00	150	883-315-22	150	884-004-00	151
882-315-18	149	883-003-00	150	883-316-18	150	884-005-00	151
882-315-22	149	883-004-00	150	883-316-22	150	884-006-00	151
882-316-18	149	883-005-00	150	883-317-18	150	884-007-00	151
882-316-22	149	883-006-00	150	883-317-22	150	884-008-00	151
882-317-18	149	883-007-00	150	883-318-18	150	884-009-00	151
882-317-22	149	883-008-00	150	883-318-22	150	884-010-00	151
882-318-18	149	883-009-00	150	883-319-18	150	884-011-00	151
882-318-22	149	883-010-00	150	883-319-22	150	884-012-00	151
882-319-18	149	883-011-00	150	883-320-18	150	884-013-00	151
882-319-22	149	883-012-00	150	883-320-22	150	884-014-00	151
882-320-18	149	883-013-00	150	883-321-18	150	884-015-00	151
882-320-22	149	883-014-00	150	883-321-22	150	884-016-00	151
882-321-18	149	883-015-00	150	883-322-18	150	884-017-00	151
882-321-22	149	883-016-00	150	883-322-22	150	884-018-00	151
882-322-18	149	883-017-00	150	883-323-18	150	884-064-00	148
882-322-22	149	883-018-00	150	883-323-22	150	884-150-00	148
882-323-18	149	883-064-00	147	883-324-18	150	884-301-18	151
882-323-22	149	883-150-00	147	883-324-22	150	884-301-22	151
882-324-18	149	883-301-18	150	883-325-18	150	884-302-18	151
882-324-22	149	883-301-22	150	883-325-22	150	884-302-22	151
882-325-18	149	883-302-18	150	883-326-18	150	884-303-18	151
882-325-22	149	883-302-22	150	883-326-22	150	884-303-22	151
882-326-18	149	883-303-18	150	883-327-18	150	884-304-18	151
882-326-22	149	883-303-22	150	883-327-22	150	884-304-22	151
882-327-18	149	883-304-18	150	883-328-18	150	884-305-18	151
882-327-22	149	883-304-22	150	883-328-22	150	884-305-22	151
882-328-18	149	883-305-18	150	883-329-18	150	884-306-18	151
882-328-22	149	883-305-22	150	883-329-22	150	884-306-22	151
882-329-18	149	883-306-18	150	883-330-18	150	884-307-18	151
882-329-22	149	883-306-22	150	883-330-22	150	884-307-22	151
882-330-18	149	883-307-18	150	883-331-18	150	884-308-18	151
882-330-22	149	883-307-22	150	883-331-22	150	884-308-22	151
882-331-18	149	883-308-18	150	883-332-18	150	884-309-18	151
882-331-22	149	883-308-22	150	883-332-22	150	884-309-22	151
882-332-18	149	883-309-18	150	883-333-18	150	884-310-18	151
882-332-22	149	883-309-22	150	883-333-22	150	884-310-22	151
882-333-18	149	883-310-18	150	883-334-18	150	884-311-18	151
882-333-22	149	883-310-22	150	883-334-22	150	884-311-22	151
882-334-18	149	883-311-18	150	883-335-18	150	884-312-18	151
882-334-22	149	883-311-22	150	883-335-22	150	884-312-22	151
882-335-18	149	883-312-18	150	883-350-18	147	884-313-18	151
882-335-22	149	883-312-22	150	883-350-22	147	884-313-22	151

Numerical index

REF	Page	REF	Page	REF	Page	REF	Page
884-314-18	151	885-001-00	152	885-315-18	152	886-003-00	153
884-314-22	151	885-002-00	152	885-315-22	152	886-004-00	153
884-315-18	151	885-003-00	152	885-316-18	152	886-005-00	153
884-315-22	151	885-004-00	152	885-316-22	152	886-006-00	153
884-316-18	151	885-005-00	152	885-317-18	152	886-007-00	153
884-316-22	151	885-006-00	152	885-317-22	152	886-008-00	153
884-317-18	151	885-007-00	152	885-318-18	152	886-009-00	153
884-317-22	151	885-008-00	152	885-318-22	152	886-010-00	153
884-318-18	151	885-009-00	152	885-319-18	152	886-011-00	153
884-318-22	151	885-010-00	152	885-319-22	152	886-012-00	153
884-319-18	151	885-011-00	152	885-320-18	152	886-013-00	153
884-319-22	151	885-012-00	152	885-320-22	152	886-014-00	153
884-320-18	151	885-013-00	152	885-321-18	152	886-015-00	153
884-320-22	151	885-014-00	152	885-321-22	152	886-016-00	153
884-321-18	151	885-015-00	152	885-322-18	152	886-017-00	153
884-321-22	151	885-016-00	152	885-322-22	152	886-018-00	153
884-322-18	151	885-017-00	152	885-323-18	152	886-019-00	153
884-322-22	151	885-018-00	152	885-323-22	152	886-020-00	153
884-323-18	151	885-064-00	148	885-324-18	152	886-021-00	153
884-323-22	151	885-150-00	148	885-324-22	152	886-022-00	153
884-324-18	151	885-301-18	152	885-325-18	152	886-023-00	153
884-324-22	151	885-301-22	152	885-325-22	152	886-024-00	153
884-325-18	151	885-302-18	152	885-326-18	152	886-025-00	153
884-325-22	151	885-302-22	152	885-326-22	152	886-026-00	153
884-326-18	151	885-303-18	152	885-327-18	152	886-027-00	153
884-326-22	151	885-303-22	152	885-327-22	152	886-028-00	153
884-327-18	151	885-304-18	152	885-328-18	152	886-029-00	153
884-327-22	151	885-304-22	152	885-328-22	152	886-030-00	153
884-328-18	151	885-305-18	152	885-329-18	152	886-031-00	153
884-328-22	151	885-305-22	152	885-329-22	152	886-064-00	147
884-329-18	151	885-306-18	152	885-330-18	152	886-150-00	147
884-329-22	151	885-306-22	152	885-330-22	152	887-001-00	154
884-330-18	151	885-307-18	152	885-331-18	152	887-002-00	154
884-330-22	151	885-307-22	152	885-331-22	152	887-003-00	154
884-331-18	151	885-308-18	152	885-332-18	152	887-004-00	154
884-331-22	151	885-308-22	152	885-332-22	152	887-005-00	154
884-332-18	151	885-309-18	152	885-333-18	152	887-006-00	154
884-332-22	151	885-309-22	152	885-333-22	152	887-007-00	154
884-333-18	151	885-310-18	152	885-334-18	152	887-008-00	154
884-333-22	151	885-310-22	152	885-334-22	152	887-009-00	154
884-334-18	151	885-311-18	152	885-335-18	152	887-010-00	154
884-334-22	151	885-311-22	152	885-335-22	152	887-011-00	154
884-335-18	151	885-312-18	152	885-350-18	148	887-012-00	154
884-335-22	151	885-312-22	152	885-350-22	148	887-013-00	154
884-350-18	148	885-313-18	152	885-364-18	148	887-014-00	154
884-350-22	148	885-313-22	152	885-364-22	148	887-015-00	154
884-364-18	148	885-314-18	152	886-001-00	153	887-016-00	154
884-364-22	148	885-314-22	152	886-002-00	153	887-017-00	154

Numerical index

REF	Page	REF	Page	REF	Page	REF	Page
887-018-00	154	888-150-00	148	890-431-00	169	903-920-40	429
887-019-00	154	889-001-00	156	890-432-00	169		
887-020-00	154	889-002-00	156	890-433-00	169		
887-021-00	154	889-003-00	156	891-180-00	68		
887-022-00	154	889-004-00	156	891-181-00	68		
887-023-00	154	889-005-00	156	891-183-00	69		
887-024-00	154	889-006-00	156	891-184-00	69		
887-025-00	154	889-007-00	156	891-188-00	70		
887-026-00	154	889-008-00	156	891-189-00	70		
887-027-00	154	889-009-00	156	891-261-00	71		
887-028-00	154	889-010-00	156	891-262-00	71		
887-029-00	154	889-011-00	156	891-295-00	68		
887-030-00	154	889-012-00	156	891-296-00	68		
887-031-00	154	889-013-00	156	891-298-00	69		
887-064-00	147	889-014-00	156	891-299-00	69		
887-150-00	147	889-015-00	156	891-333-00	78		
888-001-00	155	889-016-00	156	891-334-00	78		
888-002-00	155	889-017-00	156	891-336-00	78		
888-003-00	155	889-018-00	156	891-337-00	78		
888-004-00	155	889-019-00	156	891-339-00	79		
888-005-00	155	889-020-00	156	891-340-00	79		
888-006-00	155	889-021-00	156	891-342-00	79		
888-007-00	155	889-022-00	156	891-343-00	79		
888-008-00	155	889-023-00	156	891-375-00	70		
888-009-00	155	889-024-00	156	891-376-00	70		
888-010-00	155	889-025-00	156	891-400-00	108		
888-011-00	155	889-026-00	156	891-404-00	108		
888-012-00	155	889-027-00	156	891-406-00	46		
888-013-00	155	889-028-00	156	891-407-00	46		
888-014-00	155	889-029-00	156	891-408-00	46		
888-015-00	155	889-030-00	156	891-409-00	46		
888-016-00	155	889-031-00	156	892-000-07	304		
888-017-00	155	889-064-00	148	892-000-14	304		
888-018-00	155	889-150-00	148	892-000-17	58, 306		
888-019-00	155	890-414-00	142, 147, 148, 157, 160, 169, 170	892-000-18	109, 305		
888-020-00	155			892-000-19	58, 306		
888-021-00	155	890-414-10	170	892-000-20	305		
888-022-00	155	890-420-00	169	892-000-21	109, 305		
888-023-00	155	890-421-00	169	892-000-22	50, 305		
888-024-00	155	890-422-00	169	892-000-26	253, 293, 306		
888-025-00	155	890-423-00	169	892-000-29	27, 306		
888-026-00	155	890-424-00	169	892-001-00	303		
888-027-00	155	890-425-00	169	892-002-01	303		
888-028-00	155	890-426-00	169	892-002-02	304		
888-029-00	155	890-427-00	169	892-003-00	304		
888-030-00	155	890-428-00	169	903-920-10	429		
888-031-00	155	890-429-00	169	903-920-20	429		
888-064-00	148	890-430-00	169	903-920-30	429		

Concept/Layout

Advertising department Dentaaurum

Catalog No. 24

989-781-00

In the interest of technical improvement, we reserve the right to make improvements and changes to the design, dimensions and materials of our products.

It is prohibited to make copies of this catalog or parts thereof without our explicit consent.

Registered design.

© Copyright Dentaaurum GmbH & Co. KG

Image credits:

- © Christian Ferrari® – p. 5
- © anson tsui / stock.adobe.com – p. 370/371
- © Artco / stock.adobe.com – p. 9
- © Astrovector / stock.adobe.com – p. 439
- © Beboy / stock.adobe.com – p. 430
- © Bojan / stock.adobe.com – p. 447
- © davooda / stock.adobe.com – p. 9
- © Drobot Dean / stock.adobe.com – p. 445
- © Fons Laure / stock.adobe.com – p. 102/103
- © Jezper / stock.adobe.com – back cover
- © kwasny221 / stock.adobe.com – p. 44/45
- © master1305 / stock.adobe.com – p. 18
- © nexusby / stock.adobe.com – p. 9
- © poppystyle / stock.adobe.com – p. 431
- © S / stock.adobe.com – p. 64/65
- © Sashkin / stock.adobe.com – p. 66/67
- © spiral media / stock.adobe.com – p. 9
- © Stephanie Swartz / stock.adobe.com – p. 100/101

Dentaurum

Germany | Benelux | España | France | Italia | Switzerland | Australia | Canada | USA
and in more than 130 countries worldwide.



DENTAURUM
QUALITY
WORLDWIDE
UNIQUE



Date of information: 2024-02
Subject to modifications

D
DENTAURUM

Turnstr. 31 | 75228 Ispringen | Germany | Phone +49 72 31/803-0 | Fax +49 72 31/803-295
www.dentaurum.com | info@dentaurum.com