

Frequency of drill application

Batch / LOT

Product image	Description	ø	REF	Frequency of application ¹																			
	Marking drill	ø 1.5 mm	382-015-00	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
	Stepped drill for ø 1.8 mm	ø 1.5 mm	302-118-00	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
	Stepped drill for ø 2.2 mm	ø 1.9 mm	302-122-00	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
	Stepped drill for ø 2.5 mm	ø 2.1 mm	302-125-00	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20

¹ **Note:** With careful and correct maintenance, the drills and thread taps can be used for 15 – 20 applications thanks to the high-quality materials and workmanship (see also Surgery + Prosthetics Manual 989-360-20).

Further information is available at www.dentaurum.com.