

Chroma Concept

Farb-Zuordnungstabelle zum Reproduzieren von VITA System 3D Master* Farben.

Farbe VITA System 3D-Master*		Dentin Mischtablette ceraMotion® Me Ti Zr		Opaker oder Liner	Incisal
M1	100%	CC D 1		CC1 bleach	CC11 bleach
M2	50%	CC D 1 + CC D 2	50%	CC1 bleach	CC11 bleach
M3	100%	CC D 2		CC1 bleach	CC11 bleach
1M1	50%	D A1 + CC D 2	50%	B1	I1
1M2	60%	CC D 1 + CC D 5	40%	B2	I1
2M1	75%	CC D 3 + CC D 1	25%	B2	I2
2L1.5	67%	D A2 + D B1	33%	B2	I2
2R1.5	50%	CC D 2 + CC D 6	50%	D2	I2
2M2	90%	D A2 + D B3	10%	B2	I2
2L2.5	80%	D B3 + CC D 1	20%	B2	I2
2R2.5	55%	CC D 2 + CC D 9	45%	B2	I2
2M3	75%	CC D 5 + CC D 2	25%	B2	I2
3M1	90%	CC D 3 + CC D 5	10%	CC2	I3
3L1.5	50%	CC D 3 + CC D 4	50%	CC2	I3
3R1.5	50%	CC D 3 + CC D 6	50%	CC2	I3
3M2	50%	CC D 3 + CC D 5	50%	CC2	I3

Farbe VITA System 3D-Master*		Dentin Mischtablette ceraMotion® Me Ti Zr		Opaker oder Liner	Incisal
3L2.5	50%	CC D 4 + CC D 5	50%	CC2	I3
3R2.5	50%	CC D 5 + CC D 6	50%	CC2	I3
3M3	90%	CC D 5 + CC D 3	10%	CC2	I3
4M1	50%	CC D 7 + CC D 10	50%	CC3	I3
4L1.5	75%	CC D 7 + CC D 9	25%	CC3	I3
4R1.5	75%	CC D 10 + CC D 8	25%	CC3	I3
4M2	50%	CC D 7 + CC D 9	50%	CC3	I3
4L2.5	75%	CC D 8 + CC D 10	25%	CC3	I3
4R2.5	75%	CC D 9 + CC D 7	25%	CC3	I3
4M3	50%	CC D 8 + CC D 9	50%	CC3	I3
5M1	50%	CC D 10 + CC D 11	50%	CC4	I3
5M2	50%	CC D 10 + CC D 12	50%	CC4	I3
5M3	50%	CC D 9 + CC D 12	50%	CC4	I3

* VITA System 3D-Master ist eine eingetragene Marke der Vita Zahnfabrik H. Rauter GmbH & Co.KG, 79704 Bad Säckingen, Deutschland