

COMPLIANCE

Certificate

Test Item: REMANIUM® CS + *

Supplier: DENTAURUM J.P. Winkelstroeter KG
Turnstraße 31, 75228 Ispringen, Germany

Study performed:

The following study (BSL Project No. 961104) was carried out in accordance with ISO 10993-5, 1992: "Tests for *in vitro* cytotoxicity"

CYTOTOXICITY

The cytotoxic effect of the test material was examined *in vitro* by the cell culture test using L 929-cells (fibroblast cell line): Examination of materials for medical devices with the extraction method and protein determination.

Result: No leachable substances were released in cytotoxic concentrations (leading to a cell growth inhibition of greater than 30%) from the test item.


BSL BIOSERVICE Scientific Laboratories GmbH

Behringstraße 6/8

82152 Planegg

Germany

*originally tested as REMANIUM CS®, Lot No.: 6089


Dr. Achim Albrecht

Head of *in vivo* Pharmacology/Toxicology

Date: 02 June 2009